

NOV 12 1996

K963921

EXHIBIT 8

RESPONSE TO SMDA OF 1990

SUMMARY OF SAFETY AND EFFECTIVENESS DATA

Kenneth J. Berk
80 Oakland Street
PO Box 780
Watertown, MA 02272-0780
USA

TELEPHONE: (617) 926-6666
FAX: (617) 926-6262

DEVICE NAME: Pulpdent Seal-Rite™ UDMA

PREDICATE DEVICE: Pulpdent *Seal-Rite* and *Seal-Rite Low Viscosity*
Vivadent/Ivoclar *Helioseal* and *HelioSeal-F*
Kulzer *Estiseal*
Ultradent *Ultra-Blend*

DESCRIPTION AND INTENDED USE: *Pulpdent Seal-Rite™ UDMA* is a radiopaque, fluoride-releasing, visible light cure dental material for sealing the pits and fissures of teeth. It contains no Bisphenol A. It is supplied as a medium consistency paste, which flows readily from cusp to cusp and is easily distinguishable from enamel. *Pulpdent Seal-Rite™ UDMA* cures completely after light curing for 20 seconds. Various formulations (with and without color, with and without fluoride) will be offered to allow for the clinician's discretion and preference in different situations. *Pulpdent Seal-Rite™ UDMA* complies with ANSI/ADA Specification No. 39 for Pit and Fissure Sealants.

COMPARISON WITH PREDICATE PRODUCTS: *Pulpdent Seal-Rite™ UDMA* is substantially equivalent in composition and/or intended use to the above mentioned predicate products. Please see EXHIBIT 5 for the entire comparison.

SAFETY AND EFFECTIVENESS: According to the NIH Technology Assessment Conference on Effects and Side-Effects of Dental Restorative Materials, "General usage of these materials over about 20 years indicates a high benefit-to-risk ratio...both composites and glass ionomers are relatively trouble free. There is no evidence of short or long term risk...there is no suspicion of any problems after literally billions of procedures in the United States." Please see EXHIBIT 6 for the complete papers from the conference.