

MAR - 5 1997

K 83947

ATTACHMENT I
(Rev. 12/96)

510(K) SUMMARY

Epilaser™ Normal Mode Ruby Laser

This 510(k) Summary of safety and effectiveness for the Epilaser™ Normal Mode Ruby Laser is submitted in accordance with the requirements of SMDA 1990 and follows Office of Device Evaluation guidance concerning the organization and content of a 510(k) summary.

Applicant: Spectrum Medical Technologies, Inc.

Address: 45 Hartwell Avenue
Lexington, MA 02173

Contact Person: Mr. Sanford Lane
President, Spectrum Medical
Technologies, Inc.

Telephone: 617-676-7300
617-676-7330 (Fax)

Preparation Date: December 26, 1996
(of the Revised Summary)

Device Trade Name: Epilaser™ Normal Mode Ruby Laser

Common Name: Ruby Laser

Classification Name: Laser surgical instrument for use in
general and plastic surgery and in
dermatology (see 21 CFR 878-4810)
Product Code: GEX
Panel: 79

**Legally Marketed
Predicate Device(s):** Epilaser™ Normal Mode Ruby Laser
(K955612)
Thermolase LT-100 Nd:YAG Laser
(K950019)

**Description of the
Device:** The Epilaser™ Normal Mode Ruby Laser
operates at 694.3 nanometers with
pulse durations of 0.2 - 3.0 msec.
The energy is delivered to the
treatment site by an articulated arm.
A water-cooled handpiece (7 mm or 10
mm) is firmly held against the skin
for 2-3 seconds and the laser is then
fired.

The EpilaserTM Normal Mode Ruby Laser is capable of producing energy fluences of 10-75 J/cm².

Intended Use:

The EpilaserTM Normal Mode Ruby Laser is intended to effect hair removal of patients with skin types 1 - 4 through selective photothermolysis of hair follicles in dermatology and plastic surgery. This results in a prolonged growth delay.

Comparison:

The EpilaserTM Normal Mode Ruby Laser for hair removal has the same technological characteristics as the predicate Epilaser (K955612). It differs from the Thermo-lase LT-100 Nd:YAG Laser predicate device (K950019) in that the Thermolase applies laser energy at 1064 nanometers in 10 nsec pulses onto a cream containing carbon. The devices operate by selective photothermolysis.

Performance Data:

Animal and clinical studies were conducted to provide assurance that differences in the specifications of the EpilaserTM Normal Mode Ruby Laser and the claimed predicate device for hair removal did not result in different performance during use.

Animal Studies:

Biopsies of control skin showed no alterations the follicular epithelium or dermal collagen. Examination of irradiated skin showed damage to the follicular epithelium. Focal collagen damage in reticular dermas occurred immediately adjacent to the hair follicle.

Observed effects increased with fluence levels.

Clinical Studies:

Observations in an initial study were recorded prior to treatment and at 1, 3, and 6-months after treatment. There was no scarring in any subject.

The initial study demonstrated that selective photothermolysis targeting melanin in the human hair follicle is an effective tool for hair removal, resulting in a prolonged growth delay.

Subsequent studies, including patients treated for facial hair, confirmed the results of the initial study.

CONCLUSION:

Based on the foregoing, the EpilaserTM Normal Mode Ruby Laser is substantially equivalent to legally marketed claimed predicate devices, i.e., the EpilaserTM for dermatology and plastic surgery and the Thermolase LT-100 Nd:YAG laser for hair removal.