

DEC 12 1996

K963988

510(k) Premarket Notification  
Medi-tech Intravascular Infusion Device  
October 3, 1996

**ATTACHMENT I**

**SUMMARY OF SAFETY AND EFFECTIVENESS**

Pursuant to §513(i)(3)(A) of the Food, Drug, and Cosmetic Act, Boston Scientific Corporation is required to submit with this Premarket Notification either an "...adequate summary of any information respecting safety and effectiveness or state that such information will be made available upon request of any person." Boston Scientific Corporation choose to submit a summary of information respecting safety and effectiveness. According to §513(i)(3)(B), "Any summary under subparagraph (A) respecting a device shall contain detailed information regarding data concerning adverse health effects..."

The summary regarding the adverse health effects of the proposed Intravascular Infusion Device is as follows:

**Trade Name:** Intravascular Infusion Device

**Manufacturer:** BSC/SciMed Life Systems  
6655 Wedgewood Road  
Maple Grove, MN 55311-3646

**Device Generic Name:** Diagnostic Intravascular Catheter

**Classification:** According to Section 513 of the Federal Food, Drug and Cosmetic Act, the device classification is Class II, Performance Standards (CFR 870.1200)

**Predicate Device:** **Balt Magic Infusion Catheter**  
Manufactured by:  
Target Therapeutics, Inc.  
130 Rio Robles  
San Jose, CA 95134

**Transend Steerable Guidewire**  
Manufactured by:  
BSC/Scimed  
6655 Wedgwood Road  
Maple Grove, MN 55369

**510(k) Premarket Notification**  
**Medi-tech Intravascular Infusion Device**  
**October 3, 1996**

**Product Description:**

The Intravascular Infusion Device is a single lumen device constructed with progressively softer characteristics from proximal to distal end to aid in selective placement in the vasculature.

**Indications for Use:**

The Intravascular Infusion Device is indicated for the infusion of diagnostic agents into the general vasculature including the peripheral, coronary and neurovasculature.

**Safety and Performance:**

The following Functional Testing was performed on the proposed device:

1. Infusion Rate Testing
2. Dynamic Infusion Testing
3. Static Burst Testing
4. Tip Flexibility Testing
5. Torque Response Testing
6. Torsion Strength Testing
7. Tensile Strength Testing
8. Coatings Testing
9. *In vivo* Testing

In addition, the following Biocompatibility Testing was performed:

1. Cytotoxicity
2. Hemolysis
3. Acute Systemic Toxicity
4. Intracutaneous Toxicity
5. Sensitization
6. Pyrogenicity
7. USP Physiochemical Testing

**Conclusion:**

Based on the Indications for Use, technological characteristics and safety and performance testing, the Intravascular Infusion Device has been shown to be safe and effective for its intended use.