

DEC 18 1986

K9164250
STERLING

Diagnostic Imaging, Inc.

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

SUBMITTED BY:

Jean E. Bartlett
(302) 451-0282
Regulatory Affairs
Sterling Diagnostic Imaging, Inc.
Glasgow Business Community
Building 600/Mailbox 630
Newark, DE 19714-6101
Fax # (302) 451-0386

DEVICE NAME:

LINX™ Networking System

PREDICATE DEVICE:

LINX™ Networking System with LINX™ Laser Imager (K900935)

DEVICE DESCRIPTION AND INTENDED USE:

This Premarket Notification 510(k) is being filed to cover additions to the LINX™ Networking System with LINX™ Laser Imager (K900935) to include all capabilities that we are adding to meet market demands. The basic product concept remains - to connect various image input devices to one or more output devices. That basic concept has been expanded to meet market demands:

- new image source interfaces
- new control interfaces
- new network interfaces
- more output devices supported
- new data types managed
- image transfer mechanisms
(some output devices may require handling
in the transmission path to the output
device to compensate for differences in
devices)

DEVICE COMPARISON TO PREDICATE DEVICE:

The LINX™ Networking System is substantially equivalent to the LINX™ Networking System with LINX™ Laser Imager (K900935) as summarized below:

Characteristic	LINX Networking System	LINX Networking System (Expanded)
Image Sources		
Video	X (up to 45Mhz pixel clock speed)	X (up to 160Mhz pixel clock speed)
Digital - 405412	X	X
digitizer interface		X
computed radiography int		X
removable media		X
network interfaces		X
DICOM Devices		X
Network connections		
token - ring	X	X
ethernet		X
ATM		X
FDDI		X
protocols		
QNX	X	X
TCP/IP		X
Interface methods		
point - point	X	X
network	X	X
LAN	X	X
WAN		X
Operator Interfaces		
hardware keypad	X	X
Windows		X
remote host		X
Data Types		
image	X	X
text	X	X
patient info		X
Output devices		
LINX Laser Imager	X	X
3M 952	X	X
3M 831	X	X
3M 969	X	X
LP400	X	X
home tolerad systems		X
DICOM Devices		X
File systems		X

Characteristic	LINX Networking System	LINX Networking System -(Expanded)
Accessory Device Interfaces		X
Densitometer		X
Bar Code Entry		X
Statistics Recording	X	X
Image Data		X
Raw Data	X	X
image transfer mechanisms		
straight thru	X	X
crop		X
LUT	X	X
sharpness filter		X
interpolate		X
decimate		X
compress (lossless)		X