JAN 1 0 1997

#### 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CPR 807.92.

The assigned 510(k) number is:  $\underline{K964580}$ 

Submitted By: Jeffrey A. Dann, M.D.

President NEBL Inc.

44 Terrace Drive Worcester, MA 01609 Phone (508) 342-2224 Fax: (508) 345-0137

Date Prepared: 11/14/96

#### Device Name

Proprietary Name: NEBL, Inc. CAP-AID (Improved)

Common Name: Male Incontinence Cap

Classification Name: Penile Clamp

## Substantial Equivalence

The NEBL CAP-AID (Improved) is exactly similar in indications, design and features to NEBL's original CAP-AID device.

### Device Description

The CAP-AID (Improved) is exactly similar in indications, design and features to the original CAP-AID device. Only the materials are different. The material changes requested in the "Improved" version allow additional comfort and ease of use for the patient.

### Intended Use

The NEBL CAP-AID (Improved) is designed to effectively occlude the male meatus thereby preventing urinary leakage.

# Indications for Use

The product prevents urinary leakage in men with urinary incontinence.

# Substantial Equivalence Comparison

|                      | NEBL Inc. CAP-AID (Improved)   | NEBL Inc. CAP-AID<br>(original) K954959  |
|----------------------|--|--|
| Indication           | Prevent urine leakage in male incontinence   | Prevent urine leakage in male incontinence   |
| Mode of<br>Operation | External device which occludes the distal urethra  | External device which occludes the distal urethra  |
| Application          | 1. adheres to the glans penis 2. a shallow closed ended trough occludes the meatus 3. Patient removes device to void and re-applies a new device after each void           | 1. adheres to the glans penis 2. a shallow closed ended trough occludes the meatus 3. Patient removes device to void and re-applies a new device after each void   |
| Material             | 1. Kraton Thermoplastic elastomer G-2706 2. C-Flex Thermoplastic elastomer - CPT Technologies 3. medical grade pressure sensitive skin adhesive Dow Corning Bio PSA 7-4500 | 1. Kraton Thermoplastic elastomer G-2705 2. medical grade pressure sensitive skin adhesive Polyken Technologies #DFPE 170 or 3M Healthcare polyethylene film #1509 |
| Sterility            | Sterile/disposable   | Sterile/disposable   |

The determination of substantial equivalence is not based on an assessment of any performance date, clinical or non-clinical.

Submitter's Signature

Date