

U.S. Food and Drug Administration  
510K Document Mail Center HFZ-401  
Center for Medical Devices  
390 Piccard Drive  
Rockville, Maryland 20850

K964736

FEB 27 1997

November 15, 1996

To: Document Control Clerk

This summary of 510K safety and effectiveness for the BERGEN 610 Coagulator is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

A 0-50 digital display power setting indicator provides the operator a parallax free reading of the power setting (+/- 5 watts, full scale @ 100 ohms).

A bar graph display output current monitor provides the operator a relative measure of the output current amplitude and to denote the presence of output energy during activation.

The handswitch is powered by an independent low voltage (7 Vdc) isolated power supply.

Audio and Visual monitors are in accord with IEC 601-1-2 and ANSI HF-18 guidelines.

RF and Low Frequency leakages are well within IEC 601-1-2 and ANSI HF-18 safety guide lines lines for isolated (body floating) bipolar coagulators.

his device is substantially equivalent to existing approved devices.

Sincerely,



Roger Oosten  
Bergen Mfg  
9345 Rookery Road  
New Port Richey, Florida 34654

Date 11/15/96