K970842 Act MAY 28 1997

Section 513(j) of the Federal Food, Drug and Cosmetic Act **Summary of Safety and Effectiveness**

March 5, 1997

General Provisions 1.

Common or Usual name:

Infusion Catheter

Proprietary name:

Cordis Endeavor Infusion Catheter

Name and Address of Applicant:

Cordis Corporation

Miami Lakes Operation Center

14201 NW 60 Avenue Miami Lakes, FL 33014

Name of Predicate Devices

> SciMed Dispatch and Dispatch Gold Infusion Catheters LocalMed InfusaSleeve and InfusaSleeve II Infusion Catheters Interventional Innovations Segue Infusion Catheter

III. Classification

Infusion catheters are class II devices according to 21 CFR 870.1210.

Performance Standards 14.

> Performance standards have not been established by the FDA under Section 514 of the Food, Drug and Cosmetic Act.

8. **Intended Use and Device Description**

> The Endeavor Infusion Catheter is intended to deliver solutions to the coronary and peripheral vasculature. The device is an over-the-wire design with a distal infusion region and a proximal hub. The infusion region is indicated by a central radiopaque marker band.

\'I. **Biocompatibility**

> All appropriate biocompatibility testing was performed, and successfully passed, on the materials used for the Endeavor Infusion Catheter.

VII. In vitro Testing

A series of *in vitro* and *in vivo* tests were performed to assure that the introduction of the Endeavor Infusion Catheters does not raise new issues of safety and effectiveness. All test results met or exceeded established specifications.

VIII. Summary of Substantial Equivalence

The Endeavor Infusion Catheter is designed for the infusion of solutions to the coronary and peripheral vasculature. The Endeavor Infusion Catheters have similar intended uses, design characteristics and dimensions as the predicate devices.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

MAY 28 1997

Ms. Tamara Yount Cordis Corporation P.O. Box 025700 Miami, Florida 33014

Re: K970842

Cordis Endeavor Infusion Catheter

Regulatory Class: II (two)

Product Code: 74 KRA Dated: March 6, 1997 Received: March 7, 1997

Dear Ms. Young:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Thomas J. Callahan, Ph.D.

Thomas J. Callahon

Director

Division of Cardiovascular,

Respiratory, and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Statement of Intended Use Cordis Endeavor Infusion Catheter

The intended use statement of this product is:

The Cordis Endeavor Infusion Catheter is intended to deliver solutions, such as heparinized saline, and thrombolytic agents, such as urokinase, to the coronary and peripheral vasculature.

510(k) number

K970842

(To be assigned by FDA)

(Divisio

Division

wascural Respiratory

and i

STOR MURITION