



**Bio-Rad  
Laboratories**

ECS Division  
3726 E. Miraloma Avenue  
Anaheim, CA 92806  
Telephone (714) 630-6400  
Toll Free (800) 854-6737

K971954

June 11, 1997

## 510(k) Summary

### Submitter

Bio-Rad Laboratories, ECS Division  
3726 E. Miraloma Avenue  
Anaheim, CA 92806  
(714)630-6400  
Fax (714)666-1383

### Contact Person

Elizabeth Platt

### Date of Summary Preparation

May 23, 1997

### Device (Trade & Common Name)

Liquichek Urine Chemistry Control

### Classification Name

Class I, 75JJY  
CFR 862.1660: Quality Control Material (Assayed and Unassayed)

### Devices to Which Substantial Equivalence is Claimed

Human Urine Control  
Quantimetrix Corporation  
Redondo Beach, CA

### Statement of Intended Use

Liquichek Urine Chemistry Control is intended for use as an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in this package insert.

Description of the Device

Liquichek Urine Chemistry Control is prepared from human urine with added constituents of human and non-human origin and pure chemicals. The control is provided in liquid form for convenience. This product contains <0.1% sodium azide as a preservative.

Statement of How Technological Characteristics Compare to Substantial Equivalence Device

A table is provided below comparing the similarities between the Bio-Rad Liquichek Urine Chemistry Control and the devices to which substantial equivalence is claimed.

	Bio-Rad Laboratories Liquichek Urine Chemistry Control	Quantimetrix Corporation Human Urine Control
Intended Use	an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in this package insert.	a means for monitoring human urine assay methods to validate quantitation of patient samples.
Levels	Two	Two
Form	Liquid	Liquid
Open Vial Claim	30 Days at 2-8°C	24 Months at 2-8°C
Matrix	Human Urine	Human Urine
Storage	2-8°C	2-8°C



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JUN 11 1997

Elizabeth Platt  
• Staff Regulatory Affairs Representative  
Bio-Rad Laboratories  
3726 E. Miraloma Avenue  
Anaheim, California 92806

Re: K971954  
Liquichek Urine Chemistry Control  
Regulatory Class: I  
Product Code: JJY  
Dated: May 23, 1997  
Received: May 28, 1997

Dear Ms. Platt:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

\* If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: \_\_\_\_\_

Device Name: Liquichek Urine Chemistry Control


Indications for Use:

Liquichek Urine Chemistry Control is intended for use as an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.

(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

(Concurrence of CDRH, Office of Device Evaluation)

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number 10971954

Prescription Use \_\_\_\_\_

OR Over-The Counter Use \_\_\_\_\_