

Jl. Jend. Galot Subroto Kav. 22

Jakarta - 12930

11 th. Floor WISMA ARGO MANUNGGAL

INDUSTRIES TATA RUBBERINDO

Factory : Jl. Raya Serang Km. 13.8 : (021) 5962435 (Hunting) Phone

: (021)5962436

Fax.

Cikupa - Tangerang

K991158

Page Numbers 1 of 2

"510 (K)" SUMMARY (K991158)

(1) Name of applicant

Address

: Mr. Andy Tanaka

: PT. Tata Rubberindo Industries

Jl. Raya Serang Km. 13.8

Cikupa - Tangerang Jakarta - Indonesia

Phone No. 62-21-5962435 62-21-5962436 Fax No.

The contact persons within the firm as well as in U.S.A are given below:

Contact person in firm

: Mr. Andy Tanaka

62-21-5962436 Fax No.

Contact person in U.S.A

: Cindy Tung

: 626-913-1498 Fax No.

(2) Device details

Trade Name

: Latex Examination Gloves - Powder Free

Classification Name

: Patient Examination Gloves

Product Code

: Latex 80 LYY

(3) Equivalent device legally

marketed

: Class I Latex Examination Gloves 80 LYY Powder - Free meeting ASTM D 3578-95

(4) Intended use

: A patient examination glove is a disposable device intended for medical purpose that is worn on or finger to prevent examiner's hand

contamination between and patient examiners.

0.10 mm

0.10 mm

(5) Technological characteristic of the gloves.

Sizes	XS	S	M	L	XL
Length	240 mm	240 mm	240 mm	240 mm	240 mm
Width	80 < mm	80±10 mm	95±10 mm	111±10 mm	>111mm
Thickness 1. Cuff (min) 2. Palm (min)	0.10 m 0.10 m	m 0.10 mm m 0.10 mm		0.10 mm 0.10 mm	0.10 mm 0.10 mm

b. Physical Properties

:	After aging at 70°C 168 hrs.	
Tensile Strength : Ultimate Elongation :	21 Mpa 700 % (min.)	16 Mpa 600 % (min.)

(6) Performance data is the same as mentioned immediately above.

3. Finger Tip (min) 0.10 mm 0.10 mm 0.10 mm

- (7) Clinical date is not needed for gloves or for most devices cleared by the 510 (K) process.
- (8) Non-clinical data We certify that these gloves meet ASTM D 3578 Standard. Meets FDA pinhole requirement. Meets labeling claim.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JUL 19 1999

PT. Tata Rubberindo Industries C/O Mr. Andy Tanaka Jl. Raya Serang Km. 13.8 Cikupa - Tangerang INDONESIA

к991158 Re:

Trade Name: Latex Examination Gloves - Powder Free

Regulatory Class: Product Code: LYY Dated: May 12, 1999 Received. My 14, 1999

Dear Mr. Tanaka:

This letter corrects our substantially equivalent letter of May 28, 1999 regarding the address.

We have reviewed your Section 510(k) notification of intent to market the device ific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of yoto devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General (QS) regulation (21 CFR Part 820) and that, through periodic QS inspections, FDA will verify such assumptions.

Failure to comply with the GMP regulation may result in regulatory action. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to continue marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4639. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at their toll free number (800) 638-2041 or at (301) 443-6597 or at its internet address http://www.fda.gov/cdrh/dsmamain.htm/".

sincerely

Timothy A. Ulatowski

Director

Division of Dental, Infection Control, and General Hospital Devices

Office of Device Evaluation

Center for Devices and Radiological Health



PT TATA RUBBERINDO INDUSTRIES

11 th. Floor WISMA ARGO MANUNGGAL Jl. Jend. Gatot Subroto Kav. 22 Jakarta - 12930 Factory: JL Raya Seranç Phone: (021)5962435

Fax. : (021)5962436

Cikupa - Tange

ANNEXURE II

K991158

INDICATION FOR USE

Applicant

: Mr. Andy Tanaka

Device Name

: Latex Patient Examination Gloves Powder Free

Indication for use

A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and examiners.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Dental, Infection Control,

and General Hospital Devices 5:0(k) Number _____ QQ 1

Prescription Use

(Per 21 CFR 801.109)

OR

Over-The-Counter Use

(Optional Format 1-2-96)