

**SECTION 12.0**  
**510(k) Summary**

K991436

Company Name: Sterngold ImplaMed  
Address: 23 Frank Mossberg Drive  
Attleboro, MA 02703  
Registration #: 2921595  
Contact Person: Gordon Nelson  
Date Prepared: Thursday, April 22, 1999

Classification Name: Endosseous Implant (DZE)  
Common Name: Narrow Platform Screw Implant, Implant Fixture  
Trade Name: Sterngold ImplaMed Hex Screw Implant,  
Narrow Platform Screw Implant

**Device Description**

Device consists of titanium screws, and titanium alloy prosthetics, brass and stainless steel restorative components, and stainless steel and titanium alloy surgical instruments.

**Intended Use**

Device can be used in dental implant applications for oral rehabilitation of edentulous and partially dentate patients in the maxillae and mandible. Implant retained restorations may consist of single crowns or bridges as well as complete or partial dentures.

**Technological Characteristics**

The modified devices utilize narrow diameter titanium screws and related prosthetics and instrumentation. Materials are identical to those currently used in our predicate devices. The Narrow Platform screw products are compatible with current Sterngold ImplaMed installation instrumentation and prosthetic procedures.

**Comparative Products**

Sterngold ImplaMed currently has permission to market Titanium Screw Implants, Titanium and Titanium alloy prosthetics, and stainless steel instrumentation. Narrow Platform and diameter screw implants are currently marketed by several companies. Sterngold ImplaMed Narrow Platform implants prosthetics and instruments are substantially equivalent to these marketed devices in design, materials, performance and intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY 7 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Gordon Nelson  
Director of Operations  
Sterngold ImplaMed  
23 Frank Mossberg Drive  
Attleboro, Massachusetts 02703-0967

Re: K991436  
Trade Name: Sterngold ImplaMed Narrow Platform Dental  
Implants and Related Prosthetic and Instrumentation  
Regulatory Class: III  
Product Code: DZE  
Dated: April 22, 1999  
Received: April 26, 1999

Dear Mr. Nelson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531


Page 2 - Mr. Nelson

through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**SECTION 10.0**  
**Indications for Use**

Page 1 of 1

510(k) Number (if known): K991436

Device Name: \_\_\_\_\_

**Indications for Use:**

The Sterngold Implamed Implant can be used in dental implant applications for oral rehabilitation of edentulous and partially dentate patients in the maxilla and mandible. Implant retained restorations may consist of single crowns or bridges as well as complete or partial dentures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Runner

(Division Sign-Off)

Division of Dental, Infection Control,  
and General Hospital Devices

510(k) Number K991436

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-the -Counter Use \_\_\_\_\_