

AUG 3 | 1999

## 510K Summary

Regulatory Authority: Safe Medical Devices Act of 1990. CFR 807.87

**Company Name:**

Luma Lite, Inc.  
477 Marina Pkwy #400  
Chula Vista, CA 91910

**Company Contact:**

Joe Forehand  
Luma Lite, Inc.  
477 Marina Pkwy #400  
Chula Vista, CA 91910  
(619) 318-1251

**Device Name:**

Luma Light 2000 Cure Light

**Predicate Devices:**

Kuring and Bleaching Light	K962376	Kreativ, Inc.
Apollo 95E	K981948	DMD
Pac Light	K952333	American Dental Technologies.

**Device and indications for use:**

The Luma Light 2000 Cure Light is a high intensity filtered light source, and optical delivery system, that produces curing of photo-activated dental or photo-activated bleaching formulations at four (4) times the rate of conventional light systems. The Luma Light 2000 Cure Light cures all photo activated dental materials such as those used in direct restorations, sealants, bonding agents and primers.

**Discussion:**

Since the intended use and technical specifications of the Luma Light 2000 Cure Light are virtually identical to the predicate devices and the differences in the device only make it easier to use, more reliable and more adaptable to a variety of dental practice situations, we believe that the Luma Light 2000 Cure Light is substantially equivalent to the predicate devices and can be marketed under Section 510 (k) of the FD&C Act.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG 31 1999

Mr. Joseph M. Forehand  
Vice President of Operations  
Luma Lite, Incorporated  
477 Marina Parkway #400  
Chula Vista, California 91910

Re: K992102  
Trade Name: Luma Light 2000 Cure Light  
Regulatory Class: II  
Product Code: EBF  
Dated: June 16, 1999  
Received: June 21, 1999

Dear Mr. Forehand:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of

Page 2 - Mr. Forehand

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K992102

**Indications Statement**

510(k) Number: K992102

Device Name: Luma Light 2000 Cure Light

**Indication for Use:**

The Luma Light 2000 Cure Light is a high intensity filtered light source, and optical delivery system, that produces curing of photo-activated dental composites or photo-activated bleaching formulations at four (4) times the rate of conventional light systems. The Luma Light 2000 Cure Light cures all photo activated dental materials such as those used in direct restorations, sealants, bonding agents and primers.

*Suzanne P. [Signature]*  
Division Sign-Off  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K992102

(Please do not write below this line—continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over the Counter Use ☐

Optional Format 1-2-96