

MAR 20 2000

K994290 (P.1 of 3)

Bard Interventional Products Division
C.R. Bard, Inc.
129 Concord Road, Bldg #3
P.O. Box 7031
Billerica, MA 01821-7031
978-663-8989



VI. 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION

As required under Section 513(i)(3)(A) of the Safe Medical Device Act of 1990, an adequate summary of any information respecting the safety and effectiveness follows:

A. General Information

Name and Address of Submitter:	Bard Interventional Products Division, C. R. Bard, Inc. 129 Concord Road, Building #3 Billerica, MA 01821-7031
Contact:	Beth A. Zis, R.A.C. Manager of Regulatory Affairs Phone: (978) 262-4866 Fax: (978) 262-4878
Date of Summary:	December 6, 1999
Name of Device:	
Trade Name/Proprietary Name:	Bard® Endoscopic Suturing System
Common Usual Name:	Endoscopic Suture Device
Classification Name:	78KOG – Endoscopic Suture Device

B. Predicate Devices:

<i>Company</i>	<i>Trade Name</i>	<i>510(k)#</i>
1. U. S. Surgical Corp.	Auto Suture™ ENDO STITCH™ Suturing Device	#K934738
2. Ethicon Endo-Surgery, Inc.	ENDOPATH® Endoscopic Suturing System	#K980022

C. Description

The Bard® Endoscopic Suturing System consists of a capsule assembly with suction tubing and a fixed head knotpusher that attaches to the distal end of a

flexible endoscope or a through-the-scope knot pusher, plus a needle assembly, pusher wire, guidewire and suture cutter that all pass through an endoscope's biopsy channel. Also, included is a suture loader to facilitate loading suture tags into the needle, a syringe for flushing the suction tubing, and an detachable knot pusher handle. Only Bard® Suture Tags and the Bard® Endoscopic Handle may be used with the Bard® Endoscopic Suturing System.

D. Intended Use:

The Bard® Endoscopic Suturing System is used for endoscopic placement of suture(s) in the soft tissue of the esophagus and stomach and for the approximation of tissue for the treatment of symptomatic Gastroesophageal Reflux Disease.

E. Technological Characteristics Summary:

The technological characteristics of the Bard® Endoscopic Suturing System is the same or similar to the predicate devices, in that the materials used to manufacture these products are the same type of medical grade stainless steels and plastics. The products all share common features such as a sterile, stainless steel needle housed in a capsule or suture loading unit at the distal end of the device. They all suture soft tissue by manually actuating the needle with a handle mechanism. They all are designed to allow reloading of sutures to deliver multiple stitches under endoscopic visualization. Further, the Bard® Endoscopic Suturing System and the predicated devices have the same or similar intended use, that is to place stitches and tie suture material to approximate soft tissue under endoscopic visualization.

F. Performance Data:

Biocompatibility, in vitro bench testing, animal and clinical testing has been completed and supports the safety and effectiveness of the Bard® Endoscopic Suturing System for its intended use.

The biocompatibility test results show that the materials used in the design and manufacture of the device are non-toxic and non-reactive to biologic tissues consistent with their intended use.

Bench test results show that the materials chosen and the design utilized in manufacturing the Bard® Endoscopic Suturing System, meet the established specifications necessary for consistent performance during their intended use.

The animal studies support the safety and effectiveness of the Bard® Endoscopic Suturing System. The studies document there were no complications/adverse reactions reported, and that the system is capable of safely and consistently placing sutures in the soft tissues within the stomach and esophagus.

Finally, the summary of human clinical experience under a 64 patient prospective randomized investigation and the non-U. S. patient experience supports the safety and effectiveness of the Bard® Endoscopic Suturing System for its intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR 20 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Beth A. Zis
Regulatory Affairs Manager
C. R. Bard, Inc.
Bard Interventional Products Division
129 Concord Road, Building #3
P.O. Box 7031
Billerica, Massachusetts 01821-7031

Re: K994290
Trade Name: Endoscopic Suturing System
Regulatory Class: II
Product Code: KOG
Dated: December 20, 1999
Received: December 21, 1999

Dear Ms. Zis:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Beth A. Zis

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", followed by a small, stylized flourish or mark.

Celia M. Witten, Ph.D., M.D.
Director
Division of General and
Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K994290

510(k) Number (if known): Not Known

Device Name: Bard® Endoscopic Suturing System

Indications for Use: For endoscopic placement of suture(s) in the soft tissue of the esophagus and stomach and for approximation of tissue for the treatment of symptomatic Gastroesophageal Reflux Disease.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-the-Counter Use _____

NRD for CDRH (Optional Format 1-2-96)
(Division Sign-Off)
Division of General Restorative Devices
510(k) Number K994290