

**510(k) SUMMARY**

**Submitter's Name:** American Medical Systems, Inc.

**Address:** 10700 Bren Road West  
Minnetonka, MN 55343

**Tel:** 952-930-6120

**Fax:** 952-930-6496

**Contact Person:** Mark McIntyre

**Date of Summary Preparation:** April 19, 2002

**Device Common Name:** Surgical Mesh, Sling, Urethral Sling

**Device Trade Name:** SPARC™ Sling System

**Device Classification Name:** Surgical Mesh, polymeric

**Predicate Device:** SPARC™ Sling System – K011251, K013355, K020663

**Device Description**

The SPARC™ Sling System as currently marketed is a sterile, single use procedure kit consisting of:

- Two stainless steel, curved, 22-cm long, needle passers (also called insertion tools).
- One piece of AMS Polypropylene sling mesh with attached dilating connectors. The AMS Polypropylene sling mesh is constructed of polypropylene monofilament that is precut to 1.1 cm width x 50cm length. A fixed blue polypropylene tensioning suture runs through the middle of the sling mesh. Two plastic sheaths that overlap in the center of the sling mesh, cover the sling mesh and protect it during placement.

Dilating connectors are attached to either end of the plastic sheaths. The dilating connectors are used to attach to the vaginal ends of the SPARC™ needle passers during the procedure to facilitate sling placement.

- In K013355, two blue colored plastic cystoscopy aids are included in the kit in order to facilitate cystoscopic viewing of the bladder. The use of these cystoscopy aids is optional. The proposed device that is the subject of this 510k will not include the cystoscopy aids.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY 15 2002

Mr. Mark McIntyre  
Director, Regulatory Affairs and Biostatistics  
American Medical Systems, Inc.  
10700 Bren Road West  
Minnetonka, MN 55343

Re: K021263  
Trade/Device Name: SPARC™ Sling System  
Regulation Number: 878.3300  
Regulation Name: Surgical Mesh  
Regulatory Class: II  
Product Code: FTL  
Dated: April 19, 2002  
Received: April 22, 2002

Dear Mr. McIntyre:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
for Celia Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

INDICATIONS FOR USE ENCLOSURE

510(k) Number:

K021263

Device Name:

SPARC™ Sling System

Indications for Use:

The SPARC™ Sling System is intended for the placement of a pubourethral sling for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE  
IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
(Division Sign-off)

Division of General and Restorative Devices

510(k) Number \_\_\_\_\_

Prescription Use X  
(Per 21 CFR801.109)

OR

Over the Counter Use \_\_\_\_\_

Miriam C. Provost  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K021263