

K022186

FEB 26 2003

**510(k) Summary**  
**Scleral Plugs**  
**(per 21 CFR 807.92)**

**1. SUBMITTER NAME AND ADDRESS**

MicroVision, Inc.  
34 Folly Mill Road, Suite 200  
P.O. Box 1651  
Seabrook, NH 03874-1651

**Contact Person:** Leonard Kastrilevich, President  
Telephone: 603-474-5566  
Fax: 603-474-9336  
Email: [lkast@micro-vision.net](mailto:lkast@micro-vision.net)

**Date Prepared:** July 01, 2002  
**Date Revised:** December 20, 2002

**2. DEVICE NAME**

Trade Name: Scleral Plug  
Proprietary Name: Scleral Plugs, 19 and 20 Gauge  
Classification Name: Plug, Scleral

**3. PREDICATE DEVICE**

(K854507) Storz MVS 19 & 20 Gauge Scleral Plug  
Storz Instrument Co., St. Louis

**4. DEVICE DESCRIPTION**

Scleral plugs, 19 and 20 gauge, for use in ophthalmic surgical procedures. Both 19 and 20 gauge plugs are stainless steel. The 19 gauge plugs are gold plated.

**5. INTENDED USE**

The scleral plugs are intended to maintain patency of a previously made incision in the sclera of the eye, during ophthalmic surgical procedures.

**6. BASIS FOR SUBSTANTIAL EQUIVALENCE**

The Micro Vision Scleral Plugs are substantially equivalent to the predicate device by virtue of the following points:

6.1 No change in materials, basic components, or methods of manufacture between this device and the predicate; the raw materials used have been used in the medical industry on similar/identical products since pre-amendment, without any record of patient problems or adverse reactions:

6.2 No change in basic configuration or construction;

6.3 This device has been tested by an independent lab for biocompatibility and will be subjected to inspection and during/after manufacture and prior to release to the field.

6.4 The function and use of this product will be no different then that of the predicate device.

Predicate device is a Class II device, granted FDA marketing clearance under K854507, issued to Storz Instrument Co., St. Louis, Missouri. Substantial equivalency is being claimed to predicate device.

MicroVision, Inc. believes that based on the above comparison, the Micro Vision Scleral Plug is substantially equivalent to the above cited device, that any differences are minor and do not raise new issues of safety and effectiveness.

Signed: Leonard Rase  
Leonard Kastrilevich  
President  
MicroVision, Inc.

Dated: 12-23-02



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

FEB 26 2003

MicroVision, Inc.  
c/o Leonard Kastrilevich  
President  
34 Folly Mill Road, Suite 200  
P.O Box 1651  
Seabrook, NH 02874-1651

Re: K022186

Trade/Device Name: Scleral Plugs, 19 and 20 Gauge  
Regulatory Class: Unclassified  
Product Code: LXP  
Dated: December 20, 2002  
Received: January 6, 2003

Dear Mr. Kastrilevich:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in cursive script that reads "A. Ralph Rosenthal".

A. Ralph Rosenthal, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**510(k) Number (if known):** K022186

**Device Name:** Scleral Plug

**Indication For Use:**

The scleral plugs are intended to maintain patency of a previously made incision in the sclera of the eye, during ophthalmic surgical procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*James R. Vachner*

(Division Sign-Off)

Division of Ophthalmic Devices

510(k) K022186

Prescription Use ☒

(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐