K023291

Attachment 4 510(k) Summary

Category:	Comments	
Sponsor:	Boston Scientific Corporation	
	2710 Orchard Parkway	
	San Jose, CA 95134	
Correspondent:	Andrea Ļ. Ruth, RAC	
	Senior Associate, Regulatory Affairs	
	2710 Orchard Parkway	
	San Jose, CA 95134	
Contact Information:	E-mail: rutha@bsci.com	
	Phone: 408.895.3625	
	Fax: 408.895.2202	
Device Common Name	Electrosurgical Probe	
Device Proprietary Name	Cobra® Cooled Surgical Probe	
Device Classification	21 CFR § 878.4400, class II, product code GEI	
Predicate Device	Electrosurgical Probe	
Predicate Device Manufacturer(s)	Boston Scientific Corporation/EP	
	Technologies, Inc.	
Predicate Device Proprietary Name(s)	Cobra® Surgical Probe	
Predicate Device Classification Number	Class II	
Predicate Device Classification(s)	21 CFR § 878.4400, product code GEI	

Date Summary Was Prepared:

September 20, 2002

Description of the Device:

The Boston Scientific Corporation Surgical Probe is a sterile, single use electrosurgical device intended to be used to coagulate soft tissues. The surgical probe transmits radiofrequency energy from electrodes which are connected to an Electrosurgical unit (non-sterile; re-useable) through an Instrument Cable (sterile; re-useable).

Intended Use:

The Cobra® Cardiac Surgical Probe (Probe) is intended for the coagulation of cardiac tissue using radiofrequency (RF) energy during cardiac surgery. The Probe can be used during general surgery to coagulate soft tissues. The Probe may also be used to coagulate blood and soft tissue to produce hemostasis.

Comparison to Predicate Device:

	Predicate Device	Modified Device
510(k) Reference	K981981; K010956	Current Submission
Intended Use	Coagulation of tissue	Same
Device	Electrosurgical Probe	Same
Description		
Single Use?	Yes	Same
EO Sterilized?	Yes	Same
Manufacturer	Boston Scientific Corporation/EP	Same
	Technologies, Inc.	
Device	Class II, 21 CFR §878.4400, code	Same
Classification	GEI	

Summary of the Non-clinical Data:

Where appropriate, testing conformed to the requirements of 21 CFR Part 58 (Good Laboratory Practices (GLP)). Specifically, non-clinical tests conducted for the Device included Fluid Path Integrity, Bond Joint Tensile Strength, Bond Joint Torsional Strength, Distal Section Fatigue, Shaft to Handle Tensile Strength, Biocompatibility and, both acute and chronic, *In vivo* performance.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

OCT 17 2002

Boston Scientific Corporation Andrea L. Ruth Senior Associate, Regulatory Affairs 2710 Orchard Parkway San Jose, California 95134

Re: K023291

Trade/Device Name: Cobra® Cooled Cardiac Surgical Probe, Model 1596X

Regulation Number: 878.4400

Regulation Name: Electrosurgical probe

Regulatory Class: Class II

Product Code: GEI Dated: October 1, 2002 Received: October 2, 2002

Dear Ms. Ruth:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Andrea L. Ruth

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Célia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Attachment 2 Intended Use Statement

510(k) Number (if known): K023291

Device Name:

Cobra® Cooled Surgical Probe

Indication for Use:

The Cobra® Cardiac Surgical Probe (Probe) is intended for the coagulation of cardiac tissue using radiofrequency (RF) energy during cardiac surgery. The Probe can be used during general surgery to coagulate soft tissues. The Probe may also be used to coagulate blood and soft tissue to produce hemostasis.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)
Division of General, Restorative and Neurological Devices

510(k) Number (02329 / Over-the-Counter Use (Per 21 CFR §801.109)