

JAN 29 2003



K023659

353 Corporate Woods Parkway
Vernon Hills, IL 60061
Phone: 847-913-1113
Customer Service: 800-323-WOLF
www.richard-wolf.com

13.0 510(k) Summary of Safety and Effectiveness

Submitter:			Date of Preparation: October 31, 2002		
Company / Institution name: RICHARD WOLF MEDICAL INSTRUMENTS CORP.			FDA establishment registration number: 14 184 79		
Division name (if applicable): N.A.			Phone number (include area code): (847) 913 1113		
Street address: 353 Corporate Woods Parkway			FAX number (include area code): (847) 913 0924		
City: Vernon Hills	State/Province: Illinois	Country: USA	ZIP / Postal Code: IL 60061		
Contact name: Mr. Robert L. Casarsa					
Contact title: Quality Assurance Manager					
Product Information:					
Trade name: 1 CCD Endocam			Model number: 5520.xxx; 85520.xxx, 8934.551 85264.xxx		
Common name: Endoscopic Video Camera System			Classification name: Endoscope and/or Accessories		
Information on devices to which substantial equivalence is claimed:					
510(k) Number	Trade or proprietary or model name		Manufacturer		
1 K982965	1 1-CCD endocam 5502 and 3-CCD endocam 5507 with CF camera heads		1 Richard Wolf		
2 K983279	2 1-CCD multi-endocam 5502 with electronic CCD endoscope		2 Richard Wolf		
3 K950502	3 1-CCD endocam 5501		3 Richard Wolf		
4 K002328	4 SIOS-Interface for 3CCD Endocam 5507.752		4 Richard Wolf		
5 K964173	5 C-mount objective lenses, autoclavable		5 Richard Wolf		

1.0 Description

The 1-CCD Endocam 5520 device is a further development of previous 1-CCD Endocams with additional features. Various Camera Heads with and without integrated lenses and CCD endoscopes are available.



2.0 Intended Use

The 1 CCD ENDOCAM 5520 is designed for video endoscopy and video microscopy and can be used for diagnostic and therapeutic applications. The 1 CCD ENDOCAM 5520 is designed for connection to various camera heads and CCD endoscopes.

In conjunction with video recorders / video printers and other video equipment it can be used for recording and storing video images.

The camera heads and the CCD endoscope are used in connection with the 1 CCD ENDOCAM 5520 controller for diagnostic and therapeutic applications. Free rotation and self-alignment of the urological camera heads results in an endoscope image, which is always correctly positioned.

3.0 Technological Characteristics

The 1-CCD Endocam 5520 System has additional features such as digital zoom, mirror imaging, anti-Moiré, and interfaces for digital video output (IEEE 1394), RIWONET control, keyboard and remote control.

4.0 Substantial Equivalence

The submitted devices pose the same type of questions about safety or effectiveness as the compared devices and the new technological characteristics have not diminished safety or effectiveness. The submitted devices are substantially equivalent to existing 510(k)-devices sold by Richard Wolf.

5.0 Performance Data

The Endocam systems 5520 was tested for conformity with the specified standards UL2601-1, IEC601-2-18, CSA22.2No.601.1.


6.0 Clinical Tests

No clinical tests performed.

7.0 Conclusions Drawn

These devices are designed and tested to guarantee the safety and effectiveness, when used according to the instructions manual.

By:



Robert L. Casarsa
Quality Assurance Manager

Date:





DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Richard Wolf Medical Instruments Corporation
Robert L. Casarsa
Quality Assurance Manager
353 Corporate Woods Parkway
Vernon Hills, Illinois 60061

JAN 29 2003

Re: K023659

Trade/Device Name: 1 – CCD Endocam 5520 System
Regulation Number: 876.1500
Regulation Name: Endoscope and accessories
Regulatory Class: Class II
Product Code: KOG
Dated: October 31, 2002
Received: October 31, 2002

Dear Mr. Casarsa:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally

marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in cursive script that reads "Miriam C. Provost".

for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

5.0 INDICATIONS FOR USE

510(k) Number (if known): — K023659

Device Name: 1-CCD Endocam 5520 System

Intended use: The 1 CCD ENDOCAM 5520 is designed for video endoscopy and video microscopy and can be used for diagnostic and therapeutic applications.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

Miriam C. Provost

(Division Sign-Off)

Division of General, Restorative
and Neurological Devices

Prescription Use ☒
Per 21 CFR 801.109

510(k) Number K023659

OR

Over-The Counter ☐