

## 510(k) SUMMARY

FEB 24 2004

### AcuNav 90/10 Diagnostic Ultrasound Catheter

This summary of safety and effectiveness is provided as part of this Premarket Notification in compliance with the Safe Medical Device Act of 1990 revisions to 21 CFR, Part 807.92, Content and Format of a 510(k) Summary.

**1. Submitted by:**

Siemens Medical Solutions USA, Inc., Ultrasound Division  
1230 Shorebird Way  
Mountain View, CA 94039

**Contact Person:**

Iskra Mraković  
Manager, Regulatory Affairs

Phone: (650) 694-5004

Fax: (650) 961-6168

**2. Proprietary Name:**

AcuNav 90/10 Diagnostic Ultrasound Catheter

**Common/Usual Name:**

Intracardiac/Intravascular Ultrasound Catheter.

**Classification Name:**

Regulatory Class: II

Review Category: Tier II

	FR Number	Product Code
Diagnostic Intravascular Catheter	870.1200	74-DQO
Diagnostic Ultrasound Transducer	892.1570	90-ITX

**3. Predicate Device:**

ACUSON AcuNav® 90/10 Diagnostic Ultrasound Catheter

**4. Device Description:**

As described in the previously cleared 510(k) submissions, (K992631 on December 15, 1999, and K010950 on June 27, 2001), the AcuNav catheter is comprised of a single-use, disposable ultrasonic phased-array imaging transducer and a catheter which is 10 Fr (3.33 mm) in diameter and 90 cm in insertable length. The device is capable of imaging at multiple frequencies and obtaining blood flow data in multiple ultrasound modes. The distal portion of the catheter can be deflected in four directions in two orthogonal planes: left-right (in a plane perpendicular to the image plane) and anterior-posterior (in a

000453

plane coincident with the image plane). The range of deflection is 160° in each direction. The AcuNav catheter is comprised of three major components: (1) the catheter itself; (2) the steering mechanism; and (3) the reusable system cable.

The AcuNav 90/10 Diagnostic Ultrasound Catheter has been designed to meet the following product safety standards [as required by 21 CFR § 807.87(j)].

- UL 2601-1, Safety Requirements for Medical Equipment
- CSA C22.2 No. 601-1, Safety Requirements for Medical Equipment
- CISPR 11, Class A
- IEC 601-1-1
- IEC 601-1-2
- ISO 10993 Biocompatibility

**5. Intended Use:**

The AcuNav catheter is intended for use in the visualization of vascular anatomy and physiology, cardiac and great vessel anatomy and physiology, visualization of other devices in the heart or vasculature, as well as measurements of blood flow, cardiac and vascular anatomic dimensions.

**6. Technological Comparison to Predicate Device:**

The AcuNav 90/10 Diagnostic Ultrasound Catheter is substantially equivalent to the ACUSON AcuNav® 90/10 Diagnostic Ultrasound Catheter (K992631, 11/8/2001).

Both catheters are ultrasound-tipped catheter devices used directly within the vasculature and/or the heart for intravascular or intracardiac ultrasound imaging. Devices are specifically indicated for use in visualization of vascular anatomy, cardiac and great vessel anatomy and physiology, or other devices in the heart, as well as measurement of blood flow.

**End of 510(k) Summary.**

000454



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

FEB 24 2004

Siemens Medical Solutions USA, Inc.  
c/o Ms. Iskra Mraković  
Regulatory Affairs Manager  
1230 Shorebird Way  
P.O. Box 7393  
Mountain View CA 94043

Re: K033650

Trade Name: Acunav 90/10 Diagnostic Ultrasound Catheter  
Regulation Number: 21 CFR 870.1200 and 892.1570  
Regulation Name: Diagnostic Intravascular Catheter  
Diagnostic Ultrasound Transducer  
Regulatory Class: Class II (two)  
Product Code: DQO and ITX  
Dated: November 19, 2003  
Received: November 20, 2003

Dear Ms. Mraković:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4648. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



 Bram D. Zuckerman, M.D.  
Division Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Diagnostic Ultrasound Indications for Use Form**  
**Aspen Diagnostic Ultrasound System**

510(k) Number: K033650

Device Name: AcuNav 90/10 (IC10V5) Diagnostic Ultrasound Catheter

Indications for Use: The AcuNav™ Diagnostic Ultrasound Catheter is intended for intra-cardiac and intra-luminal visualization of cardiac and great vessel anatomy and physiology as well as visualization of other devices in the heart.

Ultrasound System: Aspen

Transducer: AcuNav

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)	Other: Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Vascular)										
Intraoperative (Neurological)										
Pediatric										
Small Organ (Specify)**										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		P*	P
Trans-esophageal										
Transrectal										
Transvaginal										
Transurethral										
Intra-Luminal		P	P	P	P	P	P		P*	P
Peripheral Vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Intra-Cardiac)		P	P	P	P	P	P		P*	P

P=Previously cleared by the FDA under premarket notification K992631.

**Additional Comments:**

\*Combinations include: B+M, B+PWD, B+CWD, B+Color Doppler, B+M+ Color Doppler, B+PWD+Color Doppler, B+CWD+Color Doppler, B+Power Doppler, B+M+Power Doppler, B+PWD+Power Doppler, B+CWD+Power Doppler

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

*Dan R. Kochner*  
 (Division Sign-Off)  
 Division of Cardiovascular Devices

510(k) Number K033650

**Diagnostic Ultrasound Indications for Use Form**  
**Cypress Diagnostic Ultrasound System**

510(k) Number: K033650

Device Name: AcuNav 90/10 (IC10V5) Diagnostic Ultrasound Catheter

Indications for Use: The AcuNav™ Diagnostic Ultrasound Catheter is intended for intra-cardiac and intra-luminal visualization of cardiac and great vessel anatomy and physiology as well as visualization of other devices in the heart.

Ultrasound System: Cypress

Transducer: AcuNav

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)	Other: Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Vascular)										
Intraoperative (Neurological)										
Pediatric										
Small Organ (Specify)**										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P			P*	P
Trans-esophageal										
Transrectal										
Transvaginal										
Transurethral										
Intra-Luminal		P	P	P	P	P			P*	P
Peripheral Vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Intra-Cardiac)		P	P	P	P	P			P*	P

P=Previously cleared by the FDA under premarket notification K992631.

**Additional Comments:**

\*The Cypress Ultrasound System does not provide combined modes where more than one scanning mode is active simultaneously.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Diana R. Vachner  
 (Division Sign-Off)  
 Division of Cardiovascular Devices

510(k) Number K033650

**Diagnostic Ultrasound Indications for Use Form**  
**Sequoia Diagnostic Ultrasound System**

510(k) Number:

K 633650

Device Name:

AcuNav 90/10 (IC10V5) Diagnostic Ultrasound Catheter

Indications for Use:

The AcuNav™ Diagnostic Ultrasound Catheter is intended for intra-cardiac and intra-luminal visualization of cardiac and great vessel anatomy and physiology as well as visualization of other devices in the heart.

Ultrasound System:

Sequoia

Transducer:

AcuNav

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)	Other: Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Vascular)										
Intraoperative (Neurological)										
Pediatric										
Small Organ (Specify)**										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		P*	P
Trans-esophageal										
Transrectal										
Transvaginal										
Transurethral										
Intra-Luminal		P	P	P	P	P	P		P*	P
Peripheral Vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Intra-Cardiac)		P	P	P	P	P	P		P*	P

P=Previously cleared by the FDA under premarket notification K992631.

**Additional Comments:**

\*Combinations include: B+M, B+PWD, B+CWD, B+Color Doppler, B+M+ Color Doppler, B+PWD+Color Doppler, B+CWD+Color Doppler, B+Power Doppler, B+M+Power Doppler, B+PWD+Power Doppler, B+CWD+Power Doppler

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Dan R. Lockner  
 (Division Sign-Off)  
 Division of Cardiovascular Devices

510(k) Number K033650