

cobas[®] HPV Test

Microbiology Devices Panel

P100020/S008

March 12, 2014



cobas[®] HPV Test

Microbiology Devices Panel

P100020/S008

March 12, 2014



CI-1

Introduction

Christoph Majewski, PhD

Lifecycle Leader, HPV and Microbiology
Roche Molecular Systems



CI-2

Agenda

Introduction	Christoph Majewski, PhD Lifecycle Leader, HPV & Microbiology, RMS
Clinical Need for HPV as Primary Screening Test	Thomas C. Wright, Jr., MD Professor Emeritus, Columbia University
ATHENA Study Objectives and Statistics	Abha Sharma, PhD Director Biostatistics, RMS
Data from ATHENA Supporting cobas® HPV Test for Primary Screening	Catherine Behrens, MD, PhD, FACOG Director, Clinical Research, RMS
Clinical Implications and Benefit-Risk	Thomas C. Wright, Jr., MD Professor Emeritus, Columbia University
Summary	Christoph Majewski, PhD Lifecycle Leader, HPV & Microbiology, RMS

CI-3

Pap Test Revolutionized Cervical Cancer Screening, but Unmet Need Still Exists

- Screening significantly reduced cervical cancer incidence
- In 2014: ~12,360 cases and ~4020 deaths in the US
- Cytology and cotesting present standard of care
- Current solutions have limitations and are highly complex
- Primary HPV screening can address some limitations

CI-4

HPV Primary Screening

Proposed New Claim for cobas® HPV Test




- Use as first-line screening test in women 25 and older to detect high-risk HPV, including HPV genotypes 16 and 18
- **Candidate Algorithm: HPV Primary Screening with HPV 16/18 and Cytology Triage**
 - Negative for HPV: Follow-up by physician's judgment
 - Positive for HPV genotype 16/18: Colposcopy
 - Positive for any of 12 high-risk HPV types: Reflex to a cytology exam to determine need for colposcopy

CI-5

cobas® HPV Test Technology Overview

Approved in 2011 for ASC-US Triage and Cotesting




4 channel design allows reporting of pooled hrHPV, as well as simultaneously providing HPV 16/18 specific identification from a single test tube

Wright TC Jr, et al. *Am J Obstet Gynecol* 2012; 206:46.e1-.e11

CI-6

Importance of Genotype

Prevalence of HPV16 and HPV18 in Cervical Cancer



Approximately 70% of all cervical cancers are associated with HPV genotype 16 or 18

¹Note distribution is for single infections only
Bosch FX, et al. *J Natl Cancer Inst Monogr* 2003; 31:3–13.

CI-7

cobas® HPV Test

Regulatory Background



ATHENA Baseline

ATHENA 3 year Follow-up

Approved

1. PMA Submission

PMA Approval

- ASC-US triage
- Adjunct testing ♀ ≥30 years
- Genotyping HPV 16 and 18

2008

2009

2010

2011

2012

2013

2014

Submitted
Additional Claim

FDA Panel

2. PMA Supplement

- HPV 1° screening

CI-8

Clinical Need for HPV as Primary Screening Test


Thomas C. Wright, Jr., MD

Professor Emeritus
Columbia University

CC-9

Cytology-based Screening

Most Widely Utilized Globally and in the U.S.



This is our Comparator Algorithm

Cervical cytology at 3 year intervals is considered an acceptable approach by both the USPSTF and American Cancer Society

CC-10

Limitations of Cervical Cytology

- Interpretation is quite subjective which results in considerable intra- and inter-laboratory variation
- Relatively low sensitivity for the detection of high-grade cervical cancer precursors
- Identifies individuals with cancer precursors but not women at risk of developing cancer precursors

CC-11

Reproducibility of Cervical Cytology

Re-read of 4948 Liquid-based Cytology Slides

		QC Reviewer's Diagnosis			
		NILM	ASC-US	LSIL	≥HSIL
Original Diagnosis	NILM	78%	19%	3%	<1%
	ASC-US	39%	43%	17%	2%
	LSIL	4%	22%	68%	6%
	≥HSIL	3%	23%	27%	47%

Stoler and Schiffman JAMA, 2001.

CC-12

Variability of Cervical Cytology **ATHENA Results**

	Lab A	Lab B	Lab C	Lab D
Number	12,294	4218	16,979	12,442
Median Age	40.9	37.9	39.3	40.1
≥ASC-US	3.8%	5.2%	8.1%	9.9%
Sensitivity of Cytology*	42.0	51.0	60.5	73.0
Sensitivity of cobas®*	90.1	88.2	88.4	88.9

*Note: for ≥CIN2
 Wright et al. *Int. J. Cancer*, 2013, Oct 7 epub
 Data not reviewed by the FDA

CC-13

Limitations of Cervical Cytology

- Interpretation is quite subjective which results in considerable intra- and inter-laboratory variation
- Relatively low sensitivity for the detection of high-grade cervical cancer precursors
- Identifies individuals with cancer precursors but not women at risk of developing cancer precursors


CC-14

Performance of Cervical Cytology Sensitivity for \geq CIN2

Author	Year	Number	Method	Sensitivity	95% CI
Petry	2003	8466	Conv	44%	(30-58)
Coste	2003	3080	Conv	65%	(50-80)
Taylor	2005	3114	LBC	71%	(58-81)
Ronco	2006	22,760	LBC	74%	(62-84)
Mayrand	2007	10,153	Conv	57%	(34-78)

CC-15

Sensitivity of HPV to Detect \geq CIN2 is Higher than Cytology in USPSTF Review



HPV testing used an HPV assay other than cobas® HPV Test
Studies performed in developed countries in women 30 years and older
Whitlock EP, et al. *Ann Intern Med*.2011; 155:687-697, W214-5.


CC-16

Limitations of Cervical Cytology

- Interpretation is quite subjective which results in considerable intra- and inter-laboratory variation
- Relatively low sensitivity for the detection of high-grade cervical cancer precursors
- Identifies individuals with cancer precursors but not women at risk of developing cancer precursors

CC-17

10 Year Cumulative Risk of \geq CIN3 in Women \geq 30 Years with NILM Cytology Kaiser Portland OR




HPV testing used an HPV assay other than cobas® HPV Test
Kahn et al. *JNCI*, 2005; 97; 1072

CC-18

Cotesting with Cytology and HPV

Used in the U.S. But Not the Predominant Method



CC-19


Large Number of Cytology Categories

NILM (negative)	AGC-EC
ASC-US	AGC-EM
ASC-H	AIS
LSIL	Other
HSIL	
Satisfactory	
Sat but limited by...	Unsatisfactory

2013 ASCCP Management Guidelines have 12 different algorithms just for cytology results

CC-20


Management of Women with Low-grade Squamous Intraepithelial Lesions (LSIL)*



*Management options may vary if the women are pregnant or ages 21 to 24 years
Massad et al. JLGTD, 2013

CC-21

Primary HPV Screening – Candidate HPV with 16/18 Genotyping and Reflex Cytology



hrHPV=high risk HPV


CC-22

At What Age Should We Initiate Primary HPV Screening

- Current U.S. screening guidelines do not recommend cotesting for women 25-29 years of age
- Transient HPV infections are common in this age group and guideline makers did not want unnecessary follow-up examinations and colposcopy
- There is a high burden of CIN3 in women 25-29 years and cytology performs poorly in young as shown by UK screening audits
- In 2013 Kaiser Permanente, N. California reviewed their registry data and decided to begin cotesting at age 25 years

CC-23


Invasive Cervical Cancer in the U.S. SEER Tumor Registry (1975-2010)



http://seer.cancer.gov/csr/1975_2010/

CC-24


HPV by Age Group ATHENA



Data not reviewed by FDA
Wright et al. *Am J Obst Gynecol*, 2011.

CC-25


HPV 16/18 by Age Group ATHENA



Data not reviewed by FDA
Wright et al. *Am J Obst Gynecol*, 2011.

CC-26


HPV 16/18 vs \geq ASC-US ATHENA



Data not reviewed by FDA
Wright et al. *Am J Obst Gynecol*, 2011.

CC-27


\geq CIN3 by Age Group ATHENA



Data not reviewed by FDA
Wright et al. *Am J Obst Gynecol*, 2011.

CC-28

≥CIN3 by Age Group ATHENA




More ≥CIN3 disease in women 25-29 years than in women ≥40 years

Data not reviewed by FDA

CC-29

Proportion of Women with ≥CIN3 Who Have Negative Cytology (NILM) ATHENA



Data not reviewed by FDA
Percentages shown are for hrHPV+ women with ≥CIN3, N=252
Huh W, et al. 27th International Papillomavirus Conference, Berlin, Germany, September 17–22, 2011, OP-229.

CC-30

Need for Primary HPV Screening Starting at 25 Years of Age

- Cervical cytology appears to have reached the point where it alone is unable to reduce cervical cancer rates further
- Current management algorithms are extremely complicated and this confusion is potentially resulting in poor clinical care
- Cytology is not a good solution for identifying the majority of high-grade disease in women 25-29 years of age

CC-31

ATHENA Study Objectives and Statistics

Abha Sharma, PhD

Director, Biostatistics
Roche Molecular Systems



CT-32

Cervical Cancer Screening Evaluation

Study Design Requirements

- Cervical cancer screening study must:
 - Be a large cross-sectional cohort with sufficient follow-up for safety
 - Have sufficient cases of \geq CIN2
 - Be representative of target population
 - Adjust for verification bias

CT-33

Study Objectives and Screening Algorithms

- **Study Objective:** Compare the performance of primary HPV screening algorithm vs algorithm using cytology as first line of screening
- **Candidate Algorithm:** Primary HPV with 16/18 genotyping and reflex cytology
- **Comparator Algorithm:** Cytology alone (\geq ASC-US to colposcopy)
- **Additional Comparator:** ATRI NM \geq 30 GT
 - ASC-US Triage in women 25-29
 - Cotesting with genotyping in women \geq 30

CT-34

Screening Algorithms

Performance Comparison



- **Endpoint:** $\geq \text{CIN2}$
- **Secondary:** $\geq \text{CIN3}$
 - Results presented for $\geq \text{CIN3}$ (better surrogate for cancer)
- **Performance Metrics:**
 - Sensitivity/specificity
 - PPV/NPV
 - Likelihood ratios PLR/NLR
- **Acceptance Criteria:** PLR/NLR
 - Higher PLR and Lower NLR indicate better performance
 - Confidence interval for the difference should exclude “0”

CT-35

Performance Estimates of Screening Algorithms

PLR, NLR




Parameter	Description	Interpretation
PLR >1	Se/(1-Sp)	1. How many times more likely women with $\geq \text{CIN3}$ are to have a positive result than women with <CIN3 2. Post-test odds = PLR \times pre-test odds
NLR <1	(1-Se)/Sp	1. How many times less likely women with $\geq \text{CIN3}$ are to have a negative result than women with <CIN3 2. Post-test odds = NLR \times pre-test odds

Note: PLR and NLR do not depend on the prevalence of disease

CT-36

Graphical Comparison of Algorithms



Kondratovich M. 2008; Biggerstaff BJ. 2000

CT-37

Statistical Methods for Baseline and Follow-up Verification Bias

- **Verification Bias:** Occurs when test results determine who is “verified” for disease status
 - Women with positive test results for HPV/cytology
 - A random subset of HPV/cytology negative patients randomized to colposcopy
- **Missing at Random (MAR)** assumption
- **VBA Calculations*:** Adjust performance statistics based on observed disease in verified group, using probability of being verified

*MS Pepe. 2002; XH Zhou et al. 2003; Begg and Greenes.1983

CT-38



ATHENA Study Objectives and Statistical Analysis Conclusions

- Study Objective – Compare:
 - HPV as primary screening (Candidate)
 - Cytology as primary screening (Comparator)
- Verification Bias Adjusted (VBA) Statistics
- Acceptance Criteria
 - NLR for Candidate < NLR for Comparator
 - PLR for Candidate > PLR for Comparator
 - With additional information: Se/Sp and PPV/NPV
- Safety of negative HPV test result established by cumulative risk from 3 year follow-up

CT-39

Data from ATHENA Supporting cobas® HPV Test for Primary Screening

Catherine Behrens, MD, PhD, FACOG

Director, Clinical Research, RMS
Roche Molecular Systems



CA-40

The Challenge for ATHENA



Can We Improve Screening Methodology and Add Medical Value by Increasing Detection of Precancer?

- ATHENA is the largest prospective cervical cancer screening study in the U.S.
- Enrolled 47,208 women ≥ 21 years undergoing routine cervical cancer screening in the U.S.
 - 61 clinical sites in 23 states and 4 clinical laboratories
- Served as the registrational study for the cobas[®] HPV Test with 16/18 genotyping and FDA approval received in 2011 for:
 - ASC-US management
 - Cotesting with cytology for screening

CA-41

ATHENA Trial Design



- Specifically designed to demonstrate the performance of HPV testing in cervical screening in the U.S.
 - The ATHENA population was representative of a U.S. screening population in demographics, cytology distribution, and HPV prevalence
 - Both cytology and HPV testing (with genotyping) performed on all women
 - Rigorous disease ascertainment was achieved
 - All women who screened positive for either test (both Pap+ and HPV+) were taken to colposcopy
 - Histology determined by consensus of expert pathologists

CA-42

Demographics of ATHENA Trial Representative of the US Population (≥25)



Characteristics	Evaluable Subjects N=40,944 % (n)	U.S. Census Figures 2012 ¹ %
Age (years)	41	
Race		
White	83.4 (34,156)	79.8
American Indian or Alaskan Native	0.6 (226)	1.0
Black or African American	13.7 (5602)	12.4
Asian	1.6 (639)	5.2
Native Hawaiian/Other Pacific Islander	0.2 (98)	0.2
Any combination ¹	0.5 (220)	1.4
Missing	<0.1 (3)	
Ethnicity		
Hispanic or Latino	18.0 (7370)	12.9

Note: Any combination refers to subjects who selected more than one race

¹Based on Annual Estimates of the Resident Population by Sex, Age, Race, and Hispanic Origin for the United States and States: April 1, 2010 to July 1, 2012 (<http://factfinder2.census.gov/faces/tableservices/jsf/pages/productview.xhtml?src=bkmk>)

CA-43

Cytology Results and Prevalence of hrHPV in ATHENA are Representative of a U.S. Screening Population



Pap Test Result	Eligible Subjects ≥25 Years N=40,944 % (n)	CAP 2010 ¹ %
NILM	93.5 (38,397)	91.5
ASC-US	4.0 (1632)	4.8
>ASC-US	2.4 (986)	3.6
LSIL	1.9	2.8
ASC-H	0.1	0.3
HSIL	0.3	0.4
Squamous Cell Carcinoma		
AGC ^{a, b}	<0.1	0.1

¹CAP Cytopathology Checklist (all ages, not adjusted for ≥25 years)

^aAGC (Atypical Glandular Cells) includes: AGC - Endocervical, AGC - Endometrial, and AGC - Not Otherwise Specified

^bAGC, Favor Neoplastic includes: AGC - Endocervical - Favor Neoplastic and AGC - Favor Neoplastic

CA-44

HPV Prevalence in ATHENA is Representative of a U.S. Population




Age Groups (Years)	HPV+		HPV16+		HPV18+	
	ATHENA %	NMHPVPR* %	ATHENA %	NMHPVPR* %	ATHENA %	NMHPVPR* %
Overall Evaluable Primary Screening Subjects	10.5	14.2	2.1	3.1	0.8	0.9
25-29	21.1	21.8	5.3	5.2	1.6	1.4
30-49	9.4	11.5	1.7	2.2	0.7	0.7
≥50	6.0	6.9	0.8	1.3	0.4	0.5

*New Mexico HPV Pap Registry; assumed that carcinogenic HPV+ (HPV 16, 18, 31, 33, 39, 45, 52, 56, 58, 59, and 68) in NMHPVPR was equivalent to HPV+ (HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 66, and 68) in ATHENA
Wheeler et al. *International J Cancer*, 2013


CA-45

ATHENA Patient Flow Cross-sectional Phase for Primary HPV Effectiveness



CA-46

ATHENA Patient Flow 3 Year Follow-up Phase for Safety Evaluation




CA-47

Management of Cytology and Histology




Bethesda Classification	Management
Negative (NILM)	Routine follow-up 3 years
ASC-US	HPV testing
ASC-US/HPV-	Routine follow-up 3 years
ASC-US/HPV+	
LSIL	Colposcopy
HSIL	
ASC-H	
AGC	Colposcopy, endocervical curettage (ECC), endometrial biopsy
ACIS	
Squamous/adenocarcinoma	Colposcopy



CA-48

3-Year Cumulative Risks for \geq CIN3 Primary Screening Population (≥ 25 Years)



VBA estimates

CA-49

Primary Screening Algorithms

Effectiveness: Baseline Data




CA-50

Choosing the Optimal Screening Strategy

- Cervical cancer screening strategies should maximize disease detection (sensitivity), while minimizing the “harms”
- Colposcopy
 - Anxiety, discomfort
- Additional harms of screening
 - False negative results
 - Precancer missed by cytology
 - False positive results
 - Over-screening, over-management of lesions likely to regress
 - Treatment
 - Procedures (LEEP, conization) that may lead to longer-term complications related to pregnancy

CA-51

ATHENA Patient Flow Cross-sectional Phase for Primary HPV Effectiveness



CA-52

ATHENA Data to be Presented

- Comparisons of the performance among 3 screening algorithms will be presented
 - Comparator (cytology alone) vs Candidate (HPV with 16/18 genotyping and reflex to cytology)
 - Additional Comparator: Cotesting
 - ATRI NM ≥ 30 GT (ASC-US Triage for women 25-29 years and cotesting for women ≥ 30 years): Current strategy supported by 2012 guidelines


CA-53

ATHENA Data to be Presented

- To demonstrate effectiveness
 - Sensitivity, specificity, predictive values (PPV and NPV) and likelihood ratios (PLR and NLR)
- To demonstrate safety
 - Negative predictive value (NPV)
 - 3 Year cumulative risks (CIRs) for a negative HPV result vs negative cytology result at Baseline was calculated
- Only data using $\geq \text{CIN}3$ endpoint will be presented since $\geq \text{CIN}3$ is considered a better surrogate for cancer when assessing screening strategies


CA-54

Comparator Screening Algorithm Cytology Alone



CA-55

Candidate Screening Algorithm HPV with 16/18 Genotyping and Reflex Cytology



hrHPV=high risk HPV

CA-56

Comparison of the Candidate vs Comparator (Cytology Alone) to Detect \geq CIN3



Description	\geq CIN3	
	Relative Sensitivity ¹ %	Relative Specificity ¹ %
Comparator	1.00	1.00
Candidate	1.37*	1.02*

Using primary screening with the Candidate algorithm increases the sensitivity of HPV testing by 37% over cytology and raises the specificity to be at least equal to cytology

¹Calculated as VBA sensitivity or specificity of Candidate/VBA sensitivity or specificity of Comparator

*Difference of VBA parameters statistically significant

CA-57

Comparison of Predictive Values and Likelihood Ratios of the Candidate and the Comparator to Detect \geq CIN3



Description	\geq CIN3			
	PPV ¹ %	NPV ¹ %	PLR ¹	NLR ¹
Comparator	6.47	99.41	7.06	0.61
Candidate	12.25	99.58	14.24	0.44
Difference	5.78* (4.72, 6.94)	0.17* (0.12, 0.23)	7.18* (5.34, 9.4)	-0.17* (-0.24, -0.2)

The Candidate nearly doubles the PPV and PLR for detection of disease when compared to Cytology

The Candidate NPV and NLR also improve, indicating a superior measure of safety over Cytology

¹ Verification bias adjusted


*Difference of VBA parameters statistically significant

CA-58

ATRI NM ≥ 30 GT (Cotesting Hybrid) **ASC-US Triage for Women 25-29 Years and Cotesting** **for Women ≥ 30 Years**



Women 25-29 years: ASC-US Triage




cobas® HPV Test, high risk HPV DNA test

CA-59

ATRI NM ≥ 30 GT (Cotesting Hybrid) **ASC-US Triage for Women 25-29 Years and Cotesting** **for Women ≥ 30 Years**



Women ≥ 30 years: Cotesting



CA-60

Comparison of the Performance of the Candidate vs Cotesting Hybrid to Detect \geq CIN3



Description	\geq CIN3	
	Relative Sensitivity ¹ %	Relative Specificity ¹ %
Comparator	1.00	1.00
Candidate	1.37	1.02
ATRI NM \geq 30 GT	1.25*	1.02

The sensitivity of Cotesting Hybrid \geq 30 years decreases due to women 25-29 years having cytology screening only

¹Calculated as VBA sensitivity or specificity of ATRI NM \geq 30 GT/VBA sensitivity or specificity of Comparator

*Difference of VBA sensitivity statistically significant

CA-61

Comparison of Predictive Values and Likelihood Ratios of the Candidate and the Cotesting Hybrid



Description	\geq CIN3			
	PPV ¹ %	NPV ¹ %	PLR ¹	NLR ¹
ATRI NM \geq 30 GT	11.04	99.52	12.66	0.49
Candidate	12.25	99.58	14.24	0.44
Difference	1.21* (0.46, 1.96)	0.06* (0.01, 0.09)	1.58* (0.62, 2.71)	-0.05* (-0.10, -0.01)


The PPV and NPV of the Candidate are superior, indicating significantly improved effectiveness and safety over the additional Comparator

¹ Verification bias adjusted

*Difference of VBA parameters statistically significant


CA-62

Comparison of Candidate and Additional Comparator with Comparator for \geq CIN3 Endpoint



CA-63

Comparison of Candidate and Additional Comparators with Comparator for \geq CIN3 Endpoint



- Comparator – Cytology Alone
- ▲ Candidate – HPV 16/18 Genotyping Test with Reflex to Pap Test
- ◆ ATRI NM \geq 30 – ASC-US Triage 25-29 and ASC-US Triage NILM 16+/18+ 30+

CA-64

Clinical Implications for Various Algorithms



CA-65

Projected Measures of Clinical Management for Disease (\geq CIN3)



Algorithm	Description	No. of Screening Tests	No. of Screening Tests Per \geq CIN3	No. of Colposcopies	No. of Colposcopies Per \geq CIN3	No. of \geq CIN3 Cases Detected ¹
Comparator	Cytology alone	40,944	239.4	2618	15.3	171
Candidate	HPV with 16/18 and reflex to cytology	44,057	189.9	1890	8.1	232
ATRI NM \geq 30 GT	ASC-US Triage 25-29 and ASC-US Triage NILM 16+/18+ 30+yr	75,574	358.2	1916	9.1	211

CA-66

Summary of ATHENA Data in Support of Effectiveness of Screening



- When compared to cytology or cotesting:
 - Candidate demonstrates the best sensitivity for detection of \geq CIN3
 - The specificity of the Candidate is at least equal to cytology when 16/18 genotyping and reflex cytology is added to HPV as the primary screen
 - The Candidate PPV and PLR are 2x that of cytology and significantly greater than Cotesting Hybrid
 - The NPV and NLR of the Candidate are improved over both cytology and Cotesting Hybrid
 - The Candidate demonstrates a better balance of clinical resource management than either cytology or cotesting

CA-67


Primary Screening Algorithms

Safety: Longitudinal Data



CA-68


Evaluating the 90% of Women Who Screen HPV(-)



The lower risk of disease of a negative hrHPV at Baseline confirms the safety of a negative hrHPV result over 3 years

CA-69

The Benefit of Cotesting Over HPV in Reducing the 3 Year Risk of ≥CIN3 is Minimal



A negative Pap result added to a negative hrHPV result at Baseline adds little benefit and increases the colposcopy rate from 4.6% to 5.4%

* Statistically significant difference

CA-70

Sensitivity of cobas® HPV Test vs Cytology to Detect Invasive Cancer



	cobas® HPV Test+	Cytology+
ATHENA n=8	8	7
UNMHPVP Registry n=18 ²	17 ¹	16
Total	25	23

- Sensitivity of cobas® HPV Test: $25/26 = 96.2\%$
- Sensitivity of cytology: $23/25 = 92\%$

¹ 1 case determined to be a poorly differentiated adenocarcinoma with origin uncertain, endocervical vs endometrial
² 1 case of endometrial cancer was found to have been sent in error after HPV testing was performed at RMS and 1 case was determined to be cobas® HPV "invalid" due to clotting of sample; these cases were excluded from the analysis

CA-71

Conclusions



- HPV-based strategies for primary screening are more sensitive for detection of high-grade disease than cytology-based strategies
 - The specificity of HPV-based strategies is increased by the addition of 16/18 genotyping and reflex testing to cytology
- Effectiveness of the Candidate algorithm is demonstrated by its superior performance compared to strategies supported by the current guidelines: Cytology alone and the Cotesting Hybrid
 - For the detection of precancer, the Candidate provides the optimal balance of benefits and harms
- Safety of the cobas® HPV Test as a primary screening test is confirmed by demonstrating that a negative HPV result at Baseline predicts a lower risk of \geq CIN3 at 3 years than a negative Pap result at Baseline

CA-72

Clinical Implications and Benefit-Risk

Thomas C. Wright, Jr., MD

Professor Emeritus
Columbia University

CR-73

Candidate Algorithm

Discussion of Clinical Implications

- Screening for other cancers and STIs will **NOT** be adversely impacted if we use HPV alone for screening
- Shifting to primary HPV screening will **NOT** put women at increased risk for invasive cervical cancer or high-grade precursor lesions

CR-74

Cervical Cytology to Screen for Other Gynecological Cancers

- The sensitivity of cervical cytology for endometrial and ovarian cancer is low (10-30% depending on study)^{1,2}
- Positive cervical cytology is associated with high-stage disease or cervical involvement^{2,3}
- Therefore, detecting endometrial or ovarian cancer by cervical cytology does not improve survival rates^{2,3}
- Cervical cytology is not considered appropriate for screening for other cancers by **USPSTF, ACS, ACOG**

¹Mitchell H. et al. *Int. J. Gyn. Pathol.* 1993; 12:34

²Nawanodi O. et al. *Arch. Gynecol. Obstet.* 2008; 277:171

³Roelofssoen T. et al. *Int. J. Gyn. Pathol.* 2013; 32:390

CR-75

Detection of Sexually Transmitted Infections with Cervical Cytology

- A number of organisms such as *T vaginalis*, *candida*, shift in flora suggestive of bacterial vaginosis, *actinomyces*, and Herpes Simplex can be identified on cervical cytology
- Sensitivity of cervical cytology is considered to be too low to be a useful screening test for infectious organisms – **CDC, ACOG, and other ID societies**
- There are other tests widely available to clinicians that are **BOTH** more sensitive and more specific for infectious organisms

CR-76

Candidate Algorithm

Discussion of Clinical Implications

- Screening for other cancers and STIs will **NOT** be adversely impacted if we use HPV alone for screening
- Shifting to primary HPV screening will **NOT** put women at increased risk for invasive cervical cancer or high-grade precursor lesions

CR-77

Sensitivity of Cytology and HPV for CIN3 and Cervical Cancer

- **NO screening test will detect ALL CIN3 or cancers**
 - Occasional sampling issues
 - Rare types of cervical cancer (mesonephric carcinoma, clear cell, etc) may not be caused by HPV
- Cervical cancer is uncommon which makes it hard to determine the sensitivity of any screening test
- The only accurate approach to evaluating the performance of screening tests for cancer is to use registry data and long-term follow-up studies

CR-78

Registry Data on Screening History of Women with Cervical Cancer

Description	No Recent Cytology* %	Cytology WNL (FN) %	Failure to Follow-up %
Kaiser ¹	56	32	13
Sweden ²	64	24	11
Netherlands ³	63	23	13
New Zealand ⁴	51	37	12

* Different definitions in the different studies

¹Leyden et al. *JNCI* 2005; 97:675

²Andrae et al. *JNCI* 2008; 100:622

³Gok et al. *BJC* 2011; 104:685

⁴Priest et al. *BJOG* 2007; 114:398

CR-79

Cytology in Cervical Cancer

Kaiser N. California Experience

- Reviewed screening histories of 965,360 women ≥30 years between 2003 and 2010

	SCC	AdenoCA
Total	198	114
NILM result	41 (20.7%)	52 (45.6%)
ASC-US / ASC-H	27 (13.6%)	14 (12.3%)
Other abnormal	130 (65.7%)	50 (43.9%)

Katki et al. *JLGTD* 2013; 17: S28

CR-80


HPV Testing Prevents Cervical Cancer

- We now have evidence that using HPV testing for screening reduces the incidence of invasive cervical cancer compared to cytology
- 4 European randomized trials conducted in Sweden, Netherlands, UK, and Italy
- Included 176,464 women 20-64 years of age
- Follow-up for a median of 6.5 years (1,214,415 person-years)

HPV testing used an HPV assay other than the cobas® HPV Test
Ronco et al. *Lancet* pub online, 2013

CR-81

HPV Testing vs Cytology for the Prevention of Cervical Cancer



HPV testing used an HPV assay other than the cobas® HPV Test
Ronco et al. *Lancet* pub online, 2013

CC-82


Risk of \geq CIN3 After a Negative Screening Test 3 Years of Follow-up

	Pap	HPV	Cotest
Dillner et al.	0.50%	0.11%	0.06%
Katki et al.	0.17%	0.06%	0.05%
Rijkaart et al.	0.26%	0.06%	0.05%
ATHENA	0.78%	0.34%	0.30%

HPV testing used an HPV assay other than the cobas® HPV test (except ATHENA data)
 Dillner et al. *BMJ* 2009;377: 21,351 women \geq 20 years; Katki et al. *Lancet Oncol.* 2011;12:663; >300,000 women \geq 30 years; Rijkaart et al. *Br. J. Cancer* 2012;106:975; >25,658 women 29-61 years; ATHENA; 41,955 women \geq 25 years

CR-83


Cumulative Incidence of \geq CIN3 in Kaiser



HPV testing used an HPV assay other than the cobas® HPV test
¹Katki HA. et al. *Lancet Oncol.* 2011; 12:663

CR-84

Cumulative Incidence of \geq CIN3 in Kaiser



HPV testing used an HPV assay other than the cobas® HPV test
¹Katki HA. et al. *Lancet Oncol.* 2011; 12:663

CR-85


Summary

Discussion of Clinical Implications

- Using HPV alone will not adversely impact women with other gynecological cancers or STIs
- No screening test with acceptable specificity will detect all cervical cancers or precursors
- HPV alone offers greater protection against CIN3 and invasive cervical cancer than cytology alone – **widely used in the U.S.**
- Provides similar protection against CIN3 and invasive cervical cancer as cotesting

CR-86

Candidate Screening Algorithm cobas® HPV with 16/18 Genotyping and Reflex Cytology



Detects significantly more disease than Comparator Algorithm
(cytology alone) or Hybrid Cotesting
Greatly simplifies screening algorithms

CR-87

Summary

Christoph Majewski, PhD

Life Cycle Leader, HPV and Microbiology
Roche Molecular Systems



CS-88



Predictive Values and Likelihood Ratios of the Candidate, Comparator and ATRI NM ≥ 30 GT

	Detection of $\geq \text{CIN3}$				Colposcopy / $\geq \text{CIN3}$
	PPV %	NPV %	PLR	NLR	
Candidate ¹	12.25	99.58	14.24	0.44	8.1
Comparator ¹	6.47*	99.41*	7.06*	0.61*	15.3*
ATRI NM ≥ 30 GT ¹	11.04*	99.52*	12.66*	0.49*	9.1*

3 year CIR of Candidate is 0.34 compared to 0.78 for cytology and 0.30 for cotesting (calculated for women ≥ 25 years)

¹Verification bias adjusted

*Difference of VBA parameters statistically significant

CS-89