



*Minimal Residual Disease (MRD) as a Surrogate Endpoint in ALL  
FDA Workshop, 18 April 2012, Silver Spring, MD*

# Minimal residual disease in acute lymphoblastic leukemia


Technical, economical, and validation – QC  
considerations for multicenter MRD assessment

J.J.M. van Dongen





# Detection of minimal residual disease (MRD) in ALL



DCLSG ALL-8

I

C

II

maintenance Rx



# Detection of minimal residual disease in acute leukemia

Technique	Applicability	Detection limit	Remark
Flow cytometry (4 to 6 colors)	BCP-ALL: 85% T-ALL: 90% AML: 60-70%	( $10^{-3}$ -) $10^{-4}$	Fast, but variable sensitivity because of similarities between normal (regenerating) cells and malignant cells
PCR of Ig/TCR genes	BCP-ALL: 95% T-ALL: 95% AML: 10-15%	$10^{-4}$ - $10^{-5}$	Time consuming and relatively expensive (junctional region sequencing), but applicable in $\geq 95\%$ of lymphoid malignancies
PCR of fusion transcripts and mutations	BCP-ALL: 40% T-ALL: 25% AML: 25-40%	$10^{-4}$ - $10^{-6}$	Limited applicability in ALL, but potentially useful in specific subgroups, e.g. BCR-ABL cases in specific protocols



# MRD diagnostics in ALL

## From research tool to surrogate endpoint in ALL treatment


1. Definitions of MRD remission, MRD recurrence, relapse
2. When should MRD be measured and which sensitivity
3. Evaluation of:
  - Induction Treatment
  - Continuous Monitoring
  - Treatment blocks (new drugs)

} Dependent on disease category and treatment protocol
4. Comparison between MRD studies: time points and quantitative range & sensitivity
5. MRD techniques: Ig/TCR PCR and/or Flow Cytometry
6. Standardization and Quality Control

Collaborative networks on standardization & quality control



# From Ig gene to Ig molecule



VH3-21 (germline)

**TGTATTACTGTGCGAGA**

TGTATTACTGT  
TGTATTACTGTGCG  
TGTATTACTGTGC  
TGTATTACTGTGCGAG  
TGTATTACTG  
TGTATTACTGTGCG  
TGTATTACTGTGCGA

insertion

AGGC  
TATCCGGA  
CCGGACTG  
CTGAGTC  
ACATCGA  
CGT  
CCGG

DH3-3 (germline)

**GTATTACGATTTTTGGAGTGGTTATTATACC**

CGATTTTTGGAGTGGTTATTATA  
TTACGATTTTTGGAGTGGTTATTATAC  
TTTTGGAGTGGTTATTATACC  
TATTACGATTTTTGGAGTGGTTAT  
CGATTTTTGGAGTGGTTATTATA  
TACGATTTTTGGAGTGGTTATTAT  
TTACGATTTTTGGAGTGGTTATTATACC

insertion

GTCCA  
CGATCG  
GGT  
CGTAGCGTA  
CGTAG  
GGCTAAGG  
CGGAGC


JH4-1 (germline)

**ACTACTTTGACTACT**

TGACTACT  
CTTTGACTACT  
ACTACTTTGACTACT  
TTTGACTACT  
ACTTTGACTACT  
TGACTACT  
TACTTTGACTACT




# Relapse free survival in I-BFM-SG study according to the combined MRD information at time points 1 and 2 (n=129)






# PCR analysis of Ig/TCR genes



→ High levels of standardization required

# Guidelines for RQ-PCR analysis of TCR/Ig gene rearrangements








**EuroMRD**

European Study Group on MRD detection

Chairman: J.J.M. van Dongen

[www.EuroMRD.org](http://www.EuroMRD.org)

Europe



Israel



Japan



Australia




Singapore

43 laboratories in 18 countries



Supported by Leukaemia & Lymphoma Research, LeukemiaNet, and EuroClonality



# MRD diagnostics in infant ALL: Interfant-99 protocol



# Current MRD technique in lymphoid malignancies




## Disadvantages of Ig/TCR-based MRD-PCR techniques:

- labor intensive (junctional regions per patient);
- require specialized laboratories;
- time consuming ( target identification: 4 to 6 weeks)

➡ Faster technique needed: 8-color flow cytometry ?




# EFS of MRD-based risk groups (FCM at day 29) in COG protocol



80% of patients  
51% of relapses




# MRD window, time points, MRD techniques and QR & sensitivity







# Precursor-B-cells in BM during ALL treatment



# Therapy-induced immunophenotypic shifts




# RQ-PCR and flowcytometric MRD in childhood ALL



Ryan J., et al. MRD detection in childhood ALL patients at multiple time-points reveals high levels of concordance between molecular and immunophenotypic approaches. *Br J Haematol* 2008 144: 107-115




# RQ-PCR and flowcytometric MRD in childhood ALL




Thörn, I. et al. MRD assessment in childhood ALL: a Swedish multi-centre study comparing real-time polymerase chain reaction and multicolour flow cytometry. *Br J Haematol* 2011 152: 743-753

# RQ-PCR and flowcytometric MRD in childhood ALL (only positive data in quantitative range)



# RQ-PCR and 4-color/6-color flow cytometric MRD in childhood ALL



V.H.J. van der Velden et al., Unpublished results

# MRD-based risk groups (day 33 and day 78) in 171 patients of the DCOG-ALL10 protocol

		RQ-PCR (Ig/TCR) based risk groups			
		HR (6%)	MR (62%)	LR (30%)	NA <sup>a</sup>
FCM-based risk group	HR (5%)	6 (4%)	2 (1%)	0 (0%)	0 (0%)
	MR (38%)	3 (2%)	60 (35%)	1 (1%)	0 (0%)
	LR (59%)	0 (0%)	44 (26%)	49 (29%)	6 (4%) <sup>a</sup>

<sup>a</sup> Not applicable: Six patients could not be classified with molecular MRD analysis (no Ig/TCR marker with at least a quantitative range of  $10^{-4}$ ).



## Risk group definition

Classical clinical  
risk groups at  
diagnosis



No treatment  
received



Flow cytometry  
MRD risk groups  
at day 8/15 (33?)



Corticosteroids  
evaluated?



Ig/TCR-based  
PCR MRD risk  
groups at day 33/  
week 12



Full induction  
evaluated?



Different composition of risk groups  
(25-40% shifts between SR and MR)



## Current position of FCM in MRD diagnostics

1. FCM has proven to be useful for MRD detection, **BUT:**
  - does FCM measure the same as PCR?
  - can FCM replace PCR?
  - can FCM supplement PCR?  
(e.g. in cases without sensitive Ig/TCR targets)
2. FCM can easily become broadly available
  - no special laboratory facilities
  - high-tech FCM equipment ( $\geq 8$  colors) for moderate costs

### !!! ATTENTION FOR:

- standardization
- improvement of specificity and sensitivity (new markers)
- broadly accepted antibody protocols
- international guidelines for data acquisition and interpretation






# Aims in Flow cytometric MRD diagnostics

## Fully standardized flow cytometric MRD detection !

1. Multicenter design, standardization, and clinical (protocol-related) evaluation of innovated flow cytometric MRD detection:
  - 8 colors: increased sensitivity
  - new markers (particularly fusion proteins/oncoproteins): increased specificity
  - new software (fast, easy, automated)
2. Evaluation of flow cytometric MRD detection in full parallel to Ig/TCR based MRD detection, using strict international guidelines for instrument settings, data acquirement, and data interpretation

# BCP-ALL panel



Mix of 3 different regenerating B cell populations (Haematogones)

BCP-ALL blast cells




**EuroFlow**

Responsible scientist: L Lhermitte




# Precursor B-cell differentiation in normal vs regenerating bone marrow

Normal BM



Normal vs **Regenerating** BM







CD19 Gated B-cells (excluding PC)

# Four BCP-ALL cases vs normal precursor B-cells







EuroFlow

## APS view 1



## APS view 2



Case 1

Case 2

Case 3

Case 4



## *Minimal Residual Disease (MRD) as a Surrogate Endpoint in ALL FDA Workshop, 18 April 2012, Silver Spring, MD*

### Conclusions

1. PCR-based MRD diagnostics (IG/TCR genes or fusion genes) is currently the gold standard in many European ALL protocols
2. Differences in MRD value between protocols is mainly caused by application of different non-standardized MRD techniques, which also differ in sensitivity.
3. PCR-based MRD diagnostics can potentially be replaced by 8-color flow cytometry (Novel developments are required)
4. Standardization, regular Quality Control, and guidelines for data interpretation and data reporting are essential for international comparability of MRD results (within and between treatment protocols).

Collaborative networks on standardization & quality control are essential