

Flow Cytometric Assessment of MRD

Brent Wood MD, PhD
Department of Laboratory
Medicine
University of Washington

Identification


Abnormal population identification

- Normal
 - Antigens expressed in consistent and reproducible patterns with maturation
- Neoplastic
 - Increased or decreased normal antigens
 - Asynchronous maturational expression
 - Aberrant antigen expression
 - Homogeneous expression


Normal B cell Maturation

	Immature Early	Immature Mid	Immature Late	Naive / Mantle	Follicle center	Marginal zone	Plasma cell
CD45	+	+	+	+	+	+	+
TdT	+	+	+	+	+	+	-
CD34	+	+	+	+	+	+	-
CD10	+	+	+	+	+	+	-
CD38	-	-	-	-	+	+	++
CD19	+	+	+	+	+	+	+
CD20	-	+	+	+	+	+	+
CD22	+	+	+	+	+	+	+
HLA-DR	+	+	+	+	+	+	+
IgD	-	-	-	+	+	+	-
IgM	-	-	-	+	+	+	-
Kappa	-	-	-	+	+	+	+
Lambda	-	-	-	+	+	+	+

Normal B cell Maturation




ALL MRD



0.1% abnormal immature B cells

Informative Immunophenotype



ALL Informative Antigens

Table 1. Useful antibody combinations for Immunophenotypic MRD detection.

<i>Antibody combinations*</i>	<i>Suitable cases (n) (%)</i>		<i>LAIP frequency in childhood ALL^o</i>
CD38/CD10/CD34/CD19	36/63	57.1%	30-50%
CD45/CD10/CD34/CD19	17/64	26.6%	30-50%
CD21/CD10/CD34/CD19	0/57	0	5-10%
CD22/CD10/CD34/CD19	0/55	0	20-30%
CD58/CD10/CD34/CD19	11/62	17.7%	40-60%
TdT/CD10/CD34/CD19	22/61	36.1%	30-50%
CD13/CD10/CD34/CD19	8/62	12.9%	10-20%
CD15/CD10/CD34/CD19	3/58	5.2%	5-10%
CD33/CD10/CD34/CD19	2/59	3.4%	5-10%
CD65/CD10/CD34/CD19	0/59	0	5-10%
CD66c/CD10/CD34/CD19	15/62	24.2%	10-20%
CD10/NG2/CD34/CD19	9/60	15.0%	3-5%
CD10/CD56/CD34/CD19	3/57	5.3%	5-10%

*FITC/PE/PerCP/APC; ^oRefs. 2,3,5,8.

Minimal MRD Reagents

- Simple two-tube panel
 - CD10 FITC / CD20 PE / CD45 PerCP / CD19 APC
 - CD34 FITC / CD9 PE / CD45 PerCP / CD19 APC
- Abnormalities detected (N=82)
 - Tube 1 = 93% of cases
 - Tube 2 = 94% of cases
 - Tubes 1 and 2 combined = 99%


Risk Adapted Therapy - AALL03B1

- 6 color flow cytometric assay

	FITC	PE	PerCP-Cy5.5	PE-Cy7	APC	APC-H7
Tube A	CD20	CD10	CD38	CD19	CD58	CD45
Tube B	CD9	CD13/3 3	CD34	CD19	CD10	CD45
Tube D	Syto16		CD3	CD19	CD71	CD45

- Collect 750,000 events


CD45/SS




Borowitz et al (1993) AJCP 100:534-40.

Steltzer et al (1993) Ann NY Acad Sci 667:265-280

Normal Blast Maturation



Normal Granulocyte Maturation



AML Informative Antigens

Table 2

Distribution of LAIP and of LAIP classes in 1400 patients with newly diagnosed and untreated AML (Laboratory for Leukemia Diagnostics, Munich, Germany)


LAIP class	LAIP	n	%
Asynchronous	Total	652	20.6
	CD11b+CD117+CD34-	156	4.9
	CD11b+CD117+CD34+	92	2.9
	CD11b+CD117-CD34+	36	1.1
	CD34+CD116+CD33+	113	3.6
	CD34+CD15+CD33+	193	6.1
	CD65+CD87+CD34+	12	0.4
	CD65+CD87-CD34+	50	1.6
Cross-lineage	Total	742	23.5
	CD34+CD13+CD19+	48	1.5
	CD34+CD2+CD33+	51	1.6
	CD34+CD56+CD33+	83	2.6
	CD34-CD13+CD19+	21	0.7
	CD34-CD2+CD33+	33	1.0
	CD34-CD56+CD33+	189	6.0
	CD4+CD13+CD14-	87	2.8
	CD7+CD33+CD34-	75	2.4
	CD7+CD33+CD34+	155	4.9

Lack of expression	Total	625	19.8
	CD15+CD13+CD33-	6	0.2
	CD15+CD13-CD33+	7	0.2
	CD34-CD135+CD117+	17	0.5
	CD38-CD133+CD34+	10	0.3
	CD4+CD13-CD14+	7	0.2
	CD9-CD34+CD33+	30	0.9
	CD9-CD34-CD33+	34	1.1
	HLA-DR+CD33-CD34+	12	0.4
	HLA-DR-CD33+CD34-	143	4.5
Overexpression	HLA-DR-CD33+CD34+	37	1.2
	MPO+LF-eCD15-	315	10.0
	MPO+LF-eCD15+	4	0.1
	MPO-LF+eCD15+	3	0.1
	Total	1139	36.1
	CD11b-CD117++CD34+	9	0.3
	CD13++CD34++	163	5.2
	CD15++CD13++CD33++	52	1.6
	CD34++CD135+CD117++	35	1.1
	CD34++CD33++	65	2.1
	CD34-7.1++CD33+	53	1.7
	CD36++CD235a++CD45(+)	25	0.8
	CD38++CD133++CD34++	16	0.5
	CD4++CD64++CD45++	144	4.6
	CD4+CD13++CD14++	19	0.6
	CD61++CD14-CD45+	5	0.2
	CD65++CD87++	162	5.1
	CD90++CD117++CD34+	23	0.7
	HLA-DR++CD33++CD34++	41	1.3
	TdT(+)+eCD33++eCD45++	327	10.4
Total		3158	100.0


Average 2.3 LAIP per patient

Population Identification

- Leukemia-Associated immunophenotype (LAIP)
 - At diagnosis
 - Immunophenotypic deviation relative to normal is identified
 - Informative reagent combinations selected (reduced or custom set)
 - Gate is created for monitoring
 - Follow-up
 - Run reduced reagent combination
 - Count events in pre-defined gate
- Deviation from normal maturation
 - At diagnosis
 - Immunophenotypic deviation relative to normal counterpart is identified
 - Uniform reagent combinations utilized
 - Follow-up
 - Identify discrete population having immunophenotype different than normal
 - Use diagnostic immunophenotype as starting point




LAIP



Focus is immunophenotype, not population


Immunophenotypic Stability

- ALL




Immunophenotypic Stability T-ALL


Figure 2




CD48



CD48 MRD





Population Identification

- Leukemia-Associated immunophenotype (LAIP)
 - Advantages
 - Conceptually simple and objective
 - Reduced reagent expense for follow up
 - Disadvantages
 - Requires pre-treatment sample to define LAIP
 - Requires immunophenotypic stability
 - Any event in pre-defined gate regarded as MRD
- Deviation from normal maturation
 - Advantages
 - Does not require pretreatment sample
 - Uniform reagent combinations utilized
 - Improved specificity through population identification
 - Less sensitive to immunophenotypic instability
 - Disadvantages
 - Requires detailed immunophenotypic knowledge (expert)
 - Subjective
 - More time consuming

Timing

- Induction nadir (day 14)
 - Reduced background populations
 - Hypoplastic with many apoptotic cells
- End of induction
 - ALL - Few immature B cells
 - AML - Active marrow regeneration, increased precursors
- End of consolidation
 - ALL - Larger number of immature B cells
 - AML - Normal marrow populations

Enumeration


Sample Acquisition

- Identification
 - Distinguish normal from abnormal
 - Degree of immunophenotypic aberrancy
 - Number and immunophenotype of background populations
 - How many events define a population?
 - 10-50 events
- Enumeration
 - Requires complete discrimination of population
 - Insufficiently informative immunophenotype
 - Maturational expression
 - Reproducibility (Poisson counting statistics)
 - $CV \sim \sqrt{N}/N$
 - 100 events gives CV of 10%
 - Sensitivity of 0.01% requires 1,000,000 events


Denominator

- Total nucleated cells
 - Most comparable to morphology
 - DNA binding dye often used (Syto16, Draq5, etc.)
 - Incomplete RBC lysis, platelet aggregates
 - Under-representation of NRBCs with lysis and washing
- White cells
 - CD45 positive cells + neoplasm
 - Variable CD45 on early NRBCs
 - Overestimation with erythroid hyperplasia
- Mononuclear cells
 - Exclude granulocytic lineage (high side scatter)
 - Most comparable to Ficoll-prepared samples
 - Early MRD literature used Ficoll
 - Reduced variability due to granulocytic degeneration
 - Shipped samples

Denominator



Utility of Denominator Tube



Unpublished data





Hemodilution

- Bone marrow is a semi-solid tissue
 - Absolute cell concentration has little meaning
- Marrow aspiration is a traumatic procedure
 - Variable amount of peripheral blood introduced
 - Increased amounts of blood with each subsequent aspiration
 - 1st aspirate should be used for MRD
- Not a major problem for many samples
 - Problem in hypocellular marrows, high PB WBC count or poor quality aspirates
- No method for accurate correction
 - One method for normalization proposed for blast counts

Sources of Variability

- Identification (false positive or negative)
 - Insufficiently informative reagents
 - Improper assay validation
 - Immunophenotypic shift
 - Inexperienced interpreters
- Quantitation
 - Too few events acquired
 - Denominator effects (2 fold)
 - Sample degeneration
 - Hemodilution


Reproducibility




Listmode file exchange

Dworak, et al (2008) Cytometry 74B:331-340

Reproducibility




Paired sample exchange




Artificial diluted samples

Reproducibility




Reproducibility of MRD detection



Tube 1
Unpublished data, courtesy Mike Borowitz

Flow MRD on AALL03B1




* day 8 M1 patients excluded

Unpublished data, courtesy Mike Borowitz


Day 29 Flow MRD on AALL0232

%
MRD +



COMPARISON OF OUTCOMES


UW vs JHH MRD .01-0.1%



Unpublished data, courtesy Mike Borowitz

COMPARISON OF OUTCOMES


UW vs JHH MRD .1-1%



Unpublished data, courtesy Mike Borowitz

COMPARISON OF OUTCOMES

UW vs JHH MRD >1%



Unpublished data, courtesy Mike Borowitz

Conclusions

- Flow cytometry
 - Capable of accurate MRD assessment
 - Interpretative assay
 - Sensitivity
 - Dependent on
 - Antibody combination (informative)
 - Number of cells evaluated
 - Time point
 - Less than PCR (10 cells vs 1)
 - 0.01% routine for ALL
 - 0.1% routine for AML
 - Reproducible
 - Correlates with clinical outcome