

In Vitro & CFD Bioequivalence Testing for Orally Inhaled Drug Products

Peter R Byron, Michael Hindle and Worth Longest
Schools of Pharmacy and Engineering, Virginia
Commonwealth University, Richmond, VA



Research support: MCV Foundation and FDA (Award # FD004570)

Based on a belief that 2 inhalers should be “BE” when their drug deposition in lung occurs in the same form, doses and locations... *“we set out to research [biorelevant in vitro and CFD] methods to partner realistically-designed airway models with representative inhalation profiles, so that ...proving drug deposition equivalence was facilitated....”*

Peter Byron, 2010.


Where are we now and what should we do to move forward?

Biorelevant Test Methods



VCU Mouth-Throat Model (RDD 2010)

Filter



- ☐ Realistic geometries
- ☐ Internal surfaces coated
- ☐ Realistic airflow profiles
- ☐ Total Lung Dose in vitro = $TLD_{in\ vitro}$ = Drug mass escaping MT

Outline

■ New method development

- VCU's 2010 model was “hypothetically medium - sized”
- “Large” and “small” models developed and paired with “simulated, realistic, inhalation profiles”
 - Models validated “geometrically”- anatomy literature
 - Results from “*in vitro* methods” compared to deposition data in literature from trained humans.. **IVIVC**

■ Choosing the inhalation profiles

- Realistic ranges for DPI inspiratory maneuvers

■ Predicting regional lung deposition based on aerosol properties of ***TLD_{in vitro}*** (use of validated CFD)

■ The future



VCU

School of Pharmacy


New Method Development and IVIVC

- Scaled MT-TB models for normal human adults



- “VCU Medium” MT model scaled by volume ± 2 SD from literature
 - $S = 27.2$, $M = 65$; $L = 107.8 \text{ cm}^3$
 - Same scaling factors used for TB
 - This normal distribution of volumes appears consistent with anatomy literature & linear scaling factors
 - Length x 0.748 = small
 - Length x 1.0 = medium
 - Length x 1.165 = large
- ...models shown at left
- www.rddonline.com/resources/tools
- MT designed to accept inhaler mouthpiece adapters

Budelin Novolizer: $TLD_{In Vivo}$ vs $TLD_{In Vitro}$



- ❑ In vivo results - gamma scintigraphy [Newman, Eur. Resp. J. 16: 178]
- ❑ IVIVC from 3 models – flow profiles simulated to match Newman
- ❑ Error bars = entire range(all cases)

Summary

- ❑ Median & range of $TLD_{in vivo}$ correlates with $TLD_{in vitro}$
 - ❑ when simulated flow extremes coupled with upper airway geometry extremes for a mixed gender, adult population.
 - ❑ Statistically significant differences between S, M and L model correlations (Budelin Novolizer)
- ❑ Median $TLD_{in vivo}$ also correlates with $TLD_{in vitro}$ in VCU_{medium} model for Handihaler (tiotropium + lactose), Aerolizer (formoterol + lactose), Easyhaler (albuterol + lactose), Turbohaler (terbutaline)
 - ❑ Delvadia et al, JAMPDD 2013, 26: 138 - 144
- ❑ Product comparisons best performed with inhaler - representative breath profiles
- ❑ Need to determine how $TLD_{in vitro}$ deposits regionally.




VCU

School of Pharmacy

Profile Analysis – toward standard profiles

- Normal profiles, across resistances, DPI trained, 20 adult volunteers
- Gray profiles = Flow rate from mouthpiece
- **Red** profiles = 10, 50 and 90 percentiles
- **Black** = sine wave curve-fit to 50% profile (breath simulator)




VCU

School of Pharmacy


Where does CFD come in?

- ❑ Coupling careful modeling with *in vitro* testing enables CFD model validation.
 - ❑ e.g. Novolizer (75 LPM for 4 s); Respimat at 37 LPM (Medium MT-TB)
- ❑ Tian et al. (2012) *Aerosol Sci. Technology* 46, 1271-1285

Novolizer




Respimat




CFD Models for Regional Distribution

❑ Based on size distribution of **TLD_{in vitro}** (drug aerosol entering lung) and validated CFD model predict regional distribution in lung.

❑ Tian et al. (2012) *Aerosol Sci. Technology* 46, 1271-1285



Respimat: SD polydisperse



Stochastic Individual Pathway model



The Future for Inhaler Comparisons

- ☐ Validated “realistic” mouth-throat models (adult human: S, M, L)
- ☐ Public database of inhalation profiles
 - ☐ Median & CIs for different airflow resistance DPIs
 - ☐ “Leaflet training” vs “personal training (Rx)” [VCU in preparation]
 - ☐ pMDIs, Gender, age, disease effects needed (*TBD*)
- ☐ Use new *in vitro* tests (with *IVIVC*) to compare values for *TLD_{in vitro}*
- ☐ Measure APSD emitted from MT... or MT-TB models with realistic profiles [use *in vitro* data from MT-TB to validate CFD model]... *TBD*
- ☐ Predict *regional* lung deposition using CFD for realistic breath profiles (noting that CFD is most reliable under the lower Reynold’s number conditions typical of lung generations 4 through 23)... *TBD*
- ☐ Accepted bio-relevant *in vitro* tests coupled with CFD predictions
 - ☐ Easier bridging
 - ☐ Easier “bioequivalence” arguments
 - ☐ Improved understanding (QbD) and ↓ testing



VCU

School of Pharmacy