



Assurance of physical equivalence (size, charge, concentration, aggregation, interactions) of complex nanomedicine substances



Hans van der Voorn, CEO
Dr. Darby Kozak, Chief Scientist
Izon Science United States

Context:

- **Nanomedicine regulation**
 - Nanoparticle based delivery of generics
 - Nanomilled APIs
 - Vaccines
- **Innovation will provide medical and cost benefits, better ways of using existing drugs**
- **Safety and certainty are still required**

Rapid development and use of complex nanomedicine products

• From the lab



• To Market




• To Generic



AZAYA THERAPEUTICS



Importance of the physiochemical properties of drug delivery particles



Concentration

- Drug delivery dosage
- Disease stage/progression
- Toxicity

Size

- Localized accumulation
Nasal, liver, kidneys, EPR effect, blood-brain barrier,
- Effluent pathway
- Toxicity

Surface Charge (ζ -potential)


- Endocytic mediated pathway
- Toxicity

*Need for accurate and sensitive analysis techniques
that measure the properties of each particle in physiological conditions*

Proposal:

- Demonstrate physical equivalence by accurate, precise and repeatable suite of particle measurements
- “Particle fingerprinting” (previous speaker’s proposal)
- This would replace present method of requiring no change to manufacturing process, which is not practical in scale up situations
- Needs to be practical and cost effective

Long established techniques have insufficient resolution for nanomedicines





Dynamic Light Scattering is the most common particle analysis technique but lacks resolution and repeatability

Particle-by-particle “**number-based**” distributions can now provide a more realistic and accurate measure of particle size distribution

New techniques offer solutions but need FDA support

TRPS vs DLS vs EM of complex particles: Carbon Nanohorns



Possible Suite of Measurements:



- Particle number (concentration)
- Accurate, calibrated size distribution based on particle by particle analysis
- Particle by particle surface charge (charge affects behaviour as much as size)
- Particle interactions, including aggregation and biomolecule interaction

Nanomedicine examples

- Liposomes, eg Doxil
- Emulsions
- Bioparticles
- Virus like particles (VLPs)
- Polymer capsules
- NP based vaccines
- Virus based vaccines

TRPS – Single Particle Analysis

When a particle enters the pore it blocks a proportion of the ions flowing through the pore – this causes a pulse in the ionic current the magnitude and duration of the pulse is proportional to the particle size and velocity, respectively.




Particle-by-particle measurements provide high resolution analysis

TRPS Principles





Assurance through calibration





Particle Rate \propto Concentration

Particle Concentration




Size-specific concentration

Measuring the properties of each particle and their concentration enables the deliverable liposome volume fraction “dosage” to be calculated



Liposome volume fraction

Particle volume
 $\frac{4}{3}\pi r^3$ (size)




Particle number
(concentration)




iZON

TRPS Size Analysis for Optimizing Liposome Synthesis




High resolution means confidence and precision in tailoring the size of all the liposomes in solution.

Caelyx vs a generic




TRPS Size Analysis for Monitoring Liposome Stability



High resolution means confidence in understanding the affects of treatment and storage conditions of long-term stability of liposome solutions.

Field based measurement may be required

Using Izon instruments to measure particle charge



- Particle velocity through the pore is due to:

Fluid Velocity (convection)

Directionally dependent and proportional to the applied pressure.


Electrophoretic Mobility

Proportional to particle charge (ζ -potential) and applied voltage.

Electro-osmosis

Proportional to pore charge (ζ -potential) and applied voltage.

High Resolution Particle-by-Particle Size & Charge Analysis




Represents a new way and ability to characterize and understand particles


Each dot represents a single particle and histograms illustrate the distribution

Size & Charge Analysis of Liposomes

Modifying the surface of liposomes is used to improve their circulation lifetime (e.g. PEGylation) and their specific cellular targeting (e.g. antibodies). These surface features can give rise to a change in the liposome size and charge.



Particle-by-Particle Size & Charge Analysis of PEGylated & 'normal' Liposomes



PEGylation – replacing a phosphocholine lipid with a glycolated phospholipid - gives rise to a negative charge.

The degree of liposome PEGylation can be monitored from the corresponding negative shift and distribution in particle zeta-potential.

Analysis for nanomedicine regulation

- **Need for particle by particle sensitive analysis which is accurate, reliable, repeatable and validated**
- **Size, charge, concentration, interactions**
- **Leverage current nanotechnology development to improve regulatory outcome and confidence**

Analysis for nanomedicine regulation

- **Comprehensive high resolution analysis to characterize a wide range of nanoparticle systems and individual particle properties**
- **Benefits in assurance for FDA**
- **Benefits in assurance, practicality and cost savings for industry**

Next Step:

- **Form project consortium for pilot scheme**
- **Input from:**
 - **FDA, EMA**
 - **Measurement institutes – NIST, NCL, NPL (UK)**
 - **Nanomedicine companies (4 identified)**
 - **Izon**



Thank You

www.izon.com