

Impact of Formulation and Process Variations on the Quality of Freeze Dried Products: How do we identify, control, and characterize critical variations?

**Michael J. Pikal
School of Pharmacy
University of Connecticut
A NIPTE Institution**


Some Quality Attributes for Freeze Drying- not all Critical

- Sterility-critical
 - Low endotoxin-critical
 - Stability-critical*
 - adequate potency
 - absence of toxic degradation products
 - Rapid and easy reconstitution
 - Cost effective process (i.e., fast)
 - Fast development process (speed to market)
 - “Elegance” - “beauty is in the eye of the beholder”
-
- Note: some are impacted by formulation, some by process, and many by both formulation & process

Storage Stability is Very Sensitive to Formulation

Human Growth Hormone Formulations

Rate Constants, $k(vt)$, at 40C (40 C \ll Tg)



* all formulations except Gly:Mann are glassy

**** Trend is same for both chemical degradation and aggregation!**

Significance of Results

- Standard, “Current most used” commercial formulation
 - $\approx 3\%$ aggregation in 2 yrs at Refrigerated Storage
(estimated from literature data)
- Simple Sucrose Based Formulation
 - ≈ 3000 years to form 3% aggregate at 25 C°
(estimated from literature data)
- Formulation Does Matter!

Control of Stability by Process Control

Use “Good Freeze Drying Practice” to guide control

- **Freezing-convert water to ice**
 - Control by control of ice nucleation temperature
 - Stability implications
- **Primary Drying-sublime the ice**
 - Control Product Temperature History
 - You freeze dry the product, not the shelf! (Felix Franks)
 - “Possible” stability implications
- **Secondary Drying-remove non-freezable water**
 - Control Product Temperature
 - Stability implications

Control of Freezing


- **Conventional freezing is uncontrolled**
 - Variable ice nucleation temperature
 - Development to production, intra- and inter-batch
 - Variable ice surface area means variable drying time and (potentially) variable stability
 - High surface area means potential for protein aggregation.
- **Solution: Controlled Ice Nucleation**
 - At least two commercial versions available for both manufacturing and development
 - Nucleate at moderate super-cooling (low surface area)
 - Same process for development as in manufacturing

hGH Aggregation Study

Compare Stability “in-process” and “in storage” for conventional uncontrolled nucleation and controlled nucleation

- Experiments

- 2 mg/mL hGH + 2 mM sodium phosphate + 6 mg/mL sucrose
- Assay hGH aggregation via HPLC (pre- & post-freeze drying and after 14 days at 50°C)



Examine Variance by Propagation of Errors[#]


Sources of variance:

Fill Volume, Pressure Variation,


Product Resistance* (ice nucleation), Heat transfer Coefficient*

* Major Contributor

Vial Heat Transfer



Product Resistance




Note: 5% RSD more typical of most experimental systems


[#] NIPTE-FDA Collaboration Projects

Impact on Primary Drying Time

Low Variances (desirable)



High Variances (typical)




Solution: Minimize variance and need to adjust “Design Space” to accommodate expected variation



Product Temperature in Secondary Drying

- Exposure to high temperature Stabilizes!

Decarboxylation in Moxalactam Na



Aspartame: sucrose (1:10) formulation



Water Contents are essentially Identical!

Conclusion: Secondary Drying Temperature matters

Conclusions

- **Formulation Matters**
- **Process Matters**
 - Even if processes are “nominally” the same
- **Intrinsic Variation in Process Matters**