

BREATH BIOMARKERS OF RADIATION EXPOSURE

Menssana Research, Inc, Newark, NJ Michael Phillips, Renee Cataneo, Anirudh Chaturvedi, Peter Kaplan, Mayur Mundada, Urvish Patel University of Louisville, Kentucky Xiang Zhang Schmitt & Associates, Newark, NJ Peter Schmitt Southwest Research Institute, San Antonio, TX Mark Libardoni MD Anderson Cancer Center, Orlando FL Naren Ramakrishna New York University – Langone Medical Center Peter Schiff North Shore Radiation Therapy, Long Island, NY Richard Byrnes Swedish Medical Center, Seattle, WA Vivek Mehta

This project has been funded in whole or in part with Federal funds from the Biomedical Advanced Research and Development Authority, Office of the Assistant Secretary for Preparedness and Response, Office of the Secretary, Department of Health and Human Services, under Contract No. HHSO100201000010C



MENSSANA
RESEARCH, INC.

Objectives

Short-term

- *Identify biomarkers of radiation in breath*
- *Determine dose-effect response*

Long-term

- *Migrate findings to a point-of-care platform*
→ *a radiation biodosimeter for first responders*

Background

- *Normal human breath* - many volatile organic compounds (VOCs)
- *Concentrations* - low! – picomolar (parts per trillion)
- *Specialized instruments* → detect ~2,000 different VOCs in a sample
- *Breath VOCs contain biomarkers*

oxidative stress




lung cancer

breast cancer

pulmonary TB

transplant rejection





Hypothesis

RADIATION

- *free radicals in tissues*
- *oxidative stress*
- *volatile biomarkers in breath*
 - alkanes*
 - alkane derivatives*
 - other metabolites?*



Study design

- *Human subjects (n=32)*
 - *Age 18-75*
 - *Treated for lung, prostate or breast cancer*
 - *Daily fractionated regimens for 4 or 5 days*
 - *Low-dose group (n=27)*
 - Standard radiotherapy*
 - Daily dose = 180-400 cGy*
 - *High-dose group (n=5)*
 - Stereotactic body radiotherapy (SBRT)*
 - Daily dose = 700-1200 cGy*
- Breath sample collections – on days 1, 2, 3, 4, and 5*
- Before treatment*
- <30 min after treatment*



Sample collection


Breath collection apparatus (BCA)



MENSSANA
RESEARCH, INC.

Sample analysis

Two-dimensional gas chromatography and time-of-flight mass spectrometry




MENSSANA
RESEARCH, INC.

Data analysis



Data pre-treatment




Display aligned data

“Breath chip” heat map

At each time point:

*Mean of all subjects
compared to
mean baseline pre-radiation level*




Increased

Decreased

Day 1 post-radiation


t zero vs 2 hrs



Day 2


Pre-radiation

t zero vs 23 hrs



Post-radiation

t zero vs 25 hrs




MENSSANA
RESEARCH, INC.

Day 3


Pre-radiation

t zero vs 45 hrs



Post-radiation

t zero vs 47 hrs




MENSSANA
RESEARCH, INC.

Day 4


Pre-radiation

t zero vs 71 hrs



Post-radiation


t zero vs 73 hrs



MENESSANA
RESEARCH, INC.

Day 5

Pre-radiation




Post-radiation: only 2 high dose samples – did not align




MENESSANA
RESEARCH, INC.


t zero vs 2 hrs




t zero vs 23 hrs




t zero vs 25 hrs




t zero vs 45 hrs




t zero vs 47 hrs




t zero vs 71 hrs



t zero vs 73 hrs



t zero vs 96 hrs



*Patterns show
changes over five
days...*

*...but can they
identify radiation
exposure?*



MENESSANA
RESEARCH, INC.



MENSSANA
RESEARCH, INC.

Is it a new biomarker? Or just a face in the clouds?

- *Costs of assays falling*
- *BUT...costs of human and animal studies rising*
- *> 1,000 candidate biomarkers, <100 subjects*
 - *high risk of over-fitting data in statistical models*
 - *false-positive IDs of biomarkers*
 - *seeing “faces in the clouds”*



A rigorous method for biomarker identification

1. What was the experimental accuracy of a candidate VOC?

Area under curve (AUC) of receiver operating characteristic (ROC curve)

Range: 0.5 (coin-flip) to 1.0 (perfect accuracy)

2. What was its accuracy by chance alone?

Multiple Monte Carlo simulations

Assign random diagnosis (radiation or no radiation)

→ determine AUC of the ROC curve


→ repeat 40 times

→ mean chance accuracy

3. Compare the outcomes




No.
candidate
biomarkers




Breath VOCs

Day 1: Baseline pre-radiation
vs day 1 post-radiation



Day 5: Baseline pre-radiation
vs final breath collection



- *Similar Monte Carlo outcomes at all time points on days 1 through 5*
- *~50 breath VOCs greater than chance at each time point (>5 sigma)*
- *15 recurring breath VOCs seen at all time points*

15 biomarkers of radiation

Consistently observed at all time points

2,4,4-Trimethyl-3-(3-methylbutyl)cyclohex-2-enone

Cyclononene

1-Methylcyclohexa-2,4-diene

2,4-Octadiene

3-Heptene, 2,2,4,6,6-pentamethyl-

Benzene, (2-methylpropyl)-

Decanoic acid, nonadecafluoro-

Decane, 2,2,3-trimethyl-

Benzene, 1-cyclopenten-1-yl-

Benzene, (2-methyl-1-propenyl)-


4,5,6,7-Tetramethylphthalide

3-Tetradecene, (Z)-

4-Tetradecene, (E)-

Nitrous Oxide

Tetradecanal



*Methylated alkane –
oxidative stress product*



2,4,4-Trimethyl-3-(3-methylbutyl)cyclohex-2-enone

Cyclononene

1-Methylcyclohexa-2,4-diene

2,4-Octadiene

3-Heptene, 2,2,4,6,6-pentamethyl-

Benzene, (2-methylpropyl)-

Decanoic acid, nonadecafluoro-

Decane, 2,2,3-trimethyl-

Benzene, 1-cyclopenten-1-yl-

Benzene, (2-methyl-1-propenyl)-

4,5,6,7-Tetramethylphthalide

3-Tetradecene, (Z)-

4-Tetradecene, (E)-

Nitrous Oxide

Tetradecanal

**Can these biomarkers
provide a radiation
biodosimeter?**



MENSSANA
RESEARCH, INC.

Dose-response curve

*Correlate abundance of 15 biomarker VOCs
with total dosage (Gy)*

third order polynomial regression




Potential application as a biodosimeter

Collect breath sample

→ assay 15 VOC biomarkers

→ derive response function

→ read off dose from curve



input
response
function

output
dosage




MENSSANA
RESEARCH, INC.

Another approach:
Dose as a binary variable


Group 1	Group 2
<1.8 Gy	>= 1.8 Gy

Advantages:


- *Robust binary algorithm*
- *can be cross-validated*




Training set




Test set



Group 1	Group 2
<5 Gy	>= 5 Gy




*Training
set*



Test set

Binary radiation biodosimeter



Cutoff at 1.8 Gy

YES/NO accuracy

= 99 %

Cutoff at 5 Gy

YES/NO accuracy

= 78 %

Conclusions

In 32 humans irradiated over 4-5 days:

- *>50 breath VOCs changed significantly (>5 sigma)*
- *15 consistent biomarker VOCs at all time points
included oxidative stress products*
- *Continuous dose-effect curve
progressive increments through dose range 1.8 - 20 Gy
changes significant at all doses ≥ 1.8 Gy*
- *Binary dose-response
Accurate YES/NO responses at 1.8 Gy and 5 Gy cutoffs*



Next steps

Currently:

*Analyzing data - completed animal study
Irradiated Gottingen minipigs*

Near future:

*Total body irradiation studies
Humans
Non-human primates*

Longer term:

*Transfer algorithms to point-of-care breath test
→ rapid cost-effective biodosimeter*



BreathLink


Point-of-care system

Current BreathLink system

- Collects and concentrates breath and air VOCs
- Separates VOCs (1D GC)
- Detects VOCs (SAW)
- Assay in <6 min
- Detects active pulmonary TB

Next-generation BreathLink

- 2D POC GC x GC
- Smaller, more portable
- Will detect biomarkers identified by GCxGC TOF MS
- Optimized for first-responders and triage



MENESSANA
RESEARCH, INC.

Simple color-coded output for first responders



→ Rapid triage following a nuclear event

Thanks for listening!

