

# A Protein Biomarker Approach to Biodosimetry

**FDA/CDRH Public Meeting**  
***Regulatory Science Considerations for Medical Countermeasure***  
***Radiation Biodosimetry Devices***  
**October 27, 2012**




**George B. Sigal**  
**Research Director**  
**Meso Scale Diagnostics, LLC.**  
[gsigal@mesoscale.com](mailto:gsigal@mesoscale.com)

# Diagnostics for Large-Scale Radiation Incidents

## Scenario:

- Improvised nuclear device (10 kT) in metropolitan area
- Potential casualties from radiation exposure > 100,000
- Large numbers of injuries from blast effects and broken glass (2 mi radius)




Washington Post, 2007

## Need for Diagnostics (from BARDA BAA-BARDA-09-36):

- High-throughput (HT) and point-of-care (POC) tools that will:
  - Identify patients requiring urgent medical treatment
  - Provide assurance to individuals with low-dose exposure
  - Improve risk assessment for delayed effects of radiation exposure
  - Provide patient tracking
  - Potentially, monitor therapy and recovery
- Specific dose assessment needs:
  - Rapid classification of patients into low (< 2 Gy) and high (> 2 Gy) exposure
  - Quantitative assessment of dose between 0.5 and 10 Gy

# Radiation Response and Biomarkers




Factors affecting biological response to radiation

- Individual susceptibility
- Heterogeneity of received dose

# Biomarkers of Radiation Injury & Response

Measure Radiation Biomarkers



# MSD Biodosimetry Program

## Protein Biomarker-Based Dosimetry

- Multi-parameter approach employs biomarkers associated with different mechanisms of radiation response
- Able to provide dose assessment over extended time period
- Provides information that is directly linked to physiological injury and treatment

## Technology and Instrumentation

- State-of-the-art platform for multiplexed biomarker measurements in both laboratory (HT) and field (POC) settings
- Instrumentation designed for broad utility in medical countermeasures and commercial diagnostics – likely to be widely deployed in the future
- Offers a realistic solution to biodosimetry in response to a large-scale radiation event – with speed, throughput and accessibility



Project funded by the Biomedical Advanced Research and Development Authority (BARDA) which is part of the Office of the Assistant Secretary for Preparedness and Response within the US Department of Health and Human Services

Collaborating Laboratories:

- Armed Forces Radiobiology Research Institute (AFRRI): Drs. Ossetrova and Blakely
- National Cancer Institute (NCI): Dr. Camphausen

# Lab & POC: Common Biomarker Approach

## Multiplexed Biodosimetry Assay Panel

High  
Throughput  
Testing



POC  
Testing






# Biodosimetry Sample Collection Kit Components



## HT Kit

1. Bar-coded patient wrist band
2. Adhesive bandage
3. Single-use safety lancet
4. Dried blood spot collection card (with bar code)
5. Sample transport bag with desiccant

Note: Assay plates/reagents provided separately



## POC Kit

6. Single-use POC test cartridge
7. Single-used capillary with plunger for blood transfer to cartridge

# Approach to Biomarker Identification & Selection

- Test known and proposed biomarkers of radiation response in animal models (mouse and NHP)
- Cover a wide range of potential physiological effects – hematopoietic system damage and recovery, tissue damage, inflammatory response, DNA damage, etc.
- Include cell surface proteins in plasma as surrogates for hematology measurements (cell counts)
- Identify radiation responsive markers, and evaluate combinations of markers in a dose-assessment algorithm
- Test dose-assessment algorithm in a blinded study (mouse model)




# Biomarker Evaluation


- Protein biomarkers were tested in mouse and NHP models of radiation response (> 50 candidates)
- For evaluation, biomarkers tested in plasma samples
- Biomarkers represented different categories of response:
  - DNA damage
  - Inflammatory response
  - Tissue damage (organ-specific)
  - Hematopoietic system injury and recovery
  - Protein surrogates for hematology measurements
- Selected biomarkers with clear radiation response, and minimal response to confounding injury (wound) in the mouse model

# Example of Biomarker Response: Flt-3 Ligand

Mouse model



NHP model




Dose  
equivalence  
for different  
species ➡

Species/Strain	LD <sub>50/30</sub> or LD <sub>50/60</sub>		Estimated Critical Threshold for ARS <sup>4</sup>
	Measured <sup>1</sup>	Ratio to Human	
Female B6D2F1/J Mice	9.3 Gy <sup>2</sup>	~2.5	~ 5 Gy
Rhesus Monkey	4.4-6.7 Gy <sup>3</sup>	~1.4	~ 3 Gy
Human	3-5 Gy <sup>3</sup>	1	2 Gy

<sup>1</sup>. X-ray or γ-ray with minimal care (mice, humans, rhesus) and mixed neutron and γ-ray with minimal or normal care (humans)  
<sup>2</sup>. Ledney et al. (2010) Health Physics 98: 145-152.  
<sup>3</sup>. DiCarlo AL et al. (2011) Disaster Med Pub Health Prep 5: S32-S44.  
<sup>4</sup>. ARS = Acute Radiation Syndrome

# Dose Assessment of Blinded Samples (Mouse Study)

- All mouse samples were blinded during testing.
- The predicted dose is shown as a function of the actual dose.
- Points for which the predicted dose exactly matches the actual dose will fall on the solid line.
- Points within the dashed lines meet our dose prediction accuracy criteria and are within 1.5 Gy of the actual dose (below 6 Gy) or within 25% of the actual dose (above 6 Gy).
- The inset shows the percentage of the predicted doses that fall within our accuracy criteria and the root mean square error in the predicted doses across the data set.




*The panel provides excellent performance in assessing dose for blinded samples.*

# Algorithm Verification: Dose Classification of Blinded Samples

- Blinded study data was analyzed for the ability of the algorithm to classify dosage above or below the critical 2 Gy dose threshold in humans (~5 Gy in mice).
- ROC curve for distinguishing doses  $\geq 6$  Gy from non-irradiated controls is in blue.
- ROC curve for distinguishing doses  $\geq 6$  Gy from doses  $\leq 3$  Gy is in green.
- Inset (histogram) shows distribution of the predicted doses for samples receiving 0 Gy, 3 Gy or 6 Gy doses and demonstrates the separation of these distributions.
- Classification performance at the optimal predicted dose thresholds is provided in the table.

*Near perfect discrimination of doses above and below critical threshold of 5 Gy (~2 Gy in humans)*

Dose Algorithm Employing 6 Biomarker Panel



Classification Performance Using the Optimal Threshold		
Classification criteria	0 Gy vs. $\geq 6$ Gy	$\leq 3$ Gy vs. $\geq 6$ Gy
True positive fraction	1.0	0.992
True negative fraction	1.0	1.0
Prediction accuracy (%)	100	99.59
Area under curve (AUC)	1.000	0.998

# Anticipated Challenges with Clinical Study

## Human clinical samples

- Accidental exposure to high-dose total body irradiation (TBI) is rare, actual dose received is not accurately known, and 1-7 day samples are usually not available.
- Most medical radiation treatments use highly localized radiation with exposure limited to 1-2% of total body. These treatments do not produce the effects seen in TBI.
- TBI is sometimes used in medical practice (e.g. bone marrow transplants in blood cancers), but always in combination with chemotherapy. Very few clinical protocols treat with radiation prior to chemotherapy.
- When used medically, TBI is always fractionated into smaller doses to avoid radiation syndromes (e.g. 2 Gy fractions twice daily, or 1.25 Gy three times daily).

*Development and validation of a biodosimetry test has to rely heavily on experimental animal studies, where the radiation dose can be carefully controlled and samples collected over a time course. FDA guidance on use of the Animal Rule for drugs/biologics can provide a starting point for designing studies.*

# A Potential Clinical/Regulatory Approach: Animal and Human Studies

## Mouse Model: Developing and Validating Dose Algorithm, Testing Key Variables

- Identification of biomarkers
- Radiation dose and time response
- Algorithm development
- Effect of combined injury
- Radiation type (photon vs. mixed photon/neutron) and dose rate
- Partial body exposure vs. total body
- Age and gender study

## NHP Model (Rhesus): Validating Biomarkers and Algorithm

- Dose studies for algorithm development and validation
- Studies on the effect of therapeutics (G-CSF) on biomarker levels

## Humans: Connecting Results from Animal Studies to Humans

- Evaluate algorithm performance in patients receiving TBI for stem cell transfer
- If available, confirm in samples from TBI in radiation accidents (extremely rare)
- Demonstrate specificity with samples from individuals representative of US population



# An Alternative Approach to Radiation Diagnostics?

## Biodosimetry (Estimation of Received Dose)

- Advantages
  - Single easily interpretable output
  - Direct connection to physical dosimetry measurements
- Disadvantages
  - Not directly linked to patient health status
  - May not accurately reflect physiological injury (e.g. partial body exposure)
  - Difficult to achieve regulatory clearance: clinical trials are not feasible

## Diagnosis of Specific Radiation Syndromes (Hematopoietic, GI, etc.)

- Advantages
  - Directly tied to patient health status
  - Provides an indication for treatment (e.g., G-CSF for hematopoietic syndrome)
  - Provides alternative approaches and patient populations for clinical trials (e.g., non-radiation related causes of myelosuppression)
- Disadvantages
  - No direct information on radiation dose
  - Different syndromes require different panels of biomarkers

## A Possible Alternative Approach: Syndrome-Specific Diagnostics

### Hematopoietic Syndrome and Other Myelosuppressive Conditions

- The primary radiation pathology in the critical range for treatment (~2 to 6 Gy) is drop in blood cell counts from loss of bone marrow activity
- Other conditions associated with severe myelosuppression:
  - Side effect of chemotherapy for cancer treatment
  - Intended effect of myeloablative chemotherapy prior to stem cell replacement for blood cancers or other disorders
- Alternative strategy for clearing diagnostics for hematopoietic syndrome: Demonstrate ability of diagnostic to detect myeloablation associated with chemotherapy as a surrogate for radiation
  - Much larger accessible patient pool for clinical studies
  - Potential commercial applications in conventional medical care

# Acknowledgements

## MSD

- PI's: George Sigal and Eli Glezer
- Assays: John Kenten, Anu Mathew, Navaratnam Manjula, Galina Nikolenko, and John Casper

- Algorithm Development: Sudeep Kumar, Srikant Vaithilingam

Armed Forces Radiobiology Research Institute (AFRRI)

- Natalia Ossetrova and William Blakely

National Cancer Institute (NCI)

- Kevin Camphausen

University of Maryland, School of Medicine

- Thomas MacVittie

Wayne State University

- Sorin Draghici (biostatistics)

## BARDA

- Donna Boston, Rodney Wallace, Lynne Wathen, Paul Eder, Nina El-Badry