

Laser Sintered Resorbable PCL Splints for Treating Tracheobronchomalacia (TBM)

**Scott J. Hollister^{1,2,3}, Colleen L. Flanagan¹, David A
Zopf⁴, Robert J. Morrison⁴, Richard G. Ohye⁵,
Glenn E. Green⁴**

¹Department of Biomedical Engineering

²Department of Mechanical Engineering

³Department of Surgery

⁴Department of Otolaryngology

⁵Department of Cardiac Surgery

The University of Michigan




Outline

- **Tracheobronchomalacia**
- **Tracheal Splint Clinical Goals & Design**
- **Laser Sintering PCL Splints**
- **Clinical Use and Outcomes**
- **Quality Control: Current & Future**

Tracheobronchomalacia (TBM)

Tracheobronchomalacia (TBM) in Humans

- Compression of airway, typically by malformed vascular structures
- Complete collapse on expiration
- Currently treated by tracheostomy/ventilators 1-2 years
- Significant complications, including death
- Need for patient specific implants due to different defect geometry (length, diameter, number)
- Stents have failed in children; FDA warning metallic tracheal stents
- Implanted splints external to airway found to give better results, but produced in an ad hoc manner





Tracheal Splint Clinical Goals and Design

Clinical Design Goals: Implanted Splint External to Airway


Mechanical Requirement: M; Biomaterial Requirement: B; Surgical Requirement: S

- The splint should provide radial compressive mechanical support to keep the airway open and patent: M/B – 0.12 MPa artery; .01 MPa exhalation
- The splint should provide this radial mechanical support for a period of 24-30 months to allow tracheal remodeling and development: M/B
- The splint should allow transverse and bending displacement, not interfering with cervical motion: M
- The splint should allow growth and expansion of the tracheobronchial complex during this 24-30 month period: M - estimated 15N growth force
- The splint should not cause adverse tissue reaction or remodeling: B/M - Biocompatible
- The splint should not interfere with the mucociliary architecture with the trachealbronchial lumen; it therefore should be placed externally: B/S
- Second surgical procedure should be avoided to remove the splint; therefore, the splint should be bioresorbable: S/M – resorbable in 3 years
- Surgical placement of the splint and attachment of the tracheobronchus into the splint should be straightforward: S; suture holes in splint to “sling” airway
- Patient Specific to account for different malacic airway diameter/length: S/M



Patient Specific Image-Based Design for Splint



3rd Patient
(Only blue
splint used)




2nd Patient
(1 of 2 splints)*



- MATLAB program to automatically generate bellow design w suture holes
- Design variables: inner diameter, open angle, spiral angle, bellow height, wall thickness, suture hole width, etc (> 3,000,000 design perturbations)
- Input parameters from CT measurements from MIMICS Digital Model
- Fit splint to patient model in MIMICS
- Perform finite element analysis: compression, bending, opening (growth)
- Complex patient specific design requires 3D printing

Laser Sintering PCL Splints

Design and Manufacture Process: Outline




Scaffold/Implant Manufacturing by 3D Printing


- Modular Image-Designed Scaffolds fabricated by laser sintering

Complete Video at

<http://www.mottchildren.org/news/archive/201403/babys-life-saved-after-3d-printed-devices-were-implanted-u>



Modular
Image-
Designed
Scaffold



PCL Laser Sintering



Final Manufactured
Scaffold

Materials and Equipment

- EOS P100 Laser Sintering System (www.eos.info/en)
- CAPA 6501 Polycaprolactone (PCL) purchased from Polysciences (www.polysciences.com) **Target Mw = 50kDa**
- Hydroxyapatite (HA) Plasma Biotal (www.plasma-biotal.com)
- Need to Cryogenically Mill Resorbable Polymers (PCL, PLA) Jet Pulverizer (www.jetpulverizer.com); Fraunhofer (<http://www.umsicht.fraunhofer.de/en.html>) ; Evonik (<http://north-america.evonik.com>); Target Particle Size Range: $25\mu\text{m} < x < 125\mu\text{m}$; Median 40-60 μm

References:

Partee et al., 2006, J Man Sci Eng, 128:531-540

Williams et al., 2005, Biomaterials, 26:4817-4827

Eshraghi, Das, 2010, Acta Biomaterialia; 6: 2467-2476

Eshraghi, Das, 2012, Acta Biomaterialia; 8: 3138-3143

Lohfield et al., 2012, Acta Biomaterialia; 8:3446-3456

Eosoly et al., 2010, Acta Biomaterialia; 6:2511-2517


PCL Laser Sintering Parameters

- Important PCL Laser Sintering Parameters: Bed Temperature, Laser Power, Laser Scanning Speed, Scan Spacing, Hatch Spacing, Beam Offset
Our parameters established at University of Michigan: Partee, Hollister, Das. (2006) J Mfg Sci Eng, 128:531-540
- Laser Power: 1 - 5.4 Watts; Typically 4 Watts (UM)
- Bed Temperature: 38 – 56°C; Typically 50-56°C (UM)
- Laser Scanning Speed: 900 – 1800 mm/s; Typically 1000-1500mm/s (UM)
- Scan Spacing: .07 - .2mm; Typically 0.15 – 0.2mm (UM)

References (see prev slide): Eshragi/Das (2010/2012); Lohfield (2012); Eosoly (2010; 2012); Partee (2006); Williams (2005)

Clinical Use and Outcomes

Design & Implantation of Patient Specific Splints




Patient 1:

Left Bronchus;
IRB Approval,
Emergency through FDA
NEJM (2013), 368:2043-2045.
31 months post-surgery



Patient 2:

Bilateral Bronchi;
IRB Approval,
Emergency through FDA
8 months post-surgery



Patient 3:

Left Bronchus;
IRB Approval,
Emergency through FDA
6 months post-surgery

Pre-Op and Post-OP Patency

Pre-Op




Post-Op




Patient 1:
Left Bronchus;
Exhalation Scans

Pre-Op




Post-Op




Patient 2:
Bilateral Bronchi;
Exhalation Scans

Pre-Op




Post-Op




Patient 3:
Bronchoscopy

Bronchial Growth in Patients




Hydraulic Diameter Measures Averaged along Bronchus in MIMICS



All Patients Pre- and Post-Op



**Patient 1 –
Pre-Op**



**Patient 1 –
2nd Birthday**



**Patient 2 – Pre-op
16 months**



**Patient 2 –
First time
sitting up**



**Patient 3 –
Pre-Op**



**Patient 3 –
Post-Op 2
months**

Quality Control: Current & Future


Quality Control Checks for Each Build

- Powder: Check particle size range; Powder Visual Inspection; Humidity Solid Hygrometer Should be 10% to 35% relative humidity
- Build: Check for errors on build log; visual inspection for part dragging; visual inspection for sintered “islands” when unpack build; stair stepping on parts
- Geometry: Caliper Measures (current); Micro-CT to assess part geometry/density (implementing)
- Mechanical Properties: Standard cylindrical test specimens for modulus; splint specimens opening, compression, bending geometric stiffness (implementing)


Geometry Quality Control

- For topology optimized (optimized for stiffness/permeability) microstructures, compare designed vs manufactured geometry by microCT (implementing for splint)

Design/Fabrication Process




Δ design Struts Areas (mm²)




■ Designed
 ■ Actual in Z dir.
 ■ Actual in X/Y plane

2uc: 14-50%;
 3uc: 4-30%;
 4uc: .3-20%

Throats Areas (mm²)



Design to Fabricated Strut/Throat Comparison



Fidelity depends on Unit Cell Size & thus
 Feature resolution; Dias et al, (2014), 36:448-457




Mechanical Testing Quality Control - General

Solid Test Cylinder Modulus

1. Affected by Laser scanning \longrightarrow
parameters: Bed Temp, Laser Power,
Scan Speed 1200-2500 mm/s
2. Anisotropic due to layering
Ex = 295.5 ± 4.4 MPa parallel to bed
Ey = 292.7 ± 9.9 MPa parallel to bed
Ez = 311.7 ± 1.2 Mpa laser direction

Optimized Microstructures


1. Topology optimized for desired stiffness/permeability
2. Compare FE idealized to laser sintered mechanically tested
3. Correlation deviates from 1 to 1 as feature sizes get smaller ($< 0.8\text{mm}$)



Mechanical Testing Quality Control - Splint

Compression:

Simulate
exhalation
loading



Opening:

Simulate growth and
inhalation loading



$$Ku = f; \Rightarrow$$

$$K = \frac{f(\text{artery} / \text{exhalation pressure} * \text{length} * 1\text{mm}; \frac{N}{\text{mm}^2} * \text{mm} * \text{mm})}{u(\text{target} \langle .1\text{comp}; .2\text{open} \rangle * \text{inner diameter}; \text{mm})}$$

$$\text{Stiffness in compression } (.12N * \text{length} / .1 * ID) > \approx 10N / \text{mm}$$

$$\text{Stiffness in opening } (.12N * \text{length} / .2 * \text{OpenAngle}) \leq \approx 2N / \text{mm}$$

Design Target:

1. Withstand arterial compression & respiration pressure
2. Allow growth




Patient 1 - Compression: 128.6 ± 11.8 N/mm; Opening: 2.77 ± 0.26 N/mm;

Patient 2 - Compression: 72.2 ± 14.6 N/mm (11mm); Opening: 1.43 ± 0.12 N/mm;
 195.8 ± 16.2 N/mm (23mm); Opening: 2.43 ± 0.15 N/mm;

Pig Preclinical: Compression: 28.5 ± 1.6 N/mm; Opening: $.43 \pm 0.05$ N/mm;
20% growth over 8 months

Fatigue & Degradation Quality Control

- For resorbable materials, need to determine affect of sintering on fatigue & degradation
- Sintering doesn't significantly change/degrade PCL molecular weight prior to implantation; ~40% loss of Mw by 18 months *in vivo* (spine cage in pig).
- Fatigue properties depend significantly on geometry; Have run spine cages to 5 million cycles in dry environment – need to test in solution



Conclusions

- Developed Laser Sintered, resorbable PCL patient specific splint for treating tracheobronchomalacia; Successful in 3 patients up to 31 months
- Fabricated topology optimized scaffolds with complex microstructure
- Splints with 0.4 to 2.8 N/mm opening stiffness allowed growth in patients and preclinical pig model; 28 to 195 N/mm compression stiffness protects malacic airway
- Laser parameters (scan speed, bed temp, scan power, particle size) significantly affect device geometry, mechanical properties (stiffness, strength, fatigue) & degradation (need to be tested)
- Ability to meet geometric and mechanical requirements depends on how close feature size is to minimum resolvable sintering feature -> closer to minimum feature size will mean larger deviation between design & actual properties

Acknowledgments

- Glenn Green, MD collaborator on tracheal splint; tracheal splint surgery
- David Zopf, MD; Robert Morrison, MD; tracheal splint
- Will Giannobile, DDS, PhD – Periodontal Scaffold
- Giulio Rasperini, DDS – Periodontal scaffold
- Richard Ohye, MD – tracheal splint (surgery)
- Marc Nelson, MD – tracheal splint(referring)
- Chia-Ying Lin, PhD cervical spine fusion
- Colleen Flanagan, MSE - Splint SLS manufacture
- Annie Mitsak, PhD; Eiji Saito PhD
- Matthew Wheeler, PhD collaborator on large animal models

ADSC

- Jonathon Mosley, Dr. Chanaka Rabel, Aaron Maki – Animal Surgeries
- Supported by the NIH and a MICHR grant