

Clinical aspects of insulin bolus calculators and how patients use them

Adam Brown

The diaTribe Foundation (www.diaTribe.org)

*FDA Public Meeting: Regulatory Science Considerations
for Software Used in Diabetes Management*

November 13, 2014



Thank you!

- Thank you so much for the opportunity to bring a patient perspective to this important meeting!
- We encourage you to think of patients as **partners** to help your team better understand the risks/benefits of all diabetes devices – we're here to help!

diaTribe®

Presentation Outline




- **Patient use of bolus calculators**
- The risks of bolus calculators
- Regulatory considerations
- Needed innovation in insulin bolus calculators



**“I believe that the bolus
calculator is one of the
biggest leaps pump
technology has ever made.”**

diaTribe®

Why are BCs such a huge advance for patients?



diaTribe®

Why are BCs such a huge advance for patients?

A “**simple**” calculation: correction only

BG=165 Target = 100 ISF: 1:35

Correction: $165 - 100 = 65 / 35 = 1.86$ units

Total dose: 1.86 units

diaTribe®

Why are BCs such a huge advance for patients?

A “**medium**” calculation: correction + carbohydrates

BG=165 Target = 100 Carbs: 45 ISF: 1:35 ICR: 1:12

Correction: $165 - 100 = 65 / 35 = 1.86$ units

Carbohydrates: $45 \text{ carbs} / 12 = 3.75$ units

Total dose: $1.86 + 3.75 = 5.61$ units

diaTribe®

Why are BCs such a huge advance for patients?

A “**hard**” calculation: correction + carbohydrates + IOB

BG=165 Target = 100 Carbs: 45 ISF: 1:35 ICR: 1:12

Last bolus taken: 3 units, 3 hours ago

Duration of insulin action (DIA): 4 hours

Correction: $165 - 100 = 65 / 35 = 1.86$ units

Carbohydrates: $45 \text{ carbs} / 12 = 3.75$ units

IOB: $3 \text{ units} / 4 \text{ hours} = 0.75$ units left

Total dose: $1.86 + 3.75 - 0.75 = 4.86$ units

diaTribe®


Why are BCs such a huge advance for patients?

They make the
math *much* easier
and less subject to
human error!


diaTribe®

However, most patients don't have access to a bolus calculator

Type 1 - Pump Penetration



Basal-Bolus Type 2s - Pump Penetration



diaTribe®

However, most patients don't have access to a bolus calculator

- Roche's Accu-Chek Aviva Expert is the only standalone US meter with a built-in bolus calculator
- Launched Sept 16, 2014



diaTribe®

Presentation Outline

- Patient use of bolus calculators
- **The risks of bolus calculators**
- Regulatory considerations
- Needed innovation in insulin bolus calculators

1. Bolus calculators *reduce* risk for patients and improve quality of life

Patients without access to a bolus calculator:

Guess/
“ballpark”

Use “mental”
math

No sense of
insulin on
board!

Feel need to eat
standardized
meals

The real-world impact of “ballparking”

	Actual Scenario	Ballpark Scenario
Current BG Target Last Bolus	165 mg/dl 100 mg/dl 6 units, 3 hours ago	165 mg/dl 100 mg/dl 6 units, 3 hours ago
Settings	ISF: 1:35 , ICR: 1:12 DIA: 4 hrs	
Carbs / insulin needed	45 g / 3.75 units	~60 g / 5 units
Correction	1.86 units	~2 units
IOB	1.5 units remaining	No IOB
Total Dose	4.11 units	7 units

2. The risks of bolus calculators often trace to poor settings

- ICR, ISF, target, Duration of Insulin action
 - Often end in 5 or 0 for easier math
 - Setup is usually based on very general formulas that don't account for individual variance
 - Rarely change over course of day, though many may benefit from this

Most of these settings are sub-optimal, and the feedback loops are slow/non-existent to improve and optimize them

3. Risky bolus calculations can be driven by patient errors – Carb Counting



3. Risky bolus calculations can be driven by patient errors – Carb Counting



3. Risky bolus calculations can be driven by patient errors – Not washing hands

“With a sample of 0.3 μl ,
1 μg of glucose (the
weight of a dust particle)
will increase the blood
glucose by 300 mg/dl.”

- Ginsberg, *J Diabetes Sci Technol* 2009

4. Bolus calculators CANNOT be made 100% fool proof – diabetes is too complicated

22+ Variables Impacting Blood Sugar!



An illustrative calculation

A “**real**” calculation

BG=165 Target = 100 Carbs: 75 ISF: 1:35 ICR: 1:12

Last bolus taken: 3 units, 3 hours ago

Duration of insulin action (DIA): 4 hours

6 hours of sleep last night, woke up multiple times

Less than normal level of activity yesterday



45 carbs includes high fat

Morning bolus

CGM trend arrow: ↑

Total dose: ????

diaTribe®



**“Insulin pumps precisely
deliver inaccurate amounts
of insulin.”**

diaTribe®

5. The risks of bolus calculators *are much lower* than other aspects of diabetes care



Presentation Outline

- Patient use of bolus calculators
- The risks of bolus calculators
- **Regulatory considerations**
- Needed innovation in insulin bolus calculators

Regulatory Focus: Presence vs. Absence

**Contains Bolus
Calculator**
Carries addressable
risk


**No Bolus
Calculator:**
Carries Clear Risk

diaTribe®

Recommendations: transparency, simple human factors

- Ideally, patients would have access to clear labeling as to how a particular bolus calculator works (one page, available online)
- Simple software validation and very basic human factors
- **Clinical trials and standardization are likely not necessary:** innovation might be hampered for large and small organizations and patients. This is at a time when the diabetes community needs innovation most.

Shared responsibility – no single party should be held liable for 100% safe operation




- Bolus calculators will never be perfect
- Industry, HCPs, patients, and the FDA should jointly **share** the responsibility for proper operation
- Zero risk is probably impossible

Presentation Outline

- Patient use of bolus calculators
- The risks of bolus calculators
- Regulatory considerations
- **Needed innovation in insulin bolus calculators**

1. More standalone meters and validated apps/software with bolus calculation

- More patients should have access to bolus calculators
- The proper balance between safety/innovation can ensure that happens



diaTribe®

2. Software to optimize bolus calculator settings

BACK **Timed Settings** **ADD**

Time	BASAL	CORRECT	CARB	TARGET BG
12:00 AM	0.75	1:50	1:15	140
3:00 AM	0.80	1:50	1:15	120
7:00 AM	0.70	1:55	1:12	100


Navigation buttons: Up arrow, Down arrow

3. Addition of CGM trend information

*RT-CGM users with Type 1 and Type 2 use rate of change arrows to make **significant changes** in their mealtime and correctional insulin doses, much larger than the 10-20% adjustments commonly recommended.*

Trend Arrow	Mean Patient Adjustment to Bolus Dose*	Example Bolus Dose*
→	---	3.0 units
↑	+111%	6.3 units
↑↑	+140%	7.2 units
↓	-41%	1.8 units
↓↓	-47%	1.5 units

4. Treat-to-range controller after meals – a patient safety net!



diaTribe®

Key Takeaway Considerations

- Bolus calculation is filled with complexity and math, even for the most engaged patients.
- Most patients do not currently have access to a pump or meter with a bolus calculator.
- Bolus calculators *reduce* risk for all types of patients – including those with literacy/numeracy issues – and are much less risky than other things patients do every day.
- There is no such thing as a perfectly safe bolus calculator – diabetes is too complicated

Key Takeaway Considerations

- Ideally, bolus calculators would have clear labeling and human factors to verify that they perform as intended.
- We need more innovation to get bolus calculators into the hands of more patients (standalone meters, software/apps)
- Ultimately, bolus calculators could become more accurate and safer with CGM trend information and treat-to-range control.

diaTribe®

Thank you!

adam.brown@diaTribe.org

www.diaTribe.org

@diaTribenews

*Committed to improving the lives of
people living with diabetes and
prediabetes and advocating for action.*