

CENTRE HOSPITALIER INTERCOMMUNAL
TOULON - LA SEYNE SUR MER



Benefits, risks, and outcomes for fully closed loop ventilation

Dr Jean-Michel Arnal
jean-michel@arnal.org

Intensive Care Unit
Home ventilation unit
Hôpital Sainte Musse
Toulon, France



Medical research manager
Hamilton Medical
Bonaduz
Switzerland





IntelliVent-ASV: passive patient

Patient

Additions

Modes

40
24 Ppeak
cmH2O

25.0
3.0
7.9 ExpMinVol
l/min

530
471 VTE
ml

35
17 fTotal
b/min

0.55 RCexp
s

 7 / 12




8.9 VT/IBW
ml/kg

0 VLeak
%

0 VLeak
ml

0.00 MVLeak
l/min

49 Oxygen
%



Adult Female
160 cm
IBW = 53 kg





Diagram of a human torso showing the lungs and heart. The lungs are depicted with a grid pattern, and the heart is shown in red.

Rinsp	Cstat	PetCO2	SpO2	pulse	HLI
14	44.3	38	97	110	1
cm H2O/A/s	ml/cm H2O	mm Hg	%	l/min	%



 Trend

 View 2/4

 IntelliVent

129 %
%MinVol

12 cmH2O
PEEP/CPAP

46 %
Oxygen

Controls

Alarms

IntelliVent-ASV: spontaneous breathing patients

Patient

Additions

Modes

40

27

Ppeak
cmH2O

20.0
4.0

7.9

ExpMinVol
l/min

750

636

VTE
ml




35

14

fTotal
b/min

0.51

RCexp
s



Trend

View 2/4

IntelliCuff

109

%

%MinVol

6

cmH2O

PEEP/CPAP

29

%

Oxygen

Controls

Alarms

Benefits of IntelliVent-ASV

Safety and efficacy

- **Safety:**
 - Number of safety events
 - Number of switch to other modes
 - % time spent in non desirable ranges
- **Efficacy:**
 - % time spent in optimal ranges


Underlying condition	Resp muscles	“Desired” values
None	Passive	$92 \leq \text{SpO}_2 < 96\%$
ALI/ARDS	Passive	$88 \leq \text{SpO}_2 < 96\%$
None or ALI/ARDS	Active	$92 \leq \text{SpO}_2 < 96\%$
Brain injury	Passive or active	$95 \leq \text{SpO}_2 < 99\%$
Chronic hyperCO ₂	Passive or active	$88 \leq \text{SpO}_2 < 94\%$
None	Passive	$30 \leq \text{EtCO}_2 < 41 \text{ mmHg}$
ALI/ARDS	Passive	$30 \leq \text{EtCO}_2 < 45 \text{ mmHg}$
None or ALI/ARDS	Active	$25 \leq \text{EtCO}_2 < 45 \text{ mmHg}$
Brain injury	Passive or active	$30 \leq \text{EtCO}_2 < 33 \text{ mmHg}$
Chronic hyperCO ₂	Passive	$40 \leq \text{EtCO}_2 < 45 \text{ mmHg}$
Chronic hyperCO ₂	Active	$40 \leq \text{EtCO}_2 < 55 \text{ mmHg}$
None	Passive	$6 \leq \text{TV} < 10 \text{ mL/Kg PBW}$
ALI/ARDS	Passive	$5 \leq \text{TV} < 8 \text{ mL/kg PBW}$
None or ALI/ARDS	Active	$6 \leq \text{TV} < 10 \text{ mL/Kg PBW}$
Brain injury	Passive or active	$6 \leq \text{TV} < 10 \text{ mL/Kg PBW}$
Chronic hyperCO ₂	Passive or active	$7 \leq \text{TV} < 10 \text{ mL/kg PBW}$
None	Passive	$6 \leq \text{Pplat} < 25 \text{ cmH}_2\text{O}$
ALI/ARDS	Passive	$6 \leq \text{Pplat} < 30 \text{ cmH}_2\text{O}$
None or ALI/ARDS	Active	$6 \leq \text{Pplat} < 25 \text{ cmH}_2\text{O}$
Brain injury	Passive or active	$6 < \text{Pplat} < 25 \text{ cmH}_2\text{O}$
Chronic hyperCO ₂	Passive or active	$6 < \text{Pplat} < 25 \text{ cmH}_2\text{O}$
None or ALI/ARDS	Active	$12 < \text{RR} < 30 \text{ b/min}$
Brain injury	Active	$15 < \text{RR} < 35 \text{ b/min}$
Chronic hyperCO ₂	Active	$15 < \text{RR} < 35 \text{ b/min}$

Benefits of IntelliVent-ASV

Safety and efficacy: passive ICU patients

Randomized cross-over study:

- 50 passive ICU patients: 19 normal lungs and 31 ARDS patients
- Two periods of 2 hours
- Control period: ASV



Results:

- No safety event, no switch to other mode

Benefits of IntelliVent-ASV

Safety and efficacy: passive ICU patients



n= 50

Parameter	ASV	IntelliVent-ASV®	<i>p</i>
MV (L/min)	7.6 (6.5–9.5)	6.8 (6.0–8.0)	<0.001
V_T /PBW (mL/kg)	8.3 (7.8–9.0)	8.1 (7.7–8.6)	0.003
RR (breath/min)	15 (14–17)	14 (13–17)	0.004
P_{INSP} (cmH ₂ O)	28 (24–33)	25 (22–29)	<0.001
P_{PLAT} (cmH ₂ O)	24 (20–29)	20 (19–25)	0.005
PEEP (cmH ₂ O)	10 (6–14)	8 (5–10)	0.011
FiO ₂ (%)	40 (30–50)	30 (30–39)	<0.001
C_{STAT} (mL/cmH ₂ O)	37 (31–48)	37 (29–44)	0.935
R_{INS} (cmH ₂ O s/L)	16 (14–18)	17 (14–19)	0.699
RC_{EXP} (s)	0.7 (0.6–0.8)	0.7 (0.6–0.8)	0.326
pH	7.3 (7.3–7.4)	7.3 (7.2–7.4)	0.104
PaO ₂ (mmHg)	92 (81–124)	84 (75–104)	0.052
PaO ₂ /FiO ₂ (mmHg)	240 (163–318)	259 (197–323)	0.117
PaCO ₂ (mmHg)	37 (34–42)	37 (33–49)	0.026
SaO ₂ (%)	97 (95–98)	96 (93–98)	0.028
Dead space (mL)	144 (99–224)	134 (85–209)	0.009



Benefits of IntelliVent-ASV

Safety and efficacy: passive ICU patients

% of breaths



% of breaths



n= 50




Benefits of IntelliVent-ASV

Safety and efficacy: spont ICU patients

Randomized cross-over study:

- 14 ICU patients in spontaneous ventilation
- Two periods of 24 hours
- Control period PS




Results:

- No safety event, no switch to other mode

Benefits of IntelliVent-ASV

Safety and efficacy: spont ICU patients



Benefits of IntelliVent-ASV

Safety and efficacy: post cardiac surgery

Randomized controlled trial:


- 60 patients after elective non complicated cardiac surgery
- Over 4 hours
- Control group: protocolized SIMV + PS

Results:

- No safety event, no switch to other mode

Benefits of IntelliVent-ASV

Safety and efficacy: post cardiac surgery





n= 60

Lellouche. Intensive care Med 2013

Benefits of IntelliVent-ASV

Safety and efficacy: post cardiac surgery




Note, conventional patients set at 10mk/kg/PBW

n= 60

Benefits of IntelliVent-ASV

Safety and efficacy studies

- ICU patients: RCT 80 patients



Bialais. Critical Care 2013 [abstract]

- Pediatric patients: sequential study

	PSV	ASV	IntelliVent	PSV
n= 10	(1 hr)	(1hr)	(1hr)	(1hr)
Time in normal breathing range(%)	82 ± 31	91 ± 10	94± 2	
number of ventilator setting modifications (n)	0 (0-0)	0 (0-0)	77 (57-120)	0 (0-0)

Jouvet . Critical Care 2012

Benefits of IntelliVent-ASV

Feasibility study in ICU

Prospective observational comparative study:

- 100 unselected ICU patients with expected duration MV > 24 hours
- Intubation to weaning or death , MV duration 3 (2-7) days

Results:

- No safety event, no switch to other mode

Benefits of IntelliVent-ASV

Feasibility study in ICU

100 patients
392 days
of ventilation
analyzed

%MinVol

Automatic

Manual

2 patients for one day
Large CO₂ gradient

Oxygen

Automatic

Manual

7 patients for one day
Poor SpO₂ signal

PEEP/CPAP

Automatic

Manual

3 patients
Increased PEEPi, PNO, P_{ESO}
Arnal. Critical Care 2013

Benefits of IntelliVent-ASV

Prevent hypoxemia

Oxygenation




Benefits of IntelliVent-ASV

Prevent hyperoxia

1 770 patients with sepsis

5 498 patients



PaO ₂ category (kPa)	n (%)	PaO ₂ (kPa, mean \pm SD)	FiO ₂	PEEP (cmH ₂ O)
<8	4,544 (3.6)	7.3 \pm 0.7	0.5 (0.4–0.61)	10 (7–12)
8.01–10	21,335 (16.8)	9.2 \pm 0.6	0.45 (0.4–0.51)	9 (5–12)
10.01–13	45,545 (35.9)	11.5 \pm 0.8	0.4 (0.4–0.5)	8 (5–10)
13.01–16	27,132 (21.4)	14.4 \pm 0.9	0.4 (0.4–0.42)	5 (5–10)
>16.01	28,222 (22.3)	20.7 \pm 6.0	0.4 (0.4–0.45)	5 (5–10)


Dahl. Acta Anaesth Scand 2015

De Graff. Intensive Care Med 2011

Benefits of IntelliVent-ASV


Prevent hyperoxia

1 770 patients with sepsis



Dahl. Acta Anaesth Scand 2015

3 322 patients



De Jonge. Crit Care 2008

Benefits of IntelliVent-ASV

Prevent hyperoxia

	1 Normal lung	2 ARDS	3 COPD	4 Others	p (ANOVA)
Number of patients	45	16	19	20	
Number of days	139	90	79	84	
PaO ₂ /FiO ₂ (mm Hg)	326 (267-380)	206 (172-252)	260 (206-328)	241 (189-304)	< 0.001
PaO ₂ (mm Hg)	100 (85-117)	76 (69-84)	86 (74-105)	82 (73-95)	< 0.001
pH	7.36 (7.30-7.41)	7.35 (7.23-7.42)	7.30 (7.25-7.35)	7.33 (7.26-7.41)	< 0.001
PaCO ₂ (mm Hg)	35 (31-40)	41 (35-49)	45 (37-53)	39 (34-46)	< 0.001

Prospective observational comparative study:

100 unselected ICU patients

Arnal. Critical Care 2013

Benefits of IntelliVent-ASV

Application of recommendations

	ARDS mild Passive	ARDS moderate Passive	ARDS severe Passive
RC_{EXP} (s)	0.45 (0.43-0.55)	0.47 (0.44-0.57)	0.40 (0.24-0.43)
C_{STAT} (L/cm H ₂ O)	42 (35-51)	43 (38-48)	19 (14-35)
V_T /PBW (mL/kg)	6.5 (5.9-8.0)	6.9 (6.3-7.5)	5.4 (5.0-5.9)
Driving pressure (cmH ₂ O)	10 (8-11)	9 (7 – 10)	10 (8 -12)
P_{PLAT} (cm H ₂ O)	20 (18-24)	20 (17-23)	24 (20-26)


n= 44

Personal data 2015

Benefits of IntelliVent-ASV


Reduce workload

Number of manual intervention / 48h




80 ICU patients

Bialais. Intensive Care Medecine 2013 [abstract]



60 post cardiac surgery patients

Lellouche. Intensive care Med 2013



128 post cardiac surgery patients

Beijers. Intensive care Med 2014

Benefits of IntelliVent-ASV



Reduce workload

60 ICU patients ventilated for more than 48 hours

IntelliVent-ASV versus VAC+PS with protocolized weaning


Manual settings

Automated settings



Benefits of IntelliVent-ASV

Outcomes





— IntelliVent-ASV
— VAC and PS

Arnal [submitted]

Benefits of IntelliVent-ASV

Easiness of use



n= 60 * p < 0.001

Arnal [submitted]

IntelliVent-ASV additional benefits

- Increase consistency of care
- Timely adjustment when mobilization
- Less sedation required
- ... Improved organization



2014-08-08
11:20:50

INTELLiVENT

ASV
Adult

Patient

Additions

Modes



Trend



View
2/4

IntelliGuff

40

18 Ppeak
cmH₂O

25.0
4.0

9.8 ExpMinVol
l/min


920

461 VTE
ml


35

13 fTotal
b/min

0.95 RCexp
s



CO₂ elimination



109

%
MinVol


10

cmH₂O
PEEP/CPAP

48

%
Oxygen

Adult Male
182 cm
IBW = 77 kg



Rinsp	Cstat	PetCO ₂	SpO ₂	pulse	HLI
7	77.2	41	96	81	2
cm H ₂ O/s	ml/cm H ₂ O	mm Hg	%	1/min	%

Oxygenation



10 / 12



7.7 V'alv
l/min

35 VDaw/VTE
%

Controls

Alarms

Monitoring

Graphics

Tools

Events

System

INT AC

IntelliVent-ASV risks

Sensor accuracy

SpO₂

E_TCO₂

n = 100 measurements

Compare SpO₂ to SaO₂

Bias = 0.2 %

Limits of agreement: -0.4 to 0.8 %


Precision= 2.4 ± 0.2 %

	n	PaCO ₂ – E _T CO ₂ (mm Hg)
All patients	100	4 (0 -9)
Passive patients	47	1 (-1 -5)
Active patients	53	7 (1 -13)
Normal lungs all	46	3 (0 -10)
Normal lungs passive	20	1 (-2 -6)
Normal lungs active	26	7 (0 -11)
COPD all	21	7 (0 -14)
COPD passive	9	0 (-5 -2)
COPD active	12	9 (7 -16)
ARDS all	34	3 (-1 -6)
ARDS passive	17	2 (-1 -5)
ARDS active	17	5 (2 -8)

IntelliVent-ASV risks

Safety features


- Maximum pressure limitation




- PEEP manual in hypercapnic and brain injury




- Ranges of PEEP



- Automatic PEEP limitation according to hemodynamic




- ASV automatic limits for V_T and RR



IntelliVent-ASV risks

Not appropriate for some patients?



- Absolute : ventilation with leaks
 - NIV
 - Bronchopleural fistula
- Particular situations:
 - Oxygenation controller and dyshemoglobin: CO poisoning, methemoglobin, monozygotic Sickle cell disease.
 - PEEP controller in case of hemodynamic instability



IntelliVent-ASV risks

Loss of situation awareness

- Monitoring is slightly different
- Alarms are important
- Require a dedicated training



Why do we make mistakes?

- I was tired
- I forgot
- I was distracted
- My attention lapsed
- I was stressed
- I was short of time
- My staff was new
- My staff was inexperienced
- The equipment was difficult to use
- I made an error of judgment




The future

- Will we have more staff?
- Will we have more experienced staff?
- Will our patient be less complex?
- Will we have more time?



Conclusion

- IntelliVent-ASV full closed loop ventilation mode
- Clinician set and adjust targets and limits
- Uses well known physiological measurements
- No safety event reported out of $\approx 30\,000$ patients
- Greatest benefits: safety, lung protection, and organization
- Greatest risk: insufficient user training
- One way to prepare our future challenges



Thank you... 😊