

Electromagnet Compatibility (EMC) for Neurostimulation Devices and Security Systems

Donald Witters

Division of Biomedical Physics
Office of Science and Engineering Laboratories
Center for Devices and Radiological Health
Food and Drug Administration

FDA EMC laboratory research helps ensure safe and effective neurostimulation devices through publication, standards, collaborations, and regulatory guidance and review.



What is EMI/EMC and Why the Concern?

Electromagnetic Interference (EMI)

Disruption of the electrically powered device function by Electromagnetic energy (EM)

- Radiated EM (e.g. radio waves)
- Conducted (e.g. AC power surges)
- Electrostatic discharge (ESD)





Electromagnetic Compatibility (EMC)

Device functions properly in its EM environment (immunity), and the device EM emissions do not cause EMI in other devices

EMI Concerns for Active Neurostimulator Devices





- Implantable-body worn neurostimulators
 - Spinal cord, Deep brain (DBS), Vagus nerve stimulator
- Hazards related to EMI
 - Over/under stimulation
 - Disturb control settings or telemetry communication
- Potential sources of EMI
 - Security systems: WTMD, HHMD, anti-theft (electronic article surveillance (EAS))
 - RFID, radio transmitters, diagnostic imaging (e.g., MRI), medical diathermy
 - Industrial/commercial environments



Leveraging Neurostimulator EMC Research

- Interagency agreement with TSA, DHS for AIT body scanner
- Collaborations with the medical device manufacturers
- Emissions measurements
- Develop test methods & perform exposure testing
- No effects seen on devices in AIT
- Findings are used in regulatory EMC guidance and reviews and national/international standards
 - ISO 14708-3
 - IEC 60601-1-2



Continuing Challenges for Medical Device EMC



- Changing electromagnetic environment, newer security systems & medical devices
- Continue to ensure devices are safe and effective for EMC
- Raise awareness of EMC/EMI
- Engage patients and caregivers
- How can EMC efforts be better engaged for patients?



Questions/Discussion



Donald Witters

Center for Devices and Radiological Health,
Food and Drug Administration

10903 New Hampshire Ave.

Building 62, Room 1130

Silver Spring, MD 20993-0002

Phone: 301/796-2483

donald.witters@fda.hhs.gov

