

1           Having said that, I would say that this is  
2           not a good device to close unoperated post-infarction  
3           VSDs and I wouldn't do it.

4           DR. AZIZ:     What about in the primary  
5           situation?.

6           DR. LOCK:     Post-infarction VSDs that have  
7           not already undergone surgery to fix their coronary  
8           artery disease, our results have not been good and I  
9           don't think this is a good device for that clinical  
10          situation.   The holes are all 10 mm, 12 mm.

11          When you put a device in the septum'  
12          continues to resorb and the infarct gets bigger and  
13          the hole gets bigger.   While you may stabilize them  
14          for 12 to 36 hours, the holes invariably have come  
15          back in the unoperated first five-day post-infarction  
16          VSDs.

17          The successes that we've had, and I don't  
18          know what the number is but it's maybe half, I think  
19          have all been post-operative, or all but one have been  
20          surgery to fix the coronaries, to fix the VSDs, the  
21          defect has recurred and that's when we have gone back  
22          and made those patients better.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. AZIZ: You couldn't see this being used  
2 as a bridge to sort of stabilizing the patient for  
3 five or six days and then going in?

4 DR. LOCK: I think there is a new device  
5 that's in development which is much larger and has  
6 partial self-centering characteristics and might, in  
7 fact, be a very successful device for stabilizing. We  
8 hope to start using that device for post-infarction  
9 VSDs but not device. I'm not going to use this device  
10 for post-infarction VSDs anymore.

11 DR. AZIZ: Thanks.

12 DR. HOPKINS: I'll echo some of the other  
13 panelists. I don't think you see a lot of surgeons  
14 fighting for these patients. I think the major  
15 outcome of significance is really the survival some  
16 six to 12 months after you've had to do something of  
17 which this is a good choice.

18 I am interested about the thoughts about the  
19 post-infarction VSD. I, too, was going to ask about  
20 that. In your indications for use, there's no  
21 specific either indication or contraindication for its  
22 use in that subset of patients.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1           If you feel strongly that it should not be  
2       used, I just wonder. I just throw it out and suggest  
3       that perhaps that should be put in as a  
4       contraindication to its use or, at least, a lack of  
5       indication.

6           A question of there were two devices. As I  
7       read through the various sections it appeared that in  
8       the pivotal series there were two devices which were  
9       explanted at surgery that were not one of the  
10      mortalities. Does anybody know the story on those two  
11      patients or why?

12           DR. JENKINS:       Two were at heart  
13      transplantation for ventricular failure. One was a  
14      failed septation that was taken out at the time of a  
15      Fontan operation. It was basically a failed  
16      procedure. The fourth explant was done in the cath  
17      lab. It was that same patient who had the four  
18      embolizations. One of the devices got taken out late  
19      and that patient ultimately went to the operating  
20      room.

21           DR. HOPKINS:       That was taken out  
22      transcatheter.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. JENKINS: Yes, it, was.

2 DR. HOPKINS: There were two.

3 DR. JENKINS: There were three. Two at  
4 transplant and one at Fontan.

5 DR. HOPKINS: Two surgical.

6 DR. JENKINS: And one at Fontan.

7 DR. HOPKINS: Okay. Thanks. In the summary  
8 of safety and effectiveness, as well as in the  
9 indications for use -- and there have been a number of  
10 references to this. Some references to poor anatomy  
11 as being a contraindication or bad anatomy or  
12 unfavorable anatomy for its use -- and sort of left it  
13 at that in terms of a qualitative sort of statement.

14 Can you provide more precise guidance for  
15 what constitutes bad anatomy for its use or should  
16 that be more specifically part of the training  
17 component? Is there some quantitative approach within  
18 2 mm of the mitral valve, etc.?

19 MS. KULIS: Certainly we can add additional  
20 detail as far as what anatomy is unfavorable.

21 Dr. Jenkins?

22 DR. JENKINS: It will primarily be with

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 relation to the valves. I don't know if one of the  
2 interventionalists could comment on anatomy where it's  
3 just not technically possible to pass a sheath or a  
4 wire through such an extraordinary pathway.

5 DR. HOPKINS: I just want you to know it's  
6 being used now in just a couple of superb centers. As  
7 it spreads out, I'm just wondering if there doesn't  
8 need to be a little bit better guidance for those.

9 MS. KULIS: I'd just like to make one point  
10 as far as you said used in a couple of centers. We  
11 have a total of 30 centers right now in the United  
12 States that do have institutional approval to perform  
13 VSD closures using this device.

14 DR. HOYER: Mark Hoyer again. As far as  
15 location of defects and difficult ones to get to,  
16 obviously I told you we have done three so I don't  
17 have an extensive experience that I'm going to be able  
18 to convince a lot of people but I can tell you that  
19 down at the apex of the heart it can be very  
20 cumbersome.

21 There's a lot of trabeculations in the right  
22 ventricular side of the septum. In fact, the device

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 won't necessarily even open completely so that it's  
2 flat on both sides but it will be darn close.        **But**  
3 you still have accomplished the task of opening the  
4 left ventricular side and then releasing the device as  
5 you open up the right ventricular side before letting  
6 go of it and is in a stable position.        Perfectly  
7 stable. That, again, is a muscular defect much closer  
8 to the apex but well away from semilunar valve or AV  
9 valve.

10                    DR. BOUCEK:    Yes.    I think you're correct  
11 that there are some locations where it is more  
12 difficult in the anterior portion of the septum  
13 sometimes it's difficult to get the sheath to go up  
14 into that portion. These are difficult procedures to  
15 begin with. I think they represent the sort of new  
16 unfortunate era, if you happen to be an interventional  
17 cardiologist,    of    where    pediatric    cardiology  
18 interventions are going.

19                    I think with experience with other types of  
20 complex interventional procedures in pediatrics, it's  
21 just a matter of a problem to be solved rather than an  
22 insurmountable problem.    It tends to be lengthy.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 Sometimes you have to try the sheath from a different  
2 approach rather than from the neck. Maybe from below.  
3 It ends up being problems that need to be surmounted  
4 rather than ones that shouldn't be attempted. They  
5 tend to be long cases. They are like some of the more  
6 complex oblation procedures or some of the more  
7 complex stent procedures that we do in terms of the  
8 duration of time that we're in the cath lab. I  
9 finally understand how much I respect the surgeons for  
10 spending eight hours in the operating room.

11 DR. HOPKINS: Well, don't misunderstand me.  
12 I'm not going to 'sign up to get trained on this  
13 device. I think that, in fact, I am on your side on  
14 this. I want this to succeed as it rolls out. I'm  
15 just concerned about the training. I think we'll  
16 probably talk about training a little bit later, but  
17 that there be a little bit more precision in the  
18 guidance of this.

19 I think, also, knowing these patients and  
20 looking at the study information and also reading  
21 between the lines, these are patients that are being  
22 managed in centers that have full cardiac surgical

1 backup.

2 In the indications for use and guidance  
3 documents, it basically says surgical support should  
4 be readily available. I think that may be more bland  
5 than it needs to be. I think this needs to be done in  
6 centers where it is truly complete support.

7 Also you talk about the' transient  
8 hemodynamic compromises. It sounds to me like the  
9 reason the mortality rate in this extraordinarily  
10 difficult group that you presented being so low is  
11 that they are managed by cardiac anesthesia,  
12 cardiology, interventionists simultaneously.

13 I wonder if there shouldn't be a little more  
14 stronger guidance about that either in the training  
15 document or in the indications for use because this is  
16 not your standard coronary stent that's going in.  
17 You're using a whole team approach here.

18 Like others, I congratulate you.

19 DR. TRACY: Dr. Zahka.

20 DR. ZAHKA: This is certainly a very diverse  
21 group of patients and a very challenging group of  
22 patients. You all deserve congratulations as well.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1           The assessment of them is not always easy as  
2           evidenced by the child with a single ventricle that  
3           was attempted to be septated, and the 12 patients who  
4           were felt to have larger VSDs and turned out to be  
5           small. Did those patients have a band on that made it  
6           impossible to really judge the VSD size, the 12  
7           patients that got enrolled but did not get implants.

8           DR. JENKINS: Had no intent of planting a  
9           device. Part of that is factual just in the way that  
10          we set up the study because we had to have the prior  
11          pier review. There was a lot of paperwork that had to  
12          be done just to have it possible to put a device in at  
13          the time of the procedure.

14          In order to have the procedure go forward in  
15          a timely fashion, we tried to anticipate cases where  
16          it might be necessary even before the hemodynamics had  
17          been done. Obviously everyone is always hoping these  
18          defects go away on their own and they sometimes do.

19          DR. ZAHKA: Does that then reflect our  
20          inability to really assess these people, these  
21          children accurately and how does that speak to the  
22          follow-up data?

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. JENKINS: I think assessment has to be  
2 made in the cath lab once the final pictures are  
3 there. For the band patients it's very difficult for  
4 the echocardiographers to always judge appropriately.  
5 Even for the nonbanded patients I think the angiograms  
6 and the hemodynamics help a lot.

7 I think in this case, though, it's partly an  
8 artifactual reflection that if there was even a small  
9 probability like 15 or 20 percent likelihood we might  
10 want to close a defect. We did peer review of the  
11 patients so then they are counted as enrolled in the  
12 study.

13 DR. ZAHKA: It's also been my sense, in  
14 fact, that infant cardiac surgery has progressed  
15 dramatically over the last 12 years. Although there's  
16 not a lot in the literature about closure of multiple  
17 muscular VSDs and that there are still problems with  
18 that, that this process has, in fact, progressed and  
19 that there are probably more children who could be  
20 done surgically as well.

21 I look at the illustration in the operator's  
22 manual of this ventricular septal defect which looks

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 like it would be good to close by intervention or by  
2 surgery. Perhaps what is the risk benefit of each at  
3 what age.

4 I think about the process you have for  
5 reviewing who should be enrolled in this approach and  
6 notice that you have a surgeon and a cardiologist  
7 review every case beforehand. Is that surgeon and  
a cardiologist also part of Boston Children's Hospital  
9 or are they kind of separated from this whole process?

10 DR. JENKINS: They are within our  
11 institution. The reason we did that was simply for  
12 expediency except for the adults enrolled in the trial  
13 where the peer reviews are done by adult cardiologists  
14 at partnership centers. The peer reviews at all the  
15 centers in the trial, that was similarly the case.

16 I think that some of it is taken as a  
17 success if the surgeons get better partly because of  
18 some of the alternatives that patients have available.  
19 I think in response to Dr. Skorton's earlier question,  
20 I did do a pretty extensive literary review looking  
21 for almost anything that was more recent than what Dr.  
22 Mayer presented.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1           What I found was a series of 11 cases in Dr.  
2       Bovey's paper that was buried between categories where  
3       they weren't really broken down by ventriculotomy.  
4       That group of 11, according to the authors of that  
5       manuscript, it does suggest that maybe some left  
6       ventriculotomies are doing a little better than they  
7       were, you know, 10 or 15 years ago.

8           There was only one other single case report  
9       from the European literature where a large  
10      ventriculotomy was presented as a good outcome short  
11      term. There was a series of letters to the editor  
12      afterwards, you know, kind of worrying about late  
13      results. That was all I found in the literature.

14           If it's true that the surgeons are doing  
15      better, it's not out there where we can review it and  
16      see the results. I'm sure there's going to be a  
17      series of defects that are always difficult to close  
18      surgically, a series that are easy to close surgically  
19      and a series in between where, you know, the cardiac  
20      surgeons will evaluate the best outcome as time  
21      passes.

22           DR. TRACY: Dr. Williams.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. WILLIAMS: Well, my questions will be  
2 related to what is the best way to transfer the  
3 experience at Boston Children's Hospital to other  
4 institutions as they become involved. And if there  
5 should be any limits on the kinds of patients that are  
6 attempted by hospitals earlier in their learning curve  
7 or who have a lower total volume experience with  
8 surgery, echo, and the other factors that are  
9 important to this process.

10 The first one was the illustration showed  
11 passage of the catheters through the simplest kind of  
12 lobe and muscular defect. Then we heard that the  
13 adverse events were more related to technical issues.

14 I have a suspicion that maybe technical  
15 issues were greater in the far interior or far  
16 posterior or apical positions. Were you able to look  
17 at those separately to see if those kinds of defects  
18 had a higher incidence of adverse events than the more  
19 favorable position?

20 DR. JENKINS: We looked at the differences  
21 in outcomes by the post-operative residual defects  
22 versus the congenital defects and we didn't really

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 find any differences in our safety or efficacy  
2 outcomes for those two groups but we never looked by  
3 the specific location in the septum where the defects  
4 were-.

5 DR. WILLIAMS: I wonder if Jim Lock, who has  
6 such large experience with this, has an impression?

7 DR. LOCK: I think Dr. Williams is correct.  
8 One can predict where the trouble will occur from  
9 choosing catheter passage. I do believe that most of  
10 the catheter induced -- most of those five patients  
11 with catheter induced heart block were posterior  
12 muscular VSDs near the tricuspid valve.

13 I do think that the patients with the  
14 catheter induced mitral regurgitation were also  
15 posterior muscular VSDs. That is the particular -- if  
16 you were going to -- I think actually the anterior  
17 septum turns out to be the easiest and the safest  
18 place to fool around.

19 I think if you were going to apical muscular  
20 VSDs, mid-muscular VSDs, intramural VSDs near the  
21 aortic valve are actually pretty safe. I think the  
22 one place where people should be more cautious really

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 in their experience is in the posterior muscular  
2 septum near the attachment to the tricuspid valve.

3 DR. WILLIAMS: Thank you. You might want to  
4 keep your seat because I've got another question  
5 coming up.

6 It seems to me that considering the  
7 difficulty sometimes intellingthe difference between  
a multiple VSDs and a patient who really truly has no  
9 septum but has bundles that are running at different  
10 angles to each other, essentially have no wall but a  
11 collection of bundles, in high referral centers by  
12 echo you often see this as a misdiagnosis from other  
13 centers. I think even in the best of hands it's  
14 possible to miss it. I think probably it was.

15 I would say that probably echo is superior  
16 to angiography in recognizing this lesion if it's done  
17 very, very carefully. I think MRI in some  
18 circumstances can also add some information.

19 My question is really what should be the  
20 experience requirements for the echo cardiography who  
21 is evaluating these patients prior to attempt or prior  
22 to talking to the family about the potential for doing

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 a device closure.

2 And if there might be some way where the  
3 mother institution could produce a teaching tape or a  
4 series of teaching evaluations to show  
5 echocardiographers how to recognize this lesion -- it  
6 ought to be done anyway -- in order to avoid this  
7 particular pitfall.

8 Or how to recognize what you would view as  
9 the higher risk defects and how to recognize that  
10 margin along the posterior -- that posterior margin of  
11 the defect, where you think the pitfalls are so they  
12 are not going to be able to recognize this with their  
13 lower volume and lower experience. Is there a way to  
14 shorten the experience, the learning curve?

15 DR. LOCK: Yes. You're exactly right. I  
16 mean, if you look carefully at the data, we made that  
17 mistake three times. We thought there were three  
18 patients that were septable that probably really  
19 weren't and they had exactly the anatomy that you  
20 describe, and that is that you could sort of talk  
21 yourself into thinking there was a septum but then  
22 when the surgeon goes in, there just is not a septum.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1 We will and have analyzed those three patients and  
2 will continue 'to do so.

3 I think you are right. Sophisticated  
4 echocardiography and probably 3-D reconstruction is a  
5 better way to assess this than angiographically which  
6 was inferior to those two techniques in deciding who  
7 is septable and who isn't. I would agree that is part  
8 of our responsibility.

9 DR. WILLIAMS: And I think that will be part  
10 of the general recommendation on my part that when you  
11 talk about what are the institutional requirements to  
12 carry this out, that it specifically states training  
13 and experience requirements for the echocardiographer  
14 and the cardiac anesthesiologist since the total  
15 outcome is so dependent upon those individuals as well  
16 as the main operator.

17 Could I just ask in the far anterior and far  
18 posterior defects, I recognize that this device is  
19 flexible and soft. It's not likely to impinge on  
20 structures so much. Has there been any indication of  
21 interference with the anterior or posterior descending  
22 coronary artery and would you have recognized it given

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 the kind of surveillance? What would you expect to  
2 have seen if you had encountered that?

3 DR. LOCK: We haven't done selective  
4 coronaries in any of- the patients. The only thing  
5 that I tried to do, and I'm not sure this is an  
6 adequate test, obviously we tried to look at  
7 ventricular performance in all of the patients and  
8 haven't recognized to my knowledge localized  
9 ventricular dysfunction.

10 There's no question that the device can sit  
11 right next to the septum and, therefore, you know, one  
12 of the anterior or posterior descending arteries. We  
13 just haven't seen it.

14 DR. WILLIAMS: Okay. So you haven't seen  
15 segmental wall motion?

16 DR. LOCK: We look pretty carefully for it  
17 because obviously it was one of the clinical concerns  
18 about ventriculotomy patients.

19 DR. JENKINS: We haven't seen signs of  
20 ischemia on the electrocardiograms or things like that  
21 on the surveillance.

22 DR. WILLIAMS: Okay. Great. Do you believe

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 the best use of this device in those patients who have  
2 complex conal truncal abnormalities or pulmonary  
3 artery bands is ultimately to do the catheter closure  
4 after you've attempted to do the surgical closure or  
5 to do the catheter closures of the more difficult  
6 defects in preparation for attempting as a stage  
7 before deciding whether to attempt a complete repair?

8 DR. LOCK: We do it both ways. I think that  
9 if the patient has a band in place, then we tend to  
10 close everything we can close safely in the cath lab.  
11 If the patient doesn't have a band in place rather  
12 than commit the patient to two cardiac operations, the  
13 surgeons decide if they think they can close most, if  
14 not all, of the defects.

15 If they think they can close most, if not  
16 all, the defects using John's requirements without a  
17 left ventriculotomy or without an extensive right  
18 ventriculotomy, then they get the first crack at those  
19 patients. It's really very patient dependent.

20 DR. WILLIAMS: Given the variation of  
21 surgical experience with these lesions, do you  
22 recommend to other institutions that they do it one

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 way before doing it another?

2 DR. LOCK: I think the safest technique now  
3 is bands for people with multiple muscular VSDs.

4 DR. WILLIAMS: But rather if you anticipate  
5 you might need to do both, which one to do first for  
6 those institutions that may have variable surgical  
7 experience?

8 DR. LOCK: I think the risk of catheter  
9 closure in banded patients is actually pretty small.

10 DR. WILLIAMS: And since the indications of  
11 the catheter closure are so closely related to the  
12 ability of the surgeon to close defects, particularly  
13 if you're going to do the surgery anyway, do you have  
14 any recommendations on the volume experience of the  
15 surgical team or the institution in terms, of surgical  
16 experience knowing that by your studies and others  
17 have been directly related to surgical outcome? I'm  
18 sorry to be asking all these questions.

19 DR. JENKINS: The wrong hat, Roberta. I'm  
20 not sure what specific volume standard, for that would  
21 be or whether a volume standard is the correct  
22 measure. I do know that through the Agency for Health

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 Care and Research that there is going to be a proposed  
2 volume standard of around 100 surgical cases a year  
3 being dictated to pediatric cardiology based on  
4 relatively little information. Whether that would  
5 apply to a specific patient with complex ventricle  
6 septum I think would be hard to say.

7 I think at this point one would need to  
8 emphasize that if the surgeon is wrong and they can't  
9 close these multiple defects safely, that the patient  
10 is likely to be very sick and the patients where we  
11 did it in the opposite direction and the VSD was left  
12 are often taken to the cath lab for a VSD closure on  
13 a fairly urgent basis.

14 I think in those cases where people were  
15 less certain about what they could do, it would be  
16 important to have really all of the alternatives  
17 available in order to get safely to the other side.  
18 It's a bit of a judgment call whether you would do the  
19 device first or the surgery first and hope for the  
20 best with the device later if the surgeon wasn't able  
21 to accomplish everything they had hoped to do.

22 DR. WILLIAMS: In the larger scheme of

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 things whether one should use device closures at all  
2 in centers that are not large volume experienced  
3 centers. I think this comes to the question of  
4 whether one should electively regionalize the sickest  
5 of the sick patients with known complex disease.

6 One easy question to end. There seems to be  
7 more fractures for the PFOs, 37 percent, than for the  
8 ASDs, 15 percent. Is that because the septum flops  
9 around more and it bends it more or is that incorrect?

10 DR. GAWREAU: We've actually noticed that  
11 larger devices are more likely to fracture. Larger  
12 devices are needed to close the PFOs and that's why  
13 you see the larger fracture rate and the higher  
14 fracture.

15 DR. WILLIAMS: Thanks.

16 DR. TRACY: Dr. White.

17 DR. WHITE: What are you planning to do  
18 about nickel- allergy?

19 DR. JENKINS: We actually have a lot to say  
20 about nickel and also nickel allusion. I think I'm  
21 going to refer that question to Carol Ryan, the  
22 engineer on the project because there are issues

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 beyond nickel allergy.

2 MS. RYAN: We've gotten that question many  
3 times and actually looked at that very early on in the  
4 design process. Significant studies were done to look  
5 at the medal ion to solution rates to be assured that  
6 they were very low. Significant literature searches  
7 have been done and discussions' with multiple  
a consultants regarding nickel allergy.

9 The one paper I tend to refer cardiologists  
10 to now when they ask that question because they have  
11 a patient with nickel allergies, a paper written by  
1 2 Katherine Merritt who actually works for the FDA. She  
13 did a nice summary on immune responses to metallic  
14 devices and their leechables.

15 Her conclusions were that -- she basically  
16 looked at all the literature that's out there as well  
17 as her own studies -- that there is no obvious  
18 relationship between a dermal response and a systemic  
19 one.

20 Her recommendation is that surgeons or  
21 clinicians should not deviate from their normal  
22 surgical practices based upon if a patient has a

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 nickel allergy or an allergy to any sort of metal ion.  
2 Devices should be designed so that the metal ion to  
3 solution rates are kept to the lowest possible amount  
4 and that was pretty much our conclusion.

5 I can think of at least 10 accounts to date  
6 where we've been approached because a patient was  
7 allergic to nickel and they've received a device and  
a we've had no adverse reports from that usage. The ion  
9 to solution rates for this device are actually  
10 extremely low. All the possible metal ions that could  
11 leech out of it were evaluated. In most cases they'  
12 were undetectable levels.

13 DR. WHITE: The second thing I have to say  
14 is a minor one. In Section 49.2 you describes the  
15 device as being 11 French and I think you've said  
16 today that it's 10.

17 MS. RYAN: It's 10.

18 DR. WHITE: You need to fix that.

19 Can you tell me, just educate me, in your  
20 tables about how well the patients did on one of the  
21 slides here, it says, "Clinical status CL by patient  
22 VSD pivotal cohort." Why did you assess the benefit

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

www.nealrgross.com



1 by median scale value? Why did you not use mean? Is  
2 there something about an ordinal scale evaluation that  
3 I don't understand?

4 DR. GAUVREAU: When you're working with an  
5 ordinal scale it's more appropriate to use medians  
6 rather than means. One reason is that the data are  
7 usually not normally distributed. The second reason  
8 is something I had mentioned earlier where the  
9 difference between a two and a three is not the same  
10 as the difference between a three and a four. It  
11 doesn't make sense to use means.

12 DR. WHITE: Fair enough. In terms of the  
13 doctor training in Section 5 you have several classes  
14 of physicians outlined. The third class is a  
15 fellowship trained doctor who you state may or may not  
16 have had a lot of experience. You were going to have  
17 your representative decide whether he needed to have  
18 Category II or Category IV training.

19 DR. JENKINS: I think that would depend on  
20 where the fellowship training was. For example, there  
21 are some people who spend an entire year in  
22 inte-rventional training fellowship.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 2344433

www.nealrgross.com

1 DR. WHITE: What I'm suggesting is that you  
2 delete the class and that you make your decision based  
3 upon whether the physician is qualified with implant  
4 or not. He's either a two or a four.

5 DR. JENKINS: Okay. I understand.

6 DR. WHITE: Take away No. 3. There's no  
7 point in that. You save the embarrassment. You save  
a your company walking up to a young doctor who thinks  
9 he knows what he's doing and you have to tell him he  
10 doesn't. It's never very pleasant.

11 The other thing is that under No. 4 you talk  
12 about proctoring doctors but 'you don't specify the  
13 number of cases that will be done. Have you given  
14 that any thought? How many cases will a proctor take  
15 an experienced physician and when is it enough?

16 MS. KULIS: Certainly, I'll ask Dr. Jenkins  
17 or one of the other clinicians to elaborate but as a  
18 company we thought that a minimum of five proctor  
19 cases would be what we would consider acceptable  
20 before we would certify the site to receive devices.

21 DR. WHITE: Given that this busy hospital  
22 did 57 in four years, how long is it going to take

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 somebody to get proctored?

2 DR. JENKINS: I think we would be very open  
3 to suggestions about how the training should be done  
4 for this project.

5 DR. WHITE: Okay. I think that it's a very  
6 complicated procedure. I don't do this procedure but  
7 it looks as if more than half your patients had  
a multiple devices placed and that more than two  
9 operators participated in 67 percent of your cases.  
10 It sounds like a little bit different than closing an  
11 ASD. I'm a little concerned about the infrequency of  
12 the procedure and then how are you going to get people  
13 trained to do this.

14 I don't want to be rude but I would  
15 challenge your primary endpoint. Everybody here seems  
16 real happy that you've done this but I'm not happy.  
17 I'm used to endpoints that say that we had a procedure  
18 success and no major complication.

19 If you subject your data to that analysis,  
20 how many of your patients were successfully closed and  
21 walked or crawled out of the cath lab without a major  
22 complication? It seems to me like so many patients

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

1 had big complications that not very many people got  
2 out of this unscathed.

3 DR. JENKINS: I guess the question would be  
4 whether you mean a manageable complication or  
5 something that would meet a definition of a serious  
6 hemodynamic impairment. .I think if you use --

7 DR. WHITE: Most of the time we don't get to  
a make excuses. I mean, you set an endpoint and you say  
9 procedure success or technical success is deployment  
10 of the device. Procedure success is successful  
11 technical deployment with no major complication. You  
12 get to pick what your major complications are. Under  
13 those criteria what would be your --

14 DR. JENKINS: In those criteria I would have  
15 personally chosen probably survival as my outcome so  
16 we might have disagreed on what was the major  
17 complication.

18 DR. WHITE: I guess what I'm saying is that  
19 your ordinal scale has its own merits or demerits but  
20 you're not balancing a successful procedure with a  
21 pretty bad complication may not be such a desirable  
22 outcome.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. JENKINS: We didn't create a composite  
2 endpoint for this study. We gave the safety data and  
3 the efficacy data in parallel without an overall  
4 measure that combined the two.

5 DR. WHITE: I don't want you to think I'm  
6 being unreasonable. I understand that you can take a  
7 band off the kid and, you know, the baby is better  
8 than he was without the band off.

9 It's just that everything else we think  
10 about has to be graded according to the risk benefit  
11 and so you don't get to claim a success if you have a  
12 major complication even if technically the procedure  
13 was effective.

14 What is a STARFlex? You had three patients  
15 crossover to STARFlex. Is that a competitive device  
16 or is that just another iteration?

17 DR. JENKINS: It's the third generation of  
18 this one that has been introduced.

19 DR. WHITE: Of this device?

20 DR. JENKINS: Yes. There's not as of yet  
21 sufficient STARFlex data to put before our panel.

22 DR. WHITE: Why did you cross patients to

**NEAL, R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 the newer device?

2 DR. JENKINS,: They weren't crossed over.  
3 The device was introduced within the time frame where  
4 the CardioSEAL was -- the CardioSEAL is still  
5 available in this study and it's the selection-.of the  
6 implanting cardiologist whether ,a CardioSEAL or a  
7 STARFlex is chosen.

8 They weren't crossed over to a STARFlex but  
9 we were just being strict that when we gave you  
10 information on all VSDs enrolled through 2/1 2000  
11 there were three that were not enrolled with  
12 CardioSEALs that were not included in this data  
13 summary. Maybe I'm not being clear. They didn't  
14 crossover into a STARFlex.

15 DR. WHITE: How did they get a STARFlex and  
16 get reported in this database?

17 DR. JENKINS: They are not reported in the  
18 database. That's the point. We gave you data through  
19 2/1 2000, all of the VSDs that were enrolled in the  
20 trial.

21 DR. WHITE: In this trial?

22 DR. JENKINS: In this trial. Everyone that

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 was enrolled through 2/1 2000 but we're only reporting  
2 on -- excuse me?

3 DR. WHITE: Where are the three STARFlex  
4 patients?

5 DR. JENKINS: The STARFlex was introduced  
6 into the study in the early part of 2000. There  
7 happened to be three patients who met that  
8 definitional criteria who had a VSD who were enrolled  
9 in the study who were included in the overall dataset.

10 But because this particular part of the data  
11 was intended to show the performance of CardioSEAL,  
12 the STARFlex patients were not included in the 57.  
13 However, just to be maybe ultra conservative in our  
14 reporting, we told you that there were three that fell  
15 within the time frame of our enrollment.

16 DR. WHITE: So have you now gone past the  
17 CardioSEAL device and' are using STARFlex for this  
18 disease?

19 DR. JENKINS: At the Children's since we  
20 -have the STARFlex device for the high risk trial on an  
21 ongoing basis, VSDs are being done with both of the  
22 devices but quite a few of the recent ones are being

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 done with the STARFlex.

2 DR. WHITE: Why did you choose not to  
3 include the catheterization complications when you  
4 reported the adverse events? You told me that out of  
5 the 222 total adverse events, there were 32 that were  
6 device related and 35 that were implantation related  
7 and 85 that were related to the cath. But when you  
8 went to look at the summary of the adverse events, you  
9 didn't include cath complications in that.

10 DR. JENKINS: They are all in the Panel Pack  
11 in exhaustive detail.

12 DR. WHITE: I mean in the --

13 DR. JENKINS: In the primary income.

14 DR. WHITE: You said you were interested in  
15 the --

16 DR. JENKINS: The reason is that we chose --  
17 the reason is that most patients would be having a  
18 catheterization anyway. That's the spirit of choosing  
19 the outcome as the specific part of the study whereby  
20 the device was placed or the implant procedure was  
21 done.

22 What we did instead is that our safety

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1 committee spent an inordinate amount of time figuring  
2 out if a specific event was due to the implant part of  
3 the procedures, or do just having a catheterization.  
4 They made that distinction.

5 With them having done that, we counted as  
6 the primary safety outcome just the device or the  
7 specific part of the procedure where the large sheaths  
8 and the wires and all that were in the heart rather  
9 than simple things that were just the result of a  
10 patient having a cath.

11 DR. WHITE: Well, the problem with that is  
12 that because you're not comparing this to anything  
13 else and the catheterization is integral to the device  
14 implantation and delivery, it's a little bit  
15 disingenuous. It makes the procedure seem safer than  
16 it might actually be.

17 If you want to know what's the risk of this  
18 baby or this child to undergo this procedure to take  
19 the cath complications out when, in fact, they were --  
20 maybe they weren't as serious but they outnumbered the  
21 number of other complications.

22 DR. JENKINS: There are a large number and

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

1 they are all listed in the Panel Pack in a lot of  
2 detail.

3 DR. WHITE: When you look at the primary  
4 safety outcome, it looks like that number may be less  
5 than it really was if you count the cath complications  
6 into it.

7 DR. JENKINS: One could have used a  
8 different definition That's true.

9 DR. WHITE: I'm really troubled by the  
10 fractures of the device. I'm really troubled by -- I  
11 mean, I know that you tell me that it hasn't called a  
12 problem but it bothers me that devices are breaking  
13 and I want to know what the company is doing about  
14 that. Are you making them so they won't break or you  
15 want me to keep putting them in to break?

16 DR. JENKINS: Again, I would like Carol Ryan  
17 to come up and talk about that.

18 MS. RYAN: We're actually -- the device, as  
19 I said, is made from MP35n and MP35n is the material  
20 that is used in pacemaker leads and pacemaker leads  
21 fracture and their fractures are unacceptable and  
22 usually have significant clinical sequelae.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1           The vendor who makes the MP35n for all of us  
2           who use MP35n wire in the medical device industry has  
3           a significant program that's ongoing to improve the  
4           quality of the raw material. We work very closely  
5           with them in evaluating each new generation of this  
6           material that comes out and implementing it into the  
7           product.

8           Kathy could probably comment to this better  
9           than I but an analysis was done of devices made from  
10          a variety of generations of this wire. We have shown  
11          that there is a statistically significant improvement  
12          in the fracture resistance of devices of the recent  
13          generation that has been incorporated.

14          We are continuing currently to evaluate  
15          future generations of the material that the vendor has  
16          provided us so we expect over time that the fracture  
17          rate will only get lower. Maybe Dr. Jenkins can  
18          comment on her analysis.

19          DR. JENKINS: We actually did do an analysis  
20          maybe three-quarters of the way through the data that  
21          showed you looking at determinants of fracture to  
22          figure out if there was specific manufacturing issues,

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 specific device design issues, or issues related to  
2 implantation that could be associated with fracture.  
3 It was a little bit of a fishing experiment. We  
4 looked at quite a few variables. We actually found  
5 three that were significantly related to fractures.

6 By far and away the most important one is  
7 device size as Kim pointed out earlier and as is shown  
8 in the fracture section of your Panel Pack whereby  
9 larger devices are more fracture prone than smaller  
10 ones. That confounder actually confounds a whole lot  
11 of other analyses that one might do looking at  
12 fractures.

13 The second one was a specific lot of devices  
14 that seemed to have an especially high fracture rate  
15 which was part of the impetus for Carol to go back and  
16 continue to look at the specific metal that's being  
17 used for manufacturing.

18 The third one was a very broad stroke  
19 variable whereby somewhere in the cath reports are  
20 follow-up letters. The procedure was described as a  
21 difficult device placement leading us to believe that  
22 pushing devices around bends in the sheath and things

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 like that may actually also be part of the determinate  
2 of fracture.

3 That was the most easily avoidable one. But  
4 we've done quite a bit to try to look into this. I  
5 think as clinicians having watched a large number of  
6 patients have fractures in the original Clamshell I  
7 cohort that we have also done extensive analyses on,  
8 and now quite a few patients experience this later.

9 We've had an increasing level of comfort  
10 around the issue that fractures really are incidental  
11 in the vast majority of cases probably because most of  
12 them are, in fact, occurring after the devices  
13 endothelialize and are completely covered. Just so  
14 you're aware, in the original Clamshell I registry  
15 series, there were seven events that were attributed  
16 to fractures in the hundreds of events that occurred  
17 in that cohort.

18 Those events were three masses that were  
19 associated with a fractured arm friction lesions,  
20 three devices that moved, and one arm that actually  
21 broke off and impeded in the free wall of the RV. I  
22 think we all wish that fractures would just go away

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 and not keep happening.

2 Even in the large number of patients in that  
3 original series who had fractures, the overall even  
4 rate was fairly small and fortunately we just haven't  
5 seen it all since 1996 despite screening extensively  
6 for them.

7 DR. WHITE: That's all.

8 DR. TRACY: Maybe this is a stupid question  
9 but why is the arm on the surface and not some place  
10 within so that it can't break lose and fly into the  
11 free wall or wherever it wants to go?

12 MS. RYAN: The predecessor, the Clamshell,  
13 where a piece of an arm migrated is somewhat of a  
14 mystery. It had to have been some sort of  
15 manufacturing defect. That device was made under a  
16 completely different processing controls than the  
17 current product.

18 The CardioSEAL device actually has each  
19 individual coil sewn to the fabric which did not  
20 happen with the Clamshell device. The nature of a  
21 fatigue fracture once one occurs in an arm, that arm  
22 really isn't under any significant stress at that

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 point. You shouldn't have a fracture at two points.  
2 With the coil sewn down there shouldn't be any  
3 migration.

4 DR. TRACY: Thanks.

5 Do any of the panel members have any  
6 additional questions they would like to ask the  
7 sponsor?

8 Dr. Williams.

9 DR. WILLIAMS: Just one very brief one.  
10 Under the contraindications, I think it would be  
11 reasonable to say the anatomy in which the CardioSEAL  
12 size required or position would interfere with  
13 intracardiac or intravascular structures because of  
14 the issue that you do select defects in which the  
15 position of the device would not interfere. I would  
16 put that specifically on the contraindications.

17 DR. TRACY: Any other members of the panel?

18 DR. LASKEY: Did I understand you correctly  
19 to say that you have not had a fracture since 1996?

20 DR. JENKINS: No. We haven't had any  
21 adverse consequences of a fracture in the entire high  
22 risk cohort.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. LASKEY: Just for my own clarification,  
2 two hours ago I asked the question who should this not  
3 be put in. I got a rather cursory answer which wasn't  
4 helpful. Now I come away hearing that there are  
5 defects where it shouldn't be approached.

6 Can you just give me a Reader's Digest  
7 summary of who this is appropriate for vis-a-vis which  
8 patients are not surgical candidates which, of course,  
9 you have in your IFU, but more specifically the  
10 anatomic subset which is not likely to do well with  
11 this procedure.

12 DR. JENKINS: That are not likely to do well  
13 with the cath procedure? Is that what you're asking?  
14 I think that the subgroup of patients that are not  
15 likely to do well with this procedure would include  
16 patients with VSDs in locations that are within 5 mm  
17 of semilunar or AV valves or valve apparatus. Or  
18 patients who are too small to have placement of 10  
19 French catheters in their vasculature.

20 DR. LASKEY: And the postero-septal defects  
21 that are perhaps a little too close to the base and to  
22 the insertion of tricuspid leaflets. I took something

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1 away from that discussion as perhaps being not as  
2 ideal a situation as other regions.

3 DR. JENKINS: I'm going to ask Dr. Lock to  
4 answer this question.

5 DR. LASKEY: Over the last couple of hours  
6 the answer to that question changed.

7 DR. WILLIAMS: My interpretation of his  
8 answer is it's harder than the other ones but it may  
9 be the only alternative. The question we have to  
10 determine is whether in hands other than Dr. Lock's it  
11 is likely to be successful.

12 DR. LASKEY: And that summarizes my concern.  
13 Dr. White, thank you for getting my adrenaline going  
14 again. You guys are experts beyond two standard  
15 deviations of the average interventional cardiologist.

16 If you expect this technology and capability  
17 to penetrate into the lower levels or the lower  
18 echelons of this profession, I don't have any desire  
19 to do this. I'm not even sure I could but if I wanted  
20 to.

21 Frankly, I'm intimidated and I've been doing  
22 intervention in sick people for 20 odd years but this

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 is a whole other order of magnitude here. I don't  
2 know if I speak for the profession or just for myself,  
3 but I get the feeling that there is a body of  
4 knowledge here and the level of expertise which  
5 desperately needs to see the light of day in order to  
6 make informed judgements about who should get this.

7 It has to be done in the context of expert  
8 surgery, expert anesthesia, a whole group of experts  
9 which is to be found only in 30 centers, did you say?

10 This all started out with my unease as the  
11 afternoon developed about, well, it's going to work  
12 better in some than in others. I think that's not  
13 clear from this material. I think users other than  
14 you need to know what to expect.

15 DR. ZAHKA: I would agree with you but  
16 disagree in the sense that the community of  
17 interventional pediatric cardiology is a very broad  
18 one. We heard about one center who has done three  
19 successful ones. I think that Dr. Lock is probably  
20 very articulate, because he is very articulate, at  
21 telling what are the tricks of the trade.

22 There's a great body of experience in

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

1 pediatric interventional cardiology that I think can  
2 be brought to bear on this so that the situation is  
3 perhaps not as bleak as it might seem from the adult  
4 world.

5 DR. LOCK: This is probably gratuitous and  
6 unnecessary, but there was a period of time 20 years  
7 ago where there was really only one place in the  
8 country that did hypoblast surgery. There was a  
9 period of time when really it was thought that only a  
10 few places could successfully perform that procedure.

11 It did take five or 10 or 15 years for that  
12 operation to become a national standard. Now, it  
13 isn't done in every center in the country but it is  
14 done in quite a few centers around the country. I  
15 expect exactly the same transition will happen with  
16 this kind of complicated intervention in children.

17 Therewillbmore complicated interventions  
18 in children like this that won't be done in two or  
19 three hundred places but will be done in 50 or 80 or  
20 30 or 20 very successfully as time goes on.

21 DR. TRACY: I think there is some difficulty  
22 because the only real concrete thing here is the death

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 rate which was about 7 percent. But you're talking  
2 about a procedure that has a 99 percent adverse event  
3 rate which anybody could go out and say, "I'm going to  
4 do a procedure now because there's almost a 100  
5 percent chance that something will go wrong."

6 I think that in the education of the  
7 physicians, all of these intangible things really have  
8 to be conveyed very clearly. Who best is this suited  
9 for? Who is this not suited for? What are the things  
10 that we have learned from our experience?

11 That kind of information has to be passed  
12 along because not even well-trained interventional  
13 cardiologists will have had that much experience doing  
14 transeptals. There's about a 1,000 pitfalls in this  
15 procedure where things can go wrong.

16 Each of those steps require some training.  
17 It's not everybody who should be taking on this type  
18 of procedure. I think that is the unease that many  
19 people feel about this procedure.

20 Dr. Williams.

21 DR. WILLIAMS: This is getting a little bit  
22 more into the domain of discussion than question so

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 I'll include Dr. Lock in this discussion point.

2 It seems to me those of us who have looked  
7 at surgical outcome relative to institutional and  
4 operator volume know that in general there is a  
5 difference between large and small but there are many,  
6 many exceptions that have to do with institutional  
7 organization accumulation of learning curve.

8 One option that we would have is to try to  
9 put some very arbitrary volume limits on this. But I  
10 wonder whether in the end more patients would be  
11 served if we put very, very heavy educational  
12 requirements' on the team and institutional record  
13 keeping. And if there were very, very careful post-  
14 market surveillance and that perhaps taking the most  
15 difficult type of VSD which would be the posterior  
16 muscular VSD and say in order to qualify to do that  
17 type of VSD, that institution would have to have both  
18 efficacy and safety record equivalent to Boston  
19 Children's Hospital. Now, that would be tough but it  
20 would be -- you know, you could earn your --

21 DR. LOCK: I intend to make it impossible.

22 DR. WILLIAMS: Of course it would. At least

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 equivalent to the average of the institutions of the  
2 group which would be a lower target. The indications  
3 to this, who is going to do it, is whether you belong  
4 to the tribe that believes in no stone unturned in a  
5 dying patient, or you belong to the tribe that says  
6 above all do no harm.

7 That's a matter of philosophy. That is also  
8 a matter of what your other alternatives are. It is  
9 an imponderable when we talk about different  
10 institutions because the resources of those  
11 institutions are different and every patient is like  
12 a snowflake. They are different.

13 I personally would feel more comfortable  
14 saying go ahead, but putting these stringent  
15 requirements on education of the team on post-market  
16 surveillance and letting that be as close as we can  
17 get to what is the right thing.

18 MR. DILLARD: Dr. Tracy, Jim Dillard. Just  
19 a point of reminder for the advisory panel is that we  
20 are sort of skirting that line and going over and  
21 coming back a little bit in terms of practice of  
22 medicine and just to remind you that we really don't

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 get involved with a lot of the practice of medicine.

2 I think Dr. Williams brought it back a  
3 little bit to say what some of those training  
4 requirements might be which is something we'll work  
5 obviously very closely with the company on.

6 A number of these in terms of who's going to  
7 do it and how many you have to do, I think, really  
8 gets in much more to the practice of medicine and  
9 something that I think their profession needs to  
10 regulate a lot more than the agency is going to. I  
11 just wanted to remind everyone.

12 DR. TRACY: I would agree with that except  
13 to the extent that this is a team approach and I think  
14 that part of the physician training -- what I would  
15 take from this as a concrete thing is part of the  
16 physician training has to include all the different  
17 pieces of the team that are going to be present or  
18 potentially present including the cardiac surgery  
19 team.

20 DR. HOPKINS: I just want to say thank you.  
21 I'm just about ready to raise that issue. We're  
22 talking about 57 patients here in four years. I mean,

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 we can get so stringent that no patient ever -- the  
2 patients are out there dying because they don't have  
3 access to this device because we've created this  
4 philosophically stringent.

5 When I first went to medical school at Duke,  
6 the only place in North Carolina that did aneurism  
7 surgery was Mr. Duke's hospital and now it's probably  
8 done in every hospital that has 50 beds or more.

9 I think we are getting way afield of  
10 labeling and indications and what is intended here  
11 which is moving a device that has been remarkably  
12 effective in a very tough set of patients from a  
13 humanitarian device to a premarket approval. All of  
14 the other stuff that sort of in the last 20 minutes  
15 has been very philosophical but I don't think has  
16 anything to do with this.

17 I agree with requesting of the company to do  
18 rigorous training but that's different than limiting  
19 access of the device to some subset of a subset of a  
20 subset.

21 DR. TRACY : I think there is one more  
22 question from the panel.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)



1 DR. WITTES: I feel like I'm in Never Never  
2 Land. I don't understand. I need to hear some  
3 numbers about what the mortality would have been.  
4 What we're hearing is this is remarkably effective,  
5 what the mortality would have been had the device not  
6 been here.

7 What would the shift have been in the  
8 clinical efficacy? I worry exactly as Dr. Laskey does  
9 about whether -- how much of this is regression to the  
10 mean, It may be none of it is but I need to hear you  
11 tell me that if I had 57 patients and I didn't give  
12 them this device, X number would die within six months  
13 and nobody would shift over in the improvement.  
14 Otherwise, I'm feeling like it's a matter of faith.

15 DR. JENKINS: I think we should have John  
16 answer that.

17 Fifty-seven patients, John. Half had failed  
18 VSD surgery elsewhere. The other half had passed a  
19 peer review whereby a surgeon, maybe yourself, maybe  
20 someone else, and a cardiologist had declared that the  
21 VSD would have been very difficult to approach in the  
22 operating room

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. MAYER: Well, I guess what was running  
2 through my mind there is to give you some context  
3 about what's the natural history. Forget surgeon,  
4 cardiologist, or anybody.

5 The natural history of patients with large  
6 ventricular septal defects, large defined as having a  
7 big left-to-right shunt is as follows. There's a  
8 large number of those patients who will die from  
9 congestive heart failure.

10 There's a huge volume load placed on the  
11 heart. There's three times as much blood going  
12 through the lungs every minute as go through the body.  
13 Those patients are highly susceptible to pulmonary  
14 infections.

15 A virus that you or I would throw off will  
16 kill those children. You know, they can't grow  
17 because they are wasting so much metabolic energy  
18 pumping all that extra blood around that they can't  
19 devote energy to getting bigger like babies are  
20 supposed to get.

21 And there are a significant number of those  
22 patients who have elevated pulmonary blood flow so a

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 lot of extra blood going through the lungs at very  
2 high pressure who will then progress to develop what  
3 is called pulmonary vascular obstructive disease. The  
4 natural history of an untreated large ventricular  
5 septal defect in children is particularly unfavorable.

6 That's why 40 years ago when cardiac surgery  
7 started, cardiologists were willing to send patients  
8 to surgery even who had surgically easily accessible  
9 VSDs because the mortality rate was 25 percent with an  
10 operation, but it was still better than what the  
11 natural history was.

12 So that's the sort of floor context. If we  
13 take the subset of patients who had a pulmonary artery  
14 band which is a palliative procedure that you can do  
15 that will limit the amount of pulmonary blood flow  
16 drops the pulmonary artery pressure down strained to  
17 the band, keeps them from getting pulmonary vascular  
18 obstructive disease, and we don't have an adjustable  
19 band.

20 What might work pretty well for a baby age  
21 six months, by the time that child is three or four  
22 years old, they're not going to have left-to-right

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 shunt. They're going to have right-to-left shunt.

2 They are going to be blue. They are going  
3 to be exercise limited. They are going to be at risk  
4 for strokes and all the things that kids with cyanotic  
5 heart disease get. That's another subset of what can  
6 happen.

7 Certainly the patients who went to surgery  
8 to have a VSD closed in whom it didn't work -- the  
9 surgeon couldn't get access to it because it was in a  
10 difficult location or whatever other reasons there  
11 might be, complicated anatomy -- those patients  
12 presumably went to surgery because there was an  
13 indication for doing an operation.

14 From my standpoint, and I guess I would  
15 hearken back to the practice of medicine question  
16 versus what is the device related issue, at least in  
17 our place this has been a pretty rigorous process  
18 **because** you have to get a surgeon and a cardiologist  
19 both to agree that this is something that is the best  
20 course of action, least risk path of treatment for  
21 this particular patient. It's really done on a **case-**  
22 **by-case** basis. That's always informed by a whole

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 variety of personal experience, literature experience,  
2 so forth.

3 I guess from my perspective, and having been  
4 a reviewer on a number of these cases as they have  
5 come along and, to be honest with you, having kicked  
6 some out saying, "I think I can close that hole," ones  
7 that came through, and some of which I actually did  
8 operate on and close the whole, I think all of those  
9 factors make it, I think, 'extremely difficult to  
10 construct a control group.

11 In the same way that there were difficulties  
12 with having what is clearly a multiple clinical  
13 presentation set of patients, and trying to figure out  
14 a scale how you deal with the banded patients who then  
15 got their device closed and then had their pulmonary  
16 band taken off, and construct a scale that is also  
17 consistent with the patient who had multiple  
18 ventricular septal defects and hadn't been banded and  
19 had one or more VSDs closed by device, I mean, it  
20 inherently is just a complicated set of patients.

21 I think that is the problem with -- I mean,  
22 I understand from a statistical standpoint why one

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 would like really to have a comparable group of  
2 patients.

3 DR. WITTES: Well, I'm not even asking that  
4 much. I'm asking for a number. I'm hearing  
5 essentially 11 percent, six-month mortality in this  
6 group is what there is. Is that right?

7 DR. JENKINS: Four patients died and one  
8 died because of the catheterization for a mortality of  
9 1.7 percent. One patient out of 58 patients died  
10 directly due to the procedure.

11 DR. WITTES: But, to me, it's still four out  
12 of 57. However --.

13 DR. LOCK: Can I interrupt for a second?  
14 The other three patients who died died from their  
15 underlying disease.

16 DR. WITTES: That's what I'm asking. What  
17 percentage of people -- if you had 57 --

18 DR. LOCK: Those were the patients who  
19 weren't made better necessarily. For all the patients  
20 who were made better, it improved their overall  
21 survival.

22 I don't know how to put this but there have

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 probably been 10 patients who were in one fashion or  
2 another didn't come for the device and I know of three  
3 who died waiting. This is a very difficult patient  
4 population to get mortality rates on. If we gave you  
5 a number, it would be arbitrary.

6 DR. WITTES: I don't care if it has a 20  
7 percent spread. I just want to know --

8 DR. JENKINS: The old-fashioned number that  
9 is widely taught to cardiologists was that there was  
10 a 20 percent of patients with this disease that didn't  
11 come off pump.

12 I actually tried to track down where that  
13 number came from because it's been widely quoted. I  
14 had trouble actually finding it so I tend not to give  
15 information I can't find.

16 It's wide quoted that the mortality rate of  
17 not coming off pump, if you take somebody with  
18 multiple VSDs to the OR and you don't close all of  
19 them, it's 20 percent.

20 DR. TRACY: Okay. I think I'm going to ask  
21 Mr. Morton and Mr. Dacey if they have any additional  
22 questions. I know there are many sort of unanswered

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 questions here but unless there is something very  
2 specific that can be answered by the sponsor, I think  
3 we need to move on to the FDA questions.

4 Can I ask the sponsor to stand back and  
5 we'll move on to the FDA questions if somebody can  
6 flash those back up.

7 The first question is dealing with the  
8 complexity of the VSD in patients entered in this  
9 registry has been defined variously as VSD not  
10 accessible to closure through an atrial or aortic  
11 approach associated with other cardiac pathology  
12 patients with single or multiple muscular septal  
13 defects or simply patients at high risk for surgery.

14 Question 1a. Based on the information  
15 provided, please discuss the description of "complex  
16 VSD" as the defining indication for use of the  
17 CardioSEAL for VSD closure.

18 I think in the indication in Section 2, I  
19 think it is, the indication is the CardioSEAL  
20 inclusion system is for use in patients with a complex  
21 ventricular septal defect of a significant size to  
22 warrant closure, but that based on location cannot be

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1 closed with standard trans-atrial or trans-arterial  
2 approaches, which is a little bit more simplistic than  
3 what Dr. Mayer detailed or than the patients that are  
4 actually included in this study.

5 I would suggest perhaps using something that  
6 is a little bit more reflective of Dr. Mayer's, I  
7 believe, sixth slide that listed the definition of  
8 high risk which included low probability of  
9 satisfactory surgical exposure, left ventriculotomy,  
10 excessive right ventriculotomy, high probability of  
11 residual VSD, failed previous VSD, multiple apical  
12 and/or anterior muscular VSDs, and posterior apical  
13 VSD covered by trabeculae.

14 I think maybe more specifically stating in  
15 the indications the actual patients that were included  
16 would be helpful.

17 MR. DILLARD: Can I ask -- excuse me. Jim  
18 Dillard. Can I ask a real quick question, which is is  
19 that all encompassing? I mean, are we even missing  
20 anything with that that may be important if we don't  
21 have the general statement. That would be my only  
22 question.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. SKORTON: I think there were a couple of  
2 other things that will be in the transcript from Dr.  
3 Lock's remarks that should be folded into there too  
4 about post-infarction VSDs and posterior versus  
5 anterior. I think the sense of what she brought up is  
6 right.

7 DR. WILLIAMS: But the indications, I think,  
8 are, as you say, are good. The contraindications may  
9 indicate the post-infarction VSD. I think defects  
10 that interfere with the valve would be in the  
11 contraindications. It happens that most of those  
12 defects are accessible so I think that is the correct  
13 -- you have the correct definition for indications.

14 DR. WHITE: I don't think we saw any data  
15 about contraindications. Did we? I mean, I think we  
16 just don't want to list it as an indication but I  
17 d o n t think we saw any data regarding the  
18 contraindication.

19 DR. TRACY: I think the contraindications  
20 are what are listed here, the obvious things on clots,  
21 etc. I believe Dr. Lock's comments have to be  
22 reflected somewhere in there. I don't know that I

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

1 would put them down as contraindications but perhaps  
2 data is less than optimal results or some type of  
3 qualitative statements could be made regarding that.

4 Question 1b. In the absence of a control  
5 group, please discuss how to evaluate the safety and  
6 effectiveness of the CardioSEAL device.

7 I think you've heard the discussion. There  
8 is no control group. It's what it is in a very high-  
9 risk patient population.

10 Question 2. Does the use of the Clinical  
11 Status Scale allow for a clinically meaningful  
12 assessment of effectiveness for the device?

13 Again, I think you've heard the discussion  
14 about that. It's difficult to get a handle on it but,  
15 again, it is the definition that was used. There are  
16 data here that are useful. Any other comments  
17 specifically on that?

18 DR. WHITE: I'm just troubled by the fact  
19 that there is no composite endpoint that should be --  
20 I mean, I'm not asking for a randomized trial there.  
21 I'm asking for a very conventional way that we assess  
22 outcomes and this didn't do that.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. WILLIAMS: I would suggest that this  
2 isn't a conventional group and that's why we can't  
3 because there really is no composite. They were asked  
4 to do it and they did the best that they could under  
5 the circumstances.

6 In truth, to mix the indications of left-to-  
7 right shunt in more complex right-to-left shunts is  
8 probably meaningless and I think they made as good an  
9 attempt as they can possible do.

10 DR. WHITE: I don't think that's true. I  
11 think given the data here I could tell you how many  
12 people had the procedure done, a technically  
13 successful procedure, and had a major complication.

14 I mean, it's just a matter of how you  
15 measure the data and whether you accept or whether you  
16 require the fact that success happened without or with  
17 a major complication and whether you're willing to let  
18 that happen.

19 DR. WILLIAMS: You could look at technical  
20 success with closing the hole but if the issue is the  
21 effect on the patient's course, then you cannot mix  
22 those two things together, I don't think.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. SKORTON: I think the answer to No. 3  
2 sort of resolves the issue of No. 2.

3 DR. TRACY: Question 3. Based on the data  
4 provided and your comments regarding questions 1 and  
5 2, please discuss whether these data provide  
6 reasonable assurance of safety and effectiveness.

7 I think that's obviously what we're  
8 struggling with. This is not a safe group of people  
9 to be working. However, it does appear to be a viable  
10 option for treatment in this very high-risk group of  
11 patients.

12 DR. WHITE: I think that is the reason for  
13 an HDE.

14 DR. TRACY: Anything else troubling? Ms.  
15 Moynahan seems troubled by that. I'm not sure why.

16 DR. JENKINS: It's kind of the pivotal  
17 question and I think a couple of the comments might  
18 help;

19 What do you think, Jim?

20 MR. DILLARD: Well, I mean, I think we heard  
21 Dr. White have perhaps a little bit different  
22 perception. There's not a right or wrong answer even

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 to the question I think that you're raising, which is  
2 how do you differentiate what is an HDE versus what is  
3 a PMA.

4 Let me try to boil it down into something  
5 pretty simple which is this product is on the market  
6 at 30 institutions because the company has  
7 demonstrated that there is reasonable assurance of  
8 safety and that there is probable benefit.

9 Now today what we're saying is the data that  
10 we're looking at today pushes over the line from  
11 reasonable assurance of safety and probable benefit to  
12 reasonable assurance of safety and reasonable  
13 assurance of effectiveness.

14 I think that is perhaps the pivotal question  
15 here today which is the data now presented here with  
16 57 patients enough to say there is reasonable  
17 assurance of effectiveness.

18 At the time we looked at the HDE a lot of  
19 that information wasn't complete. Safety seemed to be  
20 there. Is this really enough to judge effectiveness  
21 of the product for this patient population.

22 DR. HOPKINS: I would have to say for me the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealgross.com](http://www.nealgross.com)

1 answer to that is yes, that there is reasonable  
2 assurance and that one suppose. I can actually give  
3 you my answer to your question because I'm not bound  
4 by the data. As a surgeon who would have to make a  
5 decision whether to operate on these patients, I would  
6 typically quote these parents 25 to 50 percent  
7 mortality so if that gives you a figure compare.

8 DR. WITTES: Yes, that's the sort of figure  
9 I needed.

10 DR. WHITE: But the question then, Dr.  
11 Hopkins, is what has persuaded you that they need more  
12 than an HDE, you know, if this device isn't ready for  
13 prime time? I'm not arguing that 'this device  
14 shouldn't be used and I'm not arguing that you have a  
15 need for this in your patients. What I'm suggesting  
16 is I haven't been convinced that there is a need more  
17 than a HDE.

18 DR. HOPKINS: I think Jim Lock actually  
19 referred to it. The actual dynamics of 'what happens  
20 with these patients is that if you don't have such a  
21 device available, you either get pushed towards  
22 surgery or the patient sits waiting for resolution in

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 terms of referral to a center that does have this  
2 available.

3 DR. WHITE: This device is available in 32  
4 centers of which we received no data. We don't know  
5 how those people performed. One of my concerns is  
6 that this all-star group here who had significant  
7 problems is not going to be translatable to those  
8 other 32.

9 DR. HOPKINS: Yes, I share those concerns.  
10 I think in the questions to come is where I would  
11 recommend that we resolve that. That is, in the  
12 training issues and then perhaps the post-release  
13 surveillance issues rather than in the PMA.

14 DR. TRACY: So I'm going to leave the answer  
15 as being within this very small group there is some  
16 assurance of the effectiveness of this procedure as  
17 well as the safety.

18 Moving on to the training program. The  
19 summary of that is in Section 5 of the Panel Packet.

20 Question 4a. Please discuss any improvements  
21 that could be made to the training program.

22 I think it's just a very, very difficult

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)



1 thing to come up with a training program that will  
2 reflect getting trained as a superb and highly  
3 talented interventional cardiologist who has access to  
4 the world's best cardiac surgeons in the presence of  
5 a highly trained and expert group of cardiac  
6 anesthesiologists, but somehow you have to convey that  
7 all of those pieces are needed in this training  
8 program.

9 I think to reflect all of our concerns, the  
10 training has to somehow haul in all these people and  
11 get them to understand the seriousness of the clinical  
12 situation. I don't know exactly what to do with the  
13 fact that in this protocol a group decision was made  
14 between surgery and cardiology as to whether the  
15 person was a candidate for this device.

16 Is that something that we would recommend  
17 that that discussion be held on each individual  
18 patient, or is this the decision that the cardiologist  
19 is going to make and then the surgeon is going to have  
20 to live with? I don't know. I'm asking the surgeons  
21 whether they would like that.

22 DR. HOPKINS: I would have to say the nature

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 of the practice of pediatric and adult cardiac surgery  
2 is actually very different in terms of the dynamics  
3 between the cardiologist and the cardiac surgeon.

4 I would think that in every center where I  
5 have ever been and have ever visited, the decision on  
6 therapy, particularly invasive therapy for pediatric  
7 patients, is done in concert and as a group and rarely  
8 done in the same fashion that adult decisions are made  
9 where a single cardiologist makes a decision and  
10 refers the patient to a single cardiac surgeon.

11 I think the actual general dynamics of the clinical  
12 care model is so different that it takes care of that.

13 DR. WILLIAMS: I would add in terms of the  
14 training, I certainly agree with what Dr. Hopkins  
15 said. I think in terms of training requirements I  
16 would specify that there be a locus of responsibility,  
17 echo, anesthesia, surgery, and cath. They meet as a  
18 team, and that the learning curve be concentrated in  
19 those individuals because it's terribly important  
20 starting out to accumulate the learning curve under  
21 one umbrella.

22 DR. SKORTON: I have a question about that

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 from doing this a few times. It's one thing to  
2 suggest that a person who does that has no teeth  
3 whatsoever. Once the thing is marketed you can do  
4 anything you want. You can put it in the very first  
5 person you see.

6 A question for you. If we bought into Dr.  
7 William's ideas is it practical or doable to insist  
8 that before being given access to the device someone  
9 go through a particular training program? Because if  
10 it is or it isn't, that would have a big effect on  
11 whether this is a practical idea or not.

12 MR. DILLARD: Jim Dillard. I think that one  
13 of the responsibilities on the part of the agency is  
14 to certainly work directly with the sponsor to try to  
15 come up with a reasonable training program.

16 I think our first approach to that is much  
17 of what you have already discussed here which is what  
18 have the world's experts learned in terms of the  
19 initial clinical approaches as well as what the data  
20 says and how do we translate that then to the general  
21 teams that might be at the other institutions.

22 I think we are at maybe a little bit of an

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 advantage than we usually are at this stage because  
2 perhaps they have already done it 29 more times than  
3 they would have otherwise done because they have been  
4 through that training and there are other institutions  
5 based on the HDE.

6 They probably learned even a lot more than  
7 the companies who would be sitting before us here  
8 saying, "We've only trained a couple three centers  
9 that we've done the clinical study on."

10 I think actually the sponsor may have some  
11 additional comments on that, No. 1, but 'beyond that,  
12 No. 2, we would work very closely with them, we would.  
13 learn from what their experience is, and that would be  
14 part of our conditions of approval to come up with a  
15 training program that is satisfactory to the agency.

16 DR. WHITE: If we simply required that a  
17 physician be proctored for three cases, which is  
18 common in many devices and other things, you could  
19 pocket veto this PMA because there aren't enough cases  
20 out there for the physicians to be proctored for three  
21 each. I think that is one of the big issues here.  
22 Who is going to save three of these up for a proctor?

**NEAL R. GROSS.**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. TRACY: However; there are 30 centers  
2 that somehow have managed to get the device up and  
3 running so there is a way to do this.

4 I was just concerned because the  
5 verification form only deals with the interventional  
6 cardiologist. There must be something, as Mr. Dillard  
7 says, that the company and the mentors already know  
8 that have permitted this thing to expand out to a  
9 number of centers.

10 DR. WHITE: I'd be careful about what you  
11 think the 32 centers are doing. I think we haven't  
12 seen any data regarding that.

13 DR. WILLIAMS: But I'd also be careful --  
14 I'm not myself interested so much in pocket veto. I'm  
15 more interested in helping the company set out the  
16 conditions that will end up with the best result  
17 because I think this is something that should be  
18 propagated safely.

19 DR. HOPKINS: I think the sense is that the  
20 group wants some rigor in the training. Ultimately in  
21 the latter questions of the post-market evaluation we  
22 are going to deal with some of those issues.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. TRACY: Okay. 4b. More than one device  
2 was placed in 26 patients. Please discuss training  
3 issues regarding the placement of multiple devices in  
4 a single patient.

5 Obviously, the more you do the more complex  
6 it is. The more training you need, the more  
7 sponsoring you need.

8 DR. WILLIAMS: But you might not always know  
9 when you're going to have to do that so I don't know  
10 that you can necessarily in advance decide that.

11 DR. WHITE: Remember that two-thirds of  
12 these procedures had two guys working. You talked  
13 about your anesthesiologist and other people but this  
14 is somebody pulling on this wire and somebody pulling  
15 on that wire and they are a team. This isn't what one  
16 good guy can go do. This is a real tour de force, I  
17 think, to do these well.

18 DR. TRACY: Again, emphasis on the team  
19 approach.

20 Product labeling and that information is  
21 contained in Section 2.

22 5a. Please comment on the INDICATIONS FOR

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 USE section as to whether it identifies the  
2 appropriate patient populations for treatment with  
3 this device.

4 I think we already discussed that.

5 5b. Please comment on the CONTRAINDICATIONS  
6 section as to whether there are conditions under which  
7 the device should not be used because the risk of use  
8 clearly outweighs any possible benefit.

9 The only thing that I would add there is  
10 that the thrombus that's mentioned is in various  
11 vessels but if you have somebody with a clot in the  
12 left atrium, you probably shouldn't be doing this  
13 either.

14 I think that -- I had written in my notes  
15 posterior muscular defects are at higher risk. I  
16 don't know if this necessarily rises to the level of  
17 contraindication but probably comes somewhere down in  
18 the warning section to just state that.

19 DR. WILLIAMS: But position that would  
20 interfere with the function of a valve, any of the  
21 cardiac valves, would be in addition.

22 DR. TRACY: Right. One of the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1        contraindications that is already listed is anatomy in  
2        which the CardioSEAL size required would interfere  
3        with other intracardiac or intravascular structure  
4        such as valves or pulmonary veins.

5                DR. WILLIAMS:    That only says size. It  
6        doesn't say position so I would say size or position,

7                DR. TRACY:    Okay.    I think we had discussed  
8        pulmonary veins as not appropriate to this particular  
9        application so we would probably take that wording  
10       out.

11               Any other specifics on contraindications?

12               5c.        Please        comment        on        the  
13        WARNING/PRECAUTIONS    section    as    to    whether    it  
14        adequately describes how the device should be used to  
15        maximize benefits and minimize adverse events.

16               I think this would likely be where we would  
17        add those other anatomic caveats.

18               5d.        Please        comment        on        the    OPERATOR'S  
19        INSTRUCTIONS as to whether it adequately describes how  
20        the device should be used to maximize benefits and  
21        minimize adverse events.

22               I read through this and thought that it was

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)



1 quite good but a picture is worth a thousand words and  
2 I think this is where the education would come in.

3 5e. Please comment on the remainder of the  
4 device labeling as to whether it adequately describes  
5 how the device should be used to maximize benefits and  
6 minimize adverse events.

7 Any additional comments?

8 Post-market evaluation. Question 6. Based  
9 on the clinical data provided in the Panel Package, do  
10 you believe that additional follow-up data or post-  
11 market studies are necessary to evaluate the chronic  
12 effects of the implantation of the CardioSEAL device?  
13 If so, how long should patients be followed and what  
14 endpoints and adverse events should be measured?

15 This is extraordinarily hard to come up with  
16 something like that in a population that is so limited  
17 to start out with. The numbers are so small to start  
18 out with. I think to recommend in a group of patients  
19 that are going to die of their underlying cardiac  
20 condition or other conditions anyway, it's extremely  
21 difficult to come up with a concrete recommendation on  
22 this.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, DC. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 I personally think that -- I hate to use the  
2 word registry but I personally think that something of  
3 that ilk is probably the right way to do this but I  
4 don't know. Do any of the other panel members have  
5 better comments than mine?

6 Dr. Wittes.

7 DR. WITTES: Well, can we take up Dr.  
8 White's suggestion that there are 29 centers out there  
9 with presumably data. Can those data be looked at?  
10 Is that legal? I mean, that actually would be part of  
11 the training. If those centers are having trouble,  
12 there may be information in the data that is already  
13 there.

14 DR. HOPKINS: There's really two questions  
15 that are being asked here, and that is the outcome of  
16 the individual patient in which that is probably known  
17 within six months of the implantation of the device or  
18 certainly within 12 months.

19 The other is the issue of the center  
20 efficacy as opposed to the patient based efficacy.  
21 That is sort of more of a registry, I think, type of  
22 approach. Maybe the follow-up should be suggested to

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 be 12 months for the individual patient and a number  
2 picked for a center.

3 The center has to maintain appropriate  
4 records and report them to the company and ultimately  
5 thus to the FDA for 15 or 20. Just pick a number.  
6 You are really measuring two completely different  
7 things here.

8 MR. DILLARD: I might make just a real quick  
9 comment and then the sponsor may just want to address  
10 it, too. I'm not sure, Dr. Wittes, whether or not  
11 those other institutions really have "data" per se.

12 They may have information and they might  
13 come up and even say they could go so far as to say  
14 whether or not they actually have some mortality  
15 information on perhaps what I would expect to be a  
16 very small number of patients even at some of those  
17 other centers.

18 I don't know how much we will actually glean  
19 from the knowledge of what we may know up to this  
20 point in time, but I think what might be important is  
21 if you are sitting here today, and I heard some issues  
22 that came up about what might be nice to know even in

1 the future if we come back three, four, five years  
2 from now, what is going to be important to be able to  
3 say about the CardioSEAL device for VSDs, especially  
4 complicated VSDs, that U.S. clinicians might want to  
5 know about, about how the product is doing and how  
6 would we assess it in a little bit longer term. Would  
7 that then be important to the post-market period to  
8 look at.

9 DR. TRACY: Dr. Skorton.

10 DR. SKORTON: I think it would be and I  
11 wonder in the interest of efficiency when I make the  
12 motion if I could present some specific ideas how to  
13 do that in the motion.

14 DR. TRACY: Do you want to wait until we get  
15 to the --

16 DR. SKORTON: Instead of discussing it twice  
17 because I have a motion.

18 DR. TRACY : That's fine. Okay. I think  
19 that is all of the FDA questions unless the FDA has  
20 any additional questions at this time or comments.

21 MR. DILLARD: No, thank you.

22 DR. TRACY: Does the sponsor have any

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 additional comments they would like to make at this  
2 time?

3 Mr. Morton, Mr. Dacey, any additional  
4 questions or comments?

5 Okay. Dr. Skorton, would you like to  
6 make --

7 MS. MOYNAHAN: You need to do open public  
8 hearing.

9 DR. TRACY: Oh, I apologize. Is there any  
10 member of the public here present who would like to  
11 make any comments at this point at an open public  
12 hearing?

13 If not, I'll close the open public hearing.  
14 Sorry I forgot that.

15 MS. MOYNAHAN: In case any of you forgot  
16 since this morning, I'll read them again.

17 The Medical Device Amendments to the Federal  
18 Food, Drug, and Cosmetic Act as amended by the Safe  
19 Medical Devices Act of 1990 allows the FDA to obtain  
20 a recommendation from an expert advisory panel on  
21 designated medical device premarket approval  
22 applications that are filed with the agency.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1           The PMA must stand on its own merits and  
2           your recommendation must be supported by the safety  
3           and effectiveness data in the application or by  
4           applicable publicly available information.

5           Safety is defined.in the Act as reasonable  
6           assurance based on valid scientific evidence that the  
7           probable benefits to health under conditions on  
8           intended use outweigh any probable risks.

9           Effectiveness is defined as reasonable  
10          assurance that in a significant portion of the  
11          population the use of the device for its intended use  
12          as conditions of use when labeled will provide  
13          clinically significant results.

14          Your recommendation options for the vote are  
15          as follows:

16               (1) Approval if there are no conditions  
17               attached.

18               (2) Approvable with conditions. The panel  
19               may recommend that the PMA be found approvable subject  
20               to specified conditions such as physician or patient  
21               education, labeling changes, or further analysis of  
22               existing data. Prior to voting all of the conditions

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 should be discussed by the panel.

2 (3) Not approvable. The panel may recommend  
3 that the PMA is not approvable if the data do not  
4 provide a reasonable assurance that the device is safe  
5 or if a reasonable assurance has not been given that  
6 the device is effective under the conditions of use  
7 prescribed, recommended, or suggested in the proposed  
8 labeling.

9 Following the voting the chair will ask each  
10 panel member to present a brief statement outlining  
11 the reasons for their vote.

12 DR. TRACY: Right. At this point, Dr.  
13 Skorton, I'll ask if you have a motion to make  
14 regarding this application.

15 DR. SKORTON: Yes. I move that the device  
16 be approvable with conditions and then, at the  
17 appropriate time, I have four conditions to suggest.

18 DR. TRACY: Go ahead.

19 DR. SKORTON: We have to have a second first  
20 to the motion.

21 DR. WILLIAMS: Second.

22 DR. SKORTON: Okay. My first condition is

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 that I believe there should be mandatory post-market  
2 studies for five years, that the studies should be  
3 annually, that a patient should annually get  
4 fluoroscopy and echocardiography, and that the six  
5 endpoints that should be looked for are the status of  
6 the device arms where fractures have occurred,  
7 thrombosis, global and regional ventricular function,  
8 endocarditis, evidence of ventricular arrhythmias or  
9 conduction disturbances, and evidence of residual  
10 shunt.

11 DR. HOPKINS: Could I address the issue of  
12 fluoroscopy? I don't think the arm fractures as we  
13 know it are really that important late because while  
14 it sounds like a bad engineering thing to have happen,  
15 actually late the device is locked in by the fibrous  
16 ingrowth.

17 From a practical standpoint an echo can be  
18 done in multiple outpatient facilities where fluro  
19 requires bringing them in to the hospital. Adding  
20 fluro adds a real increment of difficulty in the  
21 follow-up of these patients. I'm not so sure it's as  
22 important as the other criteria that you mentioned.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)



1 DR. SKORTON: I don't feel strongly about it  
2 but I'm responding to what I heard the investigator  
3 say was the way they discovered the fractures.. Since  
4 there will be new ones put in and since I thought I  
5 heard the engineering aspect of the sponsor say there  
6 was a little bit of a moving target in terms of the  
7 materials they were made out of and the way they were  
8 constructed, I'm uncomfortable not following up in  
9 some fashion.

10 If there is something that can be done  
11 besides fluoroscopy to look for arm fractures, that's  
12 great with me but I don't think echo would be the  
13 right way to do it.

14 DR. WILLIAMS: Would it be okay just to not  
15 specify the technique but to say what is best in that  
16 institution because even fluro if it's not done by the  
17 same person might not be as adequate.

18 DR. HOPKINS: I think mandating echo  
19 annually for five years is not inappropriate.

20 DR. SKORTON: Something, however, to look at  
21 the presence and outcome of device fractures because  
22 there were 16 percent fractures. Even though I agree

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 with what you said from the data we've seen here, this  
2 is now going to be open to a much broader denominator  
3 and I'm just uncomfortable. Maybe the device fracture  
4 rate is a lot lower but I don't know that.

5 DR. HOPKINS: But even if it is, I think the  
6 point is the arms could be absorbable and the ultimate  
7 outcome once it's locked in doesn't really matter.

8 DR. WHITE: I don't think you know that it's  
9 locked in. I think you -- I mean, I worry about that.  
10 I think we wouldn't be considering any device that had  
11 a one in five chance of breaking or a one in seven  
12 chance of breaking for most other applications.

13 I think it's a little cavalier anyway. This  
14 is an opportunity if we're going to do this to at  
15 least track it and at the end of five years be able to  
16 say whether any came out or not.

17 DR. HOPKINS: I was just pointing out that  
18 fluoroscopy is much more of an impediment to the  
19 mandated follow-up that you are suggesting.

20 DR. SKORTON: Maybe it's a certain kind of  
21 x-ray. I don't know, but I would ask that the  
22 condition be discussed with the sponsor and the

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 investigators who have collected pivotal data.

2 DR. WHITE: An impediment to the patient to  
3 come back, you mean?

4 DR. HOPKINS: Yes.

5 DR. WHITE: I mean, these are kids that are  
6 looking at getting transplanted. I mean, this is  
7 serious stuff. I mean, I don't think that's a big  
8 deal.

9 DR. HOPKINS: If they're out five years,  
10 they are doing pretty well.

11 DR. TRACY: Mr. Morton.

12 MR. MORTON: Regarding the diagnostic that's  
13 used and the effect it might have on the patient,  
14 might we not ask what is the result of the fracture  
15 and maybe look for those sorts of things rather than  
16 look for the fracture itself? We examine for  
17 fractures and we leave that up to the sponsor to get  
18 back.

19 DR. TRACY: I think, though, that the point  
20 regarding that is that we don't know what the  
21 consequences of the fractures are., We don't know if  
22 that later on that this will lead to some kind of an

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1 edge that can create even a rupture in the  
2 endothelialized surface that could lead to thrombus  
3 formation. We don't know.

4 There are other devices that have had  
5 fracture type of instances with them and they are  
6 followed by cardiac fluoroscopy. It is cumbersome but  
7 we do this. I don't think it's unreasonable in a  
8 device that has a 20 percent problem rate to request  
9 that fluro be done.

10 I personally would support that. I'm not  
11 committed to saying that they have to do fluro but I  
12 do think that is something that we don't know where  
13 that's going to go.

14 DR. SKORTON: Would you be more comfortable  
15 with, say, fluoroscopy or an equivalent technique?

16 DR. TRACY: Okay. So then your condition is  
17 that --

18 DR. WHITE: I'm running through those  
19 equivalent techniques here. There's fluoroscopy and  
20 fluoroscopy and fluoroscopy.

21 DR. TRACY: Well, you could get a flat PA  
22 and lateral. If you saw a big thing sticking off of

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 it, you would know there was a big thing sticking off  
2 of it.

3 DR. SKORTON: I heard the investigator say  
4 that they discovered some of the fractures with chest  
5 x-ray and some with fluoroscopy. I would be  
6 comfortable understanding that this is only advice for  
7 the agency and for the agency to work with the  
8 sponsor.

9 DR. WHITE: I've had the experience of  
10 looking at the fractures for the valves. The York-  
11 Shileys and the chest x-ray is not of the same -- I  
12 mean, you miss the little things with the chest x-ray  
13 so it's an underestimation, whereas with the fluro,  
14 and even sometimes sine is necessary depending on the  
15 thickness of the wires in order to be able to see that  
16 break. I think that it's not the same.

17 DR. TRACY : All right. Then for this  
18 particular condition, shall we vote on this particular  
19 condition for a five-year follow-up with the details  
20 as stated by Dr. Skorton.

21 MR. DILLARD: Jim Dillard. Just one quick  
22 question. I thought I heard the answer, but I'm not

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealrgross.com](http://www.nealrgross.com)

1       sure. It sounds like you're advocating, in addition  
2       to potentially following the cohort of patients that  
3       we currently have, you're talking about new patients  
4       that otherwise would receive the device? Is that  
5       correct?

6               DR. WHITE: Yes, that's correct.

7               DR. TRACY: So for new patients also. All  
8       in favor?

9               MS. MOYNAHAN: Ten in favor.

10              DR. TRACY : So that's unanimous. Any  
11       additional conditions?

12              DR. WHITE: Yes. I don't know exactly how  
13       to state it without taking three hours to do it but  
14       all this stuff that Dr. Williams said about augmented  
15       training procedures, something that could be boiled  
16       down by the agency and the sponsor I think needs to be  
17       added as a condition.

18              DR. TRACY: Okay. I won't even attempt to  
19       summarize the three hour discussion but some type of  
20       augmented training as a condition. All those in  
21       favor?

22              MS. MOYNAHAN: That's unanimous at 10.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 2344433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. TRACY: Any additional conditions?

2 DR. SKORTON: I guess just one more, and  
3 that is as one condition all the labeling  
4 clarifications that we mentioned under indications,  
5 warnings, and so on, all those together to be made as  
6 a condition.

7 DR. TRACY: Okay. So the third condition is  
8 verification of the changes in the labeling that we've  
9 suggested. All in favor?

10 MS. MOYNAHAN: Okay. That's 10.

11 DR. TRACY: All right then. The motion has  
12 been made that this is approvable with conditions.  
13 The conditions have been stated and voted on. At this  
14 point let's vote on the major motion approvable with  
15 conditions. All in favor?

16 MS. MOYNAHAN: Is your hand up, Dr. White?

17 DR. WHITE: No, it's not.

18 MS. MOYNAHAN: Nine.

19 DR. TRACY: Opposed? Can I then ask each of  
20 the panel members to individually state what your vote  
-21 was and the basis for your vote.

22 We'll start with you, Dr. White.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. WHITE: Well, I think I was the only one  
2 who thought this was not approvable. It's not because'  
3 I don't think the device is good or doesn't have a  
4 good use and it isn't valuable, but I was not  
5 convinced that it needs to be more than an HDE.

6 The administrative inconvenience of HDE to  
7 me doesn't justify the release of this device. I  
8 think we have a lot of chance to do a lot of harm here  
9 without doing a lot of good. I think the efficacy  
10 endpoint really was not -- didn't satisfy me.

11 I think the safety is questionable. I would  
12 have a lot of concern being on record for a device  
13 that has this fracture rate and approving that.

14 DR. TRACY: Dr. Williams.

15 DR. WILLIAMS: Well, from my clinical  
16 experience, I believe this is a group that has few  
17 other options. I believe that they have demonstrated  
18 reasonable efficacy and safety relative to what I  
19 understand the natural history of this disease to be.  
20 I believe that our conditions have set forth  
21 protections for the significant multiple operator  
22 dependence for this particular type of device

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234433

[www.nealgross.com](http://www.nealgross.com)



1 placement.

2 DR. TRACY: Dr. Skorton.

3 DR. SKORTON: I voted for approval for two  
4 reasons. One is that I've had the experience of not  
5 knowing what to do with a handful of people like this.  
6 It's just been a handful and I've become more  
7 convinced today that the surgical options are quite  
8 limited.

9 Secondly, I believe, although I do agree  
10 absolutely with safety concerns, which is why I  
11 brought up one of the conditions, I think this is not  
12 going to be one of those procedures that people are  
13 going to be running to do.

14 I think it will be somewhat self-correcting  
15 because of the very difficult nature of it. I have  
16 confidence that the agency before issuing an approval,  
17 if it chooses to, will develop some sort of training  
18 and surveillance system that will make me more  
19 comfortable.

20 DR. TRACY: Dr. Zahka.

21 DR. ZAHKA: I voted for approval because I  
22 think this is a difficult group of patients who need

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 this kind of approach. After I convinced myself that  
2 there would not be slippage of an approach to patients  
3 who, in fact, would be better done surgically. That  
4 was a major concern for me. I did come away convinced  
5 that this device would, in fact, find it's way only  
6 into patients for whom surgery was not a good option.

7 DR. TRACY: Dr. Hopkins.

8 DR. HOPKINS: I voted for it for the reasons  
9 that the two folks preceding me mentioned. I actually  
10 think it will increase the efficacy or the outcomes.  
11 Also for the surgical patients because of the kinds of  
12 conversations that the clinicians will have by having  
13 this device availability will foster the team  
14 approach.

15 DR. TRACY: Dr. Aziz.

16 DR. AZIZ: Well, I voted for it because I  
17 think this may be an option for a very difficult group  
18 of patients who really don't have much else even  
19 though I think I echo Dr. White's concerns that it  
20 does have a lot of questionable issues.

21 DR. TRACY: Dr. Laskey.

22 DR. LASKEY: Well, I voted for approval as

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 well with the qualifications noted but I'm terribly  
2 uneasy because this is the first time I've certainly  
3 reviewed anything which was not rigorously controlled.  
4

5 I think that many of us were responding  
6 emotionally and overreaching and, yes, this is a  
7 desperate population and, yes, it is nice to have  
8 another option and, yes, this probably will be used  
9 correctly by a small handful.

10 I think that ultimately came down to saying  
11 yea rather than nay. I just don't see 'hundreds of  
12 people using this device. I see it centrally  
13 controlled in expert hands. I hope it is as  
14 efficacious as we all hope.

15 DR. TRACY: Dr. McDaniel.

16 DR. MCDANIEL: I voted to approve with  
17 conditions as stated for the same reasons as my  
18 colleagues. I think that it's a limited number of  
19 patients. It will offer something to some children  
20 that may be expiring in institutions without the  
21 ability to do this. It's critical that the FDA follow  
22 some of our suggestions in terms of the training, but

**NEAL R. GROSS**

**COURT** REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, DC. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 I agree that it's not going to be done in a tremendous  
2 number of patients.

3 DR. TRACY: Dr. Wittes.

4 DR. WITTES: I voted yes for much the same  
5 reason. I became convinced that this is a desperate  
6 group that needs something. I wish there had been  
7 some more control data of one kind or another.

8 DR. TRACY: Dr. Crittenden.

9 DR. CRITTENDEN: I voted for approval with  
10 conditions. Again, I share a lot of the concerns  
11 voiced by previous panel members but this is a  
12 desperate group of patients who have few options so I  
13 think we've done the right thing.

14 DR. TRACY: Mr. Morton, any comments?

15 Mr. Dillard?

16 MR. DILLARD: Yes. I would just like to  
17 thank not only the two sponsors today but certainly  
18 this group of individuals who came in mostly for this  
19 day. There will be a few that I think will be back  
20 tomorrow, but I appreciate you all coming in today and  
21 taking a look at these occluder devices with us.  
22 Appreciate it.

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 20005-3701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

1 DR. TRACY : Thank you, everybody. I'll  
2 adjourned this meeting.

3 (Whereupon, at 6:00 p.m. the meeting was  
4 adjourned.)  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22

**NEAL R. GROSS**

COURT REPORTERS AND TRANSCRIBERS  
1323 RHODE ISLAND AVE., N.W.  
WASHINGTON, D.C. 200053701

(202) 234-4433

[www.nealrgross.com](http://www.nealrgross.com)

CERTIFICATE

This is to certify that the foregoing transcript in the  
matter of                      Circulatory System Devices Panel Meeting

Before:                      DHHS/FDA/CDRH

Date:                      September 10, 2001

Place:                      Gaithersburg, MD

represents the full and complete proceedings of the  
aforementioned matter, as reported and reduced to  
typewriting.

