

VSM

1

VSM

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
OPHTHALMIC DEVICES PANEL

OPEN SESSION

PMA P010059

*This transcript has not
been edited and FDA
makes no representation
regarding its accuracy.*

Thursday, January 17, 2002

9:40 a.m.

Hilton Washington DC North/Gaithersburg
Salons A, B and C
620 Perry Parkway
Gaithersburg, Maryland

MILLER REPORTING COMPANY, INC.
735 8th Street, S.E.
Washington, D.C. 20003-2802
(202) 546-6666

P A R T I C I P A N T S

PANEL PARTICIPANTS:

Jayne S. Weiss, M.D., Acting Chair
Arthur Bradley, Ph.D., Voting Member
Michael R. Grimmett, M.D., Voting Member
Alice Y. Matoba, M.D., Voting Member

Timothy T. McMahon, O.D., Consultant, deputized
to vote

Allen C. Ho, M.D., Consultant, deputized to
vote

Anne L. Coleman, M.D., Ph.D., Consultant,
deputized to vote

Joel Sugar, M.D., Consultant, deputized to vote

Richard Casey, M.D., Consultant, deputized to
vote

Janine A. Smith, M.D., Consultant, deputized to
vote

Woodford S. Van Meter, M.D., Consultant,
deputized to vote

Glenda V. Such, M.Ed, Consumer Representative

FDA:

Sara M. Thornton, Panel Executive Secretary
A. Ralph Rosenthal, M.D. Director, Division of
of Ophthalmic Devices

Everette T. Beers, Ph.D.

James F. Saviola, O.D.

Donna R. Lochner

Bernard P. Lepri, O.D., M.S. M.Ed.

Joel P. Glover, M.S.

Eleanor McGhee

C O N T E N T S

| | |
|--|-----|
| Call to Order | 5 |
| Jayne S. Weiss, M.D., Acting Chair | |
| Introductory Remarks | 5 |
| Sara M. Thornton, Executive Secretary | |
| Conflict of Interest Statement | 9 |
| Appointment to Temporary Voting Status | 10 |
| OPEN PUBLIC HEARING | 11 |
| OPEN COMMITTEE SESSION - Jayne S. Weiss, M.D. | 11 |
| BRANCH UPDATES - Donna R. Lochner, Chief | |
| Intraocular and Corneal Implants Branch | 11 |
| PMA P010059 | 12 |
| Sponsor Presentation | 13 |
| Morcher Capsular Tension Ring [CTR] | |
| Dr. Roger F. Steinert | |
| Ophthalmic Consultants of Boston | 13 |
| Hillard W. Welch | |
| Medical Device Consulting | |
| U.S. Representative for Morcher | |
| Panel Questions for the Sponsor | 30 |
| FDA Presentation | 76 |
| Donna Lochner - Introductions | 76 |
| Joel P. Glover, Team Leader | 79 |
| Bernard L. Lepri, O.D., M.S., M.Ed. - Clinical | |
| Reviewer | 84 |
| Panel Questions for FDA | 99 |
| Additional Comments from the Sponsor | 104 |
| COMMITTEE DELIBERATIONS | 105 |
| Primary Panel Reviewers: | |
| Joel Sugar, M.D. | 105 |
| Woodford S. Van Meter, M.D. | 114 |
| PANEL DISCUSSION OF P010059 | 125 |
| OPEN PUBLIC HEARING SESSION | 172 |
| FDA - CLOSING COMMENTS | 172 |

vsm

| | |
|---------------------------------------|-----|
| | 4 |
| SPONSOR - CLOSING COMMENTS | 173 |
| Voting Options Read | 179 |
| PANEL RECOMMENDATIONS TAKEN BY VOTE | 181 |
| POLLING OF PANEL VOTES | 230 |
| COMMENTS FROM CONSUMER REPRESENTATIVE | 235 |
| FINAL PANEL COMMENTS | 235 |
| OPEN MEETING ADJOURNED | 237 |

1 P R O C E E D I N G S

2 **Call to Order**

3 DR. WEISS: I would like to call the
4 Ophthalmic Devices Panel to order, and we will have
5 introductory remarks from Sallie Thornton.

6 **Introductory Remarks**

7 MS. THORNTON: Good morning and welcome to
8 the 103rd Meeting of the Ophthalmic Devices Panel.
9 Before we proceed with today's agenda, I have a few
10 short announcements to make. I would like to
11 remind everyone to sign in on the attendance sheets
12 in the registration area just outside the meeting
13 room here.

14 I just checked, and there are very few
15 signatures, and lots of people in here. So I think
16 there are some folks that need to see Annmarie out
17 there at the registration area.

18 All handouts for today's meeting are
19 available at the registration table. Messages for
20 the panel members and FDA participants, information
21 or special needs, should be directed through Ms.
22 Annmarie Williams or Mr. Hashim Khalif, who are
23 available in the registration area.

24 The phone number for calls to the meeting
25 area is 301/977-8900. In consideration of the

1 panel, the sponsor, and the agency, we ask that
2 those of you with cell phones and pagers either
3 turn them off or put them on vibration mode while
4 in this room. We're serious about this.

5 We ask that all meeting participants speak
6 into the microphone and give your name clearly so
7 that the transcriber will have the accurate
8 recording of your comments. All available
9 information for the meeting tentatively scheduled
10 for March 14-15 will be on the FDA Advisory
11 Committee website in approximately one week.

12 Now, at this time, I would like to extend
13 a special welcome and introduce to the public the
14 panel and the FDA staff two panel consultants who
15 are with us for the first time today and our new
16 panel consumer representative.

17 Dr. Richard Casey comes to us from Los
18 Angeles--there he is--where he is an Associate
19 Professor of Ophthalmology at the Jules Stein Eye
20 Institute and the Interim Chairman of the
21 Department of Ophthalmology at the Charles Drew
22 University of Medicine and Science.

23 His clinical practice involves the
24 management of the corneal and anterior segment
25 disease, cataract and refractive surgery.

1 Dr. Janine Smith is the Deputy Clinical
2 Director of the National Eye Institute of the
3 National Institutes of Health in Bethesda,
4 Maryland. Her basic science research has been
5 immune-based disease of the ocular surface with
6 additional responsibilities for the NEI intramural
7 clinical research program.

8 And Ms. Glenda Such, the consumer
9 representative to the panel, is the Director of
10 Computer Training Programs in the Department of
11 Career Services at Lighthouse International in New
12 York City. She is a recognized expert in the field
13 of adaptive technology for visual impairments and
14 the functional implications of visual disabilities,
15 particularly low vision.

16 We very much appreciate your commitment to
17 serve and welcome you to the panel table today.

18 To continue, will the remaining panel
19 members please introduce themselves beginning with
20 Dr. Van Meter?

21 DR. VAN METER: Woodford Van Meter,
22 University of Kentucky in Lexington, Kentucky,
23 practice in corneal and external disease.

24 DR. HO: Allen Ho, Philadelphia, Thomas
25 Jefferson University, Wills Eye Hospital.

1 DR. COLEMAN: Anne Coleman, Associate
2 Professor, glaucoma specialist at UCLA, Los
3 Angeles.

4 DR. GRIMMETT: Michael Grimmett,
5 University of Miami, Bascom Palmer Eye Institute.

6 DR. WEISS: Jayne Weiss, Professor of
7 Ophthalmology and Pathology at Kresge Eye
8 Institute, Wayne State University, Detroit.

9 DR. BRADLEY: Arthur Bradley, Professor of
10 Visual Sciences, Indiana University.

11 DR. MATOBA: Alice Matoba, Associate
12 Professor at Baylor College of Medicine.

13 DR. McMAHON: Tim McMahon, Professor and
14 Director of the Contact Lens Service at the
15 University of Illinois in Chicago.

16 DR. SUGAR: Joel Sugar, University of
17 Illinois in Chicago.

18 DR. ROSENTHAL: Ralph Rosenthal, Director
19 of the Division of Ophthalmic and Ear, Nose and
20 Throat Diseases, FDA.

21 MS. THORNTON: Thank you. I would like to
22 note for the record that at the sponsor's request,
23 the panel industry representative, Mr. Ronald
24 McCarley, will not be at the table today.
25 Therefore, the change will necessitate a slight

1 correction in today's agenda. The comments of the
2 industry rep that are requested following the
3 voting will not be included. Mr. McCarley will
4 return to the table for Friday's proceedings.

5 With the chair's permission, I would now
6 like to proceed to read the Conflict of Interest
7 Statement for this meeting and the Appointment to
8 Temporary Voting Status for the Panel Consultants.

9 **Conflict of Interest Statement**

10 MS. THORNTON: The following announcement
11 addresses conflict of interest issues associated
12 with this meeting and is made part of the record to
13 preclude even the appearance of an impropriety.

14 The conflict of interest statutes prohibit
15 special government employees from participating in
16 matters that could affect their or their employers'
17 financial interests.

18 To determine if any conflict existed, the
19 agency reviewed the submitted agenda for this
20 meeting and all financial interests reported by the
21 committee participants. The agency has no
22 conflicts to report for today's agenda. In the
23 event that the discussions involve any other
24 products or firms not already on the agenda for
25 which an FDA participant has a financial interest,

1 the participant should excuse him or herself from
2 such involvement and the exclusion will be noted
3 for the record.

4 With respect to all other participants, we
5 ask in the interest of fairness that all persons
6 making statements or presentations disclose any
7 current or previous financial involvement with any
8 firm whose products they may wish to comment upon.

9 **Appointment to Temporary Voting Status**

10 MS. THORNTON: The Appointment to
11 Temporary Voting Status. Pursuant to the authority
12 granted under the Medical Devices Advisory
13 Committee charter dated October 27, 1990, and as
14 amended August 18, 1999, I appoint the following
15 individuals as voting members of the Ophthalmic
16 Devices Panel for this meeting on January 17, 2002:
17 Drs. Allen Ho; Timothy McMahon; Joel Sugar; Anne
18 Coleman; Richard Casey; Janine Smith; and Woodford
19 Van Meter.

20 In addition, I appoint Dr. Jayne Weiss to
21 serve as acting panel chair for the duration of
22 this meeting.

23 For the record, these individuals are
24 special government employees and consultants to
25 this panel or other panels under the Medical

1 Devices Advisory Committee. They have undergone
2 the customary conflict of interest review, and have
3 reviewed the material to be considered at this
4 meeting.

5 Signed, Dr. David W. Feigle, Director,
6 Center for Devices and Radiologic Health, January
7 9, 2002.

8 Thank you.

9 **OPEN PUBLIC HEARING**

10 DR. WEISS: Thank you, Sallie. This now
11 closes this portion, and we're going to continue on
12 to the Open Public Hearing. If anyone has any
13 comments to make, they need to come up to the
14 podium, identify themselves, and any financial
15 conflicts or potential conflicts that they may
16 have.

17 **OPEN COMMITTEE SESSION**

18 DR. WEISS: Seeing no one approach the
19 podium, we will close the public hearing session
20 and move on to the committee session and begin with
21 the FDA Division Update. Dr. Rosenthal. I'm told
22 that Donna Lochner, Chief of the Intraocular and
23 Corneal Implants Branch, has the update.

24 **Branch Updates**

25 MS. LOCHNER: Thank you. I have one

1 announcement of a personnel nature, and that is
2 Ashley Boam, a biomedical engineer in the
3 Intraocular and Corneal Implants Branch, has been
4 temporarily reassigned to the Office of the
5 Commissioner in FDA. She has accepted this six-
6 month assignment in the Office of Planning and
7 Legislation and is working primarily on the
8 Prescription Drug Users Fee Act.

9 We anxiously await her return in July and
10 I'll note that while she is reassigned, she will,
11 however, continue her responsibilities representing
12 FDA on the ophthalmic standards committees, perhaps
13 most notably and importantly the phakic IOL
14 standard committees.

15 Thank you.

16 PMA P010059

17 DR. WEISS: If there is no other
18 information to be updated from the agency, I would
19 like to move ahead to discuss and review the
20 sponsor's PMA P010059. We will begin with the
21 sponsor presentation. The sponsor can approach the
22 podium and there is one hour.

23 I would like each presenter for the
24 sponsor to first identify themselves at the
25 beginning of their presentation.

SPONSOR PRESENTATION

DR. STEINERT: Good morning. My name is Dr. Roger Steinert. I am not the medical monitor on this study. Dr. Howard Fine is the medical monitor. I want to just explain a few things. I have no financial interest in this product. I am not paid to be here. I have never received a cent from Morcher and I never will, as far as I know.

I am here because I was one of the investigators, and Dr. Fine could not be here today. I felt that it was extremely important that we try to focus on the clinical aspects of this implant, how it works, and the results, and when Mr. Welch asked me if I would pinch hit for Dr. Fine, I agreed.

The first time that I saw the data at all was mid-December. So I have put about 50 to 100 hours into this over the holidays and the past couple of weeks, trying to bring this into a form that made sense to me as a surgeon and as a clinician, and I want to convey that to you.

So if the format here is a little different than you might be used to, that's the reason for this.

I'd also like to take this opportunity to

1 thank Ms. Thornton and Ms. Lochner and Mr. Glover
2 and Dr. Lepri from FDA and especially Dr. Sugar and
3 Dr. Van Meter, who were the primary reviewers.

4 I know that the submission was not as
5 clean, to put it mildly, as you would like, and
6 this has been a kind of a difficult task for you,
7 and we are very appreciative of the support you
8 have given to Mr. Welch in allowing us to finally
9 get to this day of presenting to the FDA, so thank
10 you very much.

11 I think it would be helpful to start with
12 the description of the capsular tension ring
13 itself, and just give you a little bit of
14 background.

15 This device was invented a little over ten
16 years ago by Dr. Witschel in Germany. And the
17 purpose has always been in my mind one thing and
18 one thing dominantly, and that is to enhance the
19 mechanical stability of the lens capsule in the
20 presence of weak or absent zonules. That's it.
21 There's been a lot of other stuff connected to this
22 that I think is inappropriate and we're not going
23 to pursue any of those other things. This is what
24 this device is for.

25 So how does it work? Well, the basic

1 mechanical concept is recruitment of adjacent
2 zonules. When you have weak zonules or missing
3 zonules, the idea is to mechanically interconnect
4 other zonules at the equator so that the
5 neighboring zonules provide more support than they
6 would otherwise.

7 Now, I'm going to show you a brief video
8 clip here, and you'll see this is a surgical tape
9 of a phako, and you can see that there are weak
10 zonules--they're not completely absent--to the
11 right on that screen, and you can see how the
12 equator is visible out here at the edge, and you're
13 going to see a brief edited video with phakoing and
14 then implementation of the ring, and you will see
15 the shift in the position of the capsular bag as a
16 result of that.

17 This is the ring itself going in. It's a
18 very simple device, very thin piece of PMMA. This
19 is a manual insertion. There's also a shooter
20 insertion which I use regularly. It makes life a
21 lot easier. And you see as it goes around, you can
22 see how that equator is now closer to the normal
23 position. It's not perfect. This device does not
24 recreate zonules. It simply recruits mechanical
25 stability from the adjacent zonules.

1 So now the ring is in place. I think that
2 was just to demonstrate the position of it by the
3 surgeon, and now we've got a one piece PMMA lens
4 in, and you'll see a before and an after to
5 emphasize that the position of the capsular bag is
6 improved by the presence of the ring and the
7 presence of the implant.

8 So that's basically how it works. Now the
9 next slide is a Meoki [ph] presentation. This is a
10 normal cadaver human eye seen from the posterior
11 side so it's a Meoki view, and you'll see a few
12 things within this. First, this is just the
13 standard posterior view, and in just a moment,
14 we'll get a close up of what's going on in the
15 periphery.

16 You can see zonules out there attached to
17 the capsular bag running this way. Those are the
18 two little eyelets, the beginning and the end, and
19 typically there's a little bit of a space, the
20 relationship of the implant to it.

21 Now, the purpose of this is to demonstrate
22 why--the one thing you'd fear is that this thing
23 would poke through the capsular bag during
24 insertion. You can see actually those ends can be
25 pressed fairly hard against the equator without

1 puncturing, and that's the purpose of that little
2 demonstration.

3 So that's it. That is what this is all
4 about. Now, a brief history of Morcher. It was
5 founded in 1943 as a manufacturer of contact
6 lenses. It began as an IOL manufacturer in 1955 so
7 they've been in this business for a long, long
8 time. It's the Dannheim lens in 1955 and then
9 working with Binkhorst as early as 1958.

10 In 1981, to the best of my knowledge, they
11 were the first IOL manufacturer to use gamma
12 sterilization to improve the biocompatibility and
13 reduce toxicity in IOL sterilization.

14 In 1987, they developed something that
15 they called the compression forge method, which is
16 the thing that I am told allows them to create
17 these rings so that they have a high degree of
18 fracture resistance with a very flexible PMMA.

19 They distribute internationally throughout
20 the world, and as far as I know, they have an
21 excellent track record with the ring, and it is
22 consistent with the highest standards of
23 manufacturing quality.

24 In my personal opinion, there is one
25 indication for the use of the capsule tension ring,

1 and that is stabilization of the crystalline lens
2 capsule in the presence of weak or absent zonules.
3 I'm a believer in keeping things simple, and I
4 think this is what this is all about. Trying to
5 attach other things to this that cannot be
6 substantiated in any easily done clinical study I
7 think is a mistake, and I am told that Morcher and
8 the sponsor are in agreement this is the one
9 indication that we are looking for approval for
10 today.

11 I think typical conditions as guidance to
12 a clinician would be patients with
13 pseudoexfoliation, prior trauma, prior pars plana
14 vitrectomy, and Marfan's Syndrome, but it's not
15 limited to that.

16 Now, the IDE, as many of you know, and the
17 rest of you will hear several times today, occurred
18 in several phases. Phase I was the original study.
19 11 surgeons at five sites who are referred to as
20 the core group were allowed to enroll 75 eyes, and
21 there has now been a minimum two-year follow-up
22 period on most but not all of those patients.

23 There has been difficulty with follow-up
24 because a lot of these patients are referred from a
25 distance and it is impossible to extract

1 information out of the following originating
2 ophthalmologists at a distance in some cases.

3 This is the core group. Dr. Fine, as I
4 mentioned, is the medical monitor; Dr. Garbow; Dick
5 Lindstrom's group in Minneapolis; Bobby Osher's
6 group in Cincinnati; and myself.

7 Now, Phase II had two groups. One was the
8 same core group allowing additional enrollment.
9 Ultimately, in Phase II were 202 patients and 240
10 eyes. And further independent investigators,
11 ultimately totaling 43, who were begging for the
12 ability to use this device for patients who needed
13 it, were allowed to implant under the auspices of
14 this sponsor's study 204 patients and 225 eyes.

15 There is a Phase III. Now you won't be
16 hearing data on this because this is not part of
17 the submission and wasn't required, but you should
18 know that the core group of investigators has been
19 allowed to do limited ongoing implantations for
20 patients who need it and that has resulted in 54
21 further implants at four sites.

22 So, first I'd like to talk about efficacy.
23 I think the appropriate primary measures of these:
24 does it help the IOL center; is it stable in the
25 long term; and does it reduce vitreous loss at

1 surgery? These are very high risk patients by and
2 large who are at risk for vitreous loss.

3 So let's talk about each of these in turn.
4 IOL centration. Well, I asked about reportable
5 clinically detectable decentration. And let me say
6 right up front, this is a real problem if you've
7 approached this in a rigid scientific way. We
8 don't have great methodology for determining IOL
9 centration. And I don't believe any of the IOL
10 studies submitted by other sponsors just for
11 conventional IOLs where they talk about centration
12 do anything differently.

13 It's very subjective. I'd love to see a
14 practical clinical test that we could export to the
15 field that would allow us to really figure out
16 where the center of an implant is and measure it at
17 a millimeter level of accuracy, but we don't have
18 it, and it is absolutely true that as long as you
19 can't see something shifting in the pupillary zone,
20 you won't know whether it's moved.

21 So could these rings be moving one or two
22 millimeters? Could the implant be moving one or
23 two millimeters? Yes. Without detection?
24 Absolutely. And when the investigator says it's
25 one millimeter decentered, how do they determine

1 that? They don't have digitized photos. We know
2 that. This is the world we live in.

3 This is as best as we can do is to ask the
4 surgeons do they see detectable decentration? So
5 in Phase I, the core group, five out of 50 of the
6 reported patients at two years, were reported to
7 have some clinically detectable decentration. In
8 Phase II, the core group was 12 out of 157 at one
9 year, which was the requested follow-up interval
10 resulting in a rate of 7.6 percent, and for the
11 independent investigators in Phase II, it was seven
12 out of 109 at one year, or 6.4 percent, reporting
13 some clinically detectable decentration.

14 What about long-term stability? There
15 were nine reports of decentration of IOLs, and at
16 the last report, one of them was said to be two
17 millimeters and eight eyes were said to be one
18 millimeter or less. That's the total number of
19 reports of decentered IOLs. Now, I don't believe
20 that that's the total amount of decentration at any
21 level.

22 But this is what the surgeons reported and
23 so it must be what they perceive as at least
24 clinically detectable.

25 There have been no reports of extrusion of

1 the ring from the capsular bag. It stays in the
2 bag.

3 What about vitrectomy? You have to
4 remember these are high risk patients with bad
5 zonules, and we don't know. Unless somebody does a
6 prospective randomized study, nobody will know what
7 the incidence of vitrectomy will be with and
8 without the ring.

9 So the best we can say is that the
10 expected incidence approaches 100 percent. Core
11 group, Phase I, eight out of 75, or 13.3 percent.
12 Core group, Phase II, 19 out of 240, or 7.9
13 percent. And independent investigators Phase II,
14 17 out of 225, or 7.6 percent.

15 What about visual acuity issues? Well,
16 we're going to talk about visual acuity under
17 safety, but visual acuity is not an appropriate
18 efficacy outcome measure of the capsule tension
19 ring. That's not what this ring does. It's not an
20 intraocular lens. And comparison of the results to
21 the FDA grid for our results is irrelevant.

22 The cases in which the ring are implanted,
23 and I'm sure you all understand this, are selected
24 for a high degree of pre-op pathology and intra-op
25 pathology. These are high risk patients.

1 So let's talk about safety. I think the
2 primary measures were these: stability after YAG
3 laser capsulotomy; evidence of inflammation;
4 explantations; people who had best corrected visual
5 acuities less than 20/40; and other relevant post-
6 op pathology. So we'll address each of those in
7 turn.

8 Stability after YAG laser capsulotomy. We
9 went through and pulled out reports of anything
10 that was report of more decentration after YAG
11 capsulotomy than before YAG capsulotomy, and we
12 found three reports of possible new or increased
13 decentration.

14 And this is what we have to deal with for
15 reports from the doctors. One eye was reported as
16 slight pre-YAG, whatever slight means. The YAG was
17 done at four months post-op, and all the post-op
18 reports report the IOL as being two millimeters
19 decentered. It has not required reoperation.

20 One eye was reported as one millimeter
21 decentered at the first report after YAG at 10 to
22 14 weeks. All subsequent reports failed to report
23 any decentration on that eye.

24 And one eye had a very complex procedure
25 with cutting of vitreous strands with the YAG laser

1 and an anterior capsulotomy at seven weeks.
2 Decentration was reported as two millimeters at ten
3 to 14 weeks, one-half millimeter at 22 to 26 weeks,
4 no decentration 11 to 13 months, and one millimeter
5 at 23 to 25 months.

6 And those are the only reports of any
7 change in position after YAG capsulotomy. And
8 again, no cases of extrusion of the ring after
9 laser capsulotomy.

10 What about inflammation? Well, the FDA
11 has raised issues of biocompatibility of the PMMA
12 used in the ring. So we looked at possible
13 correlations with reports of inflammation, and if
14 you look at iritis, the incidence, any occurrence,
15 it was six patients, or 1.2 percent, reported as
16 having iritis at any post-op interval; zero at the
17 last reporting interval.

18 And for CME, there's 11, or 2.1 percent,
19 incidence at any time, and four, or .76 percent,
20 persisting at the last reporting interval.

21 Frankly, I'm stunned that it's that low.
22 These are very complicated cases. A lot of
23 vitrectomy is being done.

24 Now, technical problems with the ring.
25 There were 540 total implants. Three were reported

1 as having broken eyelets. As I said, none were
2 reported as extruding post-op. No surgeon felt
3 that there were complications attributable to the
4 ring. There were no infections, and there were no
5 adverse events that the surgeons felt were
6 attributable to the ring.

7 There were four cases where the ring could
8 not be fixated in the bag at the time of surgery.
9 And therefore it was not left in the eye.

10 So let's talk about ring explantations.
11 Again, 540 total implants. There were eight
12 explanations for a rate of 1.5 percent. Seven of
13 those eight were during the primary surgery. Two
14 of them were due to procedural issues. Four of
15 them I've already referenced were due to inadequate
16 capsule or zonules to support the ring, and one was
17 because the surgeon didn't feel it was the correct
18 ring size.

19 There are three slightly different sized
20 rings depending on level of myopia and the size of
21 the capsular bag.

22 And there was one post-op explanation. It
23 was a reintervention at one week post-op. The ring
24 along with the IOL was removed due to the judgment
25 that the whole capsular complex was unstable.

1 Again, I want to emphasize as somebody who
2 has used this, this doesn't manufacture zonules.
3 And there absolutely are people who have way to few
4 zonules for this ring to rescue their situation.
5 So there is surgical judgment involved, and it's a
6 learning curve, and sometimes you put it in, and it
7 doesn't work. There's always that potential.

8 Now, let's talk about retinal detachments.
9 In Phase I core, there were three RDs; in Phase II
10 core, five retinal detachments reported; and the
11 independent investigators reported no retinal
12 detachments.

13 Of those eight detachments, five were
14 present pre-op; two were detected immediately post-
15 op, at the first post-op interval. It is unclear,
16 but raises the question as to whether these were
17 also present pre-op. And there was one that
18 definitely occurred post-op at the two year post-op
19 interval.

20 Other major posterior pathology is as
21 follows: early phthisis was reported in one
22 patient. We went back and looked at that. That
23 was a patient who had one of these pre-op total
24 retinal detachments and light perception vision.
25 The lens had been removed in order to visualize the

1 retina to see whether it could be repaired. The
2 patient never did get reattachment of the retina
3 and eventually started to fade into phthisis at the
4 last report.

5 One patient had a vitreous hemorrhage that
6 was present post-op, and one patient was reported
7 as having a branch retinal vein occlusion. It was
8 not there at the one to two week report and at the
9 10 to 14 week, it was, and it was detected due to
10 the drop of best corrected acuity from 20/25 to
11 count fingers. And the surgeon did not think there
12 was any plausible connection between the branch
13 retinal vein occlusion and the presence of the
14 ring.

15 Now what about the people who lose, not
16 even necessarily lost, but let's say failed to gain
17 acuity at the level of 20/40 or better? So these
18 are people whose post-op best corrected acuities
19 were less than 20/40.

20 In yellow, we have the raw number, and
21 then white is percentages. And you can see that
22 the highest levels were for age-related macular
23 degeneration. And then other macular issues.
24 Typically that was things like traumatic
25 maculopathy and epiretinal membranes.

1 Some of these were the retinal detachment
2 patients. There are a fair number who had
3 posterior capsule opacity, but had not undergone
4 YAG capsulotomy at the time of the reporting
5 interval. Some with irregular corneas, largely
6 dryness. Allegedly, only a couple of CMEs
7 responsible for less than 20/40. One patient with
8 optic atrophy. A couple of people who were said to
9 have severe glaucoma. Two with diabetic
10 maculopathy, and then a very small number of
11 miscellaneous patients.

12 What about glaucoma? Well, glaucoma was
13 reported by two in Phase I, nine core patients in
14 Phase II, and six of the independent investigators.

15 All of the Phase I core patients reporting
16 glaucoma, it was preexisting preoperatively. And
17 Phase II, eight of the nine had preexisting
18 glaucoma. One of the nine was an acute post-op
19 pressure elevation that was treated, and the
20 glaucoma in this case, the definition was IOP
21 requiring medication. So this was reported, but it
22 was only the first post-op day and was gone
23 thereafter.

24 In the Phase II independent patients, two
25 of them it was preexisting. One it was first day

1 post-op only. Two were early post-op only, and
2 then there was one lost to follow-up--and one of
3 those two were lost to follow-up, and then one
4 we're pending longer follow-up reports on and have
5 not received it from the independent investigator.

6 These are the worldwide sales of the
7 capsule tension ring. I thought this would be
8 interesting to you to get a sense of how often or
9 more precisely how infrequently the ring is used.

10 These are sales. No one has figures on
11 actual implantation. So you could guess maybe 50
12 percent of these actually get implanted. And you
13 can see when you consider worldwide cataract
14 surgery of many millions a year, this is not a
15 large number. This is a device restricted to
16 patients who are very specific and have a very
17 unusual but very needy condition.

18 So, in conclusion, the Morcher capsule
19 tension ring has been in use for a decade
20 internationally. It's available throughout the
21 world. It enjoys consistently positive clinical
22 reports, absence of complications attributable to
23 the ring, and a track record of long-term stability
24 and biocompatibility throughout the world.

25 The U.S. clinical trials under this IDE, I

1 think, reflect the positive experience that has
2 been present worldwide with the Morcher ring. The
3 capsule tension ring in my opinion effectively
4 stabilizes the capsular bag in cases of weak or
5 partially absent zonules, and it reduces the rate
6 of serious complication such as vitreous loss,
7 dislocation of the nucleus posteriorly and
8 inability to implant a PC IOL.

9 No safety concerns about the ring have
10 arisen in the course of this trial, and there is no
11 alternative device or technique to achieve these
12 clinical objectives. Thank you very much for your
13 attention.

14 DR. WEISS: If that ends the sponsor's
15 presentation, Dr. Steinert, I'll ask you to stay at
16 the table, and we'll have 15 minutes of questions
17 from the panel for you, and then we'll have the FDA
18 presentation.

19 **Panel Questions for the Sponsor**

20 DR. WEISS: Dr. Sugar.

21 DR. SUGAR: Thank you, Jayne. This is
22 Joel Sugar. I'd like to thank Roger for his
23 candor. I have a bunch of questions and stop me if
24 these are out of the range of what I'm supposed to
25 ask now.

1 First of all, the indications for
2 implanting the lens and the numbers of patients are
3 things that have confused me throughout my review
4 of the data that's been presented and represented.

5 And Roger just presented in Phase II, I
6 guess, two independent, that there were 225 eyes
7 and 204 patients. The information presented to me
8 said that there were 241 eyes and 215 patients. I
9 can understand if the cutoff date or the date of
10 freezing of the data was changed, that the numbers
11 would increase. I can't understand the numbers
12 decreasing.

13 I assume that we work under the principle
14 that once randomized, once assigned in a study,
15 always analyzed. So I guess if you could begin
16 with that.

17 MR. WELCH: When we did the revision, in
18 order to--

19 MS. THORNTON: Excuse me. Could you
20 identify yourself, please?

21 MR. WELCH: Certainly, sorry. My name is
22 Hillard Welch. I am the U.S. representative for
23 Morcher, Stuttgart, Germany. When we did the
24 revision of the statistics and the data, we did it
25 against a different date, and the original

1 submission was a random one unfortunately. That
2 was my error in compiling it because I picked
3 different dates when I shut off various parts of
4 the tabulation.

5 And we finally settled on a date of
6 October 1, and all of those figures that you're now
7 referring to, the 225 and the 204, are based on
8 that date, and the data that had been received as
9 of that date. So that is the figure you should use
10 and not the preceding one.

11 DR. SUGAR: So the number got smaller
12 because some in the original submission didn't--

13 MR. WELCH: Yeah, they should not have
14 been included--

15 DR. SUGAR: --submit an update.

16 MR. WELCH: --in part because when the
17 tabulation was originally done, it picked up a
18 variant of the ring which is not included in this
19 study.

20 DR. SUGAR: Thank you. Can I continue?

21 DR. WEISS: Yes. Dr. Sugar.

22 DR. SUGAR: Joel Sugar. Another question
23 that comes up: the indications were never clear to
24 me. That is many patients who had
25 pseudoexfoliation. All patients supposedly had

1 cataract. Yet, in the submission--again, this data
2 was not reviewed by Roger--about 44 percent I think
3 in the core group had an acuity of 20/40 or better
4 preoperatively. Could that be explained to me?

5 DR. STEINERT: I was troubled by that as
6 well, Joel. This is Roger Steinert again. And I
7 haven't had the ability to extract all the
8 information on all those patients, but I understood
9 that a lot of those patients were, in fact,
10 patients implanted by Dr. Fine, so I specifically
11 got information from Dr. Fine. So that represents
12 a subset to be sure.

13 Almost all of those patients who were
14 20/40 or better had significant glare problems and
15 documented glare acuities in the 20/60 to 20/80
16 range generally. To the best of my knowledge,
17 there was one that was used in the course of an IOL
18 exchange, and I think there were one or two used in
19 the course of a clear lensectomy for high myopia
20 where the zonules were then judged to be
21 suboptimal.

22 But I think the vast, vast majority were
23 patients who had glare decrement in their acuity
24 and did have cataracts.

25 DR. SUGAR: Can I follow up on that?

1 DR. WEISS: Yes, Dr. Sugar.

2 DR. SUGAR: So the indications were not
3 just indications listed in the submission? That is
4 some patients had clear lens extraction for myopia
5 in this study?

6 DR. STEINERT: Apparently a few were put
7 in in patients who had clear lensectomy. That's
8 right.

9 DR. SUGAR: As you know, and I assume as
10 you experience as well, that makes it difficult to
11 assess when we're not given all the information on
12 the indications for the study.

13 DR. WEISS: I would remind if everyone can
14 identify themselves before speaking into the
15 microphone.

16 MR. WELCH: My name is Hillard Welch
17 again. I'm usually referred to as Hid, so I'm
18 going to give it to you that way each time. It
19 will simplify things.

20 DR. SUGAR: Hid often is used to mean
21 hidden.

22 MR. WELCH: Beg your pardon?

23 DR. SUGAR: I'm sorry.

24 [Laughter.]

25 MR. WELCH: I missed that. The question

1 again had to do with the indications that were
2 recorded at the pre-op.

3 DR. SUGAR: The indications for entry into
4 the study to the best of my understanding of the
5 original submission did not include myopia with
6 clear lens.

7 MR. WELCH: That's correct.

8 DR. SUGAR: But Dr. Steinert just told us
9 that some of the patients had this.

10 MR. WELCH: This was a notation made by
11 Dr. Fine on a couple of the patients. He had other
12 inclusion criteria that he used in enrolling those
13 particular patients. In the reference that Dr.
14 Steinert made to the review, there were 70 percent
15 of those patients exhibited an incidence of glare
16 and an inability to drive, an inability to read
17 small print. These were some of the additional
18 qualifications that Dr. Fine used in evaluating the
19 patient for inclusion.

20 DR. WEISS: Alice.

21 DR. MATOBA: Alice Matoba. My question
22 also is in regard to the inclusion criteria.
23 Presence of cataract is one of the inclusion
24 criteria listed, and in Volume II, page nine, the
25 sponsor states that the presence of cataract alone

1 could be an inclusion alone, that alone. If that's
2 true, I wonder how many patients were entered into
3 the study for that criterion alone and how you can
4 then say that these were all patients at risk, at
5 high risk?

6 MR. WELCH: I'm not sure I understood the
7 question. You want to know how many were--if
8 cataract alone was an inclusion criteria, how many
9 were--

10 DR. MATOBA: Cataract alone is--presence
11 of cataract is listed as one of the inclusion
12 criteria.

13 MR. WELCH: Yes.

14 DR. MATOBA: And the sponsor has stated
15 that that could stand alone as an inclusion
16 criterion to enter someone into the study. If that
17 is true, I wonder how many patients were entered
18 with just that inclusion criterion and if so how
19 can you state that these patients were all at high
20 risk?

21 MR. WELCH: I don't think I can give you--
22 my name is Hid Welch--and I don't think I can give
23 you a specific answer to that in terms of numbers,
24 but, yes, cataract was listed in the manner in
25 which it was in the original protocol.

1 I understand the question is concerning
2 did anybody get enrolled just because of a cataract
3 really? That's a different interpretation, and the
4 answer to that would be no. And I'd have to go
5 back in order to provide you with the specifics as
6 to what was the other inclusion criteria for that
7 particular patient. That is all in the database.

8 I can't pull numbers out for you right now
9 to say that there were so many that had this, that,
10 or the other thing, but the--yes. To the best of
11 my knowledge of reviewing the cases, there are no
12 instances of a single criteria for inclusion.

13 DR. BRADLEY: This is Arthur Bradley.
14 Just for clarification, then, is that an error then
15 in the report?

16 MR. WELCH: I beg your pardon?

17 DR. BRADLEY: Just following up on Alice
18 Matoba's question, is that an error then in the
19 report, because the report does state--I saw it
20 myself--that cataract alone is an inclusion
21 criteria.

22 MR. WELCH: It does?

23 DR. STEINERT: Can you refer us to exactly
24 what you're looking at?

25 DR. MATOBA: Well, let's see. Volume II,

1 page nine of 22.

2 MR. WELCH: Page one of which?

3 DR. STEINERT: I'm sorry. Dr. Matoba,
4 could you give us that number again?

5 DR. MATOBA: Page nine of 22 on Volume II
6 is what I've written down, Exhibit 8.

7 DR. STEINERT: You say nine of 22?

8 DR. BRADLEY: Page nine.

9 DR. STEINERT: Page nine on Volume II.

10 DR. MATOBA: Yes. Exhibit 8.

11 MR. WELCH: And that may be an error. It
12 is true that that is what I put in the initial
13 response, and I would have to admit that that's
14 probably an error because I don't think that is
15 correct. I believe there was always an additional
16 condition even though it does state--thank you--it
17 does state that cataract is a single inclusion
18 criteria.

19 But my memory is that that is not a
20 correct final interpretation. Hold on. I'll look
21 for that.

22 DR. WEISS: I think Dr. McMahon had a
23 comment.

24 DR. McMAHON: Tim McMahon. My
25 understanding from the submission is that there

1 were 40 patients that had two eyes in implanted
2 rings; is that correct?

3 MR. WELCH: What was that again?

4 DR. McMAHON: My understanding from the
5 submission, that there were 40 patients that had
6 two eyes where rings were implanted? Right and
7 left. Is that correct?

8 DR. STEINERT: Yeah, there bilateral
9 implants.

10 MR. WELCH: Yes.

11 DR. STEINERT: Yes, there were patients
12 who were bilaterally implanted.

13 DR. McMAHON: And that protocol was agreed
14 to by the FDA to do second eye in an investigative
15 device?

16 MR. WELCH: It was never so stated as a
17 separate condition, no. No. At no time, though,
18 there was recognition on the part that there were
19 bilateral implants.

20 DR. WEISS: Dr. Grimmer.

21 DR. GRIMMETT: Michael Grimmer. I have
22 just two questions at this time. Number one,
23 regarding one of your slides, Roger, that you had
24 up regarding best corrected visual acuity loss. In
25 looking at the Phase II core and Phase II

1 independent, just roughly eyeballing, adding up the
2 percentages of best corrected loss, worse than
3 20/40, looks like they're adding up Phase II core
4 is 15 to 17 percent or something like that.

5 In other area of the study, best corrected
6 visual acuity loss, worse than 20/40 was up near 40
7 percent in one of the data tabulations. So your
8 slide looks like it's missing 20 to 23 percent or
9 something of the causes. Do the rest of those best
10 corrected visual acuity loss remain under-
11 determined?

12 DR. STEINERT: Those, no. To the best of
13 my knowledge, what I presented to you was supposed
14 to be the total number. So I don't know. What is
15 the table that shows 40 percent being worse than
16 20/40? Can you direct us to that?

17 DR. GRIMMETT: I'll look it up. There
18 were so many different tables in the submission.

19 DR. STEINERT: I know.

20 DR. GRIMMETT: That I got confused. So
21 I'll look that up. My second question is a
22 procedural one. In Volume I, Tab Exhibit C, page
23 two, under the Operative Methodology, it states
24 that the intercapsular ring would be implanted just
25 after tearing of the capsular rexis. This

1 insertion would be prior to the hydra dissection,
2 hydra delineation and phakomulsification.

3 The video that you showed which I think
4 showed phakomulsification of the lens, removal of
5 the entire crystalline lens, and then implantation
6 of the ring, how would one insert the ring before
7 hydro dissection and hydro delineation? How is
8 that possible?

9 DR. STEINERT: Before hydro dissection and
10 hydro delineation?

11 DR. GRIMMETT: Yeah, because that's what
12 it says in Volume I.

13 DR. STEINERT: Yeah. Well, first of all,
14 that video segment is not from any of the
15 investigators. That actually was from Germany just
16 because we could get our hands on it quickly.

17 DR. GRIMMETT: Okay.

18 DR. STEINERT: But that aside, from
19 practical point of view, I know what really went
20 on, and what went on is that as we got experience
21 with the ring, it becomes apparent that the later
22 you can put it in in the case, the easier your life
23 is.

24 In some cases of extreme laxity of
25 zonules, you're lucky to get through the capsular

1 rexis, and you want, you need stability
2 immediately.

3 So the very next thing done is the
4 implantation of the ring. The ring will because of
5 its forces, will basically act like a hydra
6 dissector. It will find the equator just because
7 of its outward pressure. So you can insert it
8 under the anterior capsule prior to hydra
9 dissection, and it will nevertheless end up out at
10 the equator.

11 However, given my choice as a surgeon, I
12 always deferred it as long as I could, and
13 sometimes I'd be part way through the phako, and
14 then say, well, this is clearly starting to
15 unzipper on me; I need to put it in. But because
16 it does tend, it has the potential for trapping
17 some cortex between the ring and the equator, and
18 then making cortical stripping more difficult, it
19 is desirable to defer the implantation as far into
20 the case as possible.

21 DR. GRIMMETT: Thank you.

22 DR. WEISS: Dr. Rosenthal and then Dr.
23 Matoba.

24 DR. ROSENTHAL: Thank you. Ralph
25 Rosenthal. I just wanted to clarify to Dr. McMahon

1 that usually at the beginning of an IDE, we agree
2 to monocular implantation or monocular treatment,
3 and then as we become more comfortable with the
4 device and its performance, we will allow the
5 sponsor to move into bilateral implantation or
6 bilateral treatment.

7 So this IDE has gone on for five years,
8 and so over that five year period, we certainly--
9 I'm not sure at what point in the five year period,
10 we agreed to the second eye as being implanted, but
11 we had confidence based on the annual reports from
12 the sponsor that there were no problems with the
13 device.

14 DR. McMAHON: Thank you for clarifying
15 that.

16 DR. WEISS: Dr. Matoba.

17 DR. MATOBA: Alice Matoba. I have two
18 questions. First, in your protocol, you specify
19 the range of dates at which the follow-up visits
20 should have to occur. Did you specify when the
21 patients had to be dilated post-op?

22 MR. WELCH: No. Hid Welch. The answer is
23 no, we didn't specify dilation in the protocol at
24 any specific period.

25 DR. MATOBA: Then it seems to me that in

1 undilated people, it would be very difficult to see
2 the lens decentration, especially if it were a
3 small amount, and so I wanted to know how you could
4 have any confidence in your data that long-term
5 stability or that the decentration rate was very
6 low during the follow-up period.

7 MR. WELCH: You mean without a requirement
8 of dilation?

9 DR. MATOBA: Uh-huh. Without knowing
10 whether the patients were dilated or not during the
11 follow-up period.

12 MR. WELCH: We may have made an inaccurate
13 assumption, but there are a few instances in case
14 reports where the examination is noted as not
15 dilated. And thus, the inference is that the
16 others were under dilation at the time of the
17 report.

18 We did not, and the protocol doesn't
19 specify, that there be dilation at every exam. But
20 just with the way the original protocol--I did not
21 write the original protocol, and I think maybe that
22 should be explained. I barely got on this train
23 after it had left the station. And I picked it up
24 and ran it.

25 As a consequence, there may have been some

1 things I should have stopped and gone back and
2 redone in order such as you're asking now to be
3 more definitive, but they were not done, and we
4 continued on the track as it had originally been
5 established, and thus there was no requirement for
6 dilation or no stipulation within the protocol.

7 DR. STEINERT: This is Roger Steinert.
8 Dr. Matoba, first of all, I totally agree with you
9 that in retrospect that would have been a good
10 thing to specify because it would have improved the
11 ability to see what was going on.

12 The clinical reality, as I tried to
13 indicate, is we all know, that many of these
14 patients won't dilate beyond if you're lucky five
15 or six millimeters, and so even then we're not
16 going to pick up all levels of decentration, and
17 furthermore we have no decent truly scientific way
18 of even measuring decentration anyway on a clinical
19 basis in clinical practice.

20 So this is a deficiency. There is no
21 question about it. I think that the minimum
22 statement that you can make is that there was no
23 decentration large enough to become a clinical
24 issue or a clinical problem. That's about all you
25 can say.

1 DR. MATOBA: My second question is in
2 terms of long-term stability, is there any evidence
3 to indicate that the presence of the ring will
4 stabilize the zonules long term? Many of these
5 patients have conditions in which they're really
6 progressive weakening of zonules over time, and so
7 after many years might not the IOL, the whole
8 thing, just become destabilized?

9 Dr. Witschel, I think, had one case where
10 IOL and the ring became subluxed after six years.

11 DR. STEINERT: This is Roger Steinert
12 again. There is no question that pseudoexfoliation
13 in particular and possibly some of the other
14 conditions are associated with progressive ongoing
15 degeneration of zonular integrity over time.

16 And I think all of us who do cataract
17 surgery are seeing patients coming in, sometimes
18 years, even decades, after PC IOL implantation who
19 have lenses that are shifting, dehissing, even
20 falling back into the vitreous against the retina,
21 and that is an issue which we're all going to have
22 to deal with clinically for some time to come.

23 Whether or not the ring can affect that
24 rate is unknown and a study to prove that would
25 probably start to approach and end off the Midas

1 study in terms of difficulty in terms of number of
2 patients enrolled, not to mention the complexity of
3 a five to ten year follow-up.

4 For that reason, I feel that it is
5 inappropriate to make any claim that this ring
6 enhances long-term stability of the capsular bag
7 process. We simply don't have data to support such
8 a claim.

9 On the other hand, logically, and on a
10 clinical basis, I also cannot conceive that this
11 ring would in any way accelerate decentration, and
12 if you--understanding how it does take tension off
13 of the zonules and get some recruitment from
14 adjacent zonules mechanically, logically one would
15 think it would slow down that degenerative process,
16 but it's certainly not going to stop it.

17 DR. MATOBA: My concern is just that
18 having the lens might, as a crutch, might encourage
19 the implantation of IOLs in some patients who they
20 should not be implanted whether or not they were
21 agreed to help stabilize the IOL.

22 DR. STEINERT: Gee whiz. I mean how can
23 you legislate being smart? You know it's a
24 judgment question, and there will be errors in
25 judgment, and I agree with you. But I don't think

1 --certainly--and that's part of why I wanted to
2 show those numbers of worldwide sales. There is no
3 evidence that this thing has become, you know,
4 everybody's favorite play thing and gets implanted
5 willy-nilly in every single lens case or anything
6 close to it.

7 It slows you down. It adds cost to the
8 case, and it adds surgical time. So I think there
9 are some significant natural barriers to
10 inappropriate use of the ring.

11 DR. WEISS: We have Dr. Grimmett, Dr. Van
12 Meter, and then Dr. Smith.

13 DR. GRIMMETT: Mike Grimmett. Just in
14 follow-up to my best corrected visual acuity
15 statement to Dr. Steinert. The numbers I was
16 quoting of the 40 percent worse than 20/40 best
17 corrected vision actually came from Dr. Lepri's
18 review, page 15, his amended review, as where he
19 tabulated the numbers again. It doesn't, I don't
20 think, agree with the summary slide you had out.
21 There's about half of the patients apparently
22 missing if these numbers are correct.

23 DR. STEINERT: Okay. Certainly--this is
24 Roger Steinert--as I said, the numbers I presented
25 were the numbers I got from Mr. Welch. If there

1 are tables that disagree, they should be reconciled
2 and explained. Absolutely.

3 DR. GRIMMETT: An additional question I
4 had, I didn't locate a physician information
5 booklet typical of other PMAs, and I was just
6 curious regarding the three ring sizes, how does
7 one clinically go about measuring the appropriate
8 width of a capsule diameter to pick the appropriate
9 ring size? Just as a clinician, how do you do
10 that?

11 DR. STEINERT: Good question. Really
12 there are three ring sizes. What has tended to
13 evolve is I think the majority of people use the
14 average ring because that is an issue. It is not
15 measurable. The one suspicion many people have is
16 that high myops or perhaps extremely advanced large
17 cataracts may have larger bags. So one of the ring
18 sizes is a larger diameter.

19 The reason that is not used routinely on
20 all is that then it is too big for the average
21 capsule, and I think, although it can be inserted,
22 it makes life more difficult. So I think most
23 surgeons have gravitated toward a strategy of using
24 the middle sized ring, you know, the Mama bear, the
25 Papa bear, and the Baby bear, and they go for the

1 middle to the one that's just right. And in the
2 vast majority of cases, that works.

3 DR. GRIMMETT: Mike Grimmett again, just
4 as a final comment. If in general use, I think it
5 would be beneficial for the sponsor to have some
6 type of comments to guide the average practicing
7 ophthalmologist as to how to select the ring size
8 or something of that nature.

9 MR. WELCH: Hid Welch responding to that.
10 I've noted in my response to the FDA, which they
11 will receive, that we will look at collecting such
12 information and publishing it. Unfortunately, it's
13 of little value in the package insert. It's got to
14 be done educationally on a broad basis because
15 otherwise you get to the point of insertion and you
16 open the package and it is too late.

17 You're not going to get the information
18 you need at that point for any size determination.
19 So we will look at how we can collect such
20 information and publish it.

21 DR. STEINERT: Well, certainly I think--
22 this is Roger Steinert again--that the manufacturer
23 should do that, and all of us--in fact, there is an
24 intention among the investigators to publish not
25 only the data but a surgical procedure and what

1 we've learned along the way in terms of guidance as
2 a separate document in the peer review literature.

3 However, to the extent that the FDA wishes
4 some guidelines in the package insert, I would be
5 pleased to get that far as to be working with them
6 on that. We can do that.

7 DR. WEISS: Dr. Rosenthal.

8 DR. ROSENTHAL: Yeah. Let me address two
9 issues that the panel has raised. The first has to
10 do with package insert or labeling, and certainly
11 we would appreciate whatever recommendations the
12 panel would have concerning the labeling of the
13 device.

14 The second issue had to do with
15 inappropriate implantation. As you know--or
16 inappropriate use--we would appreciate from the
17 panel some idea as to when it would be most
18 appropriate to implant and when it would be
19 contraindicated, and you might want to choose a
20 percentage of zonular abnormality or something in
21 which the labeling would then state that this was
22 the appropriate time to use it.

23 But, of course, as you well know, the
24 practice of medicine kicks in. Once a device is
25 approved and a physician has the opportunity and

1 the right to use any approved device as they so see
2 fit. But certainly labeling and making a statement
3 about contraindications, precautions and so forth,
4 would be appropriate if you felt that that was the
5 thing to do.

6 DR. STEINERT: This is Roger Steinert. If
7 I could expand on that. You're totally right, Dr.
8 Rosenthal, and I know from talking to other
9 surgeons as well as my own experience that there is
10 almost an inevitable tendency to underestimate the
11 amount of zonule loss in trauma cases. No matter
12 what you think you see pre-op, it will be worse
13 once you get in there.

14 So, some cautionary statements about not
15 getting overly enthused and also not confusing this
16 ring with something that makes zonules grow back is
17 very, very important. There is a point where there
18 just aren't enough zonules, and I still work
19 closely with my vitreal retinal colleagues to do
20 planned pars plana lensectomies and sutured PC
21 lenses in some of the cases that are referred to me
22 for the ring, because they just don't have enough
23 zonules.

24 DR. WEISS: Dr. Rosenthal.

25 DR. ROSENTHAL: Rosenthal. I think, you

1 know, you've made a correct point, Dr. Steinert,
2 and that is there is enough option to this ring.
3 I'm not a cataract surgeon, but I think it is a
4 complex option, a pars plana vitrectomy and a
5 suturing of the lens implant. But there certainly
6 is another option, and that also would have to be
7 spelled out in the labeling so that the physician
8 would have some idea as to the appropriate time at
9 which the device would be used.

10 DR. WEISS: Jayne Weiss. Roger, what
11 percentage of zonules absent are the max that you
12 would try to implant the ring in your own practice?

13 DR. STEINERT: Personally, if I feel that
14 there are more than three to four clock hours of
15 totally absent zonules, I wouldn't go with this
16 ring.

17 Now you are probably all aware there is a
18 modification that Dr. Robert Cionni came up with
19 that involves a little loop to attach a piece of
20 suture to, and that can then hold the ring in one
21 direction. This is not the subject of this PMA
22 application. That will be an issue for a
23 supplemental application later.

24 So that ring is not under discussion here,
25 but there is that coming down the pipeline. So

1 that expands potentially the range. But for the
2 device we're talking about today, I would say three
3 to four clock hours of complete absence. The
4 bigger challenge are these people who have partial
5 absence, because you never really know what else is
6 going on elsewhere. Now are the rest of the
7 zonules nice and strong and happy, or are they all
8 damaged from this injury? It's just that they're
9 more damaged in one area. And that's the kind of a
10 thing that you don't discover until you get into
11 the case.

12 DR. WEISS: I'm going to have Dr. Van
13 Meter, Dr. Smith, and then Dr. Ho and then Dr.
14 Coleman.

15 DR. VAN METER: Thank you. I was a
16 primary reviewer for this, and I'd like your
17 reaction to the fact that I feel insulted to have
18 to review data that is as abysmal as this is.
19 There are roving denominators that change. By your
20 own admission, these are patients that are in a
21 referral practice, and it's very difficult to get
22 them back in for examination, which is not the kind
23 of study patients that you want to have.

24 The company has sold 12 to 16,000 of these
25 annually, and the comment was made that they don't

1 know how many were implanted, and there are some
2 gaping holes in this. It's almost an affront to us
3 to have to deal with this data and make some
4 meaningful conclusions on them.

5 Would you please clarify these issues?

6 Dr. Grimmett asked when do you put the ring in.

7 The protocol that we have says the ring is
8 implanted after capsulotomy. Is that incorrect?

9 DR. STEINERT: Yes. The reality is that
10 this protocol was written a long time ago, five or
11 six years ago, by neither of us at this table, and
12 that is what became the subject of the study. As
13 surgical experience evolved, it was discovered
14 that's just not a smart thing to do.

15 DR. VAN METER: All right. So that's not
16 the case.

17 DR. STEINERT: That's the earliest it
18 would be implanted, but not necessarily at that
19 point. That's correct.

20 DR. VAN METER: Okay. One of your slides
21 showed that 100 percent of these patients would
22 likely go on to vitrectomy if they didn't have the
23 ring implanted. The leading indication was for
24 pseudoexfoliation, and many patients with
25 pseudoexfoliation can successfully have a cataract

1 extraction with a lens implant and do well.

2 Can you put these two pieces together?

3 DR. STEINERT: Yes. The real question is
4 what is the surgeon's judgment? Now, all I can
5 tell you is that we can't go into the operating
6 room and pass judgment on every single case at that
7 time.

8 The intention and the discussion among the
9 investigators was that it would be cases of
10 pseudoexfoliation when there was evidence of laxity
11 of the zonules which sometimes becomes obvious
12 right away when you start your capsulotomy.
13 Sometimes it becomes obvious further into the case.

14 Not for the routine use in a patient just
15 because of the presence of pseudoexfoliation
16 material on the anterior lens capsule. So now
17 whether that was complied with, I have no ability
18 to tell you. I don't know.

19 DR. VAN METER: All right. One other
20 question. There were 25 patients that had 20/20
21 preoperative vision admitted into the study, and
22 Mr. Welch stated that 70 percent of the patients
23 that presumably Dr. Garbow did had preoperative
24 glare, means 30 percent of them do not have
25 preoperative glare.

1 So if you've got 25 patients that probably
2 didn't have any preoperative glare symptoms, were
3 these done for high myopia? Or do we have any way
4 of knowing? I mean is this--

5 MR. WELCH: Hid Welch. Not that I'm aware
6 of. I don't remember any statistic or data that
7 would show that or that recorded that. There are
8 indications given on the report that we asked Dr.
9 Fine to provide, and they range from all kinds of
10 things, not just glare.

11 And most, some of them even called it off
12 for quality of life because the patient was
13 continuously complaining of the inability to do
14 whatever it was they wanted to do in their daily
15 life. And these were listed as part of the--

16 DR. VAN METER: Okay. But since this was
17 a device that is used in complicated patients that
18 have a higher than usual risk factor, why are we
19 operating on 20/20 patients with no glare in a
20 procedure that has a higher than usual risk
21 criteria?

22 DR. STEINERT: First of all--this is Roger
23 Steinert--I understand your question and I agree.
24 I had the same reaction this past month when I came
25 across these data. And I did not personally use

1 any of these implants under those circumstances.

2 Since we didn't have the ability to track
3 down every one of those, as I said earlier, we did
4 ask Dr. Fine, who seemed to have done a large
5 number of those, to give an accounting. I'm not
6 sure about the 70 percent number. I haven't done
7 that calculation. Mr. Welch--

8 DR. VAN METER: That came from Mr. Welch.

9 DR. STEINERT: Yes. This said, but my
10 recollection is that Dr. Fine did use this on one
11 or two, and I don't think there were more than
12 that, high myops who are undergoing clear
13 lensectomy for refractive surgical reasons.

14 DR. VAN METER: Okay.

15 DR. STEINERT: And I believe there was one
16 that was an IOL exchange where the capsule seemed
17 to be unstable after the original implant was out.

18 All the rest, they either documented glare
19 or there were these functional complaints that
20 suggested glare and night vision issues, et cetera,
21 et cetera. But they might have 20/20 or 20/25 high
22 contrast acuity.

23 DR. VAN METER: Okay. Thank you. I have
24 another question that I'd like for ya'll to answer
25 if you can. If you envision this round plate as

1 being the lens capsule complex and when you put the
2 ring inside it, it stabilizes it and makes, you
3 know, a round device that supposedly uses all of
4 the zonular fibers for stability.

5 But if you're missing three clock hours,
6 say from 11 o'clock to eight o'clock, of stability,
7 and you have a lens implant and the ring in place,
8 most of the tension of, you know, mobility to the
9 eye is going to be on the 11 o'clock fibers and the
10 eight o'clock fibers. And actually by missing
11 three clock hours, you're going to have increased
12 pressure on those that are on the edge of the gap,
13 if you will.

14 And do you all have any studies to show
15 that this is stable long term? Two years is really
16 not long enough to show what's going to happen. I
17 mean my fear is that this diaphragm is going to
18 tear the 11 o'clock and the eight o'clock fibers
19 and then you have four clock hours of instability
20 from 11:30 to 7:30.

21 And from the data that you have, this is
22 really not anything that's going to show up in a
23 year or two.

24 DR. STEINERT: Well, first of all, I agree
25 with you about any claim about long-term stability.

1 I don't think we have data that can substantiate
2 such an indication, and that's why that is not
3 being asked for.

4 From a mechanistic point of view, I'm not
5 aware of any sophisticated biomechanical studies
6 one way or the other. From a kind of naive
7 conceptual point of view, I think that you're
8 right, that the zonules on the edge of the total
9 defect are going to be the ones under the most
10 pressure.

11 The concept is that in the presence of the
12 ring, though, that at least those--let's say your
13 11 o'clock zonule is getting a little help from the
14 guy at 12:30, which is getting a little help from
15 the one at 12. In the absence of the ring, that's
16 not happening. So it ought to be making the best
17 of a bad situation.

18 DR. VAN METER: Okay. I have one more
19 question, and thank you for taking the time to
20 answer some of these questions that were not
21 necessarily of your doing. But since you're the
22 one that we have to direct these to, forgive my
23 frustration.

24 The number of explantations came from
25 surgeons that put the ring in, and then thought

1 maybe it would be better not to have a ring in.
2 And it shows that the difficulty of preoperative
3 evaluation of zonular stability is one of the main
4 problems that surgeons would have in having this
5 device on hand.

6 Do you have any feel for how many patients
7 were designed to get the device from preoperative
8 planning versus what percentage had the device
9 implanted when the surgeon determined
10 intraoperatively that the ring would be helpful?

11 MR. WELCH: Hid Welch. The answer to that
12 is no, I don't have that information for you here,
13 but there are several instances in the case reports
14 of the decision being made intraoperatively to use
15 the ring.

16 DR. VAN METER: Right.

17 MR. WELCH: And we can segregate those
18 from the database to then evaluate it.

19 DR. VAN METER: Well, I was really asking
20 Dr. Steinert as a practicing cataract surgeon that
21 if he had a feel that, you know, 20 percent of them
22 you decide intraoperatively, ten percent or 50
23 percent? I just--

24 DR. STEINERT: Oh--

25 DR. VAN METER: I don't use the ring so I

1 don't know practically how it shakes down.

2 DR. STEINERT: Yeah. Offhand, I can only
3 think of a cases, Woody, where I truly didn't
4 anticipate anything until I got into surgery.
5 Every once in awhile, you know, there are ones
6 where the patient dilates so much better at
7 surgery, and all of a sudden you say, whoa, I can
8 see the edge of that. I didn't expect that.

9 And then even that, I wouldn't necessarily
10 put a ring in, but then you start manipulating and
11 everything starts moving and you get surprised.
12 But that's a real minority. What you're alluding
13 to, though, and I think is a bigger subgroup, is
14 the pseudoexfoliation group. That's the one where
15 you're most likely to not know going in whether you
16 need it. So what I do is I consent all of my
17 pseudoexfoliation patients in advance and tell them
18 that I want to be able to use this if they need it,
19 but I will only use it if they need it. It ends up
20 being about five percent probably of the
21 pseudoexfoliations.

22 DR. VAN METER: So, by and large, you
23 would order this ring for patients ahead of time?

24 DR. STEINERT: Yeah, we have--since it's
25 not like an implant with power and everything, it's

1 easy to have a stockpile, so you just have a couple
2 lying around.

3 DR. VAN METER: Okay. Thank you.

4 DR. SUGAR: Can I comment on that issue?

5 DR. WEISS: Yeah, I'd prefer to keep
6 discussion of this particular matter to be placed
7 later on in the game, and then now do you have a
8 particular question on this?

9 DR. SUGAR: Well, just that in Mr. Welch's
10 data they presented to us, they listed 133 patients
11 having the decision made intraoperatively in the
12 independent Phase II.

13 MR. WELCH: Would you repeat that, please,
14 sir?

15 DR. SUGAR: In your data that you
16 presented to us--

17 MR. WELCH: Yeah.

18 DR. SUGAR: --there were 133 patients
19 listed as having the decision made
20 intraoperatively. I was going to ask the same
21 question of how were they consented if the decision
22 was not made until the time of the operation?

23 DR. WEISS: Thank you.

24 MR. WELCH: What table is that to which
25 you're referring, sir?

1 DR. SUGAR: I will try to retrieve it.

2 DR. STEINERT: Quite frankly, Joel, part
3 of that answer, I'm sure, relates to the difference
4 between the independent investigator group and the
5 core investigator group and how they approach the
6 study and why they were asking for the ring.

7 DR. WEISS: While we're looking into that,
8 Dr. Rosenthal had a comment.

9 DR. ROSENTHAL: With regard to Dr. Van
10 Meter's comment about greater than two years, this
11 device has been under investigation since 1996.

12 MS. THORNTON: Would you speak into the
13 microphone, Dr. Rosenthal?

14 DR. ROSENTHAL: The device has been under
15 investigation since 1996, so I think there's
16 probably a lot of patients who were enrolled longer
17 than two years.

18 The other thing is as Dr. Steinert noted
19 on his slide, from the years 1992 to 1996, there
20 were several thousand that have been implanted
21 worldwide, and I would hope that we would have a
22 recommendation from the panel regarding this issue
23 of long term.

24 DR. VAN METER: Woody Van Meter. Weren't
25 those devices sold?

1 DR. ROSENTHAL: Oh, sold, sorry. Yes.

2 Well--

3 DR. HO: Allen Ho. Furthermore--

4 DR. WEISS: I think Dr. Smith was first,
5 then Dr. Ho, and then Dr. Coleman. We're going to
6 try to go back to the original.

7 DR. SMITH: Janine Smith. In Volume I,
8 Exhibit C, there is a protocol evaluation listed
9 and a dilated fundus exam is specified in the
10 evaluation process. Do we know that those were
11 performed at least on a certain interval in all of
12 the patients since somebody asked previously how
13 could you determine if there was IOL decentration
14 if you're not certain that the patient was dilated?
15 Are we certain that every patient post-operatively
16 had a dilated fundus exam because there is no place
17 on the data report form to document that?

18 DR. STEINERT: And that's exactly right.
19 It was not called out as a specification. So we
20 have no way of certifying that.

21 DR. SMITH: But it was specified in the
22 protocol that it would be done, on Exhibit C.

23 DR. STEINERT: Which page are we talking?

24 DR. SMITH: Page two of three.

25 DR. STEINERT: Yeah, that's the page I'm

1 on.

2 DR. SMITH: Second paragraph at the top
3 under Evaluation, it describes the evaluations to
4 be performed and dilating fundus exam is listed
5 there.

6 DR. STEINERT: Depending on how you read
7 that, it might be interpreted as just a pre-op
8 dilated fundus exam, I believe. It's ambiguous.
9 It's not great wording.

10 DR. SMITH: Perhaps that referred--well, I
11 don't think that can refer to just preoperative,
12 because later on in the sentence, well, there's a
13 semi-colon there, and then there is intraoperative
14 complications. So maybe that is a pre-op
15 evaluation. It was only specified to be performed
16 preoperatively.

17 DR. STEINERT: Really, it's ambiguous.

18 DR. SMITH: The second question is on
19 Exhibit 2 in Volume II which is a data report form.
20 Under the Pathology and Complications, it lists
21 inflammatory deposits on the IOL and fibrin in the
22 pupil as complications. And Dr. Steinert, on your
23 slide for documenting the rate of inflammation
24 post-operatively, you listed 1.2 percent. Were
25 these the criteria used to determine a diagnosis of

1 iritis, because that's very low, and it states that
2 this time period was one week to closure, and one
3 week after surgery in a complicated case like this,
4 it wouldn't be unusual to have an anterior segment
5 reaction. So I'm wondering if it was just
6 inflammatory deposits on the IOL and fibrin used to
7 define iritis?

8 DR. STEINERT: It, first of all, I should
9 clarify. The acute post-op phenomenon was
10 discounted because that's essentially 100 percent
11 iritis on the first day.

12 So this would be beginning at the--I'm
13 pretty sure--we did this at ten to 14 weeks.

14 DR. SMITH: So there was a specific--
15 Janine Smith--there was a specific time point for
16 iritis to be evaluated?

17 DR. STEINERT: It was a report that--yes,
18 it was part of the post-op report, but that number,
19 to the best of my knowledge, was generated by
20 people specifying there was iritis. They used that
21 word in the post-op report so the presence of a
22 precipitate on an IOL in the absence of anterior
23 chamber reaction would not have been called out as
24 iritis.

25 DR. SMITH: So then since it's not listed

1 as already on the form, they would have checked
2 under "Other" and written iritis?

3 DR. STEINERT: Yes.

4 DR. SMITH: One other question about the
5 method of insertion. In the video, you showed a
6 forceps being used to insert the ring, but you
7 commented that you personally used an injector. Do
8 we have any information on which procedures were
9 used by surgeons in this study?

10 MR. WELCH: Hid Welch. The answer to that
11 is that we do partially because what happened was
12 in the early stages of the study, the doctors found
13 that there was an injector on the market. Then an
14 injector was requested by the Morcher Company from
15 another manufacturer.

16 At that point, we also got a notice from
17 the FDA that the injector was not approved. So we
18 submitted a 510(k) for the injector, which was
19 subsequently approved, and a restriction on that
20 510(k) stipulated that it could be used by the
21 investigators only. It was not to be used for any
22 other surgeon for obvious reasons.

23 The core group, to the best of my
24 knowledge, does have the injector, and once it was
25 approved, they have used it, but they did not note

1 it on the form. We do have a number of case
2 reports where it was noted that an injector was
3 used. We have indications of the benefit of the
4 use of that injector on those case reports, but
5 that's all.

6 DR. WEISS: In the interest of time--
7 Roger, sorry. You had a comment?

8 DR. STEINERT: I was--just to clarify, Dr.
9 Smith, that I've done it both ways, and I can just
10 tell you from personal experience, it is a little
11 easier and a little faster to use the injector in
12 my personal opinion, but I have done it the other
13 way, and I have not encountered any adverse issues.
14 It's not that you break the ring or puncture the
15 capsule. It's just a little--you have to do a hand
16 over hand maneuver. So it's just a little bit more
17 complicated.

18 I also know Dr. Witschel doesn't like the
19 injector. So there's a range of opinion as in many
20 surgical things.

21 DR. WEISS: In the interest of time, we
22 are going to have another question by Dr. Ho and
23 then followed by Dr. Coleman, and then we will
24 move on to the FDA presentation.

25 DR. SUGAR: Could I add the documentation

1 for my comment earlier?

2 DR. WEISS: Yes.

3 DR. SUGAR: In Volume II, Exhibit F-2
4 revised, there are listed for Phase II independent
5 133 intraoperative zonular dehiscences. In Phase
6 II core or phase I, there were no intraoperative
7 dehiscences listed. I don't know how to interpret
8 that, but the only way that I could interpret it
9 was that these were recognized at the time of
10 surgery and not preoperatively, because none of the
11 other ones were listed that way.

12 DR. WEISS: While you're looking at that,
13 maybe Dr. Ho could ask, proceed with his question.

14 DR. HO: Allen Ho. Just a question for
15 Dr. Rosenthal, first of all. The comment was made
16 regarding long-term stability and you had suggested
17 that because this had been implanted for so many
18 years, that there might be some information
19 available to us.

20 I'm not aware of long-term stability
21 information, unless I'm mistaken. Did I miss
22 something?

23 DR. WEISS: I would just remind you the
24 sponsor is not looking for approval for long term.

25 DR. HO: Okay.

1 DR. WEISS: So I don't really think we
2 need to discuss that issue. I think it's been
3 suitably handled.

4 DR. HO: Nor am I comfortable making a
5 comment about that with that data. So that's my
6 first comment.

7 The second comment is that, you know, our
8 charge here is to advise the FDA based on data, and
9 one of the principles of a good study design is to
10 identify, first of all, the patients with whom--in
11 whom you're studying so you can make relevant
12 recommendations to those patients based on the
13 results of the study.

14 It seems to me that the core group is the
15 group of patients that might have the best follow-
16 up, the best accountability, but I'm still at a
17 loss in defining who those patients are. And I'll
18 give you a for example.

19 If you turn to Volume I, Exhibit F, as a
20 retina surgeon on this panel, for example, I am
21 struck by the relative low rate of retinal
22 detachment post-operatively that was not present
23 preoperatively, the relative low rate of CME, and
24 if those methods of ascertainment are valid and
25 reliable, then I think that's great.

1 But I can't comment to what patients they
2 apply. For example, if you go to this Exhibit
3 Chart 1 on the Phase I core group, I see a listing
4 of the diagnoses on the far left side. But if you
5 look at the Y axis, it's defined as number of
6 patients/eyes. And that's confusing to me, because
7 some of the patients had more than one eye
8 implanted. So I need to identify what the study
9 group is a little bit better before I can make
10 comments.

11 DR. WEISS: We're going to have questions
12 by Dr. Coleman next and then we'll move on to the
13 FDA presentation.

14 DR. COLEMAN: Dr. Coleman. Did you want
15 to respond?

16 MR. WELCH: No, go ahead.

17 DR. COLEMAN: My question is do you have
18 those numbers of subjects or eyes that had
19 preexisting glaucoma prior to entering the study?
20 Since, as in the core group, about 39 of those 75
21 eyes had pseudoexfoliation, you would expect there
22 to be a high incidence of preexisting glaucoma.

23 And it's important because it appears that
24 the majority of the elevated intraocular pressures
25 after surgery were on those with preexisting

1 glaucoma. So it's nice to have also the
2 denominator of the eyes that started out with
3 preexisting glaucoma.

4 DR. HO: Jayne, can I comment?

5 DR. WEISS: Yes, briefly.

6 DR. HO: It's relevant to this because it
7 again identifies that you need to clarify who those
8 patients are. They may have cataract and
9 pseudoexfoliation. Some of the patients, in fact,
10 if you look at the table in the core group, are not
11 expected to have cataracts here, because the number
12 is about 75 percent. So I'm a little bit concerned
13 on commenting when I don't know exactly whose those
14 patients are.

15 MR. WELCH: Understand. Hid Welch. I'm
16 looking at Exhibit F-1b, which is I believe the
17 table you were looking at, the chart; is that
18 correct?

19 DR. HO: Right. You described--Allen Ho--
20 you describe it as the etiology table. So I'm
21 looking for trying to identify your study
22 population here.

23 MR. WELCH: Yeah.

24 DR. HO: It's not clear to me.

25 MR. WELCH: And F-1b is simply a--

1 DR. STEINERT: Are you talking about F-1b?

2 DR. HO: F-1b.

3 MR. WELCH: He's looking at this.

4 DR. STEINERT: The chart, F-1b, or are you
5 looking at the table F-1a?

6 DR. HO: 1a.

7 DR. STEINERT: Yeah, okay.

8 MR. WELCH: He's looking at 1a.

9 DR. STEINERT: And I'm sorry. I'm still
10 not exactly following the question. What is--

11 DR. COLEMAN: This is Dr. Coleman. My
12 question was is just in terms of the number of
13 preexisting cases of glaucoma.

14 MR. WELCH: Right.

15 DR. STEINERT: But what is--

16 DR. COLEMAN: Why? Because in terms of
17 the patients that had elevated intraocular
18 pressures that required treatment, that was one of
19 your points and variables that you were following
20 afterwards. The majority of those individuals were
21 said to have preexisting glaucoma. And so it's
22 also nice to know how many preexisting glaucoma
23 individuals in the study didn't have elevated
24 intraocular pressures afterwards.

25 DR. STEINERT: In other words, did the

1 ring somehow reduce their intraocular pressure? Is
2 that the question?

3 DR. COLEMAN: No, it just gives you kind
4 of a denominator whether you know exactly what's
5 going on because if only 20 individuals in the
6 study had preexisting glaucoma prior to
7 implantation of the intraocular lens and the ring,
8 and 20 had problems with intraocular pressure after
9 the surgery, that's a little higher than you
10 usually see. And so it gives you some idea of
11 exactly what the population is.

12 It might be that 30 people or 20 people
13 had preexisting glaucoma, and you only had trouble
14 with ten of them with intraocular pressure
15 afterwards which would be expected in a population
16 like this. That's the main thing because it just
17 gives you a denominator to work with.

18 DR. STEINERT: The purpose--to be sure
19 that's an interesting question. I think, you know,
20 from our point of view, the point of that table was
21 simply to address the question as to whether there
22 was a safety issue and whether there was any
23 indication that these procedures and the use of
24 this ring in particular caused an undue or alarming
25 or concerning rate of elevated intraocular

1 pressure, and it appeared that the answer was no.
2 And that's where the analysis stopped.

3 DR. WEISS: I'd like to thank the sponsor
4 for their presentation, and we're going to have you
5 move back from the table and have the FDA come up
6 and give their presentation. Please.

7 **FDA PRESENTATION**

8 MS. LOCHNER: This is Donna Lochner. I'm
9 going to give some introductory comments to the
10 PMA. The PMA for the Morcher Capsular Tension Ring
11 was received by FDA on October 16, 2001, and was
12 accepted into the Office of Device Evaluation's
13 expedited review program. Expedited review is
14 granted for first-of-a-kind devices for which no
15 approved alternative treatment devices exists, and
16 in the case of the capsular tension ring, which may
17 potentially reduce the risk of morbidity for the
18 indicated patient population.

19 Expedited review is intended to move
20 applications to the front of FDA's review queue,
21 but does not waive clinical or scientific safety
22 and effectiveness endpoints.

23 Rather, consideration of the difference in
24 the risk-to-benefit analysis because of the lack of
25 alternatives is addressed in the design of the

1 clinical study protocol.

2 The sponsor chose to participate in OED's
3 modular PMA program which allows for review of
4 sections, or module, of the PMA as they are
5 completed. When all information that is required
6 to be included in a PMA application has been
7 submitted, the PMA may be filed and review
8 proceeds.

9 The final clinical section of this PMA was
10 submitted on October 16 and so this is the filing
11 date of the Morcher PMA.

12 Three modules preceded the submission of
13 the clinical data and contain the sterilization
14 procedures and validations, the manufacturing and
15 engineering procedures and validations, and the
16 biocompatibility data.

17 There are outstanding issues in each of
18 these scientific areas. At this time, FDA is
19 awaiting adequate responses from the sponsor. In
20 addition, completion of the bioresearch monitoring
21 inspections and scheduling and completion of the
22 good manufacturing practices inspections are also
23 outstanding.

24 Today, we are asking the panel to review
25 and make recommendations on the clinical data

1 contained in the PMA. However, the status of the
2 PMA is provided to make sure that you understand
3 the FDA awaits adequate responses to the remaining
4 scientific sections of the PMA prior to any final
5 decision on the application.

6 We appreciate the efforts of the panel and
7 particularly the primary panel reviewers, Drs.
8 Sugar and Van Meter, who reviewed the document on a
9 compressed schedule to the expedited status of the
10 PMA.

11 In balance, we felt it most efficient to
12 proceed with the panel meeting in consideration of
13 the potential benefit of this device to public
14 health.

15 I'd also like to acknowledge the
16 exceptional efforts of the FDA review team, and
17 particularly Joel Glover, the engineering and lead
18 reviewer, and Dr. Bernard Lepri, the clinical
19 reviewer.

20 I'd also like to acknowledge Dr. Kesia
21 Alexander who was the lead reviewer for most of the
22 IDE. All three of these individuals have made
23 significant efforts in consulting with the sponsor
24 over the years.

25 Now, I'd like to introduce the lead

1 reviewer, Joel Glover, who will present an overview
2 of the scientific non-clinical sections of the PMA.

3 MR. GLOVER: My name is Joel Glover. As
4 Donna mentioned, I'm the team leader for the
5 application, and as she also mentioned I'd like to
6 acknowledge Kesia Alexander who was the team leader
7 before I became involved with the application.

8 I'm going to present a brief history of
9 our experience with the PMA. The capsular tension
10 ring PMA was done under our Modular PMA Program.
11 It was actually initiated with a shell outline of
12 the intended modules back in August of 1998, and
13 followed shortly thereafter with the submission of
14 some of the preclinical modules.

15 The clinical module, if you will, that
16 actually triggered the PMA was received in October
17 of last year. The modular PMA was composed of four
18 modules. Module 1 was general information. Module
19 2, biocompatibility. Module 3, the
20 microbiology/sterilization module. And Module 4
21 contained the manufacturing. And again, the PMA
22 was to be the clinical data.

23 Module 1, general information, contained
24 general information about the device, the
25 applicant, manufacturing sites and FDA considers

1 that complete.

2 Module 2 was to address biocompatibility
3 of the device. The device is made of
4 polymethacrylate. The sponsor provided
5 cytotoxicity test data and residual levels.

6 FDA has some outstanding issues with this
7 module, in particular the nature of the specific
8 PMMA material that the sponsor is using to
9 construct their device, and some issues with
10 residual levels and identifying what those
11 residuals are.

12 The sponsor only performed essentially
13 cytotoxicity testing. So FDA has concern about the
14 lack of biocompatibility testing or having a
15 justification for omitting the testing. The issue
16 has been discussed many times with the sponsor and
17 FDA and the outcome of that is essentially that
18 the sponsor has chosen to use the clinical data
19 from the PMA to demonstrate the biocompatibility of
20 their device as opposed to conducting further tests
21 to support biocompatibility.

22 And this is an agreement that--I shouldn't
23 say agreement--but this is an argument that FDA is
24 willing to consider for this device.

25 Module 3 was the microbiology/

1 sterilization. It contains such things as how the
2 device is sterilized, validations for the
3 sterilization procedure, and a study of the shelf
4 life of the device. There are some outstanding
5 issues with this module as well.

6 Module 4 is the manufacturing section.
7 Both the Office of Compliance which looks at the
8 good manufacturing practices and quality control
9 procedures and reviews that and conducts the GMP
10 inspection has outstanding issues, as does the
11 Office of Device Evaluation, and these will need to
12 be addressed before an ultimate approval of the
13 PMA.

14 I would point out that FDA doesn't believe
15 that these outstanding preclinical issues warrant
16 delaying the panel's review of the clinical data
17 for the PMA, and that's why we brought it forward
18 to you at this early stage.

19 Finally, the PMA was submitted last year
20 and contained the results of the clinical study of
21 the device, and in a moment, Dr. Lepri will present
22 his analysis of the clinical data. And finally,
23 I'd just like to thank the panel for their review
24 and deliberations today and also the members of
25 FDA's review team for their reviews and quick

1 responses, as well as the sponsor's responses with
2 regard to this expedited PMA.

3 And if there are no questions, I'll
4 introduce Dr. Bernard Lepri, the clinical reviewer.

5 DR. WEISS: Dr. Bradley has a question.

6 DR. BRADLEY: I just wondered if there are
7 any reasons to question the use of using the
8 clinical data to ascertain biocompatibility? You
9 said there was some argument about it. I got the
10 impression that the FDA has some reservations about
11 the validity of that approach. Could you comment
12 on that?

13 MS. LOCHNER: Well, I think first of all
14 typically biocompatibility testing is done
15 preclinically to screen for potential problems
16 before a material is implanted in the eye.

17 While we didn't de facto accept the
18 sponsor's argument that this is
19 polymethacrylate and so should be allowed in
20 the eye, because we believe the sponsor had to
21 identify what type of polymethacrylate was
22 used, we felt that their worldwide experience to
23 date warranted initiation of the IDE study.

24 It was our understanding at the time that
25 they would then proceed to collect the usual

1 biocompatibility testing. When you go through the
2 usual battery of tests and look at what would the
3 clinical study not address, you come down to two
4 areas. One is the ocular implant test in rabbits.

5 Certainly the clinical data would suffice
6 as a replacement for that, but what the ocular
7 implant test does do is look at histopathology that
8 isn't provided in the clinical study. But I think
9 a reasonable argument could be made by the sponsor
10 that their outcomes wouldn't suggest problems. And
11 so not having that histopathology from the rabbit
12 study I think they could make a valid argument for
13 that.

14 The second area that the particular
15 clinical data does not address is the
16 carcinogenicity testing. And again, we believe the
17 sponsor can make an argument that even though this
18 is a PMMA that has not been used in the U.S.
19 previously, we think they could make a valid
20 argument that there is no reason to expect that
21 this PMMA would be a carcinogen.

22 So this is where we are today. However,
23 we felt it important that the panel understand that
24 unlike your usual review of clinical data, an
25 additional component is that these clinical data

1 are being used to support the biocompatibility.

2 So if there was any question that you
3 might have, not even--you know, any issue with the
4 outcomes that in the past you would have said, oh,
5 this is no problem, at least you're aware that the
6 usual biocompatibility testing was done so that
7 when you say you have no concerns with the
8 outcomes, you know, you're aware that the usual
9 biocompatibility wasn't done.

10 So we're providing all this background so
11 you understand, you know, what is atypical about
12 this document.

13 DR. WEISS: Thank you. We're going to go
14 on to Dr. Lepri's presentation.

15 DR. LEPRI: Thank you. Good morning,
16 members of the panel, sponsors, FDA, staff members,
17 and guests.

18 I'd like to make some introductory
19 comments, but I have taken the liberty of compiling
20 data and trying to present an overall picture of
21 what was given to us in this PMA, and I apologize
22 for any inaccuracies I may have, but they are
23 limited by the numbers that I was presented with,
24 and as one mentioned, we had a roaming end, but we
25 slowed down the speed of that roaming end with the

1 first deficiency letter.

2 At this time, I'd also like to thank my
3 fellow FDA members for their encouragement and
4 support: Dr. Rosenthal; Donna Lochner, the early
5 days with Dr. Kesia Alexander, and especially my
6 own personal "Lord of the Rings" hero, Joel Glover,
7 who has helped keep this entire application and
8 process moving smoothly and accurately and helping
9 me to meet these very compressed time schedules.

10 Okay. The PMA application presents
11 varying forms of the indications statement, and I'm
12 going to present those to you in the next few
13 slides for your consideration when we later on in
14 this process ask you for your labeling
15 recommendations.

16 The initial indication in the beginning of
17 the PMA and the IDE read that it is used for the
18 stabilizing the capsular bag in cataract surgery
19 with IOL implantation, in cases of
20 pseudoexfoliation syndrome, where there is
21 subluxation of the lens or zonular damage as in
22 Marfan's syndrome and in traumatic cases, and cases
23 where pars plana vitrectomy has been performed.

24 Next. The next indication statement is in
25 Exhibit 1, and all the items you see listed there

1 have been added to the indications statement in
2 this proposed labeling. I want you to particularly
3 note stabilizing the capsular bag in high myopia,
4 stabilizing operating conditions, implantation of
5 foldable IOLs, circular expansion of the capsular
6 bag, and prevention of unilateral shrinkage of
7 capsular bag, and prevention of capsular fibrosis.

8 The next is in Exhibit K where once again
9 high myopia is missing, is added, and some of the
10 same items are repeated as in Exhibit I, but
11 compressed.

12 The CTR, the capsular tension ring, is a
13 flexible, one-piece ring of PMMA that ranges from
14 ten to 12 millimeters in diameter. Utilized in
15 this trial were three types: the 14, 14A, and 14C
16 type rings, which differ in dimensions to
17 accommodate the differences in capsular bag sizes
18 of individual eyes.

19 The study was comprised in this PMA of two
20 phases: the Phase I core and Phase II which
21 included core investigators as well as independent
22 investigators.

23 And Phase II was conducted primarily to
24 provide confirmatory data. This was a prospective,
25 open label, multi-site/multi-investigator trial.

1 Demographically, the core group was
2 comprised of 27 males and 48 females. The age
3 stratification was that 24 of those 75 were between
4 the ages of 70 to 79, and 26 of them were either 80
5 years of age or older.

6 Phase II combined core and independent
7 investigators was reported as 238 males and 237
8 females. You may note that those numbers do not
9 add up to the total of 415 patients as noted
10 elsewhere in the PMA. That's another one of those
11 discrepancies.

12 The demographics of the preoperative
13 pathology show that there were some significant
14 pathologies that were majorly represented in this
15 investigation. There were a combined total of 161
16 pseudoexfoliation patients, followed most
17 frequently by trauma cases, and as you can see,
18 there are cases, 12 cases of Marfan's and 22 cases
19 of vitrectomy, that were at the time of surgery
20 when they implanted the ring.

21 Of course, the most widely represented was
22 the pseudoexfoliation. The data results for
23 effectiveness and/or safety were not stratified by
24 preoperative pathology other than for capsular
25 fibrosis and contraction and IOL decentration.

1 The only data presented on the 98 trauma
2 cases was visual acuity. I also wanted to make
3 note that we have not had any data presented to FDA
4 in the early post-operative periods where we would
5 get a very good indication in the response of these
6 patients to the implantation of the ring, and it
7 has been noted in your deliberations as well as the
8 sponsor's presentation on very low rates of iritis,
9 that those were calculated for the overall period,
10 and out to one year or two years in the case of the
11 core.

12 Accountability. The accountability, when
13 calculated by FDA's criteria, was higher than that
14 obtained by the sponsor and overall reasonably
15 good. The Office of Device Evaluation recommends a
16 minimum accountability at the time of submission of
17 a PMA to be at least 80 percent. You can see from
18 this chart that Phase I was at 88 percent at one
19 year with three lost to follow-up. And Phase II
20 was at about 81 percent at one year with
21 approximately 63 lost to follow-up.

22 The endpoints established for this
23 investigation were IOL centration both pre and post
24 YAG as a major effectiveness criterion, and the
25 safety criteria were the FDA IOL grid. These were

1 safety variables. There were no standardized
2 criteria.

3 Let me step back. There were no
4 standardized criteria, as mentioned by Dr.
5 Steinert, for measuring IOL centration and there
6 was no establishment of a criterion of what would
7 be considered significant or expected.

8 The FDA IOL grid was only used as a guide
9 for the sponsor to use in evaluating complication
10 rates of implantation of this device. In no way
11 was it intended for the sponsor to have to meet the
12 criteria established in the FDA IOL grid since this
13 is not an IOL.

14 IOL centration at one year post-op. This
15 slide presents the number of eyes in each phase as
16 identified by the individual investigators. Phase
17 I core had the highest at ten percent. Auffarth,
18 et al., in 1994, conducted post-mortem studies of
19 eyes with PXE patients and noted a higher incidence
20 of decentration in bag fixated IOLs, this resulting
21 from intraoperative zonulysis.

22 Intraoperative zonulysis ranges from 13.1
23 percent to 17.9 percent according to the
24 literature. So we can see the IOL centration
25 measured with the methods used and not firmly

1 established were well within those ranges
2 postoperatively.

3 We further analyzed these data to
4 establish how the IOL centration was rate wise even
5 though the end values were small for those
6 individuals who had YAG capsulotomies performed,
7 and we can see that a rate analysis of IOL
8 decentration greater than or equal to one
9 millimeter post-YAG produces results that are
10 comparable to those reported in the literature for
11 the amount of the zonulysis and decentration.

12 Phase I core had 12 YAGs by the last set
13 of data that I have received, and two of them
14 reported that they had decentration of greater than
15 or equal to one millimeter. And in Phase II, there
16 were seven YAGs performed. One of them reported as
17 having greater than or equal to one millimeter of
18 decentration at a rate of 14.29 percent.

19 PXE patients who were the bulk of the
20 patients in this investigation often exhibit postop
21 IOL decentration due to intraoperative zonulysis.
22 As I mentioned before, the literature rates of
23 zonulysis in PXE range from 13.1 to 17.9, and the
24 CTR post-YAG decentration rates range from 14.29 to
25 16.0.

1 Also, very common in this, these types of
2 cases, are capsular fibrosis and contraction, which
3 were issues that were mentioned in the proposed
4 labeling. Since capsular contraction results from
5 the fibrosis of the capsule, I took the liberty of
6 combining these data and presenting them all on one
7 chart.

8 The total fibrosis for all phases reported
9 was 9.5 percent and the total amount of contraction
10 reported was 3.2 percent. And the last column,
11 that's correct. There was an error last night when
12 I was preparing for this.

13 YAG rates. The rate of PCO calculated for
14 core group eyes evaluated at one year is 28
15 percent. There were 14 out of 50 eyes. I believe
16 it was Exhibit H-1 and H-2 that listed the
17 complications. And at the bottom of that list, it
18 claims that the percentages were based on the
19 number of eyes examined. So for the core group
20 that was 50.

21 And that's where the 28 percent is from.
22 I transferred the rates of the percentages of YAGs
23 and the percent of fibrosis and capsular
24 contraction on to this chart also. That way you
25 can compare them coming to this review.

1 39 of the 75 subjects enrolled in Phase I
2 had PXE, and it's reported in the literature by
3 Naumann and many others that PXE patients have
4 higher rates of PCO postoperatively.

5 The next slide, please. This slide
6 presents the numbers of postop IOP increases in
7 eyes at one year, 11 to 13 months, who did not have
8 preexisting glaucoma. These were reported by the
9 sponsor in Amendment No. 3 of this PMA. They were
10 not classified as adverse events by the sponsor.
11 The sponsor claimed that they were not classified
12 as adverse events because the patient's other
13 conditions were more serious in the investigators'
14 opinions.

15 The FDA considers all post-op IOP
16 increases as adverse events whether they are device
17 related or not, and you could see the rates that I
18 calculated based on the numbers and the level of
19 accuracy that was presented to me that there were
20 two in the core group, which gives you a rate of
21 four percent, and 14 in the independent group, and
22 at that time point at one year, it's reported in
23 Amendment 3 that 297 patients were evaluated, and
24 that's the denominator that I used to obtain a rate
25 of 4.7 percent.

1 Visual acuity. While the endocapsular
2 tension ring is not directly responsible for visual
3 acuity outcomes, their analysis is valuable in
4 representing the benefits to the subjects of
5 cataract surgery within the scope of this
6 investigation.

7 One can see that there were significant
8 numbers of eyes with better than 20/40 BSCVA
9 preoperative. 74 of the 75 core group eyes were
10 reported as having cataracts. The sponsor reports
11 BSCVA of greater than or equal to 20/40 post-op in
12 the core group at a rate of 87.87 percent, which is
13 close to the target value of the FDA IOL grid of
14 92.5 percent.

15 Phase II subjects did not fare as well
16 postoperatively with respect to best corrected
17 visual acuity. The sponsor did not provide
18 sufficient detail for the Phase II results such as
19 best case analysis or results stratified by
20 preoperative pathologies to document the cause of
21 the lower than average acuity outcomes
22 postoperatives.

23 The sponsor reported that there were 12
24 eyes of 11 subjects in Phase II with macular
25 degeneration, but this in no way accounts for the

1 percentages of BSCVA reported.

2 Also, the sponsor did not calculate the
3 rates in this table. I calculated these rates.
4 They did not present this data.

5 Explants and secondary interventions. I
6 think this is a repetition of Dr. Steinert's slide.
7 There was one secondary reintervention due to
8 capsule/IOL problems a one week.

9 There were five explants in Phase II
10 subjects' eyes; three were in the core group and
11 two in the independent group.

12 And there were seven others performed
13 during initial surgery, two of which were due to
14 procedural complications, four due to inadequate
15 capsular/zonular support, and one for an incorrect
16 ring size.

17 Inflammatory complications. The most
18 noted in these populations would be iritis,
19 synechiae, IOL lens deposits and CME.

20 The inflammatory complications in this
21 report were these and many others which were not of
22 significant numbers to mention at this time.

23 Iritis. The sponsor's presentation
24 reported six cases of iritis that are not presented
25 in Amendment No. 3. At the time of their

1 occurrence is not noted, and the PMA did not report
2 any postoperative data earlier than one year. So
3 we have no information on the critical
4 postoperative, immediate postoperative, time
5 periods.

6 Pseudoexfoliation patients according to
7 the literature exhibited an impaired blood-aqueous
8 barrier which would yield higher rates of iritis
9 postoperatively which we have not seen any of the
10 data presented in the PMA.

11 They also have increased fibrinoid
12 reactions which lead to potential posterior
13 synechiae and IOL cell deposits.

14 There were some synechiae reported in
15 Amendment No. 3 in Exhibits H-1 and H-2, and there
16 were--although the rates were low--okay--and one
17 would expect as well as hope to see these lower
18 rates at one year post-op. As I mentioned before,
19 we didn't see anything early on. The rates were
20 essentially one percent in Phase II and
21 approximately two percent in Phase I, and it was in
22 the nature of anterior synechiae.

23 Cystoid macular edema. The sponsor's
24 presentation reported--that was forwarded to FDA
25 prior to today--reported 11 cases of CME out of 524

1 patients, which I later learned that those were 524
2 implants, not patients, and these 11 cases are not
3 explained with reference to the time point of
4 occurrence, and the rate is not calculated using a
5 denominator of the eyes examined, as it should be,
6 but rather they used a denominator of the total
7 number enrolled and treated. And so we might get
8 more valuable information for purposes of labeling
9 in the performance of this device if we knew this
10 occurred and it was based on the number of patients
11 actually evaluated.

12 When we compiled--we--I compiled the one
13 year data that was presented in Amendment No. 3, I
14 found two percent rates of CME at one year for
15 Phase I core subjects and all Phase II combined
16 subjects.

17 Now I will present the questions. Some of
18 the questions I will make some reference to some of
19 the information found in our literature review,
20 just as a matter of background, and I fully
21 acknowledge your expertise and that you may already
22 know this.

23 Question No. 1: The sponsor has not
24 performed the standard battery of biocompatibility
25 testing on the device, and has proposed to use the

1 clinical data to document the biocompatibility of
2 the device. Do the adverse events and their rates
3 reported in the PMA support raise any safety
4 concerns from your clinical perspective?

5 Question No. 2: Patients with high myopia
6 were not included in the U.S. clinical study. Do
7 the data in the PMA support these proposed
8 indications for use?

9 Question No. 3: Do the clinical data
10 presented in the PMA provide sufficient evidence
11 and effectiveness of the device for the proposed
12 indications for use, taking into account the
13 revisions in response to question number two, if
14 any?

15 Question No. 4: Do you have any
16 recommendations for revisions or additions to the
17 labeling as proposed by the sponsor? Please
18 consider the following issues in your
19 deliberations:

20 Part a, high myopia, lens extraction
21 without IOL implementation;

22 Part b, progressiveness of syndromes such
23 as pseudoexfoliation and Marfan's.

24 And part c, late onset of dislocation of
25 capsular bag containing IOL and ring in

1 pseudoexfoliation syndrome.

2 And I will note that in the literature
3 review, we found Jehan, et al., at 2001, and he
4 presented the results of an eight eye/seven patient
5 study, and these were patients who had previously
6 undergone uncomplicated cataract surgery with IOL
7 implantation. All of them, 100 percent of them,
8 experienced delayed dislocation into posterior
9 chamber.

10 And the mean time for dislocation was
11 seven years. And there was one other literature
12 report that reported this occurrence as late as 12
13 years post-op.

14 And part d, the use of Type 14 rings in
15 pediatric patients, size issues and potential
16 radial tears in capsular bag. And the origin of
17 this concerned is an article published by Dietlien,
18 et al, in the year 2000, of complications in a
19 four-year-old who experienced upward displacement
20 of the bag, capsular bag, after the ring was
21 implanted interoperatively.

22 They claim in this article that the adult
23 rings were not a good choice for pediatric patients
24 for two main reasons: the proliferation of lens
25 epithelial in a growing eye and the weak zonules.

1 And these combinations led to CTR dislocation and
2 distortion. And the size of the rings may be too
3 large for some pediatric patients, particularly
4 Marfan patients or trauma victims, and has the
5 potential to cause radial tears of the capsular
6 rexis.

7 Thank you.

8 DR. WEISS: Thank you, Dr. Lepri. Now
9 we're going to open the floor for any questions
10 from the panel to Dr. Lepri or the agency.

11 **Panel Questions for FDA**

12 DR. WEISS: Dr. Van Meter.

13 DR. VAN METER: In your last page of
14 questions, you did not address the efficacy--

15 MS. THORNTON: Dr. Van Meter, could you
16 speak into the microphone, please?

17 DR. VAN METER: Yes. Woodford Van Meter.
18 On your questions on 4, part a, b and c, you did
19 not get into the demonstrated efficacy of reducing
20 capsular fibrosis or capsular contraction. Is that
21 still a concern of yours?

22 DR. LEPRI: That's still a concern. That
23 was presented; it was part of the presentation, and
24 I presented that data. We just had these
25 additional concerns that were obtained from the

1 literature and that's why they were mentioned
2 separately at the end, and I didn't want to bore
3 you to death with providing you exhaustive
4 literature.

5 DR. WEISS: Dr. Bradley.

6 DR. BRADLEY: Two questions. You gave us
7 summary data on the best corrected visual acuity
8 post-op, and the summary statistic is basically 40
9 percent end with visual acuities worse than 20/40.
10 That was my read on the table, and I just wondered
11 whether that was anticipated, and what were the
12 root cause of these poor acuities in such a large
13 percentage of these patients?

14 DR. LEPRI: I had addressed that question
15 to the sponsor in one of our deficiency letters
16 that were issued, and the explanation given to me
17 was that many of these patients had severe
18 preoperatively pathologies which would lend then to
19 not have good post-operative visual acuity
20 outcomes.

21 My contention with that is that there were
22 a large number of patients preoperatively,
23 particularly in the core group, who had BSCVAs that
24 were better than 20/40, and that if they had
25 provided a best case and worst case analysis of