

elw

JH

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
OPHTHALMIC DEVICES PANEL

*This transcript has not
been edited and FDA
makes no representation
regarding its accuracy*

OPEN SESSION

PMA P870024/S043

Friday, January 18, 2002

8:34 a.m.

Hilton Washington DC North/Gaithersburg
Salons A, B and C
620 Perry Parkway
Gaithersburg, MD

MILLER REPORTING COMPANY, INC.
735 8th Street, S.E.
Washington, D.C. 20003-2802
(202) 546-6666

P A R T I C I P A N T S

PANEL PARTICIPANTS:

Jayne S. Weiss, M.D., Acting Chair
Arthur Bradley, Ph.D., Voting Member
Michael R. Grimmett, M.D., Voting Member
Alice Y. Matoba, M.D., Voting Member

Timothy T. McMahon, O.D., Consultant, deputized
to vote

Allen C. Ho, M.D., Consultant, deputized to
vote

Anne L. Coleman, M.D., Ph.D., Consultant,
deputized to vote

Joel Sugar, M.D., Consultant, deputized to vote

Richard Casey, M.D., Consultant, deputized to
vote

Janine A. Smith, M.D., Consultant, deputized to
vote

Woodford S. Van Meter, M.D., Consultant,
deputized to vote

Glenda V. Such, M.Ed, Consumer Representative

FDA:

Sara M. Thornton, Panel Executive Secretary
A. Ralph Rosenthal, M.D. Director, Division of
of Ophthalmic Devices

Everette T. Beers, Ph.D.

James F. Saviola, O.D.

Donna R. Lochner

Bernard P. Lepri, O.D., M.S. M.Ed.

Joel P. Glover, M.S.

Eleanor McGhee

C O N T E N T S

Call to Order	
Dr. Jayne S. Weiss Acting Chair	4
Introductory Remarks	
Sara M. Thornton, Executive Secretary	4
OPEN PUBLIC HEARING	
Dr. Marjorie Rah, Assistant Professor, New England College of Optometry	11
OPEN COMMITTEE SESSION (Dr. Jayne S. Weiss, Acting Chair)	
PMA P870024/S043	
Sponsor Presentation	
Dr. Bill Meyers	18
Dr. Mark Bullimore	20
Dr. Jerry Legerton	26
Dr. Oliver Schein	35
Dr. Michael DePaolis	40
Dr. Jerry Legerton	48
Panel Questions for the Sponsor	51
FDA Presentation	
James Saviola, O.D., Introduction	90
Eleanor McGhee, Team Leader	95
Bernard L. Lepri, O.D., Clinical Reviewer	95
COMMITTEE DELIBERATIONS	
Primary Panel Reviewers:	
Timothy T. McMahon, O.D.	99
Timothy B. Edrington, O.D.	110
PANEL DISCUSSION OF P870024\S043	113
OPEN PUBLIC HEARING SESSION	195
FDA CLOSING COMMENTS	195
SPONSER - CLOSING COMMENTS	197
VOTING OPTIONS READ	197
PANEL RECOMMENDATIONS TAKEN BY VOTE	199
POLLING OF PANEL VOTES	233
COMMENTS FROM CONSUMER AND INDUSTRY REPRESENTATIVES	238
ADJOURNMENT	239

1 P R O C E E D I N G S

2 **Call to Order**

3 DR. WEISS: I would like to call this
4 meeting of the Ophthalmic Devices Panel to order,
5 and we will have introductory remarks from Sara
6 Thornton.

7 **Introductory Remarks**

8 MS. THORNTON: Good morning, and welcome
9 to the 103rd meeting of the Ophthalmic Devices
10 Panel. Before we proceed with today's agenda, I
11 have a few short announcements.

12 I'd like to remind everyone to sign in on
13 the sign-out sheets that are out on the
14 registration table. That's just outside the
15 meeting room. You probably walked by it. But
16 please sign in for us, because we appreciate
17 knowing who has attended.

18 All handouts for today's meeting are
19 available on the registration table for you.
20 Messages for panel members and FDA participants,
21 information or special needs, should be directed
22 through Ms. Ann Marie Williams, who is available in
23 the registration area. The phone number for calls
24 to the meeting area is (301) 977-8900. That phone
25 is at the registration desk.

1 In consideration of the panel, the sponsor
2 and the agency, we ask that those of you who have
3 cell phones and pagers either turn them off,
4 please, or put them on vibration mode for the
5 duration of the meeting while you're in this room.
6 And we ask that all meeting participants speak
7 clearly into the microphone, give your name, so
8 that the transcriber will have an accurate
9 recording of your comments. We had some problems
10 yesterday, and I would just like to remind everyone
11 to speak close and directly into the microphone.

12 All available information for the meeting
13 that we have tentatively scheduled for March 14th
14 and 15th will be on the FDA Advisory Committee web
15 site in approximately one week.

16 Now at this time I would like to extend a
17 special welcome, and introduce again to the public
18 and the panel and the FDA staff, two panel
19 consultants who are with us for the first time at
20 this meeting, and we have a new panel consumer
21 representative as well.

22 Dr. Richard Casey, to my right, comes to
23 us from Los Angeles, where he is an Associate
24 Professor of Ophthalmology at the Jules Stein Eye
25 Institute and the Interim Chairman of the

1 Department of Ophthalmology at the Charles Drew
2 University of Medicine and Science. His clinical
3 practice involves the management of corneal and
4 anterior segment disease, cataract and refractive
5 surgery.

6 Dr. Janine Smith, to my left, is the
7 Deputy Clinical Director at the National Eye
8 Institute of the National Institutes of Health in
9 Bethesda, Maryland. Her basic science research has
10 been immune-based diseases of the ocular surface,
11 with additional responsibilities for the NEI
12 Intramural Clinical Research Program.

13 And Ms. Glenda Such, again to my left, the
14 consumer representative to the panel, is the
15 Director of the Computer Training Programs in the
16 Department of Career Services at the Lighthouse
17 International in New York. She is a recognized
18 expert in the field of adaptive technology for
19 visual impairments and the functional implications
20 of visual disabilities, particularly low vision.

21 We very much appreciate your commitment to
22 serve, and we welcome you to the panel table today.

23 To continue, will the remaining panel
24 members please introduce themselves, beginning with
25 Dr. Harris?

elw

7

1 DR. HARRIS: Oldest first. Thank you.

2 [Laughter.]

3 MS. THORNTON: I didn't say that.

4 DR. HARRIS: Sally said that I'm the one
5 with historical perspective because I was on the
6 panel when Franklin came here with bifocals.

7 [Laughter.]

8 DR. HARRIS: Not true. I'm Michael
9 Harris. I'm the Associate Dean and Clinical
10 Professor at the University of California School of
11 Optometry, where I am also chief of the Contact
12 Lens Clinic. I'm also an attorney at law and
13 member of the California State Bar.

14 DR. EDRINGTON: Tim Edrington, Professor
15 at Southern California College of Optometry.

16 DR. McMAHON: Tim McMahon, Professor and
17 Director of the Contact Lens Service in the
18 Department of Ophthalmology at the University of
19 Illinois/Chicago.

20 DR. MATOBA: Alice Matoba, Associate
21 Professor of Ophthalmology, Baylor College of
22 Medicine.

23 DR. BRADLEY: Arthur Bradley, Professor of
24 Vision Science, Indiana University.

25 DR. WEISS: Jayne Weiss, Professor of

1 Ophthalmology and Pathology, Kresge Eye Institute,
2 Wayne State University, Detroit.

3 DR. GRIMMETT: Michael Grimmatt, Assistant
4 Professor, University of Miami School of Medicine,
5 and Medical Director of Bascom Palmer Eye Institute
6 of the Palm Beaches.

7 DR. COLEMAN: Anne Coleman, Associate
8 Professor of Ophthalmology, UCLA/Los Angeles.

9 DR. HO: Good morning. Allen Ho,
10 Associate Professor of Ophthalmology, Thomas
11 Jefferson University, Wills Eye Hospital.

12 DR. VAN METER: Woodford Van Meter,
13 Associate Professor of Ophthalmology, University of
14 Kentucky in Lexington, Kentucky.

15 MR. MCCARLEY: Rick McCarley, President of
16 the Voc Tech in Boca Raton, Florida. I'm the
17 industry representative.

18 DR. ROSENTHAL: Ralph Rosenthal, Division
19 Director of Ophthalmology and ENT, Food and Drug
20 Administration.

21 MS. THORNTON: I'd like to read now the
22 conflict of interest statement for today's meeting.

23 The following announcement addresses
24 conflict of interest issues associated with this
25 meeting, and is made part of the record to preclude

1 even the appearance of an impropriety. To
2 determine if any conflict existed, the agency
3 reviewed the submitted agenda for this meeting and
4 all financial interests reported by the committee
5 participants.

6 The conflict of interest statutes prohibit
7 special government employees from participating in
8 matters that could affect their or their employer's
9 financial interests. However, the agency has
10 determined that participation of certain members
11 and consultants, the need for whose services
12 outweighs the potential conflict of interest
13 involved, is in the best interests of the
14 government.

15 Therefore, a waiver has been granted for
16 Dr. Michael Harris for his interests in a firm that
17 could potentially be affected by the panel's
18 recommendations. The waiver allows this individual
19 to participate fully in today's deliberations.
20 Copies of this waiver may be obtained from the
21 agency's freedom of information office, Room 12A-15
22 of the Parklawn Building.

23 We would like to note for the record that
24 the agency took into consideration matters
25 regarding Drs. Arthur Bradley, Timothy Edrington,

1 Michael Harris, Allen Ho, and Timothy McMahon, who
2 reported interest in firms at issue but in matters
3 that are not related to today's agenda. The agency
4 has determined, therefore, that they may
5 participate fully in all discussions.

6 In the event that the discussions involve
7 any other products or firms not already on the
8 agenda, for which an FDA participant has a
9 financial interest, the participant should excuse
10 him or herself from such involvement, and the
11 exclusion will be noted for the record. With
12 respect to all other participants, we ask in the
13 interest of fairness that all persons making
14 statements or presentations disclose any current or
15 previous financial involvement with any firm whose
16 products they may wish to comment upon.

17 Thank you, Dr. Weiss.

18 DR. WEISS: Thank you, Sally. We now will
19 begin the open--

20 MS. THORNTON: I had one more thing.
21 Sorry. I do need to do this before we proceed, the
22 appointment to temporary voting status for those
23 who are here at the table.

24 "Pursuant to the authority granted under
25 the Medical Devices Advisory Committee charter

1 dated October 27, 1990, and as amended August 18,
2 1999, I appoint the following individuals as voting
3 members of the Ophthalmic Devices Panel for this
4 meeting on January 18, 2002: Drs. Allen Ho,
5 Timothy McMahon, Anne Coleman, Richard Casey,
6 Janine Smith, Woodford Van Meter, Timothy
7 Edrington, and Michael Harris. In addition, I
8 appoint Dr. Jayne S. Weiss to serve as Acting Panel
9 Chair for the duration of this meeting.

10 "For the record, these individuals are
11 special government employees and consultants to
12 this panel or other panels under the Medical
13 Devices Advisory Committee. They have undergone
14 the customary conflict of interest review and have
15 reviewed the material to be considered at this
16 meeting." Signed, Dr. David Feigal, Director of
17 the Center for Devices and Radiological Health,
18 January 9, 2002.

19 Thank you, Dr. Weiss.

20 **OPEN PUBLIC HEARING**

21 DR. WEISS: We will now begin the open
22 public hearing for up to a half hour. I would ask
23 that anyone who has any comments to make will
24 identify themselves and any financial or other
25 potential conflicts that they may have.

1 We have been informed that Dr. Marjorie
2 Rah has a presentation to make during this portion.
3 Dr. Rah?

4 DR. RAH: Dr. Marjorie Rah. I'm an
5 Assistant Professor at the New England College of
6 Optometry. What I'd like to present this morning
7 is results, preliminary results from an independent
8 study that I have conducted there. I do not have
9 any personal financial interest in the sponsor and
10 I am not an investigator on their PMA. I'd like to
11 thank the panel for giving me the opportunity to
12 speak this morning, and I would like to acknowledge
13 that I do have an unrestricted grant from Paragon
14 Vision Science that was provided to the Ohio State
15 University to in part fund the project.

16 The purpose of the project, "The Lenses
17 and Overnight Orthokeratology Study," was it's a
18 pilot study designed to evaluate the safety and
19 efficacy of overnight orthokeratology. It's a
20 multi-center pilot study. Data collection was
21 conducted and is ongoing at the Ohio State
22 University College of Optometry, the New England
23 College of Optometry, Southern California College
24 of Optometry, and we have one investigator who has
25 relocated to the Southern College of Optometry in

1 Memphis, but there is no data collection being
2 conducted there. The study is being conducted
3 under the FDA IDE G990205.

4 Sixty patients have been fitted in the
5 study, 25 with the Fargo 6 orthokeratology design
6 and HDS material. The remaining 35 were fitted
7 with the Paragon CRT design and HDS material.
8 Patients are being followed for one year, and what
9 I'm presenting today is a collection of the three-
10 month data.

11 Inclusion criteria included a sphere power
12 of minus 1 to minus 350; cylinder up to 2 diopters
13 at any axis but no more than 3.5 diopters at any
14 meridian. All participants must be at least 18
15 years of age and not more than 39 years of age at
16 enrollment, and have a visual acuity of at least
17 20/20 in each eye with manifest refraction.
18 Current RGP wearers were asked to abstain from lens
19 wear for at least two weeks prior to enrollment.
20 Soft lens wearers were not asked to do so.

21 Our patients were excluded if they showed
22 any sign of corneal disease or other ocular
23 disorders which would affect their vision; if they
24 had undergone previous refractive surgery or
25 previous or current orthokeratology treatment; if

1 they had any known sensitivities to the contact
2 lens solutions; use of any medications that might
3 alter the corneal curvature. We excluded pregnant
4 and lactating women and anyone who was intolerant
5 to contact lenses.

6 The examination visit schedule included a
7 baseline examination, which was a comprehensive
8 exam with the contact lens fitting; following that,
9 a dispensing visit; a one-day, one-week, one-month,
10 three-months, six-months, and one-year visit. At
11 the one-day through six-months visits, this
12 included a morning visit which was immediately
13 after awakening; they would get ready, come
14 straight in, we would remove their lenses for them.
15 The second visit was at least six hours later.

16 The testing that was conducted at each
17 visit included an unaided high and low contrast
18 Bailey-Lovie visual acuity; keratometry; corneal
19 topography; manifest refraction; biomicroscopy; a
20 contact lens comfort and self-rating of vision
21 questionnaire; and the RSVP quality of life
22 questionnaire.

23 At three months, 31 of 60 patients had
24 completed the three-month visit. Fourteen patients
25 discontinued prior to the one-month visit for

1 reasons such as poor adaptation, lack of
2 motivation, and treatment failure. Fifteen
3 additional patients discontinued between the one-
4 month and three-month visits for similar reasons.
5 And I might add that we believe that at least a
6 quarter to a third of the attrition may be
7 attributable to a learning curve on the part of the
8 investigators in fitting the lenses.

9 At baseline, the means spherical
10 equivalent refractive error was just over minus 2
11 diopters in each eye for the participants, and the
12 mean change that we noted at three months was also
13 just over 2 diopters. This graph shows a trend in
14 the manifest refraction data. You can see from
15 this graph that most of the change occurs between
16 baseline and one month, with about a quarter to a
17 half diopter of change also occurring between one
18 month and three months.

19 Unaided visual acuity at three months, 74
20 percent of the right eyes were 20/20, 61 percent of
21 the left eyes were 20/20 or better, and 93 percent
22 of right and left eyes were 20/40 or better.

23 At one month, 89 percent of the subjects
24 were within plus or minus 1 diopter of the target.
25 Those who were not, all except one were

1 undercorrected at that point. At three months, 90
2 percent of the subjects were within plus or minus 1
3 diopter of target, and again, mostly
4 undercorrected.

5 Biomicroscopy findings, no patients had
6 corneal infiltrates at any of our visits. Corneal
7 staining was noted in 77.4 percent of the patients
8 at the morning visit, but I must mention that only
9 two of those patients had a grade higher than a two
10 on that staining. By the afternoon visit this had
11 dropped to 37.9 percent, and none were higher than
12 a grade two, and only one case of staining had to
13 be treated.

14 Imprinting, which does not appear to be an
15 adverse finding, was found in 13.4 percent of the
16 patients at the three-month morning visit. There
17 was no sign of imprinting at the afternoon visit at
18 three months.

19 Microcysts were found in 38.7 percent of
20 the patients at the morning visit. All were less
21 than 10 in number, and that had reduced to 26.7
22 percent of the patients by the three-month
23 afternoon visit, again, all less than 10 in number.

24 Complications that were reported, we have
25 four documented. Two were for the same patient.

1 One patient had an insect that flew into her eye
2 while she was riding her bicycle. She was not
3 wearing her lenses at the time. However, she felt
4 better by the time she got home. She put the lens
5 back in. It aggravated the eye. Her corrected
6 visual acuity was 20/25 the morning that she came
7 in after that. It was 20/20 by the second morning.
8 She was treated with an antibiotic, and was able to
9 return to the lenses after treatment was complete.

10 A second patient reported pain in his
11 right eye in the morning after he inadvertently
12 wore his left lens on his right eye. Again, he was
13 treated prophylactically with an antibiotic, and
14 uncorrected visual acuity was 20/20. He also after
15 treatment was able to go back to wearing lenses.

16 One patient came in for a regular follow-
17 up visit and a foreign body was found in her eye.
18 The patient was unaware of the foreign body. It
19 was removed at that visit. She was treated with an
20 antibiotic, and she was able to return to lens wear
21 after treatment.

22 No permanent corneal damage was observed
23 in any of these complications. All patients were
24 able to resume lens wear after their treatment.

25 And in conclusion, I would just like to

1 say that we have found that orthokeratology can
2 produce an improvement in visual acuity and can
3 reduce myopic refractive error as long as retainer
4 lenses are worn. The changes can be maintained,
5 from our data, for at least six hours after lens
6 removal, and from our sample we have found it to be
7 a safe and efficacious procedure.

8 Thank you.

9 DR. WEISS: Thank you, Dr. Rah. I would
10 ask if any members of the panel have any questions
11 for Dr. Rah?

12 [No response.]

13 DR. WEISS: Thank you for that
14 presentation. Is there anyone else who would like
15 to make a comment during this open public hearing?

16 If not, this will close the public hearing
17 session. There are no division updates in the open
18 committee session which we will now start, so we
19 will then go on to the sponsor presentation of PMA
20 P870024/S043, and you can begin.

21 **PMA P870024/S043**

22 **SPONSOR PRESENTATION**

23 DR. MEYERS: Good morning, good morning.
24 My name is Bill Meyers. I'm the Vice President of
25 Science and Technology for Paragon Vision Sciences.

1 Today we are presenting the data on our PMA for
2 contact lens corneal refractive therapy.

3 I'd like to present the members of the
4 team who will be making that presentation:

5 Michael DePaolis, Dr. Michael DePaolis,
6 clinical investigator in this study; Dr. Mark
7 Bullimore, who is Associate Professor at Ohio State
8 University; Dr. Oliver Schein, who is a Professor
9 of Ophthalmology at Johns Hopkins University; Dr.
10 Jerry Legerton, who was the clinical monitor for
11 the Paragon Vision Sciences Study; and Dr. John N.
12 Quiring, who was the biostatistician and is
13 President of QST Consultations.

14 The device for which we are seeking
15 approval is a contact lens for corneal refractive
16 therapy. The purpose of the device is the
17 temporary reduction in myopia by the application of
18 a rigid gas permeable contact lens having on its
19 back surface a greater apical radius than the
20 pretreatment apical radius of the cornea.

21 The indication is for the temporary
22 reduction of naturally occurring myopia from minus
23 one-half to minus 6 diopter sphere in the presence
24 of cylinder up to 1.75 diopters. This is performed
25 in an overnight fitting program, differing from

1 extended wear in that in this indication one puts
2 the lenses on before going to bed and removes them
3 shortly after awakening.

4 We are seeking approval for this device in
5 both paflucocon-B and paflucocon-D. These
6 materials are already approved by FDA for seven-day
7 extended wear. We are submitting labeling for two
8 configurations of the device, the Paragon CRT
9 configuration, which is the one that was used in
10 the study exclusively, and the Quadra RG, a reverse
11 geometry design for which Paragon has already
12 received approval in the daily wear indication.

13 So, with no further ado, I'd like to
14 introduce Dr. Mark Bullimore for historical
15 considerations for this PMA.

16 DR. BULLIMORE: Good morning. My name is
17 Mark Bullimore. I'm an Associate Professor at the
18 Ohio State University. I have no direct financial
19 interest in the products being discussed this
20 morning. I am a paid consultant for Paragon Vision
21 Sciences, and I'm also a recovering panel member.

22 [Laughter.]

23 I'd like to start with a brief historical
24 review so that we can consider today's PMA in the
25 appropriate scientific perspective. There is a

1 literature on orthokeratology and the reshaping of
2 the cornea with rigid contact lenses.

3 Many of you would be familiar with the
4 Polse study, the Berkeley orthokeratology study.
5 Interestingly, this was one of the first randomized
6 clinical trials funded by the National Eye
7 Institute, conducted some 20 years ago, and this
8 involved randomization of patients to daily wear of
9 lenses. Patients were either randomized to a
10 standard fitting lens or an orthokeratology
11 modality where the lens was considerably flatter
12 than the corneal curvature.

13 Polse and his colleagues found a mean
14 reduction in myopia of around 1 diopter, an
15 improvement around about two and a half lines of
16 visual acuity, and a return to baseline after
17 discontinuation of lens wear. They also
18 established the safety profile for this modality.

19 Subsequent advancements in lens designs,
20 particularly the introduction of reverse geometry
21 lenses, resulted in the FDA approving two devices
22 for daily use. The data that are on file for these
23 two submissions demonstrate a mean reduction in
24 myopia of around about 1 and three-quarters
25 diopters, and both of these devices are approved

1 for up to 3 diopters of myopia.

2 A couple of studies I'd like to mention
3 that elucidate the mechanism underlying these
4 refractive changes. The first, by Swarbrick and
5 her colleagues, examined daily wear of
6 orthokeratology for a one-month period, and they
7 found about a 1 and three-quarter diopter reduction
8 in the myopia of their subjects. Interestingly,
9 this study was the first to establish that the
10 response was explained by a redistribution of
11 corneal tissue; that is, rather than a warping or
12 bowing of the cornea, it was a redistribution of
13 corneal tissue, principally the epithelium, that
14 resulted in these refractive changes.

15 A second study by Nichols and co-workers
16 corroborated this finding in terms of the
17 mechanism, but this time looking at an overnight
18 wear modality. Patients wearing the lenses for 60
19 days experienced a mean reduction of almost 2
20 diopters in myopia and a mean improvement of
21 uncorrected visual acuity of five to six lines on a
22 logMAR chart.

23 The final study I want to put before you
24 this morning is another study by Polse and his
25 collaborators. This was not a study of corneal

1 reshaping or orthokeratology. This was a study of
2 extended wear using traditional fitted lenses.

3 I raise this study for two reasons. One
4 is that it involves 201 subjects fitted to the same
5 FDA-approved extended wear materials that we're
6 going to discuss in this mornings PMA. The
7 subjects were randomized to one or either of these
8 materials and wore the lenses for a 12-month period
9 on a seven-day extended wear basis.

10 It's important to note as well that only
11 62 percent of these subjects completed the 12-month
12 trial. This is particularly interesting because
13 the subjects were first adapted to daily wear, then
14 went through a period of adaptation to extended
15 wear, prior to these 201 subjects being randomized
16 for these two materials.

17 I'd like to now talk about some of the
18 considerations in designing the clinical study
19 before you. As you can see, there was a multitude
20 of evaluations conducted, both pre-treatment and
21 during the post-treatment visits. The three I want
22 to draw your attention to are on the next slide.

23 These include low-contrast visual acuity.
24 This was done using a standardized protocol, using
25 Bailey-Lovie logMar charts and a bilateral scoring

1 methodology. The study also employed a
2 standardized manifest refraction technique. And
3 these two techniques are important because they
4 speak to the safety and effectiveness outcomes that
5 were the focus of this PMA. There was also careful
6 control regarding the time of examination after the
7 lenses had been removed in the morning.

8 Let me tell you the data that you're not
9 going to see today: cycloplegic refraction,
10 topography analysis, and endothelial cell count.
11 We decided not to perform cycloplegic refraction.
12 One of the reasons for doing this was, in terms of
13 data quality, cycloplegic refraction has been shown
14 to be less repeatable than manifest refraction.

15 Our other reasons for this choice speak to
16 the issue of respondent burden. We didn't want to
17 subject our subjects to periodic cycloplegic
18 refraction when we were dealing with a temporary
19 treatment for myopia, and we felt that the only
20 benefit of conducting cycloplegic refraction would
21 be to find some very rare cases of pseudo myopia or
22 other accommodative anomalies.

23 Corneal topography was performed
24 throughout the study by these investigators.
25 Topography was used as a screening tool at study

1 entry, and it was also used for quality assessment
2 of lens centration during the treatment phase of
3 the study. However, we did not subject post-
4 treatment corneal topography to any rigorous
5 statistical analysis.

6 Our reasons for not doing this are as
7 follows: We had 11 clinical sites with a multitude
8 of different corneal topographers. So in addition
9 to the problems of data reduction with these
10 devices, there was the issue of different machines
11 using different algorithms and different
12 technologies, and really there is no accepted
13 standard for the analysis of these data.

14 We also did not measure endothelial cell
15 count. I would like to point out, though, that
16 data on endothelial cell count are on record for
17 both of these materials, related to the seven-day
18 extended wear approval by the FDA of both of these
19 materials. I should also like to point out that
20 the average Dk/L, the transmissibility of the
21 lenses used in the current study, is greater than
22 those used in the extended wear methodology. I
23 would also again like to point out that we're
24 talking about overnight wear of lenses here, and
25 not extended wear.

1 One of the other key issues concerning
2 this PMA is the issue of consistent wear. The
3 patient doesn't get any benefit unless they
4 consistently wear the lens. It's a temporary
5 correction. It involves compliance on the part of
6 the subject throughout the study. They need
7 consistent wear in order to get the efficacy of the
8 device, and you'll see this when we address the
9 effectiveness data later on in this presentation.

10 I would now like to turn the podium over
11 to Dr. Jerry Legerton to discuss the device.

12 DR. LEGERTON: I am Dr. Jerry Legerton. I
13 am a paid consultant of Paragon Vision Sciences,
14 and have a financial interest in this product.

15 I want to help in our understanding here,
16 do a brief description of the device. That's
17 really made up of the design, the materials, and
18 the prescribing system.

19 As far as the design is concerned, the
20 primary interest in a lens that is fit with an
21 apical radius that's greater than the underlying
22 corneal apical radius is that of proximity control.
23 The other issue that we'll speak to is that of the
24 FDA-approved extended wear materials.

25 In lenses that are fit with a greater

1 apical radius than the underlying cornea, as the
2 lens radially--we proceed to the periphery, the
3 lens will then distance itself from the underlying
4 cornea, and a lens design then must have a
5 secondary zone that returns that lens to the
6 proximity of the cornea. It also must have a third
7 zone that is tangent to the cornea in the
8 peripheral area.

9 The materials involved are the paflucocon-
10 B and paflucocon-D with the ISO/ANSI Dk's of 40 and
11 100. Again, as Mark has indicated, these have
12 approval for seven days of extended wear.

13 The prescribing system used by the
14 investigators was one where they had a 65 lens
15 diagnostic system. They calculated the base curve
16 and the power. The base curve was determined to
17 have a radius that the cornea would need to have to
18 have emmetropia or low hyperopia, and that was
19 simply determined by calculation. But in each base
20 curve they had three depths of the return zone, and
21 could then by bracketing determine the depth that
22 would put them in proximity to the cornea. Then
23 there were nine peripheral configurations that were
24 used to find the design that would allow the lens
25 to be tangent to the peripheral cornea, actually

1 outside of a zone that topography could measure if
2 it were utilized anyway.

3 In terms of the investigational plan,
4 great effort was made to, as a sort of "first
5 through the gate" here, to make this
6 investigational plan comprehensive. We had 11
7 geographically distributed sites. We used a
8 double-masked, randomized protocol for the two
9 materials. We submitted the protocol and received
10 the FDA IDE approval for the protocol, the informed
11 consent, and the case report forms, and enrolled
12 the first adult subject and initiated treatment on
13 June 17th of 2000; the first adolescent on April
14 26, 2001; and the last subject initiated treatment
15 on August 23rd.

16 You have before you the list of
17 investigators and their locations. For the
18 inclusion criteria, we initially started with age
19 18 and older, and then that was changed to age 12
20 and older by approval of the FDA after the three-
21 month interim report was analyzed by staff and
22 approval was given.

23 The inclusion included patients with
24 between a half and six diopters of myopia with up
25 to 1 and three-quarters diopter astigmatism.

1 Patients were required to have demonstrated 12
2 months stability in their refractive status, within
3 a half diopter, and if they were rigid lens
4 wearers, were required to be out of their lenses
5 for four weeks and to demonstrate half diopter
6 keratometry stability on two serial measures. They
7 also had to be willing to return for 12 months of
8 follow-up.

9 The exclusion criteria follow those that
10 would normally be considered in a refractive
11 surgery procedure or refractive study.

12 The pretreatment evaluation again was
13 comprehensive, used to screen and determine the
14 inclusion and exclusion criteria, and had
15 particular measurements that would allow us to
16 measure the endpoint, comparing baseline and
17 follow-up measures. Those particularly were the
18 distance on corrected high-contrast logMAR, the
19 manifest refraction, the distance best corrected
20 logMAR, the keratometry, and the psychometric
21 questionnaire.

22 Hence, the follow-up evaluations have both
23 measures, as well, in your unaided and best-
24 corrected logMARs, your manual keratometry, your
25 applanation tonometry, and to your segment exam and

1 your psychometric questionnaire.

2 As far as the examination schedule, beyond
3 the baseline we did conduct examination at
4 dispensing; a one-day visit which was required to
5 be within one hour of awakening. I want to further
6 mention that in the protocol, the patient was
7 instructed to place the lens no more than 30
8 minutes before retiring, and to remove the lens
9 within 30 minutes of awakening; that no wear beyond
10 that 30-minute time was sanctioned, with the
11 exception of this one one-day visit. So there was
12 no notion that more is better. This was strictly
13 an overnight procedure.

14 On the two-week, one-, two-, three-, six-,
15 nine-month visits, the patient was to return to the
16 office within three hours of awakening, and hence
17 with no lens in place. The patients were to bring
18 their lenses in case a practitioner wanted to
19 reapply the lens, but no lens was in place at the
20 time of the visit. Subjects were required to have
21 a, we'll call it a wash-out series of examinations
22 at 8, 24, 48, and 72 hours after lens removal in at
23 least one eye.

24 And this was double-masked and randomized
25 for the two materials. The randomization schedule

1 provided for intra-investigator randomization that
2 would net an equal number of subjects in the two
3 materials by 10 subjects, and there was an inter-
4 investigator randomization utilizing six different
5 randomization schedules. This was masked from both
6 the investigator and the subject.

7 In terms of study monitoring, study
8 initiation was conducted to both review the
9 protocol and then to conduct a training in the
10 diagnostic lens evaluation. Through the course of
11 the study there were a minimum of two monitoring
12 site visits, and the utilization of telephone, e-
13 mail and written correspondence to keep a healthy
14 flow of communication between the monitor, clinical
15 monitor, and the investigators.

16 In terms of the baseline data, we are
17 directing your attention to the demography of the
18 treated and the completed subjects, and we've put
19 them on a single table here for you. You will see
20 that there is a disproportionate number or make-up
21 in gender, more female subjects than male, and if
22 you look to the right, in comparing the treated
23 with the completed, you'll see that shifts a
24 little, where we have a combination of more
25 persistence and more compliance on the part of

1 female subjects.

2 In terms of eyes, there were a near-equal
3 number. There were only two subjects that had only
4 unilateral treatment. And in terms of current
5 contact lens history at the time of baseline,
6 you'll see that the make-up has a high composition
7 of prior contact lens wearers. None of the rigid
8 lens wearers were in a former orthokeratology
9 format, nor had any of the subjects in this study
10 previously had any corneal reshaping technology
11 applied.

12 The age, if you again compare the treated
13 to the completed, the treated includes the
14 adolescent cohort or adolescent arm of the study,
15 with a mean age of just under 34. The completed
16 subjects, as we will speak in the efficacy
17 analysis, are all adult subjects with a mean age of
18 37.5.

19 Further, if we look at the make-up of the
20 treated and completed subjects relevant to manifest
21 refraction spherical equivalent, you will see the
22 percentages applied there from zero to 2 diopters,
23 26 percent, which dropped to 20.2 percent in the
24 completed group. One might speculate a reason,
25 that the need for the treatment may be less in

1 lower groups. The modal group, we'll say from 2 to
2 4, increased, and we saw a slight lesser decrease
3 in the 4 to 6 than we saw in the zero to 2 group.

4 In terms of analysis, we believe it's
5 helpful for you if you can look at the analysis
6 cohort here, and with some understanding of the
7 breakdown and how we met different groups to do
8 different analyses. The baseline evaluation was
9 applied to 218 subjects. Thirteen subjects, most
10 of which lenses were ordered, were not dispensed,
11 or they were not ordered at all.

12 The net result is that 20 subjects in fact
13 were treated, and we defined "treatment" as at
14 least one night of overnight wear. Of that, we
15 have those that were not due, that 50 subjects
16 that--correct, 50 subjects that were not due for
17 their nine-month exam at the time that the data
18 base was frozen for this analysis. That resulted
19 in 155 on the nine-month, that were due for a nine-
20 month visit. Sixty of those had discontinued,
21 which then gave us the 95 that completed the nine-
22 month visit.

23 Eleven of those had intermittent wear, and
24 as Mark mentioned, the need for consistent wear,
25 especially prior to the visit, is important to get

1 an efficacy, we'll say qualified or valid cohort.
2 And so those 11 subjects that either had lost a
3 lens just prior to the nine-month visit or had had
4 intermittent wear were not included in the efficacy
5 analysis, but were included in the safety analysis,
6 and we then net your 84 subjects that were studied
7 for the efficacy analysis.

8 Discontinuation did occur. Our
9 discontinuation rate is very much similar to the
10 Polse contact lens extended wear study using the
11 same two materials, and you will see that the
12 primary reason for discontinuation is unacceptable
13 vision. In fact--coming back to that, if you
14 would, Tim--of note here, I want to point out, so
15 that we don't raise an inconsistency concern, that
16 this report of 71 subjects discontinuing, I said 60
17 subjects discontinued, but that was 60 of those
18 that were due for the nine-month. But also we had
19 another 11 that discontinued, that were not yet due
20 for the nine-months.

21 In terms of accountability, you will see
22 that our accountability for safety exceeds 94
23 percent at all visits. Most visits it's higher
24 than that. And our accountability for efficacy
25 exceeds 93 percent at all visits, and most visits

1 again higher than that.

2 There were lens reorders allowed in the
3 clinical trial, and these lens reorders could be
4 studied in terms of two rationales. One rationale
5 would be that precision is required to get the
6 lenses to perform the way that we desire; that
7 controlling of proximity is important, and
8 centration is absolutely instrumental and essential
9 to get a good result. And hence you will see that
10 poor centration is the predominant reason for lens
11 replacement.

12 Ten lenses were replaced for loss, damage,
13 or deposit. So in a sense for the completed eyes,
14 188 completed eyes, we have 60 lenses that were
15 reordered for purpose of either refining the
16 parameters of the lens, improving centration, or
17 because of, the next leading reason would be an
18 undertreatment.

19 I now introduce Dr. Oliver Schein, who
20 will present to you the safety analysis.

21 DR. SCHEIN: Good morning. My name is
22 Oliver Schein, and I am here as a paid consultant
23 to Paragon, and in this capacity I have provided
24 input to them regarding the analysis and
25 interpretation of the data from their clinical

1 trial. I have no proprietary interest in the
2 company or the device.

3 I'm going to present, probably in the next
4 five to six minutes, data related to the safety on
5 the entire cohort. Following this, Dr. DePaolis
6 will present analogous data on the efficacy from
7 the entire cohort. This is a randomized trial.
8 Dr. Legerton will return following that
9 presentation and show you some comparison data
10 between the two devices that were compared.

11 Now, there were five safety endpoints that
12 were specified in the original protocol, and first
13 and probably the most important is loss of best
14 spectacle corrected visual acuity. So as not to
15 say that times, that will be BSCVA. The second is
16 a state of having a BSCVA of worse than a 20/40
17 threshold. Third is serious adverse events,
18 fourth, slit lamp findings, both as traditionally
19 defined in FDA premarket contact lens trials. And
20 the fifth is symptoms and complaints.

21 So starting with BSCVA, there was no
22 permanent BSCVA loss observed. By permanent we
23 mean persistent through the course of the study.
24 However, there were transient losses that were
25 observed on each visit. For example, at the one-

1 and two-month visits there were approximately 9 and
2 7 percent of the population had a BSCVA loss of at
3 least two lines. From three months and
4 subsequently, this leveled off to about 4 percent.

5 This is termed "transient" because indeed
6 these individuals at subsequent visits were not
7 found to have a BSCVA loss. I interpret this slide
8 myself as indicating that at any random moment in
9 time, at least based on the data available
10 currently, after the initial break-in one might
11 expect about 4 percent of individuals using this
12 technology to have a BSCVA loss of two or more
13 lines at a particular moment in time.

14 So these are BSCVA losses, not a BCVA
15 loss, which is a best corrected visual acuity loss.
16 In other words, the acuity could be regained with a
17 contact lens refraction at that same visit. I've
18 already explained why they've been termed
19 transient.

20 The second outcome there, then, was the
21 BSCVA status worse than 20/40 at each follow-up
22 visit. These are listed for you here. They're
23 largely in the 1 to 2 percent range. None was
24 recorded at the nine-month visit. Again, for each
25 one of these occurrences, they were in a sense

1 transient, since at subsequent visits they were not
2 measured in the same people.

3 Adverse events. There were no adverse
4 events that met the criteria for serious or
5 unanticipated. There was one instance of a corneal
6 abrasion upon removing the lens; another case
7 diagnosed as bacterial conjunctivitis; and a
8 moderate case of microcystic edema.

9 Grade 3 slit lamp findings. The
10 traditional way of presenting these is simply to
11 count the observations of a particular finding
12 across all the visits. It gives you a total number
13 of observations without regard to whether they were
14 all within the same person or not.

15 Using that kind of crude approach, there
16 were 28 observations overall of Grade 3 or higher.
17 Eighteen of them were of edema. Of note, 17 of
18 these 18 occurred or were reported from one
19 practice, and that practice is at a very high
20 altitude, and that might well be relevant. Ten
21 eyes of six subjects of these 18 observations
22 occurred on the first day, and six of the
23 observations were contributed by one subject over
24 multiple visits.

25 Nine observations of staining at Grade 3

1 level occurred in seven eyes of six subjects, and
2 there was one observation of conjunctival
3 injection. Again, no Grade 4 or what are termed
4 severe observations were noted in the study.

5 Symptoms and complaints. Discomfort was
6 quite prevalent at dispensing. About three-
7 quarters of individuals fit with the lens
8 complained of some degree of discomfort. Now,
9 unfortunately this data or this discomfort data was
10 not in any way scaled. In fact, it's apparent that
11 it's really uncertain as to whether the discomfort
12 was reported even with or without the lenses, but
13 some discomfort was certainly present and prevalent
14 in three-quarters of individuals at dispensing, and
15 down to about 20 percent by nine months.

16 The only indirect indication that I could
17 find as to the clinical significance of the
18 discomfort was to look at the reasons for dropout.
19 And as you have been shown, the most common reason
20 for dropout was inadequate level of vision. 3.4
21 percent of the population reported discomfort as
22 the reason for discontinuation.

23 So to summarize, there were no persistent
24 BSCVA losses or BSCVA worse than 20/40 at nine
25 months, nor were there serious adverse events

1 observed in this cohort. There were 28
2 observations of Grade 3 slit lamp findings. They
3 all resolved. None of the Grade 4 was observed.
4 Discomfort is prevalent initially, goes down over
5 time, and doesn't seem to be the principal reason
6 for discontinuation.

7 Dr. DePaolis will now present the efficacy
8 data.

9 DR. DePAOLIS: Thank you, Dr. Schein.
10 Thank you, members of the panel. I may draw to
11 your attention the cohort analysis sheet is part of
12 your handout. I know Dr. Legerton discussed it in
13 great detail earlier, but it may also help better
14 understand the different sample sizes at different
15 points during the analyses.

16 My name is, again, Michael DePaolis, and I
17 am a clinical investigator for Paragon Vision
18 Sciences in the device discussed here within. In
19 that capacity, I am a paid consultant. However,
20 aside from being a clinical investigator, I have no
21 proprietary interest in Paragon Vision Sciences or
22 any of their products.

23 I'm here to discuss the five primary
24 measurements that we looked at in terms of a device
25 efficacy. They included unaided visual acuity,

1 with particular attention to the percentage of
2 patients who were able to achieve 20/20 or 20/40
3 visual acuity. We looked, number two, at reduction
4 in manifest refraction spherical equivalent, or
5 MRSE. We, number three, looked at predictability,
6 with tolerance levels of a half a diopter in 1
7 diopter. We looked at, number four, stability
8 between two consecutive visits, again with changes
9 of less than a half diopter in 1 diopter
10 refractively. And then we looked at alterations in
11 corneal curvature and absolute corneal astigmatism.

12 If we looked at uncorrected visual acuity
13 at the nine-month visit--and you'll notice that
14 these are folks who were targeted for emmetropia,
15 and include a smaller cohort for those who were
16 actually able to reach 20/20 acuity predicated on
17 pretreatment 20/20 acuity levels--you'll see that
18 we were able to achieve acuities of 20/20 or
19 better, uncorrected, in 58.4 percent of the
20 population, and in 89.8 percent of the population
21 we were able to achieve uncorrected visual acuities
22 of 20/40 or better. Recognizing that this is just
23 one data point in time, we'll look at the next
24 slide, which will give us unaided visual acuity
25 over time.

1 I think we glean a couple things from this
2 chart. First and foremost, you can see most of the
3 therapeutic effect was demonstrable within the
4 first month of treatment. And, number two, we can
5 see that over time the percentage of eyes
6 stratified for different refractive errors will
7 actually be similar across the different time
8 points for 20/20 uncorrected visual acuity as well
9 as 20/40 uncorrected visual acuity.

10 The second barometer we looked at was
11 reduction in manifest refraction spherical
12 equivalent, and I did want to draw to the panel's
13 attention that the efficacy cohort for this
14 determination included patients with a mean
15 refractive sphere of almost 3 diopters of myopia
16 and a mean refractive cylinder of almost an
17 additional half diopter. As you can see, all but
18 one of our eyes at the nine-month data point was
19 able to achieve some reduction in myopia, with the
20 mean reduction for the cohort being 2.59 diopters.

21 Predictability was the third issue that we
22 wanted to take a look at from an efficacy
23 standpoint, and as this graph depicts, we can see
24 patients' refractive outcomes within plus and minus
25 a half a diopter. This is achieved in terms of

1 plus or minus a half diopter of target, and
2 achieved within a plus and minus 1 diopter of
3 target.

4 You can see throughout the various data
5 points there's a slow upward trend from one month
6 to three months, and then things stabilize pretty
7 much beyond the three-month window.

8 I'm not sure that the discrepancy between
9 plus and minus a half a diopter outcomes and plus
10 and minus 1 diopter outcomes is exclusively related
11 to the efficacy of the device, and may also speak
12 specifically to the limits of agreement in manifest
13 refraction techniques with half diopter tolerances
14 versus one diopter tolerances.

15 Another way of looking at predictability
16 is a scattergram, where you can see the limits are
17 plus and minus 1 diopter. What I glean from this
18 particular chart is three variables. One is, the
19 device seems to be fairly predictable. Number two,
20 it seems to be fairly predictable over attempted
21 correction ranges from 1 diopter of myopia up to 6
22 diopters of myopia. And last, but certainly not
23 least, particularly germane for our presbyopic
24 patients, there were no overcorrections of 2
25 diopters or greater. In fact, we only had one eye

1 with an overcorrection of a diopter.

2 Our fourth variable was stability of
3 manifest refraction spherical equivalent between
4 subsequent visits. And when you look at this
5 chart, in one respect you can look all the way over
6 to the right graph which depicts changes between
7 six and nine months that are a diopter or less in
8 manifest refraction spherical equivalent, and argue
9 that that may be a bit liberal. You could also go
10 to the other end of the chart and again, given the
11 limitations of manifest refractions, say that
12 perhaps the half diopter limit is not a good limit
13 to look at either.

14 I will draw your attention to the middle
15 set of graphs, which is those patients, the
16 percentage of patients that between the six- and
17 nine-month visit had refractive changes, spherical
18 equivalent refractive changes, of less than three-
19 quarter diopters, and demonstrate that stratified
20 across different refractive categories, we had a
21 fairly high percent of patients achieve that
22 outcome.

23 Another way of looking at refractive
24 stability is the mean of difference between
25 subsequent visits for manifest refraction spherical

1 equivalent and for average keratometry findings.
2 And, as you can see here, at data points between
3 three and six months and then again at data points
4 between six and nine months, we had fairly good or
5 fairly tight limits.

6 In fact, if we depict this in the next
7 slide, again you'll see that most of the
8 therapeutic impact occurred within the first month,
9 with a reduction in keratometry and mean refractive
10 spherical equivalent being aptly noted. And then
11 as we go forward from two months through nine
12 months, the mean of difference between visits
13 became sequentially smaller, without evidence of
14 regression, I might add.

15 The fifth and final variable we looked at
16 was change in corneal curvature. You could see
17 from the previous slide that there was in fact a
18 positive change in corneal curvature, but in the
19 spirit of full disclosure I share this data with
20 you, simply because it's a big enigmatic for me.

21 If you look at the graph, you'll see that
22 we actually had changes, absolute changes in
23 corneal cylinder, that were fairly inconsistent.
24 The percentage of eyes that had no change, the
25 percentage of eyes that had modest decreases, and

1 the percentage of eyes that had modest increases
2 clearly do not depict a trend.

3 Fortunately, I think what we glean from
4 this graph are two things: Number one, the changes
5 in absolute corneal cylinder, be they positive or
6 negative, were relatively modest. And, number two,
7 this may again implicate that our ways of
8 traditionally measuring corneal cylinder in this
9 subset of the population may not be as accurate as
10 we would like it to be.

11 In light of that last primary outcome,
12 there were two secondary outcomes we elected to
13 look at. One was a reduction in refractive
14 cylinder, and if we look at the next slide, we'll
15 see that the number of eyes that actually had a
16 reduction in refractive cylinder was approximately
17 50 percent, with a fair number of eyes having no
18 change, and fewer having an increase in refractive
19 cylinder.

20 The last secondary outcome that we looked
21 at was derived from our psychometric analyses in
22 which we determined patient satisfaction with
23 unaided vision. As you would anticipate with a
24 device of this nature, as we went from pretreatment
25 to six month to nine month, the percentage of

1 patients who reported good, very good, or excellent
2 unaided vision increased and fortunately remained
3 fairly consistent.

4 In fact, if we take a look at the nine-
5 month data point, we find the percentage of
6 patients who reported that their unaided vision as
7 a result of corneal refractive therapy was either
8 good, very good, or excellent, far surpassed their
9 pretreatment unaided visual acuity, and was
10 actually fairly comparable with the percentage of
11 patient in the cohort who reported good, very good,
12 or excellent visual acuity with their habitual
13 means of correction going into the study.

14 So if I were to summarize the efficacy
15 outcomes, particularly at the nine-month visit,
16 this is depicted on the slide and I think it pretty
17 much reiterates what we've just discussed in the
18 past few minutes. We did see a pretty consistent
19 reduction in--improvement, I should say, excuse me--
20 -in uncorrected visual acuity. We did see a
21 reduction in both manifest refraction spherical
22 equivalent and flat keratometry. We found the
23 device to be predictable and stable within
24 measurable tolerance levels of a half diopter in
25 both scenarios.

1 I will at this point turn the podium back
2 over to Dr. Jerry Legerton, who is going to give us
3 the specific analysis of material.

4 DR. LEGERTON: Jerry Legerton. The
5 analysis by material was conducted using the Rosner
6 statistical method, which recognizes that two eyes
7 within an individual are not independent, but in
8 fact correlated, and makes adjustments in the p-
9 values for that correlation.

10 We looked at the variables for efficacy of
11 uncorrected visual acuity, predictability and
12 stability, and for safety we looked at the
13 frequency of slit lamp findings or the proportion
14 of slit lamp findings and the symptoms, problems
15 and complaints. For the analysis by material--and
16 if you want to take a pen here and correct your
17 handout, thanks to our trusty statistician, we
18 found that we had made an error in the slide last
19 night by looking at a wrong proportion after the
20 materials were unmasked and switched for some
21 subjects.

22 But I want to explain how the material
23 masking and unmasking went for some subjects. That
24 of the 408 eyes, 208 by the randomization schedule
25 were put into paflucocon-B, 200 in paflucocon-D.

1 The protocol allowed the investigator, in the event
2 of concerns in outcome, to contact the monitor and
3 have subjects unmasked. Over the course of the
4 trial, 12 subjects were unmasked.

5 The reasons for unmasking would be, for
6 example in a case of edema that they believed was
7 due to hypoxia, they could unmask, they could call
8 and have an unmasking, and if it was deemed to be
9 an appropriate step, we agreed that they could do
10 that and they could reorder the lenses. Also, if
11 there were problems of indentation or adherence,
12 they could unmask and they would have the potential
13 to go either direction.

14 Again, that occurred in the trial in 12
15 subjects overall, and all of those switched from
16 the pafluvocon-B to the -D. Of those that
17 switched, one subject switched back. So the end
18 result, which will be consistent in the report that
19 you were supplied with the slit lamp findings by
20 material, your end will check out that you will
21 have 184 eyes that were only in pafluvocon-B--in
22 other words, they didn't involve those that were
23 switched--and there will be 200 that were in
24 pafluvocon-D.

25 At any rate, it was those, it was the non-

1 switchers that were studied, so we took those that
2 were only primarily in one material. And using the
3 Rosner method to these efficacy variables,
4 uncorrected visual acuity, predictability and
5 stability, there were no statistical significant
6 differences in the efficacy outcome.

7 The Rosner method was also applied to the
8 slit lamp findings, the symptoms, problems and
9 complaints. And for slit lamp findings there were
10 then 48 statistical tests of hypothesis, that being
11 the different slit lamp measures over all of the
12 visits would give 48 statistical tests, and edema
13 was the only one that approached statistical
14 significance at a 95 percent confidence level, and
15 even that did not support. I'll go back to that.

16 On symptoms, problems and complaints,
17 there were 72 statistical tests of hypothesis, and
18 halos at the two-month visit was the only one that
19 showed statistical significance at the 95 percent
20 confidence level. Given the number of statistical
21 tests over both materials, over all the findings,
22 this failed to support a hypothesis that there was
23 a difference between materials, or supported the
24 null.

25 You have before you a summary table of the

1 nine-month safety and efficacy variable. In the
2 first column it is all eyes, both materials. The
3 second two columns would be the paflucocon-B only
4 and the paflucocon-D only. In other words, those
5 second columns would not include the switchers.

6 I think what you will appreciate here is
7 that in numerical value, as well as given the lack
8 of statistical significance, but also to appreciate
9 numerical value, that you will see that the results
10 with the two materials--again, one being a high Dk
11 material, one being a super Dk material, are
12 substantially equivalent through both the safety
13 and the efficacy variables.

14 We thank you for your consideration of our
15 presentation and our submission, and we ask you to
16 consider that these data establish the safety and
17 efficacy of contact lens corneal refractive therapy
18 in paflucocon-B and -D for the overnight treatment
19 of myopia and myopia with astigmatism.

20 **Panel Questions for the Sponser**

21 DR. WEISS: The sponsor can remain at the
22 table, and we are going to proceed to questions
23 from the panel for the sponsor. Dr. Van Meter?

24 DR. VAN METER: Since there was a
25 substantial dropout because of patients that were

1 uncomfortable with their lenses or to some degree
2 dissatisfied with their vision, is there anything
3 that you would do differently if you were fitting
4 another 200 patients to rule out some of these
5 patients who you might argue in retrospect probably
6 shouldn't have been fit with the lens anyway? Did
7 you learn anything from the dropouts?

8 DR. LEGERTON: Yes. More than just life
9 is learning. It was important to us to establish a
10 good diagnostic procedure before we even started,
11 and certainly select qualified investigators, but
12 through the course of the clinical trial we
13 certainly discovered that refinements could be made
14 to the prescribing system. The real importance of
15 establishing that lens centration early on, lenses
16 that don't have adequate depth are going to
17 decenter or tip-tilt and they're going to be more
18 uncomfortable. So I think a lot of it is not so
19 much in the screening of the patient but rather in
20 the diagnostic lens procedure. The next 200 we
21 would do differently in terms of more refinement of
22 the prescribing system.

23 DR. VAN METER: Well, this is my point.
24 So it's really a fitter's learning curve rather
25 than a wearer's learning curve.

1 DR. LEGERTON: Could you repeat that?

2 DR. VAN METER: So you would say it's more
3 of a contact lens fitter's learning curve rather
4 than a contact lens wearer's learning curve?

5 DR. LEGERTON: Correct. Like a well-fit,
6 well-centering lens improves in comfort. And also
7 the whole issue of unacceptable vision, if you do
8 have lens decentration, much like a decentered
9 ablation, you're going to have an unacceptable
10 visual result. It may not be so much that you
11 didn't result the myopia but that you don't have a
12 centered applanation.

13 DR. VAN METER: All right. Thank you.
14 This is Van Meter, for the record.

15 DR. WEISS: Thank you. We're going to
16 have Dr. Grimmett, Dr. Matoba, Dr. McMahon, and
17 then Dr. Harris.

18 DR. GRIMMETT: Mike Grimmett. First a
19 comment, then a question to Dr. Bullimore. I
20 certainly appreciate the optimistic intimation of
21 hope for us panel members, given your role as a
22 recovering panel member.

23 [Laughter.]

24 And now the question. For this indication
25 of overnight orthokeratology, and given the

1 theoretical benefits of higher Dk materials, I was
2 wondering why would a clinician fit the lower Dk
3 material, the -B type, rather than the higher Dk
4 material, the -D type, for any circumstance. Is
5 there a huge cost difference between them, or some
6 other issue?

7 DR. BULLIMORE: This is Mark Bullimore.
8 As far as I'm aware of, there are no cost
9 differences between the materials. In addition to
10 the different permeability of the materials, they
11 do have subtle differences in the handling, surface
12 properties. But I would address the question in a
13 couple of ways.

14 First, I would like to point out that,
15 again, that both of the materials currently are
16 approved for seven-day extended wear. Secondly, we
17 believe our data support that there is no
18 significant difference between the performance of
19 those two materials. And, thirdly, just to
20 reiterate that we are again talking about an
21 overnight wear modality and not extended wear, so
22 one could argue that the permeability of the
23 material is not quite so critical in this modality
24 as it would be in extended wear.

25 I hope that answers your question.

1 DR. GRIMMETT: Yes, it does. Mike
2 Grimmett again. I certainly realize that the data
3 did not show a benefit or superiority of one design
4 over the other, and I recognize that. However, as
5 a commonsense approach, as a clinician it would
6 seem to me, if I were the one fitting it, I would
7 want the higher Dk material just because it has
8 theoretical advantages.

9 DR. BULLIMORE: Mark Bullimore again.
10 Yes, I mean, we would be happy to acknowledge the
11 theoretical advantages, but we really want to give
12 maximum flexibility to the people manufacturing or
13 fitting these lenses to the patient, so that they
14 could go with their own choice. I mean, I don't
15 want to call it a practice of medicine or practice
16 of optometry issue, but it's sort of at that level
17 of preference. I would certainly anticipate that
18 most of our colleagues, given the choice, might go
19 for a higher material, but we're trying to give as
20 broad an opportunity here as possible. So I would
21 acknowledge your concern. Okay?

22 DR. GRIMMETT: Thank you.

23 DR. WEISS: Dr. Matoba?

24 DR. MATOBA: I'm directing this question
25 to the patients who discontinued the study because

1 of dissatisfaction with their vision. Do you know
2 what the mean refractive error was in these
3 patients? That is, I'm wondering if they were the
4 people with a higher level of myopia, or if there
5 was something to be gleaned from the
6 characteristics of those patients pretreatment.

7 DR. LEGERTON: The make-up of the subjects
8 that discontinued for unacceptable vision is not
9 skewed for the higher pretreatment refractive
10 errors, and that is very much consistent with the
11 answer that I gave to Dr. Van Meter, that this may
12 be more an issue--unacceptable vision can also be
13 tagged to problems of actually parameter selection
14 within the lens and gaining centration and so
15 forth.

16 DR. MATOBA: My second question is in
17 regard to the patient satisfaction with vision
18 data. I know that at nine months you said that 91
19 percent had very good unaided vision. At one
20 month, what percentage of patients would have said
21 the same thing, roughly? I'm wondering how long it
22 takes to reach that level of 90 percent satisfied,
23 good to excellent.

24 DR. LEGERTON: I have to--that analysis,
25 the analysis has not been conducted yet on the one

1 month psychometric questionnaire.

2 DR. MATOBA: Those patients who had
3 suboptimal correction, and maybe not quite good
4 enough to function comfortably, were those people
5 that had just gotten to this state, or what were
6 they told, what to do? What were they told to do,
7 how to deal with that situation?

8 DR. LEGERTON: Well, first, if one were to
9 speculate on the one month, because of the
10 clustering of discontinuation, where
11 discontinuation occurred, the majority is between
12 dispensing and two months. So if one were still in
13 treatment at one month and filling out a
14 psychometric questionnaire, I think you would
15 speculate that your percentage of good, very good,
16 or excellent is lower and that that is part of the-
17 -that would pave the way to an understanding of why
18 someone would discontinue for unacceptable vision.

19 But the strength of this modality is the
20 capacity to redesign the lens, and we were rather
21 parsimonious about allowing freedom to do all kinds
22 of retreatments in the clinical trial, where in
23 clinical practice that's something that the
24 practitioner might be more free to do, to do more
25 parameter changing and refinement. So I think

1 there are two sides to that. One is when you
2 intervene to redesign a lens, and how you respond
3 to patient complaints of unacceptable vision.

4 DR. WEISS: Dr. McMahon?

5 DR. McMAHON: Tim McMahon. I have a few
6 questions that I hope you can clear up some
7 confusion for me.

8 You stated this was a randomized mass
9 trial, yet on Tab D, page 265, your opening
10 statement, it says that "This is a prospective
11 multi-center, open label, non-randomized trial."
12 In the very next statement you indicate that it is
13 randomized by material. Can you clarify to me what
14 is--is the only thing that's randomized is the two
15 materials, and everything else is--

16 DR. LEGERTON: The only thing that was
17 randomized was two materials.

18 DR. McMAHON: Okay. The second question I
19 have is more of a practical question, and correct
20 me if it's inappropriate. You indicate that the
21 CRT is a trademarked name from the sponsor, from
22 Paragon, yet you also use the term in the
23 procedure, so that CRT corneal refractive therapy
24 and CRT lens are the same thing. Do you propose
25 that the CRT be the name of the procedure for your

1 product only, or are you going to ask the community
2 to rephrase the term "orthokeratology" to "CRT" for
3 all practitioners

4 DR. MEYERS: Bill Meyers with Paragon
5 Vision. The name "CRT" is trademarked by Paragon.
6 Corneal refractive therapy is to be available to
7 anyone. We have made an attempt to trademark that
8 name, with the intention that we assure that no one
9 could lock it up. We do not insist that anyone use
10 that name. Anyone who prefers the name
11 "orthokeratology" can certainly do so. We prefer
12 not to.

13 DR. McMAHON: Thank you. The next
14 question is, thank you for reterming "retreatment"
15 to "reordering" for at least 10 of those patients.
16 There are 70 lenses that were reordered for
17 patients, or 70 patients that had reorders, 60 of
18 which you indicated were for improvements or
19 changes in centration or overcorrection or
20 undercorrection, yet I didn't see any information
21 as to the effect of the change of those lenses. Do
22 you have that data, and can you describe what the
23 effect of these retreatment lenses were?

24 DR. LEGERTON: There were 60 lenses
25 ordered, reordered for purpose of retreatment,

1 we'll say, or refinement of prescription, for 188
2 eyes. We have not conducted a specific analysis as
3 to what the effect of those retreatments were, and
4 I can only speak as, let's say, the monitor and
5 having some processing of the case report forms.

6 Some patients--well, again, these are for
7 the completed eyes, so it would be fair to conclude
8 that the retreatment was supportive to the fact
9 that the patient survived the clinical trial. But
10 what specifically it did, what it took them from
11 and to, we have not conducted that analysis.

12 DR. McMAHON: And can you clarify for me
13 when these retreatments occurred? Was it typically
14 the first month?

15 DR. LEGERTON: The protocol stipulated
16 that retreatments could occur up to that point, and
17 only thereafter with permission from the monitor,
18 from the clinical monitor.

19 DR. McMAHON: Do you know what percentage
20 of those were after that one-month point?

21 DR. LEGERTON: I would be speculating, but
22 from again being involved with it, my speculation
23 is probably fairly accurate, that about two-thirds
24 of them were within one month and about another
25 third were thereafter.

1 DR. McMAHON: Okay. Thank you. And
2 additional question: In the cohort groups that Dr.
3 DePaolis mentioned, there is the group of 125
4 patients that achieved 20/20 BCVA. Since the
5 entrance criteria for all your patients was
6 basically to have spectacle correction 20/20 VA,
7 does that suggest then there is a number of
8 patients who were unable to see 20/20 with the CRT
9 design? Is that how I should interpret that?

10 DR. LEGERTON: Yes. That was cited as one
11 of the deviations in the report. The technique
12 used for high contrast logMAR was the Bailey-Lovie
13 letter counting technique, stipulated at a minimal
14 test distance of 14 feet. All but one site was
15 testing at less than 20 feet, most clustering about
16 that 14-foot distance. Since they were recording
17 it in letter count, 50 letters at 14 feet, 45
18 letters at 14 feet, when then that is calculated to
19 logMAR value, and we allow a logMAR equivalent of
20 20/20 minus 2.04, when that calculated out to
21 logMAR, there were logMARs that were below .04;
22 could have been .06, .08.

23 The practitioner believed that their
24 patient had 20/20 because they had done a
25 refraction in an exam chair on a Snelling chart,

1 and may have had to trial frame that and move to
2 their 14-foot location to do their logMAR test with
3 the illumination control, the proper distance and
4 so forth. So, again, their belief was the patient
5 was a 20/20 patient, but when we statistically--
6 when we enter the database and do a statistical
7 analysis, apply the logMAR calculation, they
8 calculated out at below 20/20.

9 DR. McMAHON: Last question for right now
10 would be, Dr. DePaolis presented some data that I
11 hadn't seen with regard to decrease in cylinder.
12 Did you guys do vector analysis, which is really in
13 a circumstance like this, since you haven't
14 mentioned anything about what happens to cylinder
15 axis--and Dr. Bradley can actually comment on this
16 much more elegantly than I can--did you guys
17 perform vector analysis on the cylinder changes?

18 DR. LEGERTON: We did not conduct vector
19 analysis, and that does require clarification.
20 It's the indication, the difference between an
21 indication requesting a claim of treatment of or
22 full correction of astigmatism, and an indication
23 allowing for the inclusion of patients with
24 astigmatism.

25 DR. McMAHON: You missed me on that one.

1 DR. LEGERTON: So we are not claiming, we
2 are not claiming and these data don't support that
3 we are correcting astigmatism.

4 DR. McMAHON: I understand that. That
5 wasn't the point of my question.

6 DR. BULLIMORE: May I? We did consider--
7 this is Mark Bullimore--we did consider subjecting
8 our data to vector analysis, but we have declined
9 to do that thus far. Really vector analysis, as
10 you mentioned, is a very powerful technique, but
11 its ultimate value is when there is an attempted
12 astigmatic correction, when it allows you to
13 compare, for example in astigmatic PRK, the
14 attempted astigmatic correction and the actual
15 achieved astigmatic correction.

16 In this PMA, the subtle changes that we
17 observed in refractive astigmatism--and I think we
18 could all agree that they were modest at best--were
19 not an intended outcome. They were in some
20 respects a secondary outcome and a pleasant benefit
21 to the procedure. So we would submit that vector
22 analysis would be of little benefit here because
23 there was no attempt to treat a specific amount of
24 astigmatism in any given patient.

25 On average, there was only half a diopter

1 of cylinder in the entire cohort, and I think the
2 change, the reduction that we got would be in the
3 order of a sixth of a diopter of astigmatism. I
4 have to check those numbers. But really it might
5 be of intellectual interest, but in terms of
6 practical benefit here and the indication for the
7 treatment, we thought it was no real benefit.

8 DR. McMAHON: Thank you.

9 DR. WEISS: Dr. Harris?

10 DR. HARRIS: Michael Harris. Several
11 questions having to do with the submittal and the
12 presentation today. In the submittal you're asking
13 for approval of two different designs: the CRT
14 which you presented data on today, and which was
15 submitted to us in the Quadra RG. I have seen
16 nothing to support, nothing on this particular
17 design, and I'm wondering what is the basis for us
18 to consider approving this other design.

19 DR. MEYERS: Bill Meyers with Paragon.
20 Corneal refractive therapy or orthokeratology
21 requires a contact lens that has three elements.
22 Those elements are a central base curve, typically
23 spherical; a peripheral region which is designed
24 to, in one way or another, align itself to or be
25 tangent to the cornea, in order to maintain

1 centration; and a connecting zone that connects
2 those two.

3 The difference between these designs is
4 primarily in that connecting zone. The function of
5 the lens as such is similar between the two. And,
6 as I mentioned earlier, we do have an approval for
7 that design in daily wear. We felt that the issues
8 here were issues of safety, not of efficacy, since
9 the efficacy of that other design had been
10 demonstrated and is well known in the industry.

11 So we did not include it in the study, and
12 of course including it in the study would have made
13 the study much more difficult, in that there would
14 have been an additional randomization required.
15 This lens is manufactured, this design was
16 manufactured at Paragon under Paragon control, and
17 we could essentially assure that it was identical
18 one to another, and we wanted to avoid that
19 confusion.

20 DR. HARRIS: As a follow-up--Michael
21 Harris--do you have any data on the Quadra RG
22 material, either as to safety and efficacy, that
23 you would like us to consider, that hasn't been
24 submitted previously?

25 MS. THORNTON: Michael, could you repeat

1 that a little louder into the microphone, please?

2 DR. HARRIS: Yes, ma'am. Thank you.

3 MS. THORNTON: Thank you.

4 DR. HARRIS: Do you have any data
5 supporting the safety and efficacy of the Quadra RG
6 material that you would like us to consider, that
7 has not been previously submitted to the panel?

8 DR. MEYERS: I'm not exactly sure--I'll
9 let Jerry Legerton--do you want to go ahead, Mark?

10 DR. BULLIMORE: Mark Bullimore. In terms
11 of the safety of the Quadra RG, traditionally--and
12 maybe somebody could correct me here if I misstate
13 this--approval for contact lenses has been based on
14 material and indication. As far as the Quadra RG
15 lens that we're asking you to approve, or the
16 Quadra RG design we're asking you to approve as
17 part of this PMA, the material and the indications
18 are identical. We would submit, in the absence of
19 any data, that the safety profile of the Quadra RG
20 design should be considered equivalent to the
21 safety profile of the CRT design, based on the
22 materials and the general underlying principles of
23 the design.

24 As far as the efficacy is concerned, we
25 would again submit that the efficacy has already

1 been demonstrated through the daily wear approval.
2 Now, the panel and the FDA may want to consider
3 subtle differences in the labeling for the two
4 designs, and we would be happy to sort of discuss
5 that further, but I hope that addresses your
6 question.

7 DR. HARRIS: Thank you. The protocol
8 calls for following 150 patients for nine months,
9 yet only 85 completed the nine-month study. What
10 happened to the other some 65 subjects that you
11 were supposed to be completing?

12 DR. LEGERTON: Through the course of the
13 clinical trial--Jerry Legerton. Thank you.
14 Through the course of the clinical trial, we had
15 discussions with FDA staff relative to submission
16 of the PMA, and proceeded with the understanding
17 that a submission of at or about 100 subjects
18 completing a nine-month clinical trial would be
19 acceptable to proceed with the PMA. We appreciate
20 that the initial guidance given by the Ophthalmic
21 Device Panel was for that greater number.

22 DR. HARRIS: Michael Harris again. The
23 data submitted indicates that there were 11 eyes
24 fitted that were outside the eligibility criteria.
25 Were they included in the data analysis?

1 DR. LEGERTON: The 11 subjects that
2 completed but weren't in the efficacy analysis, is
3 that what you're speaking to?

4 DR. HARRIS: You indicate, on 105 of the
5 submittal it says that there are 11 eyes that were
6 outside of the eligibility criteria, departures,
7 and I'm wondering if those are included in the data
8 analysis.

9 DR. LEGERTON: Yes. We have, again, all
10 eyes treated are in the safety analysis, and we
11 have the breakout for the cohorts.

12 DR. HARRIS: And my last question for this
13 particular time: There is some concern about
14 treating adolescents, patients under 18 who are
15 minors. You have indicated that some were treated
16 and eventually not included. Can you discuss with
17 us how you wish us to evaluate the use of this lens
18 for patients under the age of 18?

19 DR. BULLIMORE: Mark Bullimore. It was
20 the intention at the commencement of the trial that
21 the population under study reflect the population
22 that the sponsor believed that the lens would be
23 used upon, and in initial discussions with the FDA
24 the sponsor's plan was to include subjects in the
25 trial of age 12 and above. Initially, the FDA

1 exhibited a preference that the initial cohort not
2 include people below the age of 18, and so the
3 initial cohort of subjects that were recruited were
4 18 and above.

5 It was only subsequent to that that the
6 FDA permitted the enrollment of subjects 12 and
7 above, and it's because of that time line we are
8 unable to provide you today with safety and
9 efficacy data on a complete adolescent cohort.
10 Now, the adolescent subjects are included in the
11 safety cohort and the analysis that Dr. Schein
12 presented, but they are, because none of them have
13 reached the nine-month period, none of them are
14 included in the efficacy cohort.

15 We would hope that the wisdom of the panel
16 and the sound clinical judgment of the members
17 might illuminate some discussion on this issue,
18 both in terms of the inclusion of adolescents in a
19 PMA trial and also, of course, in the labeling
20 indications for use of the product. But that's the
21 history of it, and that's why unfortunately you
22 don't have data to review today.

23 DR. HARRIS: Thank you.

24 DR. WEISS: We're going to have Dr.
25 Edrington, followed by Dr. Matoba and Dr. Bradley.

1 DR. EDRINGTON: Tim Edrington. Relative
2 to the materials, was there any difference in terms
3 of adherence or imprint?

4 DR. LEGERTON: The incidence of imprinting
5 and adherence was reported only in comment sections
6 and very rarely. There was no statistical analysis
7 on the frequency, but it is in fact very low, and
8 the low incidence was in both materials. We
9 actually speculated early on--even the literature
10 supports that that might be hypoxic-related as
11 well--that it may occur more in a high Dk versus a
12 super Dk material, but in fact it appeared to be
13 near equal in the two. In other words, it occurred
14 in both materials. I can't say near equal. It
15 occurred in both materials, low incidence.

16 DR. EDRINGTON: Relative to materials
17 again--Tim Edrington--are there any issues in terms
18 of fabrication, with the laboratories making the
19 lenses, when you look at the two materials?

20 DR. LEGERTON: Jerry Legerton again. Back
21 to the two material issue, and I want to respond--
22 this in part answers a former question.

23 Again, it is our experience, especially in
24 the rigid gas permeable industry, that a material
25 is approved for an indication and that some

1 parameter variance is allowed within that. And in
2 fact manufacturing methodologies and the role of
3 one who has the clearance to market that material
4 is to authorize laboratories to manufacture lenses
5 under that approval, and therein have supervision
6 over Good Manufacturing Practices and their
7 methodologies, in essence to say that we must--we
8 take that responsibility to be sure that they have
9 some design control, or they have design control
10 and that they have the process control to produce
11 what they say they produce. That really is our
12 understanding with this, and that is the purpose of
13 the Quadra RG, that that would--it would be under
14 that name, that independent laboratories would be
15 allowed to produce the very lenses they produce
16 today.

17 We fully appreciate in this particular
18 indication that the precision of the proximity
19 depths, however they achieve them, is a very
20 critical issue, and in applying the approval for
21 the Quadra RG and executing it in the real market,
22 which is where independent laboratories make their
23 lenses, there is a need to also assure that they
24 have the metrology to measure what they say they
25 make. Dr. Edrington, you are right on the mark

1 there, you know, what the panel want, understanding
2 that the world of rigid gas permeable lenses is
3 primarily a world of independent laboratories
4 manufacturing lenses.

5 DR. EDRINGTON: Relative to the CRT
6 design, is there a difference--

7 MS. THORNTON: Dr. Edrington, could you
8 speak into the microphone a little more? It's hard
9 to capture your voice. Sorry.

10 DR. EDRINGTON: Relative to the
11 independent laboratories fabricating the CRT
12 design, is there any reason they would be reluctant
13 to do it in the HDS, Hunter design, I mean
14 material, as opposed to the HDS? I'm trying to get
15 at why the HDS would be utilized at all

16 DR. MEYERS: Bill Meyers from Paragon.
17 Laboratories have preferences. Paragon has
18 demonstrated that any laboratory can manufacture
19 either of these materials. However, there are
20 process differences, and laboratories often have
21 limitations in the number of different processes
22 they like to use. So they have their preferences
23 among materials, and I'm confident there would be
24 some laboratories who would prefer one over the
25 other. But that is not because they cannot do it;

1 it is only because they choose not to for internal
2 reasons.

3 That's one of the reason that we studied
4 both materials, is because we knew there would be
5 such preference. Many of these laboratories
6 manufacture a daily wear lens today which are used
7 by professionals across the country. There would
8 be no extended wear, or rather overnight wear
9 labeling, if we do not approve this material. And
10 yet those practitioners who may already feel
11 comfortable doing that, may work off-label without
12 the benefit of labeling, and that was I think
13 primarily our reason for including.

14 DR. EDRINGTON: From Paragon's
15 perspective, then, assuming that there weren't
16 doctor preferences and laboratory preferences, you
17 would pretty much--would you recommend everybody be
18 putting the HDS 100 material in the CRT design?

19 DR. MEYERS: Not speaking as a clinician,
20 yes.

21 DR. EDRINGTON: Okay. One additional
22 question. Maybe this has nothing to do with the
23 FDA. Will practitioners be required to be trained
24 and certified in fitting of this device?

25 DR. LEGERTON: Jerry Legerton.

1 MS. THORNTON: Dr. Rosenthal?

2 DR. ROSENTHAL: That does--

3 DR. EDRINGTON: Is that outside of the--

4 DR. ROSENTHAL: No, that is not outside.

5 Rosenthal. That is not outside the purview of the
6 FDA. In fact, if the panel feels that it is an
7 appropriate recommendation, we would appreciate the
8 panel's view on that.

9 DR. LEGERTON: And that was my answer,
10 too.

11 DR. EDRINGTON: Thank you.

12 DR. ROSENTHAL: Thank you, Dr. Legerton.

13 DR. LEGERTON: That would be your purview,
14 if that's something you felt would be required in
15 the marketplace.

16 DR. EDRINGTON: Is it within our scope,
17 then, to ask what your recommendation would be or
18 what your plans are?

19 MS. THORNTON: No, not what his
20 recommendation would be.

21 DR. EDRINGTON: Okay. Your opinion? How
22 can I ask his?z

23 DR. ROSENTHAL: Well, I think it really
24 should--this is Rosenthal--I really think it should
25 be up to the panel to make that decision, not up to

1 the sponsor to suggest to the panel what they would
2 like the panel's decision to be.

3 DR. WEISS: Maybe I could ask the sponsor,
4 how difficult is it for a practitioner to learn how
5 to dispense this lens?

6 DR. LEGERTON: Jerry Legerton. There
7 definitely is a learning curve in this modality.
8 Technology can support--I think over time the panel
9 will see probably the application or the desire to
10 apply, let's say, corneal topography in a manner
11 that's different than it has been applied. There
12 might be--or the staff may see software that needs
13 to be validated, that would in fact come up with
14 the final parameters.

15 I'm convinced beyond a doubt that the
16 precision required in determining the final
17 parameters is far greater than the precision
18 required to determine the final parameter of a
19 rigid gas permeable contact lens. I believe that
20 the noise level of the return zone or the proximity
21 control of this lens or any other person's lens
22 doing the same thing is plus or minus 20 microns.
23 So the question is, how do you get within the plus
24 or minus 20 microns?

25 And fluorescein pattern reading, which has

1 been the modo dei, let's say, or the useful, the
2 most used technology for gas permeable lenses, may
3 well have a noise level far greater than 20, plus
4 or minus 20 microns. That would be a reason for
5 the reordering at times. It's not, you're not
6 doing a serial fitting where you're needing two or
7 three lenses, taking it in bits and pieces. But
8 rather the reorders would occur--because this was
9 intended to be a single lens treatment--or did
10 occur in 60 out of 188 eyes, or 60 times for 188
11 eyes, was for the purpose of refining the
12 parameters.

13 In answer to your question, specifically
14 how long would it take, or I believe there is
15 learning that a practitioner who has even been in
16 practice a good length of time and fit a lot of
17 lenses, there is some learning that needs to take
18 place, and that is both from an academic point of
19 view and from an experiential point of view.

20 DR. WEISS: Dr. Edrington, are you
21 finished?

22 Dr. Matoba?

23 DR. MATOBA: In regard to those patients
24 who had the best--transient loss of the best
25 spectacle corrected visual acuity, and particularly

1 the subjects that had worse than 20/40, were they
2 symptomatic at the time this was started? Were
3 they complaining of decrease in vision?

4 DR. LEGERTON: Jerry Legerton. I would
5 have to try to do a correlation between that visual
6 acuity and their symptoms, problems and complaints.
7 It's interested that that's best spectacle
8 corrected visual acuity.

9 DR. MATOBA: I would expect their
10 uncorrected to be worse.

11 DR. LEGERTON: Right. Not necessarily.

12 DR. MATOBA: Okay.

13 DR. LEGERTON: They may be 20/32
14 uncorrected, and put on their spectacles and still
15 be 20/32, and they may be gauging their
16 satisfaction with their vision compared to their
17 pretreatment, like "I was 20/200 before and now I'm
18 20/32."

19 DR. MATOBA: Well, if they're satisfied,
20 that's fine.

21 DR. LEGERTON: So they may not be
22 complaining at all.

23 DR. MATOBA: But if you don't know the
24 answer--

25 DR. LEGERTON: I don't. I have not done

1 that correlation.

2 DR. MATOBA: And so this was transient, so
3 at subsequent visits they no longer had that. Say
4 at six months we have 3 percent, 4 percent of
5 patients having loss of best corrected acuity,
6 that's another 4 percent of patients. So I wonder
7 why someone would be doing well and all of a sudden
8 at six months they develop a loss or decrease of
9 best corrected visual acuity, and also, were there
10 a subset of patients who had recurring episodes,
11 who kept getting every so often, every few months,
12 some--

13 DR. LEGERTON: Jerry Legerton. Yes, in
14 the clinical report we have given you a line
15 listing of eyes that lost two or more lines of
16 acuity and then categorized it, and by that it is
17 possible to look and see were there repeat eyes.
18 And if I'm not mistaken, there was one subject that
19 had repeat loss on three visits of let's say the
20 eight visits. There were a couple of others that
21 had repeats on two visits.

22 But in all cases we were able to verify
23 even those that, let's say, had a loss at a nine-
24 month visit, have been seen subsequently, so we
25 have on every patient the knowledge that on a

1 subsequent visit the loss was not present. But
2 think the answer to your question is yes, there
3 were some subjects that there were repeat measures
4 of loss of acuity.

5 DR. MATOBA: And the first part of that
6 second question was, why at some--late in the
7 study, six months into the study, when you assume
8 they are at a pretty steady state in terms of the
9 wearing time and the adaptation to contact lenses,
10 would a patient suddenly develop best corrected
11 visual acuity worse than, say, 20/40?

12 DR. LEGERTON: Jerry Legerton. I don't
13 recall whether this was Dr. Schein or Dr. Bullimore
14 that I'm quoting, but in this modality it's
15 possible to have a "bad hair day," that on a given
16 night--

17 [Laughter.]

18 DR. LEGERTON: --that on a given night a
19 lens can decenter or have less centration, less
20 than optimum centration. And so it is possible on
21 a single measure to have someone who would have a
22 loss and then would not have at another time.

23 DR. WEISS: Dr. Bradley?

24 DR. BRADLEY: A couple of questions.
25 There was the implication that the primary cause of

1 dropout was due to unsatisfactory vision quality,
2 uncorrected, and the primary reason for that was a
3 fitting failure on the part of the clinician.
4 Another way to assess failure is not in the dropped
5 out subjects but those who actually continued
6 through the nine months.

7 And I'm looking at your achieved versus
8 attempted and there is, as we might imagine, a
9 considerable scatter in the data set, and I'm
10 wondering whether the failures to achieve the
11 attempted, would you claim that that was due to a
12 fitting error, because I think the presenter
13 indicated that might be due to errors in just doing
14 refractions. It seems to me there are two quite
15 different explanations for these types of errors,
16 and I wonder if we could have some clarification on
17 that.

18 DR. BULLIMORE: Mark Bullimore. As is
19 typical in any refractive therapy, you do get some-
20 -you don't get that perfect ratio line and
21 everybody falling about that line. In the order of
22 50 percent of the patients were within half a
23 diopter of attempted correction, and in the order
24 of 90 percent were within 1 diopter of attempted
25 correction. And yes, you're subject to the whims

1 of variability of manifest refraction.

2 But even though patients may be
3 undercorrected by say a diopter or so, it's still
4 possible for that patient to be quite happy with
5 that outcome due to the reduced dependence on other
6 modes of refractive correction. I had dinner with
7 a gentleman the other night who was wearing his
8 lenses on an intermittent basis, and he was trying
9 to find that sweet spot in the middle where he had
10 adequate distant vision and, as a nascent
11 presbyope, had adequate near vision. So I'm not
12 proposing that the label include that sort of
13 intermittent wear, but you don't have to have 20/20
14 vision and emmetropia to be satisfied with this
15 modality.

16 DR. BRADLEY: I have a couple of follow-up
17 question on that. In the refitting that went on,
18 were these patients refitted because of
19 undercorrection? And were they refit to increase
20 the intended correct?

21 DR. LEGERTON: Jerry Legerton. There
22 certainly were a percentage of patients that were
23 refit for undertreatment, and the table says
24 undertreatment, and then poor centration also nets
25 undertreatment. I put back up the scattergram

1 because you were addressing that, as well.

2 DR. BRADLEY: Yes. Thanks.

3 DR. LEGERTON: And if we put sort of a
4 mean line through it, where the mean is, again at
5 nine months it was minus let's say approximately a
6 quarter. One notion is that if one can get a 5
7 diopter correction on a 5 diopter eye, why do we
8 have undertreatments on 2 diopter eyes?

9 DR. BRADLEY: Exactly the point.

10 DR. LEGERTON: Two and a half diopter
11 eyes. And we were fairly rigid in our
12 recommendations in the fitting guide to the
13 investigator, and one piece of rigidity was that we
14 suggested that they or recommended that they
15 calculate the base curve at a half diopter greater
16 than what they wanted the end result to be. If one
17 looked at that and said, "Well, gee, if you're
18 always a half under, why don't you just adjust your
19 nomogram here, your algorithm?" And that in fact
20 is a reality, that if you can get a 5 on a 5 or a
21 450 on a 5, you should be able to get a 450 on a
22 450.

23 DR. BRADLEY: So my final--you've actually
24 made my next two points for me, but the impression
25 I had looking at these data was that you could

1 certainly, the clinician could certainly do better
2 outside of this study protocol, in the sense that
3 the device has been shown to be able to achieve a
4 fairly large refractive error, and with refitting
5 presumably you can tweak the final achieved
6 refractive error. Is that your expectation?

7 DR. LEGERTON: Jerry Legerton. Yes. Yes.
8 We expect that the real world result will exceed
9 the clinic trial, these data in terms of the
10 presentation.

11 DR. BRADLEY: Now I have a completely
12 different topic, really a follow-up from Mike
13 Harris's question about seeking approval for the
14 Quadra lens. The argument seemed to be based upon
15 a couple of equivalences: One, that the underlying
16 theory behind this lens' effect was basically the
17 same as the one that was studied in this PMA. And
18 the second argument was that the device, the Quadra
19 already exists for daily wear.

20 And it seemed to me that the critical
21 comparison for us, then, was for you to convince us
22 that the overnight effects of the study lens are
23 equivalent to the daily wear effects of that same
24 lens, and by implication we could then stretch the
25 argument to say, well, perhaps we expect a similar

1 equivalence for the Quadra lens, and therefore we
2 can accept Dr. Bullimore's suggestion that these
3 lenses are equivalent.

4 DR. LEGERTON: Jerry Legerton. Yes, you
5 have summed that up well. We believe that these
6 data support the safety of this material for this
7 indication, and that parameter variance should be
8 allowed in that, similar to other rigid gas
9 permeable applications. However, since we have not
10 tested all reverse geometry designs for their
11 efficacy, we have in our labeling reduced the
12 efficacy that would be in the labeling to what was
13 established on an eight-hour open eye application
14 instead of an eight-hour closed eye application.

15 DR. BRADLEY: Thank you.

16 DR. WEISS: I had one question. Then I
17 think we'll continue with a question by Dr.
18 Grimmer, and then we'll take a 15-minute break
19 before going on to the FDA presentation.

20 You indicated that many of the cases of
21 corneal edema occurred in the practice which was in
22 high altitude, and so that might be associated with
23 the cause of corneal edema. With that in mind, do
24 you have any suggestions to practitioners as far as
25 a different fit or anything they can be doing to

1 decrease the chance of corneal edema if they're
2 fitting in a high altitude area?

3 DR. LEGERTON: Jerry Legerton. I want to
4 address that question back to a little bit more
5 elaboration on the two materials themselves, some
6 statements that weren't made. Both materials have
7 seven-day approval for extended wear.

8 One of the questions by the reviewers had
9 to do with the overall transmission or average
10 thickness of these lenses, and that's the other
11 thing I wanted to address, that these lenses are
12 all huddled about plano. They have near parallel
13 surfaces, and the thickness of the lenses is about
14 .15 over the whole profile of the lens. That's
15 what we use to support the fact that the actual
16 average transmission of these lenses is greater
17 than the same material used for minus 3 myopes and
18 plus 3 hyperopes in extended wear.

19 Further, I wanted to point out the work of
20 Hill and others in the past, that the equivalent
21 oxygen percentage requirements of humans follow a
22 normal standard distribution, and there are those
23 that have higher requirements and lower
24 requirements. And one of the values we see of
25 having both materials is that they in fact do

1 complement each other, and while a particular
2 practitioner may have a prescribing philosophy that
3 their comfort zone is, "I'm only going to use the
4 super Dk material," another may prefer the lower Dk
5 for other reasons, but they vary accordingly.

6 But I would certainly think that, for
7 instance, the doctor at the 7,000-plus altitude--
8 who happened to have been on the Device Panel for
9 17 years himself, and was involved in scaling,
10 developing the scaling requirements for slit lamp,
11 so I think he called it right--I don't think he
12 would be inclined to use much of the lower Dk
13 material, and would probably recommend to people at
14 high altitude, you know, if we have that one extra
15 compromise, and given these data, pick the high Dk.

16 DR. WEISS: Just a follow-up of that. Of
17 those, I think there were seven cases that occurred
18 at the high altitude or something along that--

19 DR. LEGERTON: There were 17 of 18 eyes
20 with Grade 3--

21 DR. WEISS: Were in the high altitude.
22 Okay, so of those 17, do you know what percentage
23 had the high Dk/A?

24 DR. LEGERTON: They were both materials.

25 DR. WEISS: So there was--

1 DR. LEGERTON: And it's in the line
2 listing. I can get it for you later and answer,
3 but it's in the line listing in your report.

4 DR. WEISS: Well, if you suggest that you
5 use a higher Dk if you're in a high altitude area,
6 and let's say half the patients who developed the
7 problem had high Dk, then it has no significance.

8 DR. LEGERTON: If I may speculate quickly,
9 and maybe one of my colleagues is tabulating, but I
10 would say it was probably--

11 DR. WEISS: Dr. Bullimore.

12 DR. LEGERTON: Mark Bullimore can answer.

13 [Laughter.]

14 DR. BULLIMORE: Mark Bullimore. Again
15 referring to Attachment 5 on page 76 of the
16 submission to the FDA, I believe I'm right in
17 saying that those 17 cases occurred in 8 patients.

18 DR. WEISS: Okay, so I said seven.

19 DR. BULLIMORE: Of those eight patients,
20 so these were the 17 instances or presentations of
21 edema at that high-altitude practice, there were
22 eight patients, of which three were fit with the
23 high Dk material and five with the low, so--

24 DR. WEISS: Okay, so we don't know.

25 DR. HARRIS: What page are we looking at?

1 DR. BULLIMORE: Page 76 in my copy.

2 DR. HARRIS: Which tab?

3 DR. BULLIMORE: I'm not working from the
4 same that you are. It's Attachment 5.

5 DR. HARRIS: Which tab?

6 DR. BULLIMORE: It's Tab 6 in mine, but
7 that may not be relevant to yours.

8 DR. LEGERTON: Mike, it's Jerry Legerton.
9 It's Attachment 5 in the Clinical Report.

10 DR. WEISS: So essentially what we're
11 saying, or what I understand you to say, is we
12 don't really know. You would suggest five patients
13 had the lower Dk but three patients had the high
14 Dk, so they wouldn't have been--we wouldn't be able
15 to offer them any suggestions to decrease the
16 chance of getting corneal edema. So perhaps that
17 will be a labeling issue, just to inform
18 clinicians.

19 Dr. Grimmett next. Oh, Dr. Edrington.

20 DR. EDRINGTON: Tim Edrington. Just a
21 follow-up on that. Were those patients
22 discontinued at that time in that practice, or did
23 the edema go away at follow-up visits, when they
24 saw the patient down the road? Were they all at
25 the first visit, or--

1 DR. LEGERTON: The majority of the reports
2 were day one, and in the majority of the cases the
3 practitioner did not intervene, and in the majority
4 of the cases the edema resolved.

5 DR. WEISS: Dr. McMahon, was it on this
6 question or something else?

7 DR. McMAHON: Related.

8 DR. WEISS: Related? Fine.

9 DR. McMAHON: I was the reviewer that
10 asked for the transmissibility data, and I still
11 haven't heard. Even though the lenses themselves
12 are in form now of being equivalently plano, I
13 haven't heard what the transmissibility--not
14 permeability, transmissibility--of these two
15 different designs are, of the materials.

16 DR. LEGERTON: Of the two different
17 designs?

18 DR. McMAHON: Two materials. Materials.
19 Do you have that?

20 DR. LEGERTON: To get transmissibility you
21 would need the average thickness profile of both
22 designs. You want to know Dk/L?

23 DR. McMAHON: Right. I asked that several
24 weeks ago.

25 DR. LEGERTON: The two materials, yes.

1 Dk/L of paflucocon-B in the Paragon CRT design is
2 40. That's a 60 Dk divided by .15. The
3 transmissibility of the paflucocon-D is 100, a Dk
4 of 150, thickness .15.

5 DR. McMAHON: Thank you.

6 DR. WEISS: Dr. Grimmett, and then we will
7 take a 15-minute break.

8 DR. GRIMMETT: Mike Grimmett. I was
9 pleased to see in the open public hearing that Dr.
10 Rah noted she used the RSVP survey, and I'd like to
11 just acknowledge the superlative work that Dr.
12 Schein has done in done in this area in his
13 development of that instrument. I have a question.
14 Was the RSVP survey used in this PMA, or any
15 quality of life survey of that nature?

16 No. Dr. Bullimore is shaking his head no.
17 Thank you.

18 DR. WEISS: I want to thank the sponsor
19 for a very clear presentation, and we will
20 reconvene in 15 minutes.

21 [Recess.]

22 **FDA PRESENTATION**

23 DR. WEISS: We are now going to begin with
24 the FDA presentation, and Dr. Saviola will begin.

25 DR. SAVIOLA: Thank you, Dr. Weiss. I'm

1 Jim Saviola. I'm the branch chief of the
2 Vitreoretinal and Extraocular Devices Branch, and
3 I'm going to introduce our review team in just a
4 moment, but first I have some brief general
5 introductory remarks about the topic of
6 orthokeratology and also some specific remarks
7 about this application.

8 Back on February 12, 1998, the topic of
9 RGP lenses for overnight orthokeratology was the
10 subject of discussion for guidance purposes at an
11 advisory panel meeting held at the Parklawn
12 Building. We also discussed 30-day extended wear
13 lenses for guidance at that meeting, as well.

14 As membership of the panel changes, most
15 of the advisory members who participated at that
16 meeting have rotated off the panel, with the
17 exception of Dr. Harris.

18 [Laughter.]

19 He will serve as a historical reference
20 today, since he was in attendance. Dr. Bullimore
21 was also present that day, but as you can see, he
22 has crossed over to the other side.

23 [Laughter.]

24 The other side of the table, that is.

25 Later that year, on September 25, 1998,

1 FDA issued a public health notification regarding
2 illegal promotion of contact lenses that addressed
3 both orthokeratology and tinted lenses, and that
4 advisory contained the following statement:

5 "FDA is not aware of any well-controlled
6 clinical studies published in the literature on the
7 overnight use of orthokeratology lenses. The
8 overnight use of lenses is not considered daily
9 wear, and is considered extended wear, since the
10 lens is worn while the user is asleep."

11 Since that time, publications have begun
12 to appear in the literature on the topic of
13 overnight ortho-k. However, the study reported
14 before you today represents the most comprehensive
15 effort undertaken and completed to date on this
16 topic. I would like to commend the sponsor for the
17 quality of the submission, not just from the
18 standpoint of making the panel's job and our job at
19 FDA easier, but for advancing the scientific
20 knowledge base on this topic by way of their
21 clinical study.

22 We have referred this application to you
23 today for your review and recommendations, since
24 this is the first application that we have received
25 that is suitable for panel review for this proposed

1 indication. Current FDA guidance for industry on
2 submission of ortho-k RGP lenses was issued back on
3 April 10th of 2000. That document recommends a
4 study of 150 to 200 completed subjects, followed
5 for 12 months, for an overnight orthokeratology
6 study.

7 We brought this PMA to the advisory panel
8 today with nine-month follow-up data because of the
9 sponsor's ability to demonstrate stability of the
10 intended effect at an early time point of the
11 investigation. In your panel pack under the
12 protocol Section 4.5.2, the original protocol did
13 stipulate that it was a nine-month study and they
14 were targeting 150 to 200 subjects to be completed.
15 Excuse me. They were targeting 150 subjects
16 complete at nine months with an enrollment of 200,
17 so they anticipated a discontinuation rate of 25
18 percent, and as you can see from the data, that was
19 overly optimistic.

20 However, in the back of your panel packs,
21 the Berkeley study that the sponsor made reference
22 to earlier in their presentation does reflect a
23 certain consistency with the outcome achieved in
24 this study, from the standpoint of successfully
25 completed patients. It seems from current

1 literature, as well as in this PMA dated today,
2 that the discontinuation rate is going to be much
3 higher than 25 percent.

4 We do not expect everybody who enters into
5 this therapy to be successful. We do expect
6 limitations with this procedure, and indeed over 50
7 of the 72 pages of our orthokeratology guidance
8 comprise sample labeling, to inform those
9 candidates and practitioners of the clinical
10 outcomes, to give them some idea of what the
11 expectations might be.

12 I would like to make an important point
13 for the public to consider. I would like to
14 caution those observers who leave this meeting
15 today with the idea that the panel recommendation
16 applies to all lenses and all RGP materials, from a
17 regulatory standpoint, FDA is considering overnight
18 orthokeratology for the materials and the lens
19 designs contained in this application only. I
20 would like to emphasize the point for the
21 professional community and the public, that today's
22 outcome is not a blanket general recommendation for
23 all orthokeratology designs and all RGP materials.

24 I would like to thank the FDA review team,
25 both the primary panel reviewers and the internal

1 reviewers. I would now like to introduce Ms.
2 Eleanor McGhee, the lead reviewer for this
3 application.

4 MS. MCGHEE: Good morning. As Dr. Saviola
5 said, my name is Eleanor McGhee, and I am the team
6 leader for this application. I would like to take
7 the opportunity to introduce the team members,
8 those individuals who played a role in getting this
9 application to the panel:

10 Dr. Edrington, Dr. McMahon, served as
11 primary reviewers for the application. Dr. Bernard
12 Lepri conducted an in-house clinical review. Dr.
13 Jean Hilmantle for statistical review. Dr. Jimmy
14 Chin for chemistry and manufacturing. Mr. Lev
15 Keely for panel coordination. Ms. Laura Mendelsohn
16 for patient labeling review. Sally Thornton,
17 executive secretary of the panel, for helping us
18 with our coordinating. And Jim Saviola, our branch
19 chief. To all of you, I thank you very much for
20 your timely responses.

21 I'd like to now turn it over to Bernard
22 Lepri for clinical review.

23 DR. LEPRI: Good morning, panel members,
24 sponsors, FDA members, and other guests. This
25 morning I am going to present to you a very brief

1 presentation of some clarification information and
2 the questions, which should all be completed within
3 under 15 minutes, and you'll notice that I'm
4 sitting this one out today because I'm suffering
5 from PTPD, post traumatic PMA disorder.

6 [Laughter.]

7 So I now present to you PMA P870024,
8 Supplement 43, Paragon Vision Sciences, Contact
9 Lens Corneal Refraction Therapy. The sponsor is
10 presenting two materials and two designs. They are
11 rigid gas permeable contact lenses for overnight
12 contact lens corneal refractive therapy of myopia
13 with or without astigmatism.

14 You know who I am. This is a refresher of
15 some of the materials, of the nature of the
16 materials presented to you today. I'd like to
17 remind the panel that both the materials are
18 approved in this PMA, have been approved for seven-
19 day extended overnight wear. The Quadra RG
20 trademark design in both materials has previously
21 been cleared by FDA for daily wear
22 orthokeratology/corneal refractive therapy, not
23 conventional daily wear of RGP.

24 As some further background information for
25 question number one, the CRT design was the only

1 one studied overnight with both materials in this
2 PMA, and the primary safety and effectiveness
3 endpoints were all met.

4 Question number one: Do the data reported
5 for the two different generic lens materials
6 evaluated during this study raise any questions of
7 safety and effectiveness?

8 Background information for question number
9 two. Both materials were evaluated in overnight
10 clinical trials, and approved. Only one of the
11 designs was evaluated in the trials reported in
12 this PMA. Efficacy of the Quadra RG design is
13 based on prior FDA daily wear clearance for
14 orthokeratology/corneal refractive therapy.

15 Also, I want to direct the panel's
16 attention to the draft labeling in this PMA, where
17 you would find the outcome results from the prior
18 clearance and review of FDA for the Quadra RG
19 design, so those data are available in the draft
20 labeling in the PMA.

21 Question number two: Do the data reported
22 for the CRT reverse geometry lens design evaluated
23 during the study raise any questions of safety and
24 effectiveness?

25 Question number three: Is the length of

1 follow-up sufficient to demonstrate the stability
2 of the intended myopic reduction with the
3 prescribed maintenance regimen?

4 Question number four: What are the
5 panel's recommendations for the proposed product
6 labeling, for example, warnings, precautions,
7 terminology to describe the procedure, etcetera?

8 Question number five: Does the panel see
9 any issues that suggest a need for a post-approval
10 follow-up of the study subjects or a post-approval
11 study?

12 And question number six: Do the data
13 presented in this PMA provide reasonable assurance
14 of safety and effectiveness for the proposed
15 indications?

16 Thank you.

17 DR. WEISS: Dr. Saviola?

18 DR. SAVIOLA: Before we retreat from the
19 table, I just want to make sure that you folks
20 don't have any questions regarding this issue of
21 the two designs that we might need to clarify for
22 you before you enter in your discussion.

23 DR. WEISS: Yes, we will entertain those
24 questions or any other questions the panel has for
25 FDA at this time.

1 There are no questions for FDA? Well, in
2 the spirit of moving ahead, it seems that because
3 this is such a well put together proposal by the
4 sponsor, we are just going to forge ahead with
5 additional comments for the sponsor. And I have
6 had some interest, and I have accepted the interest
7 on the panel's behalf, of skipping lunch and then
8 moving ahead with the committee's deliberation,
9 because it seems that this is going ahead at a good
10 clip.

11 So are there any additional comments from
12 the sponsor at this point?

13 There appear to be no additional comments
14 of the sponsor. Lunch being so skipped, we will
15 then proceed with committee deliberations and have
16 the reviews of the primary panel reviewers, first
17 beginning with Dr. Timothy McMahon.

18 **COMMITTEE DELIBERATIONS**

19 DR. McMAHON: Thank you. While we are
20 getting set up, let me preface my comments by
21 acknowledging my appreciation of the sponsor for
22 two things: one, his putting together a superior
23 proposal, certainly relative to what we experienced
24 yesterday. And the second is for being the first
25 out of the gate to deal with this aspect of

1 overnight orthokeratology in a meaningful and
2 scientific manner, and taking overnight corneal
3 reformation for refractive error from the cult
4 arena to the arena of science.

5 So this pertains to the PMA as described,
6 and I will limit my comments to issues and
7 questions that are raised in my review, and
8 actually won't go over any overall detail, but
9 specific issues.

10 As has been previously mentioned, there
11 really is no data provided with regard to the
12 Quadra design, so it really makes it difficult for
13 the panel to make any or certainly for me to make
14 any comments specific to that. At the same time,
15 it is probably reasonably similar. Certainly the
16 materials employed have been previously approved
17 for overnight wear, and the designs have in large
18 context reasonable similarities. However, there is
19 no data.

20 The issue with regard to transmissibility
21 has been provided just recently. I think a table
22 in the labeling would be useful, particularly so
23 that investigators know whether they're dealing
24 with lenses that adhere to the Holden, Mertz
25 criteria or not. The HDS-100 does. The HDS does