

sgg]

1

SG

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

GENERAL AND PLASTIC SURGERY DEVICES PANEL
OF THE MEDICAL DEVICES ADVISORY COMMITTEE

OPEN SESSION

60th Meeting

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

Monday, July 8, 2002

1:20 p.m.

Gaithersburg Holiday Inn
Two Montgomery Village Avenue
Gaithersburg, Maryland

MILLER REPORTING COMPANY, INC.
735 8th Street, S.E.
Washington, D.C. 20003-2802
(202) 546-6666

PARTICIPANTS

Thomas V. Whalen, M.D., Chairman

David Krause, Ph.D., Executive Secretary

VOTING MEMBERS:

Phyllis Chang, M.D.

Michael A. Choti, M.D.

David L. DeMets, Ph.D.

Robert L. McCauley, M.D.

Michael J. Miller, M.D.

TEMPORARY VOTING MEMBERS:

Nancy A. Dubler, LLB

Amy E. Newburger, M.D.

CONSUMER REPRESENTATIVE:

LeeLee Doyle, Ph.D.

INDUSTRY REPRESENTATIVE:

Debera M. Brown

C O N T E N T S

Conflict of Interest Statement, David Krause, Ph.D	5
Introductions	8
Panel Update, Mr. Anthony Watson	10
Classification of Silicon Sheeting for Scar Management:	
Industry Presentation:	
Mr. Mark E. Dillon, Bio Med Sciences	15
Mr. Tom Fallon, ReJuveness Pharmaceuticals	18
FDA Presentation, Sam Arepalli, Ph.D.	26
Panel Discussion	30
Classification Questionnaire and Vote	53
Reclassification of Absorbable Hemostatic Agents and Dressings:	
Industry Presentation:	
John D. Paulson, Ph.D., Johnson & Johnson	82
Ms. Ronnemoes Bobak, Ferrosan A/S	93
Ms. Judith E. O'Grady, Integra LifeSciences	100
FDA Presentation, David Krause, Ph.D.	111
Panel Discussion	121
Open Public Comment:	
Mr. Michael Ivey	152
Reclassification Questionnaire and Vote	157

P R O C E E D I N G S

Call to Order and Conflict of Interest

DR. KRAUSE: I think we have reached critical mass so we can start the open session of the panel meeting. Good afternoon, everyone. We are ready to begin the 60th meeting of the General and Plastic Surgery Devices Panel.

I am David Krause and I am the executive secretary of this panel and also a reviewer in the Plastic and Reconstructive Surgery Devices Branch, in the Division of General and Restorative and Neurological Devices.

I would like to remind everyone that you are requested to please sign in on the attendance sheets, which are available at the tables just outside the door. You may also pick up an agenda, panel meeting roster and information about today's meeting at those tables. The information includes how to find out about future meeting dates through the advisory panel phone line and how to obtain meeting minutes or transcripts. This and other panel meeting information, including panel meeting summaries and transcripts, are now also available on the worldwide web. Advisory panel meeting activities are available by clicking on the CDRH

1 home page from the FDA website, which is
2 www.FDA.gov. By clicking on premarket issues and
3 then advisory committees, the summaries,
4 transcripts and other advisory committee material
5 section may be accessed. You can then access the
6 CDRH advisory committee database.

7 Before turning this meeting over to our
8 Chairman, Dr. Whalen, I am required to read two
9 statements into the record. First I will read the
10 conflict of interest statement into the record:

11 The following announcement addresses
12 conflict of interest issues associated with this
13 meeting, and is made part of the record to preclude
14 even the appearance of an impropriety. To
15 determine if any conflict of interest existed, the
16 agency reviewed the submitted agenda and all
17 financial interests reported by the committee
18 participants. The conflict of interest statutes
19 prohibit special government employees from
20 participating in matters that could affect their or
21 their employers' financial interests. However, the
22 agency has determined that participation of certain
23 members and consultants, the need for whose
24 services outweighs the potential conflict of
25 interest involved, is in the best interest of the

1 government.

2 Therefore, waivers have been granted for
3 Drs. Michael Choti and Michael Miller for their
4 financial interests in and firms at issue that
5 could potentially be affected by the panel's
6 recommendations. The waivers allow these
7 individuals to participate fully in today's
8 deliberations. Copies of these waivers may be
9 obtained from the agency's Freedom of Information
10 Office, Room 12A-15 of the Parklawn Building.

11 We would like to note for the record that
12 the agency took into consideration certain matters
13 regarding Drs. Choti and McCauley. These panelists
14 reported current interests in firms at issue but in
15 matters that are not related to today's agenda.
16 The agency has determined, therefore, that they may
17 participate fully in all discussions.

18 In the event that the discussions involve
19 any other products or firms not already on the
20 agenda for which an FDA participant has a financial
21 interest, the participant should excuse him or
22 herself from such involvement and the exclusion
23 will be noted for the record.

24 With respect to all other participants, we
25 ask in the interest of fairness that all persons

1 making statements or presentations disclose any
2 current or previous financial involvement with any
3 firm whose products they may wish to comment upon.
4 Thank you.

5 The second statement I am going to read
6 into the record is the temporary voting memo. This
7 is a memo that is signed by Dr. Feigal who is the
8 Director of the Center for Devices and Radiological
9 Health:

10 Pursuant to the authority granted under
11 the Medical Devices Advisory Committee Charter,
12 dated October 27, 1990 and as amended August 18,
13 1999, I appoint Nancy Dubler and Amy Newburger as
14 voting members of the General and Plastic Surgery
15 Devices Panel for this meeting, on July 8 and July
16 9, 2002.

17 For the record, these individuals are
18 special government employees and consultants to
19 this panel or other panels under the Medical
20 Devices Advisory Committee. They have undergone
21 the customary conflict of interest review, and have
22 reviewed the material to be considered at this
23 meeting.

24 At this time, I would like to turn the
25 meeting over to our Chairman, Dr. Tom Whalen.

1 DR. WHALEN: Thank you, Dr. Krause. Good
2 afternoon. My name is Dr. Thomas V. Whalen. I am
3 the chairperson of the General and Plastic Surgery
4 Devices Panel.

5 Today the panel will be making
6 recommendations to the Food and Drug Administration
7 on the classification of silicone elastomer for
8 scar management devices and on the proposed
9 reclassification of absorbable hemostatic agents
10 and dressings from Class III to Class II. I would
11 like to note for the record that voting members
12 present constitute a quorum, as required by 21 CFR
13 Part 14.

14 Before we begin this meeting, I would like
15 to ask our distinguished panel members, who are
16 generously giving their time to help the FDA in the
17 matters being discussed today, and the other FDA
18 staff seated at the head table to introduce
19 themselves. I would ask that each state their
20 names, affiliations and positions and area of
21 expertise, starting to my right with Dr. Witten,
22 please.

23 **Introductions**

24 DR. WITTEN: I am Dr. Celia Witten,
25 division director of the Division of General and

1 Restorative and Neurological Devices at FDA, which
2 is the reviewing Division for these products.

3 DR. DEMETS: I am David DeMets. I am
4 professor and chair of the Department of
5 Biostatistics and Medical Informatics at the
6 University of Wisconsin, in Madison. I am a
7 statistician by degree and have been involved in
8 clinical trials for a long time.

9 DR. CHANG: I am Phyllis Chang, associate
10 professor in the Division of Plastic Surgery and
11 also in the Division of Hand and Microsurgery for
12 the Departments of Surgery and Orthopaedic Surgery
13 at the University of Iowa. I am an FDA panel
14 member.

15 DR. MILLER: I am Michael J. Miller. I am
16 an associate professor of Plastic Surgery at the
17 University of Texas, MD Anderson Cancer Center.

18 DR. NEWBURGER: I am Amy Newburger. I am
19 a dermatologist in New York, in private practice,
20 and I am an attending physician at White Plains
21 Hospital Medical Center, and I teach at St. Luke's
22 Roosevelt Medical Consortium.

23 DR. KRAUSE: I am Dave Krause.

24 DR. CHOTI: I am Michael Choti, associate
25 professor of surgery at Johns Hopkins University in

1 Baltimore, Maryland, and I am a general surgeon and
2 surgical oncologist.

3 DR. DUBLER: I am Nancy Dubler. I am
4 trained as an attorney. I direct the Division of
5 Bioethics at Montefiore Medical Center, and I am a
6 professor of epidemiology and social medicine at
7 the Albert Einstein College of Medicine.

8 DR. MCCAULEY: Robert McCauley, professor
9 of surgery and pediatrics at the University of
10 Texas Medical Branch, and chief of plastic surgery
11 services for the Shriner's Burn Hospital.

12 DR. DOYLE: I am LeeLee Doyle. I am a
13 professor of obstetrics and gynecology, and
14 associate dean for continuing medical education and
15 faculty development at the University of Arkansas
16 for Medical Sciences, College of Medicine, and I am
17 the consumer representative on the panel.

18 MS. BROWN: I am Debera Brown. I am the
19 vice president of regulatory affairs for Broncus
20 Technologies, which is a medical device company. I
21 am also the industry rep on this panel.

22 DR. WHALEN: As stated, my name is Dr.
23 Thomas Whalen. I am chief of the Division of
24 Pediatric Surgery and professor of surgery and
25 pediatrics at Robert Wood Johnson Medical School in

1 New Brunswick, New Jersey.

2 Before we continue with the classification
3 and reclassification portion of the hearing, we
4 will have Mr. Anthony Watson, acting branch chief
5 of the Plastic and Reconstructive Surgery Devices
6 Branch, provide an update on general and plastic
7 surgery device activities since the last meeting.
8 Mr. Watson?

9 **Panel Update**

10 MR. WATSON: Thank you, Dr. Whalen, and
11 good afternoon. I am Anthony Watson, the acting
12 branch chief of the Plastic and Reconstructive
13 Surgery Devices Branch at FDA. Welcome, members of
14 the panel, members of the public and manufacturers
15 to this two-day meeting of the General and Plastic
16 Surgery Panel.

17 This panel last met on July 17, 2001 and
18 recommended approval of Ortec's PMA application for
19 OrCel Bilayered Cellular Matrix for use on donor
20 sites on burn patients. The agency approved this
21 product on August 31, 2001.

22 On November 19, 2001, the agency approved
23 a PMA application for Lifecore's Intergel Adhesion
24 Prevention Solution. This application was reviewed
25 by this panel at the January, 2000 panel meeting

1 and was recommended to be not approvable. The
2 agency agreed and, after receiving a not approvable
3 decision, the sponsor requested review at the newly
4 formed Medical Device Dispute Resolution Panel.
5 This panel met on September 6, 2001 and recommended
6 that the application be approved.

7 On June 18, 2002, the agency released an
8 updated guidance document, entitled, "Guidance for
9 Resorbable Adhesion Barrier Devices for Use in
10 Abdominal and/or Pelvic Surgery."

11 Today, you will make classification
12 recommendations to the agency on two types of
13 medical devices: the silicone elastomer for scar
14 management and the absorbable hemostatic agent and
15 dressing intended for hemostasis during surgical
16 procedures. Tomorrow the panel will be presented
17 with an update of the conditions for approval for
18 the two saline-filled breast implants approved in
19 May of 2000. As a reminder, tomorrow we will not
20 be discussing silicone gel-filled breast implants,
21 and I request that panel members and members of the
22 public limit their comments to saline-filled breast
23 implants.

24 Panel members, we appreciate your
25 commitment. Members of the public who have

1 requested time to address the panel, we appreciate
2 your comments. Manufacturers, we appreciate your
3 participation in presenting the information you
4 have to the panel and answering questions that the
5 panel may have. Thank you for your attention.

6 DR. WHALEN: Thank you, Mr. Watson. We
7 will now proceed with the first open public hearing
8 session of this day. I would ask at this time that
9 any and all persons addressing the panel, please
10 come forward, speak clearly into the podium
11 microphone as the transcriptionist is dependent
12 upon this to provide an accurate record of this
13 meeting.

14 We are requesting that all persons who
15 make statements to the panel during the open public
16 hearing portion of the meeting disclose whether or
17 not they have financial interests in any medical
18 device company whatsoever. Before making your
19 presentation to the panel, in addition to stating
20 your name and affiliation, please state the nature
21 of your financial interests and if you have none,
22 please so state. Is there anyone who wishes to
23 address the panel? Please indicate by show of
24 hands.

25 Since there are no requests to speak in

1 the open public hearing, we will now proceed to the
2 open committee discussion. At this time, we will
3 begin the discussion with the classification of
4 silicone elastomer for scar management. We will
5 start with the presentation by Mr. Mark Dillon,
6 president of Bio Med Sciences. This will be
7 followed by a presentation by Mr. Carey Rehder,
8 plastic reconstruction division engineering manager
9 of PMT Corporation, who will be followed by Mr. Tom
10 Fallon and Mr. Mike O'Brien, of ReJuveness
11 Pharmaceuticals.

12 The FDA presentation and a reading of the
13 FDA questions will follow the industry
14 presentations. We will then have a general panel
15 discussion of this topic, followed by a more
16 focused panel discussion aimed at answering FDA's
17 questions.

18 Following the panel discussion, we will
19 complete the reclassification worksheet and
20 supplemental worksheet. The vote on these
21 worksheets will constitute the panel's
22 recommendation to the FDA.

23 I would like to remind public observers at
24 this meeting that while this portion of the meeting
25 is open for public observation, public attendees

1 may not participate except at the specific request
2 of a panel member. If any of the industry
3 representatives addressing the panel have copies of
4 the remarks that they are making to us today, it
5 would be greatly appreciated if they could pass
6 them to the transcriptionist so that accuracy can
7 be assured in what you are bringing to us today.
8 We will begin with Mr. Dillon's presentation.

9 **Classification of Silicon Sheeting for Scar**
10 **Management Industry Presentation**

11 MR. DILLON: Thank you very much, Dr.
12 Whalen. I am Mark Dillon, the president and
13 founder of Bio Med Sciences. We have been
14 marketing silicone-based products for scar
15 management since the early 1990's.

16 As we are all aware, these are products
17 that are used for the prevention and reduction of
18 hypertrophic scars and keloids. It is my opinion
19 that these devices have substantial importance in
20 preventing impairment of human health and present a
21 potential risk of illness or injury if misused. I
22 think it is common knowledge that some devices are
23 intended for lay use instead of use by healthcare
24 professionals.

25 I, therefore, believe that these products

1 should be classified as Class I, reserved or
2 non-exempt, therefore, requiring a 510(k)
3 notification. I have several reasons for this
4 position. First, there is a wide variety of
5 devices that are on the market. There are rigid,
6 non-adhesive silicone elastomer materials and these
7 generally require some type of tape to hold them in
8 place. There are also adhesive gel type products.
9 Some of these contain other materials as an
10 embedded mesh or some type of reinforcing
11 mechanism. There are past products that are
12 essentially massaged onto the surface of the scar.
13 There is even one product that I am aware of that
14 is a silicone gel-filled cushion that is indicated
15 for this purpose. There are also mineral oil-based
16 materials that are silicone-containing, as well as
17 products that are called tri-block copolymer
18 compounds. In addition, there are a number of
19 composite type structures such as splinting
20 materials that are lined with silicone, padding
21 type products and even textiles that are laminated
22 to silicone.

23 I think that there are likely to be new
24 designs and new products that are introduced to the
25 market, and I think it would be difficult to show

1 substantial equivalence without having some type of
2 review process involved with that.

3 Another consideration is the indications
4 for use. Some of these products are marketed
5 strictly for cosmetic purposes, but others are
6 marketed more for a professional audience, for use
7 with burn patients. Functionality and the
8 patient's health is a critical issue. Furthermore,
9 I think that there has been a wide variety of
10 claims that have appeared in the marketplace with
11 these types of products. I have seen over the
12 years products that claim to heal scars or are even
13 positioned as an alternative to surgery. Likewise,
14 I think these claims should be confirmed through
15 the 510(k) process.

16 Additionally, I believe there are some
17 risks involved with the use of these products. I
18 think patients need to be adequately warned not to
19 use these products on open wounds. Also, there is
20 a possibility of skin irritation or rash,
21 particularly with some of the products that require
22 use of adhesive tape or contain other materials,
23 other than silicone. Lastly, I think that some of
24 these products can be positioned to discourage
25 adequate professional supervision or compliance.

1 Therefore, my concern is that without the
2 premarket notification system some devices may
3 emerge in the marketplace that are not
4 substantially equivalent, are positioned with
5 inappropriate indications and claims, and may pose
6 undue risk, including the discouragement of
7 professional supervision when appropriate. Thank
8 you.

9 DR. WHALEN: Does any panel member have a
10 question for Mr. Dillon?

11 [No response]

12 Thank you, sir.

13 MR. DILLON: Thank you very much.

14 DR. WHALEN: We will now continue with Mr.
15 Rehder's presentation, if Mr. Rehder is available.

16 [Mr. Rehder is not present]

17 Very well, the final identified industry
18 speakers today jointly are Mr. Fallon and Mr.
19 O'Brien.

20 MR. FALLON: Hello, panel. Thank you for
21 letting me speak today. My name is Tom Fallon. I
22 am president of ReJuveness Pharmaceuticals. We
23 market a silicone sheeting product for hypertrophic
24 and keloid scarring.

25 I would like to address the proposed

1 regulation identification in two parts. A scar
2 management device is a silicone sheeting product
3 intended for use on uncompromised skin for scar
4 management.

5 The first part--a scar management device
6 is a silicone sheeting--we fully concur with this
7 identification. Silicone should not be treated as
8 some homogeneous category. Only silicone sheeting
9 has been demonstrated to be effective on problem
10 scarring. The oil and liquid forms of silicone
11 have never been shown, in any peer-reviewed study
12 that I know of, to be effective and are potentially
13 toxic. The difference between the two, as we see
14 it, is that the silicone sheetings give off silicon
15 when hydrolyzed. The silicone oils do not. They
16 give off just silicone.

17 We fully agree with the first part; we
18 fully disagree with the second part of the
19 identification--intended for use on the
20 uncompromised skin. The skin covering keloids and
21 hypertrophic scars seems to be compromised in every
22 way but appearance. The FDA's position is that the
23 skin is not compromised because it is visually
24 intact. We admit that it seems ironic but the
25 functional measures of the stratum corneum covering

1 these scars have been demonstrated to be
2 compromised in three ways.

3 In a study of the functional analysis of
4 the stratum corneum in scars, which I have included
5 in the package that I recommend you read, it was
6 shown that these problem scars yield the same
7 measures as open blister wounds in the categories
8 of transepidermal water loss, electric conductivity
9 and stratum corneum turnover rates. They are four
10 times higher in keloid and hypertrophic scars as
11 they are in atrophic scars and normal skin. Since
12 scar management refers to hypertrophic and keloidal
13 scars and not to atrophic scars, we would have to
14 conclude that scar management refers to compromised
15 skin.

16 I also include a couple of papers by Dr.
17 Peter Elias and his group, out at the University of
18 San Francisco Veterans Administration Hospital. It
19 demonstrates the theories emerging centering around
20 the driving function of the stratum corneum and
21 many maladies of the skin once thought to be
22 originating in the dermis.

23 In "The Mystery Widens" he applies the
24 skin-drive principle to hypertrophic scars and
25 keloids. Another paper included is the

1 "Investigation of the Keloid-Derived Keratinocytes
2 on Fibroblast Growth." It demonstrates that the
3 production of collagen in keloidal and hypertrophic
4 scars is caused by the compromised skin covering
5 them.

6 Our proposed mechanism of action is taken
7 from a paper, "Hypertrophic Scars and Keloids:
8 Immunophenotypic Features and Silicone Sheets to
9 Prevent Recurrences." In this study they took 20
10 keloidal scars, excised them and in ten of them
11 they put silicone sheeting over them; in the other
12 ten they put nothing. In nine of the ten of the
13 scars without the silicone sheeting keloid scars
14 came back. Almost all of them, except four, under
15 the silicone sheeting the scars came back.

16 They did immunophenotypic analysis and
17 they found the scavenger receptor CD36 in large
18 amounts under the silicone sheeting. These
19 scavenger receptors are essential in rebuilding the
20 stratum corneum. The most important component of
21 the stratum corneum is cholesterol. It is not
22 effective when applied topically, and it is
23 transported by these scavenger proteins from high
24 density lipid proteins to the stratum corneum. So,
25 with the proposed identification it will be quite

1 difficult to make the correct structure, function
2 and mechanism of action claims.

3 The last paper that I have included is
4 release and distribution of silicone-related
5 compounds in the skin exhibits the release of
6 silicon from silicone sheeting into a buffer
7 solution and into normal and keloidal skin. If
8 silicon is the active ingredient, then there are
9 dosage and shelf-life issues.

10 We did our own study, which is also
11 included in the packet, where we put ReJuveness
12 Spenco gel sheeting in ten-year old Cica-Care under
13 the same testing, and what we found was that the
14 rates of silicon release were different for
15 different sheetings, and in the ten-year old
16 Cica-Care there was no silicon released at all.

17 In conclusion, we feel as though the
18 silicone sheeting is completely safe and that it
19 should be a Class I but that the scars that it is
20 addressing, hypertrophic and keloid scars, are
21 composed of compromised skin and it is the driving
22 mechanism in these maladies. That is it.

23 DR. WHALEN: Are there questions for Mr.
24 Fallon?

25 MS. BROWN: I have a question. You

1 proposed Class I. Would that be with or without a
2 510(k)?

3 MR. FALLON: I would say with a 510(k).

4 MS. BROWN: Thank you.

5 DR. WHALEN: With your objection to the
6 second part to what has been proposed, is your
7 statement in the center of your third page what you
8 are proposing as an alternative wording, "scar
9 management device...?"

10 MR. FALLON: I really didn't know if we
11 were going to participate on that level to suggest
12 what the wording should be until yesterday. So,
13 yes, I think that it should be changed. I mean, if
14 you would like a suggestion from me, I just need
15 probably a day or two where I could come up with a
16 suggestion.

17 DR. WHALEN: But your viewpoint or your
18 company's viewpoint is that focally hypertrophic
19 scars and keloids are not uncompromised skin. You
20 are not saying to the panel that you think we
21 should consider applying this on open wounds, fresh
22 wounds in the operating room when we have just made
23 an incision, etc.?

24 MR. FALLON: Well, I really don't see why
25 not; I don't see why they shouldn't be applied to

1 open wounds. In the Italian study, where they put
2 it over excised keloids, they do put it on
3 excisions that were open. They did it
4 prophylactically. So, I they are safe enough, yes,
5 to put on open wounds and perhaps that would be a
6 different classification for that use.

7 DR. WHALEN: Any other questions? Dr.
8 Miller?

9 DR. MILLER: Thank you for your
10 presentation. I just want to make sure I am clear
11 about what you are calling an open wound. I mean,
12 when I think of an open wound I think of a wound
13 where the epithelium is not in contact across the
14 wound; you have exposed tissue below the epithelial
15 level that is exposed.

16 MR. FALLON: Right.

17 DR. MILLER: So, your thought is that it
18 is okay. You would suggest that we can place these
19 devices on those types of wounds?

20 MR. FALLON: If they are properly
21 sterilized, yes. I mean, I don't know what effect
22 they would have on open wounds. I know they work
23 on ulcerated wounds, and keloid and hypertrophic
24 scars are very similar to an ulcerated wound in
25 that they are microvascularly cut off. Keloid and

1 hypertrophic scars are composed of essential fats,
2 basically fats. So, the difference between the two
3 is slight.

4 DR. WHALEN: Dr. Newburger?

5 DR. NEWBURGER: Excuse me, have you see
6 any evidence of any type of foreign body reaction
7 from the silicon which is released from the gel
8 across this compromised epidermis?

9 MR. FALLON: Yes, we have had a couple of
10 reports of women using the sheeting applied to the
11 scars after breast implantation. I don't know if
12 there is something going on between them, but we
13 have had several complaints on that. But, for the
14 most part, we have sold over 100,000 of these
15 devices and we have had basically no complaints,
16 just the tape occasionally.

17 DR. NEWBURGER: I am asking specifically
18 about a foreign body reaction as opposed to an
19 irritation or a folliculitis and occlusion. In
20 other words, here is a molecule that is going
21 through the compromised epidermi, are you getting a
22 soft tissue reaction in the dermis with foreign
23 body cells?

24 MR. FALLON: I don't know that exactly.
25 Our scientific advisor is Dr. Arthur Brawer, who is

1 a noted expert on silicone. But it basically is
2 like coal miner's disease, that is, the action of
3 the silicone sheeting on the hypertrophic scar. It
4 goes down as an antigen, stimulates the scavenger
5 CD36 and marshals them to the site, it seems, and
6 then from there they are able to produce and serve
7 their many different functions, versatile functions
8 that they are able to do--transporting cholesterol,
9 essential fatty acids, as well as taking away
10 excess materials in the extracellular matrix. So,
11 that is what we think is going on.

12 DR. WHALEN: Mr. Fallon, are your remarks
13 everything or is Mr. O'Brien still going to be
14 speaking?

15 MR. FALLON: Oh, Mr. O'Brien couldn't
16 show, I am sorry.

17 DR. WHALEN: Thank you. We will continue
18 now with the FDA's presentation with Dr. Sam
19 Arepalli.

20 **FDA Presentation**

21 DR. AREPALLI: Good afternoon. We are
22 here this afternoon to seek a panel recommendation
23 to classify scar management devices indicated for
24 management or scars. My name is Sam Arepalli,
25 reviewer in Plastic and Reconstructive Surgery

1 Branch, Division of General, Restorative and
2 Neurological Devices. I will be presenting device
3 identification and health risks aspects of the
4 device. Reviewers from the Office of Surveillance
5 and Biometrics, CDRH are in the audience to clarify
6 any questions regarding Medical Device Reports.
7 After my presentation, Ms. Marjorie Shulman will
8 walk you through the classification worksheets.

9 This slide is on regulatory history. As
10 you know, medical devices are classified into three
11 classes, namely, Class I, Class II and Class III.
12 Examples of Class I exempt products include
13 hydrogels or hydrogel wound dressings and manual
14 surgical instruments. Class II devices include
15 implantable surgical meshes and sutures. Examples
16 of Class III devices are interactive wound
17 dressings and barriers.

18 At the time of the Medical Device
19 Amendments of 1976, a few medical devices were
20 unclassified. They include devices like scar
21 management devices, the one that we are going to
22 discuss today. They were unclassified. These
23 devices are currently regulated as unclassified
24 devices. The FDA has been making efforts to
25 classify and reclassify medical devices since 1976

1 into the lowest regulatory class that can
2 reasonably assure safety and effectiveness of their
3 intended use.

4 I would like to bring to your attention
5 that the same panel several years ago provided a
6 recommendation to classify non-interactive wound
7 dressings. This slide gives a brief description of
8 the proposed identification of the device: A scar
9 management device is a silicone sheeting product
10 intended for use on uncompromised skin for scar
11 management.

12 This slide gives a brief description of
13 FDA-cleared scar management devices. FDA has
14 regulated silicone sheeting intended for scar
15 management as an unclassified pre-amendment device.
16 It has been cleared for marketing under several
17 names. They are silicone sheeting, silicone
18 elastomer and silicone gel for hypertrophic and
19 keloid scar management. Also, the agency cleared a
20 hydrogel for the same intended use.

21 There are about 75 scar management devices
22 on the market. We searched the Medical Device
23 Reports database for device adverse events. Two
24 adverse events were found. The first adverse
25 event, reported in January of 1998, was a

1 significant blistering caused shortly after using
2 gel sheeting followed by full-thickness skin
3 necrosis due to secondary infection. The
4 blistering was not at the site of the gel sheeting
5 application but in the areas nearby. It was
6 determined by the reporting physician that the
7 event was unrelated to the device but we could not
8 rule out the possibility that the device was
9 involved.

10 The second adverse event, reported in
11 June, 2001, was an allergic reaction following the
12 use of silicone sheeting. Following 39 hours of
13 continuous use, the patient developed a severe red
14 rash and flaky rough skin. This was determined as
15 an isolated event and not likely that it was due to
16 the use of the device. Some possible causes for
17 the reported incident may be a reaction to the tape
18 used to hold the sheeting in place or moisture
19 created under the silicone sheeting after wearing
20 the product for such an extended period of time.

21 This slide is the questions to the panel.
22 Can I read them out?

23 DR. WHALEN: I will just interject that we
24 will not be answering the questions at this time;
25 we will at a later point in time, but please do

1 proceed, if you would.

2 DR. AREPALLI: Thank you. We have these
3 two topics for panel discussion. Following this,
4 Ms. Marjorie Shulman will walk you through the
5 classification work sheet. Here are the two
6 issues:

7 Please discuss the proposed classification
8 for the scar management device for the management
9 of hypertrophic and keloid scars. Also, discuss
10 what descriptive information and intended use
11 should be included in the proposed classification
12 identification.

13 Number two, please discuss the risk of
14 possible adverse skin reaction due to lack of
15 biocompatibility for the scar management device and
16 identify any other risks to health for these
17 devices. Thank you. Marjorie?

18 **Panel Discussion**

19 DR. WHALEN: Just a moment, does any of
20 the panel members have questions of Dr. Arepalli on
21 his presentation?

22 [No response]

23 Thank you. We will get to Ms. Shulman in
24 a moment but we are going to have a general
25 discussion first. Are there comments or questions

1 of any of the panel members about what we have just
2 heard on these devices? If I could kindle the fire
3 by asking Dr. Chang if she has any comments on the
4 subject that was raised about intact skin.

5 DR. CHANG: I have a comment and a
6 recommendation for the panel to consider which is,
7 rather than saying that this management is for
8 uninjured skin, to describe it as that it should be
9 intended for closed or intact skin, "closed, intact
10 skin" as a replacement for the wording "uninjured"
11 because, by definition, we are--

12 DR. WHALEN: Uncompromised.

13 DR. CHANG: Yes, uncompromised. We are
14 proposing this device for scar which is not the
15 same as uninjured skin.

16 The other question I had would be to ask
17 Dr. Newburger's opinion regarding whether or not a
18 skin rash could occur at a site distant from where
19 this product might be applied. In other words,
20 what is the potential for development of a rash to
21 be related to use of a gel padding in one location
22 and seeing a rash appear at a different site in the
23 body?

24 DR. NEWBURGER: To my knowledge, the issue
25 of silicones and true allergic reactions is pretty

1 well limited to foreign body reactions. We use
2 various types of silicone gel sheeting in our
3 practice and we have never seen a true allergic
4 reaction to it, and you might say, well, in a
5 private practice how substantial is this? Well, we
6 have over 30,000 patients and we are using it
7 multiple times every week and we have never seen a
8 true allergic reaction. We have seen folliculitis.
9 We have seen irritant reactions. We have seen
10 problems with tape. But when we are putting it
11 over skin which has healed over, and that is the
12 only time we use it, we have not seen allergic
13 contact dermatitis, nor have we seen distant
14 reactions.

15 We have had a number of patients who have
16 had solid silicone implants and we have seen
17 distant reactions which have been identified as
18 silicone granulomas. So, I am a little concerned
19 about this information that silicone is actually
20 released into the scar tissue. This is new to me
21 and I thought I had read rather extensively on the
22 area. I am concerned about that possibility.
23 Certainly, the identity is not the same as you are
24 going to see in the solid silicone rubber implant,
25 which is what we have seen occur, but this raises

1 more questions to me rather than fewer.

2 DR. WHALEN: Dr. McCauley?

3 DR. MCCAULEY: One of the things that we
4 actually think about when we are using these types
5 of products is that we actually talk about using
6 them in hyperproliferative scar disorders, which
7 keloids and hypertrophic scars fall under. I guess
8 that would distinguish it from some of the other
9 hyperproliferative skin disorders that occur
10 without trauma and that occur in dermatology. But
11 I would propose that actually rather than say
12 "uncompromised skin" we actually focus in on
13 hyperproliferative scar disorders which hypertrophy
14 scars and keloids represent.

15 DR. WHALEN: Dr. Witten, if I could ask a
16 question, it strikes me that if FDA were going to
17 be considering the use of this product on an open,
18 fresh wound that it would go well beyond the scope
19 of a reclassification or classification type of
20 process. Am I correct in that? Would that require
21 some other initiative on the part of a manufacturer
22 or sponsor wishing to have that indication?

23 DR. WITTEN: Well, we are really asking
24 you to classify what we have actually seen or
25 cleared, which have been these devices for scar

1 management.

2 DR. WHALEN: But within the scope of what
3 we are classifying, we are going to be defining the
4 safety and efficacy of the product with its
5 intended use in mind. So, are we at liberty in a
6 classification hearing to be considering a broad
7 scope of indications?

8 DR. WITTEN: We are only asking you to
9 consider the scope of indications for which we have
10 cleared the product.

11 DR. WHALEN: Thanks. Other issues or
12 points? Dr. Miller?

13 DR. MILLER: Yes, I would like to
14 emphasize that I think this product should probably
15 not be used on an open wound. It is a very
16 different situation than a closed wound with
17 hypertrophic scarring and I think that should be
18 emphasized. The use should be limited. I like the
19 words of an intact wound, an epithelialized wound
20 or a closed wound. I agree, uncompromised skin is
21 not very precise but, certainly, it needs to be a
22 closed wound.

23 DR. WHALEN: Any other comments
24 specifically on the semantics of the indication or
25 any of the issues that have been brought up?

1 DR. CHANG: I would like to second and
2 ditto Dr. Miller's comments that this should be
3 limited to closed, intact skin and not be placed on
4 an open wound.

5 DR. WHALEN: Dr. Witten?

6 DR. WITTEN: Yes, I just wanted to say
7 what Mr. Hurts clarified for me. Actually, I
8 should have remembered to say this, but we already
9 have a classification for open wounds. I mean,
10 there are Class I exempt wound dressings for open
11 wounds. So, there already are classifications for
12 products intended for open wounds; they would
13 already fit into a different classification.

14 DR. WHALEN: Is there a consensus on the
15 wording, that we are going to go forward with? Are
16 we are going to say closed wounds or, Dr. McCauley,
17 if you could say it again?

18 DR. MCCAULEY: Closed hyperproliferative
19 scar disorders.

20 DR. WHALEN: That implies that you are
21 talking about closed wounds.

22 DR. MCCAULEY: Right.

23 DR. WHALEN: Are there any other comments?

24 [No response]

25 Then, at this time we would like to begin

1 to focus our discussion on the FDA questions that
2 Dr. Arepalli has brought forward to us and that
3 remain projected on the screen. At this time we
4 will not refer to the reclassification
5 questionnaire. We will do that after this
6 discussion that is focused upon those questions.
7 Please consider, panel members, the silicone
8 elastomer for scar management device while
9 responding to the questions before us taken one at
10 a time.

11 The first question again, discuss the
12 proposed classification for the scar management
13 device for the management of hypertrophic and
14 keloid scars. Also, discuss what descriptive
15 information and intended use should be included in
16 the proposed classification identification.

17 Dr. McCauley, would you care to start off
18 on that one?

19 DR. MCCAULEY: First, I would like to have
20 some comments. Basically, this is related to the
21 information which has been presented to us relative
22 to this whole classification of silicone polymers.

23 Number one, they have been around for
24 quite a while and, number two, they have not, in my
25 opinion, posed a significant danger, if you will,

1 to patients. However, what bothers me is the fact
2 that there are a number of studies which have been
3 published that are, number one, anecdotal or,
4 number two, if they have been controlled,
5 randomized studies they are very small. Number
6 three, the mechanism of action for these materials
7 really has not been clarified.

8 I think that is very important in our
9 deliberations in terms of exactly how you want to
10 classify these devices. If you say that silicone
11 leaches out of these polymers into the wound and
12 affects CD36 cells, then you are really talking
13 about something that is more interactive and
14 something that may be classified as a Class III.

15 If you feel that the silicone in and of
16 itself is non-interactive but that it achieves this
17 effect just by coverage, although we know it is
18 probably not pressure that gives this effect, there
19 is some controversy in terms of whether temperature
20 really matters. Some studies by Lee suggest that
21 two degrees centigrade elevation in the temperature
22 underneath these materials causes a tremendous
23 increase in the action of collagenase, which is how
24 these effects are achieved. Other studies have not
25 shown that. Other studies have said that hydration

1 may be the mechanism by which we see improvement in
2 the wounds.

3 But I think it is very important to try to
4 decide what is the mechanism of action before we
5 can actually properly classify these compounds.

6 DR. WHALEN: Just to play the devil's
7 advocate, if these have been in use for so many
8 years and, in your opinion, you say you feel pretty
9 much that they are safe, from a pragmatic point of
10 view do you think it is that critically important
11 after all these years to delineate that mechanism
12 of action?

13 DR. MCCAULEY: I think it is important to
14 delineate that. Whether or not that is important
15 enough for classification, I think if we consider
16 the fact of this new data which was presented
17 relative to the leaching of silicone out of the
18 compound into the wound, I think that is a little
19 disturbing to me.

20 DR. WHALEN: Dr. Dubler, any comments?

21 DR. DUBLER: Dr. McCauley, in the study
22 that showed there might be some leaching of the
23 silicone into wounds---

24 DR. MCCAULEY: I am sorry, this is
25 information that was just provided to us by Mr.

1 Fallon.

2 DR. DUBLER: That is right, because there
3 are no published studies that we have reviewed that
4 have indicated that.

5 DR. MCCAULEY: Exactly.

6 DR. DUBLER: I also don't know what to do
7 with that piece of information. If that were, in
8 fact, the case then I think it would require the
9 sort of monitoring and data collection that would
10 probably only happen in Class III and, yet, the
11 published studies thus far--I can't comment,
12 obviously, on their statistical validity. They are
13 somewhat small but they didn't indicate that sort
14 of a problem so I wasn't prepared for that.

15 DR. MCCAULEY: Exactly.

16 DR. DUBLER: Therefore, based on the
17 studies that were here, it seemed to me that the
18 descriptive information would be relatively easy to
19 compile and prepare; now I am a little uncertain.

20 DR. WHALEN: Dr. Choti??

21 DR. CHOTI: I agree. I was kind of
22 expecting this to be fairly straightforward but if
23 it is really a topical application of an agent that
24 really has a direct impact on the local tissue,
25 then it muddies the water a little bit. But simply

1 based on the track record that these have been safe
2 and the reactions have been very minimal, I think
3 that the clinical safety data presented is quite
4 good. I think there is little clinical evidence to
5 suggest that there is any untoward effect of this
6 material and, therefore, I am not sure that the
7 Class III classification is warranted.

8 DR. WHALEN: Dr. Newburger?

9 DR. NEWBURGER: I concur with Dr.
10 McCauley's assessment about the need to clarify the
11 mechanism of action. Historically what we know
12 about this type of dressing, I would agree with Dr.
13 Choti.

14 DR. WHALEN: Dr. Miller?

15 DR. MILLER: Yes, I agree with all the
16 comments that have been made. I wonder, can we
17 invite our discussant back to the podium?

18 DR. WHALEN: The panel is free to ask
19 anyone a question that they wish.

20 DR. MILLER: Could Dr. Fallon come back up
21 because I too have been unaware of data which shows
22 there is a leaching of silicone into the wound that
23 is speculated to be the cause of the effect that we
24 see.

25 MR. FALLON: It is study number nine in

1 the packet.

2 DR. MILLER: Who is the author of that
3 one?/

4 MR. FALLON: That was Shigeki, Nobuoka.
5 It is a study done in Japan, published in Skin
6 Pharmacology Applied Skin Physiology.

7 MS. BROWN: I would like to ask a
8 question. Is this study relevant to silicone
9 sheeting or silicone gels?

10 MR. FALLON: Silicone sheeting, and we
11 believe the distinction between the sheeting and
12 the gel sheetings and the ointments is the release
13 of silicon, not silicone, from the sheetings. It
14 is hydrolyzed and these layers of silicone are
15 released.

16 I also have a study that shows all the
17 other proposed mechanism of actions have pretty
18 much been disproved. I really can't give it out
19 but it was done at Northwestern University by Dr.
20 Mustow.

21 DR. WHALEN: I am just perusing this for
22 the very first time, but it seems to me that these
23 are mostly in vitro skin specimens--

24 MR. FALLON: Yes.

25 DR. WHALEN: --where they are looking at

1 the distribution simply locally in a piece of skin.

2 MR. FALLON: Yes.

3 DR. WHALEN: It is not like they put this
4 somewhere in the groin and they--

5 MR. FALLON: Yes, correct, in vitro study,
6 yes.

7 DR. WHALEN: So, I am having a hard time,
8 again from first blush rapidly absorbing this,
9 saying that there is documentation of absorption
10 and systemic redistribution of silicone by this.
11 There is nothing in here that states that to me.

12 MR. FALLON: Yes, I was very surprised
13 when I saw that too. In reconsidering it with the
14 mechanism of action, as it is known, these scars
15 mostly happen to people that are from the tropics
16 and the CD36 has been connected to the prevention
17 of malaria and an antigen could have come--well, I
18 wasn't prepared for this, but I can prepare--

19 DR. WHALEN: Well, let me ask you a more
20 focused question. It is perhaps a slightly touchy
21 area. I would think from your particular vantage
22 point that you would not want to demonstrate that
23 this is systemically absorbed. Am I correct there?

24 [Laughter]

25 MR. FALLON: Yes, I am just presenting

1 what appears to be happening. I don't know, I
2 mean, I know I have a duty to my company, I own my
3 company but I am just presenting and I know it is
4 not helping and, you know, that is what it is. I
5 am presenting to the board what our findings are
6 and what we think.

7 DR. WHALEN: Are there any other
8 questions? Dr. Doyle?

9 DR. DOYLE: If we have something that has
10 been on the market this long with no untoward
11 effects, is it necessary that we know the mode of
12 action before we can approve it for classification?

13 DR. WHALEN: This is akin to the question
14 I asked Dr. McCauley a short time ago, and there
15 are certainly two answers to that question. There
16 is long-term demonstrated efficacy and, yet--it may
17 not be a related example--if we look at the
18 explosion of latex allergy, I am sure twenty years
19 ago people would have said there is millions and
20 millions of use of latex without too much of a
21 problem. So, without putting words in Dr.
22 McCauley's mouth, I think he is suggesting that if,
23 indeed, this is being now proposed as an effect we
24 certainly can't ignore it, and I would agree with
25 him. But I personally, from the perusal of this

1 letter, don't think that what we earlier had
2 suggested to us has been demonstrated by this
3 particular investigation.

4 MR. FALLON: And also, the study we did,
5 the ten-year old Cica-Care did not throw off any of
6 the silicon particles and that needs to be
7 investigated. It is thought that silicone becomes
8 toxic.

9 DR. WHALEN: Are there other questions for
10 Mr. Fallon?

11 DR. MILLER: I just have one more. Based
12 upon what you are telling us and what you have
13 learned, your recommendation remains that we
14 classify this as a Class I device?

15 MR. FALLON: Yes, on hypertrophic and
16 keloid scars, yes. Yes, definitely. I mean, I
17 don't see any safety issue. You could call Dr.
18 Brawer. His number is right there. He is an
19 expert. He is fairly articulate and you can ask
20 him directly. He is more of an expert than I am.

21 DR. WHALEN: Dr. Dubler?

22 DR. DUBLER: Therefore, you have no other
23 data but for the article that is numbered 9 that
24 would demonstrate any danger from use on scar
25 tissue?

1 MR. FALLON: Yes. No, I don't have any
2 other data, no.

3 DR. DUBLER: So, if article number 9 is
4 extinguished, then we are--

5 DR. MILLER: Yes, it also supports the
6 fact that this should be used only on intact
7 wounds.

8 DR. CHOTI: Although we still don't know
9 the mechanism of action.

10 MR. FALLON: That is the proposed
11 mechanism of action.

12 DR. CHOTI: Well, I think there is
13 hydration, there is temperature, there are other
14 modes of action. The real answer is whether you
15 discount this article or not, we don't know how
16 this works.

17 MR. FALLON: Yes.

18 DR. WHALEN: Thank you, sir. I think we
19 are getting around to Dr. Chang.

20 DR. CHANG: I would use that same analogy
21 of latex gloves, long history of use, relative
22 safety, low percentage of side effects although
23 there has been in certain populations, such as
24 those with spinal cord injury meningomyelocele,
25 certain increased risk for development of latex

1 allergies. So, that is not to say that with
2 increasing use of silicone gel products individuals
3 having this type of reaction may come forth. So, I
4 would use that analogy to say, yes, there remains
5 the potential, particularly if there is shedding or
6 potential for absorption of silicone products with
7 long-term use, that we may see this increasing
8 prevalence.

9 But I believe, looking at the data
10 presented by both FDA and industry, that there has
11 been a long record of relative safety in the face
12 of efficacy for this product. So, I would
13 emphasize that it is intended for intact, closed
14 skin and that it should be put into Class I.

15 DR. WHALEN: We haven't heard the words
16 t-test today. Dr. DeMets, any comments?

17 DR. DEMETS: I just want to second what
18 Dr. McCauley said. When I looked at the articles
19 that were in our tab, I was struck by sort of two
20 things. One, these studies are small and,
21 therefore, whatever the effectiveness is, is going
22 to be determined somewhat imprecisely, and some of
23 them were uncontrolled and those that were
24 controlled have a lot of missing data. So, I was
25 less than overwhelmed with the benefit side of the

1 equation. Because these were small, I was sort of
2 pondering the side effect side because if it is a
3 low but serious event the number of patients in
4 these studies would be way too small to detect
5 that.

6 Now, maybe there are registries or the FDA
7 database that could address that, but based on that
8 literature I reviewed, I did see that we have
9 enough numbers of patients exposed to really say
10 too much about the safety. I am coming to this
11 totally cold but that is just what I reflected when
12 I read it.

13 DR. WHALEN: Ms. Brown?

14 MS. BROWN: I would support the Class I
15 classification. As I understand, these have been
16 regulated under 510(k)s since 1976. Is that
17 correct, David?

18 DR. KRAUSE: They are considered
19 pre-amendments. So, they have been around since
20 before 1976. I don't think we had our first 510(k)
21 for them until sometime in the '80's but they were,
22 you know, identified as a pre-amendments device by
23 that submission.

24 MS. BROWN: But it sounds like there has
25 been a fair history of marketing with the product,

1 and there is a medical device reporting mechanism
2 so if there are problems they do get reported to
3 FDA. And, from what Sam Arepalli said, it sounds
4 like there have only been two. So, it sounds like
5 the risk is very minimal.

6 DR. WHALEN: Dr. Doyle, anything further
7 on the first question? No? All right, the second
8 question then that Dr. Arepalli has posed to the
9 panel is still projected. Please discuss the risk
10 of possible adverse skin reaction due to lack of
11 biocompatibility for the scar management device and
12 identify any other risks to health for these
13 devices. Dr. Dubler, any thoughts?

14 DR. DUBLER: I have one question about the
15 devices we haven't talked about thus far, which is
16 they aid in the resolution of certain complex or
17 difficult scar tissue. Would that scar tissue heal
18 on its own over time, or does this do something
19 that will create a different outcome?

20 DR. WHALEN: Yes, there might be multiple
21 answers but I think the answer is yes to your
22 question, both because it can both expedite and
23 change outcome. But the primary beneficial effect,
24 I think, has been to ameliorate the degree of
25 hypertrophy within the scarring process so outcome

1 would be changed. Any of the plastic surgeons like
2 to have another opinion on that?

3 DR. DUBLER: Since there can be a
4 beneficial outcome which would not occur but for
5 the use of this, and since there doesn't appear to
6 be any documented negative reaction but for two
7 cases, and since I don't know what to do with
8 article number 9, it seems there are no serious
9 adverse health reactions that would argue against
10 classifying as a Class I.

11 DR. WHALEN: Dr. Choti?

12 DR. CHOTI: I agree.

13 DR. WHALEN: Dr. Newburger?

14 DR. NEWBURGER: I agree as well.

15 DR. WHALEN: Dr. Miller?

16 DR. MILLER: I agree.

17 DR. WHALEN: We are on a concise streak.
18 Dr. Chang?

19 DR. CHANG: I agree. The caveat is in the
20 usage. In one of the two examples a patient had
21 the product on for over 30 hours. So, I believe
22 that in labeling patient education has to be very
23 important to limit to to 12 hours, I believe, the
24 consecutive hours that this product should be
25 applied to try to decrease the amount of skin rash

1 as a result of excessive moisture from overuse of
2 the product. But, otherwise, I agree that this
3 should be a Class I product.

4 DR. WHALEN: Dr. DeMets?

5 DR. DEMETS: I agree with the previous
6 comments.

7 DR. WHALEN: Ms. Brown?

8 MS. BROWN: I agree with the previous
9 comments.

10 DR. WHALEN: Dr. Doyle?

11 DR. DOYLE: I agree.

12 DR. WHALEN: And, Dr. McCauley?

13 DR. MCCAULEY: Same, I agree.

14 DR. WHALEN: Now that the panel has
15 discussed the FDA questions and our deliberations
16 seem complete, we have time for any final remarks.
17 Dr. Arepalli, is there any final comment from FDA,
18 or anyone else on behalf of FDA?

19 DR. WITTEN: I just want to clarify that
20 the second question was to discuss or identify any
21 other risks that you all see. I think some were
22 noted and if those are all the risks, that is fine
23 but I just wonder if there are any other risks that
24 haven't been discussed that anybody wants to
25 comment on.

1 DR. WHALEN: I think we have hit them. Is
2 there any final comment from anyone in the silicone
3 elastomer for scar management industry? If so,
4 would you please raise your hand? Yes, sir? Would
5 you please again, even though you have spoken to us
6 before, give your name and affiliation and any
7 financial interest in the devices being discussed?

8 MR. DILLON: I am Mark Dillon. I am the
9 president of Bio Med Sciences and, obviously, I
10 have a financial interest in the company.

11 I have a couple of comments. One is that
12 I am aware of one paper that was done by Dr.
13 William Monofeld where he looked for traces of
14 silicon metal in skin biopsies taken, and I believe
15 from a control source as well, underneath the
16 treated area. If I recall correctly, he saw a
17 fairly high baseline content of silicon metal in
18 the skin which he concluded could be from a number
19 of different sources--the fact that silicone is
20 often coated on capsules to make them easier to
21 swallow and on hypodermic needles and so forth--and
22 he concluded that there was no increased amount of
23 silicon in the skin treated with silicone sheeting.
24 So, I thought I would share that with the panel.

25 I would strongly agree with the idea of

1 not indicating this product for use on open wounds.
2 For one, the obvious reason is it is unknown and
3 there is a classification for that already.
4 Secondly, these products are always reusable and
5 after the first application they are no longer
6 sterile, even if they were provided sterile and
7 most of them are not. So, I think with strong
8 labeling to indicate against use on open wounds,
9 that issue would be largely put to rest.

10 I would be willing to share a theory on
11 mechanism of action if the panel would like to hear
12 some of my experience. I have noticed clinically
13 that with the use of these products on hypertrophic
14 scars, particularly in contractures over a joint,
15 you can see a benefit in range of motion occur with
16 a period of hours of use. To me, this is an
17 indicator that there is a hydration mechanism and
18 that this effect will reverse itself if the use is
19 discontinued.

20 Secondly, we see over a longer period
21 of time a remodeling of the scar, which may or may
22 not be due to hydration, but that second effect is
23 what is more permanent. DR. WHALEN: Thank you,
24 Mr. Dillon.

25 MR. DILLON: Thank you.

1 DR. WHALEN: Mr. Fallon, do you have any
2 final remarks?

3 MR. FALLON: The trace of silicon in in
4 vivo models will probably not come up because it is
5 the job of the CD36 and scavenger cells to take
6 those away. So, I really can't see how one could
7 set up an experiment, besides in vitro, to show
8 that it is getting into the skin. So, I just
9 wanted to clarify that.

10 **Classification Questionnaire and Vote**

11 DR. WHALEN: Thank you. Now we will
12 complete the classification questionnaire and
13 supplemental data sheet. Ms. Marjorie Shulman, in
14 the Office of Device Evaluation Classification,
15 Reclassification, will assist us as we go along.
16 After the formal panel discussion of each question
17 we will note the answers for each blank on the data
18 sheet as Ms. Shulman reads them out, and she will
19 record it on the overhead for all of us to see.
20 The voting members of the panel will vote then on
21 the completed questionnaire and supplemental data
22 sheet and this will then constitute the panel's
23 recommendation to the FDA. Procedurally, are there
24 any questions on what we are about to do next?

25 [No response]

1 MS. SHULMAN: Are we ready? the first
2 part on the sheet is just your panel name and you
3 can fill that out. That is administrative, and the
4 date; the generic type of device.

5 Then the first question, is the device
6 life-sustaining or life-supporting?

7 DR. WHALEN: We can just go around the
8 table, and this is for voting members. So, we can
9 start on this first question, please, with Dr.
10 McCauley.

11 DR. MCCAULEY: The answer to the first
12 question would be no.

13 DR. DUBLER: The answer to the first
14 question is no.

15 DR. CHOTI: No.

16 DR. NEWBURGER: No.

17 DR. CHANG: No.

18 DR. DEMETS: No.

19 MS. SHULMAN: The first one is no. Is the
20 device for use which is of substantial importance
21 in preventing impairment of human health?

22 DR. WHALEN: Just to stagger the way we
23 answer them, Dr. Dubler?

24 DR. DUBLER: No.

25 DR. CHOTI: No.

1 DR. NEWBURGER: No.
2 DR. CHANG: No.
3 DR. DEMETS: No.
4 DR. MILLER: No.
5 MS. SHULMAN: The second one is no.

6 Number three, does the device present a potential
7 unreasonable risk of illness or injury?

8 DR. WHALEN: Dr. Choti?
9 DR. CHOTI: No.
10 DR. WHALEN: Dr. Newburger?
11 DR. NEWBURGER: No.
12 DR. CHANG: No.
13 DR. DEMETS: No.
14 DR. MILLER: No.
15 DR. MCCAULEY: No.
16 DR. DUBLER: No.

17 MS. SHULMAN: The third one is no. We now
18 go to number four, did you answer yes to any of the
19 above three questions? That answer is no.

20 Then we go to number five, is there
21 sufficient information to determine that general
22 controls are sufficient to provide reasonable
23 assurance of safety and effectiveness?

24 DR. WHALEN: Starting with Dr. Newburger?
25 DR. NEWBURGER: Yes.

1 DR. MILLER: Yes.

2 DR. CHANG: Yes.

3 DR. WHALEN: Dr. DeMets?

4 DR. DEMETS: I will vote no.

5 DR. WHALEN: Dr. McCauley?

6 DR. MCCAULEY: Yes.

7 DR. DUBLER: Yes.

8 DR. CHOTI: Yes.

9 MS. SHULMAN: The answer to that one is
10 yes. On your sheets, you may mark whatever you
11 voted yourself. So, if the answer to that is yes,
12 it is classified into Class I.

13 So, we can skip two. We actually get to skip all
14 the way to the second page because all the rest of
15 the questions apply to Class II or Class III
16 devices.

17 Question 11 is a prescription question.
18 Can there otherwise be reasonable assurance of its
19 safety and effectiveness without restrictions on
20 its sale, distribution or use because of any
21 potentiality for harmful effect or collateral
22 measures necessary for the device? If you answer
23 yes, you are saying it is not a prescription
24 device. If you answer no, you are saying it is a
25 prescription device.

1 DR. WHALEN: Beginning with Dr. Miller.

2 DR. MILLER: No.

3 DR. CHANG: Yes.

4 DR. MILLER: No means that it requires a
5 prescription, right?

6 DR. WHALEN: Just re-explain, please.

7 MS. SHULMAN: The question is backwards.

8 If you answer yes, it is not a prescription device.

9 If you answer no, it is a prescription device.

10 DR. WHALEN: Do you still wish, Dr.

11 Miller, to vote no?

12 DR. MILLER: Yes, no, I mean--

13 [Laughter]

14 --I feel it should be a prescription
15 device.

16 DR. WHALEN: That is a no. Are you still
17 yes, Dr. Chang? Dr. DeMets?

18 DR. DEMETS: I will be a no.

19 DR. WHALEN: Dr. McCauley?

20 DR. MCCAULEY: Prescription device, so
21 that makes it a no.

22 DR. WHALEN: Dr. Dubler?

23 DR. DUBLER: Are we allowed to talk among
24 ourselves?

25 DR. WHALEN: Sure.

1 DR. DUBLER: In other words, do you think
2 it would be important for a physician to know that
3 this use was taking place and to direct and
4 supervise its use?

5 DR. WHALEN: Well, there are individual
6 state laws, if I can interject, that regulate who
7 can write prescriptions, and there are certainly
8 many places in the United States now where
9 prescriptions can be independently written by
10 non-physicians, but it would be by a licensed
11 practitioner.

12 DR. DUBLER: So, someone should be aware
13 of the use and supervise the use who has
14 specialized medical knowledge.

15 DR. WHALEN: Correct.

16 DR. DUBLER: I agree. So, that should be
17 a no.

18 DR. WHALEN: If you agree with that
19 practice, yes.

20 DR. DUBLER: Okay. No.

21 DR. WHALEN: Dr. Choti?

22 DR. CHOTI: No.

23 DR. WHALEN: Dr. Newburger?

24 DR. NEWBURGER: No.

25 MS. SHULMAN: The answer to that is no so

1 we go to 11(b), identify the needed restrictions
2 for the device. The first one is only upon the
3 written or oral authorization of a practitioner
4 licensed by law to administer or use the device.
5 The second, to use only by persons with specific
6 training or experience in it. Third, to use only
7 in a certain facility. Or, you could come up with
8 any other.

9 DR. WHALEN: Among those choices, Dr.
10 Chang, what would you suggest?

11 DR. CHANG: Well, in the State of Iowa it
12 is available in drug stores without prescription.
13 To me, I will just repeat what I had said as an
14 aside, it is about 70 percent effective overall,
15 looking at the literature. It has low danger
16 provided the label says to not wear it more than 12
17 hours; to discontinue it if there is a skin rash;
18 and it is helpful but kind of a very fancy bandage
19 over intact skin. So, I don't believe a
20 prescription is necessary. It is available already
21 in my state.

22 DR. WHALEN: So, you wouldn't want to
23 check any of these off?

24 DR. CHOTI: I don't think it is indicated.

25 DR. WHALEN: Fair enough. Dr. DeMets?

1 DR. DEMETS: I am going to stick with my
2 colleague.

3 DR. WHALEN: Dr. McCauley?

4 DR. MCCAULEY: Since I voted no, I would
5 say that only upon the written or oral
6 authorization of a practitioner licensed by law to
7 administer these. From my interpretation of this,
8 this is the least restrictive of the three that we
9 have, is that not correct?

10 MS. SHULMAN: Yes, that is correct.

11 DR. DUBLER: I want to come back to Dr.
12 Chang because I think I might have voted
13 differently on 11(a) if I had heard your comment
14 before I voted.

15 DR. CHANG: That it is available?

16 DR. DUBLER: It is available in Iowa
17 over-the-counter, and there have been no reports of
18 excessive use reactions. I mean, there is nothing
19 negative in the literature.

20 DR. CHANG: Aside from the two.

21 DR. DUBLER: Aside from those two, right.
22 Do we know if it is available in the same way in
23 other states? It is?

24 DR. WHALEN: Can I just ask in that regard
25 though, since we know that adverse events are

1 grossly under-reported by physicians, number one,
2 what mechanism exists for the lay public on OTCs--

3 DR. CHANG: To see a physician because of
4 a rash.

5 DR. WHALEN: But they may just have rashes
6 and not be doing anything about it, other than stop
7 using it and seeing if the rash goes away.

8 DR. CHANG: And it should go away.

9 DR. MILLER: Can I make a comment? I
10 mean, I think the goal of this is to treat the scar
11 and this is a tool to treat the scar. I think
12 that, you know, a physician needs to evaluate the
13 patient and decide on how to treat the scar. For
14 people just to go on their own and select this, I
15 think it doesn't make as much sense to me because
16 they may be selecting it for the wrong types of
17 scars, the wrong types of problems, and I think it
18 should be guided by a physician.

19 DR. WHALEN: Getting back to going around
20 the table, Dr. Dubler?

21 DR. DUBLER: Yes, I find this puzzling
22 because if we lived in a nation where everybody had
23 access to physicians or nurse practitioners or
24 other people who could manage their care, then I
25 would tip in one direction. But since some 42

1 million people don't and so requiring a
2 prescription will, in fact, be a barrier to access
3 for something that could be helpful in the long
4 run, I would like to change my vote on 11(a) to a
5 yes and, therefore, I don't need to choose anything
6 from 11(b). Correct?

7 DR. WHALEN: I think that is perfectly
8 acceptable. That would make the vote still 5-2 in
9 favor of no in question 11, unless there is anyone
10 else who wishes to reconsider.

11 MS. BROWN: Could I ask a question? If
12 the panel votes that this needs a prescription, is
13 the State of Iowa now going to have to take it off
14 the over-the-counter mechanism of distributing the
15 product?

16 DR. WHALEN: Dr. Witten?

17 DR. WITTEN: I am not sure what the State
18 of Iowa would do, but if we make it a prescription
19 use, then they have to interpret what that means.

20 DR. WHALEN: Keeping in mind that we are
21 an advisory panel and our advice is going to the
22 FDA to deal with this as they see fit, again, on
23 question number 11, is there anyone else who voted
24 in either direction and wishes to change their
25 vote?

1 [No response]

2 Dr. Choti, which among the options in
3 11(b) would you choose?

4 DR. CHOTI: I would say that written or
5 oral authorization is warranted in this situation,
6 the first one.

7 DR. WHALEN: The first? Very good. Dr.
8 Newburger?

9 DR. NEWBURGER: Also written or oral
10 authorization of a practitioner.

11 DR. WHALEN: Dr. Miller?

12 DR. MILLER: I agree, written and oral
13 authorization.

14 MS. SHULMAN: That is it for the general
15 device questionnaire. We will move on to the
16 supplemental data sheet. The first question for
17 your sheet, the generic type of device we have
18 covered that. You can just write that in. The
19 advisory panel is surgery, General and Plastic
20 Surgery.

21 Number three, is the device an implant?
22 No. The indications for use, here we won't have to
23 rewrite it if everybody agrees to the indication
24 that was presented during the meeting earlier.

25 DR. WHALEN: This indication was one that

1 had the wording in it about uncompromised skin.

2 So, is there anyone who wishes to modify that in
3 any way?

4 DR. DUBLER: I thought Dr. McCauley--

5 DR. MCCAULEY: I would like to modify it
6 and say intact hyperproliferative scar disorders,
7 which includes keloids and hypertrophic scars.

8 DR. WHALEN: Is there consensus on that
9 wording?

10 DR. CHANG: To clarify, would that be not
11 using for the control of hypertrophic and keloid
12 scar--I mean, if we have both the words
13 hypertrophic and keloid scar and then put in the
14 words for intact hyperproliferative skin disorder,
15 then we have it duplicated.

16 DR. MCCAULEY: I think the way it reads is
17 for intact or uncompromised skin for scar control.
18 Is that not correct?

19 DR. WHALEN: Scar management.

20 DR. MCCAULEY: Or scar management.

21 DR. WHALEN: Scar management device is a
22 silicone sheeting product intended for use on
23 uncompromised skin for scar management.

24 DR. MCCAULEY: Yes, for scar management.

25 So, intact--

1 DR. CHANG: And for proliferative---

2 DR. MCCAULEY: Scar disorders. Then in
3 parentheses you can put hypertrophic scar and
4 keloids. But I think the key basically is to
5 encompass both of them and make sure that you use
6 them for intact.

7 DR. CHANG: I agree as proposed.

8 DR. MILLER: Instead of intact could we
9 say epithelialized wounds or closed wounds? I just
10 happen to like epithelialized wounds, it is more
11 specific to me.

12 DR. MCCAULEY: I have no objection.

13 DR. CHANG: I would vote for keeping it
14 simple, and if you want to be explicit about intact
15 I would vote to say closed, intact.

16 DR. MILLER: I like that, closed, intact.

17 DR. WHALEN: Closed intact? Technically,
18 we would almost have to have a biopsy to
19 definitively declare that it is epithelialized.

20 DR. CHOTI: Well, the distinction is
21 perhaps a fresh incision and it may be semantics as
22 to whether it is a closed wound or not but,
23 clearly, it is not to be applied on a freshly
24 closed incision.

25 DR. MCCAULEY: What you are saying is that

1 it has to be a hyperproliferative problem.

2 DR. CHOTI: Yes, I think that just
3 replacing the word "uncompromised" with "intact" is
4 sufficient, and leave it just for scar management
5 without specifying hyperproliferative state.

6 DR. WHALEN: The only problem with healed,
7 as you plastic surgeons know better than I, you can
8 make an argument for a year that it is not entirely
9 healed.

10 DR. MILLER: That is true.

11 DR. WHALEN: Even though it is totally
12 epithelialized.

13 DR. MILLER: That is correct; that is
14 true.

15 DR. WHALEN: So, where are we?

16 DR. CHANG: Back to intact.

17 DR. MCCAULEY: Hyperproliferative scars.

18 DR. CHANG: Yes.

19 DR. WHALEN: So, intact skin for
20 management of hyperproliferative scars?

21 DR. MILLER: Right.

22 DR. WHALEN: Is that the way we are doing
23 it?

24 DR. CHANG: Intact skin with
25 hyperproliferative scars, parentheses, hypertrophic

1 or keloid.

2 DR. WHALEN: Okay.

3 DR. AREPALLI: Are you going to stick with
4 "management of?"

5 DR. WHALEN: No, I don't think "management
6 of," Sam. I didn't hear that. The word management
7 is the third word, the scar management device is a
8 silicone sheeting product intended for use on
9 intact hyperproliferative scars, parentheses for
10 keloids and hypertrophic scars.

11 DR. NEWBURGER: Question.

12 DR. WHALEN: Yes, ma'am?

13 DR. NEWBURGER: Is not one of the intents
14 of these dressings to be used in an area where you
15 strongly feel that there is going to be a keloid?
16 Can't you use that on a preventative basis? I
17 thought that was some of the information that we
18 got, if you have an incision that is, you know, in
19 this triangle and you have someone who
20 characteristically forms keloids, wouldn't you want
21 to use this as soon as the area has epithelialized,
22 Dr. McCauley?

23 DR. MCCAULEY: There was one study in our
24 packet that actually dealt with that issue. You
25 all may know better than I, but I didn't feel that

1 it really had a lot of good data, but what they did
2 suggest was that maybe in areas which are prone to
3 the development of hypertrophic scars it may be
4 useful in terms of prevention, but there was just
5 one paper.

6 DR. WHALEN: Phyllis?

7 DR. CHANG: I would be content to leave
8 that as an off-the-shelf use.

9 DR. WHALEN: All right.

10 MS. SHULMAN: So we agree upon the
11 wording. Number five, the identification of any
12 risks to health presented by the device. We can
13 say as covered in the panel meeting or anyone can
14 add anything they wanted to.

15 DR. WHALEN: Agreeable to say as covered
16 in the panel meeting? All right.

17 MS. SHULMAN: Number six, the recommended
18 advisory panel classification and priority--we only
19 need a classification, which is Class I, and the
20 priority we only need for Class II or three.

21 Number seven, if device is an implant or
22 is life-sustaining or life-supporting and has been
23 classified in a category other than III, explain
24 fully the reasons. We can skip this.

25 Number eight, the summary of information

1 including clinical experience or judgment upon
2 which a classification recommendation is based. If
3 you wish, we could put in there what was covered in
4 the panel meeting as a summary for the reasons.

5 DR. WHALEN: Seeing no objections, we will
6 do that.

7 MS. SHULMAN: Number nine, the
8 identification of any needed restrictions on use of
9 the device. That is a prescription question again.
10 We can just refer to question 11(a) of the general
11 device questionnaire.

12 Number ten, if the device is in Class I,
13 recommend whether FDA should exempt it from
14 registration and listing, premarket notification,
15 records and reports and good manufacturing
16 practices. It can be all, any or none.

17 DR. WHALEN: I may have lost our order
18 track but I think Dr. McCauley, if we could start
19 with you on this?

20 DR. MCCAULEY: Shall I go through each one
21 individually?

22 DR. WHALEN: Or any of those that you wish
23 to say it should be exempted from.

24 DR. MCCAULEY: None.

25 DR. WHALEN: Dr. Dubler?

1 DR. DUBLER: Then I don't understand the
2 question. I would assume, given Class I, we would
3 want to exempt it from (a) and (b). Doesn't that
4 follow from Class I?

5 MS. SHULMAN: Well, registration is where
6 you register your manufacturing facility and
7 listing is where you list the device. Number two
8 is premarket notification. Most Class I devices
9 are exempt from premarket notification, however,
10 there are about 53 reserved Class I so they do
11 require 510(k)s to come in even though they are
12 Class I.

13 DR. MCCAULEY: Can I ask for discussion,
14 particularly from Dr. DeMets? In your review of
15 the statistical data, what is your opinion in terms
16 of efficacy based on the data that you were
17 presented?

18 DR. DEMETS: Well, only what was in our
19 packet because that is all I know. The studies
20 were small. Some of them were uncontrolled. Some
21 of them, even though they had controls, had
22 substantially missing data or patients were
23 excluded from the analysis which leaves it open to
24 some potential to bias. So, I am not saying it was
25 not effective, I am just saying I am not sure how

1 effective it is based on the data that was
2 presented. There are pretty small numbers and, you
3 know, by good clinical trial design these are not
4 particularly strong studies. So, my earlier
5 remarks were not on the safety part so much as on
6 the efficacy part.

7 DR. WHALEN: Does that answer it for you?

8 DR. MCCAULEY: Yes, well, if I can get
9 some clarification again on registration and device
10 listing?

11 MS. SHULMAN: Registration is a paper
12 format where people send into our Office of
13 Compliance where their manufacturing facility is
14 located. Listing is where they list what devices
15 they are making, and that is for inspectional
16 purposes.

17 DR. MCCAULEY: Is that not part of GMP or
18 is that a separate issue?

19 MS. SHULMAN: It is separate. Few devices
20 are exempt from registration and listing, but there
21 can be some that are.

22 DR. DUBLER: And what does C mean, records
23 and reports?

24 MS. SHULMAN: That is another compliance
25 issue on their record keeping. Everyone has to

1 keep records, but the manner in which they do it,
2 would they have to follow our rules, the rules they
3 would follow under the GMPs, the good manufacturing
4 practices.

5 DR. DUBLER: So, we have had two reports
6 of adverse reactions. If there were other such
7 reports I would want them to come forward. If I
8 check any of these things, A, B, C and D, does it
9 prevent that adverse reporting or does it
10 discourage it?

11 MS. SHULMAN: No, I don't believe it does.

12 DR. MCCAULEY: It means they are not
13 required to report them.

14 MS. SHULMAN: I am sorry, I misunderstood,
15 yes, they would not be required to report.

16 DR. DUBLER: So, the first one just means
17 they have to tell us where they are.

18 MS. SHULMAN: They have to tell us where
19 they are and list what devices they are making in
20 that facility.

21 DR. WHALEN: All right.

22 DR. CHANG: And, could you clarify again
23 what is item B? What does that mean?

24 MS. SHULMAN: The second one is premarket
25 notification, also known as 510(k). So, if it is

1 exempt from that they can go to market without
2 coming in and getting a clearance from us.

3 DR. DUBLER: But they are already on
4 market.

5 DR. WITTEN: A sponsor with a new device,
6 if somebody comes in with a new device, if they
7 want to market and they need to submit a premarket
8 notification, that means they send a premarket
9 notification to us to review before they go to
10 market. If they are exempt from premarket
11 notification and, as Marjorie Shulman already
12 mentioned, most Class I devices are but there are
13 some that are reserved, then, if they are exempt,
14 they don't need to send an application. If they
15 are not exempt they need a specific clearance from
16 us prior to going to market.

17 DR. MCCAULEY: As I recall, each of the
18 industry representatives recommended Class I with a
19 510(k). Is that not correct?

20 DR. WHALEN: They did. With no offense to
21 our industry representatives, they certainly would
22 have an interest in so recommending. Any industry
23 representative would have it in their own best
24 interest to put up a potential wall for competitors
25 entering the marketplace. That is not to say that

1 their intent is not noble and scientifically
2 founded. Are we clear on what we are talking about
3 in this question? Dr. McCauley?

4 DR. MCCAULEY: I think I will stick with
5 my original vote for no exemption.

6 DR. WHALEN: And Dr. Dubler, will you
7 still go with A or B?

8 DR. DUBLER: I would exempt B but not the
9 others.

DR. WHALEN: Solely B,
10 but you would not vote for A?

11 DR. DUBLER: Solely B.

12 DR. WHALEN: Dr. Choti?

13 DR. CHOTI: I think I would not exempt any
14 of them.

15 DR. WHALEN: Dr. Newburger?

16 DR. NEWBURGER: No exemption.

17 DR. WHALEN: Dr. Miller?

18 DR. MILLER: No exemptions.

19 DR. WHALEN: Dr. Chang?

20 DR. CHANG: No exemptions.

21 DR. WHALEN: Dr. DeMets?

22 DR. DEMETS: No exemptions.

23 MS. SHULMAN: Then number 11, if there are
24 any existing standards to the device assemblies,
25 components, device materials or parts or

1 accessories that you know of that you would like us
2 to apply to these devices, then this is where you
3 can list them.

4 DR. WHALEN: I don't know that we need to
5 go around for this. Is there anyone who wishes to
6 stipulate such? I see none.

7 MS. SHULMAN: Then that is the end of the
8 sheet. You go around once and vote for these
9 sheets to be voted on as discussed as a Class I
10 reserve device, requiring 510(k).

11 DR. WHALEN: So, in effect then, we are
12 asking for a motion to accept the classification
13 worksheet as filled out, with a recommendation for
14 Class I silicone elastomer for scar management
15 intended for use on intact skin, hyperproliferative
16 scars, parentheses, keloid and hypertrophic scars.
17 Is there a motion to that effect?

18 DR. CHANG: So moved.

19 DR. WHALEN: Is there a second?

20 DR. CHOTI: Second.

21 DR. WHALEN: It has been moved and
22 seconded that silicone elastomer for scar
23 management intended for use in intact skin
24 hyperproliferative disorders, parentheses, keloid
25 and hypertrophic scars, be classified into Class I.

1 All those in favor, voting members signify by
2 raising their hands, please.

3 [Show of hands]

4 All of those opposed? Dr. Dubler
5 abstains. So, it is six yes, one abstention.

6 DR. DUBLER: Can I take just one minute,
7 Dr. Whalen?

8 DR. WHALEN: You certainly can because
9 each member has to take, maybe not one minute but
10 ten to fifteen seconds to explain why they have
11 voted in the way that they did. If we could start
12 with Dr. DeMets?

13 DR. DEMETS: I am not sure I can explain.
14 Well, I think that we have discussed the issues and
15 I can accept what we voted earlier.

16 DR. WHALEN: Dr. Chang?

17 DR. CHANG: My comment is that the science
18 is soft, as previously mentioned, but for some this
19 is efficacious. The track record over many years
20 is that it is a safe product. Side effects can be
21 prevented if it is used correctly. I would compare
22 use of nonsteroidal anti-inflammatories. They can
23 have serious side effects but they have become in
24 common use relatively safe. They can cause ulcers
25 but they are available over-the-counter and the

1 price has gone way down. So, if this is a useful
2 product, very safe, with many individuals having
3 the hyperproliferative scars, you know, I voted to
4 have it as a non-prescription Class I device.

5 DR. WHALEN: Dr. Miller?

6 DR. MILLER: I think it is a very
7 practical device that certainly appears safe.
8 Although we don't understand exactly why it works,
9 I don't think that should prevent us from making it
10 available.

11 DR. WHALEN: Dr. Newburger?

12 DR. NEWBURGER: I think this is a very
13 useful device. It has been very effective for many
14 people. By making it a prescription Class I device
15 I think we have the potential for avoiding a lot
16 more side effects, and I feel there are many more
17 than certainly have been reported. I think this
18 gives more safety to the community.

19 DR. WHALEN: Dr. Choti?

20 DR. CHOTI: I agree with the comments. I
21 think it sounds like this device is already being
22 used a lot. It sounds like it is safe and it
23 probably is effective. I think we are all a little
24 bit frustrated by the fact that we don't know how
25 it works; that we don't really have a lot of

1 records of its application. So, what can be done
2 by industry, academics and others to study a little
3 bit more the mechanisms and registry, and I think
4 we have voted, or I have voted in a way to, best as
5 we can, encourage some kind of additional record
6 keeping. But I think it sounds like it potentially
7 has a clinical role so that is how we voted.

8 DR. WHALEN: Dr. Dubler?

9 DR. DUBLER: I abstained for a very
10 particular purpose. I think that the panel's
11 discussion was very thoughtful but indicated to me
12 that this was, when used correctly for the right
13 indications, a safe application that could have a
14 real effect on someone's quality of life and on the
15 outcome of the resolution of these scars. I work
16 in the Bronx. We have a lot of people who don't
17 have health insurance, and when they have a problem
18 and they can deal with it over-the-counter they
19 have a chance of helping themselves. When they
20 have to go through a licensed practitioner they
21 don't get that help. A lot of kids have these
22 scars, a lot of mobility problems.

23 I abstained because I think we should not
24 have voted this to be a prescription item. I think
25 the single greatest ethical problem in American

1 medicine is access to care, and we have just put up
2 a barrier to what may be safe and helpful.

3 DR. WHALEN: If I can parenthetically add,
4 it is perhaps only fair, if I have in any way
5 impugned manufacturers in putting up 510(k)
6 restrictions, to state that physicians putting up
7 prescription barriers is probably not the most
8 disinterested party to do so.

9 DR. DUBLER: Here, here.

10 DR. WHALEN: Dr. McCauley?

11 DR. MCCAULEY: I think the device is safe.
12 I think that it probably is efficacious. I think
13 the data is somewhat soft, and I think the way we
14 voted probably will lend itself to really
15 determining how efficacious this product really is.

16 DR. WHALEN: Though not voting, any
17 comments, Dr. Doyle?

18 DR. DOYLE: I feel very strongly as Dr.
19 Dubler does since it is considered safe and has
20 been shown efficacious, or has not been shown not
21 to be efficacious, why are we limiting people's
22 access to it?

23 DR. WHALEN: Ms. Brown:

24 MS. BROWN: I have the same question about
25 access. If it is available currently

1 over-the-counter, it seems like we may have put up
2 another barrier to its use that wasn't here before
3 the panel meeting.

4 DR. WHALEN: Thank you. I would like to
5 announce that the recommendation of the panel, with
6 six votes for the motion and one abstention, is
7 that the silicone elastomer for scar management
8 intended for use in the management of intact
9 skin--which I still won't get right but what you
10 see up on the screen, with hyperproliferative
11 scars, parentheses, keloid and hypertrophic scars,
12 be classified into Class I.

13 In so doing, I would like to thank the
14 panel and thank Dr. Arepalli and the industry reps
15 for what they have done for us. We have another
16 item of business but we will take a ten-minute
17 break and reconvene promptly at 3:25 to being that
18 process.

19 [Brief recess]

20 **Reclassification of Absorbable Hemostatic Agents**
21 **and Dressings**

22 DR. WHALEN: I would like to call this
23 meeting back to order. Could I first ask that the
24 voting panel members pass toward the center of the
25 table, toward me, their classification

1 questionnaires on the last item of business so that
2 we can collect them for the FDA?

3 I would like to remind the public again
4 that while this portion of the meeting is open for
5 the public for their observation, public attendees
6 may not participate except at the specific request
7 of the panel.

8 We now will proceed to the open committee
9 discussion. We will begin the discussion on the
10 reclassification of absorbable hemostatic agents
11 and dressings with serial presentations from
12 industry, first by Dr. John D. Paulson, vice
13 president for quality assurance and regulatory
14 affairs, Johnson & Johnson Wound Management, a
15 Division of Ethicon, Inc., followed by Ms.
16 Ronnemoes Bobak, vice president for product
17 development, Ferrosan A/S, and then Ms. Judith E.
18 O'Grady, senior vice president, regulatory, quality
19 and clinical affairs, Integra LifeSciences
20 Corporation.

21 The FDA presentation and a reading of the
22 FDA questions will follow these presentations. We
23 will then have a general panel discussion of this
24 topic, followed by a more focused panel discussion
25 aimed at answering FDA's questions. Before we

1 complete the reclassification worksheet and
2 supplemental worksheet, we will have an open public
3 comment period. Then we will complete the
4 reclassification worksheet and supplemental
5 worksheet. The vote on these worksheets will
6 constitute the panel's recommendation to the FDA.

7 I would like to remind public observers at
8 this meeting that while this portion of the meeting
9 is open for public observation, again, public
10 attendees may not participate except at the
11 specific request of the panel. I probably should
12 seriously point out that even though I am from
13 Robert Wood Johnson Medical School in New
14 Brunswick, New Jersey, there is not a financial
15 interrelationship with Johnson & Johnson although,
16 God knows, my dean would love to have a stronger
17 one. We will begin with Dr. Paulson's
18 presentation.

19 **Industry Presentation**

20 DR. PAULSON: Dr. Whalen, Dr. Krause, Dr.
21 Witten and panel, thank you for the opportunity to
22 present here today.

23 There are several different types of
24 products in the category of absorbable hemostatic
25 agents. I am here today to present concerning

1 Surgical, oxidized regenerated cellulose,
2 representing one of these product types.

3 I would like to talk to you about globally
4 available ORC products while noting that Surgical
5 is currently the only available ORC product in the
6 United States. There was previously another
7 manufacturer making a similar product, using
8 essentially the same chemistry and manufacturing
9 process which we licensed jointly from the
10 third-party company. They have since stopped
11 making that product. But their safety record is
12 going to be discussed, I am sure, by Dr. Krause and
13 will reflect product made by the same manufacturing
14 process in essence.

15 I will talk to you briefly about the
16 manufacturing process; the mechanisms of
17 hemostasis, just very briefly; biocompatibility and
18 hemostasis data; and then provide a brief summary
19 of my conclusions. I will try to do that all in
20 less than fifteen minutes.

21 The Surgical family of absorbable
22 hemostats includes three basic product types,
23 Surgical, Surgical Nu-Knit and Surgical Fibrillar,
24 representing different physical forms of product
25 made with essentially the same process, although

1 representing different weaves and manufacturing
2 processes after the chemistry has taken place.

3 These products are used adjunctively in
4 surgical procedures for the control of capillary,
5 venous or small arterial bleeding and rapidly stop
6 the bleeding by acting as a matrix for the
7 formation of a clot, and some other mechanisms that
8 I will talk to you about a little bit later. The
9 product is often left behind in part or in whole
10 and is readily absorbed from the site of
11 implantation with minimal tissue reaction, which is
12 very important because it is frequently used in
13 cardiovascular procedures and frequently in
14 neurosurgical procedures where other methods of
15 hemostasis may not be suitable, for instance,
16 electrocautery.

17 There are other ORC products available in
18 the global market. This is Cellulostat. it is
19 available in Taiwan and China. You can see it
20 magnification at 12x having a slightly different
21 pattern of knit or weave. And, there is an ORC
22 product developed from Europe, by the name of
23 Curacel. So, I will talk a bit about those
24 products as well.

25 This is just to remind you that

1 regenerated cellulose products are more than just
2 an isolation of cellulose. They are derived from
3 wood pulp which contains about 50 percent cellulose
4 by weight, and also contains significant amounts of
5 lignin and other inter-fibrillar materials which
6 act as adhesives to kind of keep the physical
7 structure of the wood intact. There are
8 significant chemical processes in place here that
9 affect the qualities of the fabric which becomes
10 the raw material for the oxidation process. This
11 basically digests the cellulose and then
12 reconstitutes it prior to oxidation. The material
13 that we are talking about here is bright rayon, and
14 it is essentially pure cellulose.

15 This then goes into knitting and
16 purification processes and conditioning of the
17 fabric, controlled oxidation reactions which are
18 used to define chemistry and define processes,
19 involve displacement of solvents and reactants,
20 purification of the materials, dehydration and then
21 processing the material into its final product
22 forms, along with sterilization and QA testing and
23 release.

24 What I wanted to call your attention to in
25 all of this is really the complex nature of the

1 processing that is involved. This is not just an
2 isolation of cellulose as a chemical derivation,
3 and you can think of this as a biosynthetic
4 material rather than an isolated biological
5 material.

6 Cellulose is a polymer of glucose
7 basically, and oxidized regenerated cellulose in
8 its simplest form involves the oxidation at the
9 sixth position, changing it from an alcohol
10 function to a carboxylic acid function. There are
11 also other chemical byproducts, and I will call
12 your attention to the 2- and 3-ketone ORCs as well
13 as aldehydes, ketones, dialdehydes, and so on.
14 These can vary in ratios depending on the controls
15 and nature of oxidants in the oxidation process.

16 Again, we are talking about
17 cellulose-related materials. This is a reminder
18 that cellulose itself does not absorb. This is a
19 cotton suture that has been implanted for two
20 years. You can see a chronic inflammatory reaction
21 here and continued presence of the cellulose.
22 Cellulose is also well-known to surgeons from use
23 in gauze, and lint from gauze is well-known to
24 cause chronic inflammatory reactions and adhesions.
25 We don't want to end up with cellulose so

1 consistency of the process is important to achieve
2 a biocompatible and degradable material.

3 This is just a brief reminder of the
4 complex relationship between physiologic processes
5 involved in hemostasis, involving vasoconstriction,
6 platelet activation, coagulation activation,
7 conversion of fibrinogen to fibrin by thrombin. I
8 will mention briefly that Surgicel acts both in
9 terms of platelet activation and activation of
10 intrinsic and extrinsic pathways or coagulation
11 activation.

12 Surgicel at the wound site has multiple
13 mechanisms. Here it is applied to a vessel. There
14 is fluid absorption which results in a relative
15 hemoconcentration. There is hemoglobin oxidation
16 resulting in a gel formation or false clot which
17 helps achieve tamponade in conjunction with manual
18 compression.

19 You can see in the upper right-hand
20 depiction adherence of platelets to the fibrillar
21 structure of the material and ultimately clot
22 formation on the matrix of fabric.

23 I will use this swine spleen incision
24 picture to demonstrate some of the actual use of
25 the product. You can see here that in a model

1 which we used to measure hemostasis time, using a
2 controlled incision of 1.5 cm by 0.3 cm deep, that
3 you see a darkening of the blood which has to do
4 with its oxidation. You can see gel formation and
5 false clot formation, fluid absorption, and so on.
6 In this model we applied digital compression and
7 can measure time to hemostasis, as I mentioned
8 earlier. I will refer to some data later from this
9 model.

10 I am going to try to relate to you some of
11 the mechanisms of action that we talked about
12 earlier in hemostasis and some of the attributes of
13 Surgicel. I will also talk to you about how they
14 relate to controls which exist in the U.S.
15 Pharmacopeia.

16 In the first column you will see
17 mechanisms of action include tamponade due to
18 digital compression, fluid absorption, swelling and
19 gel formation. Then, at the chemical and biologic
20 level we talked earlier about protein adsorption,
21 platelet adherence, platelet aggregation and
22 platelet activation, and intrinsic and extrinsic
23 pathway activation. This relates to a variety of
24 physical and surface chemistry attributes of the
25 product that are shown in this panel. Then, if one

1 looks at USP specifications for oxidized
2 regenerated cellulose, you see that they are very
3 incomplete in addressing those items which we
4 identify as important to achieving effectiveness of
5 the product and in achieving hemostasis.

6 Similarly for biocompatibility, we have
7 highlighted just a few of the key areas of
8 biocompatibility that are important
9 here--cytotoxicity, acute inflammation,
10 biodegradation and absorption, immunogenicity and
11 neurotoxicity.

12 Surgicel properties include carboxyl
13 content and pH. Degradation, interestingly, does
14 not appear to be related to carboxylation of the
15 alcohol functions but, rather, the ketone formation
16 at C2 and C3, which is not controlled by USP
17 specifications or other recognized standards for
18 these products.

19 In terms of immunogenicity and
20 neurotoxicity, the exact determinants are not
21 clear, but it is clear that they depend on high
22 material purity and controlled chemical processes
23 and ingredients. Again, USP specifications appear
24 to be inadequate in addressing these essential
25 requirements of the product.

1 We have done some physical and chemical
2 analysis of the three product types that I showed
3 you earlier. Many of these are USP tests, for
4 instance, identification loss on drying, nitrogen
5 content and carboxyl content are USP tests. You
6 will see that Surgicel passes all of these, while
7 Curacel fails for carboxyl content; Cellulostat
8 fails on several accounts. In the right-hand
9 column you can see those USP specifications for
10 these parameters.

11 We have also assessed pH, which is
12 somewhat related to carboxyl content but
13 post-oxidation treatment can neutralize the pH and
14 add back other ions. You can see that Cellulostat
15 has a different pH, while Curacel appears to have
16 calcium added into the process.

17 Physical strength varies which, of course,
18 can affect clinical use. Water solubility varies
19 for these products and spectral identification
20 suggests that Cellulostat is not, in fact, ORC
21 despite its labeled claim to be so.

22 Time to hemostasis was measured in the
23 model that we showed you earlier. The top two bars
24 represent some historical data that we have for
25 Surgicel Nu-Knit and Surgicel. Then, due to

1 limited numbers of samples available for Curacel
2 and Cellulostat in the short time since this
3 meeting was announced, we have done some
4 head-to-head comparisons of Cellulostat and
5 Curacel, and you can see that with Surgicel the
6 mean time to hemostasis for those specific wounds
7 that we discussed earlier is approximately eight
8 minutes; for Curacel, approximately ten minutes;
9 and for Cellulostat, essentially none of them
10 achieved hemostasis in the 12 minutes that we
11 defined as the maximum time period for this model.

12 So, what do we conclude from all that?
13 Our summaries are that Surgicel Absorbable Hemostat
14 does, indeed, have a long history of safety and
15 effectiveness. I think Dr. Krause will speak to
16 you about that. Given the limited time for this
17 presentation, I haven't gone into it but it
18 certainly does have a commendable history.

19 There is complex chemistry and processing
20 required to create the unique product properties
21 here. There are multiple physiological
22 interactions required for safety and effectiveness.
23 Other ORC products are not equivalent, and USP
24 requirements do not address many critical product
25 attributes.

1 So, our conclusions are that the USP is
2 not adequate to control key product attributes, and
3 that we do not know of other standards for these
4 products which are established and accepted.

5 Finally, in the absence of recognized
6 standards, we believe that reclassification is not
7 appropriate. I would again refer back to the fact
8 that these products are often left as implantable
9 devices in critical portions of the circulatory
10 system and the central nervous system. So, I am
11 sure that as physicians and scientists you can
12 appreciate the great degree of control and
13 assurance of biocompatibility and effectiveness
14 that are needed there. Thank you.

15 DR. WHALEN: Questions for Dr. Paulson?

16 MS. BROWN: I have a question. If the FDA
17 developed a guidance document would you be in
18 support of a down-classification?

19 DR. PAULSON: I think if an adequate
20 guidance document can be created and is dictated in
21 the regulation, then that would be a reasonable
22 approach. However, at this time I am not aware of
23 what the contents of that would be so I would be
24 reluctant to say that that is the way to go in the
25 absence of a standard that we could ponder and

1 consider the adequacy of.

2 MS. BROWN: thank you.

3 DR. WHALEN: Other questions?

4 [N response]

5 Thank you. Next, we will hear from Ms.
6 Bobak.

7 MS. BOBAK: First, I would like to thank
8 you for giving me the opportunity to speak here,
9 and having the ability to have an impact during
10 this panel discussion on reclassification of
11 absorbable hemostatic agents and dressings.

12 My name is Lone Ronnemoes Bobak. I am
13 representing Ferrosan. Ferrosan is a Danish
14 medical device manufacturing company, and we have
15 given the distribution rights of our absorbable
16 gelatin sponge Surgifoam to Ethicon and they
17 distribute our absorbable gelatin sponges.

18 What I will be talking about is the
19 current regulatory status because it is different
20 in Europe and in U.S. I will talk about the
21 essential quality control elements which we have
22 implemented. I will talk about the usage of
23 Surgifoam in critical surgical procedures and give
24 my conclusion.

25 The current regulatory status in EU is

1 that it has a long product history of safety and
2 effectiveness. We have had the product on the
3 market for more than 40 years and, due to some
4 changes in regulation, it was dropped prior to '96
5 and then the CE regulation came and we got it
6 classified as CE Class III medical device in
7 December of 1996.

8 Right now there is discussion in EU about
9 reclassification borderline products into drugs,
10 meaning that products like ours might have an
11 impact and could become a drug again.

12 The current regulatory status in the
13 United States is the fact that Ferrosan, in March
14 '97, submitted an IDE, and over the next years we
15 performed extensive clinical trials on humans, and
16 they were multicenter trials. In 1999, based on
17 the outcome of these clinical trials, we submitted
18 the PMA to FDA, had the inspection in August of '99
19 and got the license to export products to the
20 United States in September of '99.

21 Since 1999 we have submitted an amendment
22 to this PMA several times. The new license has
23 been granted based upon clinical studies on humans
24 as well as clinical trials on animals and, of
25 course, by use of design controls and risk

1 management.

2 Ferrostan has achieved a commission also
3 for usage during neurological surgery. For
4 achieving that authorization, we got some
5 preclinical data as well as clinical trials and
6 filed for getting information for use in our
7 logical procedures. As you are, of course, aware
8 elevated sensitivity to toxic and infective agents
9 and the gelatin could be an infectious agent if not
10 treated right. About toxicity, I am also talking
11 about endotoxin testing which we are performing
12 both on raw materials, as well as the finished
13 products. During neurological surgery there is a
14 potential for physical damage and there are fewer
15 choices for adjunctive hemostasis.

16 The surgical product consists of gelatin,
17 water and nitrogen. Surgifoam is a very safe
18 product and we have had no MDRs in the years that
19 we have supplied product to the U.S. market. But
20 this is only caused by the additional quality
21 control measurements that we have had since the USP
22 standards that are right now in place don't fulfill
23 all the requirements that we feel must be in place.

24 When we are searching for raw materials,
25 and please recall that the gelatin is a very

1 complex biological raw material, we must make sure
2 that our pig skin gelatin is not deteriorated by
3 bovine originated gelatin or by alkaline based
4 gelatin. We have required that the animal
5 supplying the gelatin has been subjected to both
6 pre- and postmortem veterinary controls, and that
7 this is stated in the veterinary certificate
8 accompanying each batch of the gelatin. We have
9 required that the supplier certify that the gelatin
10 is in accordance to an EC standard as regards
11 manufacturing of these materials.

12 We are using a safe but also sensitive
13 method for sterilization of the finished product.
14 When I am talking about sensitive, it is dry
15 sterilization in comparison to sterilization using
16 formaldehyde. We are working on stringent hygienic
17 and very low microbial conditions during the
18 manufacturing of the sponges.

19 As we have had products on the European
20 market for more than 40 years, we feel it is safe,
21 and we felt safe about our products though we did
22 agree with the FDA on the list of requirements put
23 forward when we issued the IDE in 1997. We
24 performed extensive clinical trials in '99 and it
25 seems as if a lot has changed during these two

1 years.

2 We don't feel that the current standards
3 control all the critical elements for gelatin
4 hemostats. So, we don't feel, from a Ferrosan
5 point of view, that reclassification is appropriate
6 as long as there is no proper guidance and
7 controls. Thank you very much.

8 DR. WHALEN: Questions for Ms. Bobak? Dr.
9 Newburger?

10 DR. NEWBURGER: With all of the furor
11 about spongiform encephalopathy that seems to have
12 swept Europe, I assume with the herd and animal
13 controls that is looked at as well?

14 MS. BOBAK: Yes, definitely. But since
15 our gelatin is originated from pig skin, we don't
16 have the BSE impact on the products but, of course,
17 we must make sure that our gelatin has no contact
18 at all to bovine-originated gelatin, and the
19 guidance in Europe is to control that issue which
20 we, of course, follow very stringently.

21 DR. WHALEN: This product has been in
22 clinical use for how long?

23 MS. BOBAK: Forty-eight years.

24 DR. WHALEN: Certainly in my medical
25 school years it was nothing new. I remember seeing

1 it in the operating room and that was in the early
2 '70's. There is a multitude of manufacturers of
3 it?

4 MS. BOBAK: Yes and no. In the U.S. there
5 is manufacturing of a gelfoam product. In Europe
6 there is actually our manufacturer of a sponge,
7 Curacel. Then there are some in China and in India
8 which is sterilized and manufactured in strange
9 ways. So, I won't say many but there are some,
10 yes.

11 DR. WHALEN: And over the several decades
12 that it has been in existence, has there been any
13 scientific advance or manufacturing change that has
14 made a substantial increment in the effectiveness
15 that you perceive?

16 MS. BOBAK: We have changed our
17 manufacturing procedures several times. We have
18 had--how do you say that in English?--something
19 that changed the surface of the product. We had
20 sodium laurel sulphate in the product twenty years
21 ago. We changed that because we found out, doing
22 other precautions during our manufacturing process,
23 that it wasn't necessary to have that. Twenty-five
24 years ago we had formaldehyde as a sterilization
25 agent. We have omitted that completely. So, from

1 a Ferrosan point of view, we have definitely made
2 more stringent our manufacturing process a great
3 deal over the years.

4 DR. MILLER: So you, and I think Dr.
5 Paulson also before you, seem to be not in favor of
6 reclassifying this to a lower classification than
7 it is now. Is that what I am to understand?

8 MS. BOBAK: I would answer yes and no, if
9 I may. I can, of course, see something positive in
10 a declassification but as long as there is no
11 guidance document stating all the additional
12 quality control elements that must be in place, or
13 not having any risks to the consumer, then I am not
14 in favor of a declassification on the sponges, no.

15 DR. WHALEN: Dr. Choti?

16 DR. CHOTI: The comparisons between the
17 two products, they are both hemostatic agents but
18 prepared very differently and they are very
19 different products. Tell us a little bit about the
20 mechanism of action of the gelatin sponge versus
21 the cellulose product we have just heard about.

22 MS. BOBAK: I am sorry, I don't think I
23 would be the right person to answer that question.

24 MS. GORMAN: Hi. My name is Anne Gorman,
25 from Johnson & Johnson. Both products have very

1 similar mechanisms of action in that they act as
2 surfaces on which platelets can be bound, the
3 gelatin being more specific, Surgicel being less
4 specific. Gelatin has a specific binding site for
5 the platelets. With the Surgicel it is more of a
6 physical phenomenon and once platelets are
7 activated you have coagulation activation and clot
8 formation.

9 DR. WHALEN: Thank you, Ms. Bobak.

10 MS. BOBAK: Thank you.

11 DR. WHALEN: Now we will continue with Ms.
12 O'Grady's presentation.

13 MS. O'GRADY: Good afternoon. I am Judy
14 O'Grady, senior vice president of regulatory,
15 quality and clinical affairs for Integra
16 LifeSciences Corporation.

17 I would like to thank the Food and Drug
18 Administration and all the members of the General
19 and Plastic Surgery Devices Panel for allowing me
20 the time to speak at this public advisory committee
21 regarding the reclassification of transitional
22 Class III devices, the absorbable hemostatic agent
23 and dressing devices intended for hemostasis during
24 surgical procedures.

25 I will try to keep within the fifteen