

AT

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGIC HEALTH

NATIONAL MAMMOGRAPHY QUALITY ASSURANCE
ADVISORY COMMITTEE

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

Monday, August 26, 2002

9:00 a.m.

Holiday Inn Gaithersburg
Two Montgomery Village Avenue
Gaithersburg, Maryland

PARTICIPANTS

Maryanne Harvey, Chair
Charles Finder, M.D., Executive Secretary

MEMBERS

James F. Camburn
Nancy J. Ellingson, R.T. (R) (M)
Alisa M. Gilbert
Miles G. Harrison, Jr., M.D.
Jessica W. Henderson, Ph.D.
Debra M. Ikeda, M.D.
Andrew Karellas, Ph.D.
Amy F. Lee, M.D.
Melissa G. Martin, M.S.
Etta D. Pisano, M.D.
Linda S. Pura, R.N., M.P.
Catalina R. Ramos-Hernandez, M.D.
Amy R. Rigsby, R.T. (M)
Donald C. Young, M.D.

C O N T E N T S

	<u>Page No.</u>
Opening Remarks: Maryanne Harvey	4
Conflict of Interest Statement: Charles Finder, M.D.	4
Introductions	7
Alternative Standards Charles Finder, M.D.	10
Roger Burkhart, Ph.D.	11
Open Public Hearing Ken Crocker	23
Pamela Gormley (by Dr. Finder)	33
Overview of MQSA Inspection Findings and Current Inspection Follow-up Actions: Michael Divine, M.S.	35
Good Guidance Practices and Directions for Discussion of the MQSA Guidance under the Final Regulations: Charles Finder, M.D.	48
ACRIN Trial of Full Field Digital Mammography (FFDM): Etta Pisano, M.D.	121
Potential New Applications of FFDM: Daniel Kopans, M.D.	169
Inspection Demonstration Project Update: Charles Gunzburg, M.S.	221
Full Field Digital Mammography Accreditation and Certification Update: Kish Chakrabarti, Ph.D.	226
Priscilla Butler, M.S.	236
States as Certification Agencies Update: Roger Burkhart, Ph.D.	256
Review of Summary Minutes of August 2001 Meeting	263
Future Meetings	264

1 P R O C E E D I N G S

2 **Opening Remarks**

3 MS. HARVEY: Good morning. The National
4 Mammography Quality Assurance Advisory Committee is
5 called to order. We have a very busy schedule
6 today so it will be important for us to keep to our
7 schedule as we go. My name is Maryanne Harvey, for
a those of you who haven't met me as yet.

9 Dr. Finder will now give us our conflict
10 of interest statement.

11 **Conflict of Interest Statement**

12 DR. FINDER: The following announcement
13 addresses conflict-of-interest issues associated
14 with this meeting and is made a part of the record
15 to preclude even the appearance of any impropriety.
16 To determine if any conflict existed, the agency
17 reviewed the submitted agenda and all financial
18 interests recorded by the committee participants.

19 The confl'ict-of-interest statutes prohibit
20 special government employees from participating in
21 matters that could affect their or their employer's
22 financial interest. However, the agency has
23 determined that participation of certain members,
24 the need for whose services outweighs the potential
25 conflict of interest involved is in the best

1 interest of the government.

2 Therefore, waivers permitting full
3 participation in general matters that come before
4 the committee have been granted for certain
5 participants because of their financial involvement
6 with facilities that will be subject to FDA's
7 regulations on mammography quality standards with
8 accrediting, certifying or inspecting bodies, with
9 manufacturers of mammography equipment or with
10 their professional affiliations since these
11 organizations could be affected by the committee's
12 deliberations.

13 These individuals are James Camburn, Nancy
14 Ellingson, Alisa Gilbert, Maryanne Harvey, Melissa
15 Martin, Linda Pura, Amy Rigsby and Drs. Miles
16 Harrison, Jessica Henderson, Catalina
17 Ramos-Hernandez, Debra Ikeda, Andrew Karellas,
18 Daniel Kopans, Amy Lee, Etta Pisano and Donald
19 Young.

20 Copies of the waivers may be obtained from
21 the agency's Freedom of Information Office, Room
22 12A-15, of the Parklawn Building.

23 Several of our members also reported that
24 they received compensation for lectures they have
25 given or will give on mammography-related issues.

1 However, they have affirmed that these lectures
2 were offered because of their expertise in the
3 subject matter and not because of their membership
4 on the committee.

5 We would like to note, for the record,
6 that if any discussion of states as certifying
7 bodies was to take place in any meetings of the
8 committee, it would be a general discussion only.
9 No vote would be taken and no consensus sought.

10 In the interest of getting as many
11 viewpoints as possible, all SGEs, including state
12 employees, would be allowed to participate in the
13 general discussion so that all viewpoints could be
14 heard.

15 In the event that the discussions involve
16 any other matters not already on the agenda in
17 which an **FDA** participant has financial interest,
18 the participants should excuse him or herself from
19 such involvement and the exclusion will be noted
20 for the record.

21 With respect to all other participants, we
22 ask, in the interest of fairness, that all persons
23 making statements or presentations disclose any
24 current or previous financial involvement with
25 accreditation bodies, states doing mammography

1 inspections under contract to FDA, certifying
2 bodies, mobile units, breast implant imaging,
3 consumer complaints and mammography equipment.

4 MS. HARVEY: Thank you, Dr. Finder.

5 Introductions

6 Since we have four new members at this
7 meeting--it is nice to see our new members and to
8 see our returning members from last year--I would
9 ask that each one of us give a very short bio so
10 that we can get more acquainted with each other's
11 experience and background.

12 Dr. Pisano, could I ask you to begin?

13 DR. PISANO: I am Dr. Etta Pisano. I am
14 the Chief of Breast Imaging at the University of
15 North Carolina in Chapel Hill. I am a radiologist.

16 DR. YOUNG: I am Don Young. I am a
17 radiologist, a professor of clinical radiology at
18 the University of Iowa College of Medicine and
19 practice at the hospital and clinics where I direct
20 the breast imaging and diagnostic center.

21 DR. RAMOS-HERNANDEZ: I am Catalina Ramos
22 with the National Breast Cancer Organization. We
23 are a not-for-profit advocacy and counseling
24 services for patients with breast cancer.

25 MS. RIGSBY: I am Amy Rigsby. I am the

1 Technical Director of the Rose Breast Imaging
2 Center in Houston, Texas. I am a mammographer.

3 MS. MARTIN: I am Melissa Martin. I am a
4 medical physicist running a consulting practice in
5 Southern California.

6 DR. IKEDA: I am Debra Ikeda. I am
7 Director of Breast Imaging at Stanford University
8 Medical Center. I am a radiologist.

9 DR. KARELLAS: I am Andrew Karellas. I am
10 a medical physicist. I have been with the
11 University of Massachusetts as of two weeks ago and
12 now I have moved to join the faculty at Emory
13 University in Atlanta.

14 DR. HARRISON: I am Miles Harrison. I am
15 a surgeon by training. I am part of the Sinai
16 Surgical Associates in Baltimore, Maryland and one
17 of the designated breast surgeons at the Lapedes
18 Cancer Center which is a Hopkins affiliate.

19 DR. FINDER: Dr. Charles Finder. I am a
20 radiologist working for the Food and Drug
21 Administration and I am also the Executive
22 Secretary of this committee.

23 MS. HARVEY: I am Maryanne Harvey. I am
24 with the New York State Department of Health. I am
25 a section chief who is responsible for mammography

1 and also the Chairman of this committee.

2 MS. PURA: Good morning. I am Linda Pura.
3 I am one of the clinical nurse coordinators for the
4 Los Angeles County Breast Cancer Early Detection
5 Program under the Department of Health California,
6 Cancer Detection Section. I am also the co-founder
7 and present President of the Los Angeles County
8 Susan G. Komer Breast Cancer Foundation.

9 MR. CAMBURN: I am Jim Camburn. I am
10 Chief of the Radiation Safety Section for the State
11 of Michigan.

12 MS. ELLINGSON: I am Nancy Ellingson. I
13 am a radiologic technologist and mammographer. I
14 am with the American Society of Radiologic
15 Technologists in Albuquerque, New Mexico. We
16 represent about 100,000 radiologic technologist
17 members.

18 MS. GILBERT: I am Alisa Gilbert. I am a
19 seven-year breast cancer survivor. I work with
20 Alaska Natives and American Indians. I am the
21 Director of the National Native Cancer Survivor
22 Support Network.

23 MS. HENDERSON: I am Jessica Henderson. I
24 am an eight-year cancer survivor and I represent
25 the Oregon Breast and Cervical Cancer Coalition.

1 DR. LEE: Hi. I'm Amy Lee. I am
2 Assistant Professor of Community Medicine at
3 Northeastern Ohio University's College of Medicine
4 and Administrative Director for the Master of
5 Public Health Program also located there. I am
6 also a physician consultant for the local Breast
7 and Cervical Cancer Program and, in my former life,
8 before academia, I was OB-GYN.

9 MS. HARVEY: Thank you.

10 I will now ask Dr. Finder to talk to us
11 about alternative standards.

12 **Alternative Standards**

13 DR. FINDER: I just want to give a little
14 bit of background on approval for alternative
15 standards. For those not familiar with this
16 section of the regulations, FDA may approve an
17 alternative to a quality standard that currently
18 exists under Section 900.12 when the agency
19 determines that, one, the proposed alternative
20 standard will be at least as effective in assuring
21 quality mammography as the standard it proposes to
22 replace and, two, the proposed alternative is too
23 limited in its applicability to justify an
24 amendment to the standard or it offers an
25 unexpected benefit to human health that is so great

1 that the time required for amending the standard
2 would represent an unjustifiable risk to human
3 health and also that the granting of the
4 alternative is in keeping with the purposes of the
5 statute.

6 Since our last meeting, the division has
7 approved several alternative standards and these
8 will be discussed by Dr. Roger Burkhart.

9 DR. BURKHART: I might begin by referring
10 you to one of the documents that you were given in
11 preparation for the meeting, the Modifications and
12 Additions to the Policy Help Guidance System, No.
13 5. You will find the new alternative standards
14 included within this document.

15 The first one, the first new one which we
16 approved last September, is found on Page 57 of the
17 document. It is entitled The Manufacturer's
18 Software Modification of the Automatic Exposure
19 Control, but really what it applies to is the
20 testing which has to take place after such
21 modifications occur.

22 Software upgrades or modifications are
23 defined by FDA as being major changes in the system
24 which means that, after they are to take place,
25 there has to be a mammography equipment evaluation

1 conducted of the equipment and any problems that
2 are found have to be corrected before the equipment
3 is used on patients.

4 **Also**, the regulations require that this
5 evaluation and the testing which is involved be
6 done by the physicist on site. But the applicant
7 for this particular alternative requirement made a
8 convincing case that, in this particular situation,
9 we can assure mammography quality if on-site
10 testing is done under conditions of medical
11 physicist oversight.

12 By medical physicist oversight, what we
13 mean is that the physicist has to be consulted, but
14 it is his or her decision as to whether they
15 actually have to come on site to do the testing or
16 whether somebody else can do the testing and send
17 the results to them for evaluation.

18 **As** I said, the applicant for this standard
19 made a convincing case that medical physicist
20 oversight would assure quality in this particular
21 case. **So**, for this specific software modification
22 has given an alternative standard, and when it is
23 applied to the units and the models which are
24 listed in the standard, medical physicist oversight
25 is an option for the facility.

1 The second newly approved alternative
2 standard was reapproved last May as found on page
3 58, and it, too, refers to the testing, the testing
4 conditions or how the testing is done after the
5 modification takes place.

6 Like the last one I just mentioned, it
7 also started out as a request related to a specific
8 software upgrade, but in this case, the
9 justification was that the testing which would take
10 place after this particular upgrade, was the same
11 type of testing which is done routinely by the
12 quality control technologist.

13 *So*, it was felt that if the quality
14 control technologist is qualified to do this as
15 part of the routine QA program, then, they should
16 be able to do it after the software modification,
17 and so it is to be done in conditions of medical
18 physicist oversight.

19 But we got to thinking if it applies to
20 this particular upgrade, it could also apply to any
21 with that particular qualification, so we do have
22 the authority to expand on requests, and we did
23 that.

24 We applied this alternative standard to
25 any upgrade or modification of the computer

1 software in which the testing afterwards is
2 normally done by the quality control technologist.

3 *So*, if the manufacturer feels that the
4 modification meets this particular standard, they
5 need to consult with *us*, confirm with *us* that we
6 agree with this, and if we do, then they can go
7 ahead and this testing after this modification can
8 be done under conditions of medical physicist
9 oversight.

10 The third new alternative standard, which
11 begins on page 60, is different in nature. It
12 deals with the quality assurance program including
13 the quality control testing, full field digital
14 mammography units, again to be more specific, it
15 deals with the time period for which the corrective
16 actions can be taken if the testing reveals that
17 there are problems with the system.

18 Now, in the case of screen-film systems,
19 we were able, at the time of the regulations, to
20 define two classes of test failures. There are
21 those test failures which are significant enough
22 that the problem causing them has to be corrected
23 before the piece of equipment is used again on
24 patients.

25 Then, there are those which for various

1 reasons mostly dealing with the fact that there are
2 compensating methods, we can allow more time for
3 the correction of the test failures for the problem
4 causing the test failures, and the time which was
5 set for that was 30 days.

6 For screen-film systems, as I say, we can
7 make that distinction, but for full field digital
8 mammography systems, which were still in the
9 research stage at the time the regulations were
10 being developed, we couldn't make such a
11 distinction.

12 So, for **full** field digital mammography
13 systems, and any other new modality that might
14 appear in the coming years, we took a conservative
15 public health safety position. We said that any
16 quality control test that has failed, the problem
17 causing the failure has to be corrected before the
18 equipment can be used on patients.

19 So, this **is** what the regulations say, and,
20 of course, the expectation is that eventually, if a
21 technology like full field digital mammography,
22 once it becomes fully established, eventually, we
23 would be able to rewrite the regulations and make a
24 similar distinction as we did with screen-film.

25 But to begin with then, we took this

1 position, but at the same time, we tried to make it
2 clear to the manufacturers and also to any facility
3 that has a full field digital mammography unit that
4 if they feel that there are some of the tests that
5 they can make a case that a 30-day correction
6 period could be allowed for test failures of that
7 test, they could always apply for an alternative
8 requirement.

9 About two months ago, General Electric
10 made such an application for their Senographe
11 2000D, full field digital mammography unit, and we
12 approved that alternative on July 14th.

13 Basically, what this alternative does is
14 divide the quality control tests of the 2000D into
15 three general groups, and an important one as far
16 as terms of the changes which are involved is the
17 third group or the group which is labeled with the
18 letter C.

19 These were tests, quality control tests of
20 the 2000D, which were equivalent to quality control
21 tests of screen-film systems for which the 30-day
22 correction period was already allowed.

23 Some of them involved the exact same
24 testing methods using the exact same action limits
25 and even testing components which were identical to

1 the full field digital mammography system as in
2 some of the screen-film systems, so they were than
3 equivalent in that case, they were virtually
4 identical.

5 The other tests were not quite as
6 identical, but basically, for most of them, the
7 only difference is that the measurements are done
8 off of digital images rather than off of films as
9 it would be in the case in screen-film systems.

10 So, it seemed logical then that if these
11 tests already were allowed a 30-day correction
12 period when problems are found, then, it would be
13 logical to allow the same 30-day correction period
14 for them with full field digital mammography unit,
15 with the 2000D.

16 During this 30-day correction
17 period--perhaps I didn't make this quite clear
18 earlier--during this 30-day correction period, the
19 facility continued to use the part of the system
20 which failed the test, they wouldn't have to take
21 it out of service. a

22 The other two groups are the remaining
23 tests, and these are still tests which, if they are
24 failed, the component which failed the test has to
25 be taken out of service until the problem is

1 corrected.

2 The only reason really for dividing them
3 into two groups is to emphasize something which
4 also exists in the screen-film area, and that is,
5 it is only the unit or the part of the unit which
6 fails the test which has to be taken out of
7 service.

8 That may mean the whole system in some
9 cases, but it may mean only part **of** others. In the
10 case of the 2000D, the A Group tests are tests that
11 the image acquisition part of the system, and so if
12 these tests are failed, the facility may have to
13 stop acquiring new images until the problem is
14 corrected, but as long as the B Group of tests are
15 passed which relate to the interpretation of
16 images, they can continue to interpret old images.

17 The opposite is true, if a test in the B
18 Group is failed, the facility may have to stop with
19 the interpretation of images, but as long as the A
20 Group of tests are passed, they can continue to
21 acquire and store images for each interpretation.

22 **So**, then this alternative requirement
23 applies only to the Senographe 2000D. As I
24 mentioned, other manufacturers and facilities which
25 own full field digital mammography units or the

1 manufacturers, have the option of applying for a
2 similar alternative requirement and we will
3 consider it and determine whether or not it can be
4 accepted.

5 These, then, are the three new alternative
6 requirements that we have accepted since the last
7 NMQAAC meeting. I would be happy to answer or try
8 to answer any questions that you might have on
9 them.

10 DR. PISANO: What do the other
11 manufacturers have to do, do they have to apply for
12 the same thing **to** get, because obviously, this only
13 applied to **GE**?

14 DR. BURKHART: It only applies to **GE** as
15 you would expect as **GE** was the applicant. They
16 would have to go through a similar process, and
17 their lists **of** tests, depending upon their system,
18 will be different. They might have different
19 tests, more or less, and are different, that could
20 be given a 30-day correction period.

21 DR. PISANO: Just as a comment, some of
22 the things listed under Item C were probably added
23 to the QC manuals **of** **GE** at least, and the other
24 companies, because of the MQSA requirements for
25 film. **So**, my prediction, I actually have two other

1 units myself, I have **GE** Fischer and a Fuji system
2 at UNC, is that all these tests listed under Item C
3 are also required by the other manufacturers.

4 So, it might expedite things a bit if a
5 more general statement could be made about these
6 particular tests rather than every company, I don't
7 know what the process is, but we would be grateful,
8 those of us who are using the equipment, if there
9 could be a more expedited process for this.

10 I hate to shut down a machine just because
11 one of these things that wouldn't shut down my **GE**
12 unit, you know, if my Fischer unit had one of these
13 problems, you know, before Fischer had the chance
14 to go in and apply for it.

15 So, it would be nice if the FDA could make
16 a more general statement about this, just as a user
17 of the device.

18 MS. HARVEY: One more reminder. Before we
19 speak, let us give our names for the record.

20 DR. PISANO: Oh, sorry. I am Etta Pisano
21 from UNC/Chapel Hill. It is just a suggestion of
22 making this more general.

23 DR. BURKHART: Roger Burkhart again. We
24 thought about that at the time we were looking at
25 the **GE** application. As I mentioned, we did expand

1 the second alternative requirement to be broader
2 than the original.

3 We decided at this stage, though, that
4 since there is still I guess not much of a
5 consensus in general on the testing with the
6 different models, that it would be best to take
7 each one in turn, but once the ground is broken,
8 obviously, you know, it makes it much easier for
9 the people coming along behind if the tests are
10 really the same, it would only be for any unique
11 things that they might have that we have to look at
12 further.

13 The process I might mention, the process
14 for looking at alternative requirements, actually
15 also is described in a general way, and, in fact,
16 it is right after the listing of the three
17 alternative requirements, the new ones, on page 62,
18 and it basically involves when a request comes in,
19 a staff member is assigned to evaluate it.

20 Sometimes we see right upfront that more
21 information is needed, **so** the staff member will ask
22 for it, but if it looks fairly complete, then, we
23 form a committee to look at it and evaluate it.

24 The committee tries to come up with a
25 consensus, a recommendation being either to approve

1 or disapprove or to ask more information, and then
2 it goes to the branch chief of the Accreditation
3 and Certification Branch, which is the branch
4 responsible for this, to agree or disagree, and
5 then it goes to the division director for final
6 decision.

7 So, the process is not really complicated
8 and, as I said, once the ground has been broken in
9 that area, it can go fairly fast.

10 The other point I guess I should mention,
11 too, is that a facility can apply, as well as a
12 manufacturer, so you do not have to wait for the
13 manufacturer to take action to make the case.

14 MS. HARVEY: Any comments from our health
15 physicists? Ms. Martin or Dr. Karellas.

16 MS. MARTIN: These standards are fine with
17 me.

18 DR. KARELLAS: I agree with Melissa. The
19 only thing I am a little apprehensive is that
20 physicists just have to watch closer since
21 different manufacturers have different
22 requirements, so which is okay, that they can do
23 it. Physicists practicing out there are very much
24 into that. It just will add a little bit on their
25 time for it. That is my only concern, but I am sure

1 they can do it.

2 MS. HARVEY: Thank you.

3 DR. BURKHART: Thank you.

4 **Open Public Hearing**

5 MS. HARVEY: Now, we are moving into the
6 open public hearing aspect of our meeting today.

7 We have comments on quality control for
8 full field digital mammography from Ken Crocker,
9 who is Director of Marketing, Product Planning for
10 the Fischer Imaging Corporation. Welcome.

11 MR. CROCKER: Thank you. My name is Ken
12 Crocker and I very much appreciate this opportunity
13 to address the committee on what I think is
14 becoming a more important topic as time goes on.

15 As you all are probably becoming more
16 aware, digital mammography has reached a greater
17 degree of acceptance throughout the U.S., and I
18 think is becoming more of a standard operating
19 practice with probably well over 300 systems
20 already in place throughout the country.

21 [Slide.1

22 There is I think a few issues developing
23 as this acceptance has increased, and I thought it
24 would be important to bring to the attention of
25 this committee some of the issues that I think

1 apply to not only the manufacturers, but to the
2 accrediting bodies, the **FDA**, and, of course, most
3 importantly, the actual users of the equipment, as
4 well.

5 We are in a situation right now that is
6 rather unique, because when the regulations were
7 really originally developed, of course, digital
8 mammography was basically only a gleam in the eye
9 of most people, but, in fact, we have reached that
10 point that things need to happen now for it.

11 Because of that, the original regulation
12 only stated that the Operator's Manual from the
13 manufacturer should be followed as the appropriate
14 quality control procedure. That is going to result
15 in delay in providing oversight to users of the
16 full field digital mammography system.

17 While it may be true that the **FDA**, as part
18 of the **PMA** process, does review the Operator's
19 Manual, and, of course, I think they do a fine job
20 of that because they are the same people that are
21 reviewing the proposed quality standards that
22 certifying and accrediting bodies would propose, so
23 I think we should be confident that they are doing
24 a good job in that area, but nonetheless, it
25 provides only a limited amount of oversight.

1 Once the **PMA** has been approved, the manual
2 is in use, FDA will look at quality control charts
3 after six months of use from a facility, but that
4 is kind of the end of the process right there, and
5 I think we want to get to a point where we have
6 uniformity and standards.

7 So, basically, I will show you what the
8 proposal is, but the issues today, there is a lack
9 of uniformity because you do have each manufacturer
10 proposing their quality standard, and there is
11 limited oversight because it is primarily that
12 review that happens as part of the **PMA** process.

13 Dr. Burkhart described the approval of
14 alternate standards process. Certainly, you know,
15 I think it has its place, but overall it would only
16 be a stopgap measure in this particular instance
17 since we are looking at a complete new set of
18 standards for digital mammography, and I don't
19 think we could rely strictly on that to address all
20 of the needs.

21 [Slide.]

22 So, why does this issue linger?
23 Obviously, hopefully, it hasn't lingered too long,
24 but right now for us, as manufacturers, there is
25 not a tremendous amount of incentive to

1 standardize.

2 Once we get through that process, which is
3 challenging unto itself, we feel pretty confident
4 that we have produced a reasonable quality control
5 approach, and unless there is some really undue
6 needs, we would prefer to just keep things running,
7 because we want to be able to meet all of our
8 customers' needs.

9 The MQSA Reauthorization Act of 2002 is
10 not going to substantially change the landscape, I
11 don't believe, at least in the last versions I have
12 seen of it, that it addressed any of the issues
13 that we are talking about here related to full
14 field digital mammography.

15 The approval of alternative standards
16 doesn't address the needs of accrediting bodies, as
17 well. The accrediting bodies need to be able to
18 get more involved with this process of controlling
19 digital mammography.

20 [Slide.]

21 The proposal is to charge the FDA and
22 accrediting bodies with development of these
23 uniform standards, and to encourage their
24 cooperation. I know we will be hearing from Dr.
25 Chakrabarti this afternoon on that, and I know

1 there have been starts into this area, but I think
2 there needs to be more urgency applied given how
3 rapidly the acceptance is taking place with digital
4 mammography.

5 We need to charge the FDA and the
6 accrediting bodies to seek guidance from industry.
7 Certainly, individual manufacturers are interested
8 and willing to participate, and then NEMA, which,
9 of course, is an industry group, also has started
10 efforts in that area, as well, and would be more
11 than happy to participate in doing that.

12 Our third point would be that FDA should
13 allow accrediting bodies to accept the
14 Manufacturers Manual for accreditation on an
15 interim basis, and this would allow the transfer of
16 responsibility over to the accrediting bodies, so
17 that they could start becoming involved in this
18 process. As of today, you know, they don't really
19 have any sort of regulatory capabilities in that
20 area.

21 [Slide.]

22 Charging the accrediting bodies with
23 simplification and uniformity as a longer term
24 goal, would also then follow with that, but at
25 least we would have that immediate knowledge that

1 the Operators Manual is available and accepted by
2 the accrediting bodies and the **FDA**.

3 Eliminate the requirement to maintain a
4 film-screen system. **As I** mentioned, the review of
5 the manuals is already carried out by a very
6 respected group of individuals at the **FDA**, that are
7 also responsible for the oversight of film-screen
8 certifying and accrediting efforts.

9 If they are capable of doing that, I would
10 expect they would be applying the same degree of
11 rigor to digital mammography, and believe, in fact,
12 that they have, so this requirement of maintaining
13 a film-screen system, I think imposes an
14 unnecessary burden both on users of the systems, as
15 well as the entire community.

16 Lastly, getting to Dr. Pisano's point
17 about why don't we just do this for all
18 manufacturers, this synchronizing of tests that
19 allow 30 days for corrective action, certainly, we
20 agree and support that position.

21 Tests that are as common as repeat
22 analysis, I know from our standpoint we tried to
23 make things as simple, straightforward, and
24 consistent as possible for users of our system so
25 that they wouldn't have to completely rethink the

1 process of how they do a repeat analysis, for
2 example

3 I think that the FDA should be charged
4 with developing that uniform standard rather than
5 waiting for either a user or a manufacturer to come
6 to them. We actually did go to the FDA over a year
7 ago with a similar request to what had been made
8 for the GE system, but it is a rather challenging
9 process to get through.

10 I mean I am not blaming anybody for what
11 happened there, but I would just say that from my
12 standpoint, I think it would be better for all of
13 us if the FDA could take a proactive stance rather
14 than a reactive stance to requests for these kinds
15 of changes.

16 That is all I have prepared. I just want
17 to thank you again for the opportunity to address
18 the committee and appreciate your efforts to
19 provide the best quality mammography for our
20 community.

21 MS. HARVEY: Thank you, Mr. Crocker.

22 Any questions? Yes.

23 MR. CAMBURN: Jim Camburn from the State
24 of Michigan. I think I have one question for you
25 related to one of the things you commented on,

1 eliminating the requirement to maintain a
2 film-screen system.

3 Are you suggesting, then, that the
4 facility would not have any film-screen unit at the
5 site where they would have a **full** field digital
6 machine?

7 MR. CROCKER: Yes. Let me explain a
8 little bit about what the requirement is today.
9 The requirement today is that there is at least one
10 film-screen system within a particular FDA
11 jurisdiction that is under the supervision of a
12 particular radiologist who has responsibility.

13 They don't have to be at physically the
14 same location, you could have one at a hospital and
15 then you could have a digital at an off-site
16 facility, and as long as there was one film-screen
17 present at one of those two locations under the
18 jurisdiction of a particular supervising
19 radiologist, that would be acceptable.

20 But we see it all the time now that there
21 are situations where different groups of physicians
22 want to become involved with digital mammography,
23 they have the experience, they are willing to do
24 the quality control that has been approved and
25 recommended by the FDA, but they do not want to

1 invest in having a film-screen system, as well.

2 In fact, they have no intentions of using
3 the film-screen system, but because of the way the
4 regulations are today, they will go out, they will
5 buy a film-screen system, they will do the absolute
6 minimum to maintain the accreditation of that
7 system or certification of that system, and
8 therefore, I don't think it is really accomplishing
9 much.

10 I think we are better off letting them
11 focus on the quality control of the digital system
12 that they really intend to use.

13 MR. CAMBURN: We see this from maybe a
14 slightly different perspective because we have a
15 number of facilities that have one digital unit and
16 one film-screen unit, and they seem to use them
17 differently, at least some facilities do. The
18 digital full field mammography machine has a
19 relatively small image receptor compared to the
20 larger film size that you can get with film-screen
21 imaging.

22 What they do, they will--average size
23 patients might fit fine with the digital image
24 receptor, but larger patients would require two
25 exposures for each projection, and it kind of

1 doubles the patient dose in the area that the x-ray
2 beam overlaps.

3 **So**, from a radiation dose point of view,
4 isn't it better to have the ability to do both
5 types of imaging?

6 MR. CROCKER: I certainly appreciate and
7 agree with what you are saying. With the full
8 field digital mammography from Fischer, it has a
9 larger field of view, and therefore, in fact, the
10 larger field of view is 21 by 29 cm, so the percent
11 of the population that would require a multiple
12 stitching together of images is no greater than
13 what would be required under a film-screen system.

14 **So**, for our particular equipment, we don't
15 see that problem, but I certainly can understand
16 where you might be concerned about that from a
17 radiation dose standpoint with some other systems
18 that are available in the marketplace.

19 MS. HARVEY: Any other questions?

20 DR. PISANO: I just have a follow-up
21 comment. I actually think from a public health
22 viewpoint, in terms of getting digital out to
23 remote areas where the images could be beamed back
24 to a central site for interpretation, it makes a
25 lot of sense to not require the film mammography,

1 because you are probably limiting access to remote
2 areas, if they are going to use digital, they also
3 have'to have a film mammogram unit.

4 **So**, I would agree with his comments that
5 he made, that we would like to move this process
6 along, I would like to see it moved along, so that
7 it is more standard and that the film mammography
8 isn't required.

9 MS. HARVEY: Thank you.

10 Dr. Finder.

11 DR. **FINDER**: We have one comment that came
12 in, and the person who submitted it would like me
13 to read it into the record. It is a written
14 statement from Pamela Gormley, who is a mammography
15 supervisor at Epic Imaging in Oregon.

16 Her statement is as follows:

17 The following is a mammography item that I
18 believe the FDA needs to expedite the changes on.
19 We have had two of the new **FFDM** GE 2000D digital
20 mammography units since October 2000. However, the
21 FDA says we still have to have a film-screen unit
22 on the premises, plugged in and ready to use, even
23 though a film-screen is outdated technology.

24 This is approaching two years. This is
25 wasting both their resources and space for

1 mammography. We would have replaced that unit with
2 newer technology if the FDA allowed it. We also
3 have to maintain a film processor that we don't
4 want or need.

5 All of the quality control tests that we
6 do on the digital units show that they have much
7 better detail on the phantom image and on patients,
8 with one-third less radiation per view than the
9 former state-of-the-art film-screen system that we
10 have.

11 Our film-screen combo is the detailed Fuji
12 AD-M fine screen with Fuji AD-M film, dedicated
13 Kodak M-35A processor with White Mountain
14 chemicals, 135-second processing at 95 degrees,
15 using GE Senographe DMR bi-metal tube mammogram
16 machine maintained by **GE** service, but it cannot
17 begin to compare with what we see with the digital
18 system.

19 Please get this changed immediately, so
20 that we can provide the best medical care to our
21 patients without wasting money.

22 MS. HARVEY: Thank you.

23 DR. FINDER: I would like to add that we
24 are going to have some more talk about this entire
25 issue later on in the afternoon.

1 MS. HARVEY: We are a little ahead of
2 schedule, so if Michael Divine is prepared, we will
3 move on to the open committee discussion.

4 Michael is going to talk to us on Overview
5 of MQSA Inspection Findings and Current Inspection
6 Follow-up Actions.

7 **Overview of MQSA Inspection Findings and**
8 **Current Inspection Follow-up Actions**

9 MR. DIVINE: My name is Mike Divine and I
10 work in the Inspection and Compliance Branch in the
11 Division of Mammography Quality and Radiation
12 Programs.

13 [Slide.]

14 My talk today is, appropriately enough, on
15 inspections and compliance.

16 [Slide.]

17 I will be going over a summary of problems
18 that we have found during our annual inspections
19 and also an overview of the various actions FDA
20 might take when facilities have serious problems or
21 failed to correct these problems.

22 [Slide.]

23 For the inspection data, my talk will
24 cover the last two complete fiscal years for FDA
25 data plus most of the current fiscal year which

1 will end on September 30th. While the data for
2 this year is not complete, I think we have enough
3 data for comparison purposes.

4 [Slide.]

5 While most people here today are probably
6 familiar with our inspection levels, I thought a
7 slide was needed for those who might not be
8 familiar with them.

9 Level 1 is the most serious and could
10 result in **FDA** action if not corrected.

11 Level 2 is less serious, but still
12 significant enough that a facility is required to
13 respond to **FDA** with their corrective action.

14 Level 3 findings are considered minor.

15 [Slide.]

16 As you can see from this first slide,
17 facilities continue to improve and the overall rate
18 of problems has been declining, which is very good
19 news.

20 'While this slide only shows two full
21 fiscal years plus most of a third, if we extended
22 these data back to **1995**, when we started
23 inspections, the trend would be even more
24 pronounced.

25 [Slide.]

1 This slide shows Level 1 problems with
2 personnel. While the chart shows a jump for some
3 categories in 2001, the small numbers compared to
4 the overall percentage of inspections doesn't
5 indicate that this is a real problem. For the
6 medical physicist, the number of violations has all
7 but vanished.

8 I would mention at this point that this
9 data represents inspections **of** approximately 9,500
10 facilities.

11 [Slide.]

12 Processor QC problems continue to be a
13 source of problems, but these numbers are also
14 going down. The same is true for missing phantom
15 QC data.

16 [Slide.]

17 On this slide, as opposed to the previous
18 slides which showed data for the facility QC
19 testing, these data come from our inspector
20 testing. The number of violations for phantom
21 image is very **small**, as are data for processor
22 speed. Fog values are somewhat higher, although
23 these numbers have been declining.

24 [Slide.]

25 For the medical physicist surveys and

1 mammography equipment evaluations, the most common
2 problems are overdue surveys, surveys missing
3 specific tests or data, and failures to do
4 evaluations on x-ray units and processors. These
5 numbers are also declining.

6 [Slide.]

7 For interpreting physicians, the number of
8 facilities cited for initial training or experience
9 remains low. For continuing education and
10 experience, the numbers are greater, but are
11 declining.

12 [Slide.]

13 For radiologic technologists, we see a
14 similar trend. The jump in continuing experience
15 is probably due to this first thing checked in
16 2001. We expect these numbers will decline in the
17 years to come.

18 [Slide.]

19 For medical physicists, the small number
20 of facilities cited exaggerates the difference
21 between the years. The missing bars for 2002,
22 under the initial requirements, is due to no
23 facilities being cited, as with the technologists,
24 checking on the continuing experience is a
25 relatively new assessment.

1 [Slide.3]

2 For medical records, a lack of an
3 appropriate assessment category on mammography
4 reports dominates the problems, however, we have
5 seen a substantial drop in the numbers in just
6 three years.

7 [Slide.3]

8 This chart shows some other requirements
9 we check during inspections. The problems with
10 x-ray units has dropped to almost nothing. For our
11 first inspections with complaint and infection
12 control procedures, the drop in the number of
13 facilities with these problems has dropped
14 dramatically.

15 [Slide.3]

16 For the medical outcomes audit, only a
17 small number of facilities still have problems.
18 The last three sets **of** bars here reflect
19 requirements only being checked in the last two
20 years. **As** with some of the other cases like this,
21 we expect these numbers will **go** down with time.

22 [Slide.3]

23 This last slide from our inspection data
24 shows a number **of** facilities that had at least one
25 problem during their inspection for not having

1 complete documentation for their personnel.

2 [Slide.]

3 Moving away from the inspection, the next
4 few slides will focus on the various options **FDA**
5 has when facilities have continuing problems
6 complying with our regulations.

7 [Slide.]

8 These types of actions include a follow-up
9 inspection, additional mammography review, patient
10 and physician notification, which is actually a
11 follow-up action in case the additional mammography
12 review shows problems, a directed plan of
13 correction, civil money penalties, suspension or
14 revocation of a facility's certificate, an
15 injunction, which is actually a court order that
16 would shut the facility down.

17 [Slide.]

18 When facilities fail to meet specific
19 requirements, we may need to reinspect the facility
20 to see if it has corrective problems. Most of the
21 time, these inspections only focus on areas where
22 the facility has failed in the past.

23 [Slide.]

24 Additional mammography review is a review
25 of mammograms and/or mammography reports to

1 investigate previous or ongoing clinical problems
2 at the facility. The purpose of the **AMR** is to look
3 for serious problems where patients and physicians
4 need to be notified. If there was a serious risk,
5 there could be a possible patient and physician
6 notification.

7 [Slide.]

8 For additional mammography review, we
9 generally select certain types of issues that we
10 think we want to do an **AMR**. One we do which is the
11 most common although it has been significantly
12 declining the last few years is we find a phantom
13 image problem that is at Level 1 during an
14 inspection, we will do an **AMR**.

15 We could do one for an interpreting
16 physician that would fail to be qualified.
17 Clinical image quality problems would be an obvious
18 one. If there was an overall failure in the
19 quality assurance program at the facility, that
20 could trigger one, and we have done a few for
21 fraudulent recordkeeping situations.

22 [Slide.]

23 The extent of an **AMR** could range from a
24 few films to a larger sample. Our most common
25 reason for **AMR**, as I mentioned, is Level 1 phantom

1 failure. A larger sample is usually needed if a
2 smaller AMR shows serious problems or the problems
3 at the facility to make a smaller review
4 inappropriate.

5 [Slide.]

6 When an AMR shows serious problems, FDA
7 would send the facility a letter requiring the
8 patient and physician notification. These letters
9 outline options referring physicians and patients
10 have, such as getting their mammograms reread by
11 another interpreting physician or getting a new
12 mammogram. The letters are written in plain
13 language, avoids using complicated jargon with
14 patients.

15 [Slide.]

16 A directed plan of correction is a
17 regulatory action FDA may take that imposes
18 additional requirements on the facility. The goal
19 of the DPC is to force the facility to perform
20 mammography in compliance and allow FDA to easily
21 monitor this performance.

22 Under a **DPC**, the facility is usually
23 required to send FDA copies of records on a monthly
24 basis and are subject to additional inspections to
25 check on their performance.

1 [Slide.]

2 For more serious problems, FDA may suspend
3 a facility's certificate. Once a certificate has
4 been suspended, the facility can **no** longer perform
5 mammography. In most cases, facilities are usually
6 given a hearing prior to the suspension, however,
7 FDA may suspend prior to a hearing if there is a
8 serious risk to human health or other substantial
9 violations.

10 [Slide.]

11 A last list of the remaining options that
12 FDA has is rather than shutting a facility down,
13 FDA may opt for charging a facility civil money
14 penalties, and this could be up to \$10,000 per
15 violation or per day.

16 We could also revoke a facility's
17 certificate, which is equivalent to suspension,
18 however, once a certificate has been revoked, the
19 owner or operator of the facility cannot own a
20 mammography facility or operate a mammography
21 facility for at least two years after the
22 revocation.

23 [Slide.]

24 Lastly, if everything else fails and we
25 feel that we have to go **to** court, we have the

1 ability to use an injunction which actually closes
2 the facility down through a court order.

3 [Slide.]

4 In closing, this table shows the number of
5 times FDA has taken specific regulatory actions
6 based on problems occurring at facilities.

7 That concludes my talk.

8 MS. HARVEY: Thank you.

9 DR. PISANO: I have a question.

10 MS. HARVEY: Yes, Dr. Pisano.

11 DR. PISANO: Is this in the history of
12 enforcement of MQSA or is this one year?

13 MR. DIVINE: This is the entire program.

14 DR. PISANO: And what is the denominator
15 like 900-plus facilities per year times 8 years,
16 something like that?

17 MR. DIVINE: Well, we inspect about 9,500
18 facilities.

19 DR. PISANO: That is what I meant, 9,500
20 times about 8 or 10 years?

21 MR. DIVINE: Yes.

22 DR. PISANO: So, it is 95,000 facilities,
23 and these are the numbers, something like that, is
24 that right?

25 MR. DIVINE: Yes.

1 DR. FINDER: I would want to add that this
2 represents actions taken by FDA. This does not
3 include actions taken by the State, and in several
4 cases or many cases, the State has taken action
5 before we have, and in that case, we don't pursue
6 it any further, so it is not the total number of
7 facilities that ran into problems.

8 MS. HARVEY: Dr. Karellas.

9 DR. KARELLAS: You mentioned something
10 about that the equipment-associated problems are
11 something like very few or next to nothing, which
12 is very encouraging, but I would like to comment
13 for the public and for the lay press, because often
14 we read about that there are no problems with the
15 equipment or it has nothing to do with the
16 equipment.

17 The reasons that inspectors find very few
18 problems with the equipment is that equipment is
19 very well maintained. We find problems with the
20 equipment all the time routinely. Almost on a
21 weekly or monthly basis, technologists will walk in
22 and will find problems with a processor, on
23 occasion with the automatic exposure control, they
24 typically call service or physicist depending on
25 the situation, and the problems are taken care of.

1 So, this is why you don't see the
2 problems. I am sure you know that, but the public
3 perhaps doesn't understand that.

4 MR. DIVINE: That is a good point. We
5 only go in once a year to do the inspection, and
6 when we look at the equipment, basically, we find
7 that there was a problem, but it has been fixed.
8 It certainly is not something that shows up during
9 the inspection.

10 I would also point out that as the years
11 have gone by, a lot of equipment that had problems
12 and couldn't be maintained has been replaced or
13 repaired to where it can meet the requirements.

14 MS. HARVEY: Dr. Lee.

15 DR. LEE: Amy Lee. I was wondering if you
16 ever analyzed your data for specific trends, like
17 geographical areas that tended to have more
18 violations or specific kinds of equipment, and if
19 you have, have you noted any kinds of trends or
20 clusters.

21 MR. DIVINE: I am not aware if we have
22 done any geographic types of analyses.

23 MS. HARVEY: When we look at the Level 1
24 phantom image violations, would you talk to us a
25 little bit about the scores that that might

1 represent, what is passing and what scores would be
2 considered to be a serious violation, is it
3 triggered immediately after below 10?

4 MR. DIVINE: The criteria we use for
5 phantom image, we have two, Level 1 and Level 2.
6 Level 2 is where it fails at the accreditation
7 body's limit, which all the accreditation bodies
8 use the same values, which are 4-5ers, 3 speck
9 groups or 3 masses. If any of the objects go below
10 any of those, it's at least a Level 2.

11 Now, our criteria for Level 1 is if it
12 goes below 3, 2, or 2, which is one unit below the
13 criteria. So, we do have a certain number of Level
14 2 phantom failures, and those are higher than a
15 Level 1, but even those are not very high.

16 MS. HARVEY: Equipment has become much
17 better at resolution over the years. Is there a
18 debate about raising the image score?

19 MR. DIVINE: I am not aware of one. I
20 have heard some people mention that, but there has
21 been no urge for us or, as far as I know, the
22 accreditation bodies to raise the values, but it is
23 possible that there has been, I am not aware of it.

24 MS. HARVEY: Dr. Young.

25 DR. YOUNG: Don Young. Have you compared

1 your data with the States that are accrediting
2 bodies and certifying bodies that had the data
3 required relative to the inspections?

4 MR. DIVINE: Not that I am aware of.

5 MS. HARVEY: Dr. Young, do you have any
6 data?

7 DR. YOUNG: No, I don't personally.

8 MS. HARVEY: Any other questions,
9 comments?

10 Thank you.

11 MR. DIVINE: Thank you.

12 MS. HARVEY: It's time for a break. It's
13 about 5 minutes of 10:00, perhaps 15 minutes, back
14 at 10 minutes after 10:00. Thank you.

15 [Break.]

16 MS. HARVEY: Dr. Finder will provide us
17 with information on Good Guidance Practices and
18 Directions for Discussions on MQSA Guidance under
19 the Final Regulations.

20 Dr. Finder.

21 **Good Guidance Practices and Directions for**
22 **Discussion of the MQSA Guidance under the**
23 **Final Regulations**

24 DR. FINDER: Before we begin our
25 discussion of final regulation guidance, I would

1 like to briefly explain the procedures that **FDA** is
2 following as it develops new guidance.

3 In response to public comment regarding
4 the use of guidance documents, FDA held an open
5 meeting on April 26, **1996**, and on February 27,
6 **1997**, they published a federal notice outlining the
7 steps the agency needed to take prior to issuing
8 guidance.

9 In brief, it stated the following.
10 Guidance has to be developed in an open manner that
11 permitted input from the general public and the
12 regulated industry. In most cases, new or
13 controversial guidance had to allow for such input
14 prior to its implementation.

15 While the statutes and their associated
16 regulations were binding and enforceable, guidance
17 was to represent a way or ways of meeting the
18 regulations, but other ways would be acceptable as
19 long as they met the requirements of the underlying
20 regulations or statute.

21 Before we begin our discussions, I would
22 like to emphasize the following. We are here to
23 discuss the proposed guidance, not the underlying
24 regulations. The regulations have already gone
25 through their own extensive approval process and

1 while they are subject to future change, the
2 purpose of today's meeting is to address the
3 proposed guidance.

4 When you hear or see words like shall
5 require or must, they refer to the underlying
6 regulation, whereas, the words should, may, or
7 recommend refer to the guidance. I also want to
8 add that since the beginning of the program, we
9 have issued a large amount of guidance to help
10 facilities meet the underlying regulations.

11 This material, this guidance has been
12 compiled into what we call the "policy guidance
13 help system," which is a computerized search engine
14 that is now available on the Internet to aid
15 facilities in their compliance with the
16 regulations.

17 There is probably about anywhere from 5-
18 to 700 pages worth of guidance encapsulated in that
19 search engine and what we are in the process of
20 doing right now is going through all that guidance
21 to update and revise it.

22 One of the documents that you have, which
23 is Modification Document No. 5, is the first in
24 that series where we are actually going page by
25 page through all the previously issued guidance to

1 update and add material as appropriate.

2 So, with that said, I think that probably
3 the first item that we would like to talk about is
4 the issue of the agency automatic exposure control.

5 With that said, I guess we are done with
6 AEC.

7 [Laughter.]

8 DR. FINDER: Let me give a little bit of
9 background. We did send out a letter to the
10 committee for them to look at prior to the meeting,
11 and basically, this raised several issues about
12 testing of the automatic exposure controls in some
13 of the newer equipment that have multiple different
14 configurations and submodes.

15 If anybody would like to start the
16 discussion on that, I would appreciate it,
17 otherwise, we are going to have a lot of time
18 between now and lunch.

19 MS. HARVEY: Yes, Dr. Karellas.

20 DR. KARELLAS: At least I would like to
21 start in one area of the AEC issues. There are
22 certain systems that they may have various modes
23 and medical physicists may be evaluating modes that
24 they may not be actually used by the facility.

25 My own view is that there should be no

1 need to test every available mode of a complex AEC
2 system if the facility does not intend to put it to
3 use, and a facility should decide as to what they
4 use, and that should be tested.

5 Now, I understand that in real life, a
6 facility will start with something and perhaps a
7 month later, they will decide that they need to use
8 another mode, and that will happen. Although it
9 may be not a problem for a physicist to test these
10 two or three modes and have that, but if they are
11 far more complex than that, and there are too many
12 combinations, it may be unrealistic to be testing
13 all these modes.

14 Then, the physicist could come back and
15 reevaluate the system a few months later if that
16 had to be, but I am not suggesting that the
17 physicist should evaluate the AEC every time every
18 minor modification is made, and the way it is used
19 or some very minor repair.

20 I am saying that it should be tested only
21 if there **is** a very substantial departure from what
22 the system was initially tested.

23 MS. HARVEY: Ms. Martin.

24 MS. MARTIN: This is Melissa Martin.

25 I am the other medical physicist on this

1 panel. Obviously, what we are discussing affects
2 what Andrew and I do the most. As a consulting
3 physicist, just to put this in perspective, I
4 obviously provide the medical physics services for,
5 at this point, around 150 facilities, which covers
6 around 250 mammography units on an annual basis.

7 We have many of the high-level,
8 multi-mode, multi-target, multi-filter units in our
9 practice. We have made great strides to test what
10 I have considered all the clinical modes used for
11 each one of these units when we go on site the
12 first time.

13 In Southern California, I cover an area I
14 call Southern California. I cover sites that are
15 about 300 miles away. It is to my benefit and the
16 facility's benefit to make these measurements when
17 I go out there initially. That is why, as Andrew
18 said, I try to cover all what is going to be called
19 clinically useful or possible clinically useful
20 combinations.

21 For those that use, as an example, the GE
22 DMRs, I don't have any facilities that use the dose
23 node on a 2-cm breast, so it makes no sense to
24 require the physicist to test the dose mode for a
25 2-cm breast.

1 i think we do need to set an understanding
2 here of what is clinically useful and what the
3 physicist would be expected to test. Some of the
4 newer units, the low radium 4's, it is not even
5 possible to test at the low kVp's or the high kVp's
6 for thin breasts, because the grid doesn't have
7 time to move.

a So, if you technically looked at some of
9 what is proposed, it can't be. The other factor is
10 how much time are we taking a facility down. We do
11 impact the access to people to get mammography. If
12 we go in to perform our measurements, we are
13 typically in a room somewhere between 4 and 6
14 hours.

15 That is 4 to 6 hours that room is out of
16 service and available for serving patients, and I
17 think we have to be aware and very careful not to
18 set measurement criteria that is not clinically
19 relevant, but which will also add cost to the
20 facility and decrease the amount of time that the
21 patients can be examined.

22 **DR. KARELLAS:** Melissa put that very
23 nicely. I think the vast majority of medical
24 physicists feel that way.

25 **MS. HARVEY:** What percentage of the time

1 that you are testing the equipment do you think the
2 AEC testing would involve? Is it a major part of
3 the testing?

4 MS. MARTIN: At least 25 percent
5 currently.

6 MS. HARVEY: And there are no surrogates,
7 there aren't any simple tests that we can go to
8 that would be representative of larger--

9 MS. MARTIN: I am saying it's 25 percent
10 to do what I have been considering the clinical
11 modes. If one of the discussions was pursued here
12 that I had to test every kVp in every mode, you are
13 adding at least 2 to 3 hours of testing, so you are
14 roughly adding somewhere between 3- and \$500 of
15 additional cost and another 2 to 3 hours of time
16 out of the room.

17 MS. HARVEY: Plus your time.

18 MS. MARTIN: Yes.

19 MS. HARVEY: Any other comments about AEC
20 testing?

21 DR. FINDER: This is Dr. Finder again. I
22 think you framed the issue. Now, we have got to
23 get down to some of the specifics, and we did ask a
24 couple of questions in the document that went out.
25 We would kind of like some guidance on what you

1 think that reasonable testing is under the various
2 scenarios.

3 The first issue that we talked about here
4 is, as we have said, the current requirement is
5 that a mode or configuration needs to be tested
6 prior to clinical use. There are two different
7 areas where that could occur. One is at the
8 initial evaluation called the mammography equipment
9 evaluation, and the other is during the annual
10 physics survey.

11 The requirements in the regulations are
12 slightly different for those two types of testing.
13 One of the issues that has been brought up is, is
14 it enough to test just all the modes and
15 configurations at the initial mammography equipment
16 evaluation, and then do something lesser of not
17 possibly all the clinical modes at the annual
18 survey.

19 This is one of the issues that we would
20 appreciate the committee's guidance on.

21 MS. MARTIN: I will just respond and again
22 this is where I am coming from. I basically use
23 Alternative Test **No. 2**, which conforms to the ACR's
24 suggested forms that are available in the latest QC
25 Manual for the physicists.

1 That test I perform annually on all the
2 modes with every machine. I personally have not
3 skipped any **of** them, so I don't have any feel for
4 what percentage of people typically only test in
5 the contact mode. I always test the mag mode.

6 DR. KARELLAS: I always test the mag mode
7 and most people I know test the mag mode, because
8 we **go** pretty much by the ACR guidelines, so we use
9 that as a guide, and we may make one or two
10 additional measurements for other things that we
11 feel might be necessary. That is what we go by.

12 MS. HARVEY: Does that help you, Dr.
13 Finder? No?

14 DR. FINDER: Well, it partially addresses
15 some of the issues, but we also have the concept of
16 these units that have multiple different AECs. As
17 has been brought up, some of these AECs are not
18 used over the entire 2 to 6 cm range.

19 The regulations, however, say that the
20 AECs have to be tested over that range, and we have
21 some situations where a facility may say, well, we
22 never use the 2 cm range for this type of submode
23 of AEC, but we do use it at 4 and 6. Well, how is
24 that going to be tested?

25 We also have the issue **of** a facility that

1 says they are going to use one submode at 2,
2 another submode at 4, another submode at 6, what is
3 the appropriate testing under those types of
4 scenarios, do you just look at those three
5 individual submodes at those levels, or do you
6 require each one of those submodes to be tested at
7 2, 4, and 6?

8 These are some of the questions that have
9 come up, and how do you deal with those kind of
10 situations.

11 MS. MARTIN: What we have made the
12 choices, if they are clinically using it, in other
13 words, again, go back to the AOP contrast, the GE
14 DMR has three different modes, well, actually
15 several different modes, but three automatic modes
16 AOP, which is contrast, standard, and dose.

17 Typically, what Dr. Finder is saying is a
18 2- or 4-cm breast would be examined in the contrast
19 mode, 6 could be either contrast or standard. I
20 typically test 2 and 4, I don't test 6 under
21 contrast unless the facility says that is what they
22 use. A standard, I typically do 4 and 6, I don't
23 do 2 unless the facility says that is what they
24 use.

25 I really do think you have to look at what

1 the facility is using, but if those modes are set
2 up and cross with each other, that basically is
3 your test. I think it is crucial that each mode be
4 tested at least on some thickness to verify that
5 you are tracking between modes. I don't find it
6 mandatory that you test every thickness on every
7 node.

8 MS. HARVEY: So, if you are testing a unit
9 and you look at either the 2, 4, or 6, and it is
10 working properly, can you presume that it would be
11 working properly at the other two thicknesses that
12 you are not testing, or do you actually have to do
13 each thickness?

14 MS. MARTIN: I have not found a problem.
15 If it is working properly in the mode for which it
16 is pretty much designed to work in, I have not
17 found a problem with it tracking between the other
18 nodes, but certainly I am not the only one
19 performing these measurements.

20 Dr. Karellas.

21 DR. KARELLAS: We actually test at higher
22 thicknesses than 6, but it is not a requirement,
23 but we just do it because we want to see how the
24 machine works, and we just draw a line there, so if
25 it deviates, and sometimes they do, that then we

1 know at least that we have bracketed for the
2 requirement.

3 Usually, the deviation may be very
4 nominally below or above what the requirement is,
5 but we know that the deviation is not really
6 against any regulations, state or federal, but we
7 are aware of it. If we see something and we need
8 to adjust it, at least we know.

9 MS. MARTIN: One of the more crucial items
10 as far as time testing is for those units that are
11 now out there, that have individual detectors,
12 there is one manufacturer that has seven or eight
13 separate detectors, independent detectors.

14 I think the more crucial time thing is I
15 have found it and what I have been doing is doing
16 the full set of tests on one of the detectors and
17 then cross-checking all the other detectors for the
18 4 cm breast. I have found that to be sufficient.
19 I have not been testing all eight detectors for
20 every target filter for every one.

21 I think again, as long as you are making a
22 reasonable attempt to verify at typically the 4 cm
23 thickness, that your detectors cross with each
24 other, that should be considered an acceptable
25 test.

1 **MS. HARVEY:** Dr. Pisano.

2 **DR. PISANO:** Isn't it true--this is a
3 question for the physicists--isn't it true that if
4 there were a problem with one of the **AECs**, for
5 example, in that eight system, it would be obvious
6 in the clinical images?

7 Wouldn't it be that it would be either too
8 light or too dark over a certain region of the
9 breast, **so** that it is not something that is likely
10 to create real clinical problems? In other words,
11 the radiologists, the readers, and the local
12 physicist, if there is one, would be able to spot
13 it very quickly?

14 **MS. MARTIN:** You are going to spot it very
15 quickly, right.

16 **DR. KARELLAS:** It is true that eventually,
17 it will become obvious when an astute radiologist
18 or technologist will discover it, however, the
19 concern is that there will be certain studies that
20 will be done, and the patient will be gone, and
21 images may be suboptimal, so eventually, it will be
22 found, it will not go for very long, but it is just
23 that there may be something compromised, perhaps
24 not of great significance, but certainly
25 mammography may not be done at the state-of-the-art

1 level.

2 MS. HARVEY: I have another question. Do
3 you generally find a problem on the initial
4 testing, or is this a problem that occurs over
5 time, that you see that the detectors go out of
6 whack on your annual when you come back, do you see
7 problems at the annual testing, or do you see more
8 problems at the initial testing?

9 MS. MARTIN: The initial testing is
10 usually good. I think it is absolutely crucial,
11 though, that it be tested, at least sampled, again
12 cross-check some way for each one of those
13 initially.

14 My experience is the installations are
15 usually done very well now. I wouldn't say that
16 was necessarily the case two years ago, being the
17 installers have gotten much better, and I think
18 part of that is because they are trying to make the
19 criteria that has been set.

20 MS. HARVEY: Dr. Karellas.

21 DR. KARELLAS: The question that I have is
22 what we should be doing when we have a brand-new
23 machine and it has all these multiple modes, and we
24 have not discovered how we are going to be using
25 that machine, and the site needs to accept it.

1 If you accept only two or three or four
2 modes, and then six months later, you discover that
3 you want to use all the others, and you find
4 something that doesn't work very well, then, it is
5 more difficult to go back to say that that was not
6 done properly in the beginning, especially post
7 warranty.

8 So, it is somewhat of an issue as to
9 whether we need to test absolutely everything upon
10 installation, but that can be a very frustrating
11 experience because you are testing something that
12 may be so far out of the real application.

13 So, I still maintain that upon
14 installation, the site should define certain modes
15 of use and perhaps one or two or three above and
16 beyond that based on how this machine should be
17 used, and perhaps after the first or second year,
18 you could perhaps narrow it down a little narrower
19 to say that we are never going to use these modes,
20 we are only restricted at just to these three
21 modes.

22 MS. HARVEY: Dr. Pisano.

23 DR. PISANO: I have another question for
24 the physicists. Isn't it the case also that your
25 phantom testing on a weekly basis would suggest

1 :here is an AEC problem, you would see your OD
2 :hanging over time even with the same settings?

3 I am trying to get a feel for how
4 dangerous this is for patients. My sense of this
5 is that it is not very because of these two things,
6 these clinical images will change and the phantom
7 imaging will change, so even if you don't check
8 every mode, even if they use it one or two, you are
9 going to find it somehow. But I would like to hear
10 your comments on that.

11 MS. MARTIN: Again, as long as you
12 cross-check the modes, I don't think you are going
13 to have a problem. I have not found it a necessity
14 to check the complete thickness for every single
15 mode, and, yes, the idea of if you wanted to
16 extrapolate it for those instruments that have
17 eight detectors, do you want eight phantom images
18 every week to verify that you have consistency.
19 You could take it to the nth degree and make the
20 same requirement, and I don't think any of us want
21 to go there.

22 MS. HARVEY: Dr. Karellas.

23 DR. KARELLAS: I agree with Dr. Pisano.
24 First of all, I do not feel it is dangerous and I
25 feel very strongly that weekly phantom tests are

1 very good, and technologists, in fact, in most
2 situations, it is technologists, that they call us
3 about problems, and they are very vigilant about
4 image quality.

5 On the down side is that if something ever
6 might happen that shows okay on the 4.3 cm phantom,
7 but it doesn't track very well when the thickness
8 is 6 cm, that **is** possible, but at least in our
9 experience, most of the time technologists call us
10 and they say there is something wrong with my **AEC**,
11 and interestingly, you may have tested it three
12 months ago and everything was fine, and they
13 alerted us.

14 **So**, what radiologists and technologists
15 see every day is extremely critical.

16 **DR. FINDER:** One point that I would want
17 to bring up in terms of the phantom testing, the
18 phantom is one image taken under the clinical
19 conditions for that thickness of breast, so you
20 will not be checking or necessarily have any idea
21 how the other submodes might operate in terms of
22 the **AEC**.

23 In fact, depending on how a facility sets
24 up its protocols, if it does, their standard
25 patients, say, in the typical **AEC** mode, and not the

1 iull auto mode, you won't have any idea what is
2 going on with the full auto mode in terms of
3 looking at the phantom.

4 But these are the issues that we would
5 like discussed and try to come to some kind of
6 consensus.

7 MS. HARVEY: Mr. Camburn.

8 MR. CAMBURN: Maybe I can just relay some
9 of the information that we get from our inspectors
10 from time to time about this, especially in terms
11 of your fourth question that you ask here about
12 testing submodes 1, 2, and 3, when submode 1 might
13 only be used with a 2 cm thick breast and submode 2
14 with a 4 mm, and so on.

15 What we find from time to time is the
16 technologists will inadvertently use the wrong
17 submode especially if you have a number of
18 technologists working with the same equipment, they
19 don't all seem to be on the same page at the same
20 time, plus there are patients whose size fall
21 between 2, 4, and 6, and the technologist makes a
22 judgment and may sometimes judge to use a different
23 submode than what may have been initially assigned
24 for that thickness.

25 So, we kind of like to see, although a

1 reasonable amount of testing, that the submodes are
2 all tested.

3 MS. MARTIN: That is why I made the
4 comment earlier. I think each of the submodes
5 needs to be tested at least at some thickness to
6 show that it tracks, but if all the submodes are
7 tested--and I would add to Andrew's comment, we
8 always test an 8 cm breast because at least in
9 Southern California, we have several women that
10 fall in that category--so, I think the 8 cm breast
11 is absolutely critical to be tested.

12 Now, what I have found is for many of the
13 8 cm breasts, the 8 cm phantoms, it is necessary to
14 adjust the density on some units to achieve the
15 optimum density, and that option is nice to have,
16 and I think that is part of the physicist's
17 responsibility to give the facility a technique
18 that will bring their large breasts into the same
19 density range as their average breast. That is
20 part of working with the facility.

21 As long as I have tested those modes over
22 some part of the thickness and they all meet the
23 tracking, I have not found a problem with testing
24 every mode for every thickness.

25 MS. HARVEY: Dr. Pisano.

1 DR. PISANO: I guess I want to reiterate
2 Melissa Martin's point about the cost to facilities
3 of having machines down for longer periods of time
4 than are needed.

5 I think the reason why I kept going back
6 to the point of patient safety is I think that
7 maybe one could make a case that at acceptance, all
8 the modes and all the thicknesses should be tested
9 and then maybe after acceptance, then, only the
10 ones that are clinically used should be tested
11 routinely.

12 In that way, I really think you are
13 probably doing the maximum at the beginning and
14 then you are not going to hurt patients or I want
15 to also echo Andrew Karellas' point that the
16 technologists are really right on top of those,
17 when the AEC drifts out of calibration, we know
18 about it pretty quickly, so I don't see that there
19 is a practical real problem.

20 It is more because the regulations say it,
21 you have to figure out how to do it problem. In
22 reality, we are on top of this AEC, and we know
23 when it is not working properly in clinically
24 relevant modes. So, I feel like we should probably
25 try to make it as supportive of the regulations,

1 but not hugely time-consuming for the facilities.

2 DR. KARELLAS: I think we should also
3 realize that the way we are testing it, at the
4 various modes and the various thicknesses, we are
5 using Lucite, which is really a structureless
6 material as far as x-rays are concerned, and the
7 sensors, that is not what they see. They see very
8 inhomogeneous density.

9 So, we can be testing some of these things
10 forever and never really reaching perfection as far
11 as matching it to the anatomy, and we must realize
12 that. There is a point that when we go above and
13 beyond, we get diminishing returns. We just do not
14 get much better image quality.

15 I do not want to de-emphasize the
16 importance of the proper exposure. There is no
17 question that with film-screen, the correct
18 exposure is one of the most critical aspects of a
19 good mammogram, but I think it can be done without
20 going far above what we are doing today on testing
21 the AEC.

22 DR. FINDER: I just have an attempt at
23 clarification here. Suppose a situation occurs
24 where in the beginning, a full testing of all the
25 equipment modes as best as possible was done for

1 the evaluation.

2 Then, the facility decides that they are
3 only going to use, let's say, two of the submodes,
4 contrast and let's say dose, whatever, and then you
5 do your annual survey testing those submodes, but
6 sometime after your survey they decide that they
7 want to use a third submode.

8 In your opinion, would that require you
9 coming back to retest before they could use it on
10 patients?

11 MS. MARTIN: Not if you could actually
12 give a cross-check, and I think that is where you
13 would fall into the medical physicist oversight.
14 What my advice to a facility would be is have the
15 technologist on site shoot a phantom in both modes,
16 and if they cross-check with each other, so you
17 could calculate a dosage, and the dosages are
18 reasonable and the technologist is trained to read
19 out that phantom image, and so is the radiologist,
20 if I get feedback that that is acceptable, that
21 would be fine with me.

22 Frankly, as far as testing all those
23 modes, what you find is for 6 and 8 cm, standard
24
25 only have to test contrast and standard, because

1 after about 6 or 8 cm, they are all going to pull
2 the high kV and high filters anyway, there is not
3 that much difference in them.

4 MS. HARVEY: Dr. Karellas.

5 DR. KARELLAS: This is one point that I am
6 not sure I am in total agreement with Melissa
7 Martin. Perhaps I don't understand or perhaps she
8 has conducted some experiment on cross-modes, and
9 there are some data that we should look at.

10 If, for example, we have tested something
11 all in the contrast auto mode and somebody all of a
12 sudden switches to the dose mode, we do have a very
13 different situation in the equipment, and I am not
14 certain that the system would behave the way that
15 we would want to.

16 Now, I wouldn't be surprised if it does
17 happen with certain machines and when you test it
18 across all modes and somehow everything just clears
19 'through and everything is fine, but if I get a call
20 as a medical physicist, and they ask me is it okay
21 to do that, then, I would have to ask them to
22 conduct the measurement.

23 But then they are performing a measurement
24 that I should be performing, and we run into a gray
25 area although theoretically, it is possible that

1 somebody can send you three, four images under
2 these conditions that you prescribe and measure the
3 optical density and generate a report and say that,
4 yes, it actually does conform versus going there on
5 a visit and checking it **on** your own.

6 MS. HARVEY: We have Question No. 2. I
7 refer you to the document on AEC testing, the
8 second page, at the bottom, which has to do with,
9 "Since some or all of the AEC configurations may
10 share key components or algorithms, is it
11 reasonable to assume that the failure of one
12 configuration immediately makes the other suspect
13 unless the cause of failure in one configuration
14 can be isolated as unique to that mode. In that
15 case, only the manual mode could be used as back-up
16 until repairs have been made.

17 "An example of an isolated configuration
18 failure would be a system that incorporates
19 separate AEC detectors for different image receptor
20 sizes. If one detector fails and can be identified
21 as the cause of failure, then the continued use **of**
22 the AEC with other image receptor would be
23 appropriate."

24 MS. MARTIN: I think that comes back to
25 Dr. Pisano's point a while ago. If you have one of

1 these instruments with eight detectors, and
2 suddenly you find out that one of these is off, it
3 doesn't necessarily mean all of them are off if you
4 can verify that you move it to a different unit and
5 it performs fine, then, obviously, you make a note
6 and post something that says one unit is not
7 usable.

8 The same think would come down to the
9 Siemens unit. You could very well have the large
10 bucky fail or the small bucky fail, but you
11 wouldn't necessarily fail both of them. Obviously,
12 that can happen.

13 Frankly, I don't have any facilities that,
14 as the physicist, I allow them to use the manual
15 mode of exposure. If my AEC failed, they are down.
16 There is no way we use a manual technique for
17 anything, and I think that is your bigger--the idea
18 that you are going to allow screen-film mammography
19 in today's world to be performed with manual
20 techniques is out of date.

21 MS. HARVEY: Dr. Karellas.

22 DR. KARELLAS: I totally agree with Ms.
23 Martin. There is no way that I could think of
24 mammography going on in a facility, going on in an
25 manual mode. If the AEC fails, they are down,

1 period.

2 I: would be interested to know whether
3 anybody really would do that, but we don't, and I
4 think most of my colleagues would not allow that to
5 continue, and I think most technologists would
6 stop.

7 MS,. HARVEY: Dr. Pisano.

8 DR. PISANO: Just to clarify, because I
9 don't understand, may have misunderstood, you are
10 not saying, however, like for the Siemens unit,
11 where there are two separate AECs, if one of them
12 was down, the big image detector, but the smaller,
13 you would still go ahead and allow imaging with the
14 smaller detector?

15 MS. MARTIN: Correct.

16 DR. PISANO: Okay, because that is the way
17 we do it at our place, and it seems appropriate to
18 me.

19 MS. HARVEY: Is this a frequent problem,
20 that sensors go down?

21 MS. MARTIN: No, I don't hear it that
22 often,

23 MS. HARVEY: Dr. Karellas.

24 DR. KARELLAS: We have had several
25 problems with AEC on several units. In some cases,

1 we identify them, and in some cases, the
2 technologists would call us. When we say
3 "frequent," it is not too frequent, but if you have
4 10, 20 or 25 mammography units, one of them, in two
5 or three or four years, something will come up.

6 MS. HARVEY: Rather infrequent.

7 Dr. Pisano.

8 DR. PISANO: This is a question for the
9 physicists. In my experience, the AECs tend to
10 drift a little as opposed to totally failing, so
11 you notice on your phantom that the OD is changing,
12 either going up or going down. That is what you
13 notice as opposed to it's not working at all and
14 the clinical images are really terrible, it's
15 really just a very gradual, is that correct, is
16 that what generally happens?

17 MS. MARTIN: That has certainly been my
18 experience, and that is why you do QC every week,
19 and that is why you track those phantom images, and
20 that is why you have PMIs on the machines, is to
21 bring them back into your desired range.

22 MS. HARVEY: Any other comments on
23 Question No. 2?

24 Actually, I think we have sort of answered
25 No. 3 since we have talked about manual mode.

1 "In the event of AEC failure, the manual
2 node may be used for up to 30 days while the AEC is
3 being repaired." I sense that that is not what the
4 panel is recommending.

5 MS. BUTLER: Could I ask a question from
6 the floor?

7 MS. HARVEY: Certainly. I recognize Ms.
8 Butler from the floor.

9 MS. BUTLER: This is Penny Butler from
10 ACR.

11 I would just like to ask for clarification
12 on the document that was provided. I think it may
13 assist the discussion that is going on.

14 AEC failure, what exactly does that mean?
15 Does that mean that it fails the physicist test or
16 does it mean that it just doesn't work, because I
17 think a clarification on that point may sort of
18 influence how the discussion goes.

19 The other question is what is the
20 definition of manual mode, because in the current
21 guidance that is out there, there was a discussion
22 of if the AEC performance fails one of the
23 performance tests in full auto, it would be
24 appropriate to temporarily use the fixed kVp AEC
25 mode in order to continue operating.

1 I would like FDA's interpretation of how
2 this plays into this discussion.

3 MS. HARVEY: Thank you.

4 DR. FINDER: I guess it plays into the
5 discussion in the sense of as we are trying to get
6 a handle on the fact that these units have multiple
7 AEC modes, and if you can figure out that only one
8 or two or three of these modes are affected by
9 whatever problems is causing it to fail a test or
10 to cause problems, but the other remaining AEC
11 modes are not, then, obviously, you could continue
12 to use those other AEC modes.

13 If all those fail, the regulations do
14 allow a manual technique for up to 30 days. Again,
15 that is the way the regulations were written,
16 taking into account the guidance that were received
17 at the time the regs were written.

18 But the facility certainly has flexibility
19 in terms of if they have a functioning AEC that is
20 within the limits, and they can have confidence
21 that it is, they can use that if their other AEC
22 modes fail, for example, the full auto mode fails
23 in some manner, they could use the fixed kVp AEC
24 mode and continue on that basis.

25 But you do raise a good question of how do

1 you know what it fails, what is the definition of
2 failure. Obviously, there are multiple definitions
3 here, and I think that they each raise their own
4 issues.

5 There is the failure that occurs during
6 the physics testing, and then there is the failure
7 that occurs clinically when somebody suspects that
8 there is a problem and what do they do in that
9 case.

10 Generally, what obviously we would
11 recommend is if they believe that there is a
12 problem, they get their physicist and take a look
13 and see what really is going on, so that they do
14 have a better understanding of what is failing and
15 what isn't.

16 MS. HARVEY: Dr. Karellas.

17 DR. KARELLAS: In our experience, failure
18 may be gradual, as Dr. Pisano described, you see
19 some drifting and the technologists may catch that
20 before anybody else.

21 The other mode of failure is when, on the
22 annual testing, that it does not track with
23 thickness, and we will notice that it is slightly
24 off, and the other mode of failure is when you get
25 a call from the technologists and they tell you

1 that something is just very erratic.

2 This **is** not unusual, that they will tell
3 you that it works the first 30 minutes in the
4 morning, or if I take an exposure 15 minutes after
5 I turn the machine on, it doesn't work very well,
6 and then it sort of behaves somewhat better. That
7 is somewhat of an erratic mode, and the
8 technologists pick it up.

9 On the other part about the AEC, this is
10 an automatic mode, so switching to more
11 conventional AEC, as Ms. Butler indicated, fixed
12 kVp, automatic exposure control, this is a form of
13 an automatic mode versus going all manual.

14 MS. HARVEY: Ms. Martin.

15 MS. MARTIN: I would totally agree with
16 what has been said, that would be the first option,
17 if one of the full auto, auto mode fails, you would
18 go to the next level down, which is the manual
19 section of the kV and target and filter.

20 Again, I would come back to that is why
21 initially, we do check, at least cross-check for
22 the phantom with all the modes and make sure all
23 the modes are functioning properly, which certainly
24 allows the facility that option, and if they lose
25 the auto-auto mode, they can very well use the

1 phototimed mode, and that has already been checked
2 out by the physicist and it is ready to go, and it
3 won't take them down.

4 I don't consider the single auto mode as a
5 manual technique. I was thinking of manual as
6 totally manual where that technologist is setting
7 the exposure.

8 DR. FINDER: That is the correct
9 interpretation of that. AEC mode, the fixed kVp is
10 an AEC mode, it is not a manual mode.

11 MS. MARTIN: Yes.

12 DR. FINDER: Let me also ask this, follow
13 up with this. You do the cross-testing both on the
14 initial equipment evaluation and during each of the
15 annual physics surveys?

16 MS. MARTIN: Yes, I do.

17 DR. FINDER: If a physicist didn't do
18 that, would you say that if they hadn't tested it
19 during the annual physics survey, at least at the 4
20 cm cross-check level, that if the facility wanted
21 to switch to one of those modes and it hadn't been
22 tested, the physicist would have to come back out
23 and do that testing before it could be used
24 clinically?

25 MS. MARTIN: If they don't have the data

1 vailable to calculate a dose and image quality,
2 some type **of** check I would think has to be made. I
3 would think you would have to at least shoot a
4 phantom.

5 Not all facilities have the 2, 4, 6 cm
6 Lucite to test, **so** you are either going to have to
7 have a physicist on-site or some acceptable
8 procedure previously outlined that the physicist is
9 willing to accept.

10 That could be the medical physicist
11 oversight is what I am coming back to. If you have
12 checked it out, **so** that you say in the auto-auto
13 mode typically pull a 26 kV, and your mAs is 143,
14 and you shoot it in phototime at 26, and it shoots
15 145, that should be perfectly acceptable, but you
16 do need that cross-check done before you are going
17 to use it on a patient.

18 DR. FINDER: Just to clarify things, right
19 now the regulations, as interpreted, as written,
20 require that the physicist come on-site. What is
21 being mentioned here is the possibility **of** doing
22 this kind of remotely through physicist oversight,
23 which would be a modification of what we have right
24 now.

25 That is what you are proposing or

1 recommending or suggesting?

2 MS. MARTIN: Yes, I am suggesting.

3 DR. KARELLAS: I don't think this is
4 unreasonable for providing a set of data for the
5 physicists under specific conditions if it is
6 needed remotely to advise the facility on something
7 like that, on cross-checking, however, it raises
8 the question whether the physicist should be doing
9 one back-up mode on the annual inspection, because
10 there are many physicists that they will do only
11 one.

12 They will ask the facility, what do you
13 use, and they will use contrast auto all the time,
14 and they will evaluate the contrast auto, and
15 that's it.

16 If the contrast auto does not work, then,
17 you don't have any data on the other mode, so you
18 don't have a back mode, so the option is to either
19 have an evaluation done there, and you can tell
20 them go ahead, you can switch to the other mode, we
21 have the data and your other modes would work.

22 But the question is how do you know now
23 that the other modes would work? If one mode
24 doesn't, how can you assume that the other modes
25 would work. In all fairness to the patient, we do

1 not know. So, somebody has to do something at that
2 point.

3 Now, in the more real world, a
4 technologist is going to tell me on the other side,
5 I have been using contrast auto all the time, I am
6 not going to switch now to select the kVp. With
7 three or four technologists doing that, they are
8 going to get all confused, so chances are they are
9 going to tell me I am calling Service right now and
10 we are stopping.

11 I think they would be very unwilling to
12 just go and do all kinds of things because they
13 would be afraid that they would be doing the wrong
14 thing.

15 MS. MARTIN: I guess maybe I have got
16 technologists that would have no problems with
17 that. I think it strictly depends on the facility,
18 and I think that has to be part of the medical
19 physicist's understanding and agreement with that
20 facility is when they are going to be called and
21 what they allow the technologists to do.

22 MS. HARVEY: Have we completed?

23 DR. FINDER: I just want to clarify, in
24 the fourth one, where I think we had already gone
25 over this, about the testing of the 2, 4, and 6 cm.

1 I just want to clarify in my own mind the
2 consensus or at least some of the comments were
3 that you would test the 2, 4, and 6, but only at
4 the submodes that were used at those levels. Is
5 that correct?

6 MS. MARTIN: That would be my
7 interpretation. I think that is the suggestion,
8 that is certainly the training the technologists
9 are given when they are given their clinical
10 training, that it is never suggested that they use
11 the dose mode for a 2 cm fatty breast. That is
12 part of any technologist's understanding is the
13 appropriate mode to select for the type of breast
14 being examined.

15 DR. FINDER: Just again to clarify, let's
16 say the contrast is used at the 2, and the standard
17 was used at the 4, and then all of a sudden they
18 wanted to use the standard at the 2, any additional
19 testing required or no?

20 MS. MARTIN: It would depend on whether it
21 is a new unit or a reevaluation unit. If I have
22 checked it at 4, and it crosses its contrast at 2
23 and 4, and the standard crosses at 4, if they want
24 to shoot standard at 2 and the techniques are
25 reasonable, it is probably going to work.

1 It hasn't been my experience that that
2 ould necessarily be a problem.

3 MS. HARVEY: Any other comments?

4 DR. FINDER: I just wanted to check.

5 Anybody from the FDA side have any questions that
6 they would like asked, or any items? Speak now or
7 forever hold your peace.

8 MS. HARVEY: Or any other members of the
9 audience?

10 I refer you to the Modification Document
11 No. 5, Guidance for Industry and FDA. There have
12 been quite a few changes in this document, which
13 has been out for a while now for use by
14 individuals.

15 DR. FINDER: I just want to again bring
16 what this document represents. It basically took
17 the guidance that we had already issued on these
18 items, and what we are trying to do is update and
19 modify what needs to be changed. So, that is why
20 you are looking at a lot of issues that deal with
21 the same type of topics, such as accreditation and
22 certification.

23 If anybody has any comments about the
24 changes, these actually have already been published
25 and are out to the public. They are up on our web

1 site. If anybody has any items that they would
2 like to discuss, now is a good time to do it.

3 I would also mention that there is quite a
4 long, new section dealing with full field digital
5 units and their certification, and explaining what
6 we are doing presently with those types of units.

7 MS. HARVEY: Does anyone see anything?
8 Yes, Dr. Pisano.

9 DR. PISANO: I just wanted to talk a
10 little bit about the digital requirement, the
11 digital pages 32 through 37 or so, 40, I guess.
12 This ties in with what I am going to talk about a
13 little bit this afternoon, so I don't know if you
14 want me to wait and talk then or you want me to
15 talk now.

16 MS. HARVEY: Well, give me a hint.

17 DR. PISANO: There is a currently active
18 clinical trial going on for which I am the PI,
19 called the American College of Radiology Imaging or
20 Digital Mammographic Imaging Screening Trial,
21 otherwise known as D-MIST, and we actually have a
22 fair amount of data at this point about what tests,
23 you know, we have been doing the manufacturers'
24 recommendation as is required under MQSA, under
25 this law, plus we have data on other tests, plus we

1 have been doing it for quite a while. At this
2 point, we have 19 open centers and the trial has
3 been open since October of 2001.

4 So, we have a lot of information, and I
5 would encourage the FDA to move forward on kind of
6 looking at the data that exists and trying to
7 perhaps pare down the requirements over what is in
8 the manufacturers' documents or their user manual,
9 whatever it is, whatever you want to call it.

10 I am concerned about spending a lot of
11 time doing tests because we have always done them
12 on film, they may not be appropriate for digital,
13 and it is just that time is money and if we don't
14 need to do it, we probably shouldn't have to do it.

15 The reason I bring this up, it is going to
16 be presented publicly at RSNA, the Radiologic Site
17 of North American meeting this November, and Martin
18 Yaffe, out of the University of Toronto, is the PI
19 of the quality control piece of the trial.

20 I feel that once this presentation takes
21 place, there is going to be more pressure on FDA to
22 kind of respond and maybe cut down the
23 requirements, so I would like to see us kind of be
24 proactive. I want to echo what Ken Crocker said
25 from Fischer this morning in his public

1 announcement .

2 I would like to see us kind of have more
3 specifics kind of and detailed recommendations for
4)C for digital as opposed to just what the
5 manufacturers recommend. I think there are many
6 reasons. It is in the interest of patients, I
7 think, just because of the amount of time that we
8 spend doing it does cost the facilities money. It
9 is obviously in the interests of the facilities and
10 people like me who run facilities to kind of try to
11 keep the requirements to a minimum.

12 Understandably, they have kind of
13 mushroomed into a big set of requirements because
14 the companies just didn't want to leave anything
15 out that the **FDA** might want them to put in, but I
16 think we now have really kind of--I am not prepared
17 to discuss today in my talk what things should be
18 cut out, but I know that we will have, you know,
19 this thing hardly ever drifts or once in a blue
20 moon drifts or never.

21 So, I just want to encourage the **FDA** to
22 kind of perhaps talk to ACR. I am talking about
23 the American College of Radiology Imaging Network,
24 not the ACR mammography presentation program, now
25 there is two separate entities, about what is

1 available, perhaps even to hear what Martin has to
2 say before the talk in November.

3 DR. FINDER: Yes, we would appreciate any
4 access we could get to that information as soon as
5 possible.

6 DR. PISANO: There is actually a public
7 meeting in Washington, D.C. Actually, it is in
8 Arlington, Virginia, at the Ritz Carlton in
9 Pentagon City in October of the American College of
10 Radiology Imaging Network where I am sure we will
11 all get a glimpse of the information. I can't say
12 who is going to give a polished RSNA presentation
13 at that meeting, but I am sure we are going to hear
14 about this at that meeting.

15 So, I will inquire of the ACR Imaging
16 Network folks if you all can be invited to listen.
17 It is a public meeting, so you certainly are
18 welcome to come to the public sessions and perhaps
19 meet with Martin there and the other physicists.

20 All the physicists for all the sites will
21 be present or at least they are all supposed to be
22 present at that meeting, so I would expect a very
23 useful amount of information, you know, you will be
24 able to get a lot of interesting and useful
25 information at that setting, so I would encourage

1 ou to attend.

2 DR. FINDER: I do want to kind of put this
3 into perspective. The regulations dealing with
4 full field digital, when they were written,
5 obviously, there weren't any digital units that
6 have been approved yet. We took as has been stated
7 before the conservative approach and said that
8 without any data, we would rely on the manufacturer
9 of the equipment to establish a quality control
10 system that would be adequate for their unit.

11 I think the idea has always been that as
12 more information became available and the ability
13 to kind of standardize the quality control for
14 these units was developed, that that is what would
15 happen, but until we get enough data available, it
16 is going to be difficult and as we just heard,
17 there is some data that is going to become
18 available soon and as soon as it is, we are going
19 to certainly want to take a look at it and see if
20 we can progress along that frontier.

21 Another issue that has to be kept in mind
22 is that some of these digital units are quite
23 different from each other and that the quality
24 controls that might apply to one unit may not apply
25 to another. That also is an issue and the ACRIN

1 study will be dealing with a lot of different units
2 and hopefully will be able to provide us with
3 enough information, so that we can start
4 formulating the ideas for a standardized quality
5 control system.

6 I am sure that this committee is going to
7 be directly involved when that information becomes
8 available and guiding us in terms of what we would
9 require.

10 DR. PISANO: Just to follow up on that
11 point, we will have data from D-MIST, as you
12 mentioned, for manufacturers, so we will be able to
13 compare the need for different tests for each
14 manufacturer, a very rigorous quality control
15 program centrally monitored also, which is one of
16 the strong features of it.

17 It is a little stronger than what MQSA
18 does because FDA's inspections are an annual
19 snapshot of what happens. This is literally being
20 monitored by central physicists every week, so we
21 can watch for it because it is so important in the
22 trial to be sure we have the highest quality
23 images, because we don't want people to question
24 our results at the end as has happened in other
25 clinical trials.

1 We want to really be sure, and so this is
2 being monitored very, very carefully by physicists
3 every week, so I feel that we are going to have
4 about as comprehensive data on this topic as you
5 can get.

6 MS. HARVEY: Very good. Thank you.

7 Dr. Karellas.

8 DR. KARELLAS: There is no question that
9 they did some mammography units very different from
10 each other, and it will require different modules
11 for physicists and technologists who are involved
12 in that.

13 I find it somewhat more difficult when I
14 deal with the manual that comes from the
15 manufacturers although it is welcome that they have
16 that, but I believe medical physicists will find it
17 a lot easier if they have to deal with modules from
18 an accrediting body, such as the American College
19 of Radiology, and I believe that it is a lot easier
20 to communicate, a lot easier to ask a question.

21 We can always call ACR for what we want or
22 hopefully, they feel an obligation that they have
23 to do it. Frankly, the companies are excellent and
24 they will give you what you want.

25 I find it somewhat more difficult to have

1 to call each company to find each person, if I have
2 a question or the particular book or the version,
3 the accrediting body provides more of a centralized
4 area that I can deal with, who feel absolutely
5 obligated that they have to provide that
6 information that I need.

7 So, if we work in that direction, and if
8 manufacturers cooperate to provide the information
9 they feel that is very important, that their
10 systems need to be tested on, and that can be
11 worked into a document that parallels the existing
12 ACR guidelines that we have.

13 Let's not forget that quality in
14 mammography, at least from the technical point of
15 view, has been achieved to a large degree because
16 of the uniformity that we are able to achieve in
17 quality assurance through these manuals that
18 radiologists, technologists, and physicists have,
19 and we refer to them all the time.

20 MS. HARVEY: Thank you.

21 Any other areas for comments?

22 DR. FINDER: If nobody else has any
23 questions about the guidance document per se, there
24 are some other guidance issues that have come up, I
25 just wanted to bring them to the committee's

1 attention.

2 MS. HARVEY: Hold one second.

3 Ms. Gilbert.

4 MS. GILBERT: I do have a question. This
5 is Alisa Gilbert.

6 In the documentation here, under Breast
7 Implants, I know that there is not a whole lot
8 written in this, and for a lot of the native
9 patients that I work with, this is a new procedure
10 that is being offered now, but I also know that
11 there is not a whole lot of information when it
12 comes to following screening after that on the
13 remaining breast or even doing mammography with
14 implant or not even implant, but other types of
15 procedures have been done on tram flaps or some
16 information like that.

17 I would like to see some documentation or
18 additional information in this body written on
19 that, as well, on the procedure and protocol that
20 should be developed for that, just to include it.
21 I know that it wasn't even included in this, it
22 just said implant, but there is other procedures
23 that are now being taken that haven't been.

24 Walking into the procedure and asking for
25 follow-up, it is completely an unknown, and I just

1 think that that might be something that might be
2 brought to the attention.

3 DR. PISANO: Can I ask a question? I am a
4 little confused about what you are talking about,
5 and I just want a clarification. Are you talking
6 about women who have had implants placed
7 post-mastectomy for reconstruction purposes?

8 MS. GILBERT: Yes, and the follow-up for
9 that. I know that in here, it just says, on page
10 33, there is a question, "Is there a specific
11 amount of training or number of mammograms of
12 breast implant patients that the technologist must
13 perform under direction supervision prior to
14 performing these studies?"

15 It is requiring that 40 hours initial
16 training for that procedure.

17 DR. PISANO: I just have another question.
18 I am not sure what extra procedure you may be
19 referring to for patients who are
20 post-reconstruction. I am not exactly sure. You
21 said there are new procedures being offered, and I
22 am not sure what you are referring to.

23 MS. GILBERT: I just see the notation here
24 alternative requirements or breast implants.
25 Implants aren't the only procedures that are being

1 lone for women that have received mastectomies, and
2 I was just kind of interested to know if that is
3 going to be something that is going to be another
4 addition to this, for women that do have just like
5 one--maybe I am not posing my question or my
6 concern clearly.

7 DR. FINDER: I think at the time that
8 these regulations were written, there were two
9 major issues that were being discussed, and I think
10 that the idea of the implants here was more the
11 cosmetic implant use.

12 I think you bring up a very good point
13 about additional training for patients that have
14 undergone surgery and the correct procedures on how
15 to do those examinations.

16 That is an issue we may be able to deal
17 with in some way through guidance as a
18 recommendation, but in terms of the regulations,
19 what we have got is what we have got, but I do
20 think that we do have the potential to expound a
21 little bit on the guidance and deal with some of
22 these other issues that you do bring up, and I
23 think that is a possibility, and if not directly in
24 our guidance, then, referring people to other
25 sources where they can get the correct procedures

1 o do some of these patients, so I think that is a
2 ery good idea.

3 MR. GOODE: I am Claude Goode from the
4 tate of California.

5 I would like to relate a horror story. I
6 ad a patient that complained of a breast
7 ompression that had an implant, and the
8 ompression ruptured the implant. The horror story
9 as that the woman literally was screaming at the
10 echnologist to release her, and she would not
11 elease the compression. Evidently there was
12 ignificant compression applied although the
13 achine was within normal operating standards.

14 This patient has been left to suffer.
15 There were no basic regulations that we could
16 nforce, and I do believe that there is a need for
17 DA or someone to at least come forth with some
18 standards for the mammography of patients with
19 breast implants. This needs to be discussed at
20 length and in depth, and presented somehow for the
21 echnologists, the physician who is not present,
22 and this does present a major problem.

23 I would just like to relate that horror
24 story.

25 MS. HARVEY: Dr. Harrison.

1 DR. HARRISON: Miles Harrison. May I ask
2 a similar question? This is someone with breast
3 implants for cosmetic purposes, this is not
4 post-mastectomy?

5 MR. GOODE: That is correct.

6 MS. HARVEY: Dr. Ikeda.

7 DR. IKEDA: This is Debra Ikeda from
8 Stanford University.

9 I think we are talking about two different
10 things here, I would like to clarify. That is a
11 terrible story. I would like to first address the
12 question about the post-mastectomy patient and the
13 tram flap patient.

14 I think if I understand correctly, you
15 were discussing patients who have undergone a
16 mastectomy and have either had latissimus dorsi
17 flap reconstruction or a tram flap reconstruction
18 with autologous material from the abdomen, or
19 patients, for example, who have a mastectomy and
20 may have a small amount of residual tissue, and the
21 recommendations for imaging that.

22 There are various amounts of scientific
23 literature for stating that either it is not
24 recommended, for example, there are some articles
25 in which patients have had mastectomies and then

1 looked, they have tried to find out if there is a
2 breast cancer recurrence, because that is what
3 everybody is concerned about, and there is varying
4 data on that.

5 Many places state that you should not be
6 doing those patients routinely for screening, but
7 if there is a lump, then, special views are often
8 used for that, and every patient is so different,
9 that I am not really sure that if there is a
10 problem, that you can actually say that this is the
11 right view to do or that is. Oftentimes, we have
12 to come up with special views to address that
13 specific patient's problem.

14 So, it is important when the technologist
15 is initially trained in their 40 hours, for
16 example, I think guidance is--correct me if I am
17 incorrect--but guidance states that the
18 technologist must learn all of the views, as well
19 as patients who have implants, and they must be
20 doing patients with implants under supervision, so
21 they do them correctly.

22 I think we are talking about two different
23 things.

24 DR. FINDER: Let me just correct one
25 thing. The regulations do require that

1 technologists who qualify under the final
2 regulations have some training in doing patients
3 with breast implants. There is no requirement that
4 they have to do a number of breast implant
5 patients, though.

6 DR. HARRISON: I am in complete agreement
7 with Dr. Ikeda, we are clearly talking about two
8 separate populations here, and actually, the
9 standard of care in our setting is such that we
10 don't do imaging routinely of people post
11 mastectomy with either autologous reconstruction or
12 implants.

13 I guess I need to ask the question, are
14 you referring to the training of the technologist
15 to, in fact, be able to do mammography on women who
16 have had cosmetic breast implants, are you
17 addressing the training?

18 MS. GILBERT: I guess I am addressing
19 both. I know in the native population, like Alaska
20 Native specifically, that follow-up isn't
21 recommended after that.

22 DR. HARRISON: Follow-up is not
23 recommended?

24 MS. GILBERT: It's just unknown, it is
25 just one of those unknown procedures, and because