

COPY

1

FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
OPHTHALMIC DEVICES PANEL

**This transcript has not
been edited and FDA
makes no representation
regarding its accuracy**

Thursday,
August 1, 2002

Salons A-C
Hilton Hotel Gaithersburg
620 Perry Parkway
Gaithersburg, Maryland

FRIEDMAN & ASSOCIATES, COURT REPORTERS
11923 Parklawn Drive, Suite 203
Rockville, MD 20852
(301) 881-8132

IN ATTENDANCE:

Jayne S. Weiss, M.D., Chair

Arthur Bradley, Ph.D., Voting Member

Michael R. Grimmett, M.D., Voting Member

Alice Y. Matoba, M.D., Voting Member

Karen Bandeen-Roche, Ph.D., Consultant, deputized to vote

Stephen A. Burns, Ph.D., Consultant

Mark A. Bullimore, MCOptom, Ph.D.,
Consultant, deputized to vote

Andrew J. Huang, M.D., Consultant, deputized to vote

Leo J. Maguire, M.D., Consultant, deputized to vote

Cynthia Owsley, Ph.D., Consultant, deputized to vote

William H. Swanson, Ph.D., Consultant, deputized to vote

Glenda V. Such, M.Ed., Consumer Representative

Ronald E. McCarley, Industry Representative

C O N T E N T S

PAGE

Call to Order

Jayne S. Weiss, M.D.
Chair

7

Introductory Remarks, Conflict of Interest
Statement, and Appointment to Temporary
Voting Status

Sara M. Thornton
Executive Secretary

7

Open Public Hearing

Ronald J. Link
Surgical Eyes

13

Letter from LASIK Recipient
(Read by Sara M. Thornton)

21

David L. Shell
Arlington, VA

24

Open Committee Session

Jayne S. Weiss, M.D.

27

Division Update

David M. Whipple
Deputy Director
Division of Ophthalmic and Ear, Nose, and
Throat Devices

27

C O N T E N T S

PAGE

Branch Updates

Vitreoretinal and Extraocular Devices Branch

James F. Saviola, O.D.
Chief

28

Intraocular and Corneal Implants Branch

Donna R. Lochner
Chief

33

Diagnostic and Surgical Devices Branch

Everette T. Beers, Ph.D.
Chief

34

PMA P970043/S010, CustomCornea Myopic
LASIK With the LADARVision 4000 System

Sponsor Presentation by Alcon

Kathleen Chester
Director, Regulatory Affairs

37

George Pettit, M.D., Ph.D.
Vice President and Chief Scientist,
Clinical Outcomes Research

38, 61, 72

Daniel Durrie, M.D.
Investigator
Hunkeler Eye Clinic, Inc.
Kansas City, MO

50

Omar Hakim, M.D.
Investigator
TLC Canada
Waterloo, Canada

55

Stephen Brint, M.D.
Investigator
Brint Vision Correction Center
Metairie, LA

65

Questions for the Sponsor

79

C O N T E N T S

	PAGE
FDA Presentation	
Everette T. Beers, Ph.D.	127
Jan C. Callaway Team Leader Diagnostic and Surgical Devices Branch	127
Bruce A. Drum, Ph.D. Vision Scientist and Physicist Diagnostic and Surgical Devices Branch	129
Questions for FDA	134
Additional Comments from the Sponsor	135
Committee Deliberations	
Primary Panel Reviewers	
Andrew J. Huang, M.D.	137
Arthur Bradley, Ph.D.	142
Panel Discussion of PMA P870043/S010 and FDA Questions	156
Open Public Hearing	
Ronald J. Link	245
Closing Comments by FDA	247
Closing Comments by the Sponsor	248

C O N T E N T S

	PAGE
Voting Options	250
Panel Discussion and Recommendation Taken by Vote	252
Polling of Panel Votes	289
Comments from Consumer and Industry Representatives	
Glenda V. Such, M.Ed.	292
Ronald E. McCarley	293
Closing Comments	
Jayne S. Weiss, M.D.	293

1 P R O C E E D I N G S (8:35 a.m.)

2 DR. WEISS: I'd like to call this meeting of
3 the Ophthalmic Devices Panel to order and we'll have
4 introductory remarks from Sarah Thornton.

5 MS. THORNTON: Good morning and welcome to the
6 104th meeting of the Ophthalmic Devices Panel.

7 Before we proceed with today's agenda, I have a
8 few short announcements as usual that I'd like to make.
9 I'd like to remind everyone out there as well as the panel
10 and the FDA folks to sign in on the attendance sheet in the
11 registration area just outside the meeting room here.
12 Messages for panel members and FDA participants and
13 information or special needs should be directed through Ms.
14 Annmarie Williams or Ms. Jennifer Weber who are available
15 in the registration area.

16 The phone number for calls to the meeting area
17 is (301) 977-8900. In consideration of the panel, the
18 sponsor and the agency, we ask that those of you with cell
19 phones and pagers either turn them off or put them on
20 vibration mode while in this room.

21 Lastly, will all meeting participants please
22 speak clearly into the microphone, give your name clearly,
23 until I get a signal from the transcriber that he no longer
24 needs your name, so that we will have an accurate recording
25 of your comments, please.

1 Now, at this time, I'd like to announce the
2 confirmation of the new Ophthalmic Devices Panel Chair, Dr.
3 Jayne Weiss. We also have three newly appointed voting
4 members, Drs. Anne Coleman, Allen Ho, and Timothy McMahon,
5 who are regrettably unable to be with us today. However,
6 we look forward to their attendance at future meetings.

7 I'd also like to extend a special welcome and
8 introduce to the public and panel and FDA staff three panel
9 consultants who are with us for the first time today. Dr.
10 Stephen Burns. Dr. Burns comes to us from Boston,
11 Massachusetts, where he is a senior scientist at the
12 Schepens Eye Research Institute and associate professor at
13 the Harvard University. Dr. Cynthia Owsley is from
14 Birmingham, Alabama, where she is the Professor of
15 Ophthalmology at the School of Medicine and Co-Director of
16 the Center for Research on Applied Gerontology at the
17 University of Alabama. And Dr. William Swanson is a senior
18 research scientist in the Department of Clinical Sciences
19 at the State University of New York, College of Optometry,
20 in New York City.

21 Welcome to you all. Hope you enjoy your day
22 with us.

23 Will the remaining panel members take the time
24 now to introduce themselves, and I'd like to begin with our
25 industry rep.

1 MR. McCARLEY: My name's Rick McCarley. I'm
2 the industry rep. I'm the President and CEO of Ophtec in
3 Boca Raton, Florida.

4 DR. BANDEEN-ROCHE: I'm Karen Bandeen-Roche,
5 Associate Professor of Biostatistics at Johns Hopkins
6 University.

7 DR. BULLIMORE: Mark Bullimore, Associate
8 Professor, Ohio State University.

9 MS. SUCH: I'm Glenda Such, consumer
10 representative, Director of Computer Training Programs at
11 Lighthouse International, New York City.

12 DR. MATOBA: I'm Alice Matoba, Associate
13 Professor of Ophthalmology, Baylor College of Medicine.

14 DR. GRIMMETT: Michael Grimmett, Assistant
15 Professor of Ophthalmology at the University of Miami,
16 School of Medicine.

17 DR. WEISS: Jayne Weiss, Professor of
18 Ophthalmology and Pathology, Kresge Eye Institute, Wayne
19 State University, Detroit.

20 DR. BRADLEY: Arthur Bradley, Professor of
21 Visual Science, Indiana University.

22 DR. HUANG: Andrew Huang, Associate Professor
23 of Ophthalmology, University of Minnesota.

24 DR. MAGUIRE: Leo Maguire, Associate Professor
25 of Ophthalmology, Mayo Clinic.

1 MR. WHIPPLE: And I'm Dave Whipple, Deputy
2 Director of the Division of Ophthalmic, Ear, Nose and
3 Throat Devices.

4 MS. THORNTON: Thank you very much.

5 At this time, I'd like to read the conflict of
6 interest statement for the meeting today.

7 "The following announcement addresses conflict
8 of interest issues associated with this meeting and is made
9 part of the record to preclude even the appearance of an
10 ~~impropriety.~~

11 "To determine if any conflict existed, the
12 agency reviewed the submitted agenda for this meeting and
13 all financial interests reported by the committee
14 participants. The conflict of interest statutes prohibit
15 special government employees from participating in matters
16 that could affect their or their employer's financial
17 interests. However, the agency has determined that
18 participation of certain members and consultants, the need
19 for whose services outweigh the potential conflict of
20 interest involved, is in the best interests of the
21 government.

22 "Therefore, waivers have been granted to Drs.
23 Mark Bullimore and Stephen Burns for their interest in
24 firms that could potentially be affected by the panel's
25 recommendation. Dr. Bullimore's waiver allowing him to

1 participate fully in today's deliberations involves a
2 consulting arrangement with a competing technology firm.
3 For this unrelated consulting service, he receives less
4 than \$10,000 a year. Dr. Burns' limited waiver allows him
5 to participate in the panel discussion but excludes him
6 from voting. His interest involves a grant to his employer
7 with a competing firm funded for less than \$100,000 per
8 year for which he has involvement in data collection and
9 interpretation.

10 "Copies of these waivers may be obtained from
11 the agency's Freedom of Information Office, Room 12A-15 of
12 the Parklawn Building.

13 "We would like to note for the record that the
14 agency took into consideration other matters regarding Drs.
15 Arthur Bradley, Michael Grimmett, and Jayne Weiss. They
16 reported interest in firms at issue but in matters not
17 related to today's agenda. The agency has determined
18 therefore that they may participate fully in all
19 discussions.

20 "In the event that the discussions involve any
21 other products or firms not already on the agenda for which
22 an FDA participant has a financial interest, the
23 participant should excuse him or herself from such
24 involvement and the exclusion will be noted for the record.

25 "With respect to all other participants, we ask

1 in the interest of fairness that all persons making
2 statements or presentations disclose any current or
3 previous financial involvement with any firm whose products
4 they may wish to comment upon."

5 I will read now the appointment to temporary
6 voting status.

7 "Pursuant to the authority granted under the
8 Medical Devices Advisory Committee Charter, dated October
9 27th, 1990, and as amended August 18th, 1999, I appoint the
10 following individuals as voting members of the Ophthalmic
11 Devices Panel for this meeting on August 1st, 2002: Dr.
12 Karen Bandeen-Roche, Dr. Mark Bullimore, Dr. Andrew Huang,
13 Dr. Leo Maguire, Dr. Cynthia Owsley, Dr. William Swanson.
14 For the record, these individuals are special government
15 employees and consultants to this panel or other panels
16 under the Medical Devices Advisory Committee. They have
17 undergone the customary conflict of interest review and
18 have reviewed the material to be considered at this
19 meeting." Signed Dr. David W. Feigal, Jr., Director of the
20 Center for Devices and Radiological Health, dated July
21 19th, 2002.

22 Thank you, Dr. Weiss.

23 DR. WEISS: Thank you, Sally.

24 We will now start the open public hearing.

25 There are three individuals who have requested to speak

1 before us. I would appreciate when they approach the
2 podium, they should identify themselves and any financial
3 conflicts or potential conflicts, and we'll start with Mr.
4 Ron Link. If you could come to the podium and read your
5 statement, please?

6 MR. LINK: Good morning.

7 My name's Ron Link. I'm Executive Director of
8 Surgical Eyes, Tampa, Florida. I have no conflict of
9 interests with regard to this meeting.

10 Well, good morning, Ophthalmic Devices Panel
11 members and members of the audience. I'm here today to
12 advocate on behalf of thousands of people who have
13 longstanding complications of refractive surgery.
14 Together, we can confront the challenge of rehabilitation
15 of as many of these people as possible. Surgical Eyes
16 supports the advancement of wavefront technology in that it
17 may hold promise for the visual rehabilitation of these
18 patients who live with complications of LASIK and other
19 refractive surgeries. No less important is our mutual
20 obligation that fewer patients with complications be
21 created.

22 I want to reference this slide of a recent
23 survey that we did at Surgical Eyes, only ran for a week,
24 but I think the results are compelling. This first slide
25 were LASIK successes in the sense that they had 20/40 or

1 better UCVA. You can see significant quality of vision
2 issues and you expect much more in the excellent
3 categories, but we see that's not the case across different
4 lighting conditions. Now, attempts to improve that vision
5 with contact lenses or glasses yielded slight improvement.
6 So what this tells me is that we need better technologies
7 to rehabilitate these folks. Now, people who were not
8 corrected to 20/40 or better had even less good result with
9 the attempts at correction.

10 These people and many others at Surgical Eyes
11 may very well benefit from laser and contact lens wavefront
12 technology. If the panel votes for approval on the PMA for
13 this device, we ask that it do so with the following
14 conditions.

15 Number 1. Controlled studies on post-
16 refractive eyes. Clinical studies at multiple sites across
17 the United States on post-refractive eyes with a minimum
18 one-month follow-up. Anecdotal reports from the estimated
19 200 global cases of wavefront-guided treatments on patients
20 with complications reveal that results are often immediate.
21 Clinical studies should not only include those with under-
22 corrections or smaller aberrations. Decentration, central
23 islands, disparity between the affected optical zone and
24 pupil size, these and other complications result in higher-
25 order aberrations. We believe wavefront treatment on such

1 patients should first be performed in controlled
2 circumstances by surgeons who have the requisite skill and
3 technical experience to perform therapeutic studies on
4 post-refractive eyes.

5 Controlled studies are necessary to prevent a
6 rush by patients and doctors alike to avail themselves of a
7 new device in off-label use without the necessary specific
8 clinical protocols and data analysis to evaluate the safety
9 and effectiveness of a device on post-refractive eyes.
10 Should you vote for approval, such a condition attached to
11 today's PMA would be a win-win for all concerned.

12 Second. Professional use information and
13 patient information booklets. Preexisting dry eye should
14 be listed as a contraindication warning in the professional
15 use information and patient information booklets of any
16 laser approved for LASIK brought before the FDA. The FDA's
17 LASIK website has a section entitled "When is LASIK Not for
18 Me?" Under the Other Risks section, it states, "Dry eye.
19 LASIK surgery tends to aggravate this condition." Under
20 another section entitled "What Are the Risks and How Can I
21 Find the Right Doctor For Me?," it states, "Some people may
22 develop severe dry eye syndrome. As a result of surgery,
23 your eye may not be able to produce enough tears to keep
24 the eye moist and comfortable. This condition may be
25 permanent." This last slide shows that dry eye is a

1 significant factor of the 100 percent of LASIK patients
2 that are reported in this survey.

3 We believe that the information published by
4 the FDA with regard to lasers should be consistent with the
5 information presented in laser device professional use
6 information and patient information booklets. For the last
7 approved LASIK device, dry eye syndrome was listed under
8 exclusion criteria for the PMA study data. If dry eyes are
9 excluded from PMA data, the same warning should be extended
10 to the public.

11 The American Academy of Cataract and Refractive
12 Surgery recently issued LASIK Screening Guidelines for
13 Patients on June 4th, 2002. Under the less than ideal
14 LASIK candidates listed have a history of dry eyes as they
15 may find that the condition worsens following surgery. We
16 recommend that these guidelines be consulted to update the
17 information presented in both professional use information
18 and patient information booklets.

19 A similar argument can be made with regard to
20 the issue of pupil size. In any professional use
21 information and patient information booklets that we have
22 read for this same laser device, the language is not
23 consistent. By way of example, with regard to the last
24 FDA-approved laser, it states in the professional use
25 information that visual performance could possibly be

1 worsened by large pupil sizes or decentered pupils. This
2 sentence is not included in the Precautions Section of the
3 patient information booklet. The FDA website under "When
4 Is LASIK Not For Me?," states under Other Risk Factors,
5 "Your doctor should screen you for the following conditions
6 or indicators of risk: large pupils. Make sure this
7 evaluation is done in a dark room." Jumping ahead a bit,
8 "This can cause symptoms, such as glare, halos, starbursts,
9 and ghost images, double vision after surgery. In some
10 patients, these symptoms may be debilitating. For example,
11 a patient may no longer be able to drive a car at night or
12 in certain weather conditions, such as fog."

13 The overriding point is that the information
14 put forth by the FDA on its website and what it requires of
15 manufacturers in professional use information and patient
16 information booklets should be consistent with regard to
17 dry eye, pupil size, and any other pre- or postoperative
18 information provided to the public. Degree of success for
19 higher myopes being yet just another example.

20 Also, we suggest that pictures be provided in
21 patient information booklets so that patients understand
22 the visual manifestations of both lower and higher-order
23 aberrations in various lighting conditions, particularly at
24 night.

25 And lastly, post-approval studies. If required

1 for contact lenses, by way of example PMA for the CIBA
2 Extended Wear Contact Lens, it is logical to expect that
3 the same for LASIK, PRK, LASIK, or any other laser device
4 used to perform these surgeries. Lastly, questionnaires
5 should be included to account for any potential adverse
6 effects on quality of life.

7 Surgical Eyes is cautiously optimistic that
8 post-refractive patients with complications may indeed
9 benefit from controlled studies of the wavefront device
10 being presented here today. We support such advancement.
11 No less important is the identification and ready
12 disclosure of all pre- and postoperative risk factors to
13 patients and doctors alike.

14 Members of the panel and the audience, thank
15 you for your time.

16 DR. WEISS: Thank you very much, Mr. Link, for
17 your thoughtful comments.

18 I'd ask the panel at this point if anyone has
19 any questions for Mr. Link. Dr. Bullimore?

20 DR. BULLIMORE: This is Mark Bullimore. I have
21 a question for the FDA.

22 The PMA does not cover therapeutic use of the
23 device. It just covers primary LASIK, is that correct?
24 I'd like to add my thanks to the chair's. I do appreciate
25 your efforts at advocacy.

1 MR. LINK: Thank you.

2 DR. WEISS: Dr. Bradley?

3 DR. BRADLEY: Question for the speaker.

4 You are suggesting that the FDA mandate a one-
5 month follow-up post-surgical clinical study.

6 MR. LINK: On post-refractive eyes, yes.

7 DR. BRADLEY: Are you aware that those that
8 we're studying today has one-month, three-month, and six-
9 month, and that's typically what we see? I'm wondering.
10 Are you suggesting something different from that?

11 MR. LINK: Well, I think from the anecdotal
12 reports that we've heard of global wavefront treatment, the
13 results, the efficaciousness of the procedure, is noticed
14 immediately. So at a one-month period, from the doctors
15 that I've spoken to who've done the procedures, say that
16 you do see whether or not it's going to work. For
17 instance, on a decentered ablation or spherical
18 aberrations. What we want to avoid is people rushing to
19 have it done and then there not being the necessary
20 controls to see if it actually works, and if the FDA were
21 to vote for approval with conditions, we would have data
22 that could be shared openly and I think it would be a good
23 thing for all of us.

24 DR. WEISS: If there are no other questions,
25 thank you, Mr. Link.

1 MR. LINK: Thank you.

2 MR. WHIPPLE: Jayne?

3 DR. WEISS: Yes, Dr. Whipple?

4 MR. WHIPPLE: Very quickly, for Mr. Link.

5 MR. LINK: Yes.

6 MR. WHIPPLE: You've recommended some specific
7 changes to the website for updates on our FDA website.

8 MR. LINK: Well, I think by virtue of the fact
9 that the patient information booklets and the professional
10 use information, it's my understanding that what is present
11 there is mandated by the FDA, and the patient information
12 booklets are to end up in the hands of the patients. The
13 information that is on the FDA website is actually more
14 stringent than what goes in the patient information
15 booklets, and we believe it should be consistent throughout
16 all three forums.

17 MR. WHIPPLE: Periodically, obviously, we
18 update our website and we'll take into consideration what
19 you've said and see if there's any merit to making any
20 changes.

21 MR. LINK: Yes. If the information that's on
22 the site ends up in the booklets and the professional use
23 information, that would be superb.

24 MR. WHIPPLE: Thank you.

25 MR. LINK: Thank you.

1 DR. WEISS: Thank you very much.

2 We have a letter from an additional person who
3 was not able to appear here today that Sarah Thornton will
4 read to us.

5 MS. THORNTON: "Dear FDA Panel Members:
6 Unfortunately I was not able to attend today's meeting in
7 person. However, I would like to request that this letter
8 be read aloud on my behalf.

9 "Two years ago, I had LASIK surgery. My
10 initial uncorrected visual acuity was 20/50, 20/20. I was
11 therefore considered a LASIK success. The reality is that
12 as a result of having large pupils, I have debilitating
13 loss of night vision, ghosting, haloes, starbursts and loss
14 of contrast sensitivity. My eyes are extremely dry and
15 burn constantly. Over time, I experienced a loss of
16 approximately 35 percent of the surgical effect. Rigid gas
17 permeable contact lenses are my best hope for visual
18 rehabilitation. Unfortunately, I am intolerant to hard
19 lenses, probably due to my dry eyes.

20 "I wish I could describe to you the totality of
21 what LASIK surgery has done to my life. Driving at night
22 is extremely difficult and dangerous for me. I can no
23 longer enjoy things that I used to take for granted, such
24 as going to the movies or dining in a dimly-lit restaurant.
25 Now, I see two or three smeared moons in the sky at night.

1 "Today, you will be considering an application
2 for wavefront customized LASIK. I completely support the
3 advancement of this technology that has the potential to
4 treat post-surgical eyes for the correction of induced
5 higher-order aberrations such as those I suffer from.

6 "More importantly, however, I am asking the
7 panel to condition this PMA approval to include preexisting
8 dry eyes in the contraindications and to limit the approval
9 of this device based on pupil size. The scoptic pupil
10 measurement should not exceed the effective optical zone.
11 The transition zone is not receiving the full refractive
12 treatment and therefore should not be included as part of
13 the optical zone for this purpose.

14 "Wavefront custom ablation is still LASIK
15 surgery. In the words of Dr. I. Howard Fine, past
16 President of the American Society of Cataract and
17 Refractive Surgery, 'As we all know, LASIK transects the
18 corneal nerves, therefore inducing dry eyes in most
19 patients.' The size of the optical zone in wavefront
20 treatments is limited, just as in traditional LASIK, by
21 corneal thickness. Limiting the optical zone results in
22 the induction of spherical aberrations in patients whose
23 pupils dilate larger than the optical zone in low light.
24 On Table 16 in the Summary of Safety and Effectiveness Data
25 for the approval of the LADARVision laser, 29.5 percent of

1 the patients in the study reported halos that were worse or
2 significantly worse postoperatively. I believe these
3 patients reporting worse halos were primarily the patients
4 with large pupils.

5 "I know another woman with 8mm pupils who was
6 treated on the LADARVision system by one of the
7 investigators for this PMA. She, too, sees massive
8 starbursts and halos at night, has dry eyes, regressed, and
9 is battling recurrent corneal erosions. She struggles with
10 uncomfortable hard contact lenses at times just to see her
11 infant son's face clearly. The surgeon was fully aware of
12 her large pupils before he treated her. Now, she is
13 practically disabled at night.

14 "Patients are making their decision to have
15 this surgery based on a barrage of advertising that doesn't
16 disclose the risks or contraindications. My informed
17 consent did not mention dry eyes or large pupils. Dr.
18 Fine's statement shows that the industry is fully aware of
19 the magnitude of the dry eye problem. Why is it not taking
20 preventive action? The fact that 29.5 percent of patients
21 were complaining of worse halos proves the industry does
22 not even consider night vision problems as a complication,
23 even though it can be incapacitating. This is a medically
24 unnecessary, elective procedure and therefore should be
25 held to higher standards.

1 "It is not enough to simply put warnings in the
2 patient labeling. Surgeons do not always give patients a
3 copy of the patient information booklet. I was not given
4 one. To protect patients, the FDA must limit the device
5 based on pupil size and dry eyes must be listed in the
6 contraindications.

7 "Ladies and gentlemen on the panel, I believe
8 that some good can come from my terrible experience.
9 Please help me make a difference for future refractive
10 surgery patients who don't know that dry eyes, glare and
11 halos are not simply minor side effects. They are life-
12 altering complications.

13 "Thank you very much."

14 DR. WEISS: Thank you, Sally.

15 We will also have Mr. David Shell approach the
16 podium and he would also like to read a statement.

17 MR. SHELL: Members of the committee, I am
18 David Shell, mechanical engineer from Arlington, Virginia.
19 I appreciate the opportunity to speak here. The focus of
20 my testimony will be on the inclusion of an incident of
21 LASIK-induced dry eye statistic and to the patient
22 information booklet for the medical device before us today.

23 Four years ago, I underwent LASIK. Since my
24 surgery, I live in daily misery from burning and stinging
25 eyes induced by LASIK. Artificial tears don't bring much

1 relief. Eyedrops give only temporary relief and cause
2 greater visual distortion when used. My LASIK dry eye is
3 not a minor problem, as downplayed by some
4 ophthalmologists. It's a disability. I estimate that I am
5 blind approximately 10 percent of the time due to my eyes
6 being closed because of the pain. At the time of my
7 surgery, I was told only a small number of patients
8 experience a complication from this procedure.

9 There is substantial evidence that shows this
10 crippling side effect to be relatively common. For
11 example, an article in EyeWorld stated that 100 percent of
12 patients have dry eye after LASIK. While most patients
13 improve, many do not. Numerous articles in industry
14 magazines and journals talk about how to manage LASIK dry
15 eye. Internet websites, such as www.surgicaleyes.org,
16 discuss this issue frequently.

17 I know now that I did not have the information
18 that would have assisted me in making a fully informed
19 decision. No one really knows the risk of getting this
20 debilitating condition in terms of percentage or the
21 information is just not getting out. Therefore, a person
22 is unable to make an informed decision about having this
23 procedure. Should not this type of data be available to
24 the public? This type of data is no where to be found in
25 the patient information booklet.

1 My recommendations are as follows. Premarket
2 approval for this medical device should be contingent upon
3 manufacturer conducting clinical studies on the incident of
4 LASIK-induced dry eye; data to be listed in the
5 manufacturer's patient information booklet in terms of
6 percentage, not just a casual mention that one could get
7 dry eyes from this procedure. We need a percentage.

8 Adoption of these recommendations will help
9 increase public awareness about this serious overlooked
10 complication. I believe these recommendations are fair and
11 reasonable, easy to administer, and do not impose an undue
12 burden on the industry.

13 Before I conclude, I want to remind everyone
14 that our eyes are very precious. The standards for safety
15 and effectiveness need to be very high for an elective
16 procedure on one's eyes. Personally, I don't think they're
17 high enough. I didn't need this surgery and ended up with
18 inheriting a lifetime of misery and pain. I'm asking the
19 committee to make certain that any device that purports to
20 correct a relatively minor problem does not create
21 crippling visual defects as a result.

22 Members of the committee, this concludes my
23 testimony. Thank you very much.

24 Also, I'd like to add, if anyone wants to try
25 some Dry Eyes, during the break, I have a number of unused

1 vials up here, and I'd be happy to give you one.

2 DR. WEISS: Thank you.

3 Does anyone have any questions for Mr. Shell?

4 (No response.)

5 DR. WEISS: Seeing no questions, thank you very
6 much, Mr. Shell.

7 Are there any other participants who would like
8 to come forward during the open public hearing?

9 (No response.)

10 DR. WEISS: If not, that will end the open
11 public hearing, and we will move on to the open committee
12 session and start with the FDA division update.

13 Dr. Whipple?

14 MR. WHIPPLE: Thank you. It's Mr. Whipple, by
15 the way. Thank you for the promotion.

16 It's been some time since I've had a chance to
17 address this panel. I'm usually doing my thing behind the
18 scenes and occasionally going up to Ralph after the meeting
19 or during the meeting and whispering things in his ear,
20 that kind of thing. But Ralph can't be with us today and
21 he's asked me to set in for him and I'm glad to do that.
22 He does send his regards to everybody here and he wants to
23 make sure you know that he will be at the next panel
24 meeting.

25 Now, as for division updates, the branch chiefs

1 really have all the important information to provide to
2 you, so I won't steal their thunder or step on their toes.
3 But I do have one personnel piece of information that I do
4 want to present. The office director for Device
5 Evaluation, Dr. Bernie Statland, who is Ralph's boss, will
6 be leaving the center at the end of August, and on behalf
7 of the division, I want to thank Dr. Statland for his
8 guidance and support in the past two years that he's been
9 our office director. We sure will miss his kind and gentle
10 demeanor and we do wish him well as he goes on to pursue a
11 law degree at the University of Minnesota, and as soon as
12 we know who Ralph's new boss is going to be, we'll let you
13 know.

14 So that's all I have for now, and we can go on
15 to the branch updates.

16 DR. WEISS: Thank you.

17 I think Dr. Saviola will be starting.

18 DR. SAVIOLA: Good morning, panel members.

19 I'd like to update you today on one 510(k)
20 clearance and two PMA approvals. We recently cleared an
21 application on May 9th, 20002, for the ChromaGen Reading
22 Aid Soft Contact Lens manufactured by Cantor and Silver
23 Limited of England. The indication includes the correction
24 of refractive ametropia as it would for a standard contact
25 lens. It also has a statement. "In addition, the lenses

1 may also be prescribed as a colored filter for individuals
2 who experience reading discomfort not related to binocular
3 vision problems or uncorrected refractive error."

4 This lens had previously obtained a marketing
5 clearance in October of 2000 as an optical aid for people
6 with red-green color deficiencies. There was a small
7 clinical study conducted that supported the use of the lens
8 as a colored filter to aid individuals who experience
9 reading discomfort. The 510(k) Summary of Safety and
10 Effectiveness available on our website provides a
11 description of the clinical study, and for those
12 interested, the K number is K012132.

13 Current literature studies report inconsistent
14 results concerning the effect of colored filters on reading
15 rate and comprehension and symptoms of reading discomfort,
16 along with rate of reading and reading comprehension, can
17 be strongly influenced by psychological factors. The
18 clinical data in our view did not support the use of the
19 ChromaGen lens in treating dyslexia or improving the
20 general reading speed. Dyslexia is a poorly defined
21 anomaly with some controversy as to how it is identified.
22 Therefore, the lens is not cleared with the indication to
23 treat dyslexia or improve reading speed and the precaution
24 statement was added to the labeling to inform patients that
25 results are variable due to the subjective nature of visual

1 discomfort and that not all patients will experience
2 success.

3 In our view, there is minimal risk associated
4 with this device as there's no indication to aid reading
5 discomfort. The risk is comparable to other tinted soft
6 contact lens and while effectiveness is expected to be
7 variable, the lens may be beneficial for some people. We
8 had issued a homework assignment to a panel member on this
9 project, and I want to thank the panel for assistance in
10 review of this submission.

11 The first PMA approval I want to inform the
12 panel about is the Paragon CRT and Quadra RG Lenses for
13 overnight orthokeratology which we discussed at the January
14 2002 panel meeting. These were approved on June 13th.
15 Following the panel meeting, ODE issued an approvable
16 letter which the firm responded to. Following review of
17 the additional clinical data provided and interactive
18 review of the draft product labeling, final review was
19 granted. The lenses are manufactured by Paragon Vision
20 Sciences of Mesa, Arizona, and the Paragon CRT designs are
21 indicated for overnight wear in a corneal refractive
22 therapy fitting program for the temporary reduction of
23 myopia up to six diopters in eyes with astigmatism of up to
24 1.75 diopters. The Paragon RG designs have essentially the
25 same indications as the CRT lenses, except the pretreatment

1 of myopia is up to three diopters in eyes with astigmatism
2 up to 1.5.

3 In order to address effectiveness concerns of
4 the Quadra RG design used overnight, a further analysis of
5 existing data was conducted by comparing the Quadra daily
6 wear effectiveness data with the CRT overnight data for a
7 three-month time interval which was the duration of the
8 data we were studying. There are no statistically
9 significant differences in reduction of pretreatment
10 refractive error, accuracy stability or uncorrected visual
11 acuity.

12 As to the age issue that was discussed during
13 the panel meeting, among the many recommendations by the
14 advisory panel was the limitation to those 18 years of age
15 or older since limited data on this age group were
16 presented during the panel meeting. There was not an age
17 restriction included in the final approval by FDA since the
18 company provided additional data on adolescents between
19 ages 12 to 17 who completed the study, and they accounted
20 for 11 percent of the completed dataset.

21 The primary effectiveness concern for this age
22 group is their expected progression of myopia. Although
23 daily wear contact lens wear can reshape the cornea or
24 known as orthokeratology has been practiced since the
25 1960s, the long-term safety effects of overnight wear for

1 reshaping the cornea are not known for any age population
2 and there's not been shown any real long-term safety issues
3 for daily wear orthokeratology.

4 And the last PMA update is for the Menicon Z,
5 which is an RGP lens, that was just approved in July, on
6 July 12th, 2002, and this is a supplement for P990018 for
7 extended wear of an RGP up to 30 days of wear. This had a
8 prior daily wear approval in a variety of designs for
9 indications for myopia, hyperopia and presbyopia in both
10 aphakic and non-aphakic persons. In the extended wear
11 version, it's approved in the spherical, aspheric and non-
12 prism ballast toric and non-prism ballast multifocal lenses
13 for again myopia, hyperopia and presbyopia, but it's only
14 for non-aphakic persons, and there's limitation on the
15 power range of up to +8 diopters for hyperopes.

16 Although this was not discussed at a panel
17 meeting, we did solicit homework assignments from three
18 panel members in order to corroborate the internal FDA
19 clinical review. The panel reviews did not raise any
20 additional clinical issues that were unique to this device
21 or different from those identified in the internal review.
22 All three homework assignments recommended approval of the
23 device for extended wear up to 30 days. The post-approval
24 condition of conducting a clinical study was placed upon
25 this approval in the same manner as the previous two

1 silicone hydrogel lenses from CIBA and from Bausch & Lomb.

2 That concludes my remarks. Does anybody have
3 any questions?

4 (No response.)

5 DR. WEISS: Thank you, Dr. Saviola.

6 Dr. Lochner will update us.

7 MS. LOCHNER: Again, thank you for that
8 honorary medical degree.

9 DR. WEISS: I'm giving out medical degrees left
10 and right here today. Step right up, Ms. Lochner.

11 MS. LOCHNER: All right. At the January 2002
12 meeting, the panel reviewed P010059, the Morcher capsular
13 tension ring, and recommended that the PMA was approvable
14 pending additional analyses of the clinical data. I would
15 like to advise you that this document is still being
16 reviewed by FDA. The issues that the panel discussed were
17 related in a letter to the sponsor and we are currently
18 awaiting their responses.

19 Secondly, I'd like to advise you that on March
20 26th, 2002, we cleared a new glaucoma shunt from Optonol,
21 Limited, K012852, the Ex-Press Miniature Glaucoma Implant,
22 Models R-30 and R-50, and this device is different from
23 previously cleared shunts in that it is a stainless steel
24 tube with a blunt needle-shaped penetrating tip at one end
25 and a flat angled flange at the distal end. It functions

1 similarly in shunting aqueous fluid from the anterior
2 chamber into a conjunctival bleb and is intended to reduce
3 intraocular pressure in patients with glaucoma where
4 medical and conventional surgical treatments have failed.

5 That concludes my comments.

6 DR. WEISS: Thank you.

7 Are there any questions?

8 (No response.)

9 DR. WEISS: Thank you.

10 We'll go on to Dr. Beers.

11 DR. BEERS: Thank you.

12 We've approved three devices since the last
13 panel. A couple of months ago, we approved the Bausch &
14 Lomb PMA P990027 Supplement 2 for the Technolas 217A for
15 high myopia and that's up to MRSC of less than -12 with a
16 sphere of less than 11. We also on April 11th, 2002, were
17 approved the Refractec PMA, P010018, for the ViewPoint CK
18 or Conductive Keratoplasty System. That was reviewed by
19 this panel which recommended that it was approvable on
20 November 30th of 2001.

21 Based on the panel's recommendations, the
22 indication for this device is for the temporary reduction
23 of spherical hyperopia in patients who have .75 diopters to
24 3.25 diopters of psychoplegic spherical hyperopia, and also
25 added to the indications for use is the statement that the

1 magnitude of correction with this treatment diminishes over
2 time in some patients, with some patients retaining some or
3 all of their intended refractive correction.

4 The other device that we approved was approved
5 on December 19th, 2002, for the VISX Humanitarian Device
6 Exemption, or HDE, H000002 for the Customized Contoured
7 Ablation Pattern Method for the treatment of certain
8 patients, and the indication is important here, for the
9 treatment of certain patients with symptomatic decentered
10 ablations from previous laser surgery as viewed on the
11 Zeiss Humphrey topography unit.

12 Now, I'm guessing that many of you don't know
13 what an HDE is. An HDE is an application that's similar to
14 a PMA but is exempt from the effectiveness requirements of
15 PMA. An approved HDE authorizes marketing of a
16 humanitarian use device. The humanitarian use device is
17 intended to benefit patients in the treatment and diagnosis
18 of diseases or conditions that affect fewer than 4,000
19 individuals per year in the U.S. So obviously given such a
20 small patient base, it's difficult for these types of
21 devices to gain significant clinical trials to support
22 safety and effectiveness or certainly to support the
23 effectiveness. So this is a little bit different route for
24 some of these.

25 It's important, though, to remember with these

1 devices that the use of the device at each institution is
2 overseen by the IRB of that institution and the IRB may
3 make decisions on whether to use the device on a case-by-
4 case basis. So there are certain severe -- well, I
5 wouldn't say severe but there are certain limitations on
6 the HDE that you don't see with a PMA.

7 That concludes my presentation. Are there any
8 questions?

9 DR. WEISS: Dr. Bullimore?

10 DR. BULLIMORE: Mark Bullimore. I have a
11 question for Dr. Beers. It is Dr. Beers.

12 You say the approval is limited to 4,000 cases
13 per year. Is that the actual approval or is that just --

14 DR. BEERS: That's in the Act. I mean, that's
15 a limit.

16 DR. BULLIMORE: The FDA has no role or
17 responsibility to monitor the number of procedures that are
18 done after approval?

19 DR. BEERS: They are monitored. The sponsor
20 has to keep track of that.

21 DR. BULLIMORE: Okay. Thank you.

22 DR. WEISS: If there are no other questions,
23 I'd like to thank Dr. Saviola, Ms. Lochner and Dr. Beers.

24 If there's no other information to be updated
25 from the agency at this point, we're going to move ahead to

1 the discussion and review of PMA Number P970043/S010.

2 I passed muster with Sally, so I must be doing
3 okay.

4 We'd like to inform the sponsor they have one
5 hour, and I would like each presenter when they come
6 forward to identify themselves at the beginning of the
7 presentation, also to inform us of any financial conflicts
8 or potential conflicts.

9 MS. CHESTER: Good morning. I'm Kathleen
10 Chester, Director of Regulatory Affairs for Alcon's
11 Refractive Products, and today, we'll be presenting the
12 clinical results from the CustomCornea Myopic LASIK
13 Clinical Trial involving the commercially available
14 LADARVision Laser System.

15 The agenda for our presentation includes the
16 following: I will give a brief introduction. Dr. George
17 Pettit from Alcon will provide an overview of wavefront
18 technology. Drs. Daniel Durrie and Omar Hakim, two of our
19 clinical investigators, will present a summary of our
20 safety and effectiveness results. Dr. Pettit will then
21 discuss our wavefront aberration clinical outcomes followed
22 by Dr. Stephen Brint who will discuss the clinical
23 implications of wavefront correction based on our clinical
24 results. And finally, if time permits, we will respond to
25 a number of questions we received in advance from the

1 panel.

2 This PMA supplement application requests an
3 expansion of the existing indications to include wavefront-
4 guided custom cornea LASIK correction of myopia of up to -7
5 diopters of sphere and less than -.5 diopters of cylinder
6 at the spectacle plane in the subjects who are 21 years of
7 age or older and with a documented stability of refraction.

8 Alcon is pursuing approval of spherical myopia
9 only at this time. There are no safety issues related to
10 this decision. The decision is based on the intent to
11 provide the most effective astigmatic outcomes possible
12 with this new technology with minor adjustments in the
13 algorithm before seeking approval for the myopic
14 astigmatism indication.

15 The study population consists of a safety
16 cohort of 426 eyes in the range of up to -7 diopters of
17 sphere and up to -4 diopters of astigmatism. The primary
18 effectiveness cohort is comprised of a subset of those
19 eyes, of which there are 139, in the range of up to -7
20 diopters of sphere and less than .5 diopters astigmatism.

21 Now, I'd like to introduce Dr. George Pettit
22 who will give you an overview of wavefront technology.

23 DR. PETTIT: Good morning. I'm the chief
24 scientist at the Alcon Orlando Technology Center and
25 therefore I do have a financial interest in this

1 technology.

2 I'd like to start with a very simple
3 introduction to what is wavefront-guided ablation. What
4 are we talking about here? When we talk about wavefront-
5 guided customized treatment, our definition includes,
6 first, a quantitative measurement of both the lower and
7 higher-order aberrations -- i.e., those aberrations beyond
8 simple sphere and cylinder -- that are present in the eye
9 and transfer of that detailed aberration data to an excimer
10 laser which positions the treatment profile correctly on
11 the eye and calculates and delivers a specific ablation
12 pattern unique to each patient based on the preop
13 aberrations. So it's important to note that the ablation
14 pattern is unique and is based on the preop aberration
15 measurement for each eye.

16 DR. WEISS: Dr. Bullimore, you're obscuring the
17 view.

18 DR. BULLIMORE: I apologize.

19 DR. PETTIT: So the technology requirements in
20 order to effect this type of treatment, there's two
21 components. First, we have to have a wavefront system
22 that's capable of measuring the higher-order aberrations
23 obviously accurately, quantifying the wavefront aberrations
24 in the patient, and we also need the wavefront device to
25 accurately register where exactly on the eye those

1 wavefront aberrations came from.

2 The treatment laser we use employs an active
3 high-tracking system to stabilize the eye during the
4 surgery, compensate for patient eye movement, and allow us
5 to deliver the customized ablation profile as accurately as
6 we can. We use a small Gaussian excimer beam to precisely
7 scope the subtle contours in the corneal surface and we use
8 fairly sophisticated software algorithms to convert the
9 wavefront data into the appropriate treatment profile.

10 So what is wavefront sensing? I'd like to give
11 just a very simple introduction so we all know what we're
12 talking about here. Simply put, wavefront sensing is a
13 measurement of how the eye operates as a integrated optical
14 system, and the wavefront device gives you a detailed
15 refractive map over the pupil of the eye. We think in
16 simple terms of how a theoretical perfect eye sees the
17 world. When a simple perfect eye looks at a far-off
18 target, the light from each point in that target enters the
19 eye as a bundle of parallel rays. The wavefront is the
20 surface perpendicular to each of those rays. So in the
21 case of a perfect eye, the wavefront entering the eye from
22 a distant target is a flat wavefront and that flat
23 wavefront is well focused down to a very small spot back on
24 the fovea.

25 Now, in the case of myopia, which is the

1 indication we're considering this morning, when a myope
2 looks at a distant object, those parallel light rays are
3 not well focused on the retina. They're focused somewhere
4 in front of the retina, but by the time the light reaches
5 the retina, they're blurred out. A myope can see clearly
6 if the target is up close to the eye. So you imagine a
7 nearby point source, the light rays from that nearby source
8 enter the eye as a diverging bundle. The wavefront in that
9 case is part of a spherical surface and that spherical
10 wavefront is then well focused down to a small spot on the
11 retina.

12 When we perform classical vision testing, we're
13 in a sense doing a primitive form of wavefront sensing.
14 We're asking what combination of spherocylindrical lenses
15 do we need to put in front of the eye so that that flat
16 wavefront from a distant source is as best focused as
17 possible back on the retina? The limitation of that, of
18 course, is that there are higher-order aberrations that are
19 known to exist in the eye and these cannot be characterized
20 with simple spherocylindrical lenses. This is an example
21 of spherical aberration whereby the periphery of the eye is
22 more refracting, has more myopic power than the central
23 part, and the rays are blurred out in the retina. Another
24 common example is coma and coma can be simply thought of as
25 one side of the pupil being slightly more myopic than the

1 average and the other side being more hyperopic than the
2 average. Again, it causes blurring, puts the light back on
3 the retina, and this can't be characterized well with
4 simple lenses.

5 When we perform wavefront sensing using the
6 Shack-Hartmann approach, we take advantage of the fact that
7 light is reversible. If we want to study how light gets
8 from Point A to Point B, being refracted at various
9 surfaces within the eye, we can instead measure how light
10 travels from Point B back to Point A. So what we do is we
11 have the patient look into the device and then we shine a
12 narrow eyesafe probe beam into the eye and illuminate a
13 small patch back on the fovea. Some of that light is
14 scattered back out of the eye just like when you do flash
15 photography and you get a red eye effect and outside the
16 eye, we now have a reemitted wavefront which is just the
17 time-reverse process of how the myopic patient sees the
18 world. So now we have that same spherical surface but it's
19 traveling out of the eye rather than into it.

20 So how do we measure what that wavefront looks
21 like? Well ,inside our wavefront sensing device, there's a
22 group of relay optics, so that this plane over here is
23 imaged over here right at the entrance space of the actual
24 wavefront sensor. So whatever the wavefront is doing as it
25 exits the eye, it does again over here at the wavefront

1 sensor itself. Now, let's zoom in on the wavefront sensor.
2 This slide shows a simple myopic wavefront impinging on the
3 wavefront sensor itself and I've isolated one lens lit up
4 here. You're seeing part of the wavefront go through that
5 lens lit. There's an array of microlenses with a CCD
6 camera screen sitting at a fixed distance behind it and
7 through one of those lenslets, this piece of the wavefront
8 is being focused to this camera screen at this point here.

9 Now, if the wavefront was perfectly flat, that
10 light would have traveled straight through and then hit the
11 screen over here. We want to measure or actually describe
12 what the wavefront looks like in some mathematical term, W
13 of X and Y , where X and Y are the transverse pupil
14 coordinates, and let's consider that single lenslet sitting
15 at the location of Y not. We know that the wavefront piece
16 that went through that lenslet traveled this distance to
17 get back to the CCD camera and it's laterally displayed
18 from its ideal location by ΔY . From that information,
19 we can calculate the slope of the incident wavefront at
20 that lenslet and by doing this at a large number of points
21 across this lenslet array, we're actually able to rebuild
22 the shape of the original wavefront.

23 This is an example of the CCD camera screen
24 showing you a picture of the focused light dots. The
25 software in the wavefront device goes in and finds the

1 centers of all those dots, figures out how they're related
2 to one another, and more importantly which lenslet each
3 one of them came through them, and from this picture and
4 the processing of the information, we calculate the shape
5 of the original wavefront.

6 Now, we have to have a set of mathematical
7 tools to describe what does the wavefront look like. It's
8 going to be a complicated surface and we need to describe
9 it. We use what's called Zernike polynomials, which are a
10 convenient mathematical basis set for describing visual
11 aberrations. There's an infinite number of these. They
12 come in orders which are shown by the different layers in
13 this pyramid. So there's an infinite number of these going
14 off to the bottom there and you can see as you work your
15 way down the pyramid, the orders become progressively
16 higher and the shapes become progressively more complex.
17 These second-level, second order aberrations are closely
18 associated with the conventional spherocylindrical errors
19 in the eye, and then the higher-order aberrations
20 correspond to the lower layers in that pyramid.

21 Now, from the wavefront information, we can
22 estimate the optical performance of the eye. We can take
23 the wavefront and calculate what's called the point spread
24 function which is an optical analysis of what a distant
25 point source would look like on the retina. So in the case

1 of a perfectly flat wavefront over a big pupil, that
2 wavefront is focused on the very tight spot, and this is a
3 simulated optical image of what the retinal image should
4 look like when the patient looks at approximately 20/16
5 through through 20/10 lines on an eye chart. They're
6 slightly blurred out due to the effects of defraction, but
7 this is about as well as the optics can do for a 6.5mm
8 pupil.

9 I'll just show you a couple more examples.
10 This is myopia, that simple spherical shape. The way to
11 think about it, the wavefront of the eye is down here and
12 the wavefront's propagating up out of the eye. Myopia
13 causes, as you can see, blurring of the optical image on
14 the retina. I should point out that this does not take
15 into account retinal effects or neuroprocessing. We don't
16 know exactly how well an eye cognitively could see the eye
17 chart, but this is a simulation of what the optics of the
18 eye produce on the retina.

19 I'll just show you one more. This is vertical
20 coma, and you can see that in this case, you actually can
21 read all these letters, but there's this sort of comet-like
22 tail extending in the downward direction due to the coma
23 present in the eye.

24 In addition to being able to provide a detailed
25 mathematical description, we also want to have a simple

1 single parameter to characterize, well, just how distorted
2 is this wavefront? So, we use what's called the RMS error.
3 You can hear this mentioned several times this morning.
4 The RMS error is simply the standard deviation of the
5 wavefront relative to that idealized flat profile. So, if
6 the wavefront is in fact perfectly flat, the RMS error is
7 going to be zero, and as the wavefront becomes more and
8 more distorted, the RMS error becomes progressively more
9 positive. That's a very simple overview of wavefront
10 sensing, and there's a lot of details that I simply don't
11 have time to go into in this one hour.

12 So in a clinical wavefront sense, we need a
13 little bit more equipment to actually be able to measure a
14 patient accurately, and so this slide shows the principle
15 optical components in our wavefront-sensing device. The
16 eye is sitting here looking into the instrument. We
17 obviously have to give the patient a target so they know
18 where to look, and these are myopic patients, so they don't
19 see very well before surgery. So our target has an
20 adjustable focus mechanism to correct for the preoperative
21 myopia and actually to fog the eye slightly to try to
22 minimize any accommodation effect as the patient looks in
23 there.

24 We have a video camera that's staring out at
25 the eye that helps us position the eye for the measurement

1 and equally important, it helps us record exactly how the
2 eye was positioned and how it was rotated at the instant
3 that the wavefront data was taken. We've already talked
4 about the probe beam and the wavefront-sensing pathways.

5 So how do we perform a surgical wavefront
6 measurement? When a patient comes into the clinic on the
7 day of surgery, the first thing that actually happens
8 before they have any dilation applied to the eye is we have
9 them sit down at the wavefront sensor and we take a video
10 snapshot, so they look into the device and we use that
11 video path to capture a frozen picture of their eye under
12 daytime illumination conditions, and we do this because we
13 want to record where their daytime natural pupil center
14 sits relative to their limbus. We do this by asking the
15 clinician to align two software reticles in the frozen
16 video image, one to the limbus and one to the pupil, and
17 having done this, our software now knows where the daytime
18 pupil center is relative to the limbus and that's going to
19 be our anchor point later on.

20 The patient then goes off and has the eye
21 dilated and in this trial, we used a combination of
22 tropicamide and phenylephrine. Immediately before the
23 wavefront measurement, the surgeon applies two ink marks
24 using a standard eye-marking pen to the sclera just outside
25 the limbus. The patient then sits down at the wavefront

1 sensor. They're positioned appropriately. They view the
2 target. We fog the eye, and then we take five repeat
3 wavefront measurements in relatively quick succession. At
4 the instant that each of those wavefront measurements is
5 taken, the video image is frozen. So we have a frozen
6 video image that's synchronized with the wavefront capture,
7 and we ask the technician in that frozen image to align two
8 reticles, one to the limbus, so that's an elliptical
9 reticle, and we also have a linear reticle that they're
10 supposed to draw through these applied ink marks, and with
11 this information and the wavefront data, we now know
12 exactly where the wavefront came from on the eye and the
13 cyclotorsional angle of the eye at the measurement time.

14 The five measurements are then automatically
15 analyzed and the two outliers are rejected based on a
16 statistical analysis of the RMS errors. The remaining
17 three are then compared for consistency and then averaged
18 together to make a final composition wavefront and this is
19 what we actually base the surgery on. As a final sort of
20 sanity check, we can back calculate the effect of clinical
21 prescription from the wavefront data and compare that to
22 what was measured at the foropter, and in this trial, it's
23 worth noting that we require both the sphere and cylinder
24 calculated from the wavefront had to agree within 1 diopter
25 with what was measured at the foropter.

1 The wavefront and the geometry information are
2 then transferred electronically to the treatment laser.
3 Our treatment device again employs a blind spot, relatively
4 small, 193nm excimer laser, uses an active eye-tracking
5 system to stabilize the eye during the treatment, and it's
6 currently approved for all conventional treatment types of
7 refractive error.

8 The LADARVision treatment device actually takes
9 the wavefront data and calculates the appropriate ablation
10 profile. The patient lies down, is prepared for the
11 surgery, and in the tracked image, once the patient and the
12 doctor are ready, in the tracked image of the eye, a single
13 linked reticle, a combination of the elliptical limbus
14 reticle and the linear cyclotorsional reticle, comes up in
15 the tracked image screen and the clinician then aligns
16 these to the anatomical features on the eye and that's how
17 we register the ablation profile correctly both in position
18 and cyclotorsional angle.

19 The device requirements on the wavefront-
20 sensing unit itself then. First, it must record the
21 natural pupil limbus geometry. It must measure wavefronts
22 up to at least the fourth Zernike order which is what we
23 used in this study, must be able to measure pupils in
24 excess of 7mm in diameter, obviously has to have a
25 validated accuracy in wavefront measurement performance,

1 must record the geometry of the wavefront data relative to
2 the limbus and the cyclotorsion features, these are these
3 applied ink marks that I referred to, and it must be able
4 to obviously export the wavefront and the geometry data in
5 a format compatible with the LADARVision system.

6 I'd now like to turn the podium over to Dr. Dan
7 Durrie who's going to summarize our safety data.

8 DR. DURRIE: Thank you very much.

9 It's a pleasure being here, and I'm Dan Durrie,
10 and I'm one of the investigators in this clinical trial,
11 and I'm a consultant for Alcon and I'm also a paid
12 consultant for a competing technology.

13 I'd like to review the safety criteria of this
14 particular study that's under question. First, I'd like to
15 clarify a little bit of the two groups that we'll be
16 talking about. The safety cohort includes 426 eyes which
17 includes the astigmatism cases that we're not asking for
18 approval for today but are included in the whole safety
19 cohort. The primary efficacy cohort as has been shown
20 before is 139 eyes and were the ones without significant
21 astigmatism. This is based off the manifest refraction,
22 and as I go through this, I'll be showing any differences.

23 First, accountability is always important with
24 any clinical trial and this was excellent. It's always
25 great to see a 100 percent down at the bottom of the

1 accountability chart. Unfortunately, there were two eyes
2 that were lost because of the death of a patient to colon
3 cancer, but it was a 100 percent of all patients who were
4 available were accountable at all visits. So therefore, we
5 can have at every visit 139 eyes in the spherical cohort
6 and in the safety cohort of all the eyes is all the eyes
7 available except the two that were lost at the six-month
8 visit.

9 As far as the demographics, and I will compare
10 the two groups, they are very similar. Between the two
11 groups, the points of interest are the fact that this was
12 primarily a male study and as with most excimer laser
13 studies, it's predominantly Caucasian and also we note that
14 most of these patients were soft contact lens wearers. The
15 age was in the upper 30s and similar between the two
16 groups. Also as far as the amount of correction that was
17 attempted, the only difference between the two groups
18 really was the fact that there was the spherical group did
19 not have the cylinder as previously described but other
20 than that, the average amount of myopia and the range was
21 similar in the two groups.

22 Now, again, I'm going to talk about safety and
23 this is the total group with a 100 percent follow-up of
24 those available. We're all familiar with the guidance
25 documents. We've reviewed these studies before, and if

1 you're using the criteria in the lines lost of best-
2 corrected vision, we can see that this easily meets all the
3 previously discussed guidelines. There was one eye that
4 did have vision that was less than 20/25 that was 20/25
5 preop and the three eyes that had loss of two lines of
6 best-corrected vision but no eyes that lost more than two
7 lines.

8 Looking at the best-corrected data in a
9 different way and overall lines lost or gained as a
10 clinician, we like to see the graph leaning to this side
11 from preop to postop which means that there is more lines
12 gained than lost, and as you can see at the six-month
13 follow-up, 37 percent gained one line or more and only 9.4
14 percent of patients lost any vision. For most studies,
15 these are fairly even and now with these newer
16 technologies, we're seeing the graph moving in that
17 direction.

18 It also showed in looking at the best-corrected
19 vision comparing preop to postop and looking at the high
20 level of vision correction of 20/12.5, that we had a
21 doubling of that in the study from preop to six-months
22 postop and 20/16 increased from preop to postop. So the
23 best-corrected vision overall was increased with this
24 study.

25 In terms of complications and the way these are

1 reported, they're reported at any visit, any time. So
2 that, the typical things we'd see with LASIK surgery, small
3 amount of DLK, epithelial ingrowth, ghosting images, and
4 some corneal edema, and there were other findings that were
5 listed that were below the .2 rate that are listed at the
6 bottom here. Nothing out of the ordinary here for a LASIK
7 trial or LASIK clinical.

8 Now, what happened to those patients who
9 reported complications? All but two eyes that had best-
10 corrected vision of 20/20 or better and a 100 percent had
11 20/32 or better at the last visit. Also, all complications
12 resolved, except for four eyes, and these were one patient
13 with epithelial ingrowth and three patients with ghosting
14 images. All the reported DLK and epithelial ingrowth were
15 Grade 1 or less at all of the reporting visits.

16 In regards to adverse reactions, those related
17 to the device, there was recalcitrant DLK associated with
18 blepharitis in two eyes of one patient and there was one
19 miscreated flap. The patient was exited from the study and
20 had successful recut LASIK surgery with a conventional
21 laser. There were also some unrelated to the procedure.
22 There was one patient, I told you before, that died of
23 colon cancer and one patient that developed multiple
24 sclerosis during the procedure. There was one retinal
25 horseshoe tear which was felt to be unrelated to the

1 procedure by the retinal specialist. What happened to
2 those patients, of all those patients who had adverse
3 reactions, a 100 percent of them were 20/16 or better at
4 the last reported visit, and other than the multiple
5 sclerosis, all of the other adverse reactions resolved.
6 In regards to intraocular pressure, corneal haze or other
7 slit lamp findings, there was nothing unusual in this study
8 and there was nothing that was out of the ordinary that we
9 would expect. So the overall safety of this study was
10 extremely good.

11 Going to the spherical cohort, this is 139
12 eyes, just quickly showing you that if we look at that
13 group, there was all zeroes on the parameters for the FDA
14 guidance and no eyes had worse than 20/20 vision. In
15 regard to complications, it's the same distribution but
16 slightly less in this group, and there was no adverse
17 reactions in the spherical myopia group that we'll be
18 discussing for efficacy.

19 Therefore, the safety criteria in this study
20 meets or exceeds the guidelines for loss of best-corrected
21 vision, best-corrected vision worse than 20/40 and induced
22 cylinder and the incidence of adverse reactions were
23 overall, and there was no demonstrated significant safety
24 concerns.

25 I'd like to introduce Dr. Omar Hakim who will

1 be talking about efficacy and this again to define is 139
2 eyes that are in the efficacy group.

3 Thank you.

4 DR. HAKIM: Thanks, Dan.

5 Hi. My name is Dr. Omar Hakim. I'm Medical
6 Director for TLC Laser Eye Centers in Canada. I've been
7 performing laser refractive surgery since 1994 and custom
8 ablation surgery since May of 2000, using a variety of
9 different platforms, and actually I had my own vision
10 corrected with LASIK in 1998. I am a consultant on Alcon's
11 Refractive Medical Advisory Board and travel expenses for
12 this meeting were paid for by Alcon.

13 I've been asked to present the effectiveness
14 outcomes for the 139 eyes in the study. Preoperably, these
15 139 eyes had up to 7 diopters of myopia and less than half
16 the diopter of astigmatism. First, we'll review the
17 manifest refraction spherical equivalent results and we'll
18 see that following surgery here on the left, that 83.5
19 percent of the eyes had an MRSE within half a diopter of
20 emmetropia at one-month postop and 74.8 percent at six
21 months following surgery. Fully 97.1 percent of the eyes
22 were within 1 diopter of emmetropia at one month and almost
23 96 percent at six months. Both of these groups clearly
24 exceed the FDA Guidance Document guidelines calling for 50
25 percent of patients to be within half the diopter of

1 emmetropia and 75 percent within 1 diopter of emmetropia.

2 This graph shows the attempted versus achieved
3 correction at six months of all 139 eyes and it really
4 demonstrates the vast majority of eyes fall within a 1
5 diopter bracket on each side of emmetropia shown by the
6 dashed line here. At the higher ranges of correction,
7 there are three eyes that fall outside the 1 diopter
8 bracket. However, even at these higher levels of myopia,
9 excellent results were still reported with uncorrected
10 visual acuity of 20/25 or better in 92 percent of the eyes
11 between 5 and 5.99 diopters and 75 percent of the eyes
12 above 6 diopters or above. In fact, in this higher myopic
13 group, 75 percent of the eyes were still within half a
14 diopter of emmetropia. However, overall, there was a small
15 amount of undercorrection as shown by this slide. When we
16 reexamined the results around this mean MRSE line of -0.37
17 diopters, we see that actually more than 90 percent of eyes
18 were within half a diopter of this mean value, reflecting
19 really a very high level of precision and reproducibility
20 of result.

21 This chart looks at the MRSE over time and
22 shows excellent stability from one-month postop with
23 refractive MRSE of -0.27 diopters, -0.35 at three months,
24 and -0.37 diopters at six months. Again, the notable
25 precision of the MRSE is shown in the standard deviation

1 values which range from 0.34 diopters to 0.42 diopters of
2 standard deviation. Of course, surgeons routinely make
3 adjustments in treatment based on environmental factors,
4 such as temperature and humidity, and surgeon-specific or
5 site-specific factors because we know these influence the
6 accuracy of our outcome in every-day surgery. The accuracy
7 of these results reflected in the fact that we had a mean
8 undercorrection of -0.37 diopters were limited by the study
9 protocol because it restricted clinicians from making these
10 usual treatment adjustments. This undercorrection could be
11 dealt with by usual nomogram, adjustments by the surgeon or
12 could be incorporated into the software and as commented on
13 by Dr. Eydelman in her medical officer's review, Alcon has
14 already initiated a clinical evaluation of the minor
15 software adjustment to address this.

16 In terms of stability of MRSE, then we see that
17 100 percent of eyes between the one- and three-month visits
18 and three- and six-month visits had less than 1 diopter
19 change in MRSE, surpassing the FDA Guidance Document
20 guidelines. In fact, the mean change was only .07 diopters
21 in the one-to-three-month period and -0.3 diopters in the
22 three-to-six-month period, which translates into a mean
23 change per month of -0.035 diopters and in the three-to-
24 six-month group of -0.01 diopters, really incredibly good
25 stability.

1 This chart then shows the percentage of eyes
2 achieving uncorrected visual acuities of 20/20 on this side
3 and 20/40 on the right side. We see that 86.3 percent of
4 eyes at one month had 20/20 uncorrected acuity and 79.9
5 percent at six months had 20/20 acuity. Looking at 20/40,
6 we see that 99.3 percent of eyes had that acuity level at
7 one month and 98.6 percent, almost 99 percent, at the six-
8 month visit. Again, these results exceed the FDA guidance
9 calling for 85 percent of eyes to have uncorrected acuity
10 of 20/40 or greater.

11 It really should be noted that these excellent
12 uncorrected acuity results were obtained despite the mean
13 undercorrection of -0.37 diopters that we noted previously
14 and a reduction in this undercorrection along with the
15 excellent precision of effect as shown by the MRSE results
16 should provide even better uncorrected visual acuity
17 results with nomogram adjustment, and in fact, at one
18 month, 59 percent of eyes had uncorrected visual acuity
19 equal to or better than their preoperative best spectacle-
20 corrected acuity and at six months, this figure was still
21 52.5 percent.

22 In summary then, this study of CustomCornea
23 Wavefront-Guided Ablation has demonstrated uncorrected
24 visual acuity results and accuracy and stability of MRSE
25 results that exceed those called for by the FDA Guidance

1 Document. Patients were also asked to grade any ocular or
2 visual symptoms. With regards to ocular symptoms, only 2.2
3 percent of patients noted significantly worse dryness of
4 their eyes and only 1.5 percent noted significantly worse
5 burning or gritty sensation with their eyes at six-month
6 postoperatively, and in fact, many patients actually noted
7 they had significant decreases in symptoms, including
8 significant, for example, you look at dryness, 8.1 percent
9 of patients said their dryness of their eyes was actually
10 significantly better following surgery.

11 On evaluation of the visual symptoms, 2.9
12 percent of patients noted they had significantly worse
13 blurring of vision and 0.7 percent of patients noted they
14 had significantly worse night driving difficulty, double
15 vision or fluctuation of vision. What's interesting is
16 that if we looked at the mean MRSE on patients who said
17 they were better or significantly better, that ranged from
18 .26 diopters to -0.36 diopters, and in the group who noted
19 that their symptoms were worse or significantly worse, that
20 range was -0.46 to -0.7 diopters. So that, again, further
21 improvements in undercorrection may further improve upon
22 this already-low level of symptoms.

23 As a surgeon, all these visual symptoms are
24 important, but the most concerning to me, you know, really
25 the top three, glare, halos, and night driving

1 difficulties, and, you know, here we see that although the
2 numbers are small, more patients actually said their
3 symptoms were significantly better than significantly worse
4 following surgery.

5 At three months following surgery then, 85
6 percent of patients said that they were satisfied or
7 extremely satisfied overall and at six months, this was 79
8 percent. Again, the MRSE in patients who were satisfied or
9 extremely satisfied was -0.3 diopters and in the
10 unsatisfied or extremely unsatisfied group, that was -0.93
11 diopters, and again a reduction in this undercorrection
12 should shift the cohort towards even higher patient
13 satisfaction rates. Almost 90 percent of patients
14 described their quality of their vision as being equal,
15 better or significantly better following surgery at both
16 the three- and six-month intervals, and over 95 percent of
17 patients at three months and 94 percent of patients at six
18 months had no need for distance correction of any kind.

19 So again, in summary, the study of CustomCornea
20 Wavefront-Guided Treatment has clearly exceeded the
21 performance guidelines laid out in the FDA Guidance
22 Document in terms of uncorrected visual acuity and accuracy
23 and stability of postoperative manifest refractive outcome.

24 Dr. George Pettit's now going to return and
25 discuss the wavefront and higher-order aberration outcomes

1 for the cohort.

2 Thanks, George.

3 DR. PETTIT: Now we're going to talk about the
4 higher-order aberration changes that we saw in this myopic
5 cohort.

6 This slide summarizes the changes in the
7 various higher-order parameters. So the third- and fourth-
8 order aberrations are considered along the horizontal axis
9 here. We're looking at the total higher-order aberration
10 content and then the individual content from the various
11 third- and fourth-order aberrations. The vertical axis
12 indicates the magnitude of the different wavefront
13 parameters. The blue bar indicates what they were
14 preoperatively, so that's the starting baseline level, and
15 then the green, yellow, and red bars indicate the one-,
16 three- and six-month postop measurements on this cohort.

17 Now, I'd like to just note there's 137 patients
18 in this table. The entire cohort was 139 eyes, but there
19 were two eyes that missed one of these wavefront measures
20 at some interval. So, there's 137 eyes considered here and
21 that's why the N is slightly smaller. The little asterisks
22 indicate those changes that were statistically significant
23 between the preop and the six-month postop interval and you
24 can see that for almost all of the aberration parameters,
25 they're actually slightly higher by a statistically

1 significant amount after surgery than before. The trifoil
2 is actually less after treatment but it's not a
3 statistically significant difference.

4 Now, I should just also mention this is based
5 on a 6.5mm wavefront analysis diameter. This isn't a
6 surprising finding. It's well known that LASIK tends to
7 increase the higher-order aberrations. So an important
8 question to ask is how does this compare to conventional
9 surgery? We have a comparative conventional cohort. In
10 the early phase of this, our clinical trials, we ran a
11 bilateral study where we had a contralateral control arm.
12 What I mean by that is that patients would come in, if they
13 met all of the entry criteria to be enrolled in the study,
14 they were randomly selected, so that one eye received
15 conventional treatment with our system and the other eye
16 received customized treatment. Again, the eyes were
17 randomized.

18 Of that comparative arm, 50 eyes actually meet
19 the criteria of being myopic with less than half a diopter
20 of cylinder. So we went back and looked at all of the
21 conventional patients that we treated in the early phases
22 of the study and found out that 50 eyes met the criteria
23 for our current conditions of approval, and you can see
24 that the refractive parameters for that conventional
25 comparative group match up very, very well with the primary

1 Custom cohort. Those patients also had a treatment optical
2 zone of 6.5mm in diameter and we have wavefront data
3 measured in the same way available preoperatively as well
4 as the one- and three- and six-month postop intervals.

5 If we look at the aberrations in these two
6 groups, the Custom shown by the blue bars and the 50-eye
7 conventional comparative cohort shown by the red bars,
8 preoperatively, the aberration content's relatively similar
9 between these two groups. There was a small but
10 statistically significant difference in the spherical
11 aberration term but all other parameters were well matched.
12 That's not the case when we look six months after surgery.
13 Again, these asterisks indicate anything that exceeded the
14 P value for statistical significance of being less than .05
15 and you can see that the total higher order, coma, trifoil,
16 spherical aberration, and tetrafoil, are all significantly
17 lower six months after wavefront-guided treatment than six
18 months after conventional surgery. Secondary astigmatism
19 term was also lower in the Custom eyes but that wasn't a
20 statistically significant difference.

21 Now, we've tried to and we've worked with the
22 agency and tried to come up with a way of describing what's
23 the optical impact of the magnitude of these differences in
24 the higher-order aberrations, and this again is an optical
25 simulation of how a patient might see the eye chart under

1 best-corrected vision where the lower-order aberrations are
2 removed. On the left, they're left with the postop
3 aberration mean for wavefront-guided treatment, and on the
4 right, they're left with the mean values for conventional
5 surgery, and you can see there's a modest but definite
6 difference with the optical quality being better in the
7 wavefront-guided approach.

8 Now, that's all based on mean values. We also
9 looked at on an individual patient basis what percentage of
10 patients exhibited an actual decrease in the higher-order
11 aberration parameter after surgery as compared to before,
12 and so the different lines in this table indicate the
13 various ways of looking at the higher-order aberrations for
14 third and fourth order, and this middle column indicates
15 the percentage of wavefront-guided eyes that showed a
16 reduction in that particular parameter, and on the right-
17 most column, we're looking at the conventional eyes, and
18 you can see that for most of these parameters, much higher
19 percentage of wavefront-guided patients actually showed a
20 decrease six months after treatment as before surgery and
21 that's not true in the conventional eye. The percentages
22 are much lower.

23 I'd like to now invite Dr. Steve Brint to come
24 up and talk about the clinical implications of the
25 wavefront correction.

1 DR. BRINT: Thank you, Dr. Pettit.

2 I'm Steve Brint from New Orleans, Louisiana.
3 I'm in private practice and on the faculty of Tulane
4 University. I likewise am a member of the Alcon Medical
5 Advisory Board. I've been performing LASIK since 1991 and
6 am a LASIK patient myself, and I'm also the Medical Monitor
7 of a competing laser technology and my expenses for this
8 trip were also compensated by the sponsor.

9 As clinicians, we know that prior studies of
10 conventional LASIK in general have shown that higher-order
11 aberrations, particularly spherical aberration which is
12 frequently linked to poor night vision, occasionally may be
13 increased after conventional LASIK. These increased
14 higher-order aberrations after conventional LASIK are
15 pupil-size dependent with larger pupils showing decreased
16 retinal image quality as measured by point spread function
17 and modulation transfer function and visual performance as
18 measured by the clinically useful contrast sensitivity
19 testing and low-contrast visual acuity testing.

20 David Williams' group at the University of
21 Rochester has done work in this area and has shown that
22 correction of these higher-order aberrations using an
23 adaptive optics system is able to improve the visual acuity
24 and especially the contrast sensitivity.

25 As Dr. Pettit just mentioned, we do have this

1 comparative conventional cohort of 50 spherical myopic eyes
2 that was derived as he mentioned which has comparable
3 demographics and virtually identical preoperative
4 refractive error. Quality of vision in this group was used
5 just as we did in the wavefront-guided eyes using tests to
6 determine the visual performance under these low-contrast
7 and mesopic situations. I think we all realize as
8 clinicians that we've done a very good job of getting good
9 quantity at vision and now, as has been mentioned
10 throughout the morning, the goal is not only to improve the
11 quantity but improve the quality of vision for our
12 patients.

13 So the contrast sensitivity testing was done
14 using the VectorVision Chart, the CSV1000, measured at 3,
15 12 and 18 cycles per degree. This was done in both every-
16 day full illumination to simulate a photopic situation as
17 well as in a room with total darkness, other than the light
18 coming from the eye chart, to simulate a mesopic light
19 situation, and a neutral density filter was placed in front
20 of the eye which only transmitted 3.16 percent of the
21 light. We know that greater higher-order aberrations are
22 seen in these larger dark-adapted pupils as opposed to the
23 smaller light pupils.

24 Previous FDA studies have used a contrast
25 sensitivity definition as a clinically significant change

1 of greater than 0.3 log units at two levels from
2 preoperative at two or more spatial frequencies, and this
3 is used to distinguish measurement noise from actual true
4 sensitivity change. What we saw in our spherical cohort as
5 regards to percentage of eyes with this clinically
6 significant change of photopic contrast sensitivity, here
7 in the Custom group of spherical eyes, we see two to three
8 times gain as opposed to loss of contrast sensitivity and
9 clinically significant contrast sensitivity, and in the
10 conventionally treated eyes, we see actually at both the
11 three- and six-month interval no gain and some loss of this
12 clinically significant photopic contrast sensitivity.

13 The full eye larger cohort is nice in that it
14 confirms what we saw before with a tendency towards gain as
15 opposed to loss in the wavefront-treated eyes and loss as
16 opposed to gain as treated in the conventionally treated
17 eyes.

18 Looking at photopic contrast sensitivity
19 another way, at individual spatial frequencies, we see
20 conventionally treated eyes, no mean log change in the
21 lower spatial frequencies and significant loss in the
22 higher spatial frequencies. In the Custom group, however,
23 we see gain across the board at all log changes at all
24 spatial frequencies and this is statistically significant
25 at the higher spatial frequencies of 12 cycles per degree

1 and 18 cycles per degree, and this is confirmed looking at
2 the larger full eye group. In the conventional group, we
3 see across the board at all spatial frequencies a trend
4 towards loss of photopic contrast sensitivity as measured
5 in log units and somewhat modest gain in the Custom-treated
6 eyes.

7 As regards the clinically important mesopic
8 contrast sensitivity, the large pupil at night time
9 contrast sensitivity, important to our patients in tasks
10 such as driving at night, we see, also, two to three times
11 the number of eyes gaining clinically significant contrast
12 sensitivity in the Custom group, both early and at the late
13 testing intervals, while in the conventionally treated
14 eyes, initially there's more loss. However, there is
15 recovery in the later interval which simulates what we see
16 in the six-month conventional group, what we see in the
17 Custom group at three months, suggesting that these
18 customized enjoy this improved mesopic contrast sensitivity
19 throughout their recovery period without going through the
20 decreased period that our conventional eyes have. This is
21 confirmed once again in the all-eye group with improved
22 mesopic contrast sensitivity and decreased at three months,
23 recovers at six months, but once again is maintained
24 throughout the entire postoperative period in our Custom
25 wavefront-treated eyes.

1 Looking at mesopic contrast sensitivity in our
2 conventional eyes, we see a modest gain at all spatial
3 frequencies at six months. However, in the Custom-treated
4 eyes, we see a much larger gain in the Custom-treated eyes
5 as measured in log units over this period of time at six
6 months at all spatial frequencies, and this once again is
7 confirmed when we look at the larger cohort of all eyes with
8 these more difficult toric prescriptions.

9 Low contrast visual acuity was measured using
10 the ETDRS eye chart. We're all familiar with the standard
11 high contrast eye chart. You can all see how extremely
12 difficult it is to see this 10-percent low contrast eye
13 chart that was actually viewed in a room with ambient very
14 dim light, and it should be noted that even in our
15 preoperative best spectacle-corrected vision patients, only
16 8.6 percent of patients were able to read the 20/20 line
17 preoperatively.

18 Looking at the change in the low contrast best-
19 corrected vision in our spherical group, we see that there
20 is more gain than loss at both the early as well as the
21 later time interval, more so than we see in the
22 conventional eye group, but more importantly, statistically
23 significant, we see that there is significantly less loss
24 in the Custom eye group at three months as compared to the
25 conventional eye group, less loss of one or more lines from

1 preoperative value at 22 percent as opposed to 36 percent,
2 and this approaches statistical significance. If we look
3 at all the eye group, we see that this is even more
4 statistically significant at both the early as well as
5 later intervals, less loss of one or more lines of low
6 contrast vision as compared to what we're used to seeing in
7 our conventional eyes.

8 So in summary, of the Custom spherical eyes
9 evaluated at the six-month time gate, as looking at
10 photopic contrast sensitivity, we had 2.2 percent gain as
11 opposed to 0.7 percent loss with a mean gain at all spatial
12 frequencies. Looking at clinically important mesopic
13 contrast sensitivity, we had 15.3 percent gain as opposed
14 to 5 percent loss, and once again mean gain at all spatial
15 frequencies, and with this extremely difficult low contrast
16 best-corrected vision, we had one or more lines gained in
17 38.8 percent of the patients as compared to 20.9 percent
18 loss.

19 So in conclusion, the Alcon CustomCornea System
20 is unique in that it's capable of measuring each of these
21 aberrations measured by the aberrometer and taking them and
22 registering them to each other so that we have a very
23 accurate composite aberrometry reading of both the low- and
24 high-order aberrations of the entire optical system which
25 we're then able to match and link and transfer to the

1 excimer laser, register it perfectly to the treatment of
2 the excimer laser, which then calculates and delivers a
3 specific ablation pattern unique for each individual eye.
4 The ablation pattern is uniquely determined from these
5 preoperative aberrations present in each individual eye.

6 Dr. Durrie has shown that wavefront-guided
7 CustomCornea treatment easily meets all the FDA guidance
8 criteria for safety with exceptional improvement of best
9 spectacle-corrected vision, especially at the extremely
10 high levels of acuity, 20/12, 20/16. Dr. Hakim has shown
11 that the CustomCornea treatment exceeds all the
12 effectiveness criteria as established by the FDA with a
13 very, very precise type standard deviation around the mean.
14 The Custom eyes have shown a consistent trend for more eyes
15 to have a clinically significant gain as opposed to loss of
16 both mesopic and especially photopic contrast sensitivity
17 and more eyes have shown a gain of one line or more of low
18 contrast best-corrected vision as opposed to loss.

19 Compared to these conventional eyes, the Custom
20 eyes have a statistically significantly better mean
21 photopic contrast sensitivity and as has been shown, we're
22 able to preserve in the Custom eyes this mesopic contrast
23 sensitivity at three months which is lost in the
24 conventional eyes, although it does recover, but it allows
25 the Custom eyes to enjoy excellent mesopic contrast

1 sensitivity throughout their recovery period and there is a
2 statistically significant lower loss of low contrast best-
3 corrected vision of one line or more.

4 We believe that wavefront-guided LASIK produces
5 an eye that's optically superior to conventional LASIK, and
6 for our patients, this means significantly less
7 postoperative aberrations, as has just been shown in Dr.
8 Pettit's presentation, and in these Custom eyes, we've seen
9 a significantly greater reduction in the higher-order
10 aberrations of virtually all of the specific types from
11 preop as compared to what was seen in the conventionally
12 treated eyes.

13 So I think from a clinical point of view, this
14 is something that, as has been discussed this morning, is
15 extremely important for improving the visual quality of our
16 patients and in the future for perhaps going back and
17 taking care of some of the problems as Mr. Link has
18 previously discussed.

19 I think we have a little bit of time left to
20 answer some of the questions that were specifically
21 addressed by the panel. I'll turn the podium back over to
22 Dr. Pettit.

23 DR. PETTIT: Thank you, Dr. Brint.

24 By my watch, I have eight minutes, and I'd like
25 to just touch on some of the questions that the FDA and you

1 panel members have raised in reviewing some of this
2 material.

3 I'd like to just start with Dr. Huang's review,
4 and he noted that at three and six months, after a
5 CustomCornea, only 78 and 82 percent of patients had a low
6 contrast UCVA of 20/40 or better and only 5.8 percent
7 achieved a low contrast UCVA of 20/20 or better and
8 expressed some concern about that. We went back and looked
9 at that data and found that that's actually a fairly
10 unremarkable -- this is a relatively difficult test, and
11 what I mean by that is if we look at the eyes
12 preoperatively, best-corrected visual acuity of 20/20 or
13 better, only 8.6 percent of the spherical eyes were able to
14 see the 20/20 line on the eye chart. Three months after
15 surgery, 5.8 percent of our patients were able to see that
16 line, and at six months, it's 7.9 percent low contrast UCVA
17 approaching their preop BCVA. Low contrast UCVA of 20/40
18 or better is actually slightly higher in the wavefront-
19 guided-treated eyes at both three and six months than it is
20 in our comparative conventional group.

21 In Dr. Bradley's review, he raised a very
22 interesting question, and I'm paraphrasing it slightly
23 here. He's basically asking are we correcting the
24 aberrations that were in the eye before surgery? Are we
25 compensating for the treatment-induced aberrations or are

1 we doing both? Well, in all honesty, we are attempting to
2 do both, so that the postop aberrations are as small as
3 possible.

4 When we began the wavefront development effort,
5 our initial aim was simply to treat the preop aberrations.
6 However, in our early trials, looking at the pre- and
7 postop wavefront data, it became very clear that the
8 aberrations induced by the surgery were also very important
9 and some, not all, but some of these surgical effects were
10 predictable, and therefore our Custom ablation algorithm
11 evolved from one that simply calculated the profile
12 directly from the wavefront data into one that took the
13 wavefront data but then made some adjustments to the
14 wavefront-based profile to compensate for predictable
15 surgical effects. I just want to stress that all 139 eyes
16 in this cohort were treated at the end of this process with
17 a consistent algorithm.

18 Dr. Bradley went on to offer one analysis
19 possibility that we could perform. How do we know if we're
20 treating the preop aberrations? He suggested we look at
21 the correlation between the aberrations before and after
22 surgery, and if successful, the wavefront-guided postop
23 aberrations should be uncorrelated with the preop eye-only
24 aberrations. We did that analysis, and the correlation
25 coefficients are shown here. This is looking at the

1 individual higher-order third- and fourth-order aberrations
2 between preop and three months after surgery, preop and six
3 months after surgery, and you can see down here for two of
4 these fourth-order aberrations, there is a modest positive
5 correlation coefficient, somewhere in the .4 to .45 range,
6 but in general, these numbers are pretty small. Postop
7 aberrations are not well correlated to the preop
8 aberrations. That's not ironclad proof that we're treating
9 the preop aberrations. I'm going to come back to that in
10 just a second and say a little bit more.

11 How do we know that we're treating the
12 surgically-induced aberrations? Well, we looked to see if
13 there were any significant correlations between the lower-
14 order aberration changes, the changes in the myopia, and
15 changes in the higher-order aberrations, and we found there
16 were no significant correlations between the lower- and the
17 higher-order aberration changes. We also looked to see if
18 there's any correlation between the clinical refraction
19 changes and changes in the higher-order aberrations, and
20 again we found no correlation between the higher-order
21 aberrations and the refractive treatment effect. Neither
22 of these findings is true for conventional surgery, and
23 we've submitted a large body of data to the FDA over time
24 showing the trends we see for conventional treatment. We
25 don't see such coupling, so-called, effects here.

1 What I'm really trying to say is best
2 summarized by this slide which I showed earlier. This is
3 looking at the higher-order aberrations six months after
4 surgery. The blue bars are the Custom outcomes, the red
5 bars are the conventional 50-eye outcome. The surgical
6 predictable effects that we include in the ablation
7 algorithm involve only the rotationally symmetric terms.
8 So on this chart, that only involved the spherical
9 aberration term here. So the fact that we're having
10 significant success in limiting the spherical aberration, I
11 believe, is due to the fact that we're compensating for the
12 surgical induction that would otherwise occur, but these
13 other higher-order aberrations are not rotationally
14 symmetric and there's no compensation mechanism folded into
15 the treatment profile to deliberately counter these. So
16 the fact that we see lower levels of these non-rotationally
17 symmetric terms to me, I think, is the best evidence that
18 by including them in the preop profile, we are effectively
19 treating them, although not eliminating them.

20 One of the questions that the FDA is posing to
21 the panel is are any of the differences between Custom and
22 conventional outcomes clinically and/or functionally
23 significant? I think Dr. Brint very nicely summarized what
24 we think are the significant differences. Compared to
25 conventional surgery, wavefront-guided-treated eyes have

1 significantly lower postop higher-order aberrations,
2 significantly higher percentage of eyes with an actual
3 reduction in various higher-order aberration parameters
4 relative to preop. They have a statistically significantly
5 better mean photopic contrast sensitivity and preservation
6 of mesopic contrast sensitivity at three months where we
7 see a consistent dip in the conventional treatments, and
8 they have a statistically significant lower loss of low
9 contrast BCVA defined as one or more lines.

10 Part of Question Number 3 for the panel asked
11 what information about the measurement, analysis and
12 correction of higher-order aberration is needed to
13 accurately inform physicians and prospective patients about
14 the safety and effectiveness? We've worked with the agency
15 to try to come up with a simple means of describing the
16 potential optical benefit of doing this type of surgery.
17 We're going back to that optical simulation of what a high
18 contrast chart would look like, and we come up with a
19 difference between wavefront-guided and conventional
20 treatment that corresponds approximately to about 2
21 diopters effective defocus blur. Under low contrast
22 conditions, this would be a different outcome certainly.

23 DR. GRIMMETT: That would be .2.

24 DR. WEISS: Yes, I think you misspoke. I think
25 you meant .2 rather than 2 diopters.

1 DR. PETTIT: I'm sorry. Point 2 diopters.

2 Absolutely. That would be very nice.

3 (Laughter.)

4 DR. PETTIT: We're not there yet.

5 Question Number 4 for the panel talks about the
6 refractive effects of correcting the higher-order
7 aberrations and states that these are smaller than the
8 effects of correcting the lower-order aberrations,
9 suggesting that relatively modest instabilities of sphere
10 and cylinder correction could disrupt the higher-order
11 corrections.

12 There's two points I just want to touch on
13 briefly here. Number 1. The wavefront-treated outcomes
14 have the same refractive stability as conventional surgery.
15 The higher-order aberrations in our wavefront population
16 are at least as stable as the aberrations in conventionally
17 treated eyes, and therefore we believe that modest versus
18 large amounts of these higher-order aberrations should be
19 beneficial to patients in the presence of refractive
20 instability in the postoperative course.

21 Another final point is we looked at are the
22 refractive instabilities somehow linked to higher-order
23 instabilities? In other words, if the patient's myopia is
24 changing in the postoperative interval, are the higher-
25 order aberrations changing in any corresponding fashion?

1 No. That should be defined as a correlation analysis
2 looking at the higher-order aberrations between three and
3 six months, and we saw no significant correlation, no
4 correlation larger than .18, between the refractive changes
5 and the higher-order changes. We also looked at changes in
6 the lower-order aberrations as measured by wavefront
7 device, compared those to the higher-order aberration
8 changes, and again saw no significant correlations.

9 That actually concludes our presentation.

10 Thank you very much for your attention.

11 DR. WEISS: I'd like to thank the sponsor for
12 their presentation.

13 We will take a 10-minute break. I'd ask
14 everyone to be back here promptly so we can start exactly
15 in 10 minutes, and we'll proceed with the meeting then.

16 Thank you.

17 (Recess.)

18 DR. WEISS: We will be starting now, if
19 everyone is now seated. We're going to proceed with, for
20 the next half hour, with the panel questions for the
21 sponsor, and then we'll have the FDA presentation.

22 I would first like to ask the sponsor two
23 questions. One. There was a cohort in which one eye had
24 conventional treatment and one eye had customized
25 treatment. Since a question the patients will ask is can I

1 notice any difference, aside from the numbers that we see,
2 can the patients notice any difference, were those patients
3 asked which eye they preferred, the customized corneal
4 treatment eye or the conventional eye?

5 DR. PETTIT: We'll see if we can -- I'm not
6 sure. Right. That's right. There are 50 eyes in the
7 conventional arm of the study that meet the spherical
8 definition. There are actually only 19 eyes that were
9 treated conventionally in the spherical group. There were
10 19 patients that were treated one eye conventional and one
11 eye with the Custom algorithm. So it's not 50 eyes. It's
12 not 50 patients. It's only 19 where they actually were
13 treated one eye one way and one eye the other, and we'll
14 see if we can find information on those particular 19 eyes
15 for you.

16 DR. WEISS: Okay. While we're waiting for
17 that, I had a second question on one of the higher-order
18 aberrations, the tetrafoil. Does this behave different
19 than the other aberrations? For example, the conventional
20 group actually had a higher percentage of eyes, 28 percent
21 of eyes in the conventional group had a reduction in this
22 particular higher-order aberration as opposed to 22 percent
23 in the Custom group. So the conventional group had a
24 greater percentage of people with reduction in this
25 particular type of aberration, and in addition, on Table 4,

1 there is a continued reduction in this aberration which is
2 statistically significant between three and six months. So
3 it still is changing after the "eye" has stabilized
4 refractively.

5 So I was wondering why. Is this different than
6 the others, and does this treatment not treat this
7 particular aberration, and if so, why not?

8 DR. PETTIT: Well, right. It's important to
9 note that the average value, looking across the groups, the
10 average values for the tetrafoil aberration were smaller in
11 the wavefront-guided populations than they were in the
12 conventional, but there was a slightly higher percentage of
13 actually conventionally treated eyes that showed a
14 reduction. So the mean value was less looking across the
15 entire group. Slightly higher percentage of eyes, though,
16 you're right, had a reduction in aberration.

17 The simple answer is the higher-order
18 aberrations, they all have slightly different optical
19 effects on image quality, and I think Dr. Burns or Dr.
20 Bradley can very objectively speak to which aberrations
21 potentially are the most detrimental. If I had to pick one
22 that I wouldn't worry as much about, it actually would be
23 the tetrafoil as opposed to some of the others. But
24 certainly the panel members can speak to that very
25 eloquently.

1 DR. WEISS: Thank you.

2 Dr. Huang?

3 DR. HUANG: I have a follow-up question on the
4 first question Dr. Weiss raised.

5 There were only 50 patients had a treatment in
6 one conventional and one eye with the Custom treatment, but
7 on the table presented by Dr. Brint earlier this morning,
8 there were several tables indicating that in all eyes
9 group, there were N equal to about 420 some odd eyes in the
10 Custom group and then in the conventional treatment, there
11 were N equal to about 130 some odd eyes.

12 Was there a mistake in terms of the tabulation
13 or was that --

14 DR. PETTIT: No, there was no mistake.

15 Again, we're seeking approval just for the
16 spherical cohort. So those are eyes with less than half a
17 diopter refractive cylinder. We've treated a much larger
18 population of patients and contrast sensitivity. We have
19 data on the entire eye cohort which includes myopes and
20 astigmats. So if we include astigmatism, we have data on
21 something like 426 wavefront eyes and 139 or something like
22 conventional eyes.

23 For contrast sensitivity, which has
24 historically been a safety parameter and is done under
25 best-corrected conditions, we presented that information as

1 supportive of the trends that we saw in the primary cohort.
2 So we have the primary cohort which is 139 Custom spherical
3 eyes as compared to 50 conventional spherical eyes, but
4 then supportive data on the larger body which includes
5 astigmats and the Ns are larger for that reason.

6 DR. HUANG: Thank you.

7 DR. WEISS: Dr. Matoba?

8 DR. MATOBA: In looking at your protocol, I
9 didn't see any reference to pupil sizes as either exclusion
10 criteria or inclusion criteria. Is that taken into account
11 when you were entering these patients?

12 DR. PETTIT: I have to defer.

13 DR. MATOBA: And the reason I ask is because
14 the 50 patients who had the conventional ablation were done
15 earlier in the study than the Custom ablation patients in
16 general, correct? And I wondered if there might be a
17 difference in the average pupil size between the two groups
18 and that may affect patient satisfaction or other visual
19 tests.

20 DR. PETTIT: Okay. The entry criteria did not
21 change over the course of the study.

22 DR. MATOBA: If you didn't take pupil size into
23 account, then how would you know what they may have been in
24 those patients that were entered?

25 DR. PETTIT: We can look at pupil size. Do you

1 want to speak? Okay. Let me make sure I have the question
2 down exactly right. For this conventional comparative
3 group, did we control for pupil size or did we have an
4 analysis baseline?

5 DR. MATOBA: I'm sorry. Go ahead.

6 DR. PETTIT: Well, I want to make sure I get it
7 right. Did we analyze the data based on pupil size to see
8 if there was any statistically significant difference in
9 pupil size between those eyes and the wavefront group?

10 DR. MATOBA: Well, my main question is could
11 there have been a significant difference in the average
12 pupil size between the patients who had Custom ablation
13 versus people who had conventional ablation?

14 DR. PETTIT: Okay. So you're concerned that
15 there might be a significant difference in pupil size
16 between the two groups?

17 DR. DURRIE: I can just comment from an
18 investigator standpoint, is the inclusion criteria were the
19 same throughout the study, and the conventional eyes were
20 done fairly contemporary because initially we did a group
21 with the same algorithm. We did a group that had one eye
22 with conventional and the one eye with Custom. The reason
23 that the number is small is when we drop it down just to
24 spherical cohort, then it gets us down to only 50 eyes. So
25 these were done at relatively the same time period. It's

1 not like one group was done three years ago and one group
2 was done a year ago, and the inclusion criteria regarding
3 pupil size was the same.

4 I think it is interesting, and I think it'd be
5 a good thing to look at, is to actually analyze, since the
6 aberrometer itself records pupil size, so that data is
7 available, and it would be a good thing to continue to look
8 at. I think with all these things, we're coming up with
9 new things we can look at because now we have a digital
10 instrument that can actually give us data that we haven't
11 had before.

12 DR. WEISS: I'm going to ask members of the
13 panel, for the purpose of the transcription, if they can
14 identify themselves before they speak. The sponsor.
15 Excuse me.

16 Dr. Bandeen-Roche?

17 DR. BANDEEN-ROCHE: Dr. Bandeen-Roche.

18 My question follows up on Dr. Weiss's first
19 question and Dr. Huang's question. I think it's coming up
20 repeatedly. If you could just clarify very explicitly this
21 conventional cohort, you know, so they were selected early
22 on in the study. To what extent was randomization involved
23 in their selection versus the Custom eyes? Were they
24 treated by all of the same physicians who treated the
25 Custom eyes? Were there any differences with respect to

1 this cohort that could be expected to influence results,
2 including practice effects, as the study went on? I know
3 there were some differences in temperature and humidity.
4 Pupil size has been raised, et cetera.

5 DR. PETTIT: Okay. I'm going to ask Dr.
6 Christy Stevens, our Clinical Affairs Director, to come
7 forward to talk about the inclusion criteria.

8 MS. THORNTON: Excuse me. I'd like to
9 emphasize to the sponsor, the transcriber and summary
10 writer have indicated that they would like very much to
11 have help from the sponsor group and give your
12 identification before you speak.

13 DR. STEVENS: Christy Stevens, Alcon.
14 The study started with a contralateral design
15 with one eye Custom, one eye conventional, and it was
16 randomized as to which eye would be the Custom treatment.

17 DR. WEISS: I wonder. Would you be able to get
18 a little closer to the microphone and speak a little
19 slower? Thank you very much.

20 DR. STEVENS: Do you need me to repeat what I
21 just said?

22 DR. WEISS: I think it would be best. Yes,
23 please.

24 DR. STEVENS: Okay. It was a contralateral
25 study design when it began. One eye received Custom, one

1 eye received conventional, and it was randomized as to
2 which eye received the Custom treatment.

3 When we started the study, we modified the
4 Custom algorithm over the course of the early part of the
5 study. So in our Custom cohort that you've seen presented,
6 it contains the final algorithm, only contains the last
7 algorithm, but we included all conventional eyes that were
8 in the beginning of the study because they were all treated
9 the same way.

10 DR. BANDEEN-ROCHE: And I thought I heard a
11 comment that only 19 of the eyes were treated one eye one
12 way and one eye the other. So who were the others?

13 DR. STEVENS: There were 50 spherical
14 conventional eyes total, of which 19 had an algorithm or
15 current algorithm, the last algorithm in Custom, that were
16 also a spherical eye, and so yes, they were treated by the
17 same centers with the same study protocol.

18 DR. BANDEEN-ROCHE: Thank you.

19 DR. WEISS: Thank you.

20 Mr. McCarley?

21 MR. MCCARLEY: I had just a question about what
22 effects have you considered or do you expect for a patient
23 following cataract surgery? In other words, the mean on
24 these patients was 36 years, I think, and assuming they're
25 all phakic eyes, what happens when they're 70 or 75 and

1 have cataract surgery? Would they be expected to go back
2 and have a reablation to address what had been corrected in
3 their system, now that one of the components is missing?

4 DR. PETTIT: Yes, that's a good question, and
5 the honest answer is we don't have any clinical data on
6 patients that meet that criteria.

7 There is evidence in a young patient that the
8 corneal aberrations are somewhat balancing compared to the
9 internal aberrations and that, you know, obviously the
10 situation's going to change very significantly when you go
11 in and do cataract surgery. I think it's important to
12 mention, though, that with our treatment, we are keeping
13 the aberration magnitude comparable to what it was before
14 surgery. We're not grossly changing the magnitude of the
15 aberrations that were present in the eye beforehand. We're
16 keeping them more like they were before treatment. So we
17 don't anticipate that we're suddenly going to have all
18 these new problems when the patients come back for cataract
19 surgery.

20 Now, to get the best possible optical quality
21 after cataract surgery, potentially sure, they might
22 benefit by some kind of customized correction on top of
23 that, but we don't have any clinical data that that's
24 actually true.

25 DR. WEISS: Dr. Owsley and then Dr. Huang.

1 DR. OWSLEY: Thank you.

2 I just wanted to make sure I'm understanding
3 your low contrast acuity data properly. It appears that
4 whether we look at the sphere analysis or the all eyes
5 analysis, 20 percent or one in five patients experienced a
6 loss in low contrast acuity. I know that's different,
7 lower than the rate in the conventional surgery, but I just
8 want to make sure I understand. Twenty percent of the
9 patients, one in five, experienced at least a one line or
10 greater loss?

11 DR. PETTIT: That is correct.

12 DR. OWSLEY: Thank you.

13 DR. WEISS: Dr. Huang, please.

14 DR. HUANG: I just want to ask. Since we set
15 out to try to correct the higher-order aberrations by this
16 application, but the end result shows that there were
17 general increase of the higher-order aberrations, and I
18 don't know if the clinicians or the sponsor have any kind
19 of comment regarding the outcome.

20 DR. PETTIT: It is true that even after our
21 wavefront-guided surgery, that the higher-order aberrations
22 are generally higher. They're higher by an amount that's
23 significantly much less than what we get with our
24 conventional surgery. We believe that's beneficial to the
25 patient.

1 Our theoretical endpoint is to make them zero,
2 and we clearly are not achieving that yet, but by shooting
3 for that as the theoretical target, we are limiting what
4 happens to them and that's where we are with the current
5 state of the technology.

6 DR. HUANG: But my point is, instead of
7 reducing, now we are actually increasing. So what's the
8 future direction in that regard?

9 DR. DURRIE: Dan Durrie.

10 From a clinical standpoint, this is a step
11 along the way because before we weren't even measuring the
12 patient's preoperative aberrations, other than sphere and
13 cylinder. Now we're finding other things that we find
14 clinically significant in the population now that the
15 aberrometer can measure.

16 As George said in his presentation about the
17 progress, we found out then that there was some surgically-
18 induced aberrations and some of them were predictable. I
19 think as time goes on, we will learn more about the
20 surgically-induced aberrations and then may have to make
21 compensations. I think it's going to be important for all
22 of us to start thinking about how are we going to
23 accomplish that from a regulatory standpoint when you come
24 up with the next new iteration, so it isn't so onerous that
25 the companies cannot pursue that, and it isn't too onerous

1 from a regulatory standpoint, and I think it's something
2 that I know that you're having a meeting tomorrow to talk
3 about phakic eye welds, but I'd certainly like to have us
4 continue to have a discussion between the sponsors and the
5 clinicians and the agency about once we get better, what do
6 we do then? Because I think this is a step along the way,
7 but we still would like to make that zero, and we're going
8 to have to continue to evaluate data in order to make that
9 happen.

10 DR. WEISS: Dr. Burns, and then Dr. Bandeen-
11 Roche.

12 DR. BURNS: Yes. Your sample had a very low
13 percentage of Asians in it, yet it's a high-refractive
14 error group, and I just wondered if you had a comment.

15 DR. PETTIT: I think the race distribution in
16 the study was comparable to that we've seen in prior
17 studies. Do we have any further comment? I mean, there
18 was no attempt to include or specifically recruit certain
19 patient populations or not. This is just the patients that
20 came in, were interested in being in the trial and met all
21 of our inclusion criteria.

22 DR. WEISS: Dr. Roche?

23 DR. BANDEEN-ROCHE: Dr. Bandeen-Roche.

24 I have a question that goes to the fact that
25 relatively few sites participated in the study. So you

1 provided the site distribution of all eyes and it ranged
2 from 36 percent in the provider who did the most to about
3 10 percent in the provider who did the least.

4 Do you have the same distribution for the
5 spherical eyes and also the distribution of the
6 conventional, the 50 conventional eyes by site?

7 DR. PETTIT: We'll see what we can dig up in
8 that regard.

9 DR. BANDEEN-ROCHE: And finally, can you
10 describe how the sites were selected and what training,
11 just very briefly, the extent of training that the
12 physicians received?

13 DR. PETTIT: The site selection criteria, we
14 obviously were interested in trying to get innovators in
15 this field, high-volume/high-profile refractive surgeons
16 that were knowledgeable about new advances in technology.
17 I don't know that there was anything beyond that. It was
18 as simple as these doctors seemed to be very well qualified
19 and were interested in participating in the study, and we
20 wanted to work with them.

21 The training that they received, they obviously
22 received some training in how to use the wavefront device.
23 The treatment aspects are very similar to what they were
24 already using for their conventional LADARVision
25 treatments. There was a slight difference in the fact that

1 they had to mark all patients before treatment as opposed
2 to just spherical patients in their conventional surgery,
3 those little eye marks they put on the eye, and then the
4 reticle that came up during treatment, the software image
5 projected into the LADARVision tracked image screen was
6 slightly different, but the other aspects of the treatment
7 were really identical to what they'd seen before.

8 We did spend some time going through them, the
9 meaning of the wavefront measurements and, you know, when
10 their clinician brought them data, what did that mean for
11 that patient in terms of relative to the foropter, for
12 example?

13 DR. WEISS: Dr. Grimmiett, and then Dr. Matoba.

14 DR. GRIMMETT: Dr. Michael Grimmiett.

15 I just had an observation and would like to
16 hear if you have a comment. You may have none. There may
17 be no answer.

18 I found it curious that despite a very
19 comprehensive analysis and sophisticated technology, that
20 the patients that were unsatisfied or extremely unsatisfied
21 approximated 9 percent. It's notable that the PERK study
22 by comparison, using bear skins and stone knives, had an
23 11-percent dissatisfaction rate, and I found it curious
24 that one in 10 patients are unsatisfied, despite a
25 phenomenal amount of technology and analysis, and I would

1 like to commend you on a superb analysis and presentation.

2 Do you have a comment why it's still one in 10
3 despite all the sophisticated technology or is there no
4 answer to that, sir?

5 DR. PETTIT: No. Well, I don't know
6 everything. I didn't personally speak to these patients.
7 I think one factor was that they tended to be a little bit
8 undercorrected. The patients that were undercorrected on
9 average were less satisfied than the patients that were
10 right on, and again there's no latitude for the surgeon
11 trying to optimize the refractive outcome. We wanted
12 everything done exactly the same way and that led to a
13 slight undercorrection, and the patients where that
14 undercorrection was more than the mean, they ended up more
15 myopic than the mean, they tended to be less satisfied.

16 I think, you know, in all honesty, in addition,
17 there's a lot of hype surrounding this procedure, and I
18 think their expectation levels in some cases was pretty
19 high, but, you know, that's not scientific. That's just an
20 opinion.

21 DR. HAKIM: If I could just add to Dr. Pettit's
22 comments, I mean, I agree that --

23 DR. WEISS: Could you introduce yourself,
24 please?

25 DR. HAKIM: Omar Hakim. Sorry. Omar Hakim.

1 I just want to add and really reinforce
2 George's comments, you know, about the undercorrection
3 aspect of this. We weren't able to make any adjustments as
4 we normally do when we do surgery, and clearly there was a
5 difference in the patients who were satisfied versus
6 dissatisfied with their surgery, based on their residual
7 refractive error, and the expectation issue, I think, you
8 know, is a very big issue. I literally had patients coming
9 back who were now seeing 20/16 uncorrected acuity and
10 wanted enhancements. So their expectations of the surgery
11 clearly were raised as well as in the minds of their own
12 physicians who had referred them. They were talking about
13 supervision and the Popular Science article last March, you
14 know, talking about the ability to give people, you know,
15 better than 20/10 or 20/8 or 20/6 vision.

16 Clearly, what we really want to do is avoid
17 problems, you know, like Ron Link talked about this
18 morning, is try to create better quality of vision, and as
19 Dan Durrie was saying, this is really a process in
20 evolution, but if I could have my surgery done again today
21 and avoid the induction of these higher-order aberrations
22 that we create whenever we do conventional surgery, that's
23 what I would choose for myself and all my patients.

24 DR. PETTIT: Just to follow up a little bit on
25 an earlier question. This is George Pettit from Alcon.

1 I think Dr. Matoba asked the question about the
2 pupil sizes and were the patients informed. Given this
3 high-level expectation, it's important to note that the
4 optical zone was 6.5mm and we informed all patients
5 considering being in the trial that if their natural pupil
6 was larger than 6.5mm, even with this new technology, there
7 was a potential risk for them to have night vision
8 symptoms. So we tried to bring their expectations more in
9 line.

10 DR. WEISS: Dr. Grimmer has a follow-up
11 question.

12 DR. GRIMMETT: I have just a simple operational
13 question.

14 Is this software that's going to be retrofitted
15 to existing product base, the LADARVision 4000s out there?
16 Does this require a whole investment in brand-new
17 technology?

18 DR. PETTIT: No, from the LADARVision side,
19 it's a simple software upgrade.

20 DR. WEISS: Dr. Matoba?

21 DR. MATOBA: Alice Matoba.

22 I actually was going to refer to the same chart
23 for the table on Patient Satisfaction and also the previous
24 page on Patient Symptoms. Do you have these same data for
25 the people who were treated with the conventional laser?

1 DR. PETTIT: So the question is do we have the
2 patient satisfaction questionnaire-type data for the
3 patients treated conventionally?

4 DR. MATOBA: And also symptoms at six months.

5 DR. PETTIT: And do we have data on the
6 conventional patient symptoms at six months? We'll see
7 what we have in that regard.

8 DR. WEISS: Dr. Bradley was next, then Mr.
9 McCarley, Dr. Swanson, and Dr. Owsley.

10 DR. BRADLEY: Dr. Bradley.

11 I'm just curious about the apparent huge
12 discrepancy between what Ron Link presented earlier today
13 and the data on the symptoms presented by the sponsor. For
14 example, Ron Link indicated that dryness and double vision
15 are huge problems, and I think we have a couple of other
16 people indicating that, and I look at the data you just
17 presented on dryness where we have slightly more patients
18 indicating worse dryness than those indicating better, and
19 we have a very small number indicating increased double
20 vision, about the same as those who are indicating
21 decreased double vision.

22 So from the sponsor's dataset, it appears that
23 we don't have this very large and disturbing incidence of
24 dry eye and optical problems, such as double vision,
25 whereas Ron Link and a couple of the other presenters

1 indicated that these are very serious problems, and I
2 wondered therefore if we could clarify perhaps some
3 inclusion criteria from the sponsor because one wonders if
4 Mr. Link's dataset is rather biased to those who have the
5 problems and somehow you have been able to effectively
6 filter these people out of your datasets. Yours are biased
7 the other way.

8 I think it's very important to get a sense of
9 that, particularly for those people who are going to
10 utilize this technology, and if you have effectively
11 avoided these problems by your patient selection criteria,
12 then this clearly should be included in the final labeling
13 for this device.

14 DR. PETTIT: Dr. Durrie, would you like to
15 comment?

16 DR. DURRIE: Yes. I'd really like to address
17 that, and this is Dan Durrie.

18 Ron's website, which is where he gets his data,
19 are for people who've had surgery and by its own definition
20 and its goal, it's for people who have problems with
21 refractive surgery, and I really appreciate the work that
22 he's done on helping us define of those patients who have
23 problems with refractive surgery, what are their problems,
24 and obviously 25 percent of those problems are persistent
25 dry eyes. But this is a very selected group not only

1 that's had refractive surgery but is self-defined that they
2 have problems and they're logging into the website. So I
3 think that's a defined group on that side.

4 On our side, I think that the only criteria
5 that I think is significant from my clinical experience is
6 the average age of this group was 38 years old, and we know
7 that the patients, if I did LASIK on an average age of 55-
8 year-olds, they'd have more problems with dry eyes. So I
9 think that if there's a self-selection in this, there
10 certainly wasn't anything in the screening from the
11 standpoint of we had healthy eyes, there wasn't any tear
12 film screening or any special testing, but I think that you
13 do have a healthy eye group that's screened for a clinical
14 trial that certainly is a healthy eye group and on the
15 other side, in the surgical eyes group, you have the group
16 that basically is having problems, and I think both those
17 datasets are important.

18 I'd like to also, because I did look this up
19 during the discussion, is if we take the 426 eyes that are
20 available for analysis and run that same grid, that total
21 grid of symptoms, the numbers are essentially identical.
22 So here, you have -- which was requested really by the
23 public presenters -- a very good dataset with a 100-percent
24 follow-up on 426 eyes that gives you an array of symptoms
25 on how many patients were the same, better or worse, and I

1 think that could be a dataset for labeling that could give
2 you some good information with peer data because it does
3 have 100-percent follow-up and it was done under a
4 controlled fashion.

5 DR. WEISS: Dr. Bradley has a follow-up.

6 DR. BRADLEY: From your reply, you seem to be
7 saying you've done nothing special to avoid these dry eye
8 or night vision, double vision problems. Did I understand
9 that correctly?

10 DR. DURRIE: In patient selection.

11 DR. BRADLEY: Second question, and this is
12 really to --

13 DR. PETTIT: Dr. Bradley?

14 DR. BRADLEY: Yes?

15 DR. PETTIT: Could I just perhaps follow up?
16 This is again fairly much just conjecture on my part, but
17 the patients that participated in this study had to be
18 willing to come back for many, many follow-up visits. So
19 perhaps we, without attempting to, screened for a more
20 educated or, you know, patients that really wanted this
21 type of procedure and knew what the risks were ahead of
22 time, I don't know, but they weren't just your standard
23 patient coming in off the street that weren't going to have
24 to go through all these tests for six months.

25 DR. BRINT: I think my comment is similar to