

1 some surgeons will report just the cataracts that are
2 creating visual problems for the patient, and some
3 opacities that are not decreasing the visual acuity would
4 not be reported.

5 Also, there is such a variation in the age of
6 the patients included, and also such a variation in the
7 follow-up considered, and in some studies they would say
8 that the follow-up is from three months to two years, but
9 in fact maybe just two patients had the follow-up of two
10 years and the great majority were followed up for six
11 months only. And in some studies -- for example, for the
12 ICL -- during the same study different models of the same
13 design were implanted, and you know with different vaulting
14 characteristics, so the effect for cataract formation would
15 not be the same.

16 So there is a great need of standardization of
17 these studies evaluating cataract formation. First of all,
18 the parameters used for the YAG iridotomies should be
19 described, because eventually this is cataractogenic. All
20 trauma to the anterior capsule during this surgery should
21 be noted for future reference. A follow-up period should
22 be at least two years because in fact, according to the
23 literature, the majority of cataracts appear between one
24 and two years after the procedure. And of course, we have
25 to evaluate very nicely the relationship of the phakic lens

1 with anatomic structures because the contact is one of the
2 most important factors.

3 Also, of course, we need a very accurate method
4 to choose the IOL size. To use just the white to white to
5 choose the overall size of the lens that's going to be
6 implanted is not accurate at all, and that's why we're
7 having so many complications, and this could really be
8 avoided.

9 There is a need of evaluation of subclinical
10 inflammation with laser flare meters because in some cases,
11 there was a very good vaulting lens and a cataract appeared
12 anyway.

13 There is a need of evaluation of explanted
14 phakic lenses with histopathological analysis of adjacent
15 tissues, but if you ask me right now if just with
16 histopathologic studies alone we will be able to
17 differentiate cataracts caused by the surgeon and cataracts
18 caused by the lenses, the answer is that we don't know yet.
19 We have just one specimen and we would like to look for
20 more specimens to have an impression about that.

21 And of course, we need to describe the
22 evolution of the anterior subcapsular opacity.

23 So let's talk about the possibility of a
24 classification for cataract formation after phakic IOL
25 implantation. So as I mentioned, some surgeons would say

1 that all the opacities they saw are peripheral and they are
2 non-progressive, but we have this paper indicating that the
3 peripheral superior opacity progressed, involving the
4 optical zone. So if you have a classification, we should
5 maybe classify the opacity in each visit to have an
6 impression about the progression of the problem.

7 So as you know, there are systems for the
8 classification of cataracts and they are all based on
9 standard retroillumination photographs and the total area
10 of the opacity. These are three well-known systems, the
11 LOCS system, Wilmer system, and Oxford system.

12 Here, we have some pictures showing how to
13 classify nuclear, cortical, and posterior subcapsular
14 cataracts with the LOCS system, all based on these standard
15 photographs. This could eventually be applied to cataract
16 formation after phakic IOL implantation. This is how we
17 grade the cortical opacities according to the Wilmer system
18 and this is the way we score anterior subcapsular and
19 posterior subcapsular opacities according to the Oxford
20 system.

21 Also, there are very sophisticated systems
22 combining high-resolution digital retroillumination imaging
23 with image analysis systems, allowing objective and
24 quantitative measurement, for example, of posterior
25 capsular opacification, which eventually would be very

1 useful in this clinical situation.

2 So in a way, if there is a way to do a
3 classification, this classification should indicate the
4 location of the opacity. For example, peripheral or
5 paracentral or central opacity. Also, maybe there is a
6 possibility to have an index for the intensity of the
7 opacity, and of course, we would have to score the area of
8 the opacity, and by doing that in each visit, we will have
9 an impression about the evolution and the progression of
10 the opacity and we would really understand better the
11 phenomenon.

12 Thank you very much for your attention. Thank
13 you very much again for the opportunity.

14 DR. WEISS: Thank you very much for an
15 excellent presentation. Would you be able to take a seat
16 at the table, and we'll open up to the panel for some
17 directed questions.

18 Dr. Bandeen-Roche, then Dr. Matoba, and then
19 we'll continue around. Dr. Bandeen-Roche?

20 DR. BANDEEN-ROCHE: Yes, I want to thank you
21 for your presentation -- very clear -- and your careful
22 recommendations about the data to collect involving the
23 surgery I think was great.

24 Just one brief question. You talked about the
25 importance of tracking the evolution of the cataract and

1 discussed grading. What implications do you think there
2 are for the frequency of evaluation and do you think that
3 that should be by passive surveillance or by active
4 surveillance? Just elaborate a little bit more on how to
5 track the evolution of cataract.

6 DR. WERNER: I don't have a precise idea about
7 when all these gradings should be done. Of course,
8 immediately after the postoperative period, within one
9 month to anything that's very fast developed, and maybe six
10 months, one year, two years, because this is what we see in
11 the literature.

12 But this is not only to just have the score.
13 It's also for us to understand the phenomenon because still
14 there are many surgeons who believe that the cataract is
15 really not a problem because they have just opacities in
16 the periphery that never progress, and we need to
17 understand if they are not really progressing or he is just
18 not observing.

19 DR. BANDEEN-ROCHE: Thank you.

20 DR. WEISS: Dr. Matoba?

21 DR. MATOBA: In regard to the development of
22 cataracts in areas of contact between the IOL and the
23 crystalline lens, how much do you think the lens material
24 or the nature of the lens material contributes or is it
25 mostly, you think, a mechanical effect?

1 DR. WERNER: Well, I don't know if I can answer
2 this question because we only have results about the ICL,
3 which we know the material. There is another lens made of
4 silicone and I'm not aware of their results, so we cannot
5 really compare if there is a material effect.

6 Also, as I mentioned, there are papers showing
7 that cataract is formed only in the area of contact. Other
8 papers would say that it was formed in a different area,
9 but maybe the follow-up was not enough. So there are still
10 many questions about that.

11 DR. WEISS: Dr. Bradley, Dr. Huang, and then
12 Dr. Mathers.

13 DR. BRADLEY: Just a clarification on the
14 peripheral cataracts you described. How peripheral are
15 they and, for example, would they become visually
16 significant under nighttime viewing conditions where the
17 pupil would dilate?

18 DR. WERNER: Well, when you look at the
19 literature, it is really not described, and sometimes, in
20 very early papers, they would describe some peripheral
21 opacities that would cause some glare in light conditions
22 of evening or something like that. But talking to
23 surgeons, they would say that the peripheral opacities
24 which are barely visible in the pupil dilation, they would
25 not cause any problem.

1 DR. WEISS: Dr. Huang?

2 DR. HUANG: Two questions. The first question
3 was regarding your earlier presentation, you seem to have
4 implied there were two types of cells that were induced in
5 the two different locations of the cataract. One is A-
6 cells induced in the anterior capsule, and then E-cells
7 induced in the posterior capsule. Is there any vital stain
8 that can help you to distinguish what type of cells are
9 responsible for the evolution of this cataract?

10 DR. WERNER: When you perform histopathological
11 studies, in fact what you see is that both cells are always
12 involved in everything, but there is always a predominate
13 type. For example, even for PESU, you have a fibrotic form
14 of PESU and you have a firm form of PESU, and when you
15 perform normal stains, you can see even morphologically
16 they are very different because E-cells always have the
17 tendency to be bloated, and the other are elongated
18 fibrotic cells, fibrotic-like cells.

19 DR. HUANG: In the LOCS III grading system,
20 it's really a numerical system and there is a highly
21 individual variation. Do you have any suggestion how to
22 standardize if that system were to be chosen for the
23 cataract characterization? And also, I believe the LOCS
24 III does not have any geographical information about a
25 cataract, and so do you have any suggestion how to modify

1 the system?

2 DR. WERNER: I believe there should be a
3 geographic thing because it's very important for the areas
4 of contact who are outside the contact of the lens. So in
5 this case, it's very important.

6 But with regards to your first question, I
7 think we should start by collecting many pictures from
8 surgeons to create standard pictures, as they have in the
9 LOCS system, and we have some, but we need more pictures to
10 have really different grading if we use such a system based
11 on standard photographs.

12 DR. WEISS: Dr. Mathers?

13 DR. MATHERS: Do you think that the
14 photographs, the retroillumination of the opacification,
15 can the photographs detect smaller anterior subcapsular
16 cataract formation than the slit lamp can detect it? What
17 do you think is actually the finest, highest resolving
18 method?

19 DR. WERNER: We have experience with these
20 systems based on retroillumination photographs for
21 posterior capsular opacification, and this is the best we
22 can have for that.

23 DR. MATHERS: Do you think that's higher
24 resolution than the human eye achieves with the slit lamp?

25 DR. WERNER: Eventually, yes.

1 DR. MATHERS: Yes. And have you any experience
2 using confocal microscopy, which could actually focus onto
3 the anterior capsule and have even higher resolution? Do
4 you think that that's possible?

5 DR. WERNER: Well, my own experience with
6 confocal microscopy regards the cornea, and I know that
7 there are some objectives which you could switch in some
8 devices and have an imaging of the anterior surface of the
9 crystalline lens. I have no experience with that and I
10 don't know if there is any data available published about
11 that.

12 DR. MATHERS: Sizing is clearly an important
13 process here. Do you think that the high-resolution
14 ultrasound will give the best sizing data and can that be
15 used clinically to determine which size to put in?

16 DR. WERNER: We are evaluating this right now
17 with these cadaver eyes. The results have been very
18 interesting. So I don't know exactly the status of
19 development of the technique, when it's going to be
20 available -- maybe this year -- but apparently it's the
21 best we can get for the moment.

22 DR. MATHERS: Do you think that the issue of
23 visual significance could be assessed best with glare
24 testing or what would you suggest as the most significant,
25 highest-resolving method to test the small amounts of

1 visual impairment you might get from an early cataract?

2 DR. WERNER: Yes, it is a very good question,
3 because in some cases discussed with the surgeon, the
4 lenses were explanted because of glare, and not really
5 decreasing visual acuity. So glare is very important.

6 DR. MATHERS: You think it's better than
7 contrast sensitivity testing as it's normally performed?

8 DR. WERNER: Maybe both.

9 DR. MATHERS: Just your opinion.

10 DR. WERNER: Maybe both should be associated in
11 this case, yes.

12 DR. WEISS: Mr. McCarley, did you have a
13 question as well?

14 MR. MCCARLEY: One question just quickly. In
15 your experience, and maybe one of the clinical
16 ophthalmologists can answer this maybe even better, a cell
17 flare meter is used to determine subclinical inflammation.
18 Is that different in a posterior chamber lens than it would
19 be, for instance, in an anterior chamber lens?

20 DR. WERNER: Well, what we saw in the
21 literature is that there are also increased values of flare
22 cell meter with anterior chamber lenses, and apparently the
23 values are even higher and they also stabilize above the
24 preoperative values, contrary to cataract surgery, where
25 you have higher values, but these have a tendency to come

1 back to preoperative values after one year or so.

2 DR. WEISS: If there are no other questions, I
3 want to thank you, Dr. Werner, for your excellent
4 presentation.

5 We can move on to the open public hearing
6 session. Is there anyone who wanted to make any
7 statements?

8 (No response.)

9 DR. WEISS: If not, I'm going to put a question
10 to the panel. If we have perhaps a 15-minute coffee break
11 and skip lunch, we might be able to catch earlier flights,
12 as I know is in the interest of some. Are any of you
13 interested in doing that, taking a 15-minute break, rather
14 than -- so we have two hands up for Dr. Bullimore and Dr.
15 Matoba, and Dr. Mathers for sure. I would say that passes
16 without a formal vote.

17 So we'll take a 15-minute coffee break and
18 we'll see you back here in 15 minutes.

19 (Recess.)

20 DR. WEISS: We will now start the FDA
21 presentation and Donna Lochner will introduce the questions
22 for panel discussion.

23 MS. LOCHNER: Yes. I'm just going to go
24 through the questions and give a little bit of background
25 to where we were coming from with each question.

1 We're going to begin today with the endothelial
2 cell density study. After I've stepped through this
3 question, I'll turn the floor over to Dr. Grimmert, and
4 then the panel will discuss the endothelial cell issue
5 before going on to the next questions.

6 First, "Please comment upon the inclusion
7 criteria recommendations found in Table 1." This topic is
8 still being actively discussed, and particularly with the
9 ANSI Standards Committee, and so we believe any comments
10 will be very timely.

11 Table 1, which is just shown right here,
12 provides the recommended minimum endothelial cell
13 densities, and these values for minimum endothelial cell
14 density are generally being used in current U.S. phakic IOL
15 studies. These values were determined over the course of
16 several meetings over the years with input from FDA,
17 industry, and ophthalmologists that attend these standards
18 meetings. Allow me to hopefully clarify how the minimum
19 densities per age category were determined.

20 This slide is included as Attachment B in the
21 handout. First, the approximate initial cell density for a
22 21-year-old, as shown in the second column in this table,
23 was taken from the 1997 Moller-Pedersen article and the
24 citation for this article is provided in the handout, but
25 not on this slide.

1 For the 35- and 46-age categories, the cell
2 density at time of implant was approximated by assuming .6
3 percent yearly cell loss due to normal aging, with the .6
4 percent figure taken from the 1997 Bourne article, as
5 referenced earlier by Drs. Edelhauser and McCarey. This
6 was done to provide a check of whether the minimum
7 inclusion criteria per age group were reasonable.

8 The third column, the estimated rate of cell
9 loss per year, represents potential rates of loss due to
10 the phakic IOL. In other words, 1.5 and 2 percent assumed
11 loss from the phakic IOLs were used as examples to then
12 calculate the age when the cell density would be less than
13 1,200 cells per millimeter squared and less than 1,000.
14 These ages, shown in the fourth and fifth columns, assume a
15 surgical loss of 10 percent and compound the 1.5 and 2
16 percent loss annually.

17 Finally, in order to determine the minimum cell
18 density inclusion criteria, we looked at the starting
19 densities that would ensure greater than 1,000 cells at age
20 70 for the 21- to 25-age range, and at 75 for 26 and older.

21 So this table verifies that the minimum
22 inclusion criteria, as shown on Table 1, would be
23 sufficient in a worst case situation to allow for adequate
24 cell density for the health of the cornea for roughly the
25 life of the patient, assuming a 2 percent annual loss from

1 the phakic IOL and that patients would have at least 1,000
2 per millimeter squared at ages 70 to 75.

3 As you can see, there are various assumptions
4 inherent in the inclusion criteria and because of the
5 iterations this has undergone because of the committee
6 process, we are very much interested today in the panel's
7 comments on this inclusion criteria.

8 Our statistical calculations suggest that 200
9 subjects should be sufficient to detect a 2 percent loss
10 using measurements at multiple visits in order to establish
11 linearity of the loss. The measurements, as currently
12 proposed, would taken at the three- or six-month visit, the
13 12-, 24, and 36-month visits.

14 Further, although the statistics suggest that
15 200 subjects would be sufficient, we recommend that
16 specular microscopy be performed on all subjects enrolled
17 in the study to ensure that 200 analyzable photographs are
18 obtained.

19 Last, we recommend that multiple images be
20 captured at each visit and the mean endothelial density
21 from those multiple images be used in the analysis.

22 We are asking for panel comments on these
23 criteria as well.

24 Now, I'd like to turn the floor over to Dr.
25 Grimmatt for his review.

1 DR. GRIMMETT: Thank you, Donna.

2 This is Michael Grimmett. I've prepared some
3 comments in outline form which at least the panel members
4 should have on their table. It's a 10-page outline, and I
5 promise to go very quickly through it. There are four
6 tables as well.

7 Regarding the questions, I think I'll go
8 through my outline first. I was doing the bulk of this
9 thought process and review prior to having the questions.
10 I was doing the work, I'd gone on vacation out in Santa Fe,
11 and I'm relieved and grateful that in the open public
12 hearing session, my wife did not comment on the timing of
13 that review.

14 (Laughter.)

15 DR. GRIMMETT: Notwithstanding that, I think
16 I'll do this review first, and then we'll come back to the
17 questions if some of the issues are not resolved.

18 Just in general, and some of this has been
19 alluded to by the previous speakers, the peer-reviewed
20 literature on phakic IOLs has numerous limitations, and
21 it's important to recognize that when reviewing any data
22 that's reported in the literature. Mostly, the data are
23 retrospective in design, they're non-randomized case
24 series, they have extremely low numbers of eyes reported,
25 there is poor accountability for the longer follow-up

1 intervals, and in general morphometric endothelial analyses
2 are not generally reported. So coefficient of variation
3 and percent hexagonality are not there.

4 Additionally, my review is not exhaustive or
5 comprehensive.

6 I think it's instructive to look at phakic IOL
7 types because we can separate them into three different
8 animals that have different implications for the corneal
9 endothelium. I've separated them into anterior chamber
10 type and posterior chamber type. Of the anterior chamber
11 type, there are two, angle-supported and iris-fixated.

12 Looking at the literature on the angle-
13 supported lens, I'll just give you a smattering of some
14 articles to go over what are reported endothelial cell
15 losses and what's known about design parameters that would
16 impact our recommendations regarding future phakic IOL
17 studies.

18 To start with, angle-supported lenses, a first-
19 generation lens, the Baikoff ZB lens, had a distance
20 between the IOL edge and the endothelium of only 1.16
21 millimeters, and that was the key factor. It was
22 determined that there was a high endothelial cell loss
23 secondary to excessive contact between the IOL optic edge
24 and the endothelium. One report by Jimenez-Alfaro reported
25 a 16 to 18.8 percent loss at one year and a 20 to 28

1 percent loss at two years. Another report reported up to
2 56 percent loss with a shorter follow-up. Of the second
3 report, 37.5 percent had morphologic or cell density
4 changes.

5 Quotes in these articles include "We have
6 stopped doing this surgical procedure" and "Further
7 implantation of this IOL is unacceptable."

8 The second-generation lens, the Baikoff ZB5M
9 lens, was manufactured until 1997. The major difference
10 here is that they increased the distance between the IOL
11 edge and the endothelium, increasing that to 1.56
12 millimeters. Endothelial cell loss in one study was
13 reported at 4.5 to 5 percent cell loss at one year, 5.6 to
14 6.8 at two years, and 5.5 to 7.5 loss at three years.

15 A separate study by Perez -- forgive me for my
16 pronunciation -- Perez-Santonja found an endothelial cell
17 loss of 12.33 percent at one year and remaining relatively
18 stable at two years.

19 Then the larger study was recently reported in
20 Ophthalmology found that at one year there was a 5.53
21 percent loss, and interestingly, after year 2, while the
22 numbers fell down significantly in terms of total eyes
23 examined, the overall loss approximated normal aging
24 losses. That is, preop to year one was at 5.5 percent
25 loss; year 1 to year 2, 1.37 percent loss; year 2 to 3, .72

1 percent loss; year 3 to 4, there was a .28 percent loss;
2 year 4 to 5, there was a .55 percent loss; year 5 to 6, .37
3 percent; and year 6 to 7, .56 percent.

4 So based on this study, and the reason I find
5 this instructive, is that they give us some information
6 regarding what is a reasonable duration to follow these
7 particular lenses before they come to the panel. If we
8 take this data, which is one of the largest subsets that I
9 could find in the published literature, in year 3 they
10 still had 157 eyes, albeit they're mixing lens times.
11 They're either started or stabilized, and they found
12 stabilization from the two- to three-year period, and you
13 can tell three years is the appropriate duration of a
14 study.

15 A fourth-generation lens, ZSAL-4 by Morcher,
16 had a distance from the IOL edge to the peripheral cornea
17 of 1.65 millimeters and they reported an endothelial cell
18 loss of 3.5 percent at one year and 4.2 percent at two
19 years, albeit the numbers are low. There are only 18 eyes.
20 They commented that there is a fifth-generation lens, but I
21 could not locate it in the published data regarding this
22 product.

23 Another anterior chamber angle-supported lens
24 in the literature is the Nuvita lens by Bausch & Lomb, and
25 it's a single-piece PMMA IOL. In one study of 21 eyes by

1 Allemann in Ophthalmology, they reported a 10 percent cell
2 loss in the first year, another 4.3 percent in the second
3 year, and, given these large numbers, certainly if that
4 were to come to the panel, I would want to see more data.

5 The next category of anterior chamber lenses
6 are iris-fixation lenses. A Worst-Fechner lens developed
7 in the late '80s, one particular study just reported two
8 eyes of greater than 50 percent cell loss, but they didn't
9 provide any mean cell density analyses or morphometric
10 analyses. A separate study in 1996 by Perez-Santonja
11 reported a 13 percent loss at one year, another 4.6 percent
12 at two years, indicating to me that the central endothelial
13 cell loss did not stabilize over a two-year period. The N
14 was only 30 in this particular study.

15 The Artisan or Iris-Claw lens, the prior
16 nomenclature which is unfortunately called the Worst Iris-
17 Claw lens --

18 (Laughter.)

19 DR. GRIMMETT: Those are synonymous as far as I
20 understand. There may be some design changes I'm not aware
21 of, but I believe they're in the same family.

22 It had reported some distances from the corneal
23 endothelium in various publications. For example, a -15
24 diopter lens with a 3.2 millimeter anterior chamber leaves
25 1.97 millimeters from the corneal endothelium. Because

1 these lens are fixated at the iris itself and not vaulted
2 above the iris, certainly you would expect more distance,
3 and this is what we're seeing. There's more distance from
4 the endothelium. So the angle-supported lenses, in my
5 view, given their proximity to the cornea, have a higher
6 risk to the corneal endothelium than, let's say, an Iris-
7 Claw lens, assuming that the anterior chamber (inaudible)
8 is constant, whatever it happens to be.

9 Looking at some data regarding endothelial cell
10 loss on the Artisan lens, a study by Menezo showed a 6.6
11 percent loss at 12 months for 109 eyes and 9.22 percent at
12 two years, but I think it's instructive just to look at the
13 differences. It's 6.6 percent preop to year 1; 2.63 year 1
14 to 2; 2.5 percent year 2 to 3; 1.74 year 3 to 4.

15 They found that the cell loss correlated to
16 increased power of the lens, the thicker lens -- that is,
17 perhaps closer to the endothelium -- and shallower anterior
18 chamber depth, which would also tend to bring the lens
19 closer to the corneal endothelium.

20 This particular study reported morphometric
21 measurements regarding percent hexagonality and coefficient
22 of variation.

23 They also find some changes that I found
24 instructive. They found that there were statistically
25 significant decreases or changes at six months and then

1 return at approximately two years with gradual resolution
2 back to near preop factors in their four-year study.

3 Then the final category of lenses would be
4 posterior chamber lenses, which we've heard a great deal
5 about already here today by the speakers. Dr. Werner's
6 excellent presentation showed us pictures of the Fyodorov
7 posterior chamber lens, which is not in general clinical
8 use, with cell losses of 10 percent at 12 months.

9 Interestingly, the Chiron Adatomed lens that
10 Dr. Werner indicated was pulled because of cataract
11 formation, at least the two studies I pulled, did not
12 report any endothelial cell data whatsoever, reports by
13 Fechner and Marinho. I didn't locate any other studies on
14 those lenses.

15 The Staar Surgical ICL that we heard about
16 today already, one study by Zaldivar of 124 eyes didn't
17 report endothelial cell counts, and in another study by
18 Arne, cell loss was actually remarkably low, 2 percent at
19 12 months and 2 percent at 24 months, and no eye had an
20 endothelial cell loss greater than 3.8 percent at one year
21 in 58 eyes. Since these lenses are in Phase III trials,
22 I'm certain there's some data targeted and reviewed by the
23 FDA in a confidential fashion regarding larger sample
24 sizes.

25 Knowing what the literature has to say about

1 endothelial cell loss and the differences between these
2 three types of lenses would help us give guidance regarding
3 future manufacturers' studies. It's important, as Dr.
4 Edelhauser alluded to, to look at normative data, so that
5 we know what the natural rate of loss is for corneal
6 endothelium in order to gauge our comments.

7 Rates of normal endothelial cell loss range
8 between .3 and 1 percent per year, depending on where you
9 look. I've quoted the various studies. The Bourne article
10 in 1997 with the .6 percent rate is frequently quoted, but
11 there are slight differences in the literature in that
12 regard. It's in somewhere in that range.

13 As far as surgical procedures, we all know that
14 operative procedures can create both a direct surgical
15 instantaneous hit to the endothelium as well as possibly
16 change the annualized cell loss rate. Cataract surgery as
17 far as Bourne's article in 1994 reported a mean 2.5 percent
18 cell loss per year over a 10-year period. It's important
19 to realize that these were intracapsular and extracapsular
20 surgeries with iris-sutured lenses, transiridectomy clip
21 lenses, and a few posterior chamber lenses. There are some
22 other articles that indicate that cataract surgery causes a
23 1.1 percent cell loss rate per year, and a Werblin article
24 in '93 said that one year is an 8.8 percent loss, and
25 you'll note that in the FDA table previously shown by Ms.

1 Lochner and in some of the assumptions I'll later make, we
2 just round that figure to a 10 percent surgical loss at the
3 time of the procedure, which obviously is different from
4 the future cell loss rate that may be increased over
5 normal.

6 As a point of interest, penetrating
7 keratoplasty has a 7.8 cell loss per year.

8 I went ahead and quoted the age-stratified
9 normal endothelial cell density values, both from the Yee
10 article in 1985 as well as a pathologic analysis. The
11 difference between the studies, one is specular microscopy,
12 the other is pathologic. I found they were similar in
13 terms of mean values. The lower age bracket, age 20 to 29,
14 they start at 2,900 cells and by the time you get up to age
15 80-89, there are 2,300 mean cells. The main difference is
16 in the standard deviations.

17 Dr. Edelhauser earlier indicated that the non-
18 contact robo data agrees with the Yee data, and I quite
19 frankly find the standard deviation values in the path
20 study to be huge, plus or minus 500, plus or minus 690. In
21 talking with Dr. Edelhauser on the break, he's going to
22 look further into the Moller-Pedersen article regarding the
23 standard deviations, but I tend to gravitate toward the Yee
24 article with the standard deviations.

25 The peripheral cornea, if measured, is known to

1 have increased cell density. Per some data that Dr.
2 Edelhauser provided, there's a 5.8 percent increase in cell
3 density in the paracentral region and a 9.8 percent
4 increase in the peripheral cornea.

5 The next category that I did, knowing the peer-
6 reviewed literature and now knowing normative data, is I
7 tried to determine what would be some thresholds for
8 unacceptable rates of endothelial cell loss. There are two
9 different questions, I believe.

10 One, as a panel member, when an application
11 comes to panel, we all be concerned about what is an
12 acceptable cell rate when we get that application with a
13 particular observed cell loss rate. The reason it's
14 important to get some judgement about acceptable cell loss
15 rates now is it will help us give guidance regarding what
16 or how low should thresholds be that we're actually trying
17 to screen for to make the determination how big the sample
18 sizes must be. So we have to have a sense for what are the
19 maximum and minimum ranges for thresholds that we're even
20 looking at, so we can give some type of guidance regarding
21 sample sizes.

22 I certainly don't expect anyone right now to
23 define an acceptable cell loss rate. That will certainly
24 be a hotly debated topic once an application's received,
25 but I think it's instructive that we go through the

1 exercise to see what our limits are, maximum and minimum
2 unacceptable rates of cell loss.

3 I've approached this argument by first looking
4 at life expectancy data. The RP-2000 mortality table is
5 based on a study of the mortality experience of pension
6 plans conducted by the Society of Actuaries and was in
7 response to pension legislation that directed the Secretary
8 of Treasury to promulgate the use of updated mortality
9 tables for various pension calculation purposes.

10 According to that table, the life expectancy
11 for a 21-year-old male is 58 future years or an age of
12 death of 79. The life expectancy for a 21-year-old female
13 is 62 future years, so an age of death of 83. Those are
14 United States data.

15 Realize that depending on your entry date,
16 you'll have change, obviously, to your age of death. If
17 you enter at age of 80, you don't have an age of death of
18 79.

19 (Laughter.)

20 DR. GRIMMETT: But I used it as a fixed value
21 for this particular analysis, so as to not get confused
22 with multiple iterations of the tables. Suffice it to say
23 that when you enter at 20, 30, or 40, it may only differ by
24 a few years in terms of your age of death.

25 The minimal acceptable corneal endothelial cell

1 density value is critical, and Ms. Lochner presented a
2 table regarding a 1,200 threshold and a 1,000 threshold
3 near death. I picked values on either side of the range to
4 show the range of acceptable cell loss for sake of
5 argument.

6 In one of Dr. McCarey's earlier versions of a
7 slide, he quoted a minimal acceptable rate of 1,500 cells.
8 Whether we call that number 1,500 or 1,400 or 1,300 or
9 whatever it's picked as is not really the crucial value,
10 but I think a larger number has the argument that if these
11 patients get into any kind of trouble with their phakic
12 IOL, number one, it allows you to do a surgical procedure
13 to possibly correct that, whether it be explantation of the
14 IOL or manipulation of the IOL.

15 Secondly, we all know as clinicians that as
16 patients enter their early 70s, they have a much higher
17 likelihood of having us see them with cataracts just from
18 age-related phenomena. Whether or not it's increased with
19 phakic IOLs remains to be determined, but this larger
20 target value near or at the age of death will certainly
21 allow a future intraocular surgical procedure, such as
22 cataract surgery. If we run all of our calculations and
23 run them right down to the wire and leave them the bare
24 minimum and leave half of the years to the time of death,
25 you are not giving that patient the opportunity to have any

1 intraocular intervention of any kind. So that's why I
2 picked the number of 1,500, and Dr. McCarey in an earlier
3 version happened to put that number there, so I went with
4 it.

5 The next number I used is a potential corneal
6 edema of 800 cells per square millimeter. Acceptance of
7 this particular target will not allow a future intraocular
8 procedure, in my opinion, and certainly if I had a patient
9 of advanced age with immense nuclear sclerotic cataract
10 requiring higher phaco times, I would not be comfortable
11 performing phacoemulsification with entry cell count of
12 800. I would certainly advise the patient clinically that
13 they would have a higher chance of having postoperative
14 corneal edema.

15 In the literature, there's been a quote of 500
16 cells for imminent corneal decompensation, but I think it
17 depends on the actual function of the remaining cells. So
18 the exact, precise figure is not locked down, but in my
19 opinion, clinically it's somewhere in the area of 800.

20 The assumptions I made for my threshold
21 analysis is that the endothelial cell loss, as calculated,
22 was an instantaneous, exponential endothelial cell loss
23 rate, that they lose 10 percent at the time of the surgical
24 procedure, and then they have a continuous stable
25 annualized exponential cell lose rate after that time. I

1 used a round down analysis of any remaining fraction. I
2 used .99 rounded down because partial cells do not survive.

3 I did not alter the life expectancy target
4 values, as I had previously mentioned.

5 Table 1. I'm not going to go through this,
6 obviously, in detail, but all I did is create an Excel
7 spreadsheet with a formula. I set the target at the end at
8 1,500, I set the percent drop per year, and I was back-
9 calculating the cells you would need to enter the study.
10 That's all it is. So I was trying to find out what cell
11 count would you need at age 21 in order to end up with a
12 cell count of 1,500 at the time of death for male or
13 female, and then supplying the percent drop per year.

14 For that particular assumption, 1,500 cells at
15 the time of death, 2 percent cell loss per year, I found,
16 for example, age 21, you need 5,900 cells if you're female
17 and you need 5,400 male. Not possible. We looked at the
18 normative data and those numbers exceed the normative data.

19 So the conclusion from that scenario is that if
20 a 1,500 cell target is desired at death, a 2 percent annual
21 cell loss rate is not acceptable. It's not possible to
22 enter with a high enough cell count for all but the very
23 older age ranges. So everyone will fall below 1,500
24 because you can't enter with a high enough cell count.

25 I did that same type of analysis trying to

1 bracket what would be inclusive of all ages in order to
2 allow everyone to enter the study. A 1.5 percent cell loss
3 rate, which is Table 2, they're still entering with some
4 very high cell counts. At age 21, for example, 4,300 for a
5 woman and 4,100 cells for a male.

6 The conclusion from that scenario with 1,500
7 target and 1.5 percent cell loss is that no one could enter
8 with a high enough cell density, except those older than
9 approximately 50. So that cell loss rate is unacceptable,
10 the 1.5, if you're targeting for 1,500, that target, of
11 course, having the advantage of allowing someone in the
12 future an intraocular surgical procedure.

13 It turns out a .9 percent cell loss rate per
14 year allows all ages to enter with a reasonable cell count
15 that could be achieved based on the normative data. The .9
16 percent cell loss per year on Table 3 shows all the entry
17 cell requirements, and they reasonably match or were below
18 the normative data for all ages, telling us that the cell
19 density values will approximate or exceed 1,500 at the time
20 of death for all patients entering the study.

21 So the .9, assuming you'd want 1,500 to be your
22 target, is the maximum allowable rate that's inclusive of
23 all ages, and that rate is approximately 50 percent higher
24 than the normal of .6 percent cell loss rate per year.
25 We'll later hear from our statisticians regarding the

1 feasibility of trying to read for that low level of cell
2 loss rate in the face of the variability and precision
3 issues that Dr. Edelhauser brought up, because it's my
4 opinion, based on review of some initial data, that the
5 sample sizes would have to be unreasonably large.

6 Just as an example from the literature, an
7 angle-supported phakic IOL from Alio, 1999, has a cell loss
8 rate of .72 percent from years 2 to 7 after experiencing a
9 6.83 percent loss from preop to year 2. So based at least
10 on something in the literature that the numbers are not
11 huge, it's doable.

12 Looking at it a different, looking at the
13 target of 800, so running right up to the edge of corneal
14 edema, I can give you all the tables, but the number that I
15 am giving you is 1.9 percent cell loss rate per year. That
16 number of 1.9 percent cell loss rate per year would allow
17 all patients to enter based on the normative data.

18 So if you desire an 800-cell target value, a
19 1.9 percent rate of annual endothelial cell loss is the
20 maximum allowable rate of loss that's inclusive of all
21 ages, and that's about three-fold higher than the normal .6
22 percent rate per year, and at least based on the Menezo
23 1998 Artisan lens data, I have an average cell loss rate of
24 2.28 percent from years 1 to 4.

25 Based on all these analyses and review of the

1 literature, I went ahead and just prepared a bunch of
2 different issues that I would want to see or issues that I
3 would want considered in phakic IOL studies regarding the
4 endothelium, and then we can back to the specific
5 endothelial questions. Most of the question issues I
6 believe will be covered.

7 General issues that I would have for all phakic
8 IOL studies is that certainly endothelial cell density
9 measurements are mandatory. We saw in some of the
10 published literature it did not report endothelial cell
11 densities whatsoever, which I think is unacceptable for a
12 new product of this design because it's a critical issue
13 for the survival of corneal health.

14 Certainly, a central count is mandatory. A
15 peripheral count will be important, especially if it's an
16 anterior chamber lens. You'd want a cell density in the
17 region near the IOL edge or in the area of minimum distance
18 between the IOL and endothelium. That was learned from
19 early lens design. So especially with an angle-supported
20 lens, I would be highly interested in reviewing the
21 peripheral cell count.

22 It is my belief that morphometric analyses are
23 mandatory. We know that analysis of cell shape and size
24 provides a more sensitive indication of endothelial cell
25 damage than cell density alone. I understand the

1 limitations of the issues, as Dr. McCarey and Dr.
2 Edelhauser have outlined, regarding the coefficient of
3 variation and the reliability, especially with the
4 algorithms of the non-contact robo. Notwithstanding those
5 limitations, I still feel it's an important variable to
6 consider and may give us some additional information.

7 Just as a matter of course, corneal pachymetry,
8 corneal functional analysis, most certainly would be
9 suggested.

10 The duration of the study that I would be
11 interested in prior to that study coming to panel would be
12 three years. Based on the literature, it seems like there
13 is an initial larger decrease and subsequent next decrease,
14 and then stabilization after a year or perhaps two to
15 three. Depending on what is seen in panel, it would be my
16 guess that there would be discussion of postmarket
17 surveillance of endothelial cell data for possibly a year
18 or two more.

19 There are some data that it may take four years
20 to see the morphometric data return to baseline levels and
21 ensure stability, but I do not believe, based on what I've
22 seen in the peer-reviewed literature and based on my
23 concern about the corneal endothelium, especially for
24 anterior chamber lenses, that I would be comfortable
25 approving a lens with only two years of data. So for me,

1 three years would be the minimum I would favor.

2 I would favor a longer endothelial study
3 duration, perhaps postmarket-mandated, for higher risk
4 factors, such as angle-supported lenses, which have a
5 higher risk to culture the endothelium than iris-fixated
6 lenses, which are closer than posterior chamber lenses.

7 The same issue. Thicker lenses are closer to
8 the endothelium than thinner lenses. So we would want
9 longer follow-up for those.

10 A shallower anterior chamber depth, such as
11 hyperopic patients, would have a higher risk than deeper
12 anterior chamber depth, perhaps if the IOL is closer to the
13 endothelium, and certainly chronic anterior chamber
14 inflammation, as previously mentioned, would be a higher
15 risk factor for endothelial loss than a quiet anterior
16 chamber.

17 Interestingly, in one early, approximately 50-
18 page paper by Drews in 1991, he went over study parameters
19 versus the FDA grid regarding what he would expect for
20 phakic IOLs, and he recommended a five-year study duration.
21 I certainly don't disagree with that, given the importance
22 of corneal endothelium.

23 With respect to anterior chamber phakic IOLs,
24 we know that the optic-endothelium distance plays an
25 important role in potential endothelial damage. Therefore,

1 high-resolution ultrasound, as mentioned by Dr. Werner, I
2 believe would be mandatory to disclose the optic-
3 endothelial distance and distances between other eye
4 structures and the IOLs, such as the crystalline lens.

5 The peripheral endothelial cell density and
6 morphometric measurements would be mandatory in my opinion
7 in the region of the IOL optic edge, in addition to the
8 central endothelial analysis, because if the examination is
9 limited to the central cornea, it may fail to detect
10 significant endothelial injuries, and then specular images
11 can show significant morphologic changes over the edge of
12 the IOL in the absence of central cell density changes.

13 I would caution that a preop history of eye
14 rubbing may be a contraindication for entry into the study,
15 especially when we're talking about such low tolerances
16 between the distance between the IOL edge and the corneal
17 endothelium.

18 Depending on the IOL design, such as an angle-
19 supported lens, a manufacturer must specify a minimum
20 anterior chamber depth that contraindicates IOL insertion.

21 Because chronic inflammation is a known factor
22 in endothelial damage and because some studies have
23 disclosed chronic anterior segment inflammation after
24 phakic IOLs using a laser flare cell meter, I would
25 recommend laser flare fluorophotometry to evaluate chronic

1 anterior chamber inflammation, and just raise the question
2 would iris fluorescein angiography be of any help if a lens
3 were either rubbing the iris or directly affixed to the
4 iris? So that would be an issue that perhaps in a small
5 subset of patients in early studies may be relevant.

6 Posterior chamber phakic IOLs certainly would
7 have advantages for corneal endothelium. They would have a
8 maximum distance between the IOL and the endothelium. They
9 would avoid optic-endothelial contact, but of course, their
10 location behind the iris would have a higher risk for
11 pigment dispersion and introduction of cataract. If these
12 lens induce chronic anterior segment inflammation, they
13 certainly may have secondary effects on the endothelium, so
14 I do believe that endothelial studies are also important
15 for posterior chamber lenses, irrespective of the fact that
16 they're the furthest lens away from the endothelium.

17 As far as that table -- now getting back to the
18 two questions, I added these last two portions right at the
19 end. Getting back to the table, what would be the
20 recommended minimum endothelial cell density to enter these
21 studies? And we saw that table put forth by the FDA.

22 I see two approaches to attack the problem.
23 One is standard deviations. You could look at accepted
24 normative data and say you don't want anyone entering the
25 study outside of, let's say, lower than, let's say, two

1 standard deviations.

2 Let's talk about two standard deviations.

3 Assuming the standard deviations are reasonable and your
4 studies are reproducible, you would want to exclude anyone
5 out in that 2.5 percent tail, and I went ahead and just
6 listed the values as they would range from age 20 at about
7 2,700 going down to 2,400 or so by age 60. That's sort of
8 one way of looking at it.

9 The other way to look at it would be the way
10 that the FDA has approached it. Look at do you have enough
11 cells to make it near the age of death? And that's the
12 other way to look at it.

13 In reviewing the FDA Attachment B, the
14 threshold values selected don't exactly -- they do not
15 likely allow a patient to undergo a secondary ocular
16 procedure such as cataract surgery or further phakic IOL
17 manipulation during the life of the patient. I would have
18 concerns as a clinician if a patient had a cell count of
19 1,000 with a moderately dense cataract doing a phaco
20 procedure. I have to be worried that I might tip them over
21 into clinical corneal edema.

22 So I think we have to understand that the chart
23 values selected, and I selected 1,500 and the FDA selected
24 1,200 and 1,000, are variable depending upon what our
25 expectations are for these patients to have elective

1 procedures down the road, and it's impossible to precisely
2 determine the minimum entry requirements without knowing
3 the exact rate of endothelial cell loss per year at the
4 various sites that run these. So when we discuss that
5 particular table, I'm certain all those issues will come
6 into play.

7 The issue that I haven't addressed, and I'll
8 leave it to the statisticians, is regarding sample size
9 analysis. Drs. Edelhauser and Bernie McCarey talked about
10 that in the best case scenario, the best precision they can
11 get on endothelial cell measurements is 2 percent and real-
12 world data is at 9 percent according to Dr. McCarey. I
13 would like to know how that translates into the standard
14 deviation values that our statisticians have calculated
15 regarding sample sizes, if that's known, to see what kind
16 of numbers we would actually need if we wanted to screen
17 for the lower rates of annualized loss.

18 That would conclude my introductory comments at
19 this time.

20 DR. WEISS: Thank you very much for a very
21 thorough, as usual, presentation.

22 DR. GRIMMETT: Thank you.

23 DR. WEISS: Dr. Matoba has a question. Then
24 we'll go around.

25 DR. MATOBA: Mike, I might have missed this,

1 but when you did your calculations, did you assume any
2 amount of cell loss from the actual surgery?

3 DR. GRIMMETT: Yes. I assumed a 10 percent
4 instantaneous loss at the time of the surgery based on the
5 (inaudible) of having 8.8 percent with a phaco procedure.

6 DR. WEISS: Mr. McCarley, did you have a
7 question?

8 MR. MCCARLEY: I just had a couple of
9 questions. The 9 percent that you're referring to is what
10 Dr. Edelhauser and Dr. McCarey were talking about? Weren't
11 those the controlled laboratory percentages and not actual
12 real data?

13 DR. GRIMMETT: I believe it was real data from
14 a real study, but we could probably have Dr. McCarey answer
15 it, but Dr. Edelhauser's number of 2 percent precision was
16 sort of in the best of hands what is the best precision.

17 MR. MCCARLEY: Right, exactly.

18 DR. GRIMMETT: And Dr. McCarey's number was in
19 order to -- it might be better to have Dr. McCarey answer,
20 but I believe it was real data. I'm not sure how many eyes
21 were in the study.

22 DR. WEISS: You can over there. As I recall,
23 it was different centers, with the variation 9 percent
24 between different people reading the same thing.

25 DR. MCCAREY: Yes, it was real data from a

1 clinical trial that Medennium's working. It's the control
2 data. There were 58 patients, seven clinical sites. I
3 just started to collect the data and asked the question and
4 came up with that answer.

5 MR. McCARLEY: Okay. I wonder is the FDA able
6 to provide the panel, especially the voting members and
7 reviewers, with data that they currently have on phakic
8 intraocular lenses. In other words, the endothelial cell
9 data that has come already out of the studies? There must
10 be over 2,000 implants so far over --

11 MR. WHIPPLE: You mean the ones that are under
12 IDE?

13 MR. McCARLEY: Pardon?

14 MR. WHIPPLE: Dave Whipple. You mean the ones
15 that are already under IDE?

16 MR. McCARLEY: Yes, correct.

17 MR. WHIPPLE: We can summarize it.

18 MR. McCARLEY: In other words, we're reviewing
19 the literature and making decisions based on the
20 literature, and a lot of these have been recognized as
21 being small studies. We don't know what methodology was
22 used to validate -- the validity of the methodology and so
23 forth. A lot of reference goes back to original studies
24 that are even from the 1970s, from Bourne and group, and of
25 course, he did a later study, but he did longitudinal study

1 showing I think 10 years' follow-up.

2 But I'd be curious, one, is can the FDA provide
3 the panel members with the data that they have now to show
4 them what the standard deviation is right now?

5 MR. WHIPPLE: Yes. They're special government
6 employees. As long as they kept it amongst themselves,
7 they could use it, yes.

8 MR. McCARLEY: Right. I would just say before
9 a recommendation is made on limits and so forth, I would
10 say you have real-time data, real data on real phakic
11 intraocular lenses, that is up to date that you may want to
12 consider before you make a final recommendation.

13 DR. WEISS: Well, actually, the advantage of
14 this sort of meeting is that everyone's entitled to their
15 opinion, there is no final vote, and there has to be
16 consensus. So basically, the FDA will use all the
17 information gathered here today, and including any other
18 information we think would be helpful to obtain in the
19 future, to come up with some final recommendations on their
20 own.

21 Dr. Grimmett?

22 DR. GRIMMETT: Yes. Michael Grimmett. I
23 alluded to it in my introductory comments that while it
24 would be wonderful to have the Phase III endothelial cell
25 loss data versus what is published in the literature, we're

1 not trying to set a threshold rate of what is acceptable to
2 this panel or what is not acceptable. That's not even the
3 purpose of this. But I approached the analysis in that
4 fashion to show what the edges of the approach would be, so
5 that we can get some data on how many patients would we
6 need to screen in that range. That was my goal.

7 Understandably, the literature is not giving us
8 that much guidance regarding actual cell loss rates because
9 most of the lens designs have been just started and we're
10 now on the newer lens designs. So we actually don't know
11 what the newer loss rates are.

12 DR. WEISS: Dr. Mathers?

13 DR. MATHERS: Yes, I want to compliment Dr.
14 Grimmett on his excellent analysis, and it may sound like
15 it is rather strict, but I don't think that using the
16 actuarial tables as you have is in any way offbase. In
17 fact, it's probably ultraconservative because, as we have
18 seen, as medical advances continue, it isn't impossible to
19 look at ages beyond what you are saying, and we all know a
20 lot of 80-year-olds who do not feel like dying right now.

21 (Laughter.)

22 DR. MATHERS: And who are probably going to
23 live longer as these issues become more relevant. So 10
24 years from now, with one advance in cardiac pathology, this
25 would become not strict enough.

1 I would also like to ask who you think should
2 do the reading and the entrance qualification to get into a
3 study like this, because clearly the endothelium, which is
4 actually quite accessible as the cell to study, can be best
5 read perhaps by a central office, and should the entrance
6 requirements be that a central group does the initial
7 reading to get in? What would you think of that?

8 DR. GRIMMETT: Michael Grimmert. Based on the
9 comments of Dr. Edelhauser, given all the variability that
10 exists if the technician is not trained properly and all
11 the parameters to be analyzed, I obviously am in favor of
12 the highest trained, highest precision measurer because of
13 the critical nature of endothelial cell loss over time as
14 our population ages.

15 DR. MATHERS: It wouldn't be difficult to do
16 this with digital capture. You can transfer these images
17 instantaneously and you can then decide whether you have
18 good images, because it would seem that a lot of this
19 depends upon having good images to start, so that you know
20 how to get good data afterwards, and that's certainly
21 possible on an instantaneous basis.

22 DR. WEISS: The way the document reads at the
23 present time is that "The use of a reading center is
24 strongly recommended. If the use of a reading center is
25 not possible, the sponsor should establish a protocol for

1 collection and analysis of images to be used by each
2 participating site."

3 Would you then change it from strongly
4 recommended to required or you'd leave it as strongly
5 recommended?

6 DR. GRIMMETT: This is Michael Grimmatt again.
7 I think there are always multiple ways to skin a cat, and
8 if a particular study or a sponsor can demonstrate the
9 reliability that they have internal mechanisms to validate
10 precision and validity and it appears to be equivalent to a
11 standardized reading center, I think there is always
12 flexibility in that regard.

13 DR. WEISS: Dr. Huang?

14 DR. HUANG: My first comment is that regarding
15 Mr. McCarley's comment earlier, using the FDA existing data
16 is almost like using the soccer player to be the referee.
17 You know, that we are judging the safety of the data and
18 then using the data to be the reference for its own safety.
19 I think it's questionable.

20 The second comment I would like to make is also
21 that Dr. Grimmatt used a very nice life table actuary to
22 analyze this, but I think the endpoint is a little bit
23 strict because we don't do cataract surgery at the time of
24 the death. We do that cataract surgery maybe hopefully
25 five or 10 years before the patient's life expectancy in

1 order to improve their quality of life.

2 So maybe we have two points that we can use.

3 You know, maybe one at one point is 10 years before death
4 and then at the death point, and then to find a reasonable
5 middle ground for the starting point.

6 DR. GRIMMETT: This is Michael Grimmett again.
7 Just a quick response. I first attempted to do the
8 cataract surgical procedure 10 years before the age of
9 death with an increased rate of annual cell loss. The
10 table became so complicated in the formula that I would
11 have to do to change the rate of cell loss midstream and
12 back-calculate to the entry cell data that I couldn't get
13 the spreadsheet to work in that regard.

14 I took took the give-up approach. You know, I
15 was on vacation.

16 (Laughter.)

17 DR. WEISS: Well, maybe in that case, the
18 guidance could be to the FDA to recalculate this with the
19 average age of cataract surgery as the final.

20 Yes, Dr. Huang?

21 DR. HUANG: I'm just joking. I say, you know,
22 maybe he can take another vacation to calculate that.

23 (Laughter.)

24 DR. WEISS: Dr. Burns?

25 DR. BURNS: Yes, but I think all things being

1 equal, if you move the cataract surgery earlier, you're
2 going to actually end up with stricter numbers because
3 you'll have the loss from the surgery and then you'll have
4 an increased growth and you're going to come lower on the
5 curve. So I think that will only make things stricter.

6 DR. HUANG: Andrew Huang again. I think this
7 is just a recommendation. Not all the patients eventually
8 are going to need cataract surgery, but I certainly
9 understand that there will be additional loss, and I don't
10 know which one is greater, 2 percent annual loss versus the
11 10 percent initial loss. So statistics will tell us.

12 DR. WEISS: I think it would be up to the FDA
13 to crunch the numbers both ways to see what the differences
14 are and if they're clinically relevant or not.

15 Dr. Grimmer?

16 DR. GRIMMETT: I would just point out on the
17 tables, you realize that it is a 10 percent initial loss
18 for the phakic IOL surgery. Please realize there are some
19 data on some of the more current studies that say that that
20 initial loss may at that rate, 5 to 6 percent. The earlier
21 studies did have a higher rate.

22 So I used a cataract surgery
23 phacoemulsification initial rate that may not exactly be
24 true for phakic IOLs. I used the most conservative
25 approach to make sure that we leave people with enough

1 cells at the end.

2 DR. WEISS: While we're on the endothelial cell
3 topic, I just want to pose one other question relating to
4 this, and then perhaps Ms. Lochner could come up and then
5 we can go on with our other presenters and other
6 discussions.

7 The other question that I would pose to the
8 panel, as I already did to our experts, is what would you
9 think would be the reasonable number of months to tell a
10 patient they would have to be out of contact lenses? If
11 you have an opinion.

12 Dr. Bullimore?

13 DR. BULLIMORE: None. I mean, I think keeping
14 it --

15 DR. WEISS: Not even a week?

16 DR. BULLIMORE: Well, let me say less than a
17 month.

18 DR. WEISS: Less than a month.

19 DR. BULLIMORE: I mean, I think trying to sort
20 of take them out of contact lenses with the expectation
21 that the endothelium's going to change in any meaningful
22 fashion, based on my experience and what I heard from the
23 experts today, is futile.

24 DR. WEISS: So you would keep them out long
25 enough to get a proper keratometry or corneal topography?

1 DR. BULLIMORE: I mean, I'd use the same
2 guidelines that exist for corneal refractive surgery.

3 DR. WEISS: Gee, I think we just eliminated two
4 questions with that comment.

5 Dr. Mathers?

6 DR. MATHERS: Yes, I strongly agree with that.
7 Aside from the fact that it would be really difficult to
8 have patients go through that period of time, since the
9 polymethism afterwards doesn't really evolve very quickly
10 and we don't know exactly the significance of it, I think
11 that it's not terribly relevant to have them out of their
12 lenses a long period of time, except to establish their
13 refractive error issue.

14 DR. WEISS: Dr. Bullimore?

15 DR. BULLIMORE: And I think we need to be
16 generalizable, and these products and these procedures are
17 going to be done on people who, to a large extent, are
18 long-term contact lens wearers. We saw yesterday that
19 something like 80 to 90 percent of the patients enrolled in
20 a study for low myopia or low to moderate myopia were
21 contact lens wearers. So in the high group, it's going to
22 be probably even higher.

23 DR. WEISS: Another question, which I warn you
24 in advance I'm going to limit the discussion on because it
25 will come up again, is with this in mind, would you only

1 want one eye done at a time or one eye done and use the
2 other eye as a control for the endothelial cell study?

3 Dr. Bullimore?

4 DR. BULLIMORE: Again, I think you don't want
5 to burden the patients too much, and enrolling them in,
6 say, a three-year study where they can only have one eye
7 done with a device and the other eye has to wait I think is
8 an unreasonable burden to be placed on the patient. We
9 have a lot of historical control data. It seems to me we
10 have some historical data on endothelial cell count as a
11 function of age. We can use that, and what we're looking
12 for, I guess, are the extreme or the worst cases where
13 people really do loss a lot of endothelial cells within a
14 relatively short amount of time, and I don't think we need
15 a control group to necessarily look at those event rates.

16 DR. WEISS: Dr. Grimmer?

17 DR. GRIMMETT: Dr. Grimmer. My concern with a
18 unilateral one-eye study would be for quality of life for
19 the patient. These patients would be typically those that
20 don't qualify for other refractive procedures, and hence
21 they have a higher range of myopia. There is significant
22 (inaudible) of emmetropia. So in a three-year duration
23 study, I think that would be unwieldy and probably not very
24 reasonable for the patient. So I'm in agreement with Dr.
25 Bullimore.

1 DR. WEISS: Dr. Mathers?

2 DR. MATHERS: Yes, I'm also highly in
3 agreement. I think we'll need to discuss again the time
4 delay between the first and second operation, but three
5 years is too long.

6 DR. WEISS: Then I would sort of conclude --
7 Dr. Bandeen-Roche?

8 DR. BANDEEN-ROCHE: Yes, just very briefly.
9 You know, statistically, we'd obviously like to have a
10 contralateral, but I am absolutely swayed by the quality of
11 life considerations as long as we have good quality control
12 data. So that would include both a precise estimate of the
13 rate of loss, but also of the variability in rates.

14 DR. WEISS: Okay. So then I would conclude
15 that the fact that the contact lens issue in terms of the
16 change and shape of the cells and number of the cells might
17 still be evolving after the implant was placed in will be a
18 confounding variable, but not objectionable by the panel.
19 Fine.

20 Ms. Lochner, perhaps we can go on.

21 MS. LOCHNER: For analysis of the crystalline
22 lens for lens opacity, we currently recommend that "The
23 natural lens should be evaluated preoperatively and at each
24 of the postoperative intervals. The level of evaluation
25 should be commensurate with the risk of cataractogenesis or

1 lens changes identified by the risk analysis performed by
2 the manufacturer. For phakic IOLs where the design or
3 surgical procedure may lead to lens changes, a grading
4 system or quantitative method should be used to evaluate
5 lens changes over time. For IOLs for which lens changes
6 are not an identified risk, qualitative observation may be
7 adequate."

8 The analyses should include the number of
9 patients with any change in the appearance of the lens
10 stratified by the type of change and the number of patients
11 with clinically significant lens opacities, and the term
12 "clinically significant" is as yet undefined.

13 We are asking for panel comments on whether you
14 believe evaluation of lens changes should be requested of
15 all sponsors or whether this evaluation should only be
16 performed if the sponsor's risk analysis warrants and
17 whether you have any specific recommendations for defining
18 the term "clinically significant" lens opacities.

19 We're also asking for your recommendations
20 about the use of quantitative methods for measurements of
21 lens changes versus the use of grading systems, and finally
22 we'd like your thoughts on the duration of the study and
23 request that you specifically discuss the length of follow-
24 up you believe would be adequate for panel review of this
25 cataractogenesis outcome.

1 Now, Dr. Mathers will provide his review.

2 DR. MATHERS: Thank you.

3 Regarding the first question we are asked to
4 address, the evaluation of the lens changes, it is this
5 reviewer's opinion that all phakic IOLs need to be
6 evaluated for any cataract or changes in the lens.
7 Anything that happens essentially inside the anterior
8 chamber in front of the lens is an issue here, and even
9 something that would not touch the lens or is known to be
10 touching the lens would still be a problem. Certainly,
11 central touch is not the only issue nor is just peripheral
12 touch. There's anterior chamber inflammation, and all of
13 this can affect the lens changes over time. So any
14 perturbation would be of interest and all cataract
15 processes need to be assessed.

16 In this direction, I think it's going to be
17 important not just to look at the lens itself, but to look
18 at the source of the possible problem, such as we have
19 heard regarding flare assessment or anterior chamber
20 inflammation, and I think that it would useful to measure
21 flare in these patients with perhaps the laser flare
22 system, and also, regarding the evaluation of the cataract
23 process, that the sizing and the structure of the anterior
24 chamber is going to be a key issue. So the use of high-
25 resolution ultrasound I think is going to be extremely

1 important in evaluating this.

2 Regarding the clinical significance of the lens
3 opacities, this is a more subjective area and we don't have
4 very good data on this, but this of course refers to vision
5 changes, and in this regard we need I think to be as
6 precise as possible because when one measures vision under
7 various circumstances, you get very, very different
8 results.

9 So if you have a central opacification of the
10 anterior subcapsular area, it will not show up in vision
11 testing in a dim room. You're going to have to do this
12 with a small pupil under conditions that will induce some
13 glare, and in fact, the more glare, the better. This is
14 not a non-real-life situation. I mean, when people are in
15 a high-light environment, which they often are, this comes
16 into play.

17 So I think glare testing is the most relevant,
18 but the only way that this should be assessed. For this,
19 we need to have careful measurements of the patients as
20 they enter the study and a standardized method of
21 evaluating the glare, and that high-glare settings should
22 be used.

23 Now, everyone will have some decrease in visual
24 acuity with a high setting of a glare. So we have to
25 decide whether we think one, two, or more lines of loss

1 compared with the preoperative evaluation would be
2 significant. My recommendation would be two or more lines
3 compared with what they had in the loss before.

4 I think that all lens changes should be
5 reported, not just the anterior subcapsular fibrosis, but
6 we also have anterior cortical changes which may result
7 from these anterior subcapsular processes.

8 The posterior subcapsular cataract is also an
9 issue because the anterior lens cells are not the only ones
10 going to be involved. As you build lenses that do not ride
11 or touch the central lens, they may ride in the periphery
12 of the lens, and they are going to get closer to those
13 cells on the outside which can migrate posteriorally, and
14 certainly may do so. So that becomes an issue, and the
15 possibility of inflammation that occurs with the anterior
16 chamber lens of this design may affect the development of
17 nuclear sclerosis, and this is going to be something that
18 is going to be harder to assess, but I think that we need
19 to at least monitor this.

20 In Part C, we're asked to comment on the use of
21 quantitative measures for the measurement of lens changes
22 versus the more semi-quantitative grading system, and by
23 quantitative here, we mean the assessment or the
24 visualization, the optical visualization, of the changes
25 underneath the anterior capsule.

1 There is no standard way to do this and the
2 examination of this is very light-dependent. If the
3 lighting system is a little bit off, you get a little
4 different view of this and it becomes harder to see, and
5 Dr. Werner's presentation on this was excellent, but even
6 she could not really tell us exactly how much better, say,
7 a photographic system is from a slit lamp system, which of
8 course is going to be subjective. But I believe that the
9 backlighting and retroillumination of the lens with high-
10 resolution color photography probably offers us the most
11 objective and reliable way to follow this over time.

12 The development of a scale to do this has
13 already been done. I don't think this needs to be
14 reinvented, but can be perhaps modified slightly. The LSCS
15 system can grade these opacifications and can be used for
16 all of this -- the anterior subcapsular, the posterior
17 subcapsular, and the nuclear -- and I think something like
18 that would be appropriate.

19 I would strongly recommend the use of digital
20 photography to perform this process and I also think it is
21 possible, as technology improves, that we will have a
22 better understanding through other means of visualization,
23 perhaps confocal microscopy of these cells now that this
24 becomes an issue and it becomes relevant to look at this.
25 But there are no standards for this now, and that would not

1 necessarily be used in an early study like this.

2 In Part D, we're asked to comment on phakic
3 IOLs and the length of time that it might be useful to
4 evaluate this for, and here, as we've heard with the
5 endothelium, when we're studying the endothelium, we have a
6 fairly clear endpoint. We have endothelial cell loss and
7 we can follow this fairly objectively.

8 With lens changes, it's much less objective,
9 and one study, noted recently, showed that there was some
10 change in light transmittance with these lenses, not
11 necessarily based on cataract, but I think that with the
12 long time span that we're talking about and the possibility
13 of chronic inflammation associated with this that is
14 subclinical, that at least three years would be necessary
15 to evaluate it. I think that the monitoring of this
16 perhaps should go on longer than that, but I think the
17 three years is probably enough to give us an idea of what
18 is happening.

19 The capsular process of cataract formation is
20 tied to a number of different other processes. Not just
21 lens touch, but, as I said, inflammation, and as the lens
22 is redesigned to minimize the cataract process, the other
23 issues, such as iris touch and development of pigment
24 dispersion and glaucoma, become more of an issue, and I
25 think that the industry or lens manufacturers will be

1 tempted to avoid the obvious problems of lens touch by then
2 shifting the burden to the back of the iris. I think that
3 this kind of monitoring is also going to be important. I
4 know it's not part of this particular issue, but I think
5 it's relevant because there are tradeoffs here. There is
6 not much space in this area and as you design the lens to
7 perform in one way, you then create other issues of
8 significance.

9 That is my summary.

10 DR. WEISS: Thank you very much.

11 I just have just two questions on things which
12 you've already mentioned, but I just want to clarify. What
13 would you call a clinically significant cataract?

14 DR. MATHERS: Clinically significant cataract
15 refers to a loss of a number of lines, but it's highly
16 dependent upon that is assessed, and I think that needs to
17 be assessed not just with standard -- well, our assessment
18 of vision can be done in the standard way in dimly lit room
19 to optimize vision, but it needs to go beyond that. We
20 need to have glare testing as well because our standard
21 measures of vision will not be adequate to pick up the kind
22 of changes we're going to see with capsular opacification.

23 DR. WEISS: So then to just restate that, with
24 glare testing, it depends on sensitive you want to be to a
25 very early cataract, and that's the question I'm asking, is

1 how sensitive do you want to be? So if we said let's it
2 bring it to the glare testing realm and say, okay, we're
3 going to determine that by a loss of X number of lines by
4 glare testing, is there an idea you have?

5 DR. MATHERS: Certainly, it has to be more than
6 one line, so I would say two.

7 DR. WEISS: So would you prefer, if I was going
8 to make you quantify it, would you then say a loss of two
9 lines by glare testing, rather than just a straight loss of
10 two lines without the glare testing or would you want to
11 say something else or it's totally unknown?

12 DR. MATHERS: I think without the glare testing
13 loss of one line would be important, would be significant.
14 With glare testing, it's going to be two, because the glare
15 testing is much more sensitive.

16 DR. WEISS: And would you require contrast
17 sensitivity data or that would be optional?

18 DR. MATHERS: I think contrast sensitivity data
19 also should be included.

20 DR. WEISS: Okay. So you'd like to have both,
21 but at least at this point of the discussion, a loss of two
22 lines at glare testing would be considered significant.

23 Dr. Bradley?

24 DR. BRADLEY: Yes, just a comment on the means
25 by which one does glare testing and how that interacts with

1 the underlying spatial distribution of the cataract. We
2 might think of two scenarios, one that was just mentioned
3 of a central cataract that is anterior or posterior. The
4 idea of most glare testing is you employ a bright light
5 source and, in so doing, the pupil constricts, and
6 therefore the cataract fills a larger proportion of the
7 pupil, and therefore the scattered light becomes a larger
8 proportion of the retinal image, and that leads to an
9 increased visual effect.

10 Obviously, the converse is true. If you have a
11 peripheral or marginal cataract, as the pupil is
12 constricted, a smaller and smaller proportion of the pupil
13 is covered by the cataract, and therefore a smaller
14 proportion of the retinal image is scattered light, and
15 thus the visual effects are decreased under those
16 circumstances.

17 So I'm not sure there is a single way one can
18 do a glare test that would sensitize the tester to the
19 visual impact of a cataract, and it may be necessary to
20 employ more than one approach. Just an off-the-cuff
21 suggestion would be to employ the one suggested, which is
22 the standard approach, perhaps, where the pupil constricts
23 in the presence of the bright light and it's highly
24 sensitive for picking up the visual impact of a central
25 cataract. I might also suggest performing the same test

1 under cycloplegic pupil dilation to emphasize the impact of
2 a marginal or peripheral cataract that, of course, would be
3 visual manifest for the patient under, for example, night
4 driving circumstances, which arguably are the most
5 important ones.

6 DR. WEISS: The question, Arthur, has anyone
7 done cycloplegic glare testing? Do we have any data to
8 know what the results are in the normals?

9 DR. BRADLEY: I don't know.

10 DR. WEISS: Dr. Mathers?

11 DR. MATHERS: I would strongly agree with what
12 you have suggested. I meant that we should do standard
13 testing, contrast sensitivity testing, and glare testing,
14 and I do not have any experience doing glare testing in a
15 dilated pupil, but it is similar to night driving, but we
16 just don't have any data on that, and I don't think that
17 the visual assessment is going to give us all of the
18 answer. I think we're going to see a lot more with the
19 objective and quantitative than we do simply with the
20 vision changes.

21 DR. WEISS: Well, perhaps if we don't have that
22 data, we could request a subset. If panel thought that was
23 helpful and so did the agency, we could request a subset of
24 patients when they have their initial entry dilated exam to
25 be glare tested while dilated and not dilated, so not to be

1 too onerous to any sponsors.

2 Dr. Grimmett, then Dr. Bullimore.

3 DR. GRIMMETT: Dr. Grimmett. Just a quick
4 comment. If my memory serves me correctly, I think I've
5 seen some dilated glare testing data in an ARVO abstract by
6 Arthur Ginsberg out of California, San Ramon. He does
7 functional driving tests and other activities and has some
8 data on that kind of stuff. He's a contrast sensitivity
9 guru. I think I've seen that data before, so I think it
10 could exist.

11 DR. WEISS: Dr. Bullimore?

12 DR. BULLIMORE: My impression is we're going to
13 come to contrast sensitivity in a minute and glare testing,
14 but since it's on the table, this is a real sticky area,
15 and anybody who a few years ago was involved in the Eye
16 Care Technology Forum knows that agreement was quickly
17 reached on some areas, like measurement of intraocular
18 pressure and visual acuity and visual fields, but contrast
19 sensitivity and glare testing became a thorn in the side of
20 the organizers of that meeting, and Morris Waxler was the
21 point person on that and Arthur was involved in the panel
22 as well. It was very difficult.

23 As far as assessment of cataract, I think one
24 of our speakers put forward a number of mechanisms by which
25 a standardized grading system could be used. I think the

1 panel should consider whether the FDA should strongly
2 recommend a reading center for cataracts or for lens
3 opacities in the same way they're recommending it for
4 endothelial cell density.

5 It becomes difficult because in the case of
6 endothelial cell density, you have a standardized
7 instrument that can capture the image. Photographing lens
8 opacity, since you want to capture the different features
9 of the lens, is a much more difficult and sophisticated
10 procedure.

11 So I think using standardized grading systems
12 is appropriate, paying attention to the kind of opacities
13 that are likely to occur. Anterior and posterior capsular
14 opacities I think are appropriate, and cortical opacities
15 maybe, but I'm not in favor of requiring a reading center
16 in the same way that it's currently recommended for the
17 endothelium.

18 DR. WEISS: If there are no other comments on
19 this section, Ms. Lochner, if we could proceed to the third
20 and last question.

21 MS. LOCHNER: At the most recent ANSI meetings,
22 a consensus appeared to be reached on the general
23 parameters of the contrast sensitivity substudy that's
24 outlined in Section 8.3, and Dr. Bullimore, you might be
25 interested to know that we had Dr. Ginsberg, who's the

1 contrast sensitivity guru, and we had Dr. Jack Holiday,
2 who's advocated contrast acuity testing, actually agreeing
3 on this point. So it was a red letter day.

4 (Laughter.)

5 DR. BULLIMORE: For the record, I never called
6 him a contrast sensitivity guru. That was Dr. Grimmer. I
7 want to strictly go on the record on that.

8 (Laughter.)

9 MS. LOCHNER: And I think that the use of the
10 contrast sensitivity systems, rather than contrast acuity,
11 was recommended because of contrast sensitivity's ability
12 to capture the full range of spatial frequencies and
13 contrasts, and it was felt that the contrast acuity charts
14 would potentially miss significant contrast losses because
15 of the unpredictability of the spatial frequency at which
16 these losses may be seen. Of course, we will have letter
17 recognition performance under low light conditions assessed
18 by the best-corrected visual acuity testing.

19 The contrast sensitivity testing, as proposed,
20 includes mesopic and mesopic with glare conditions. Please
21 comment on the clinically significant decrease being set at
22 .3 log units, and also on whether this decrease should be
23 at one or two or more spatial frequencies to be considered
24 significant.

25 Next, should charts with the minimum contrast

1 at each spatial frequency be used to minimize the problem
2 of missing data, and perhaps first of all, are these charts
3 commercially available?

4 Please also comment on the recommended analyses
5 of these data, including how missing data should be
6 handled, and by missing data, we mean when the patient is
7 unable to see the target at a particular spatial frequency
8 at any of the available contrast levels.

9 Last, please provide any additional comments,
10 particularly any recommendations you may have to improve
11 the quality of the data generated from this testing.

12 Dr. Bullimore?

13 DR. BULLIMORE: Can you go back to the first
14 slide? I want to take these questions one at a time with
15 panel input, if that's okay, Madam Chairman.

16 DR. WEISS: Anything you want.

17 DR. BULLIMORE: Wow. I guess lunch is not on
18 the table.

19 (Laughter.)

20 DR. WEISS: That's true.

21 DR. BULLIMORE: Let me paraphrase my comments
22 by saying that I have a long and distinguished record of
23 being a fan of letter charts over grating, so anything I
24 say should be taken in that context.

25 That notwithstanding, I think first of all, the

1 statement about .3 log units being a clinically significant
2 decrease in contrast sensitivity is reasonable. Just to
3 put that in context, we've come to accept two lines of
4 visual acuity on a logMAR chart as being a meaningful
5 decrease and representing a complication or an adverse
6 event or, to put it more broadly, an unsatisfactory outcome
7 of a refractive procedure. Here, we're talking in the
8 contrast domain of an equivalent of three lines, and I
9 think this is reasonable and conservative.

10 I think saying that the drop should be at two
11 or more spatial frequencies, again, we get into the martial
12 end of contrast sensitivity testing pretty quickly. One of
13 the limitations and reservations that some people have
14 about these tests is that unlike letter testing, for
15 example, where the patient has to name a letter, the
16 patients is asked either whether they can actually see the
17 grating or not or is asked to say is the grating on the top
18 part of the chart or on the bottom part of the chart? So
19 the opportunity for a bias based on shifts in criteria if
20 you use the first approach or the opportunity to sort of
21 guess correctly when it's just a one in two chance compound
22 the analysis of some of these data.

23 But again, in the interests of the goodwill
24 exhibited between Dr. Holiday and Dr. Ginsberg, I think
25 this again is a reasonable, practical approach, and with a

1 letter chart, of course, you might not have the problem
2 with then having to say, well, is it one or two spatial
3 frequencies, but really I think the panel at this stage
4 should be presented with the data when it's available and
5 let the panel decide, so to speak.

6 Does anyone want to comment on that first
7 thing? I'm sure Arthur would. I'd appreciate Dr. Owsley's
8 input anytime she wants to say something.

9 DR. WEISS: Dr. Bradley? Or Owsley. Whoever.

10 DR. OWSLEY: Why don't keep going and then I
11 can probably just make a few comments at the end?

12 DR. BULLIMORE: Okay. Arthur, can I keep going
13 for you as well or do you want to interject?

14 DR. BRADLEY: I really appreciate the
15 opportunity to interject.

16 (Laughter.)

17 DR. WEISS: So why don't you?

18 DR. BRADLEY: Mark raises an interesting sort
19 of comparative analysis to try and decide, well, is .3 log
20 units, obviously a factor of two changing contrast
21 sensitivity, clinically significant? And he draws the
22 parallel between what people decide is clinically
23 significant in terms of logMAR for visual acuity change,
24 and two lines, obviously that's .2 log units on a logMAR
25 chart, and I question, Mark, that that is somehow

1 equivalent to .3 log units in contrast sensitivity. Was
2 the equivalence based upon some sort of Z score change,
3 Mark, or --

4 DR. BULLIMORE: I didn't mean to imply that
5 they were equivalent, but I said it sort of parallels the
6 change.

7 Now, you could say, well, if we're doing .2 log
8 units of visual acuity, we should use .2 log units for
9 contrast sensitivity. Unfortunately, most of the
10 commercially available tests for contrast sensitivity go in
11 steps of .15.

12 Again, we're going to have the data. We're
13 going to be able to look at the number of patients that
14 have lost .3 or more log units at one or two spatial
15 frequencies. We'll have mean contrast sensitivity data for
16 each spatial frequency and for each lighting condition.
17 There will be a colossal amount of data that we'll be able
18 to sort of chew over in depth when the opportunity arises.

19 DR. BRADLEY: Yes, I think we'll have the data,
20 but we'll still be left with the question about what's
21 going to be significant.

22 Just as a suggestion, then, I would make that
23 if there is consensus that a .2 log unit change of visual
24 acuity is clinically significant, would you think it
25 reasonable that we convert that into some sort of acuity Z

1 score change -- say, two, three, four standard deviations,
2 whatever it is -- and propose that an equivalent Z score
3 change in contrast sensitivity be considered significant.
4 Does that make any sense at all?

5 DR. BULLIMORE: It makes sense, but I don't
6 think it's an approach that I would advocate at this stage.

7 DR. WEISS: Dr. Owsley?

8 DR. OWSLEY: This is Cynthia Owsley. I think
9 one of the --

10 (Telephone rings.)

11 DR. OWSLEY: Could somebody get that?

12 DR. BULLIMORE: I think it's Art's.

13 DR. WEISS: It's probably Pizza Hut returning
14 someone's surreptitious call.

15 (Laughter.)

16 DR. OWSLEY: I mean, I think both Mark's
17 approach and Arthur's approach are reasonable approaches.
18 The problem for me, when I think about this, is when you
19 decide how much decline on a visual function test is bad in
20 some sense, clinically significant in a bad sense, you have
21 to ask yourself what is it you're trying to prevent.

22 Answering these questions in vacuo without
23 looking to see how much of a loss you need in contrast
24 sensitivity or acuity or glare or whatever the visual field
25 causes a problem in functional performance, without looking

1 at it in that way, I just don't see how you could answer
2 the question.

3 I'm not suggesting that we propose all kinds of
4 -- whatever they're called -- substudies to answer that,
5 but I think it's an important dilemma, because when it's
6 considered on its own, it's abstract. For the patient,
7 it's in terms of their everyday life. What implications
8 will a .2 loss in logMAR acuity mean? What implications
9 for their everyday will be a .3 loss in contrast
10 sensitivity?

11 So I haven't proposed any answers to this, but
12 I see it as a very sticky dilemma that we might just have
13 to kind of go with something that feels like it has some
14 face validity sort of on a clinical level.

15 DR. WEISS: Mr. McCarley, Dr. Mathers, Dr.
16 Bandeen-Roche, and then Dr. Bradley.

17 MR. MCCARLEY: Yes, just one question. Is this
18 the first time contrast sensitivity testing has been
19 required by the panel? My understanding is that
20 manufacturers of refractive lasers also collected this
21 data. I mean, the question behind that is is there a
22 standard for contrast sensitivity or are we now making the
23 standard?

24 DR. WEISS: Well, we're making a standard for
25 phakic IOLs, so even if it wasn't required for anything

1 else, it really --

2 MR. MCCARLEY: In my understanding, it was.
3 Maybe I'm wrong.

4 DR. WEISS: Well, I'll defer to Mr. Whipple.

5 MR. WHIPPLE: Yes, I believe we have required
6 contrast sensitivity for LASIK studies.

7 DR. BULLIMORE: Yes, and the current guidance
8 document says that either you have to measure it or
9 basically to say that you didn't measure it.

10 MR. WHIPPLE: Right, and I'll defer also if
11 Donna and Malvina --

12 PARTICIPANT: That sounds definitive.

13 DR. BULLIMORE: I'm paraphrasing a little bit,
14 but that's the general spirit of it.

15 DR. WEISS: Donna?

16 MS. LOCHNER: Well, it's not the first time
17 it's been required. As you point out, the LASIK example,
18 and of course the panel has reviewed extensive contrast
19 sensitivity data in the multifocal IOL example.

20 I think it's also important to point out this
21 question about the clinical significance because I think
22 what Dr. Bullimore was saying about, you know, you're going
23 to have to look at the data and make some judgements when
24 it's received, but this question is backing up to the
25 sample size calculations. It's not backing up that if it's

1 .3, the device is approved and if it's .29, it isn't. It's
2 really more is this a reasonable effect to be looking at
3 the sample size calculations, et cetera.

4 I think also, as Dr. Owsley said, the meaning
5 of this is really not going to be clear to many of the
6 panel members, as well as to the patients, unless this were
7 there were this daily living-type testing required, which
8 we really have not -- we're not suggesting that, but I
9 think it's a point well taken. So really, if you think of
10 this in the context of the sample size.

11 DR. WEISS: I think Dr. Bullimore --

12 DR. BULLIMORE: Yes, I mean, you can look at a
13 .3 loss as being clinically significant two ways, whether
14 it occurs in a subset of individual patients or in one
15 individual patient or whether it occurs in the population
16 as a whole, and in terms of sample size, you're interested
17 in the population as a whole. In terms of the safety of a
18 device, you might say, well, if a two-line loss of visual
19 acuity sets off alarm bells and appears on some summary
20 statistics for safety, then should there be an equivalent
21 here?

22 I'm not comfortable doing the later. You know,
23 we heard from a speaker this morning who has a lot more
24 experience with phakic IOLs than the rest of us, and he was
25 advocating that the IOLs be held to the same standard

1 visually as LASIK, and we have some historical precedent
2 here with these tests being done on everybody and I think
3 we should carry on with that.

4 In terms of the sample size for the entire
5 cohort, I think it's much more likely to be driven by
6 endothelial cell density considerations than anything
7 visual.

8 MS. LOCHNER: So it's your recommendation that
9 this testing be performed on all individuals?

10 DR. BULLIMORE: Yes. That would be my
11 recommendation.

12 DR. WEISS: I think we're going to go along
13 this time even if it's out of order I originally said.

14 DR. BULLIMORE: I'm not asking that everybody
15 flies to Iowa for a driving simulation. I mean, this a
16 test that should be reasonably easy to incorporate into a
17 protocol and I would like to see data on as many patients
18 as possible.

19 MR. CALOGERO: Don Calogero. Right now, it's
20 set up as a substudy, and I believe the sample size is
21 somewhere between 60 and 100 patients, and that gives us
22 the ability to detect down to .15 log units. So you're
23 saying that in spite of perhaps having sufficient precision
24 to the study, you'd like to see it on the entire
25 population?

1 DR. BULLIMORE: Well, I think where Dr. Mathers
2 was going with this is that he would like to see data on
3 contrast sensitivity with and without glare that would give
4 you clues, maybe not definitive decisions, as to whether
5 significant cataractogenesis had occurred in these
6 patients.

7 MR. CALOGERO: Okay. So you're also using it
8 for that purpose, then.

9 DR. BULLIMORE: Exactly. Now, if my memory
10 serves me correctly, I mean, yesterday's presentation for
11 an intraocular --

12 DR. WEISS: Wavefront.

13 DR. BULLIMORE: Wavefront. That was it.

14 DR. WEISS: How quickly we forget.

15 DR. BULLIMORE: A good dinner.

16 For that, we had data on all the patients. Am
17 I correct?

18 MR. CALOGERO: I believe so, yes.

19 DR. BULLIMORE: Yes. So I would have thought,
20 with a newer technology, which, let's face it, these phakic
21 IOLs are, and I think increased or elevated safety concerns
22 in terms of lens opacities, I'd want to have that data on
23 as many patients as possible.

24 Now, I realize we may be into latter phases of
25 some PMA investigations, and that's fine, but I think it's

1 going to inform at this stage whether the patients have
2 indeed developed anything that may cause concern in the
3 lens opacity department.

4 MR. CALOGERO: Thanks for clarifying that.

5 DR. WEISS: Dr. Mathers, then Dr. Bradley, and
6 we'll move our way around, and then come back to Dr.
7 Owsley.

8 Dr. Mathers?

9 DR. MATHERS: Yes, I agree with Dr. Owsley that
10 it's difficult to assess how this impacts the patient, but
11 what we need to do is have reasonably high resolution of
12 the data that we're picking up. When we're talking about
13 whether this is actually a good thing, you have to contrast
14 that with the struggles the patients have with enormous
15 myopia and their tremendous problem and the alternative of
16 clear lens extraction and other significant issues, but we
17 want to have reasonable resolution and an ample sample size
18 so that we can tell what's going on, and it may be that two
19 or three lines is a reasonable expectation considering the
20 other struggles that they have. But we would determine
21 that later. We need to have the data now, and I think the
22 contrast sensitivity not only gives us something about the
23 visual function of the system, but also the
24 cataractogenesis process.

25 DR. WEISS: I think that probably most people

1 here are in agreement that we need the data and whether or
2 not someone postulates that they know what the data means
3 upfront versus whether they don't know what the data means
4 upfront really won't affect the FDA. They can just tell
5 you after the fact whether you were right or not.

6 So if anyone has any comments not related to
7 their opinion on that particular topic, we can proceed with
8 them.

9 Dr. Bradley?

10 DR. BRADLEY: Yes, I think just because this is
11 a new product doesn't necessarily mean we have to measure
12 contrast sensitivity. We really need to have sort of
13 underlying theoretical rationale for why contrast
14 sensitivity measurement might or might not be useful in
15 this particular case.

16 I think the primary concern we have here is
17 degradation of optical quality, whether it be in the cornea
18 due to endothelial problems or primarily in the lens due to
19 cataract development. The likely optical and visual
20 consequences of that are related to scattered light, and
21 they have fairly predictable optical effects. Indeed, they
22 will have some effect on our ability to see fine detail,
23 which ought to be revealed by some visual resolution task.

24 In addition to that, they will have impact on
25 the image quality for larger targets, and the primary

1 impact will be on contrast reduction. One method for
2 assessing this is to examine contrast sensitivity.

3 That said, as Mark alluded earlier, there are
4 devotees of letter contrast sensitivity testing and there
5 are devotees of grating contrast sensitivity testing. One
6 of the arguments in favor of grating testing is that one
7 can test at many spatial frequencies, which academically
8 might be quite interesting, but unless one can come up with
9 a reasonable theoretical argument for why a measurement at
10 a single low spatial frequency might not provide the same
11 information, I think one is wasting one's time measuring at
12 lots of different spatial frequencies.

13 Therefore, I wonder about most grating tests
14 for that reason, and I think Dr. Owsley might be able to
15 comment on that.

16 DR. OWSLEY: Yes. I very much agree with the
17 perspective or the question you're raising. I would not
18 describe myself as a devotee or a guru of any test.

19 The reason I favor letter tests in situations
20 like this is that I know of no convincing evidence in the
21 peer-reviewed literature that shows that letter tests are
22 worse than grating tests or grating tests are better than
23 letter tests, however you want to say it.

24 I think that there are a lot of us who do
25 clinical studies, clinical intervention evaluation studies,

1 including bodies like the National Eye Institute, who go
2 the way of letter tests not only based on the evidence that
3 they're not worse than the grating tests, but because of
4 two things. One, if you're getting a measure of high-
5 resolution visual acuity, and then you do a contrast
6 sensitivity letter test, say like the Pelly-Robson, which
7 has large letters, you're basically getting information
8 about the entire shape of the function.

9 Then the second reason is that we're talking
10 about -- I don't know which specific test the sponsor would
11 be using or any sponsor would be using, but if you're doing
12 spatial frequency testing in the kinds of things we're
13 talking about in this kind of study, you're doing those
14 spatial frequencies for topically, mesopically, with glare
15 and without glare, before and after surgery. This is a
16 long testing time, okay?

17 And then I think -- what was my last point on
18 this? Well, I'll just leave it like that, but I know that
19 this issue has been visited by this panel before and I've
20 heard about it. I've never been to the panel to hear it
21 argued as a panel member, but I think that it's an
22 important point that sponsors should hear, the public
23 should hear, and the FDA should hear that there really are
24 really sound arguments for not, in every case, doing all
25 the spatial frequencies. You need to ask yourself, what

1 are you doing with that information? What's it really
2 providing for you that high-contrast acuity in a single
3 measure on a large letter contrast sensitivity chart would
4 not provide?

5 DR. WEISS: Dr. Swanson?

6 DR. SWANSON: I agree with both sets of
7 comments and wanted to add one other thing. I've used
8 letters and spatial frequencies. It depends on what the
9 task is.

10 The point that Dr. Bullimore raised becomes
11 important for sample size calculations as well as for the
12 amount of time that goes on, particularly if we're looking
13 at this question of glare, but as Dr. Bradley mentioned, in
14 general, these effects should be in a range of spatial
15 frequencies.

16 If you are using a grating and you have a two-
17 alternative forced choice, basic researchers who use those
18 things do lots of trials because that's the only way you
19 can reduce test/retest variability. For letter testing,
20 where there's going to be between 10 and 26 options that
21 the person's guessing amongst, a much smaller number of
22 trials is needed.

23 So in order to have a significant change, you
24 need to have a test where the test/retest variability is
25 low. For a small number of trials, which needed to run all

1 these conditions in all these people, a multiple
2 alternative, 10- or 20-alternative forced choice test, is
3 going to be much better than a two-alternative forced
4 choice, and there aren't any commercially available 10-
5 alternative forced choice grating tests and it's going to
6 be very difficult to create on given the response demands
7 on the patients. So a letter type of acuity test is going
8 to be much more suitable in terms of gathering a fair
9 amount of useful data in a short period of time.

10 Then the questions that come up, such as, well,
11 do we have to have two spatial frequencies down or three
12 spatial frequencies or one spatial frequencies down, those
13 become quite complicated by the high test/retest
14 variability of a two-alternative forced choice with just a
15 single endpoint.

16 DR. WEISS: So let me get to a bottom line.
17 What do you want?

18 DR. SWANSON: I'm just trying to hammer home
19 all the points they made about letter charts being superior
20 for this purpose. I understand there was agreement before.

21 DR. WEISS: So you would like a letter chart
22 for this purpose.

23 DR. SWANSON: For the purpose, a letter chart
24 is going to allow a much smaller sample size and a much
25 larger --

1 DR. WEISS: And Dr. Owsley, you agree, and Dr.
2 Bradley?

3 DR. BRADLEY: Yes, I agree, a letter chart. I
4 think Donna asked us if such charts were available. There
5 certainly is a letter contrast sensitivity chart, as I
6 think it's been referred to. Sometimes it's called a
7 Pelly-Robson chart.

8 One thing I would alert the FDA to is if you
9 are to request use of that chart, it doesn't have to be
10 used at the standard distance at which it was originally
11 designed, and there are reasons to pick your distance,
12 depending upon the scale or size of target for which you
13 wish to study.

14 DR. WEISS: Dr. Bullimore, your comment?

15 DR. BULLIMORE: I mean, the spirit in sort of
16 guidance documents before, whether for this issue or
17 others, was that the sponsor or a potential sponsor would
18 be able to speak with the FDA and the FDA, in terms of the
19 guidance document, wouldn't specify a given test, but there
20 would be some scope. I was just a little surprised that
21 this was so specific, even naming spatial frequencies. I
22 mean, you come pretty close to naming a test or naming a
23 couple of tests, and I'd just like to see some statement
24 that other tests could be considered if the sponsor in
25 consultation with the agency would be able to agree that

1 these were acceptable.

2 DR. WEISS: Don Calogero?

3 MR. CALOGERO: I can provide a little
4 background. At the last ANSI meeting, we sort of discussed
5 this issue, and there was a presentation and the basic
6 summary of the presentation was that it's really not
7 possible for us to predict where the largest effect is
8 going to be in terms of the spatial frequencies. We went
9 through the literature containing some information --

10 DR. BULLIMORE: But Don, the panel's telling
11 you that that's not an appropriate -- well, not necessarily
12 a widely held view of the world.

13 MR. CALOGERO: Well, this was from the
14 literature and this is from sort of the internal data we
15 have. Some devices, depending on the type, had the largest
16 effect at 1.5 cycles per degree, others 12 cycles per
17 degree, and the committee felt that if we had recommended
18 the letter charts, based on sort of the spatial frequency
19 domain that they evaluate, you're essentially just simply
20 getting the highest spatial frequencies, maybe 12 cycles
21 and above, whereas we could potentially miss large drops at
22 other spatial frequencies, and with our current to predict
23 where the effect is, we felt it would be judicious to
24 recommend that the entire contrast sensitivity function be
25 defined.

1 DR. OWSLEY: Can I make a comment?

2 DR. WEISS: Yes. Make a comment, Cynthia, and
3 then also if we could direct it at, from what I understand
4 from the vision scientists on this panel, they don't agree,
5 and is that the case or is that not the case?

6 DR. OWSLEY: The vision scientists on this
7 panel I think do agree the letters would be better.

8 DR. WEISS: No, you agree that letters would be
9 better, but are you agreeing with what's been put forward
10 by the FDA?

11 DR. OWSLEY: Well, let me put it like this.

12 DR. WEISS: That's what I'm asking.

13 DR. OWSLEY: And that's exactly what my comment
14 is on. I've been reading this literature for 20, 25 years,
15 like Arthur -- Mark's a little younger -- and Steve.

16 DR. WEISS: He's thinking a lot younger.

17 DR. OWSLEY: I would be happy as a consultant
18 to the panel to look at those peer-reviewed papers. Of
19 course, if you have internal data you can or cannot share
20 with me, that's your decision, but I have not seen any
21 convincing evidence that measuring all these different
22 spatial frequencies would matter in any of the decisions
23 that you would be faced with at the FDA regarding these
24 types of devices. But I'm openminded. I just haven't seen
25 those papers.

1 DR. WEISS: So I see sort of agreement by Dr.
2 Swanson.

3 Dr. Burns has a comment.

4 DR. BULLIMORE: I mean, my overall view is that
5 these --

6 DR. WEISS: Dr. Burns has a comment.

7 (Laughter.)

8 DR. BULLIMORE: This --

9 DR. WEISS: Well, I guess Dr. Bullimore has a
10 comment.

11 (Laughter.)

12 DR. BULLIMORE: Since it's my --

13 DR. WEISS: Dr. Burns is deferring to you,
14 Mark, so why don't you give your comment?

15 DR. BULLIMORE: I have no problem with gratings
16 being used. It was just having a sentence in there that
17 other avenues could be pursued if --

18 DR. WEISS: So you don't like the rigidity of
19 it, but you don't have a --

20 DR. BULLIMORE: It's more the rigidity, and I
21 mean, as I said before, you're almost sort of saying you've
22 got to use this test or this test, and those of us who work
23 in other arenas, and we've already spoken about the
24 limitations of certain types of tests, we just find that a
25 little offensive, I think.

1 DR. WEISS: Dr. Burns, and then Karen Bandeen-
2 Roche and Mr. McCarley.

3 DR. BURNS: Yes, I just want to chime in that I
4 also believe that the combination of the high contrast and
5 the low contrast letter chart will give us enough
6 sensitivity both for some of the safety issues, such as
7 contrast degradation from cataracts combined with glare,
8 and remember also, this is a refractive procedure. So we
9 do want to get general decrement information. But I do
10 believe that a contrast letter chart combined with a high
11 contrast letter chart will be acceptable.

12 DR. WEISS: Mr. Whipple?

13 MR. WHIPPLE: Yes, Dave Whipple. I just wanted
14 to address Mark's comment about the rigidity of the
15 guidance. We may recommend certain tests or prefer certain
16 things, but a guidance is exactly that. It's guidance. It
17 always carries with it the option of using other tests,
18 other test methods, and making arguments why the
19 recommended tests aren't appropriate for that particular
20 device. So that's inherently built into the guidance
21 document development.

22 DR. WEISS: We'll go on to the other two
23 comments if they're not related or they don't change what
24 my next statement is going to be. It's from why I
25 understand, the members of the panel prefer a letter chart

1 and Dr. Bullimore would also prefer that the wording
2 doesn't sound so restrictive, even if it actually isn't.

3 Having said that, do you have anything else to
4 add to that, Dr. Bandeen-Roche or Mr. McCarley?

5 DR. BANDEEN-ROCHE: My comment was a brief one
6 on a different topic, so should we finish this?

7 DR. WEISS: If it's related to this question --

8 DR. BANDEEN-ROCHE: It's related to clinical
9 significance in terms of functional performance.

10 DR. WEISS: Well, why don't you just proceed?

11 DR. BANDEEN-ROCHE: I just wanted to mention
12 that there is detailed data from the Salisbury Eye
13 Evaluation Project on glare, contrast, et cetera, with many
14 measures of functional performance. So I don't think we're
15 completely in the dark. I think that there are some data
16 that can inform that question, albeit it is a sample of
17 older adults.

18 DR. WEISS: Good. Thank you.

19 Mr. McCarley?

20 MR. MCCARLEY: Yes, just very quickly. I'd
21 like to understand or get clarified for me the intent. My
22 understanding, at least from the ANSI side of it when we've
23 been discussing this for the last several years, was that
24 doing contrast sensitivity was to see if the combination of
25 optical components in the eye degraded, similar to what the

1 intent was when they did LASIK studies. Now, it seems to
2 be that's what's being used to determine whether or not the
3 patient's undergoing a degradation as a result of cataract,
4 for instance. So one would be a sample size and one would
5 be all patients, I think.

6 DR. WEISS: Dr. Bullimore, did you want to
7 address that?

8 DR. BULLIMORE: No.

9 DR. WEISS: No one wants to address that. I
10 guess we have no takers?

11 DR. BURNS: If you're measuring a visual
12 performance like this, you're tapping the whole system, and
13 the degradation can come about from optical imperfections
14 or from what are actually just very high-order optical
15 imperfections of scattering and tissue turbidity. So it
16 probes the whole system.

17 DR. WEISS: I think we've beat this one into
18 the ground and it's no longer even gasping. So why don't
19 we go on to the next one?

20 DR. BULLIMORE: The next one I think is very
21 easy to deal with.

22 DR. WEISS: We'll hold you to that.

23 DR. BULLIMORE: Again, we're asking people to
24 be very specific about the tests, and I think the FDA have
25 enough savvy and flexibility to deal with this, but if

1 somebody can't read or can't see the highest contrast on
2 the chart, I think their contrast sensitivity has to be
3 scored as zero. One would hope that preoperatively the
4 test was sufficiently intelligently designed such that they
5 could read well above that bottom level.

6 Certainly, recalling some of the data we were
7 presented with yesterday, we were dealing with contrast
8 sensitivities of one or above most spatial frequencies. So
9 as long as we have something down at the 40 percent or 60
10 percent, which of course is .3 or so contrast sensitivity,
11 we should be okay.

12 But if they can't see it, it shouldn't be
13 counted as missing data. We should assume that their
14 contrast sensitivity is zero or some other intelligent
15 value based on the range of the test.

16 DR. WEISS: So can the transcript reflect there
17 are about six heads bobbing up and down?

18 Next question.

19 DR. BULLIMORE: And the next question, please,
20 Donna?

21 I think this is kind of tied in with the last
22 one. If I can answer, I think tests are available and if a
23 patient can't see a targeted spatial frequency at any
24 contrast, it should be scored as a contrast sensitivity of
25 zero.

1 MR. WHIPPLE: That is information.

2 DR. WEISS: Dr. Bradley, and then Dr. Swanson.

3 DR. BRADLEY: I'm having a bit of trouble with
4 the question. I mean, we're not anticipating these
5 patients are going to have horrible vision, are we?

6 DR. WEISS: Before or after the implant?

7 (Laughter.)

8 DR. WEISS: Sorry.

9 DR. BRADLEY: I mean, we're talking about some
10 serious methodology here. We're trying to ascertain
11 something fairly subtle, and I don't think in any of these
12 cases we're going to have the problem that people can't see
13 the target at all.

14 MR. CALOGERO: I think why we're asking the
15 question is under the mesopic test conditions. Under the
16 mesopic test conditions, we do in fact have cases where at
17 the higher spatial frequencies, the patients are unable to
18 perform anything, and so you might have 20 or 30 percent of
19 the population that you're testing that essentially has a
20 zero, let's say, 12 or 18 cycles per degree. So we wanted
21 some sort of standardized way of handling all of those zero
22 data points.

23 DR. BRADLEY: Well, I think the panel has given
24 you our suggestion on that, which is that you should not do
25 these grating tests. Particularly, don't try to find out

1 if people can see fine detail in the dark. We know they
2 can't.

3 DR. OWSLEY: This is Cynthia Owsley. It's very
4 unusual to find a patient who can't see the first triplet
5 of the Pelly-Robson chart, and I see Karen Bandeen-Roche
6 from the SEE Project nodding.

7 DR. WEISS: But in that case, I think we have
8 the answer to this question. We may not like the question,
9 but I think Dr. Bullimore has already given an answer to
10 it, which I assume you would agree to, but you don't think
11 it's going to even come up.

12 Fine. Next?

13 DR. BULLIMORE: Beaten to death. Beaten to
14 death.

15 I really have nothing to add on this one at the
16 moment. I'll read it again. "Please provide any
17 additional comments on the contrast sensitivity substudy"
18 -- now it's a substudy -- "including any other guidance
19 that could be provided to enhance the quality of the data
20 that are generated from this testing."

21 Use a good test would be my recommendation.

22 (Laughter.)

23 DR. WEISS: That's why we call the experts to
24 figure this out.

25 DR. BULLIMORE: And use it well.

1 DR. WEISS: So is that the end of your
2 questions?

3 DR. BULLIMORE: I'd like to hear what Dr.
4 Owsley has to say.

5 DR. WEISS: Dr. Owsley?

6 DR. OWSLEY: I think the test should have the
7 best test/retest reliability as possible and should be
8 brief.

9 DR. WEISS: And brief is always dear to my
10 heart.

11 Yes?

12 DR. BULLIMORE: This is Dr. Bullimore again,
13 not wanting to shut up. I think it would be ideal and
14 preferable and maybe mandatory that for any test that's
15 going to be used that normal data are available on a wide
16 range of age ranges and for the measured conditions.

17 DR. WEISS: Dr. Bradley?

18 DR. BRADLEY: I think just to reiterate what
19 everybody has said here and maybe to put the comment that
20 Bill Swanson made earlier, what we're arguing for regarding
21 the quality of a contrast sensitivity test can be thought
22 of by thinking of an analogy with a visual acuity test.

23 Imagine you had a visual acuity chart, one
24 letter per line, and it was either A or B. Would we
25 consider that an appropriate test? And clearly not, and

1 all we're saying is let's hold the contrast sensitivity
2 test to a similar standard psychophysically in terms of
3 rigor and test/retest reliability that we are quite
4 comfortable demanding from our visual acuity tests, and it
5 turns out at the moment the only convenient one out there
6 that we know of is a letter chart.

7 DR. WEISS: Dr. Bradley, we get the idea you
8 like the letter chart.

9 DR. BULLIMORE: He is so proud that he can
10 read.

11 DR. WEISS: He likes that letter chart.

12 DR. BULLIMORE: He is so proud.

13 (Laughter.)

14 DR. WEISS: If there are no other comments on
15 this, what I'd like to do then is move on to panel
16 discussion on remaining issues, and I have three issues,
17 and then we'll introduce if the agency has any additional
18 issues or if any members of the panel want to introduce any
19 additional issues.

20 One issue, which has been alluded to and
21 discussed, but I would like to see if we can give as many
22 opinions as possible, is duration of study before it's
23 presented to the panel. Somewhere between two and three
24 years has been mentioned. Perhaps any of those of you who
25 would like to give your opinions can raise your hand and

1 just give me a number and tell me why.

2 Dr. Burns?

3 DR. BURNS: From what I see of the endothelial
4 cell measurements and their likelihood of accuracy in this
5 population, I think three years is a minimum.

6 DR. WEISS: So Dr. Burns feels three years.

7 Dr. Bullimore, had you -- three years. We have
8 a sign. We're going to signs now. Okay. So we have Dr.
9 Matoba at three, Dr. Grimmett at three. Dr. Bradley likes
10 the number chart.

11 (Laughter.)

12 DR. WEISS: Dr. Huang, three. Dr. Mathers is
13 three. Dr. Owsley is three. Dr. Swanson is three. So I
14 think we have actually, even though we're not looking for a
15 vote, essentially it's almost unanimous, if not unanimous,
16 for three years.

17 DR. BANDEEN-ROCHE: Madam Chair?

18 DR. WEISS: Yes?

19 DR. BANDEEN-ROCHE: May I make a very brief
20 comment?

21 DR. WEISS: Yes.

22 DR. BANDEEN-ROCHE: And that is that I believe
23 it's important, based on Dr. Grimmett's comments, that at
24 least one, preferably two, evaluations be scheduled between
25 two and three years because he suggested that linearity was

1 not beginning to be achieved until then.

2 DR. WEISS: Okay. So one to two endothelial
3 cell counts should be performed between the two- and three-
4 year mark.

5 Dr. Grimmett?

6 DR. GRIMMETT: Dr. Grimmett. I couldn't agree
7 more with Dr. Bandeen-Roche that the data that is present
8 in the peer-reviewed literature do not show a linear
9 approach. We don't have enough data long enough. So that
10 is correct. In an intermediate point.

11 DR. WEISS: Dr. Mathers?

12 DR. MATHERS: I would not end the data
13 collection necessarily with three years because I think
14 that ongoing data could be important.

15 DR. WEISS: That's an important question in
16 terms of a postmarket study. It could be brought to panel
17 at three years, at least from the recommendations from this
18 panel today. After that, a postmarket study of what
19 duration of time would you like to see, Dr. Mathers?

20 DR. MATHERS: Another two years.

21 DR. WEISS: Another two years.

22 Dr. Grimmett?

23 DR. GRIMMETT: Dr. Grimmett. I would agree
24 with Dr. Mathers' sentiments, but that decision would be
25 made more accurately with the data in hand with the PMA at

1 three years obviously, but I would feel certainly
2 comfortable with probably a postmarket study.

3 DR. MATHERS: Yes.

4 DR. WEISS: Dr. Bullimore?

5 DR. BULLIMORE: This is endothelium only?

6 DR. WEISS: It can be whatever you want.

7 DR. MATHERS: Cataract.

8 DR. WEISS: And you don't have to define it
9 now. Basically, what the agency would like is your
10 opinions.

11 DR. BULLIMORE: My opinions will be more well-
12 informed once I have some data, but I have considerable
13 concerns about the long-term safety of these intraocular
14 devices. As the lens continues to grow, as it gets to
15 cohabit with the phakic IOL for a long period of time, I
16 don't think any of us can predict with any accuracy or
17 precision what that's going to do to the crystalline lens
18 on a long-term basis.

19 You know, these devices we can expect to be
20 better than some of the first-generation phakic intraocular
21 lenses, but they haven't been around long enough for us to
22 tell, so I could foresee a substantial -- i.e., five-year
23 -- postmarket study to track these patients and see, for
24 example, how many people develop significant lens opacity,
25 whatever that may be, how many people have to have a

1 cataract extraction and conventional IOL implantation, how
2 many of those patients have complications from that
3 procedure, and it can be nothing more than fairly
4 straightforward head counting, but I think tracking these
5 patients long-term, given what we know happened with, for
6 example, extended wear contact lenses in the '80s and
7 anterior chamber intraocular lenses, I think it would be
8 prudent.

9 DR. WEISS: So the impression I get -- and if
10 there's something you would like to add, then we'll add it
11 -- is that the panel members would like to see postmarket
12 studies from two years on up, maybe even to five years, but
13 pending what the data shows.

14 DR. BULLIMORE: Once we have a PMA before us,
15 we can make a more informed decision, but certainly, given
16 the fact that this is an intraocular device that's
17 cohabiting with a lens that's continuing to grow and
18 function, I would --

19 DR. WEISS: You want to see a postmarket study.

20 DR. BULLIMORE: Well, I would anticipate the
21 possibility of this requirement.

22 DR. WEISS: Dr. Mathers?

23 DR. MATHERS: We have not focused on pigment
24 dispersion, but in the literature there are a number of
25 papers that do. Not so much that they've had a lot of

1 glaucoma yet, but they often have pigment dispersion in the
2 range of 30 percent, and therefore gonioscopy and pressure
3 measurement.

4 DR. WEISS: Actually, that's an excellent
5 point. I don't recall if that's in the document,
6 gonioscopy, and then how often --

7 DR. MATHERS: It is.

8 DR. WEISS: And how often should postop.

9 DR. MATHERS: It's in there now and it hasn't
10 been focused on, but it's in the literature, and I think
11 it's an important issue.

12 DR. WEISS: So we haven't defined here, but
13 there is concern in terms of follow-up as far as not only
14 the cataract and the endothelium, but also gonioscopy and
15 how frequently that would be done and what the postmarket
16 would be would be defined in terms of when some of the data
17 comes in.

18 Another question is what is an independent
19 entity? Is the subject an independent entity or is the eye
20 an independent entity?

21 For example, classically, for intraocular
22 lenses, each patient has been an independent entity. Each
23 subject has been an independent entity. So if the FDA
24 statistically said that 300 data points were required, then
25 each patient would be considered separately an entity

1 whether or not they had one or both eyes treated.

2 On the other hand, for LASIK, each eye has been
3 considered an independent entity in the same patient, so
4 that if 300 data points were required, then you might only
5 need 150 patients.

6 So this is a very important point for sponsors
7 to know whether they need 300 patients, and we're going to
8 be getting to sample size as well, but if the sample size
9 is determined to be necessary to be 300, should that be 300
10 patients or can it be 150 patients with both eyes treated?

11 Dr. Bandeen-Roche?

12 DR. BANDEEN-ROCHE: Yes, I mean, statistically
13 speaking, from a strict point of view, sites are the
14 independent entities. Now, obviously, we're not going to
15 base sample size on that, but methods are available to
16 account for correlation of eyes within individuals, and in
17 turn, of individuals within sites if necessary in computing
18 power.

19 So one doesn't necessarily need to go to the
20 extreme of saying, for instance, we're going to count
21 people as entity. One can still account for the
22 information provided by two eyes within an individual, yet
23 accounting for the correlation in that information that two
24 eyes shouldn't count as much as two separate people, and
25 indeed that perhaps two people within a site shouldn't

1 count as much as two people in different sites.

2 I actually have a figure that I brought on
3 this, but if you prefer to save it, we can.

4 DR. WEISS: No, sure.

5 DR. BANDEEN-ROCHE: I need the overhead. So
6 can I just quickly set that up?

7 DR. BANDEEN-ROCHE: Sure. Maybe while you're
8 setting that up, we can go on to a comment by Dr. Burns,
9 and then to Dr. Matoba.

10 DR. BURNS: Yes, I just want to support that
11 strongly because we're talking here a lot about biological
12 reactions to a foreign body. So there will be a
13 significant amount of correlation, I assume, between eyes,
14 and so the studies should explain how they're going to
15 handle it statistically in the design, and not in post-hoc
16 analyses.

17 DR. WEISS: So I want to understand. Would you
18 feel that an eye should be an independent entity or a
19 subject?

20 DR. BURNS: I think Karen is going to explain
21 the statistical way to handle it, but it's important that
22 it be handled that way from the outset, that you're
23 accounting for the correlation between eyes.

24 DR. WEISS: Now, we're going to go along with
25 other comments, but we also need to look at this not only

1 statistically, but from a safety standpoint because this is
2 a new device.

3 Dr. Matoba?

4 DR. MATOBA: Yes, in regard to the safety
5 standpoint, since we've already decided that because of the
6 issues of quality of life we're going to allow patients to
7 have both their eyes done instead of limiting it to one eye
8 over a three-year period of time, we might as well get data
9 from both eyes, and so fewer people would have to be
10 involved in the initial study.

11 DR. WEISS: Dr. Mathers?

12 DR. MATHERS: Yes, I strongly agree with that
13 comment. I think we should have both eyes, but we should
14 account for the correlation, and that allows us to get over
15 this quality of life issue, which is going to be very, very
16 important to the patient because only doing one eye is not
17 going to work.

18 DR. WEISS: Well, we're not talking about only
19 doing one eye. We're talking about whether you need 300
20 subjects who can both have both eyes done or 150 subjects
21 who can still both have both eyes done.

22 DR. MATHERS: I think we shouldn't throw away
23 the data on the second eye, but we should use it with an
24 appropriate correlation.

25 DR. WEISS: So you basically want to count it

1 and not require another 150 patients.

2 DR. MATHERS: We're about to hear how we can
3 account for the correlation.

4 DR. WEISS: Okay.

5 DR. MATHERS: In estimate.

6 DR. BANDEEN-ROCHE: Yes, and so certainly one
7 can just state what is the expected design, including eyes
8 within patients, patients within sites, et cetera, and then
9 do power and all other statistical analyses accounting for
10 that correlation. So what I had diagrammed here was
11 thinking of people within sites, but it could equally well
12 apply to eyes within people. There might be several
13 levels.

14 So here, M is referring to the larger level.
15 So let me just stick to people within sites because that's
16 the way that I did it. So M is the number of sites and N
17 is the number of individuals within sites.

18 So for instance, if it were people and eyes, it
19 would be M people, and for two eyes, N would be two eyes
20 per person, and so what I'm showing you here is the
21 standard deviation of the mean -- say, a rate, say, a
22 safety rate -- or I guess to simplify things, this was
23 meant to be a continuous measure. So it would be something
24 more like an acuity measure that's measured on a continuous
25 scale.

1 So you can see that the standard deviation
2 depends both on the overall number of units, either people
3 times eyes or sites times people, and the symbol rho, which
4 I'll point to here, that's the correlation between outcomes
5 within a person -- so two eyes within a people -- or
6 between visual acuities within a site, and then sigma is
7 just the standard deviation of the measurement that's being
8 taken in the population.

9 So one thing to point out is that if N equals
10 1, then this top term goes away, and you just get the usual
11 standard deviation of the mean.

12 So what I'm showing you here is how the
13 standard deviation varies for the -- I assumed the FDA
14 sample size of 300, but what's varying here is the number
15 of clusters, and so here's where this is really more
16 relevant to sites and people within sites, because the FDA
17 guidance for the number of sites says something like at
18 least 20 per site and no more than a quarter being
19 accounted for by any one site, and so if we had four sites,
20 then that would just be meeting no more than a quarter by
21 any one site. Then here, up at 300, you only have one
22 patient per site, so that's the outer limit.

23 So the bottom line here, the straight line, is
24 what you get if there's no correlation. That's just the
25 standard deviation of the mean. If there's no correlation,