

R32

TEST FOR MUTAGENICITYin BACTERIA STRAINS in the ABSENCE and PRESENCE
OF A LIVER PREPARATION

The substance ..M.03...77...(=.OCTOPIROX).....
has been investigated with the mutagenicity test originally
described by Bruce N. AMES (see references 1 - 7).
A lot of organic chemicals causing cancer are reactive
chemical species or they must be converted to such by
metabolic processes. They will readily react with nucleophile
If reactive metabolites are generated within mammalian cells
they will attack among other macromolecules, the nucleic acid
Such reactions may form the basic steps in cancer induction.
For production of reactive chemical species can be detected
biologically by studying the mutations, which are induced
in bacteria consequent on the chemicals reaction with the
hereditary material (DNA). Therefore bacteria are used as
Salmonella typhimurium strains which have been specifically
developed for detecting reactive metabolites. For their
growth these strains require the addition of histidine. If
they are treated with agents reacting with the bacterial DNA
a certain number of bacteria in the population will be mutated
so that they no longer require the amino acid for growth.
Such mutants will grow (on an agar plate) as single colonies
on a background lawn of non-mutated bacteria. The number of
colonies can be counted, agents which cause mutation will
increase the number of colonies per agar plate. In any
bacterial population there will be a certain number of
spontaneous mutants and this will differ from strain to
strain. For TA 98⁺ and TA 100⁺⁺ the spontaneous mutation
frequency is usually higher than for TA 1535⁺⁺ and TA 1537⁺.

⁺ (strain with frameshift mutation)

⁺⁺(strain with base-pair substitution)

The first two strains are more sensitive to mutagens than the latter.

Different bacterial strains are used because each strain responds to only a particular class of chemical. A chemical which is mutagenic for one strain is not necessarily mutagenic for the other three. For better understanding the genotype of the tester strains is demonstrated in the following table.

Genotype of the TA strains used for Mutagen Testing⁺

Histidine mutation		Additional mutations		
hisG46	hisD3052	LPS	Repair	R fact..
TA1535	TA1537	rfa	Δ uvrB	-
TA100	TA98	rfa	Δ uvrB	+R

⁺All strains were originally derived from *S. typhimurium* LT₂. The deletion (Δ) through uvr B also includes the nitrate reductase (chl) and biotin (bio) genes.

MATERIAL AND METHODS

The methods used are published in detail by AMES et al. (see reference 2). In brief the procedure is the following one:

The tester strains TA 98, 100, 1535, and 1537 are maintained at -60° C and recultivated for the test in fresh nutrition brot. The S-9 fraction is prepared from the microsomes of a liver homogenate as described by AMES also. Sprague-Dawley rats are injected once with 500 mg/kg of polychlorinated biphenyl (PCB) = Aroclor 1254 (Monsanto, St. Louis, Mo. USA) in corn oil (200 mg/ml) i.p. five days before killing the animals.

The livers are prepared in cooled sterile glass ware using cooled solutions by homogenizing the livers in a Potter Elvehjem Apparatus. This homogenate is centrifuged by 9000 g for 10 minutes (Sorvall RC2-B, in SS-34 head 8700 rev/min.). 1 ml of the supernatant contains microsomes of around 25C mg of liver. This preparation is the S-9 fraction which is stored in deep frozen stage.

Samples of it are thawed immediately before each test and the S-9 mix is prepared. It contains per ml:

S-9 fraction	(0.04 - 0.1 ml)
MgCl ₂	(8 μmoles)
KCl	(33 μmoles)
Glucose-6-phosphat	(5 μmoles)
NADP	(4 μmoles)
Natriumbosphat, pH 7,4	(100 μmoles)

The exact dosage of every charge of S-9 fraction is evaluated in pretests with known negative and positive substances. The tester strains are checked for their reactivity also. Sterility controls are also routinely performed in each experiment. All dilutions of the substances in test are prepared with Dimethylsulfoxid (Schwarz-Mann, spectrophotometric grade), or with different dilutants such as ethanol or Aq. dest. This is reported in each special case in the protocols.

The dosages range from 0.2 to 500 γ, and the tests are performed without and with activation with S-9. The bacterial plates are incubated for 48 hours at 37° C and then the numbers of reverted bacterial colonies counted. Control plates for determining the spontaneous mutation rate are investigated in the same way.

The results are listed in the attached figures.

M 03 - 77

.....

tested with 4 bacterial strains showed No mutagenicity
neither in the presence nor absence of rat liver preparation
in the dosage range of 500 y to 0,2 y.

Frankfurt/Main, 22.8.1977

KREBSFORSCHUNGSLABOR

(Prof. Dr. Gericke)

ABT. FÜR CHEMOTHERAPIE

Wagner
(Prof. W.-H. Wagner)

y = γ = 10^{-6} g/ml = ppm = μ g/ml

REFERENCES

1. AMES, B.N., DURSTON, W.E., YAMASAKI, E. and LEE, F.D.: Carcinogen are mutagens: a simple test system combining liver homogenates for activation and bacteria for detection. Proc.Nat.Acad.Sci. USA 70:2281-2285 (1973)
2. AMES, B.N., McCANN, J. and YAMASAKI, E.: Methods for detecting carcinogens and mutagens with the Salmonella/Mammalian Microsome Mutagenicity Test. Mutation Res. 31:347-364 (1975)
3. McCANN, J. and AMES, B.N.: A simple method for detecting environmental carcinogens as mutagens. Ann. N.Y. Acad.Sci. 269:21-25 (1975)
4. McCANN, J., CHOI, E., YAMASAKI, E. and AMES, B.N.: Detection of carcinogens as mutagens in the Salmonella/microsome test: Assay of 300 chemicals. Proc.Nat.Acad.Sci. USA 72:5135-5139 (1975)
5. McCANN, J. and AMES, B.N.: Detection of carcinogens as mutagens in the Salmonella/microsome test: Assay of 300 chemicals: discussion. Proc.Nat.Acad.Sci. USA 73:950-954 (1976)
6. CAMUS, A., BERTRAM, B., KRÜGER, F.W., MALAVEILLE, C. and BARTSCH, H.: Mutagenicity of S-oxidized N,N-Di-n-Propylnitrosamine Derivatives in *S. typhimurium* Mediated by Rat and Hamster Tissues. Z. Krebsforsch. 86:293-302 (1976)
7. HERVOLD, S. and BUSELMAIER, W.: Induction of point mutations by different chemical mechanisms in the liver microsomal assay. Mutation Res. 40:73-84 (1976)

20.7.1977

52/77

AMES-Test..... Protocol-No.

Preparation in test..... M 03 - 77.....

Microsome fraction S9..... Rat C..... Induction by Aroclor 1254...

Tester strains: mutants of S.typhimurium TA 98
Dilutant:.....^{H2O}.....

XXXXQQ

XXXXSS

XXXXSS

Time of incubation.... 48 hrs... 37° C

Tester strains TA 98	Preparation in test mutants/petri dish 0.1 ml	Spontaneous mutants petri dish 0.1 ml	Evaluation
500 y	Ø - Ø	12 - 19 - 19	
500 y + S9	Ø - Ø	14 - 14 - 16	
250 y	Ø - Ø	12 - 19 - 19	
250 y + S9	Ø - Ø	14 - 14 - 16	
100 y	4 - 5	12 - 19 - 19	
100 y + S9	10 - 11	14 - 14 - 16	
20 y	5 - 8	12 - 19 - 19	
20 y + S9	13 - 15	14 - 14 - 16	
2 y	10 - 11	12 - 19 - 19	
2 y + S9	15 - 15	14 - 14 - 16	
Ø = tested, but no bacterial growth			

20.7.1977

52/77

AMES-Test..... Protocol-No.....

M 03 - 77

Preparation in test.....

Microsome fraction S9..... Rat¹..... Induction by Aroclor 1254....Tester strains: mutants of *S.typhimurium* ~~XXXXXX~~

TA 100

Dilutant:..... ^{H2O} ~~TA XXXX~~~~TA XXXX~~~~TA XXXX~~

Time of incubation... 48 hrs... 37° C

Tester strains TA 100	Preparation in test mutants/petri dish 0.1 ml	Spontaneous mutants petri dish 0.1 ml	Evaluation
500 y	Ø - Ø	65 - 70 - 71	
500 y + S9	Ø - Ø	62 - 74 - 77	
250 y	Ø - Ø	65 - 70 - 71	
250 y + S9	Ø - Ø	62 - 74 - 77	
100 y	32 - 33	65 - 70 - 71	
100 y + S9	33 - 35	62 - 74 - 77	
20 y	44 - 58	65 - 70 - 71	
20 y + S9	58 - 57	62 - 74 - 77	
2 y	54 - 70	65 - 70 - 71	
2 y + S9	44 - 45	62 - 74 - 77	
Ø = tested, but no bacterial growth			

20.7.1977

Protocol-No. 52/77

AMES-Test

M 03 - 77

Preparation in test.....

Microsome fraction S9..... Rat [♂] Induction by..... Aroclor 1254Tester strains: mutants of S.typhimurium XXXXXS
XXXXXXODilutant:..... H₂O..... TA 1535
XXXXXXXX

Time of incubation.... 48 hrs... 37° C

Tester strains TA 1535	Preparation in test mutants/petri dish 0.1 ml	Spontaneous mutants petri dish 0.1 ml	Evaluation
500 y	Ø - Ø	4 - 9 - 13	
500 y + S9	Ø - Ø	8 - 9 - 18	
250 y	Ø - Ø	4 - 9 - 13	
250 y + S9	Ø - Ø	8 - 9 - 18	
100 y	8 - 14	4 - 9 - 13	
100 y + S9	6 - 7	8 - 9 - 18	
20 y	13 - 19	4 - 9 - 13	
20 y + S9	11 - 14	8 - 9 - 18	
2 y	21 - 24	4 - 9 - 13	
2 y + S9	14 - 19	8 - 9 - 18	
Ø = tested, but no bacterial growth			

2.8.1977

52 a/77

AMES-Test..... Protocol-No.

M 03 - 77

Preparation in test.....

Rat C Aroclor 1254
Microsome fraction S9.....Induction by.....

Tester strains: mutants of *S.typhimurium* XAXXXg8

TA 100

Dilutant: DMSO

Dilutant.....
48 hrs 37° C. - ~~22AXXDSG7~~

Time of incubation..... 48 hrs. 37 C .

Time of incubation: _____

AMES-Test..... 2.8.1977 Protocol-No..... 52 a/77

M 03 - 77
Preparation in test.....

Microsome fraction S9..... Rat C' Induction by Aroclor 1254

Tester strains: mutants of *S.typhimurium* TA 98

Dilutant:..... DMSO

Time of incubation 48 hrs. 37° C.

Time of incubation..... 48 hrs. 37 C

20.7.1977
AMES-Test..... 52/77.....

Preparation in test..... M 03 ~ 77

Microsome fraction S9..... Rat ♂ Induction by..... Aroclor 1254

Tester strains: mutants of *S.typhimurium* XXXX88

XXXX89

Dilutant:....^{H2O} 2..... XXXX88

TA 1537

Time of incubation....48.hrs...37° C

Tester strains	Preparation in test mutants/petri dish 0.1 ml	Spontaneous mutants petri dish 0.1 ml	Evaluation
TA 1537			
500 y	Ø - Ø	3 - 4 - 5	
500 y + S9	Ø - Ø	2 - 5 - 5	
250 y	Ø - Ø	3 - 4 - 5	
250 y + S9	Ø - Ø	2 - 5 - 5	
100 y	2 - 3	3 - 4 - 5	
100 y + S9	3 - 3	2 - 5 - 5	
20 y	3 - 5	3 - 4 - 5	
20 y + S9	2 - 3	2 - 5 - 5	
2 y	2 - 4	3 - 4 - 5	
2 y + S9	2 - 4	2 - 5 - 5	
Ø = tested, but no bacterial growth			

2.8.1977

52 a/77

AMES-Test..... 2.8.1977 Protocol-No.

Preparation in test M 03 - 77

Preparation in test.....
Microsome fraction S9..... Rat C¹ Induction by..... Aroclor 1254

Tester strains: mutants of *S.typhimurium* ~~X~~XXXXX98

XXAXXOXXXX

Dilutant: DMSO TA 1535

10 hrs 37° S MAX 53°

Time of incubation 48 hrs. 37° C.

Time of incubation.....



Abteilung: Krebsforschungslabor
Verfasser: Prof.D.Gericke

Datum: 20.7.1977
Bericht-Nr.: 52/77
Seite: 9

Compound	0.1ml/plate	S-9 Mix	Revertant colonies / plate			
			TA98	TA100	TA1535	TA1537
Positive controls						
2-acetylaminofluorene	10, [/] ug	+	402	415		
Benzo(a)pyrene	5, [/] ug	+		680	730	
β-Naphthylamine	10, [/] ug	+		275	310	
Neutralred	100, [/] ug	+			309	325
4-Nitro-o-phenylenediamine	50, [/] ug	-	2700	580		80
Streptozotozin	1, [/] ug		3000	642		85
					720	766