

Report Title: Test of the Potential of Octopirox to Induce Mutation in the Skin of Muta™ Mouse *In Vivo*.

Test Type: Genotoxicity Study

Conducting Laboratory and Location: P&G Miami Valley Laboratories, Biological Testing Facility, Cincinnati, OH

Test Substance(s): G0539.06 -- Octopirox in ethanol

Species: Muta Mouse *In Vivo*

of Animals: 3 mice per group

Test Conditions: The potential of Octopirox to induce mutations in skin was evaluated using Muta Mouse. This transgenic strain was developed to allow the detection of mutation in any tissue *in vivo*. A single application of 0.1ml of a 7.5% solution of OP in ethanol was topically applied to Muta Mouse skin. This dose was previously determined to be a maximum tolerated dose.

Results: The potential of Octopirox to induce mutations in skin was evaluated using Muta Mouse. This transgenic strain was developed to allow the detection of mutation in any tissue *in vivo*. The maximum tolerated dose was not mutagenic though in previous studies it had been shown to inhibit DNA synthesis in the Muta Mouse. See study #B91-0153.

Study #: B91-0227

Report Date: 4/27/92

QA report/GLP compliance: Yes

**Test of the Potential of Octopirox to Induce Mutation in the
Skin of MutaTMMouse *In vivo*
B91-0227**

Robert L. Binder, Audrey A. Erickson, Roman E. Frank and E. D. Thompson

Human and Environmental Safety Division
Miami Valley Laboratories
The Procter & Gamble Company
Cincinnati, OH 45239

Report Date 4/27/92

Procter & Gamble

The Procter & Gamble Company
Miami Valley Laboratories
P. O. Box 39870, Cincinnati, Ohio 45239-8700

QUALITY ASSURANCE STATEMENT

STUDY NUMBER B91-0227

TEST FACILITY The Procter & Gamble Company
Miami Valley Laboratories
Cincinnati, Ohio 45239

TYPE OF STUDY Test of the Potential of Octopirox to Induce
Mutation in the Skin of MutaMouse In Vivo

DIVISIONAL REQUEST DOCUMENT None

TSIN G0539.06 R0363.01, R0362.01

DATA LOCATION YE-1402

<u>PORTION(S) OF STUDY AUDITED</u>	<u>AUDITOR</u>	<u>DATE AUDITED</u>	<u>DATE REPORTED TO STUDY DIRECTOR</u>	<u>DATE REPORTED TO MAN- AGEMENT</u>
Dosing	L. K. Klahn	6/28/91	7/18/91	8/27/91
4-Day Sacrifice	L. K. Klahn	7/12/91	7/18/91	8/27/91
Study Data	L. K. Klahn	7/18/91	7/18/91	8/27/91

In compliance with the Good Laboratory Practices regulations this study has been audited by the Quality Assurance Unit and the results of those audits have been reported to the appropriate management. The protocol was audited for GLP required elements. The study data accurately reflects the procedures described in the protocol. The reported results accurately reflect the raw data of the study.

Stephen K. Klahn 4/29/91
Quality Assurance Unit - Date

SUMMARY SHEET

Study No.	B91-0227
Animal Activity No.	AA91-0081
Testing Facility:	Biological Testing Facility Miami Valley Laboratories The Procter and Gamble Co. P.O. Box 398707 Cincinnati, OH 45239
Test Substance(s):	Octopirox (TSIN G0539.06) Hydroxyurea (TSIN R0363.01) 7,12-Dimethylbenzanthracene (TSIN R0362.01)
Storage Conditions:	Octopirox (room temp.), hydroxyurea (refrigerated), 7,12-dimethylbenz- anthracene (approx. -20°C in dark)
DRD:	HESE 360
Date Study Started:	6/26/91 (Mice first shaved)
Date Study completed (in-life):	7/12/91
Report Date:	4/27/92
Study Director:	Robert L. Binder
Study Technicians:	Audrey A. Erickson Roman E. Frank
Genetic Toxicologist:	Edward D. Thompson
Notebook:	YB-1402
Archived at:	Miami Valley Laboratories

CONTENTS

I. PURPOSE.....	4
II. SUMMARY.....	4
III. METHODS	4
IV. RESULTS	7
V. DISCUSSION.....	8
VI. REFERENCES	9
TABLE 1. Assay Results from Mice Treated with 100 μ g DMBA	10
TABLE 2. Results from Mice Treated with Octopirox and Hydroxyurea	11
APPENDIX I. Hazleton Washington Report	12

I. PURPOSE

The purpose of this study was to determine whether a single maximum tolerated dose of Octopirox, applied topically, is mutagenic in mouse skin.

II. SUMMARY

The potential of Octopirox to induce mutations in skin was evaluated using MutaTMMouse. This transgenic strain was developed to allow the detection of mutation in any tissue *in vivo* (1,2). Octopirox was topically applied to MutaTMMouse skin in ethanol at a dose of 7.5 mg (~1100 $\mu\text{g}/\text{cm}^2$), which previously was determined to be a maximum tolerated dose (3). The potential of hydroxyurea (HU) to induce mutations in skin was also evaluated for comparison to Octopirox, and the established skin mutagen, 7,12-dimethylbenzanthracene (DMBA), was used as a positive control.

Since the optimal time after dosing to sample skin for mutation analysis was not known, skin samples from mice killed 4, 7 and 10 days after dosing with 100 μg of DMBA were analyzed for mutations. At 4 days after dosing with DMBA the mean mutant frequency was about twice that in controls. In contrast, at 7 and 10 days after dosing the mutant frequencies were about 13-times the control level. Therefore, both 7 and 10 days after dosing were appropriate times to collect skin for mutation analysis. At these times the mutant frequencies in skin treated with Octopirox were not greater than those in solvent controls. Also, 7 days after dosing, hydroxyurea (0.5 mg/g body weight i.p.) did not induce an increase in mutant frequency. Therefore, under the conditions of testing both Octopirox and hydroxyurea were not mutagenic *in vivo*.

Both Octopirox and hydroxyurea had been shown previously to cause a > 90% inhibition of epidermal DNA synthesis in MutaTMMouse skin at the doses shown here not to be mutagenic. Therefore, the transient inhibition of epidermal DNA synthesis by a single high dose of either Octopirox or hydroxyurea is not associated with mutagenesis in MutaTMMouse.

III. METHODS

Materials

Octopirox (TSIN G0539.06) was obtained from Beauty Care Product Development. Hydroxyurea (TSIN R0363.01) was obtained from the Sigma Chemical Co., and DMBA (TSIN R0362.01) was from the Aldrich Chemical Co. Other chemicals were of reagent grade or higher quality and their sources are indicated in the study notebook.

Animals

Male mice of the MutaTMMouse strain were received from Hazleton Research Laboratories at approximately 7 weeks of age, and were housed 5/shoebox cage on hardwood chip bedding. A 12 hr light/dark cycle (7:00 am to 7:00 pm) was maintained in the animal rooms (with the exception noted below), and Purina Lab Chow and water were available *ad libitum*. Room temperature and humidity were maintained to Biological Testing Facility (BTF) standards (BTF SOP: ENV 3.4). The mice were carefully shaved when they were 8

to 9 weeks old using a small animal clipper, and only mice in the resting phase of the hair cycle (i.e. animals without obvious hair regrowth within two days of shaving) were used. Mice were individually housed after shaving and treatments did not begin until at least 2 days after shaving. Initially all mice were maintained in room L-42 of the BTF. However, mice to be treated with DMBA or HU were transferred to room L-43 in the Carcinogen Area after shaving, and remained there until the end of the experiment.

Dosing

The overall experimental schedule and identity of groups are indicated below. Mice were uniquely identified with the group numbers (1-4) and letters (A-C) and animal numbers (1-3), which were written on their tails with a permanent marker before dosing. The MutaTMMouse strain has 3 coat colors: black, brown and golden brown. Mice of different colors were distributed as uniformly as possible among the various treatment groups.

Group	Animals	Treatment	Days Between Treatment and Killing
1A	1-3	0.1 ml ethanol (EtOH)	4
1E	1-3	0.1 ml EtOH	7
1C	1-3	0.1 ml EtOH	10
2A	1-3	7.5 mg Octopirox (Octc)	4
2E	1-3	7.5 mg Octc	7
2C	1-3	7.5 mg Octc	10
3A	1-3	0.5 mg/g Hydroxyurea (HU)	4
3E	1-3	0.5 mg/g HU	7
3C	1-3	0.5 mg/g HU	10
4A	1-3	100 μ g DMBA	4
4E	1-3	100 μ g DMBA	7
4C	1-3	100 μ g DMBA	10

Group 1 was treated with 0.1 ml of ethanol and group 2 received a single application of 0.1 ml of a 75 mg/ml (7.5%) solution of Octopirox in ethanol. The ethanol or Octopirox solution was dripped over the shaved area using a micropipettor to achieve uniform coverage, while avoiding the border so that the dose was not wicked into the surrounding hair.

Mice in group 3 received 0.5 mg hydroxyurea/g body weight dissolved in isotonic saline (150 mg hydroxyurea/ml) by sterile i.p. injection.

Mice in group 4 received a single topical application of 100 μ g of DMBA in 0.1 ml of acetone (1 mg/ml). The dose was dripped on the skin to achieve uniform coverage as indicated above for ethanol or the Octopirox solution. During dosing with DMBA and for approximately the next 1.5 hr the mice were maintained under yellow lights, then the white room lights were turned back on. Because of an error in setting the light switches in the room (L-43), the timer was bypassed and the lights remained on overnight. The

problem was corrected the following day. Groups 1 and 2 were housed in a separate room from groups 3 and 4 and were not affected.

Since the altered light cycle potentially could affect the circadian rhythm of epidermal DNA synthesis, additional mice were dosed with DMBA and HU as indicated below. Again dosing with DMBA was done under yellow light.

Group	Animals	Treatment	Days Between Treatment and Killing
3D	1-3	0.5 mg/g HU	7
4D	1,2	100 µg DMBA	7

All necessary precautions were taken to avoid any possible cross contamination of mice in the different treatment groups. Groups 3 and 4, while maintained in the same room, were kept physically separated on different cage racks in independent cage ventilation chambers.

Tissue Collection

At the times after dosing indicated above, mice were killed by CO₂ asphyxiation. The treated skin was excised, leaving a border to ensure only treated skin was sampled. Dissection instruments were carefully washed after each mouse to ensure there was no cross contamination of samples. Also, the mice were killed in the order of groups 1 -> 4 which minimized the significance of any inadvertent cross contamination. The skin samples were sealed in labelled plastic bags, frozen in liquid nitrogen, and stored at -80°C until sent to Hazleton Research Laboratories for mutation analysis.

Safety Considerations

Handling of all test materials was consistent with personnel protection and the SOP for the Carcinogen Laboratory (Standard Method MCM-26). All cages, cage tops and water bottles from mice treated with DMBA were disposed of as carcinogenic waste. DMBA-treated mice were dissected in the chemical fume hood in L-43 (BTF). Dissection instruments and carcasses were discarded as carcinogenic waste.

Analysis of Mutations in Skin

Analysis of mutations in the skin samples was done at Hazleton Laboratories. The samples analyzed were: 1B(1-3), 1C(1-3), 2B(1-3), 2C(1-3), 3D(1-3), 4A(1-3), 4B2, 4C(1-3), 4D(1,2). Samples 1A (1-3) and 2A (1-3), collected 4 days after dosing, were not analyzed because the positive control, DMBA, induced a much lower mutant frequency at 4 days after dosing compared to the effects at 7 and 10 days after dosing.

The methods for analysis of mutations in the skin samples are described in detail in Appendix I, the report from Hazleton Laboratories.

V. RESULTS

Results of mutation assays are described in detail in Appendix I and summarized below.

DMBA

Some crusting and focal ulceration of treated skin was noted 7 days after dosing mice with 100 μ g DMBA.

The skin mutation assay results from the DMBA-treated mice are shown in Table 1. Approximately 300,000 plaques were scored, since previous experience indicated that this was adequate to detect the strong response induced by this potent mutagen and carcinogen. DMBA was dosed in acetone, and a specific vehicle control for DMBA was not included in the experimental design. However the mutant frequencies in skin treated with ethanol, the vehicle for Octopirox, are shown in Table 2. These values ranged from 23.8×10^{-6} to 55.9×10^{-6} with an overall mean value obtained at 7 and 10 days after dosing of 40.1×10^{-6} . Similar mutant frequencies have been reported in the skin of untreated and acetone treated mice (2).

Since the optimal time after dosing to sample skin for mutation analysis was not known, skin samples from mice killed 4, 7 and 10 days after dosing with DMBA were analyzed for mutations. At 4 days after dosing the mutant frequency in DMBA-treated mice was about twice that in ethanol controls. However, at 7 and 10 days the mutant frequencies were about 13-times the mean level in ethanol controls, indicating that both 7 and 10 days after dosing were appropriate times to collect skin for mutation analysis.

Mouse 4D1, sampled at 7 days after dosing, had a mutant frequency of about 7000×10^{-6} , which is more than 10-fold greater than that determined in the skin of the other DMBA-treated mice. As discussed in Appendix I, such a high mutant frequency is probably the result of a heritable mutation in the target *lac Z* gene.

As noted in the Methods section most of the DMBA-treated mice were subject to an altered light cycle during the day of dosing. Mouse 4D2 was an additional mouse subsequently dosed and maintained under the normal light cycle and killed 7 days after dosing. The mutant frequency observed in 4D2 was similar to that in 4B2, which was also killed 7 days after dosing. Furthermore, the mutant frequencies determined in all the mice killed 10 days after dosing were similar to the value obtained with 4D2. Therefore, it is unlikely that the alteration in the light cycle affected the induction of mutations in skin by DMBA.


The strong mutation response induced by the positive control, DMBA, indicates that the MutaTM Mouse skin mutation assay as performed in this study was capable of detecting a known mutagen and skin carcinogen.

Octopirox and Hydroxyurea

Results from MutaTMMouse skin mutation assays on Octopirox and hydroxyurea are shown in Table 2. Based on the relatively low response to DMBA 4 days after dosing, only skin samples from the mice killed 7 and 10 days after dosing were evaluated for mutation. Approximately, 600,000 plaques were scored, since preliminary validation studies indicated this was an adequate number to assess mutant frequency. At both 7 and 10 days after dosing, the mutant frequencies in skin from mice treated topically with Octopirox were not greater than those in solvent controls. Also, 7 days after dosing hydroxyurea did not induce an increase in mutant frequency. Therefore, under the conditions of testing both Octopirox and hydroxyurea were not mutagenic *in vivo*.

IV. DISCUSSION

The dose of Octopirox evaluated here for mutagenic potential (7.5 mg or $\sim 1100 \mu\text{g}/\text{cm}^2$) (4) was shown previously to be a maximum tolerated dose based on the severity of dermal irritation induced by higher doses (3). This dose was not mutagenic, but in an earlier study caused $> 90\%$ inhibition of epidermal DNA synthesis in MutaTMMouse (4). Similarly, the dose of hydroxyurea tested (0.5 mg/g body weight i.p.), caused a 98% inhibition of epidermal DNA synthesis in MutaTMMouse (5); yet was not mutagenic *in vivo*. Together these results indicate that the transient inhibition of epidermal DNA synthesis by a single high dose of either Octopirox or hydroxyurea is not associated with mutagenesis *in vivo*.


Robert L. Binder
4/5 7/92
Date


Audrey A. Erickson


Roman E. Frank


Edward D. Thompson

VI. REFERENCES

1. Gossen, J.A., deLeeuw, W.J.F., Tan, C.H.T., Zwarthoff, E.C., Berends, F., Lohman, P.H.M., Knook, D. L., and Vijg, J. (1989) Efficient rescue of integrated shuttle vectors from transgenic mice: A model for studying mutations *in vivo*. Proc. Natl. Acad. Sci. USA **86**: 7971-7975
2. Myhr, B. C. (1991) Validation studies with MutaTMMouse: A transgenic mouse model for detecting mutations *in vivo*.
3. Binder, R.L., Erickson, A.A., Reer, F.J. and Frank, R.E. (1992) Evaluation of skin irritation potential of acetic acid and Octopirox in the MutaTMMouse. The Procter & Gamble Company, Study B91-0116 Report.
4. Binder, R.L., Erickson, A.A. and Frank, R.E. (1992) Inhibition of epidermal DNA synthesis in MutaTMMouse by Octopirox: time-course and dose-response. The Procter & Gamble Company, Study B91-0153 Report.
5. Binder, R.L., Erickson, A.A. and Frank, R.E. (1992) Inhibition of epidermal DNA synthesis in MutaTMMouse by hydroxyurea. The Procter & Gamble Company, Study B91-0210 Notebook reference YE-1402, pg 85.

TABLE 1

MutaTMMouse Skin Mutation Assay Results from Mice Treated with 100 μ g DMBA

Animal	Total Plaques	Mutants	M.F. X 10 ⁻⁶ *
<i>4 Days After Dosing</i>			
4A1	316,065	23	73
4A2	288,384	31	128
4A3	292,683	14	<u>48</u>
		AVG \pm SD	83 \pm 41
<i>7 Days After Dosing</i>			
4B2	312,690	137	438
4D1	139,725	977	Spontaneous Mutant
4D2	312,908	179	<u>572</u>
		AVG \pm SD	505 \pm 95 (N = 2)
<i>10 Days After Dosing</i>			
4C1	286,152	148	517
4C2	300,426	193	642
4C3	307,510	147	<u>459</u>
		AVG \pm SD	539 \pm 94

* M.F. indicates mutant frequency.

TABLE 2.

MutaTMMouse Skin Mutation Assay Results from Mice Treated with
Octopirox and Hydroxyurea

Animal	Total Plaques	Mutants	M.F. X 10 ⁻⁶ *
<i>7 Days After Dosing</i>			
<i>Ethanol Control</i>			
1B1	588,018	14	23.8
1B2	619,448	20	32.3
1B3	591,809	22	<u>37.2</u>
		AVG ± SD	31.1 ± 6.8
<i>Octopirox (7.5 mg)</i>			
2B1	602,998	10	16.6
2B2	668,547	26	38.9
2B3	672,949	23	<u>34.2</u>
		AVG ± SD	29.9 ± 11.8
<i>Hydroxyurea (0.5 mg/g body wt)</i>			
3D1	627,636	14	22.3
3D2	628,494	18	23.9
3D3	574,604	24	<u>41.8</u>
		AVG ± SD	29.3 ± 10.8
<i>10 Days After Dosing</i>			
<i>Ethanol Control</i>			
1C1	590,667	33	55.9
1C2	611,948	28	45.8
1C3	611,523	28	<u>45.8</u>
		AVG ± SD	49.2 ± 5.8
<i>Octopirox (7.5 mg)</i>			
2C1	600,587	17	28.3
2C2	707,802	21	29.7
2C3	615,523	26	<u>42.2</u>
		AVG ± SD	33.4 ± 7.7

* M.F. indicates mutant frequency.

APPENDIX

HAZLETON WASHINGTON REPORT

ANALYSIS OF MUTAMOUSE TISSUES FOR LACZ MUTANT FREQUENCY

SKIN SAMPLES FOR STUDY B91-0227



ANALYSIS OF MUTAMOUSE TISSUES FOR
LACZ MUTANT FREQUENCY

SKIN SAMPLES FOR STUDY B91-0227

FINAL REPORT

AUTHOR

BRIAN C. MYHR, Ph.D.

PERFORMING LABORATORY

HAZLETON WASHINGTON, INC.
5516 NICHOLSON LANE
KENSINGTON, MARYLAND 20895
and
9200 LEESBURG PIKE
VIENNA, VIRGINIA 22182

LABORATORY PROJECT ID

HWA STUDY NOS.: 14654-0-300
14654-1-300
14654-2-300

P&G REFERENCE NO.: 91000672

SUBMITTED TO

THE PROCTER & GAMBLE COMPANY
MIAMI VALLEY LABORATORIES
P.O. BOX 398707
CINCINNATI, OHIO 45239-8707

STUDY COMPLETION DATE

March 27, 1992

14654-0-300
14654-1-300
14654-2-300

QUALITY ASSURANCE STATEMENT

PROJECT TITLE: ANALYSIS OF MUTAMOUSE TISSUES FOR LACZ MUTANT FREQUENCY

PROJECT NO.: 25501ASSAY NO.: 14654PROTOCOL NO.: 300 EDITION NO.: 1

Quality Assurance inspections of the study and/or review of the final report of the above referenced project were conducted according to the Standard Operating Procedures of the Quality Assurance Unit and according to the general requirements of the appropriate Good Laboratory Practice regulations. Findings from the inspections and final report review were reported to management and to the study director on the following dates:

Inspection/DateFindings ReportedAuditor

Plaque Scoring/ 8-14-91

8-14-91

W. Yee

Final Report Review/ 3-26-92

3-27-92

D. Wallace



Quality Assurance Unit3-27-92
Date Released

14654-0-300

14654-1-300

14654-2-300




COMPLIANCE AND CERTIFICATION STATEMENT

The described study was conducted in compliance with the Good Laboratory Practice regulations as set forth in the Code of Federal Regulations (21 CFR 58, 40 CFR 792, and 40 CFR 160). To the best of the signers' knowledge, there were no significant deviations from the aforementioned regulations or the signed protocol that would affect the integrity of the study or the interpretation of the test results. The raw data have been reviewed by the Study Director, who certifies that the evaluation of the test article as presented herein represents an appropriate conclusion within the context of the study design and evaluation criteria.

All test and control results presented in this report and the supporting raw data are maintained in the archive files of Hazleton Washington, Inc., Vienna, Virginia. Copies of the raw data will be supplied to the sponsor upon request.

SUBMITTED BY:

LABORATORY SUPERVISOR:




Hoda Khouri, M.S.
Research Associate
Department of Genetic and
Cellular Toxicology

3/27/92

Date

STUDY DIRECTOR:



Brian C. Myhr, Ph.D.
Associate Director
Department of Genetic and
Cellular Toxicology

3/27/92

Study Completion
Date

14654-0-300
14654-1-300
14654-2-300

TABLE OF CONTENTS

QUALITY ASSURANCE STATEMENT	2
COMPLIANCE AND CERTIFICATION STATEMENT	3
SUMMARY	5
I. STUDY IDENTIFICATION	6
A. Sponsor	
B. Test Article	
C. Type of Assay	
D. Protocol Number	
E. Study Dates	
F. Supervisory Personnel	
II. OBJECTIVE	7
III. RATIONALE	7
IV. MATERIALS	7
A. Animal Tissue	
B. Bacterial Host	
V. EXPERIMENTAL DESIGN	8
A. Sample Selection and Size	
B. Mutation Assay	
C. Data Presentation	
D. Assay Acceptance Criteria	
E. Response Evaluation Criteria	
VI. RESULTS AND DISCUSSION	14
VII. CONCLUSION	15
VIII. REFERENCE	15
IX. EXPERIMENTAL DATA TABLES	16

SUMMARY

A CD2 transgenic mouse strain, designated MutaMouse, which contains multiple copies of a lambda gt10lacZ construct integrated into the genome of each cell, was used to evaluate the mutagenic activities of two test articles, G0539.06 and R0363.01. Twenty-four (24) coded samples of skin tissue were provided to Hazleton Washington, Inc., to determine the frequency of lacZ mutants. DNA was extracted from each skin sample, and the lacZ genes were packaged in vitro into lambda pre-heads. The resultant phage were plated onto E. coli C (lacZ⁻) agarose cultures in the presence of a chromogenic substrate (Xgal) for the lacZ gene product, β -galactosidase. The plaques that subsequently developed were colored blue for the nonmutant, lacZ⁺ genes, whereas mutations in the lacZ gene were indicated by colorless or light blue plaques. All plaques suspected of being mutants were replated to confirm the mutant phenotype. The ratio of confirmed mutant plaques to the lacZ⁺ plaques was calculated as the lacZ mutant frequency in each skin sample.

After analysis of the skin samples for the lacZ mutant frequency, the samples were decoded. Six negative control samples, harvested in triplicate at 7 days and 10 days after study initiation, showed mutant frequencies in the range of 23.8×10^{-6} to 55.9×10^{-6} , which were consistent with the Hazleton historical data base for skin. Other samples were exposed to 7,12-dimethylbenz(a)anthracene (DMBA), as a positive control, and harvested in triplicate at 4, 7 and 10 days after treatment. The samples obtained at 7 and 10 days showed large increases in mutant frequency (505×10^{-6} and 539×10^{-6} , respectively, as averages). Test article G0539.06 showed no increases in mutant frequency for triplicate samples harvested at 7 and 10 days; the individual samples ranged from 16.6×10^{-6} to 42.2×10^{-6} . Likewise, test article R0363.01 showed no increase in mutant frequency at a single harvest time of 7 days. These results demonstrate the lack of mutagenic activity in skin by G0539.06 and R036.01 under test conditions that led to large increases in mutant frequency after treatment with DMBA.

ANALYSIS OF MUTAMOUSE TISSUES FOR LACZ MUTANT FREQUENCY

I. STUDY IDENTIFICATION

A. Sponsor: The Procter & Gamble Company

B. Test Articles:

1. Sponsor Test Substance Identification Numbers (TSIN):

G0539.06 and R0363.01

supplied within the following coded skin samples (24):

1B1, 1B2, 1B3, 1C1, 1C2, 1C3, 2B1, 2B2, 2B3, 2C1, 2C2, 2C3,
3D1, 3D2, 3D3, 4A1, 4A2, 4A3, 4B2, 4C1, 4C2, 4C3, 4D1, 4D2

Coded skin samples (6, provided but not analyzed in this report: 1A1, 1A2, 1A3, 2A1, 2A2, 2A3.

2. Sponsor Study Number: B91-0227

3. Physical Description of the Test Samples: Frozen skin samples individually sealed in labeled bags.

4. Dates Received: July 10, July 16, and October 10, 1991

5. Genetics Assay Number: 14654

C. Type of Assay: Analysis of MutaMouse Tissues for LacZ Mutant Frequency

D. Protocol Number: 300, Edition 1

E. Study Dates:

1. Study Initiation Date: July 11, 1991

2. Experimental Start Date: July 23, 1991

3. Experimental Termination Date: December 12, 1991

F. Supervisory Personnel:

1. Laboratory Supervisor: Hoda Khouri, M.S.

2. Study Director: Brian Myhr, Ph.D.



II. OBJECTIVE

The objective of this study was to evaluate the *lacZ* mutant frequency in MutaMouse skin samples provided under code identification by Procter & Gamble.

III. RATIONALE

MutaMouse is a transgenic mouse which has bacterial *lacZ* genes stably integrated into the genome of every cell. This gene is unexpressed in vivo but serves as a potential target for mutagenesis by chemicals, radiation or endogenous processes. The particular advantage of the transgene is that the *lacZ* gene is inserted within a lambda gt10 vector which can be rescued enzymatically from the mouse DNA. DNA isolated from any mouse tissue of choice is reacted with a lambda packaging extract in order to produce viable phage. Each lambda phage contains one *lacZ* target gene obtained from the mouse DNA. The phage are adsorbed to a large excess of *E. coli* C bacteria and allowed to multiply and form plaques on an agar surface. By including a chromogenic substrate in the agar for the *lacZ* gene produce, β -galactosidase, a dark blue color is developed in the wildtype (nonmutant) plaques. However, mutations that cause a loss or reduction in β -gal activity result in colorless or light blue plaques. The ratio of altered plaques to blue plaques is the *lacZ* mutant frequency in the mouse DNA being analyzed. This assay is sensitive to mutagenic agents that cause base-pair changes or deletions in regions of the *lacZ* target required for enzymatic activity. The construction of MutaMouse and its characteristics were described as strain 40.6 by Gossen et al., 1989.

IV. MATERIALS

A. Animal Tissue

Skin tissue samples for analysis were obtained by Procter & Gamble from MutaMouse transgenic mice, designated CD₂-*lacZ*80/HazfBR. Frozen tissue samples packaged in dry ice were sent to HWA, and the tissues were stored at -70°C or colder until removed for the extraction of DNA.

B. Bacterial Host

The bacterial host for the development of plaques was *E. coli* C (*lac*⁻, *tet*^r, *amp*^r) obtained from Dr. Jan Gossen, Mediscand Ingeny, The Netherlands. Stocks were stored at -20°C or colder in 50% glycerol in Luria Bertani (LB) medium. Overnight cultures were subjected periodically to growth in LB medium containing ampicillin or tetracycline to guard against adventitious bacterial contamination. Daily cultures were prepared from overnight cultures in LB medium containing 0.2% maltose (LB/M medium). The daily cultures were grown in LB/M medium to a density of $OD_{600} = 0.5$ to 1.0 prior to collection for use in the assay.

EXPERIMENTAL DESIGN

A. Sample Selection and Size

In this study, skin samples were the only tissue samples provided. The sponsor initially requested that 300,000 blue plaques (lacZ⁺ gene population) be analyzed per sample. A number of samples were then selected for further analysis up to 600,000 lacZ⁺ genes. Six coded skin samples were placed on hold and were not analyzed for this report (see Section I.B.)

E. Mutation Assay

1. Tissue preparation

Each sample of frozen skin was removed from -70°C storage and cut into small pieces with scissors. The pieces were then dispersed into lysis buffer (10 mM Tris·HCl, pH 8.0, 150 mM NaCl, 20 mM EDTA). Each dispersed tissue was digested for 2-3 hours at approximately 50°C in a shaking waterbath using 0.9% sodium dodecylsulfate and 0.9 mg/ml proteinase K in lysis buffer.

2. DNA preparation

The tissue digest was mixed with an equal volume of equilibrated phenol:chloroform (1:1), using gentle inversion to produce a homogeneous suspension. The water phase (containing the DNA) was separated from the organic phase by centrifugation at approximately 1500 x g for 15 minutes. The water phase was then transferred to a clean tube, while being careful to avoid insoluble material at the interface.

A one-fifth volume of 8M potassium acetate was added to the water phase and mixed gently. Then one volume of chloroform was added and mixed by slow inversions. After centrifugation at approximately 2500 x g for 30 minutes at approx. 4°C, the upper aqueous phase was transferred to a clean tube. DNA was precipitated by the addition of 2 to 2.5 volumes of cold absolute ethanol.

The collected DNA was dissolved in TE-4 buffer (10 mM Tris·HCl, pH 7.5, 4 mM EDTA) and refrigerated indefinitely at approx. 4°C. The DNA concentration was determined from the ultraviolet absorbance at 260 nm, using the relationship of 70 µg/ml of calf thymus DNA for an absorbance of 1.0 (Rodriguez and Tait, 1983).

The molecular weight distribution of each DNA preparation was determined by electrophoresis in 0.5% agarose. The migration of the DNA and DNA molecular weight markers were visualized by ethidium bromide staining and photography of induced fluorescence. DNA preparations of good quality for the packaging reaction, defined as DNA remaining primarily between the origin and the 25,000 base-pair position, were used for the analysis of mutant frequency.

3. DNA packaging

Sequences of λ gt10lacZ were removed from the mouse DNA and packaged into empty lambda preheads by use of lambda packaging extracts. The extracts were purchased commercially and consisted of two components (freeze/thaw extract and sonicate). The components were stored at -70°C or colder and were quickly thawed just prior to use.

Each DNA preparation was diluted with TE-4 buffer to a concentration of $1.5\text{ }\mu\text{g DNA}/\mu\text{l}$. A volume of $5\text{ }\mu\text{l}$ was carefully pipeted into $10\text{ }\mu\text{l}$ of freeze/thaw extract and stirred. The sonicate ($15\text{ }\mu\text{l}$) was then immediately added and stirred. The reaction was incubated for approximately 3 hours at $37 \pm 2^{\circ}\text{C}$, then terminated by the addition of 0.5 to 1 ml SM buffer (50 mM Tris-HCl, pH 7.5, 10 mM MgSO_4 , 100 mM NaCl, 0.01% gelatin) and 30 μl of chloroform as preservative. The reaction mixture was pulse-spun to pellet the chloroform, and the preparation was stored in the dark under refrigeration (-4°C) until analyzed.

4. Lambda plaque production

The number of reliable lambda phage in each packaged DNA preparation was estimated by a titering. A measured volume was removed and incubated for 20-30 minutes at room temperature with 1 ml of *E. coli* C prepared as described below. LB/MM agarose medium was prepared, and 14 ml were added to the tube. After mixing, pouring the contents into a 150 mm plate, and allowing solidification, the culture was incubated overnight at approximately 37°C to allow plaques to develop. The plaques were counted and used to calculate the volume of packaged DNA preparation needed to produce approximately 1500 plaques.

The packaged DNA preparations from each tissue sample were analyzed by producing 100,000 or more plaques at a nominal (precalculated) density of approximately 1500 plaques per 150 mm culture plate. Usually, several packaging reactions and platings (on different days) were used to achieve the total population desired for each tissue sample. The plating procedure is described in the last paragraph of this section.

LB/MM agarose medium was prepared on the day of use. The contents were 0.75% purified agarose in LB medium to which 10 mM $MgSO_4$, 0.2% maltose, and 0.35 mg/ml X-gal were added aseptically to molten medium (approx. 55°C) just prior to use. The X-gal solution was previously dissolved in dimethylformamide (DMF) and then diluted into the agarose solution such that the DMF concentration was less than 0.5%.

E. coli C bacteria were obtained from daily cultures as described in Section IV.B. The cultures were centrifuged and resuspended in fresh LB/MM medium at OD_{600} of approximately 0.40-0.60, which corresponds to approximately 2 to 4×10^8 bacteria/ml. This suspension was held on ice and used within 1-2 hours of preparation.

A series of tubes was prepared for the adsorption of phage to the bacteria. Normally, three plates were poured from one adsorption, so the volume of packaged DNA preparation used corresponded to a titer of approximately 4500 phage. A 3 ml aliquot of bacterial suspension was added, and the suspension was mixed and incubated for 20-30 minutes at room temperature. LB/MM agarose (42 ml) was added and mixed. The contents were then quickly pipetted into three empty 150 mm culture plates at 15 ml/plate, spread evenly, and allowed to harden. The cultures were incubated overnight at approximately 37°C to allow plaques to develop.

5. Plate scoring

All cultures were scored on the day of removal from the incubator. The total population was determined from blue plaque counts on four representative plates. An area equivalent to one-tenth the plate area was marked on the plate surface and used as the count area. The average plaque count and standard deviation were determined. The total population was calculated as the average blue plaque count $\times 10 \times$ the number of plates scored for mutant plaques.

Each plate was examined carefully for presumptive mutant plaques, identified as colorless plaques or plaques with a reduced color relative to the majority of blue, wildtype plaques. The positions were marked and the plate given to a second person to locate any additional presumptive mutants.

All presumptive mutant plaques were picked by collecting a core through the plaque with a clean Pasteur pipet. The contents were transferred to a microfuge tube containing 300 μ l SM buffer and 50 μ l chloroform. After mixing well, the phage preparations were allowed to sit for several hours or overnight in a refrigerator. The phage preparations can be held

indefinitely at approximately 4°C in the dark but were assayed for mutant content within 7 days to avoid possible loss of infectivity.

6. Mutant confirmation

Presumptive mutants were tested for mutant phenotypes by replating onto *E. coli* C cultures containing λ -gal. Bacterial suspensions in LB/MM agarose were prepared as already described in Section V.B.4, except that phage were not included. Either 150 mm plates or the much larger 24 cm x 24 cm bottom agar plates were used. In the latter case, 30 ml of LB/MM agarose were mixed with 3 ml of bacteria suspension and poured for each plate. After hardening, the plates were positioned on a grid pattern of approximately 15 mm x 15 mm. A 10 μ l aliquot of each presumptive mutant phage preparation was applied to one grid area. After the spots had adsorbed into the agarose, the plates were incubated overnight at approximately 37°C.

A presumptive mutant phage sample was scored as mutant if any colorless or light blue plaques were found in the spotted area. An all-blue response of wildtype intensity in the spotted area was scored as negative for mutation.

If the spotted area was completely lysed so that individual plaques were not distinguishable, a second confirmation test with a smaller titer was performed. Conversely, if no plaques were observed, a second trial with a larger titer was used for that phage preparation.

7. Mutant frequency calculation

The mutant frequency (MF) in a DNA sample was defined as the ratio of confirmed mutant plaques to the total plaques analyzed. The MF was expressed in units of 10^6 plaques. To increase the population of analyzed *lacZ* genes for a given skin DNA sample, the results from several packaging reactions and phage platings were combined and used for the MF calculation.

C. Data Presentation

For each skin DNA sample analyzed, the following data are reported: the total number of plaques scored, the number of colorless mutant plaques found, the number of color-mutant plaques, the total number of mutant plaques, and the calculated mutant frequency.

D. Assay Acceptance Criteria:

The assay was considered acceptable for evaluation of the results by meeting the following criteria:

- The negative control tissue yielded a mutant frequency consistent with historical experience for that tissue. The average mutant frequency among the 6 negative control samples was $(40 \pm 11) \times 10^{-6}$, which was not significantly different from the HWA historical negative control data provided below:

HWA HISTORICAL NEGATIVE CONTROL MUTANT FREQUENCY IN SKIN

ANIMAL*	MUTANT LACZ GENES	WILDTYPE LACZ GENES	MUTANT FREQ. $\times 10^6$
1	29	1,059,759	27.4
2	28	1,013,281	27.6
3	37	1,062,595	34.8
4	37	1,092,815	29.3
5	24	1,156,650	20.7
6	28	1,204,389	23.2
7	34	1,050,676	32.4

* Male, 8-10 weeks

Avg. MF = $(27.9 \pm 4.9) \times 10^{-6}$

- The plaque density on the cultures scored for mutants did not exceed 2000 plaques (to 5% precision) per 150 mm plate.
- The efficiency of the packaging reaction for each DNA sample did not fall below 10,000 pfu/ μ g DNA. Because most DNA samples were packaged in several individual reactions, an occasional reaction in the 9000 pfu/ μ g DNA range was acceptable if the majority of plaques were obtained at higher efficiency.

E. Response Evaluation Criteria

All tissue samples were analyzed at plaque populations greater than 100,000 and in accordance with directives from the sponsor. Because the concurrent negative controls were comparable to the HWA historical negative control data, the concurrent negative controls were used to evaluate the responses in the treated tissues.

An appropriate statistical method for data analysis is not yet developed. Therefore, a guideline of at least a 2-fold increase in mutant frequency over the concurrent negative control mean was used as a minimum response to indicate mutagenic activity. In this study, that mutant frequency was 80×10^{-6} or greater. In addition, a consistently elevated response among the replicate samples was necessary to consider a treatment to be mutagenic.

VI. RESULTS AND DISCUSSION

All of the skin samples evaluated in this study were coded. Upon completion of the study, the code was broken by the sponsor, and the data were arranged into the appropriate treatment groups for analysis, as shown in Table 1.

The results in Table 1 show that the negative control skin samples, harvested 7 days after initiation of the study, yielded mutant frequencies within the HWA historical distribution of $(28 \pm 5) \times 10^{-6}$. Higher mutant frequency values were obtained at 10 days after study initiation. It is not clear why this second set of controls should yield apparently higher mutant frequencies, but if all six control animals are considered together, the resultant mean and standard deviation, $(40 \pm 11) \times 10^{-6}$, was not significantly different from the HWA historical distribution obtained from approximately 1,000,000 plaques per animal. Hence, the concurrent negative control, consisting of all 6 samples, was considered to be appropriate for the evaluations of the treated skin samples.

7,12-Dimethylbenz(a)anthracene (DMBA) was used as a positive control test substance. As shown in Table 1, the DMBA treatment caused only a minimal increase in mutant frequency by 4 days after application. However, by 7 days after treatment, a large increase in mutant frequency was observed. This response increased somewhat further to an average mutant frequency of 539×10^{-6} at 10 days. (One of the animals, coded 4D1, showed a mutant frequency that was indicative of a mutant offspring in the animal breeding program, so this animal was excluded from the analysis of DMBA-induced mutagenesis.) These results demonstrated the mutagenicity of a well-known chemical mutagen and showed the importance of harvest time in detecting this activity. The time course for the maximum accumulation of induced mutants in skin is not established for any chemical, but 7-10 days' harvest time appears to approach the time necessary for maximum assay sensitivity.

Compound G0539.06 did not induce any significant increases in mutant frequency when assayed 7 and 10 days after treatment. The mutant frequencies in the skin samples from individual animals ranged from 16.6×10^{-6} to 42.2×10^{-6} , which fell within the negative control range of frequencies. The lack of any response at 7 or 10 days after treatment indicated that G0539.06 did not possess any significant mutagenic activity in mouse skin.

Compound R0363.01 was analyzed 7 days after treatment. Again, no departure from the range of mutant frequencies typical of negative control skins were observed. Thus, no evidence for mutagenic activity was obtained for this test substance.

VII. CONCLUSION

MutaMouse skin samples assayed 7 and 10 days after the application of 7,12-dimethylbenz(a)anthracene showed large increases in the lacZ mutant frequency. In contrast, test compound G0539.06 did not induce any increases in mutant frequency at these harvest times. Similarly, test compound R0363.01 caused no change in mutant frequency at 7 days following treatment. Thus, no evidence for mutagenic activity in skin was obtained for either test compound in this study.

VIII. REFERENCES

Gossen, J.A., deLeeuw, W.J.F., Tan, C.H.T., Zwarthoff, E.C., Berends, F., Lohman, P.H.M., Knook, D.L. and Vijg, J.: Efficient rescue of integrated shuttle vectors from transgenic mice: a model for studying mutations in vivo, Proc. Natl. Acad. Sci. USA, 86:7971-7975, 1989.

Rodriguez, R.L. and Tait, R.C.: Recombinant DNA Techniques: An Introduction, The Benjamin/Cummings Publishing Company, Inc., Menlo Park, CA, 1983, pp 42-43.

IX. EXPERIMENTAL DATA TABLES

14654-0-300
14654-1-300
14654-2-300

TABLE 1

LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

 P&G Study No.: B91-0227
 P&G Reference No.: 91000672
 WMA Assay No.: 14654

Tissue Sample Identity*	Packaging Reaction No. ^b	Reaction Titer (pfu/ μ g DNA)	Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF $\times 10^4$
1B1	1633-1	36,027	1163	87,225	1	0	1	
Negative	1633-1	36,027	1543	46,290	1	0	1	
Control,	1633-1	36,027	1215	21,870	2	1	3	
7 days	1633-2	18,827	1580	142,200	1	1	2	
	1633-3	20,507	1680	151,200	3	2	5	
	1633-4	22,320	1365	<u>139,230</u>	2	0	<u>2</u>	
				588,015			14	23.8
1B2	1633-1	21,413	1363	122,670	3	0	3	
Negative	1633-2	12,960	1515	99,990	1	0	1	
Control,	1633-3	14,747	1535	96,705	4	0	4	
7 days	1633-4	16,880	1428	114,240	5	1	6	
	1633-5	17,787	1595	130,790	5	0	5	
	1633-6	21,493	1835	<u>55,050</u>	1	0	<u>1</u>	
				619,445			20	32.3
1B3	1633-1	23,573	1323	119,070	2	1	3	
Negative	1633-2	16,907	1420	113,600	3	1	4	
Control,	1633-3	18,347	1385	70,635	1	0	1	
7 days	1633-4	32,640	1448	221,544	10	2	12	
	1633-5	11,440	1240	<u>66,960</u>	2	0	<u>2</u>	
				591,809			22	37.2

 14654-C-30C
 14654-I-30C
 14654-2-30C

TABLE 1 (Continued)

LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	Packaging Reaction No. ^b	Titer (pfu/ μ g DNA)	Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF $\times 10^6$
1C1	1670-1	18,853	1383	121,704	9	0	9	
Negative	1670-2	20,000	1293	116,370	5	5	10	
Control,	1670-3	14,560	2005	138,345	3	0	3	
10 days	1670-4	13,893	1588	104,808	7	0	7	
	1670-5	14,773	1575	66,150	1	1	2	
	1670-5	14,773	1443	43,290	1	1	2	
				590,667			33	55.9
1C2	1670-1	23,307	1435	159,285	5	5	10	
Negative	1670-2	13,467	1505	96,320	4	1	5	
Control,	1670-3	19,573	1578	63,120	1	0	1	
10 days	1670-3	19,573	1383	53,937	1	1	2	
	1670-4	11,093	1370	69,870	1	0	1	
	1670-5	25,573	1448	169,416	0	0	9	
				611,948			28	45.8
1C3	1670-1	32,613	1188	163,944	1	4	5	
Negative	1670-2	23,547	1093	65,580	2	1	3	
Control,	1670-2	23,547	1118	60,372	1	2	3	
10 days	1670-3	11,813	1555	88,635	3	0	3	
	1670-4	13,627	1438	99,222	4	4	8	
	1670-5	10,107	1790	85,920	3	1	4	
	1670-6	9,787	1450	47,850	0	2	2	
				611,523			28	45.8

TABLE 1 (Continued)

LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	<u>Packaging Reaction</u>		Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF x 10 ⁶
No. ^b	Titer (pfu/μg DNA)							
4A1 DMBA, 4 days	1670-1	20,800	1165	118,830	6	0	6	
	1670-2	10,667	1350	68,850	5	2	7	
	1670-3	16,293	1585	128,385	8	2	10	
				316,065			23	72.8
4A2 DMBA, 4 days	1670-1	19,520	1375	132,000	12	1	13	
	1670-2	12,907	1423	89,649	12	2	14	
	1670-3	9,707	1483	66,735	8	2	10	
				288,384			37	128.3
4A3 DMBA, 4 days	1670-1	22,267	1228	136,308	5	0	5	
	1670-2	12,640	1375	86,625	3	1	4	
	1670-3	9,867	1550	69,750	3	2	5	
				292,683			14	47.8

14654-0-30C
14654-1-30C
14654-2-30C

14654-0-30C
14654-1-30C
14654-2-30C

TABLE 1 (Continued)
LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	Packaging Reaction		Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF x 10 ⁶
	No. ^b	Titer (pfu/μg DNA)						
4B2 DMBA, 7 days	1633-1	34,587	1295	194,250	65	4	69	438.1
	1633-2	24,747	1410	<u>118,440</u>	59	9	<u>68</u>	
				312,690			137	
4D1 DMBA, 7 days	1633-1	28,720	1035	139,725	973	4	977	6992
4D2 DMBA, 7 days	1633-1	19,360	1513	139,196	72	15	87	572.1
	1633-2	17,280	1533	128,772	58	8	66	
	1633-3	31,813	1498	<u>44,940</u>	22	4	<u>26</u>	
				312,908			179	

14654-0-30C
14654-1-30C
14654-2-30C

TABLE 1 (Continued)
LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity*	Packaging Reaction No. ^b	Titer (pfu/ μ g DNA)	Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF $\times 10^6$
4C1 DMRA, 10 days	1670-1	33,307	1203	191,277	83	6	89	517.2
	1670-2	27,440	1375	<u>94,875</u>	45	14	<u>59</u>	
				286,152			148	
4C2 DMRA, 10 days	1670-1	24,880	1268	148,356	82	23	105	642.4
	1670-2	23,360	1370	<u>152,070</u>	75	13	<u>88</u>	
				300,426			193	
4C3 DMRA, 10 days	1670-1	22,453	1195	118,305	60	11	71	458.5
	1670-2	13,200	1435	94,710	27	9	36	
	1670-3	18,400	1210	22,990	6	2	8	
	1670-3	18,400	1355	36,585	14	4	18	
	1670-4	17,493	1455	<u>34,920</u>	6	2	<u>8</u>	
				307,510			141	

14654-0-30C
14654-1-30C
14654-2-30C

TABLE 1 (Continued)
LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	Packaging Reaction No. ^b	Titer (pfu/ μ g DNA)	Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF $\times 10^4$
2B1 Compound G0539.06, 7 days	1633-1	18,160	1190	107,100	2	0	2	16.6
	1633-2	23,787	1353	154,242	2	0	2	
	1633-3	14,053	1578	47,340	2	0	2	
	1633-3	14,053	1598	43,146	0	1	1	
	1633-4	12,480	953	57,180	0	0	0	
	1633-5	13,147	1795	114,880	1	0	1	
	1633-6	15,920	1515	36,360	0	0	0	
	1633-7	9,093	1425	<u>42,750</u>	2	0	<u>2</u>	
				602,998			10	
2B2 Compound G0539.06, 7 days	1633-1	36,560	1100	99,000	1	0	1	38.9
	1633-1	36,560	1283	107,772	6	1	7	
	1633-2	18,053	1435	93,275	6	0	6	
	1633-3	26,747	1700	212,500	6	1	7	
	1633-4	23,173	1625	<u>156,000</u>	5	0	<u>5</u>	
				668,547			26	
2B3 Compound G0539.06, 7 days	1633-1	21,040	1253	72,674	3	1	4	34.2
	1633-2	17,813	1515	136,350	4	0	4	
	1633-3	51,600	1430	102,960	4	2	6	
	1633-3	51,600	1733	181,965	3	1	4	
	1633-4	21,440	1790	<u>179,000</u>	3	2	<u>5</u>	
				672,949			23	

TABLE 1 (Continued)

LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	Packaging Reaction		Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF x 10 ⁶
	No. ^b	Titer (pfu/μg DNA)						
2C1 Compound G0539.06, 10 days	1670-1	20,960	1178	123,690	2	0	2	28.3
	1670-2	19,173	1333	123,969	3	1	4	
	1670-3	10,000	1433	61,619	1	1	2	
	1670-4	10,000	1695	81,360	0	1	1	
	1670-5	11,093	1950	101,400	2	1	3	
	1670-6	13,280	1723	<u>108,549</u>	2	3	<u>5</u>	
				600,587			17	
2C2 Compound G0539.06, 10 days	1670-1	19,413	863	80,250	1	1	2	29.7
	1670-2	21,920	1463	149,226	2	0	2	
	1670-3	13,147	1845	105,165	2	0	2	
	1670-4	21,067	1918	195,636	7	4	11	
	1670-5	15,333	1458	109,350	0	1	1	
	1670-6	13,253	1623	<u>68,166</u>	3	0	<u>3</u>	
				707,802			21	
2C3 Compound G0539.06, 10 days	1670-1	22,133	1335	128,160	4	0	4	42.2
	1670-2	10,000	1763	84,624	3	0	3	
	1670-3	17,520	968	78,408	5	0	5	
	1670-4	22,347	1533	133,371	5	2	7	
	1670-5	26,453	1540	<u>190,960</u>	5	2	<u>7</u>	
				615,523			26	

 14654-C-30C
 14654-1-30C
 14654-2-30C

14654-0-30C
14654-1-30C
14654-2-30C

TABLE 1 (Continued)
LacZ MUTANT FREQUENCY IN MUTAMOUSE SKIN SAMPLES

Tissue Sample Identity ^a	Packaging Reaction No. ^b	Titer (pfu/ μ g DNA)	Plating Density (LacZ ⁺ plaques/plate)	Total LacZ ⁺ Population	Colorless LacZ Mutants	Colored LacZ Mutants	Total LacZ Mutants	MF $\times 10^6$
3D1 Compound R0363.01, 7 days	1633-1	31,653	1375	206,250	2	0	2	
	1633-2	13,813	1320	89,760	2	3	5	
	1633-3	12,933	1663	89,802	1	1	2	
	1633-4	10,693	1848	116,424	1	0	1	
	1633-5	17,893	1425	<u>125,400</u>	4	0	<u>4</u>	
				627,636			14	22.3
3D2 Compound R0363.01, 7 days	1633-1	31,067	1225	180,075	6	0	6	
	1633-2	14,933	1568	114,464	3	2	5	
	1633-3	24,480	1870	220,660	2	0	2	
	1633-4	14,000	1743	<u>113,295</u>	2	0	<u>2</u>	
				628,494			15	23.9
3D3 Compound R0363.01, 7 days	1633-1	15,733	1290	96,750	7	0	7	
	1633-2	12,907	1438	89,156	1	1	2	
	1633-3	17,627	1653	138,852	5	0	5	
	1633-4	14,000	1888	130,272	1	0	1	
	1633-5	16,000	1533	<u>119,574</u>	9	0	<u>9</u>	
				574,604			24	41.8

^a The test sample code, experimental identity revealed after analysis, and the harvest times (days after dosing) are given.

^b Packaging reaction number: Gigapack II Gold lot number (nnnn) and the individual reactions performed (-n) for each tissue sample.

pfu = plaque forming unit MF = LacZ Mutant Frequency

DMBA = 7,12-Dimethylbenz(a)anthracene, 100 μ g per animal, used as positive control.