



# SmartCare®

October 27, 1999

1023 91 OCT 28 1999

Docket No. 99D-2335  
Dockets Management Branch  
Division of Management Systems and Policy  
Office of Human Resources and Management Services  
Food and Drug Administration  
5630 Fishers Lane, Room 1061 (HFA-305)  
Rockville, MD 20852

To Whom It May Concern:

We applaud the FDA's apparent overall goal to lower the protein content of natural rubber latex gloves and thereby reduce the symptom severity of latex allergic individuals and potentially reduce the rate of sensitization in high risk groups. However, we have several serious concerns detailed below regarding the proposed regulations.

**1. Protein Labeling (*Proposed 21 CFR 801.440*)**

The FDA proposes to amend the existing latex cautionary label ("*Caution: This product contains natural rubber latex which may cause allergic reactions.*") by adding an additional statement as to latex protein content: "*The FDA recommends that this product contain no more than 1200  $\mu$ g extractable protein per glove. This product contains \_\_\_\_\_  $\mu$ g extractable protein per glove.*"

First, the units used in this statement are not consistent with the recently revised 1999 ASTM rubber glove standards (ASTM D3577 and D3578) or the Malaysian Rubber Board Standards, which use mg per decimeter<sup>2</sup>. Because the latter unit reflects surface area rather than weight, it better addresses the issue of surface protein level and does not apply a different standard to a large glove versus an extra-small glove. Either these units or the ASTM units should be modified to be consistent with one another.

Based on our recent analysis of approximately 20 brands of powdered and powder-free latex gloves, most of the powder-free brands would comply with this labeling standard whereas a large portion of the powdered gloves would not. If this standard is not intended to apply to existing gloves (but only new 510k applicants), it will have little impact on the majority of gloves currently available in the United States and will, therefore, fail to accomplish the intended goal of lowering latex glove protein content for the benefit of end users.

Furthermore, the mixing of protein content labeling and cautionary statements sends a badly mixed message to the end user and is particularly confusing for gloves with lower protein levels. We frequently hear from latex allergic consumers who believe it is safe for them to use a "low protein" or "hypoallergenic" glove. Despite the current regulations, there are still latex glove products using these marketing claims. These new changes further muddy the consumers' vision of the

98N-0313

C23

issues. Content labeling should not be mixed with caution statements. An analogy that illustrates the problems encountered by this confusing approach can be found in cigarette marketing claims. If "Caution: The Surgeon General has determined smoking can be hazardous to your health," was followed by "The FDA recommended the nicotine/carcinogen level be below\_\_\_\_," the emphasis and credibility of the warning label would be severely undermined.

Finally, the additional or alternate caution statement proposed for those gloves with a protein content less than 50 µg per gram (or 300 µg per glove) should be applied to all latex gloves, not just to those below the "detection" limit for the inherently problematic Lowry assay. The safe use of any latex glove by NRL-sensitized individuals is simply an oxymoron. As noted in the medical literature, while Lowry protein levels generally correlate with symptom responses, Lowry total protein level is not an absolute predictor of antigenic protein content. Furthermore, there are no established total or antigenic protein concentration thresholds for symptom elicitation or sensitization.

In general, the Lowry assay is subject to significant variability (as noted in the 1999 ASTM standards) and should not be misinterpreted until these technical issues can be resolved. Our recent analysis of nearly 40 brands of latex gloves analyzed by the Lowry method in two different laboratories showed as much as a 20 µg/g difference between results for powder free latex brands (those most likely to bear a different label). In addition, antigenic (LEAP) protein content varied considerably on these same samples, and not all latex gloves with Lowry contents less than 50 µg per gram had similarly low LEAP values. As the FDA is no doubt aware, subtle changes in Lowry assay procedures between laboratories can influence extraction efficiency, artificially lower NRL protein content (Koch, 1997), and affect assay sensitivity (Yip, 1997). Chemicals commonly added during the manufacture of NRL can also interfere with the Lowry assay (Chen *et al.*, 1997). We have also observed LEAP and Lowry values about the detection limits on nitrile and vinyl gloves. Clearly, neither the current Lowry or LEAP assay methods are robust enough to discriminate reliably and dictate safety thresholds.

## **2. Powder labeling (*Proposed 21 CFR 801.440*)**

We believe this proposed caution statement is inappropriately overstated and strong: "*Caution: Glove powder is associated with adverse reactions. FDA recommends that this product contain...*" It seems ludicrous that natural rubber latex (associated by the FDA with several fatalities and serious occupational health problem) and powder produces only *allergic* reactions, while glove powder on synthetic non-latex gloves produces *adverse* reactions. Why is a similar statement not included on natural rubber latex glove proposed labeling where the majority of health hazard documentation has been generated?

With regard to the 120 mg per glove limit proposed, we are concerned with the determination basis for this threshold. Historically, the use of powdered gloves during surgeries has been associated with an increased incidence of granulomas and peritoneal adhesions as summarized by the FDA Glove Powder Report in September 1998. However, while the studies are well documented, they often failed to distinguish the type of powder (cornstarch, talc, etc.) or glove material (NRL or other) used. Many factors are suspected to play a role in this problem, such as the method of sterilization, allergic reactions to cornstarch (or NRL, or other contaminants), endotoxins and

contamination with microorganisms. Recent published studies demonstrate that glove powder on latex gloves contributes significantly to symptom severity and IgE levels in latex allergic workers, and may facilitate sensitization of non-allergic individuals. However, we are not aware of any peer-reviewed published study in which 120 mg of powder per glove (or any other level) has been demonstrated as a safety threshold for either symptom elicitation or sensitization or even for post-surgical granuloma/adhesion development for latex or synthetic non-latex gloves. This labeling comment is completely inappropriate when applied only to synthetic non-latex gloves as proposed of which there is little data to support adverse reaction claims at the proposed level.

Similar to our comments above concerning latex content values, 120 mg of powder per glove represents a much different powder level on a small glove than on a large glove and should be expressed in relation to surface area. Furthermore, the current analytical method for determination of powder on a powdered glove may not be sufficiently rigorous to provide this arbitrary level with much accuracy.

As with the above comments regarding latex content labels, it is not appropriate to mix caution statements with content labels as this often implies to the user that products below the recommended levels are safe for use.

**3. Expiration date labeling (*Proposed 21 CFR 801.440(d)*)**

With regard to surgeon's gloves, it seems inappropriate to require real time or accelerated aging studies to support shelf-life data that is not applicable to the product. As the maintenance of sterility or packaging integrity generally supercedes any expiration date based upon physical and mechanical integrity, it seems unwarranted to make these potentially useless determinations. Furthermore, it is unclear from the proposed language exactly how the expiration date will be timed. Specifically, does expiration begin from the time of production/manufacture, from the time of chlorination, from the time of packaging or (in the case of sterile product) from the time of sterilization?

How does the Agency plan to address the potential fire hazard issue with their recommendation for shrink wrapping exam gloves that may be of particular concern with the increase of chlorinated gloves? Does the Agency plan additional recommendations for shelf life and physical property assessment of chlorinated products?

**4. Quality System Regulation (21 CFR Part 820) Page 2-9**

It is clear that medical grade examination gloves are being used for a wide variety of applications and procedures within healthcare. Does the Agency plan to more clearly define the "intended use" definition to address performance in use (i.e., chemical exposure, concentration, duration of exposure or wear)? (Page 2-10) What is the Agency's role in determining adherence or compliance with quality system requirements post Pre-Market 510(k) acceptance?

**5. Page 3-5 Orthopedic Gloves**

Does the Agency plan to address the need for additional performance in use requirements (i.e.,

HEALTH AND HUMAN SERVICES  
FOOD AND DRUG ADMINISTRATION  
**CROSS REFERENCE SHEET**

Docket Number/Item Code: 99D-2335/C5

See Docket Number/Item Code: 98N-0313/C23

**fedEx** USA Airbill

FedEx  
Tracking  
Number

8146 4651 0880

Form  
I.D. No.

0215

From **This portion can be removed for Recipient's records** 814646510880  
Date **10-21-77** FedEx Tracking Number

Sender's Name **Carl Hermann** Phone **602 225-1070**

Company **SMART PRACTICE**

Address **3400 E MCDOWELL RD**

City **PHOENIX** State **AZ** ZIP **85008**

Your Internal Billing Reference

To Recipient's Name **Dr. L. L. Dear Alvin** Phone **301 514-5819**

Company **SECRET # 190-4335**

Address **2014 Airport Express, Phoenix, AZ**

We cannot deliver to PO boxes or PO ZIP codes

City **PHOENIX** State **AZ** ZIP **85022**

To "HOLD" at FedEx location,  
print FedEx address here.

City **PHOENIX** State **AZ** ZIP **85022**



**4a Express Package Service**

☐ FedEx Priority Overnight  
Next business morning

☒ FedEx Standard Overnight  
Next business afternoon

☐ FedEx First Overnight  
Earliest next business morning  
delivery to select locations

☐ FedEx 2Day\*  
Second business day

☐ FedEx Express Saver\*  
Third business day

\* FedEx Letter Rate not available  
Minimum charge: One-pound rate

**4b Express Freight Service**

☐ FedEx 1Day Freight\*  
Next business day

☐ FedEx 2Day Freight  
Second business day

☐ FedEx 3Day Freight  
Third business day

\* Call for Confirmation

Packages up to 150 lbs.  
Delivery commitment may be later in some areas.

Packages over 150 lbs.  
Delivery commitment may be later in some areas.

**5 Packaging**

☒ FedEx Letter\*

☐ FedEx Pak\*

☐ Other Pkg.  
Includes FedEx Box, FedEx  
Tube, and customer pkg.

\* Declared value limit \$500

**6 Special Handling**

☐ Saturday Delivery  
Available for FedEx Priority  
Overnight and FedEx 2Day  
to select ZIP codes

☐ Sunday Delivery  
Available for FedEx Priority  
Overnight to select ZIP codes

☐ HOLD Weekday  
at FedEx Location  
Not available with  
FedEx First Overnight

☐ HOLD Saturday  
at FedEx Location  
Available for FedEx Priority  
Overnight and FedEx 2Day  
to select locations

Does this shipment contain dangerous goods?

☒ No

☐ Yes

As per attached  
Shipper's Declaration

☐ Yes

Shipper's Declaration  
not required

☐ Dry Ice

Dry Ice, 9, UN 1845

x

kg

Dangerous Goods cannot be shipped in FedEx packaging.

☐ Cargo Aircraft Only

**7 Payment Bill to:**

☒ Sender  
Acct. No. in Section  
1 will be billed.

☐ Recipient

Enter FedEx Acct. No. or Credit Card No. below

☐ Third Party

☐ Credit Card

☐ Obtain Recip.  
Acct. No.

☐ Cash/Check

Total Packages

Total Weight

Total Charges

Credit Card Auth.

\*Our liability is limited to \$100 unless you declare a higher value. See the FedEx Service Guide for details.

**8 Release Signature**

Sign to authorize delivery without obtaining signature

By signing you authorize us to deliver this shipment without obtaining a signature  
and agree to indemnify and hold us harmless from any resulting claims.

Questions? Call 1-800-Go-FedEx (800-463-3339)

Visit our Web site at [www.fedex.com](http://www.fedex.com)

SRP 799 • Rev. Date 11/98 • Part #1548135 • © 1994-98 FedEx • PRINTED IN U.S.A.

359