



# Hip Implant Corrosion Mechanisms and Effects:

Mechanically Assisted Corrosion, Crevices and Voltage Effects

Jeremy L. Gilbert, Ph.D., FBSE
Professor of Biomaterials
Syracuse Biomaterials Institute
Department of Biomedical and Chemical Engineering
Syracuse University
Syracuse, NY 13244

# **Potential Conflicts**

- Consultant to and/or research contracts with
  - Medtronic
  - Stryker Orthopedics
  - Depuy Orthopedics
  - Biomet, Inc.
- Other Potential Conflicts
  - Editor-in-Chief, Journal of Biomedical Materials
     Research Part B: Applied Biomaterials
  - Past President, Society for Biomaterials

# Mechanically Assisted Corrosion

- Tribocorrosion (wear and corrosion)
- Fretting Corrosion (fretting and corrosion often in presence of a crevice)
  - Fretting = small scale cyclic motion (< 100 um) between two opposing surfaces</li>
- Stress Corrosion Cracking
- Stress Enhanced Corrosion
- Each is present in hip replacements where metals are stressed, abraded, worn
  or fretted
- Crevice-like environments add to the severity and complexity of the corrosion
- Fretting INITIATED Crevice Corrosion Fretting can produce the conditions for run-away corrosion (where fretting (or wear) is no longer required)
  - Fretting = Match
  - Crevice Corrosion = Fire

## **Some Basic Facts of Importance**


- You CANNOT HAVE WEAR of Metal Alloy Hip Implants WITHOUT CORROSION!
  - (But, you can have corrosion without wear)
- WEAR and CORROSION ARE COUPLED and NON-LINEAR
  - Wear -> Oxide Film Abrasion -> Corrosion -> Voltage Drop ->
    - -> Altered Oxide -> Altered Wear
- Oxide Films on Alloy Surfaces are Critical to Mechanically Assisted Corrosion Behavior and Corrosion Resistance
- With Corrosion: Voltage drops (cathodic excursions) and currents are generated which can affect the implant
  - depends on area abraded,
  - crevice geometries
  - and area available for reduction reactions
- Cathodic excursions have significant effects on the local biological reaction. It's not just about the metal ions and debris!

# Hip Prosthesis Corrosion Sources

**Examples of Current Modular Designs** 



Not indicative of anything corrosion related



Most severe corrosion observed at modular tapers Corrosion at one location affects all others (e.g., voltage drops)


## In-Vivo Corrosion of Modular Tapers: Head-Neck Tapers: 1993


Gilbert, et al., JBJS 1993

We have known since the late 1980's about modular implant corrosion (Svensson et al, 1988, JBJS(A), Fulminant pseudotumor with Co-Cr/CoCr modular taper)

All currently-used alloy combinations are known to be susceptible to attack (Ti, CoCrMo, 316L SS)

Passive oxide films on the surface are central to mechanism of attack.





Ti-6Al-4V/Co-Cr-Mo

#### Modular Body Interfaces: Ti6Al4V/Ti6Al4V Interfaces: 2009





Company A

Infection-related

22 months

Company B► Infection-related► Unknown -1?

Company C


→ Pain


→ 27 months


# Ti/Ti Modular Tapers:

Fretting Crevice Corrosion is STILL a problem










Pitting Corrosion and Hydrogen Embrittlement have been observed,




Fretting INITIATED Crevice Corrosion


#### Modular Neck Corrosion: 2012

In-Vivo Oxide-Induced Stress Corrosion Cracking in Ti-6Al-4V Modular Neck-stem Tapers

Medial proximal neck-stem taper cracking



# In-Vitro Fretting Corrosion Testing 316L SS Stem/CoCrMo Head Couples





# **Electrochemistry of Fretting Interfaces**

Metal ions, phosphate, chloride, hydrogen ions accumulate in crevice

Fretting currents consist of dissolution and repassivation reactions

Oxidation debris accumulates in crevice






2.  $M + H_2O \rightarrow MO + 2H^+ + 2e^-$  Excess electrons decrease voltage of implant


# Crevices and Spatially Separated Half-Cell Reactions


- Crevice and fretting localize oxidation
  - Oxide film repassivation
  - Metal ion dissolution
- Outside Crevice
  - Reduction (multiple, biologically based species available)
- Note: in crevice
  - Oxygen depletion
  - pH drops
  - Cl<sup>-</sup>, PO<sub>4</sub><sup>3-</sup> increases to balance M<sup>+</sup>
- Note: Fretting disrupts oxide in crevice and accelerates oxidation by 6 orders of magnitude.
- Build up of electrons drops voltage of surface




## Oxide Abrasion: High Speed Scratch Testing Ti- and CoCr Alloy

- 2.5 μA for 700 μm² results in 0.35
   A/cm² current density at scratch site
- 10<sup>6</sup> higher current densities than at passive oxides.







$$M \rightarrow M^{n+} + ne^{-}$$
  
 $mM + nH_2O \rightarrow M_mO_n + 2nH^+ + 2ne^{-}$ 

$$i_{total} = i_{film} + i_{dissolution} = \frac{\rho n F \upsilon}{M_{w}} \frac{d\theta}{dt} + i_{o} A e^{\frac{\eta}{\beta}} (1 - \theta)$$


Current response is voltage and alloy dependent Buckley

Buckley and Gilbert, 1994 Goldberg and Gilbert, 1997, 2004 Gilbert and Jacobs, ASTM, 1997


#### **Asperity-Based Model of Fretting Corrosion Currents**

$$I_{film} = 2 \frac{\rho nF}{M_w} \frac{A_{nom}}{\Delta} m(E - E^{onset}) \frac{d\delta}{dt}$$


Fretting currents resulting from oxide abrasion and repassivation




Note: Fretting 1 mm<sup>2</sup> will result in over 10 μA of current



Swaminathan and Gilbert, 2012, Biomaterials






Fretting Initiated Crevice Corrosion


## Reduction Reactions and Cell Viability (MC3T3-E1)

1 uA/cm<sup>2</sup> of reduction reaction kills cells

CoCrMo 24 hr at voltage, Haeri et al., 2012, Biomaterials








200

Reduction reactions at metal surfaces Kill osteoblast-like cells

Previously unknown mechanism of biological interaction


May play a role in hip prostheses in-vivo

It's not just about the metal ions



#### **Animations of Cell Behavior Under Voltage Control**





Vinculin Fluorescence Thresholded and Fourier Transform Filtered. (real time in upper left of each image)

#### Corrosion – Ion Level Estimate:

#### Can one estimate ion levels from corrosion currents?


Assume Steady State: Rate of generation of ions=rate of excretion via Urine

$$I = \sum_{i} \frac{n_{i} F}{Mw_{i}} \frac{dm_{i}}{dt}$$

$$I = \frac{dm}{dt} F \sum_{i} \frac{n_{i} \widehat{m}_{i}}{Mw_{i}}$$

$$\dot{m} = kI$$

$$k = \frac{1}{F \sum_{i} \frac{n_{i} \hat{m} i}{Mw_{i}}} = 2.18 \frac{\mu gs}{\mu C day}$$



$$\rho_m = \frac{kI}{V_u}$$

 $\rho_m$  is metal ion levels in urine, I is the corrosion current,  $V_u$  is the volume of urine excreted per day (2 L)

$$Co \rightarrow Co^{2+} + 2e^{-}$$
  
 $Cr \rightarrow Cr^{3+} + 3e^{-}$ 

This estimate (and its only an estimate) gives a sense of the relationship between amount of corrosion required to result in a systemic ion level

# Conclusions

- Mechanically assisted corrosion (MAC) continues to be a serious concern for metallic biomaterials in all applications
- Wear and Corrosion are coupled and interactive
- Fretting can INITIATE crevice corrosion
- All current alloy systems are susceptible to MAC
- Negative voltage excursions result from MAC and may lead to adverse biological responses (apoptosis in-vitro)

# Acknowledgements

#### Supported in part by:

- Medtronic, Stryker, Depuy
- To date, no research grants were provided to the author by the federal government for any of the work presented.

#### Collaborators

- Joshua J. Jacob, M.D., Robert Urban, Rush Medical University
- Eugene P. Lautenschlager, Ph.D., Northwestern University
- George Langford, Ph.D., Syracuse University
- Torsten Wollert, Ph.D., Syracuse University

#### Students:

- Christine Buckley, Ph.D.
- Jay Goldberg, Ph.D.
- Spiro Megremis, Ph.D.
- Mark Ehrensberger, Ph.D.
- Shiril Sivan, Ph.D. (expected 2013)
- Morteza Haeri, Ph.D. (expected 2012)
- Viswanathan Swaminathan, Ph.D. (expected 2012)