



**Impax Generics  
Epinephrine Injection, USP Auto-Injector  
Patient Inspection Instructions**

*PATIENT INSPECTION INSTRUCTIONS  
Epinephrine Injection, USP auto-injector  
0.3 mg, 0.15 mg*

*This is to remind you to look at the epinephrine solution through the viewing window of the auto-injector before use. If the solution is discolored (pinkish or brown color), cloudy, or contains particles, do not use the auto-injector and return it to your pharmacy for replacement. Please contact Impax Generics at 1-800-934-6729 if you have questions about this information.*