



Before you get started:

- Your sample must be dropped off on the same day you activate and collect it.
- **Make sure there is a collection today.** Refer to your program details if you're not sure where to return your sample.
- **Don't miss the last pickup.** Refer to your program details for pickup times.
- Make sure you wash your hands before and after collection.

Your kit includes:



Collection device
(tube and swab)



Biohazard bag with
Absorbent Pad



Barcode
card

Checklist before returning your kit:

- Did you activate your kit?
- Did you collect a sample from both nostrils?

Visit vaulthealth.com/covid/consumer#faq to learn more.

Questions?

Visit vaulthealth.com/covid/consumer#faq or contact us at covid@vaulthealth.com or (212) 880-5494.

This product (collection kit in combination with the authorized test) has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

RhinoStics Nasal Swab Collection Kit

for Nexus High Throughput
SARS-CoV-2 Assay

On-site Instructions

For In Vitro Diagnostic Use, For Prescription Use Only, For Use Under Emergency Use Authorization Only
For use by individuals 18 years and older.

Version 1.0

1

Don't forget:
Activate your kit online
at pcr.vlt.co.

You need your unique barcode,
which you'll find on your collection
tube.



Warning: Your sample will not be
tested if the kit is not activated.

2


Before collecting your sample,
wash your hands thoroughly for
20 seconds and dry them well.

3



Unscrew the cap from the
collection tube.
You can discard the cap.

4



Remove the swab from
its bag.
Do not touch the ridged tip
with your hands or lay it
down on any surface.


5



Rotate swab tip in
first nostril, 3 times.
Insert the swab into one nostril, until you
feel resistance. Using medium pressure
against the inside of your nostril the whole
time, rub the swab in a large circular
motion at least 3 times.


Note: Swab should not go further than 1 inch.

6




Gently slide the swab up and
down against the inside of your
nostril 3 times.
Keep medium pressure on the swab against
the inside of your nostril.

7



Hold the swab against
the inside of your nostril.
Keep it there for about
10 seconds.


8



Repeat steps 5 to 7 in the
other nostril.
Rotate at least 3 times, slide the
swab up and down, then hold
the swab against your nostril for
10 seconds.

Note: You will use the same end
of the swab for both nostrils.

9



Put swab into the
tube portion of the
collection device.
The tip of the swab that went
into your nose should go into
the tube first.

10



Screw the cap closed tightly.
Make sure the top is carefully
screwed closed.
Note: After screwing the collection
device closed, wash your hands
thoroughly for 20 seconds and dry
them well.

11



Put collection device
into the biohazard bag.
Seal the bag closed.


12

Bring the biohazard bag to your
designated test collection point
on the day you collected your
sample.
Refer to your program details if you're
not sure where to return your sample.
Note: Don't miss the last pickup time.


Before you get started:

- Your sample must be dropped off on the same day you activate and collect it.
- **Make sure there is a collection today.** Refer to your program details if you're not sure where to return your sample.
- **Don't miss the last pickup.** Refer to your program details for pickup times.
- Make sure you wash your hands before and after collection.

Your kit includes:



Collection device
(tube and swab)



Biohazard bag with
Absorbent Pad



Barcode
card

Checklist before returning your kit:

- Did you activate your kit?
- Did you collect a sample from both nostrils?

Visit color.com/sample-return to learn more.

Questions?

Visit support.color.com or contact us at mycovidtest@color.com or (844) 352-6567.

This product (collection kit in combination with the authorized test) has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. §360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

RhinoStics Nasal Swab Collection Kit

for Nexus High Throughput SARS-CoV-2 Assay

On-site Instructions

For In Vitro Diagnostic Use, For Prescription Use Only, For Use Under Emergency Use Authorization Only
For use by individuals 18 years and older.

Version 1.0

1

Don't forget:
Activate your kit online at
color.com/covid/activate.
You need your unique barcode,
which you'll find on your collection
tube.



Warning: Your sample will not be
tested if the kit is not activated.

2


Before collecting your sample,
wash your hands thoroughly for
20 seconds and dry them well.

3



Unscrew the cap from the
collection tube.
You can discard the cap.

4



Remove the swab from
its bag.
Do not touch the ridged tip
with your hands or lay it
down on any surface.


5



Rotate swab tip in
first nostril, 3 times.
Insert the swab into one nostril, until you
feel resistance. Using medium pressure
against the inside of your nostril the whole
time, rub the swab in a large circular
motion at least 3 times.


Note: Swab should not go further than 1 inch.

6




Gently slide the swab up and
down against the inside of your
nostril 3 times.
Keep medium pressure on the swab against
the inside of your nostril.

7



Hold the swab against
the inside of your nostril.
Keep it there for about
10 seconds.


8



Repeat steps 5 to 7 in the
other nostril.
Rotate at least 3 times, slide the
swab up and down, then hold
the swab against your nostril for
10 seconds.

Note: You will use the same end
of the swab for both nostrils.

9



Put swab into the
tube portion of the
collection device.
The tip of the swab that went
into your nose should go into
the tube first.

10



Screw the cap closed tightly.
Make sure the top is carefully
screwed closed.
Note: After screwing the collection
device closed, wash your hands
thoroughly for 20 seconds and dry
them well.

11



Put collection device
into the biohazard bag.
Seal the bag closed.



12

Bring the biohazard bag to your
designated test collection point
on the day you collected your
sample.
Refer to your program details if you're
not sure where to return your sample.
Note: Don't miss the last pickup time.

Before you get started:

- Your sample must be dropped off on the same day you activate and collect it.
- Refer to your program details if you're not sure where to return your sample.
- Don't miss the last pickup. Refer to your program details for pickup times.
- Make sure you wash your hands before and after collection.

Your kit includes:



Checklist before returning your kit:

- Did you activate your kit?
- Did you collect a sample from both nostrils?

Visit vaulthealth.com/covid/consumer#faq to learn more.

Questions?

Visit vaulthealth.com/covid/consumer#faq or contact us at covid@vaulthealth.com or (212) 880-5494.

This product (collection kit in combination with the authorized test) has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Rhinostics Nasal Swab Collection Kit

for Nexus High Throughput SARS-CoV-2 Assay

Mail-In Instructions

For In Vitro Diagnostic Use, For Prescription Use Only, For Use Under Emergency Use Authorization Only
For use by individuals 18 years and older.

Version 1.0

1

Don't forget:
Activate your kit online at pcr.vlt.co.

You need your unique barcode, which you'll find on your collection tube.



Warning: Your sample will not be tested if the kit is not activated.

2

Before collecting your sample, wash your hands thoroughly for 20 seconds and dry them well.



3

Unscrew the cap of the collection tube.

You can discard the cap.

4

Remove the swab from its bag.

Do not touch the soft tip with your hands or lay it down on any surface.



5

Rotate swab tip in first nostril, 3 times.

Insert the swab into one nostril, until you feel resistance. Using medium pressure against the inside of your nostril the whole time, rub the swab in a large circular motion at least 3 times.

Note: Swab should not go further than ½ inch.



6

Gently slide the swab up and down against the inside of your nostril 3 times.

Keep medium pressure on the swab against the inside of your nostril.



7

Hold the swab against the inside of your nostril.

Keep it there for about 10 seconds.



8

Repeat steps 5 to 7 in the other nostril.

Rotate at least 3 times, slide the swab up and down, then hold the swab against your nostril for 10 seconds.

Note: You will use the same end of the swab for both nostrils.



9

Put swab into the tube portion of the collection device and screw cap closed tightly.

The tip of the swab that went into your nose should go into the tube first. **Make sure the cap is closed tightly.**

Note: After screwing the collection device closed, wash your hands thoroughly for 20 seconds and dry them well.



10

Put the collection tube into the biohazard bag.

Seal the biohazard bag closed.

Note: The collection tube should go within the rigid plastic insert in the biohazard bag.




11

Put biohazard bag into the shipping box.

Seal the shipping box closed.

Note: Use the same box that your kit came in.



12

Put box into the return shipping mailer.

Seal the return shipping mailer closed. Your return shipping has already been paid.

Note: This packaging is specially designed for shipping samples suspected of carrying infectious substances. Do NOT replace it.



13

Return kit.



Your sample must be returned on the same day you collect it.

Note: If you aren't able to return your sample on the same day, contact covid@vaulthealth.com.

Before you get started:

- Your sample must be dropped off on the same day you activate and collect it.
- Refer to your program details if you're not sure where to return your sample.
- Don't miss the last pickup. Refer to your program details for pickup times.
- Make sure you wash your hands before and after collection.

Your kit includes:



Checklist before returning your kit:

- Did you activate your kit?
- Did you collect a sample from both nostrils?

Visit color.com/sample-return to learn more.

Questions?

Visit support.color.com or contact us at mycovidtest@color.com or (844) 352-6567.

This product (collection kit in combination with the authorized test) has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by the authorized laboratory.

This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Rhinostics Nasal Swab Collection Kit

for Nexus High Throughput SARS-CoV-2 Assay

Mail-In Instructions

For In Vitro Diagnostic Use, For Prescription Use Only, For Use Under Emergency Use Authorization Only
For use by individuals 18 years and older.

Version 1.0

1

Don't forget:
Activate your kit online at
color.com/covid/activate.

You need your unique barcode, which you'll find on your collection tube.



Warning: Your sample will not be tested if the kit is not activated.

2

Before collecting your sample, wash your hands thoroughly for 20 seconds and dry them well.



3

Unscrew the cap of the collection tube.

You can discard the cap.

4

Remove the swab from its bag.

Do not touch the soft tip with your hands or lay it down on any surface.




5

Rotate swab tip in first nostril, 3 times.

Insert the swab into one nostril, until you feel resistance. Using medium pressure against the inside of your nostril the whole time, rub the swab in a large circular motion at least 3 times.

Note: Swab should not go further than 1/2 inch.



6

Gently slide the swab up and down against the inside of your nostril 3 times.

Keep medium pressure on the swab against the inside of your nostril.



7

Hold the swab against the inside of your nostril.

Keep it there for about 10 seconds.




8

Repeat steps 5 to 7 in the other nostril.

Rotate at least 3 times, slide the swab up and down, then hold the swab against your nostril for 10 seconds.

Note: You will use the same end of the swab for both nostrils.



9

Put swab into the tube portion of the collection device and screw cap closed tightly.

The tip of the swab that went into your nose should go into the tube first. **Make sure the cap is closed tightly.**

Note: After screwing the collection device closed, wash your hands thoroughly for 20 seconds and dry them well.



10

Put the collection tube into the biohazard bag.

Seal the biohazard bag closed.

Note: The collection tube should go within the rigid plastic insert in the biohazard bag.




11

Put biohazard bag into the shipping box.

Seal the shipping box closed.

Note: Use the same box that your kit came in.



12

Put box into the return shipping mailer.

Seal the return shipping mailer closed. Your return shipping has already been paid.

Note: This packaging is specially designed for shipping samples suspected of carrying infectious substances. Do NOT replace it.



13

Return kit.

Your sample must be returned on the same day you collect it.

Note: If you aren't able to return your sample on the same day, contact mycovidtest@color.com