

The SleepWeaver Advance Pediatric PAP Mask

Industry Partner: Circadiance, LLC

The SleepWeaver Advance Pediatric PAP Mask was developed by Circadiance, LLC for use in pediatric patients requiring positive airway pressure.

This Class 2 device provides an interface for noninvasive Positive Airway Pressure (PAP) ventilation therapy; and is intended for single-patient use in the home and single-patient use in the hospital/ institutional environment. The device's cloth material was designed to protect a child's skin from developing pressure ulcers from prolonged contact with traditionally rigid masks.

