



BIONETICS

Final report-Mutagenicity Evaluation of FDA 75-68 (Calcium Glycerophosphate)
5/77

MUTAGENICITY EVALUATION
OF

FDA 75-68
CALCIUM GLYCEROPHOSPHATE
000126-95-4

FINAL REPORT

MIS

5546 Nicholson Lane
Kensington, Maryland
20755

MUTAGENICITY EVALUATION
OF
FDA 75-68
CALCIUM GLYCEROPHOSPHATE
000126-95-4
FINAL REPORT

SUBMITTED TO

DEPARTMENT OF HEALTH, EDUCATION AND WELFARE
FOOD AND DRUG ADMINISTRATION
NEGOTIATED CONTRACTS BRANCH
5600 FISHERS LANE, HFA-510, 5B-37
ROCKVILLE, MARYLAND 20852

SUBMITTED BY

LITTON BIONETICS, INC.
5516 NICHOLSON LANE
KENSINGTON, MARYLAND 20795

LBI PROJECT NO. 2672

MAY, 1977



BIONETICS

TABLE OF CONTENTS

	Page No.
EVALUATION SUMMARY.....	1
I. <u>OBJECTIVE</u>	2
II. <u>MATERIALS</u>	2
A. Test Compound.....	2
B. Indicator Microorganisms.....	2
C. Reaction Mixture.....	2
D. Tissue Homogenates and Supernatants.....	3
E. Positive Control Compounds.....	3
III. <u>METHODS</u>	3
A. Toxicity.....	3
B. Plate Tests.....	4
C. Suspension Tests.....	4
D. Preparation of Tissue Homogenates and 9,000 x g Cell Fractions.....	5
E. Data Recording and Reporting.....	5
IV. <u>RESULTS SECTION</u>	
A. Solubility Properties of the Test Compound.....	6
B. Toxicity and Dosage Determinations for the Test Compound.....	6
C. Plate Assay Results.....	7
D. Suspension Assay Results.....	7
V. <u>INTERPRETATION OF RESULTS AND CONCLUSIONS</u>	15
VI. <u>EXPLANATION OF EVALUATION PROCEDURES FOR PLATE ASSAYS</u>	17
VII. <u>EXPLANATION OF EVALUATION PROCEDURES FOR SUSPENSION ASSAYS</u>	19
APPENDIX - Tabulation of Data.....	A-1



EVALUATION SUMMARY

The test compound, FDA 75-68, Calcium Glycerophosphate, 000126-95-4, did not exhibit mutagenic activity in any of the assays employed in these studies.



BIONETICS

DATE: May 10, 1977

SPONSOR: U.S. Food and Drug Administration

SUBJECT: Evaluation of Test Compound FDA 75-68 Calcium Glycerophosphate
000126-95-4

I. OBJECTIVE

The objective of this study was to evaluate the test compound for genetic activity in microbial assays with and without the addition of mammalian metabolic activation preparations.

II. MATERIALS

A. Test Compound

1. Date Received: October 29, 1976
2. Description: Fine white crystalline powder

B. Indicator Microorganisms

The following strains of indicator microorganisms were used in the evaluation:

Yeast Strain: Saccharomyces cerevisiae, strain D4

Bacteria Strains: Salmonella typhimurium, strains TA-1535
TA-1537
TA-1538
TA-98
TA-100

C. Reaction Mixture

The following reaction mixture was employed in the activation tests:

<u>Component</u>	<u>Final Concentration/ml</u>
1. TPN (sodium salt)	4 μ moles
2. Glucose-6-phosphate	5 μ moles
3. Sodium phosphate (dibasic)	100 μ moles
4. $MgCl_2$	8 μ moles
5. KCl	33 μ moles
6. Homogenate fraction equivalent to 25 mg of wet tissue.	

D. Tissue Homogenates and Supernatants

The tissue homogenates and 9,000 x g supernatants were prepared from tissues of the following mammalian species: Mouse - ICR random bred adult males; rat - Sprague-Dawley adult males; and monkey - Macaca mulatta adult males.

E. Positive Control Compounds

Table 1 lists chemicals for positive controls in the direct and activation assays.

TABLE 1
POSITIVE CONTROLS USED IN DIRECT AND ACTIVATION ASSAYS

<u>Assay</u>	<u>Chemical^a</u>	<u>Solvent</u>	<u>Probable Mutagenic Specificity</u>
Nonactivation	Methylnitrosoguanidine	Water or saline	BPS ^b
	Ethylmethanesulfonate	Water or saline	BPS ^b
	2-Nitrofluorene	Dimethylsulfoxide ^c	FS ^b
	Quinacrine mustard	Water or saline	FS ^b
Activation	Dimethylnitrosamine	Water or saline	BPS ^b
	2-Acetylaminofluorene	Dimethylsulfoxide ^c	FS ^b
	8-Aminoquinoline	Dimethylsulfoxide ^c	FS ^b
	2-Aminoanthracene	Dimethylsulfoxide ^c	BPS ^b

^a Concentrations given in the Results Section

^b BPS = base-pair substitution; FS = frameshift

^c Previously shown to be non-mutagenic

III. METHODS

A. Toxicity

The solubility, toxicity and doses for the test chemical were determined prior to screening.

The test chemical was tested for toxicity against specific indicator strains over a range of doses to determine the 50% survival dose. Bacteria were tested in phosphate buffer, pH 7.4, for one hour at 37°C on a shaker. Yeasts were tested in phosphate buffer, pH 7.4, for four hours at 30°C on a shaker. The 50% survival concentrations and the 1/4 and 1/2 50% doses calculated.

If no toxicity was obtained for the chemical with a given strain, then a maximum dose of 5% (w/v) was used.

Unless otherwise specified, the doses calculated for the tests in buffer were applied to the activation tests. The solubility of the test chemical under treatment conditions is stated in the Results Section.

B. Plate Tests (Overlay Method)

Approximately 10^8 cells from an overnight culture of each indicator strain were added to test tubes containing 2.0 ml of molten agar supplemented with biotin and a trace of histidine. For nonactivation tests, the three dose levels of the test compound were added to the contents of the appropriate tubes and poured over the surfaces of selective agar plates. In activation tests 0.5 ml of a $9,000 \times g$ tissue supernatant and required cofactors (core reaction mixture) were added to the overlay tubes. Three dose levels of the test chemical were added to the appropriate tubes, which were then mixed and the contents poured over the surface of a minimal agar (selective medium) plate and allowed to solidify. The plates were incubated for 48 to 72 hours at 37°C , and scored for the number of colonies growing on each plate. The concentrations of all chemicals are given in the Results Section. Positive and solvent controls using positive compounds that are active directly and those that require metabolic activation were run with each assay.

C. Suspension Tests

1. Nonactivation

Bacteria and yeast cultures of the indicator organisms were grown in complete broth, washed and resuspended in 0.9% saline to densities of 1×10^{10} cells/ml and 5×10^9 cells/ml, respectively. This constituted the working stock for tests of a group of test chemicals and their respective controls. Tests were conducted in plastic, 24-well tissue culture plates (Linbro). Cells plus appropriate volume(s) of the test chemical were added to the wells to give a final volume of 1.5 ml. The solvent replaced the test chemical in the negative controls. Treatment was at 30°C for four hours for yeast tests and at 37°C for one hour for bacterial tests. All flasks were shaken during treatment. Following treatment, the plates were set on ice. Aliquots of cells were removed, diluted in sterile saline (4°C) and plated on the appropriate complete media. Undiluted samples from flasks containing the bacteria were plated on minimal selective medium in reversion experiments. Samples from a 10^{-1} dilution of treated cells were plated on the selected media for enumeration of gene conversion with strain D4. Bacterial plates were scored after incubation for 48 hours at 37°C . The yeast plates were incubated at 30°C for 3-5 days before scoring.

2. Activation

Bacteria and yeast cells were grown and prepared as described in the nonactivation tests. Measured amounts of the test and control chemicals plus 0.25 ml of the stock-cell suspension were added to wells of the Linbro plate containing the appropriate tissue fraction and reaction mixture. All flasks (bacteria and yeast) were incubated at 37°C with shaking. The treatment times as well as the dilutions, plating procedures and scoring of the plates were the same as described for nonactivation tests.

D. Preparation of Tissue Homogenates and 9,000 x g Cell Fractions

Male animals (except monkeys) sufficient to provide the necessary quantities of tissues were killed by cranial blow, decapitated and bled. Monkey tissues were obtained from freshly killed and bled male rhesus monkeys. Organs were immediately dissected from the animals using aseptic techniques and placed in ice-cold 0.15M KCl. Upon collection of the desired quantity of organs, they were washed twice with fresh KCl and completely homogenized with a motor-driven homogenizing unit at 4°C. The whole organ homogenate obtained from this step was divided into two samples. One sample was frozen at -80°C and the other was centrifuged for 20 minutes at 9,000 x g in a refrigerated centrifuge. The supernatant from the centrifuged sample was retained and frozen at -80°C. These two frozen samples were used for the activation studies. Protein and P-448 determinations were made for each lot of homogenate.

E. Data Recording and Reporting

1. Plate test assays

The numbers of colonies on each plate were counted and recorded on printed forms. These raw data were entered into a computer program designed to print out all data by test. The data are presented as revertants per plate for each indicator strain employed in the assay. The positive and solvent controls are provided as reference points.

2. Suspension assays

Following the specified incubation periods all population plates were scored by an automatic colony counter and the results from each plate of a set were recorded, in ink, on data processing forms. All minimal or other types of selective media plates were hand scored and the results recorded along with the respective population data. Other relevant experimental data were recorded on experimental definition forms. For bacteria strains the number of colonies recorded from either the population or selective plates represents that number in 1 ml of test suspension plated. The numbers recorded for the yeast strain D4 represent the number in 0.5 ml of test suspension plated. The data were then processed and printed from a computer program. All raw data sheets are dated and signed by the responsible technician.



IV. RESULTS SECTION

A. Solubility Properties of the Test Compound

1. Name or code designation of the test compound: 000126-95-4
Calcium Glycerophosphate
FDA 75-68
2. Test solvent: Saline
3. Solubility of the test compound under treatment conditions: Soluble
4. Additional comments: Fine white crystalline powder

B. Toxicity and Dosage Determinations for the Test Compound

1. Test date for toxicity determination: November 22, 1976
2. The 50% survival level was determined for bacteria and yeast indicator organisms by conducting survival curves with the test compound at the following concentrations:

Percent Concentration (w/v or v/v)

5.0
0.5
0.05
0.005
0.0005

3. Concentrations of the test compound used in the mutagenicity tests:

<u>Test Doses</u>	<u>Percent Concentration</u>	
	<u>Bacteria</u>	<u>Yeast</u>
1/4 50% Survival	1.25	0.3575
1/2 50% Survival	2.50	0.7150
50% Survival	5.00	1.4300

C. Plate Test Results

The plate test results are summarized in the following table. The values presented in this table are the number of revertants per plate.

D. Suspension Assay Results

The suspension test results for the test compound are summarized in the tables following the plate test summary. The values presented in these tables are the calculated mutation frequencies for each control and experimental test point. The first table of the suspension set presents the results for the nonactivation assays, and the second table through the fourth table of the suspension set presents the results for the activation assays. A listing of computer codes and abbreviations is included for reference. Tabulation of all raw data is provided in the Appendix.

SUMMARY OF TEST RESULTS

PLATE TESTS

A. NAME OR CODE DESIGNATION OF THE TEST COMPOUND: 000126954

B. TEST DATE: DEC. 16, 1976

TEST	SPECIES	ISSUE	BEVERIANIS PER PLATE									
			TA-1535		TA-1537		TA-1538		TA-98		TA-100	
			1	2	1	2	1	2	1	2	1	2
1. NON-ACTIVATION												
SOLVENT CONTROL*	---	---	24	34	17	15	18	21	53	47	63	141
POSITIVE CONTROL**	---	---	>1000	>1000	>1000	>1000	>1000	>1000	>1000	>1000	>1000	>1000
TEST 5.00000 %	---	---	14	17	28	20	31	15	40	36	61	52
2.50000 %	---	---	16	11	19	22	24	29	39	41	60	73
1.25000 %	---	---	28	13	47	33	35	28	50	38	66	77
2. ACTIVATION												
SOLVENT CONTROL*	MOUSE	LIVER	32	28	25	30	34	32	47	49	160	163
	RAT	LIVER	30	29	29	34	23	19	38	47	219	195
	MONKEY	LIVER	26	32	18	23	29	40	42	45	176	199
POSITIVE CONTROL***	MOUSE	LIVER	793	851	252	297	>1000	>1000	>1000	>1000	>1000	>1000
	RAT	LIVER	699	590	264	357	719	894	>1000	>1000	>1000	>1000
	MONKEY	LIVER	539	446	155	181	514	525	>1000	940	>1000	>1000
TEST 5.00000 %	MOUSE	LIVER	23	44	17	21	17	9	41	35	55	44
2.50000 %	MOUSE	LIVER	33	47	31	22	12	15	53	39	50	81
1.25000 %	MOUSE	LIVER	15	23	14	19	21	15	37	44	112	142
5.00000 %	RAT	LIVER	19	20	28	22	9	13	36	27	89	125
2.50000 %	RAT	LIVER	30	29	35	26	10	13	25	37	86	115
1.25000 %	RAT	LIVER	23	31	24	14	36	43	38	53	166	232
5.00000 %	MONKEY	LIVER	27	36	20	18	15	13	47	40	64	109
2.50000 %	MONKEY	LIVER	32	28	26	17	10	8	44	40	63	97
1.25000 %	MONKEY	LIVER	16	27	24	21	17	15	35	50	107	171

* NON-ACTIVATION ASSAYS CONSIST OF THE CELLS PLUS THE TEST COMPOUND VEHICLE (SOLVENT). FOR ACTIVATION ASSAYS, THE OVERLAY CONTAINS THE ACTIVATION SYSTEM PLUS THE TEST COMPOUND VEHICLE.

** TA-1535 MNNG 2 UG/PLATE
TA-1537 QM 20 UG/PLATE
TA-1538 NF 100 UG/PLATE
TA-98 NF 100 UG/PLATE
TA-100 MNNG 2 UG/PLATE

*** TA-1535 ANTH 100 UG/PLATE
TA-1537 AMQ 100 UG/PLATE
TA-1538 AAF 100 UG/PLATE
TA-98 AAF 100 UG/PLATE
TA-100 ANTH 100 UG/PLATE

NOTE: CONCENTRATIONS ARE GIVEN IN MICROLITERS(UL) OR MICROGRAMS(UG) PER PLATE.

LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
REPORT EXR34

COMPOUND FREQUENCY SUMMARY REPORT 03/14/77

NONACTIVATION COMPOUND 000126954

TEST	ORG	TA100 HIS EX-8	TA1535 HIS EX-8	TA1537 HIS EX-8	TA1538 HIS EX-9	TA98 HIS EX-8	TA98 HIS EX-8	000004 ADE EX-5	000004 TRY EX-5	
NAN		56.09	5.34	20.06	10.96	3.79	13.05	19.75	5.37	CONTROLS
NAP		164.41	507.36	137.92	187.12	141.25		101.54	67.35	
<hr/>										
NA1		11.94	8.37	6.53	5.83	3.53	17.30	15.46	3.84	TEST DATA
NA2		21.97	2.85	11.75	8.06	2.68		7.60	6.86	
NA3		21.78	11.54	19.47	10.33	2.50		11.82	5.01	

LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
REPORT EXR34

COMPOUND FREQUENCY SUMMARY REPORT 03/14/77

SPECIES ICRFLO/MOUSE COMPOUND 000126954

TEST	ORG	TA100 HIS EX-8	TA1535 HIS EX-8	TA1537 HIS EX-8	TA1538 HIS EX-8	TA1538 HIS EX-8	TA98 HIS EX-8	0000D4 ADE EX-5	0000D4 TRY EX-5	
ACT	A+C	184.08	9.01	5.44		7.54	2.55	14.04	11.23	NEGATIVE CONTROLS
ACT	A-C	70.56	8.33	4.77		10.87	2.87	16.96	13.67	
ACT	ALI	109.86	11.11	4.68	3.79	9.73	5.80	19.69	12.92	
ACT	ALU	78.31	8.72	5.66		10.04	4.99	20.41	15.61	
<hr/>										
ACT	PLI	251.33	179.78	141.85		222.79	58.00	76.42	76.94	POSITIVE CONTROLS
ACT	PLU	94.27	11.78	7.13		67.40	24.95	18.36	11.01	
<hr/>										
ACT	L11	51.42	9.95	5.22		11.17	3.58	11.60	8.46	TEST COMPOUND
ACT	L12	36.17	10.23	7.20		13.43	9.11	10.79	7.80	
ACT	L13	34.97	6.65	6.21	4.40	00.00	14.73	11.76	6.31	
ACT	LU1	32.34	8.78	4.62		17.39	9.24	15.12	7.72	
ACT	LU2	71.15	4.85	5.89		14.90	6.84	10.53	4.34	
ACT	LU3	48.11	9.00	6.99		15.94	5.32	13.98	9.13	

LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
REPORT EXR34

COMPOUND FREQUENCY SUMMARY REPORT 03/14/77

SPECIES SPRDAW/RAT

COMPOUND 000126954

TEST	ORG	TA100 HIS EX-8	TA1535 HIS EX-8	TA1535 HIS EX-8	TA1537 HIS EX-8	TA1538 HIS EX-8	TA98 HIS EX-8	0000D4 ADE EX-5	0000D4 TRY EX-5	
ACT	A+C	161.00	7.00		3.77	46.02	2.66	13.64	10.19	NEGATIVE CONTROLS
ACT	A-C	128.57	7.78		1.60	9.88	6.37	16.00	11.09	
ACT	ALI	128.30	2.66		2.05	13.54	22.50	18.91	9.31	
ACT	ALU	111.58	3.85	8.11	2.85	12.30	9.64	17.03	8.67	
ACT	PLT	173.97	138.63		89.01	128.54	299.56	97.06	72.40	POSITIVE CONTROLS
ACT	PLU	112.05	4.45		8.17	32.58	49.31	10.15	3.16	
ACT	L11	48.61	4.94		4.34	14.34	11.48	12.96	7.83	TEST COMPOUND
ACT	L12	71.38	4.86		4.10	10.91	9.01	10.12	8.63	
ACT	L13	48.25	2.86		2.50	10.04	9.69	10.00	6.52	
ACT	LU1	50.44	8.28	4.05	3.37	8.93	10.71	11.74	7.25	
ACT	LU2	60.69	5.66	11.33	2.44	6.64	4.50	11.69	5.70	
ACT	LU3	51.78	2.67	12.17	3.51	10.89	6.12	14.24	6.78	

LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
REPORT EXR34

COMPOUND FREQUENCY SUMMARY REPORT 03/14/77

SPECIES RHESUS/MONKEY COMPOUND 000126954

TEST	ORG	TA100 HIS EX-8	TA1535 HIS EX-8	TA1535 HIS EX-8	TA1537 HIS EX-8	TA1538 HIS EX-8	TA98 HIS EX-8	0000D4 ADE EX-5	0000D4 TRY EX-5	
ACT	A+C	53.27	2.58		2.45	9.78	5.79	19.41	10.65	NEGATIVE CONTROLS
ACT	A-C	57.26	2.21		2.44	13.84	5.50	18.48	9.49	
ACT	ALI	82.47	3.39	8.45	2.07	12.29	9.56	16.49	10.91	
ACT	ALU	72.94	2.78		1.10	11.95	10.63	15.33	10.17	
<hr/>										
ACT	PLI	195.14	36.71		81.50	137.05	123.31	93.48	73.43	POSITIVE CONTROLS
ACT	PLU	101.96	2.37		4.66	15.42	17.43	25.84	17.05	
<hr/>										
ACT	L11	90.75	3.51		1.73	10.42	4.93	17.13	7.01	TEST COMPOUND
ACT	L12	133.03	3.07		2.48	5.23	6.60	17.64	9.87	
ACT	L13	130.53	1.85		2.18	11.38	12.54	17.56	6.54	
ACT	LU1	126.61	1.70		1.37	7.49	1.78	18.69	9.90	
ACT	LU2	109.25	1.07	4.90	1.49	16.20	3.22	18.49	8.22	
ACT	LU3	72.75	4.73		0.92	10.13	3.46	15.64	12.42	

DATA TABLE TERMS AND ABBREVIATIONS

ABBREVIATION OR TERM	DEFINITION OR EXPLANATION
COMPOUND	Client designated compound number appears in this column.
TEST CODES	<p> NAN = Nonactivation: Solvent Control NAP = Nonactivation: Positive Control NA1 = Nonactivation: Test Compound Dose 1 NA2, etc. = Reflects the other dose level(s) </p> <p> A+C = Negative Chemical Control for ACP A-C = Activation: Solvent Control ALI or A+T = Activation: Homogenate Control (Liver) ALU = Activation: Homogenate Control (Lung) ACP = Activation: Positive Control ACT = Activation Test </p> <p> LI = Liver Tissue Activation Fraction LU = Lung Tissue Activation Fraction KI = Kidney Tissue Activation Fraction TE = Testes Tissue Activation Fraction 1,2, etc. = Dose Levels </p>
CONCENTRATION	<p>All test compound dose levels are expressed as a whole number followed by an exponent (negative) identified by the appropriate units.</p> <p>Example: 0025-2PCT = 0.25 percent concentration</p>
POPU	Total number of viable cells in the plating sample raised to some exponent printed directly below the abbreviation (i.e., EP + 6 = $x \times 10^6$).
MUT 1	Total number of mutants or convertants obtained from the sample plated raised to some exponent printed directly below the abbreviation (i.e., EP + 0 = 10^0). For strain D4, MUT 1 represents the number of ADE+ convertants.
MUT 2	Only used for strain D4 and represents the number of TRY+ convertants in the plated sample.
FREQ 1	The calculated mutation or gene conversion frequency times the negative exponent written directly below. For strain D4, FREQ 1 represents the ADE+ value.
FREQ 2	Only used for strain D4 and represents the TRY+ conversion frequency.
CONTAM	Presence of contamination on any plates.

DATA TABLE TERMS AND ABBREVIATIONS (continued)

ABBREVIATION OR TERM	DEFINITION OR EXPLANATION
AAF	2-Acetylaminofluorene
DMSO	Dimethylsulfoxide
DMN	Dimethylnitrosamine
EMS	Ethylmethanesulfonate
QM	Quinacrine Mustard
NF	Nitrofluorene
ANTH	2-Amino Anthracene
AMQ	8-Amino Quinoline
SPECIES	Animal Strains
SPRDAW	Sprague Dawley Rats
ICRFLO	Flow ICR Random Bred Mice
RHESUS	Rhesus Monkey (<u>Macaca mulatta</u>)
MIXEDB	Dog, Mixed Breed
NEWZEA	New Zealand White Rabbit
UG	Microgram
UM	Micromole
ADE	Adenine
TRY	Tryptophan



BIONETICS

Litton

V. INTERPRETATION OF RESULTS AND CONCLUSIONS

The test compound, FDA 75-68, Calcium Glycerophosphate, 000126-95-4, was evaluated for genetic activity in a series of in vitro microbial assays with and without metabolic activation. The following results were obtained:

A. Salmonella typhimurium

1. Plate tests

The results of these tests were negative.

2. Nonactivation suspension tests

The results of these tests were negative. The following tests were repeated: NA1 dose with TA-98 because low population count was observed with this dose in the initial experiment (less than 25% of the solvent control).

3. Activation suspension tests

The results of these tests were negative. The following tests were repeated: (a) LI3 dose with TA-1538 using mouse tissue because no colonies were observed on the bacterial plates which is attributed to experimental error and not due to toxicity; (b) complete test with TA-1535 using rat lung homogenate because low population was observed at LU1 and LU2 doses in the first experiment; and (c) LU2 dose with TA-1535 using monkey tissue because low revertant count was obtained (less than 10) at this dose in the initial experiment which could be due to plating error.

B. Saccharomyces cerevisiae

1. Nonactivation suspension tests

The results of these tests were negative.


2. Activation suspension tests

The results of these tests were negative.

C. Conclusions

The test compound, 000126-95-4, did not exhibit mutagenic activity in any of the assays employed in these studies.


Submitted by:



David J. Brusick, Ph.D.
Director
Department of Genetics

5/11/77
Date

Reviewed by:



Robert J. Weir, Ph.D.
Vice President

5/11/77
Date



BIONETICS

VI. EXPLANATION OF EVALUATION PROCEDURES FOR PLATE ASSAYS

Plate test data consist of direct revertant colony counts obtained from a set of selective agar plates seeded with populations of mutant cells suspended in a semisolid overlay. Because the test chemical and cells are incubated in the overlay for 2-3 days, and a few cell divisions occur during the incubation period, the test is semiquantitative in nature. Although these features of the assay reduce the quantitation of results, they provide certain advantages not contained in a quantitative suspension test.

- The small number of cell divisions permits potential mutagens to act on replicating DNA which is often more sensitive than non-replicating DNA.
- The combined incubation of the compound and the cells in the overlay permit constant exposure of the indicator cells for 2-3 days.

A. Surviving Populations

Plate test procedures do not permit exact quantitation of the number of cells surviving chemical treatment. At low concentrations of the test chemical, the surviving population on the treatment plates is essentially the same as the negative control plate. At high concentrations, the surviving population is usually reduced by some fraction. Our protocol normally employs dose levels that are selected such that the highest dose will show slight toxicity (as determined by subjective criteria) and several doses ranging down 1 to 2 logs lower.

B. Dose Response Phenomena

The demonstration of dose-related increases in mutant counts is an important criterion in establishing mutagenicity. Factors which may modify dose response results for a mutagen would be the selection of doses that are too low (usually mutagenicity and toxicity are related). If the highest dose is far lower than a toxic concentration, no increases may be observed over the dose range selected. Conversely, if the lowest dose employed is highly cytotoxic, the test chemical may kill any mutants that are induced and the compound will not appear to be mutagenic.

C. Control Tests

Positive and negative control assays are conducted with each experiment and consist of direct acting mutagens for nonactivation assays and mutagens that require metabolic biotransformation in activation assays. Negative controls consist of the test compound solvent in the overlay agar with the other essential components. The negative control plate for each strain gives a reference point to which the test data are compared. The positive control assay is conducted to demonstrate that the test systems are functional with known mutagens.



D. Evaluation Criteria for Ames Assay

Because the procedures used to evaluate the mutagenicity of the test chemical are semiquantitative, the criteria used to determine positive effects are inherently subjective and are based primarily on a historical data base. Most data sets are evaluated using the following criteria:

1. Strains TA-1535, TA-1537, and TA-1538

If the solvent control value is within the normal range, a chemical that produces a positive dose response over three concentrations with the lowest increase equal to twice the solvent control value is considered to be mutagenic.

2. Strains TA-98, TA-100, and D4

If the solvent control value is within the normal range, a chemical that produces a positive dose response over three concentrations with the highest increase equal to twice the solvent control value for TA-100 and two to three times the solvent control value for strains TA-98 and D4 is considered to be mutagenic. For these strains, the dose response increase should start at approximately the solvent control value.

3. Pattern

Because TA-1535 and TA-100 were both derived from the same parental strain (G-46) and because TA-1538 and TA-98 were both derived from the same parental strain (D3052), there is a built-in redundancy in the microbial assay. In general the two strains of a set respond to the same mutagen and such a pattern is sought. It is also anticipated that if a given strain, e.g. TA-1537, responds to a mutagen in nonactivation tests it will generally do so in activation tests. (The converse of this relationship is not expected.) While similar response patterns are not required for all mutagens, they can be used to enhance the reliability of an evaluation decision.

4. Reproducibility

If a chemical produces a response in a single test that cannot be reproduced in one or more additional runs, the initial positive test data loses significance.

The preceding criteria are not absolute and other extenuating factors may enter into a final evaluation decision. However, these criteria are applied to the majority of situations and are presented to aid those individuals not familiar with this procedure. As the data base is increased, the criteria for evaluation can be more firmly established.

VII. EXPLANATION OF EVALUATION PROCEDURES FOR SUSPENSION ASSAYS

Data obtained from mutagenicity tests are evaluated on a test by test basis followed by an examination of the total response pattern using all the data. To facilitate this type of evaluation, we have prepared two separate formats in which data are processed. The first is the Compound Summary Backup Detail Sheet, which details the essential raw data from each experiment showing surviving population counts, total mutant or revertant counts, as well as, calculated mutation frequencies. This format permits close examination of each set of test data. The following considerations are part of any assessment.

A. Surviving Population Counts

A certain level of chemically-induced toxicity is anticipated, but occasionally isolated tests or groups of tests show very low (<25%) survival compared to the tissue controls. Such isolated decreases may result from improper dilution procedures or defective growth media and decrease confidence in the calculated mutation frequencies especially if the total mutant counts appear unaffected. Data of this type are generally unacceptable and these experiments are routinely repeated at a lower dose level to reduce killing and increase confidence in the nature of the response.

B. Total Mutant Counts

For nonmutagens, the mutant/surviving population ratio should be roughly equivalent for each test point in a given experiment. If the cell number drops in response to killing, the mutant number should decrease proportionately. A mutagenic chemical, however, will produce an altered mutant/surviving population ratio. Mutant numbers as well as calculated frequencies are compared to the negative control data. In certain instances, the mutant frequencies will increase with little or no change in the absolute number of mutants especially where the test chemical is toxic. Data of this type, although not necessarily aberrant, or even rare, must be viewed with special care to ensure that the increased frequencies were not the result of selective toxicity of the test chemical for the his⁻ cells. This phenomenon, referred to as selection, can lead to erroneous conclusions. Thus we attempt to keep the surviving population of cells high and look for positive responses that show increases in both numbers of mutants and mutation frequencies. Again, occasional isolated fluctuations in mutant counts are found that can be attributed to improper pipetting or media contamination. These fluctuations are usually easy to identify by inspection of the other data points in the experiment which will be negative.

C. Dose Response Phenomena

Dose-related increases in mutants and mutation frequencies are the most convincing data to have in assessing mutagenic activity of chemicals. In some cases, however, dose-related increases are not observed for mutagens. This depends considerably on the dose levels selected. The figure on the following page illustrates how one might obtain various types of dose-related responses by a mutagen based solely on dose selection. It also emphasizes the need to keep dose levels within a relatively low range of toxicity so that data are consistently on the uphill side of the hypothetical curve.

D. Control Tests


Positive and negative control tests are conducted with each experiment and consist of direct acting positive agents for nonactivation assays and chemicals that require metabolic transformation for activation assays. In nonactivation assays, the NAN control contain the test chemical solvent plus cells, but no chemical, and is used as a reference to assess the level of response obtained in the various tests. It is not possible at this time to put precise cut-off points where negative responses become positive responses. A statistical component for our computer program is under development and will be included when available. Positive controls are only used as relative reference points and to demonstrate that the system is functioning with known mutagens. In activation assays, three types of negative controls are run: (1) A solvent control minus the chemical and minus the activation system (A-C); (2) a control plus the positive control chemical minus the activation system (A+C); and (3) a control containing the activation system and the test chemical solvent (ALI or ALU). All three controls are used collectively to assess the level of response in the various activation tests. A chemical may appear positive when compared to an A-C control but not when compared to an A+T control. The value of each of the above controls with respect to their weight in evaluation is ALI or ALU > A-C > A+C.

The other data format is the Compound Frequency Summary Report sheet in which all the calculated frequencies obtained for a given compound are displayed in a table. This format permits an overview of all data. The points form a matrix of information that should present a consistent pattern. Nonmutagens should produce a matrix with data frequencies clustered around the negative control values. Occasional random high or low fluctuations are not uncommon and seldom indicate true genetic activity. Mutagenic chemicals should, on the other hand, produce a set of consistent responses that demonstrate a logical pattern. The patterns depend on the mutagenic specificity of the chemical but can be easily recognized in the Compound Frequency Summary Report format.

These mutagenicity assays are designed to optimize the probability of recognizing mutagens from nonmutagens and, in most cases, they work well. Occasionally, the data points are such that a definitive conclusion cannot be made without additional data.



HYPOTHETICAL MUTATION AND TOXICITY KINETICS



HYPOTHETICAL EXPERIMENT

- (1) Dose levels 1, 2 & 3 were used
- (2) Dose levels 2, 3 & 4 were used
- (3) Dose levels 3, 4 & 5 were used

OBSERVED DOSE RESPONSE

A typical positive dose response set of data would be obtained.

The intermediate dose level shows a higher mutation frequency than both the low dose and the high dose.

Here an inverted dose response would be observed with the highest dose level showing the lowest response.

APPENDIX
Tabulation of Data



BIONETICS

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		EXPERIMENT 633704		DETECTOR 1A100		SPECIES		PROJECT 2672		DATE - 03/14/77	
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM				
		NAN	SOLVENT	0797	0447	56.09	0				
		NAP	FMS 0.066%	0885	1455	164.41	0				
000126954	NA1		0005-0 PCT.	0469	0056	11.94	0				
000126954	NA2		0025-1 PCT.	0519	0114	21.97	0				
000126954	NA3		0125-2 PCT.	0473	0103	21.78	0				

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		EXPERIMENT 633402		DETECTOR TA1535		SPECIES		PROJECT 2672		DATE - 03/14/77	
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM				
	NAN		SOLVENT	0543	0029	5.34	0				
	NAP		EMS 0.2%	0633	3718	587.36	0				
000126954	NA1		0005-0 PCT.	0239	0020	8.37	0				
000126954	NA2		0025-1 PCT.	0491	0014	2.85	0				
000126954	NA3		0125-2 PCT.	0156	0018	11.54	0				

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		PROJECT 2672		DATE - 03/14/77			
EXPERIMENT 633404		DETECTOR TA1537		SPECIES /			
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	NAN		SOLVENT	0643	0129	20.06	0
	NAP		QM 13 UG/ML	0240	0331	137.92	0
000126954	NA1		0005-0 PCT.	0521	0034	6.53	0
000126954	NA2		0025-1 PCT.	0366	0043	11.75	0
000126954	NA3		0125-2 PCT.	0375	0073	19.47	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		PROJECT 2672		DATE - 03/14/77			
EXPERIMENT 633403	DETECTOR TA1538	SPECIES /					
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	NAN		SOLVENT	0456	0050	10.96	0
	NAP		NF 667 UG/ML	0458	0857	187.12	0
000126954	NA1		0005-0 PCT.	0824	0048	5.83	0
000126954	NA2		0025-1 PCT.	0558	0045	8.06	0
000126954	NA3		0125-2 PCT.	0455	0047	10.33	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

EXPERIMENT	CONTRACT	223-76-2102	DETECTOR	TA98	SPECIES	PROJECT	2672	DATE	- 03/14/77
	634120					/			
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8		CONTAM	
	NAN		SOLVENT	1397	0053	3.79		0	
	NAP		NF 667 UG/ML	1275	1801	141.25		0	
000126954	NA1		0005-0 PCT.	0425	0015	3.53		0	
000126954	NA2		0025-1 PCT.	0635	0017	2.68		0	
000126954	NA3		0125-2 PCT.	1242	0031	2.50		0	

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT	223-76-2102	PROJECT	2672	DATE -	03/14/77
EXPERIMENT	701303	DETECTOR	TA98	SPECIES	/
COMPOUND	TEST	OPG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0
	NAN		SOLVENT	0498	0065
000126954	NA1		0005-0 PCT.	0370	0064
					FREQ1 EP-8
					13.05
					17.30
					CONTAM
					0
					0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634902 DETECTOR 000004 / SPECIES / DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+4	MUT1 EP+1	MUT2 EP+1	FREQ1 EP-5	FREQ2 EP-5	CONTAM
		NAN	SOLVENT	0633	0125	0034	19.75	5.37	0
		NAP	FMS 1.0 %	0585	0594	0394	101.54	67.35	0
000126954	NA1		0143-2 PCT.	0964	0149	0037	15.46	3.84	0
000126954	NA2		0715-3 PCT.	0947	0072	0065	7.60	6.86	0
000126954	NA3		3575-4 PCT.	0939	0111	0047	11.82	5.01	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

EXPERIMENT 633802		CONTRACT 223-76-2102 DETECTOR TA100		SPECIES ICRFLO/MOUSE		PROJECT 2672	DATE - 03/14/77
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0201	0370	184.08	0
	A-C		SOLVENT	0591	0417	70.56	0
	ALI		TISSUE	0507	0557	109.86	0
	ALU		TISSUE	0747	0585	78.31	0
	ACP	LI	DMN 90 UM/ML	0450	1131	251.33	0
	ACP	LU	DMN 90 UM/ML	0524	0494	94.27	0
000126954	ACT	LI1	0005-0 PCT.	0702	0361	51.42	0
000126954	ACT	LI2	0025-1 PCT.	0998	0361	36.17	0
000126954	ACT	LI3	0125-2 PCT.	0915	0320	34.97	0
000126954	ACT	LU1	0005-0 PCT.	0807	0261	32.34	0
000126954	ACT	LU2	0025-1 PCT.	0468	0333	71.15	0
000126954	ACT	LU3	0125-2 PCT.	0609	0293	48.11	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 633501 DETECTOR TA1535 SPECIES ICRFLO/MOUSE DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0433	0039	9.01	0
	A-C		SOLVENT	0444	0037	8.33	0
	ALI		TISSUE	0459	0051	11.11	0
	ALU		TISSUE	0447	0039	8.72	0
	ACP	LI	DMN 90 UM/ML	0445	0800	179.78	0
	ACP	LU	DMN 90 UM/ML	0416	0049	11.78	0
000126954	ACT	LI1	0005-0 PCT.	0432	0043	9.95	0
000126954	ACT	LI2	0025-1 PCT.	0430	0044	10.23	0
000126954	ACT	LI3	0125-2 PCT.	0406	0027	6.65	0
000126954	ACT	LU1	0005-0 PCT.	0524	0046	8.78	0
000126954	ACT	LU2	0025-1 PCT.	0392	0019	4.85	0
000126954	ACT	LU3	0125-2 PCT.	0411	0037	9.00	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 701805 DETECTOR TA1537 SPECIES ICRFLO/MOUSE

DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		AMQ 333 UG/ML	0533	0029	5.44	0
	A-C		SOLVENT	0608	0029	4.77	0
	ALI		TISSUE	0449	0021	4.68	0
	ALU		TISSUE	0583	0033	5.66	0
	ACP	LI	AMQ 333 UG/ML	0466	0661	141.85	0
	ACP	LU	AMQ 333 UG/ML	0561	0040	7.13	0
000126954	ACT	LI1	0005-0 PCT.	0556	0029	5.22	0
000126954	ACT	LI2	0025-1 PCT.	0542	0039	7.20	0
000126954	ACT	LI3	0125-2 PCT.	0580	0036	6.21	0
000126954	ACT	LU1	0005-0 PCT.	0519	0024	4.62	0
000126954	ACT	LU2	0025-1 PCT.	0543	0032	5.89	0
000126954	ACT	LU3	0125-2 PCT.	0687	0048	6.99	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 701820 DETECTOR TA1538 SPECIES ICRFLO/MOUSE DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	0716	0054	7.54	0
	A-C		SOLVENT	0561	0061	10.87	0
	ALI		TISSUE	0514	0050	9.73	0
	ALU		TISSUE	0488	0049	10.04	0
	ACP	LI	ANTH 67 UG/ML	0531	1183	222.79	0
	ACP	LU	ANTH 67 UG/ML	0543	0366	67.40	0
000126954	ACT	LI1	0005-0 PCT.	0591	0066	11.17	0
000126954	ACT	LI2	0025-1 PCT.	0499	0067	13.43	0
000126954	ACT	LI3	0125-2 PCT.	0000	0000	00.00	0
000126954	ACT	LU1	0005-0 PCT.	0460	0080	17.39	0
000126954	ACT	LU2	0025-1 PCT.	0557	0083	14.90	0
000126954	ACT	LU3	0125-2 PCT.	0483	0077	15.94	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

		CONTRACT 223-76-2102		PROJECT 2672		
EXPERIMENT 701206		DETECTOR TA1538		SPECIES ICRFLO/MOUSE		DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	ALI		TISSUE	0317	0012	3.79	0
000126954	ACT	LT3	0125-2 PCT.	0364	0016	4.40	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

EXPERIMENT 634104		CONTRACT 223-76-2102 DETECTOR TA98		SPECIES ICRFLO/MOUSE		PROJECT 2672 DATE - 03/14/77	
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	1605	0041	2.55	0
	A-C		SOLVENT	1531	0044	2.87	0
	ALI		TISSUE	1051	0061	5.80	0
	ALU		TISSUE	0941	0047	4.99	0
	ACP	LI	ANTH 67 UG/ML	1181	0685	58.00	0
	ACP	LU	ANTH 67 UG/ML	1010	0252	24.95	0
000126954	ACT	LI1	0005-0 PCT.	0838	0030	3.58	0
000126954	ACT	LI2	0025-1 PCT.	0560	0051	9.11	0
000126954	ACT	LI3	0125-2 PCT.	0258	0038	14.73	0
000126954	ACT	LU1	0005-0 PCT.	0314	0029	9.24	0
000126954	ACT	LU2	0025-1 PCT.	0570	0039	6.84	0
000126954	ACT	LU3	0125-2 PCT.	0846	0045	5.32	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634805 DETECTOR 000004 SPECIES ICRFLO/MOUSE DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP-4	MUT1 EP-1	MUT2 EP-1	FREQ1 EP-5	FREQ2 EP-5	CONTAM
	A+C		DMN 90 UM/ML	0748	0105	0084	14.04	11.23	0
	A-C		SOLVENT	0578	0098	0079	16.96	13.67	0
	ALI		TISSUE	0650	0128	0084	19.69	12.92	0
	ALU		TISSUE	0583	0119	0091	20.41	15.61	0
	ACP	LI	DMN 90 UM/ML	0386	0295	0297	76.42	76.94	0
	ACP	LU	DMN 90 UM/ML	0681	0125	0075	18.36	11.01	0
000126954	ACT	LI1	0143-2 PCT.	0638	0074	0054	11.60	8.46	0
000126954	ACT	LI2	0715-3 PCT.	0769	0083	0060	10.79	7.80	0
000126954	ACT	LI3	3575-4 PCT.	0697	0082	0044	11.76	6.31	0
000126954	ACT	LU1	0143-2 PCT.	0635	0096	0049	15.12	7.72	0
000126954	ACT	LU2	0715-3 PCT.	0760	0080	0033	10.53	4.34	0
000126954	ACT	LU3	3575-4 PCT.	0701	0098	0064	13.98	9.13	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		PROJECT 2672					
EXPERIMENT 634208	DETECTOR TA100	SPECIES SPRDAW/RAT				DATE - 03/14/77	
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0441	0710	161.00	0
	A-C		SOLVENT	0644	0828	129.57	0
	ALI		TISSUE	0643	0825	128.30	0
	ALU		TISSUE	0717	0800	111.58	0
	ACP	LI	DMN 90 UM/ML	0365	0635	173.97	0
	ACP	LU	DMN 90 UM/ML	0581	0651	112.05	0
000126954	ACT	LI1	0005-0 PCT.	0574	0279	48.61	0
000126954	ACT	LI2	0025-1 PCT.	0601	0429	71.38	0
000126954	ACT	LI3	0125-2 PCT.	0686	0331	48.25	0
000126954	ACT	LU1	0005-0 PCT.	0452	0228	50.44	0
000126954	ACT	LU2	0025-1 PCT.	0697	0423	60.69	0
000126954	ACT	LU3	0125-2 PCT.	0645	0334	51.78	0

REPORT EXR33 LITTON RIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 703901 DETECTOR TA1535 SPECIES SPRDAM/RAT DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0343	0024	7.00	0
	A-C		SOLVENT	0180	0014	7.78	0
	ALI		TISSUE	0527	0014	2.66	0
	ALU		TISSUE	0545	0021	3.85	0
	ACP	LI	DMN 90 UM/ML	0422	0585	138.63	0
	ACP	LU	DMN 90 UM/ML	0472	0021	4.45	0
000126954	ACT	L11	0005-0 PCT.	0547	0027	4.94	0
000126954	ACT	L12	0025-1 PCT.	0761	0037	4.86	0
000126954	ACT	L13	0125-2 PCT.	1153	0033	2.86	0
000126954	ACT	LU1	0005-0 PCT.	0290	0024	8.28	0
000126954	ACT	LU2	0025-1 PCT.	0265	0015	5.66	0
000126954	ACT	LU3	0125-2 PCT.	1163	0031	2.67	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 704503 DETECTOR TA1535 SPECIES SPRDAW/RAT DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
		ALU	TISSUE	0493	0040	8.11	0
000126954	ACT	LU1	0005-0 PCT.	0346	0014	4.05	0
000126954	ACT	LU2	0025-1 PCT.	0415	0047	11.33	0
000126954	ACT	LU3	0125-2 PCT.	0460	0056	12.17	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634206 DETECTOR TA1537 SPECIES SPRDAM/RAT DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		AMQ 333 UG/ML	1460	0055	3.77	0
	A-C		SOLVENT	1438	0023	1.60	0
	ALI		TISSUE	1075	0022	2.05	0
	ALU		TISSUE	1157	0033	2.85	0
	ACP	LI	AMQ 333 UG/ML	0382	0340	89.01	0
	ACP	LU	AMQ 333 UG/ML	1248	0102	8.17	0
000126954	ACT	L11	0005-0 PCT.	1060	0046	4.34	0
000126954	ACT	L12	0025-1 PCT.	0780	0032	4.10	0
000126954	ACT	L13	0125-2 PCT.	1160	0029	2.50	0
000126954	ACT	LU1	0005-0 PCT.	0802	0027	3.37	0
000126954	ACT	LU2	0025-1 PCT.	1231	0030	2.44	0
000126954	ACT	LU3	0125-2 PCT.	0969	0034	3.51	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 635201 DETECTOR TA1538 SPECIES SPRDAW/RAT DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	0528	0243	46.02	0
	A-C		SOLVENT	0587	0058	9.88	0
	ALI		TISSUE	0539	0073	13.54	0
	ALU		TISSUE	0626	0077	12.30	0
	ACP	LI	ANTH 67 UG/ML	0515	0662	128.54	0
	ACP	LU	ANTH 67 UG/ML	0488	0159	32.58	0
000126954	ACT	L11	0005-0 PCT.	0523	0075	14.34	0
000126954	ACT	L12	0025-1 PCT.	0550	0060	10.91	0
000126954	ACT	L13	0125-2 PCT.	0498	0050	10.04	0
000126954	ACT	LU1	0005-0 PCT.	0515	0046	8.93	0
000126954	ACT	LU2	0025-1 PCT.	0633	0042	6.64	0
000126954	ACT	LU3	0125-2 PCT.	0689	0075	10.89	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		PROJECT 2672				DATE - 03/14/77	
EXPERIMENT 634207		DETECTOR TA98		SPECIES SPRDAW/RAT			
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	1729	0046	2.66	0
	A-C		SOLVENT	0973	0062	6.37	0
	ALI		TISSUE	0360	0081	22.50	0
	ALU		TISSUE	0643	0062	9.64	0
	ACP	LI	ANTH 67 UG/ML	0227	0680	299.56	0
	ACP	LU	ANTH 67 UG/ML	0728	0359	49.31	0
000126954	ACT	LI1	0005-0 PCT.	0488	0056	11.48	0
000126954	ACT	LI2	0025-1 PCT.	0433	0039	9.01	0
000126954	ACT	LI3	0125-2 PCT.	0413	0040	9.69	0
000126954	ACT	LU1	0005-0 PCT.	0616	0066	10.71	0
000126954	ACT	LU2	0025-1 PCT.	0711	0032	4.50	0
000126954	ACT	LU3	0125-2 PCT.	0784	0048	6.12	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102		PROJECT 2672		DATE - 03/14/77					
EXPERIMENT 634205		DETECTOR 000004		SPECIES SPRDAW/RAT					
COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+4	MUT1 EP+1	MUT2 EP+1	FREQ1 EP-5	FREQ2 EP-5	CONTAM
	A+C		DMN 90 UM/ML	0638	0087	0065	13.64	10.19	0
	A-C		SOLVENT	0550	0088	0061	16.00	11.09	0
	ALI		TISSUE	0677	0128	0063	18.91	9.31	0
	ALU		TISSUE	0646	0110	0056	17.03	8.67	0
	ACP	LI	DMN 90 UM/ML	0442	0429	0320	97.06	72.40	0
	ACP	LU	DMN 90 UM/ML	0601	0061	0019	10.15	3.16	0
000126954	ACT	L11	0143-2 PCT.	0702	0091	0055	12.96	7.83	0
000126954	ACT	L12	0715-3 PCT.	0672	0068	0058	10.12	8.63	0
000126954	ACT	L13	3575-4 PCT.	0690	0069	0045	10.00	6.52	0
000126954	ACT	LU1	0143-2 PCT.	0690	0081	0050	11.74	7.25	0
000126954	ACT	LU2	0715-3 PCT.	0667	0078	0038	11.69	5.70	0
000126954	ACT	LU3	3575-4 PCT.	0590	0084	0040	14.24	6.78	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634209 DETECTOR TA100 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0336	0179	53.27	0
	A-C		SOLVENT	0468	0268	57.26	0
	ALI		TISSUE	0485	0400	82.47	0
	ALU		TISSUE	0436	0318	72.94	0
	ACP	LI	DMN 90 UM/ML	0350	0683	195.14	0
	ACP	LU	DMN 90 UM/ML	0306	0312	101.96	0
000126954	ACT	LI1	0005-0 PCT.	0400	0363	90.75	0
000126954	ACT	LI2	0025-1 PCT.	0333	0443	133.03	0
000126954	ACT	LI3	0125-2 PCT.	0285	0372	130.53	0
000126954	ACT	LU1	0005-0 PCT.	0233	0295	126.61	0
000126954	ACT	LU2	0025-1 PCT.	0400	0437	109.25	0
000126954	ACT	LU3	0125-2 PCT.	0356	0259	72.75	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634802 DETECTOR TAI535 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		DMN 90 UM/ML	0854	0022	2.58	0
	A-C		SOLVENT	0905	0020	2.21	0
	ALI		TISSUE	1003	0034	3.39	0
	ALU		TISSUE	0792	0022	2.78	0
	ACP	LI	DMN 90 UM/ML	0760	0279	36.71	0
	ACP	LU	DMN 90 UM/ML	0930	0022	2.37	0
000126954	ACT	L11	0005-0 PCT.	0882	0031	3.51	0
000126954	ACT	L12	0025-1 PCT.	0783	0024	3.07	0
000126954	ACT	L13	0125-2 PCT.	0757	0014	1.85	0
000126954	ACT	LU1	0005-0 PCT.	1118	0019	1.70	0
000126954	ACT	LU2	0025-1 PCT.	0655	0007	1.07	0
000126954	ACT	LU3	0125-2 PCT.	0655	0031	4.73	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

		CONTRACT 223-76-2102		PROJECT 2672		
EXPERIMENT 701207		DETECTOR TA1535		SPECIES RHESUS/MONKEY		DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
		ALI	TISSUE	0580	0049	8.45	0
000126954	ACT	LU2	0025-1 PCT.	0572	0028	4.90	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 634210 DETECTOR TA1537 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		AMQ 333 UG/ML	2044	0050	2.45	0
	A-C		SOLVENT	1436	0035	2.44	0
	ALI		TISSUE	1784	0037	2.07	0
	ALU		TISSUE	2176	0024	1.10	0
	ACP	LI	AMQ 333 UG/ML	0454	0370	81.50	0
	ACP	LU	AMQ 333 UG/ML	1803	0084	4.66	0
000126954	ACT	LI1	0005-0 PCT.	1271	0022	1.73	0
000126954	ACT	LI2	0025-1 PCT.	0966	0024	2.48	0
000126954	ACT	LI3	0125-2 PCT.	1239	0027	2.18	0
000126954	ACT	LU1	0005-0 PCT.	1385	0019	1.37	0
000126954	ACT	LU2	0025-1 PCT.	0736	0011	1.49	0
000126954	ACT	LU3	0125-2 PCT.	1191	0011	0.92	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 635102 DETECTOR TA1538 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	0501	0049	9.78	0
	A-C		SOLVENT	0549	0076	13.84	0
	ALI		TISSUE	0602	0074	12.29	0
	ALU		TISSUE	0569	0068	11.95	0
	ACP	LI	ANTH 67 UG/ML	0529	0725	137.05	0
	ACP	LU	ANTH 67 UG/ML	0493	0076	15.42	0
000126954	ACT	LI1	0005-0 PCT.	0144	0015	10.42	0
000126954	ACT	LI2	0025-1 PCT.	0459	0024	5.23	0
000126954	ACT	LI3	0125-2 PCT.	0457	0052	11.38	0
000126954	ACT	LU1	0005-0 PCT.	0187	0014	7.49	0
000126954	ACT	LU2	0025-1 PCT.	0432	0070	16.20	0
000126954	ACT	LU3	0125-2 PCT.	0474	0048	10.13	0

REPORT EXR33 LITTON BIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672.
EXPERIMENT 634302 DETECTOR TA98 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+6	MUT1 EP+0	FREQ1 EP-8	CONTAM
	A+C		ANTH 67 UG/ML	2229	0129	5.79	0
	A-C		SOLVENT	1691	0093	5.50	0
	ALI		TISSUE	0973	0093	9.56	0
	ALU		TISSUE	1355	0144	10.63	0
	ACP	LI	ANTH 67 UG/ML	1137	1402	123.31	0
	ACP	LU	ANTH 67 UG/ML	1308	0228	17.43	0
000126954	ACT	LI1	0005-0 PCT.	1439	0071	4.93	0
000126954	ACT	LI2	0025-1 PCT.	1393	0092	6.60	0
000126954	ACT	LI3	0125-2 PCT.	1300	0163	12.54	0
000126954	ACT	LU1	0005-0 PCT.	1915	0034	1.78	0
000126954	ACT	LU2	0025-1 PCT.	1892	0061	3.22	0
000126954	ACT	LU3	0125-2 PCT.	2252	0078	3.46	0

REPORT EXR33 LITTON RIONETICS MUTAGENIC ACTIVITY SYSTEM
COMPOUND SUMMARY BACKUP DETAIL

CONTRACT 223-76-2102 PROJECT 2672
EXPERIMENT 701806 DETECTOR 000004 SPECIES RHESUS/MONKEY DATE - 03/14/77

COMPOUND	TEST	ORG ID	CONCENTRATION	POPU EP+4	MUT1 EP+1	MUT2 EP+1	FREQ1 EP-5	FREQ2 EP-5	CONTAM
	A+C		DMN 90 UM/ML	0742	0144	0079	19.41	10.65	0
	A-C		SOLVENT	0790	0146	0075	18.48	9.49	0
	ALI		TISSUE	0770	0127	0084	16.49	10.91	0
	ALU		TISSUE	0698	0107	0071	15.33	10.17	0
	ACP	LI	DMN 90 UM/ML	0414	0387	0304	93.48	73.43	0
	ACP	LU	DMN 90 UM/ML	0387	0100	0066	25.84	17.05	0
000126954	ACT	L11	0143-2 PCT.	0613	0105	0043	17.13	7.01	0
000126954	ACT	L12	0715-3 PCT.	0618	0109	0061	17.64	9.87	0
000126954	ACT	L13	3575-4 PCT.	0581	0102	0038	17.56	6.54	0
000126954	ACT	LU1	0143-2 PCT.	0626	0117	0062	18.69	9.90	0
000126954	ACT	LU2	0715-3 PCT.	0730	0135	0060	18.49	8.22	0
000126954	ACT	LU3	3575-4 PCT.	0652	0102	0081	15.64	12.42	0