




# InterStrand®

**InterStrand®** consists of 10 **InterSource®** seeds threaded onto a monofilament absorbable suture. The seeds are spaced 1 cm apart measured center to center. **InterStrand®** easily penetrates a standard bone wax plug, minimizes seed migration and reduces prep time.

**InterStrand® Special** offers the same robust features as **InterStrand®** with the added flexibility of specifying the number and placement of seeds along the strand. This allows improved coverage of the base and apex of the gland without risk of increased seed migration while reducing the number of needles required in the procedure.

InterSource® hollow seed  
I-125 or Pd-103

Seeds spaced at 1 cm  
distance center to center



## **Benefits**

- Easily Penetrates Bone Wax
- Maintains Seed Spacing
- No Adverse Effects from Fluids
- Reduces Preparation Time
- Minimizes Seed Migration